DynaScope 7000 Series

Bedside Monitor



Ver.08

Operation Manual

《 Monitoring Operation 》



- Before using this device, read this operation manual thoroughly.
- Keep this manual near the device for future reference.



This operation manual is for the DS-7200 System Version 08.

CAUTION FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

CAUTION:

- This device for sale by or on the order of a physician.
- The company and product names used in this manual are trademarks or registered trademarks.
- · If this manual has pages missing or out of order, contact Fukuda Denshi for replacement.
- · Only physician or persons instructed by physicians are allowed to use the equipment.
- The information contained in this document is subject to change without notice due to improvement in the equipment.

Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or automated sphygmomanometers.

Copyright © 2013 by Fukuda Denshi Co., Ltd.

No part of this document may be copied or transmitted in any form without the prior written permission of Fukuda Denshi Co., Ltd.

Printed in Japan

Preface

Thank you for purchasing this product. Before using this product, read the following precautions to make sure the product is used correctly and safely.

	ition Manualii
	iii
	nitiii
Measurement Unit for Ea	ch Parameter ····· vi
	······viii
Precautions for Safe Ope	ration of Medical Electrical
• •	xi
-	eration of Medical Telemetry xii
	aintenance·····xiii
	cemakerxiii
	xiv
	xiv
.	xiv
	etic Resonance Imaging xv
	ections to Peripheral Devicesxv
	se ····· xv
	I Accessories xv
	S-7200 System ·····xvi
	red Network System (DS-LAN
II/DS-LANIII)·····	
	reless Network System ·····xxx
Precautions for Use of the	
	entilator Monitoring xxxiii
	OO2 Sensor ······xxxv
	Model: DS-7210M······xxxv
	BP Cuffxxxv
Disposing of Equipment,	Accessories, or Components
	portationxxxv
	or Data Backup ······xxxv
	thium-Ion Battery Pack xxxvi
	cy Use xxxvii
	ilityxxviii
	eration under Electromagnetic
	xxxviii
	Electromagnetic Emissions
	······································
	Electromagnetic Immunity (1)
	······································
	Electromagnetic Immunity (2)xl
	paration Distances between
	e RF Communications
Equipment and the	DS-7200 System ······xli

Composition of This Operation Manual

The DS-7200 System Operation Manual is composed of the following 3 sections.

≪General Description≫

This section is composed of the chapters stating the general description of the device and basic operation procedure.

1. General Description	:	Describes the outline of this equipment.
2. Basic Operation	:	Describes the basic operation for monitoring.
3. Vital Application	:	Describes the procedure for vital application, etc.

: Describes the procedure for vital application, etc.

≪Monitoring Operation≫

This section is composed of the chapters explaining the detailed monitoring procedures and setup procedures.

4. Monitoring Setup	:	Describes the procedures to set the monitor according to the monitoring purpose.
5. Admit / Discharge of a Patient	:	Describes the procedure to admit or discharge a patient.
6. Parameter Setup	:	Describes the procedure to set the measurement condition, size, scale, etc. for each parameter.
7. Function	:	Describes about the functions such as arrhythmia analysis, trend, recall, etc.
8. System Configuration	:	Describes about the system configuration such as night mode, alarm mode, display mode, etc.

«Maintenance»

This section is composed of the chapters describing the installation procedure, maintenance, technical information, accessories, etc.

9. Installation 10. Maintenance	:	Describes about the environment for use, wireless system, etc. Describes about the maintenance, troubleshooting of this equipment.
11. Technical Information		Lists the specification, default settings, pin assignments of external connector, etc.
12. Accessories	:	Lists the accessories and optional accessories for this equipment.

Safety Precautions

- Read the "Safety Precautions" thoroughly before use to ensure correct and safe use of the product.
- Be sure to follow the precautions indicated below, as these are important messages related to safety.

▲ DANGER	Failure to follow this message may cause immediate threat of death or serious injury, or complete failure of the equipment.
≜ WARNING	Failure to follow this message may result in death or serious injury, or complete failure of the equipment.
A CAUTION	Failure to follow this message may cause injury or failure to the equipment.
NOTE	A note is not related to product safety, but provides information about the correct use and operating procedures to prevent incorrect operation and malfunction of the equipment.

Labels Attached to the Unit

Make sure to read the warning labels attached to the unit and comply with these requirements while operating the unit.

	Do not damage or erase the warning labels attached to the unit.
A CAUTION	These warning labels contain descriptions important for handling and operating the unit properly and safely. A damaged label may compromise safe operation.

DS-7200 System

A DANGER

Risk of explosion if used in the presence of flammable anesthetics.

▲ CAUTION

Before connecting, read instruction manual.

▲ CAUTION

To reduce the risk of electric shock, do not remove the cover. Refer servicing to qualified service personnel.



HU-71/HU-72/HU-73 Option Unit

Risk of explosion if used in the presence of flammable anesthetics.

Before connecting, read instruction manual.

To reduce the risk of electric shock, do not remove the cover. Refer servicing to qualified service personnel.



OAO-12B Battery Pack



形 式 / Type : OAO-12B	8 2	/ Capacity : 6600mA·h	
		·	
電 圧 / Voltage : 14.8V DC	製造番号	/ Lot No.	
▲ 危険		∆ DANGER	
電池/やクを機能に扱称するときは、コネクタの向きを確かめぐ確実に装着して下さい この電池/や/うな構設には接用しないで下さい。 電池の汚電は指定の機器本体のみで行って下さい。 電性であるのショントやな違い。資料をしないで下さい。 少りンプにジンとごを置いめのと、44にを買しないで下さい。 電池になき命があります。34化した電池/やクは使用しないで下さい。 電池/やって活体化した/いンマーのいした18点の1代か。845+745mL/Lの.キズをつけ 電池/やって活動したり、24でのしたり自らっけたりしないで下さい。	When instaling the battery pack to the equipment, ensure the connect Do not use the battery pack with an equipment other than spacefied. The battery must be charged on specified equipment. Do not short the electrode or termming, or mendel/biassemble the bat Do not throw into the first, heat, or leave/charge the battery under high Do not thore the battery with metal such as do or pin. The battery detectrates with time. Do not use the detectrated battery Do not there a rail in, this harmer, steps on the battery pack, or pr Do not apply storing impact or throw the battery pack.	ttery. temperature. pack.	
・・ 限りあるこれらの資源の有効活用のため、リサイクルにご協力ください。 リサイクル処理は弊社販売店・代理店に提出いただくか、各自治体の処理方法に従って	ください。		

Measurement Unit for Each Parameter

Detail	Parameter	Display	Unit	Default
	ECG	HR	bpm	
Heart Pote / Dulce	Invasive Blood Pressure	PR_IBP	bpm	
Heart Rate / Pulse Rate	SpO ₂	PR_SpO ₂	bpm	
Nate	Non-Invasive Blood Pressure	PR_NIBP	bpm	
ST Level	ECG	ST	mm, mv	mm
VPC	ECG	VPC	bpm	
	Impedance Respiration	RR IMP	Bpm	
Respiration Rate		RR_CO ₂	Bpm	
-	Ventilator	RR_VENT	Bpm	
	Impedance Respiration	APNEA	s (second)	
Apnea	CO ₂	APNEA	s (second)	
	Ventilator	APNEA	s (second)	
Invasive Blood Pressure	Invasive Blood Pressure	BP	mmHg, kPa cmH ₂ O (CVP only)	mmHg
Non-Invasive Blood Pressure	Non-Invasive Blood Pressure	NIBP	mmHg, kPa	mmHg
Arterial Oxygen	SpO ₂	SpO ₂	%	
Saturation	Perfusion Index	PI	%	
Temperature	Temperature	TEMP	°C / °F	°C
End-Tidal CO ₂ Concentration	CO ₂	EtCO ₂	mmHg, kPa, %	mmHg
Inspiratory CO ₂ Concentration	CO ₂	InspCO ₂	mmHg, kPa, %	mmHg
	Cardiac Output	CO	L/minute	
Cardiac Output	Cardiac Index	CI	L/minute/m ²	
Blood Temperature	Blood Temperature	Tb	°C / °F	°C
Injectate Temperature	Injectate Temperature	Ti	°C / °F	°C
Airway Flow	Airway Flow	AWF	L/minute	
Airway Pressure	Airway Pressure	AWP	cmH ₂ O	
	Expiratory Tidal Volume	E_TV	mL	
Tidal Volume	Inspiratory Tidal Volume	I_TV	mL	
	Tidal Volume	TV	mL	
	Inspiratory/Expiratory Ratio	I:E	(none)	
Respiratory Minute	Minute Volume	MV	L/minute	
Volume	Spontaneous Minute Volume	SMV	L/minute	
.	Compliance	COMP	mL/cmH ₂ O	
Compliance	Static Compliance	S_COMP	mL/cmH ₂ O	
	Dynamic Compliance	D_COMP	mL/cmH ₂ O	
	Expiratory Resistance	E_RES	cmH ₂ O/L/Sec	
Airway Resistance	Inspiratory Resistance	I_RES	cmH ₂ O/L/Sec	
	Static Airway Resistance	S_RES	cmH ₂ O/L/Sec	
	Dynamic Airway Resistance	D_RES	cmH ₂ O/L/Sec	
	Mean Airway Pressure	MEAN	cmH₂O	
Airway Pressure	Maximum Airway Pressure	PEAK	cmH ₂ O	
An way FIESSULE	Pause Airway Pressure	PAUSE	cmH ₂ O	
	Minimum Airway Pressure	P_Min	cmH ₂ O	

The measurement units for this equipment are as follows.

bpm: beats per minute Bpm: breaths per minute

Detail	Parameter	Display	Unit	Default
Spontaneous	Spontaneous Respiration	S_RR	Bpm	
Respiration		•		
Peak End Expiratory Pressure	Peak End Expiratory Pressure	PEEP	cmH₂O	
Fraction of	Fraction of Inspiratory			
Inspiratory Oxygen	Oxygen	FIO ₂	%	
mophatory exygen	Mixed Venous Oxygen		0/	
	Saturation	SvO ₂	%	
	Central Venous Oxygen	ScvO ₂	%	
	Saturation	-		
	Arterial Oxygen Saturation	SaO ₂	%	
	Oxygen Uptake Index	O ₂ EI	%	
	Oxygen Transport	DO ₂	mL/minute	
	Oxygen Consumption	VO ₂	mL/minute	
	Stroke Volume	SV	mL	
	Stroke Volume (STAT Mode)	SV_STAT	mL	
	Stroke Volume Index	SVI	mL/m ²	
	Stroke Volume Index (STAT Mode)	SVI_STAT	mL/m ²	
	Heart Rate	HR	bpm (beats per minute)	
	Mean Arterial Pressure	MAP	mmHg	
	Central Venous Pressure	CVP	mmHg	
Vigilance Data	Continuous Cardiac Output	CCO	L/minute	
 Vigilance Vigilance CEDV 	Continuous Cardiac Output (STAT Mode)	CCO_STAT	L/minute	
 VigilanceII 	Continuous Cardiac Index	CCI	L/minute/m ²	
• Vigileo	Continuous Cardiac Index (STAT Mode)	CCI_STAT	L/minute/m ²	
	Systemic Vascular Resistance	SVR	dynes-sec/cm ⁵	
	Systemic Vascular Resistance Index	SVRI	dynes-sec/cm⁵	
	Blood Temperature	BT	°C	
	Ejection Fraction	EF	%	
	Ejection Fraction (STAT Mode)	EF_STAT	%	
	End-Diastolic Volume	EDV	mL	
	End-Diastolic Volume (STAT Mode)	EDV_STAT	mL	
	End-Diastolic Volume Index	EDVI	mL/m ²	
	End-Diastolic Volume Index (STAT Mode)	EDVI_STAT	mL/m ²	
	End-Systolic Volume	ESV	mL	
	End-Systolic Volume Index	ESVI	mL	
	Stroke Volume Variance	SVV	%	
	Bispectral Index	BIS	(no unit)	
BIS Monitor Data	Signal Quality Index	SQI	%	
	Electromyograph	EMG	dB	
	Suppression Ratio	SR	%	

Graphic Symbols

Refer following for the meaning of the symbols indicated on the equipment.

Symbol	Description
	Caution; refer to accompanying documents Indicates the need to refer to related accompanying documents before operation.
Å	Equipotential Terminal Indicates the terminal to equalize the potential difference when interconnecting the devices.
\otimes	Inhibition The operation is inhibited. Refer to the instruction.
	Protective Earth Indicates the protective earth inside the equipment.
\sim	Alternating Current (Main Power Input Indicator)
	Direct Current
<u> </u>	Battery Charge (Battery Charge Indicator)
Ò	"OFF" for a Part of an Equipment Indicates the "OFF" condition for a part of an equipment.
\odot	"ON" for a Part of an Equipment Indicates the "ON" condition for a part of an equipment.
Ara A	Electrostatic Sensitive Part Directly touching this connector part with hands should be avoided.
ł	Type CF Applied Part with Defibrillation-Proof Indicates the degree of protection against electric shock is Type CF Applied Part with defibrillation-proof.
۱ ۲ ۲	Type BF Applied Part with Defibrillation-Proof Indicates the degree of protection against electric shock is Type BF Applied Part with defibrillation-proof.
Ŕ	Type BF Applied Part Indicates the degree of protection against electric shock is Type BF Applied Part.
\ominus	Signal Output Part
	GAS Output Part

Symbols indicated on the equipment

Symbol	Description
-	Signal Input Part
	Manufactured Date
- 2 -2-	TCP/IP Network Connector Connects to TCP/IP network.
	RS-232C Connector Connects the related device.
	Eject Indicates the switch to remove the recorder paper cassette.

Symbols displayed on the screen

Symbol	Description
	Battery Mark During battery operation, battery status will be displayed.
×	Alarm OFF Indicates the alarm is OFF.
•	Heart Rate Synchronization Mark This mark flashes synchronizing to the heartbeat.
N	Respiration Synchronization Mark This mark flashes synchronizing to the inspiration.
æ	Event Key This mark will be displayed when an alarm generates. Whether or not to display this icon can be selected on the monitor setup menu.
	Device Configuration Icon This mark will be displayed when device configuration has changed. Whether or not to display this icon can be selected on the monitor setup menu.
θ	Message Icon This mark will be displayed inside the parameter key when an alarm message is present for that parameter. Whether or not to display this icon can be selected on the monitor setup menu.
𝓲 𝔄 𝓲 𝓲×	TCON Antenna Mark Indicates the receiving condition of the Bidirectional Wireless Communication Module (HTC-702).
000000	SEC Alarm Display Indicates the SEC alarm status.
▶◀▲▼ ⋫⋞⋨¥	Scroll Keys These keys will allow to scroll the screen.
ð	Laser Printer This mark will be displayed when a laser printer connected to the TCP/IP network is used.
	Laser Printer Output Indicates the current printing progress.

Precautions for Safe Operation of Medical Electrical Equipment

▲ CAUTION	 Read the following precautions thoroughly to correctly operate the device. Users should have a thorough knowledge of the operation before using this system. Pay attention to the following when installing and storing the equipment. Do not install or store in an area where the equipment will be subject to splashing water. Do not install or store in an area where the environmental conditions, such as atmospheric pressure, temperature, humidity, ventilation, sunlight, dust, sodium, sulfur, will adversely affect the system. Place the equipment on a stable surface where there is no inclination, vibration, or shock (including during transportation). Do not install or store in an area where there are chemical or gasses stored. Verify the power frequency, voltage and allowable current (or power consumption). Ensure the grounding is proper by connecting the accompanying power cable to the hospital grade outlet. Do not install the equipment in a location where it is difficult to unplug the power cable. Before operating the system, verify the following items. Verify the power voltage. Check the cable connection and polarity to ensure proper operation of the equipment. Make sure the power system has adequate earth ground. Ensure that all cables are firmly and safely connected. Pay special attention when the device is used in conjunction with other equipment as it may cause erroneous judgment and danger. Ensure all patient connections are proper and secure. During operation of the system, verify the following items. Always observe the system and patient to ensure safe operation of the equipment. Do not all the cables from the patient to ensure safe operation of the safest way for the patient. Do not allow the patient. Do no tall allowable come in contact with the device. After using the system, verify the following items.
▲ CAUTION	 Pay special attention when the device is used in conjunction with other equipment as it may cause erroneous judgment and danger. Ensure all patient connections are proper and secure. During operation of the system, verify the following items. Always observe the system and patient to ensure safe operation of the equipment. If any abnormality is found on the equipment or patient, take appropriate measures such as ceasing operation of the equipment in the safest way for the patient. Do not allow the patient to come in contact with the device. After using the system, verify the following items. Unplug all the cables from the patient before turning off the power. When unplugging the cables, do not apply excessive force by pulling on the cord. Pull by the connector part of the cable. Clean the accessories and cables, and keep them together in one

Precautions for Safe Operation of Medical Telemetry

	Precautions for Safe Operation of Medical Telemetry
	To operate the device correctly, read the following precautions carefully.
	 The medical institution (hereinafter referred as "Institution") must decide the telemetry installation plan for the medical institution in order to prevent interference and interference between transmitters (telemetry based on destination country's radio law).
	• When using telemetry which requires zone location, the institution is to set up the zones as an operation unit for each transmitter to prevent electronic interference between telemetry throughout the medical institution.
	 When using telemetry which requires zone location, display and identify each prepared zone in the equipment.
	• When laying receiver antenna for each transmitter, the institution has to be examined so as not to generate electronic interference.
	 Based on the above examination result, the institution places each receiver antenna as required.
	In managing, be sure to follow the precautions below.
	• The institution appoints a person to manage the wireless channels for the whole medical institution. And when using telemetry which requires zone location, the institution nominates a person to manage the wireless channels in each zone (a "Zone Manager"). However, when using such telemetry in a local medical institution, one person can perform both functions.
A CAUTION	 Select a telemetry manager who understands the characteristics and functionality of telemetry systems, and is skilled in operating telemetry.
	 When installing telemetry, the Overall Manager and the Zone Manager have to understand the precautions for use of the telemetry in advance.
	 The Overall Manager takes responsibility of wireless channel management and transmitter storage for the whole medical institution by giving proper instruction.
	• The Overall Manager creates a management log, list of wireless channels, management status for the whole medical institution (hereinafter referred to as the "management log"). When changing a wireless channel, register it in the log and give proper instructions to the zone manager or to the user.
	 The Zone Manager assumes responsibility for managing the wireless channels, storing, and managing telemetry.
	 The Zone Manager assigns the transmitter to the user, and provides enough education for use inside the zone.
	• The telemetry user verifies operation of the transmitter/receiver before use.
	 The telemetry user, if using the telemetry in a zone location, follows the instructions of the zone manager for the zone and gives instructions to the patient if required.
	• When interference or breakdown occurs in telemetry communication, the user is required to inform the zone manager and the overall manager of the problems. The zone manager and overall manager are to deal with the problem properly and/or contact their nearest Fukuda Denshi representative for service.

Precautions about the Maintenance

Safety Inspection and Maintenance

For safe operation of the equipment, regular inspection and maintenance is required. Once a year, check all cables, devices, and accessories for damage, earth impedance, earth and leakage currents, and all alarm functions. Also, ensure that all safety labels are legible. Maintain a record of these safety inspections.

Immediate maintenance has to be carried out if ;

- the equipment was subjected to extreme mechanical stress, e.g. after a heavy fall.
- the equipment was subjected to liquid spill.
- the monitoring function is interrupted or disturbed.
- parts of the equipment enclosure are cracked, removed, or lost.
- any connector or cable shows signs of deterioration.

 Reference
 Refer to "10. Maintenance" for details.

 Maintenance
 Never open the housing while the equipment is in operation or connected to hospital grade outlet as it may result in electric shock.

Maintenance, Modifications, and Repairs

Fukuda Denshi is liable for the safety, reliability, and performance of its equipment only if;

- Maintenance, modifications, and repairs are carried out by authorized personnel.
- Components are used in accordance with Fukuda Denshi operating instructions.

A full technical description of the DS-7200 system is available from your local Fukuda Denshi representative.

Precautions about the Pacemaker

₩ARNING	 Minute ventilation rate-adaptive implantable pacemakers can occasionally interact with certain cardiac monitoring and diagnostic equipment, causing the pacemakers to pace at their maximum programmed rate. The cardiac monitoring and diagnostic equipment may possibly send wrong information. If such event occurs, please disconnect the cardiac monitoring and diagnostic equipment, or follow the procedures described in the operation manual of the pacemaker. (For more details, contact FUKUDA DENSHI personnel, your institution's professionals, or your pacemaker distributors.) Reference "Minute Ventilation Rate-Adaptive Pacemakers" FDA alerts health professionals that minute ventilation rate-adaptive implantable pacemakers can occasionally interact with certain cardiac monitoring and diagnostic equipment, causing pacemakers to pace at their maximum programmed rate. [Based on a safety bulletin issued by FDA Center for Devices and Radiological Health on October 14, 1998] Rate meters may continue to count the pacemaker rate during occurrences of cardiac arrest or some arrhythmias. Do not rely entirely upon rate meter alarms. Keep pacemaker patients under close surveillance. See this manual for disclosure of the pacemaker pulse

Non-Explosion Proof

▲ DANGER	Never operate the equipment in the presence of flammable anesthetics, high concentration of oxygen, or inside hyperbaric chamber. Also, do not operate the equipment in an environment in which there is a risk of explosion. Explosion or fire may result.
-----------------	--

Defibrillation Safety

▲WARNING	 When defibrillating, keep away from the electrodes or medicament applied to the patient chest. If this is not possible, remove the electrodes or medicament before defibrillating. If the defibrillator paddles directly contact the electrodes or medicament, electrical shock may result by the discharged energy.
	 When defibrillating, make sure that the electrodes, sensor cables, or relay cables are firmly connected to the device. Contacting the metal part of the disconnected cable may result in electrical shock by the discharged energy.
	 When defibrillating, do not touch the patient and the metal part of the device or cables. Electric shock may result by the discharged energy. This equipment will return to standard operating mode within 10 seconds. The stored data will not be affected. The measurement accuracy will temporarily decrease during defibrillation, but it will not compromise the safety of patient and the equipment.

Electrosurgery Safety

[
▲WARNING	The monitoring system contains protection against interference generated by electrosurgical instruments. However, operating conditions, surgery site with respect to the location of ECG electrodes, or the type of instrument used, may cause noise on the ECG. The noise is generated at the tip of an electrical knife and is difficult to completely eliminate because of the frequency components of the ECG. To reduce electrosurgical interference, take the following precautions: Location Locate the electrosurgical unit as far as possible from this unit and the patient cable. This will help reduce interference on the ECG through the monitor or cables. Power Supply Connect the electrosurgical unit to a power supply that is different from that of the monitor. This will help prevent interference through the power cable. Electrode Placement The amount of interference is considerably different depending on the electrode position and surgery site. Place the ECG electrodes as far away as possible from the surgery site and the ground plate. If the electrodes are placed in this path, the amount of interference will be quite large. Position (+) and (-) electrodes as close as possible to each other. Ground Plate When using electrosurgical instruments, make sure the contact between the patient and the ground plate is secure. If the connection is incomplete, the patient and the ground plate is secure.

Precautions about Magnetic Resonance Imaging

	MARNING	The local heating caused by the induced electromotive force may cause burn injury to the patient (subject). For details, refer to the operation
--	---------	---

Precautions about Connections to Peripheral Devices

In the interest of safe and sufficient performance of this equipment, the connection of other manufacturers' equipment to the monitor is not authorized, unless the connection is explicitly approved by Fukuda Denshi. It is the user's responsibility to contact Fukuda Denshi to determine the compatibility and warranty status of any connection made to another manufacturer's equipment.

Awarning	For the connector with \bigwedge mark, only the peripheral devices specified by Fukuda Denshi should be connected with the given procedure. Use of an unspecified device may cause electric shock to the patient and/or operator due to excessive leakage current.
▲ CAUTION	All the peripheral device connectors on the DS-7200 system are isolated from the power supply. The connecting peripheral devices should comply with IEC 60601-1 or should be isolated with the isolation transformer in compliance with IEC 60601-1. To prevent danger of electric shock, always position the peripheral devices away from the patient.

When connecting peripheral devices to DS-7200 system, it is the user's responsibility to verify that the overall system complies with IEC 60601-1-1, "Collateral Standard: Safety Requirements for Medical Electrical Systems".

Precautions about the Fuse

▲ DANGER	If the fuse blows, contact Fukuda Denshi Service Representative. Do not continue using it as internal damage to the equipment may be considered.
-----------------	--

Accessories and Optional Accessories

AWRNING	Use only the cables specified by Fukuda Denshi.
	Not only the DS-7200 cannot deliver its maximum performance but may also result in increase in emission or decrease in immunity.

Precautions about the DS-7200 System

A DANGER	When connecting to other device, contact Fukuda Denshi service representative.
MOER	Danger such as electric shock may result to the patient and operator.
1	
▲ WARNING	 The DS-7200 system is not a life-support equipment. The DS-7200 system is not intended for use during patient transport outside a healthcare facility, and is not considered as mobile equipment. Do not connect unit or cable not authorized by Fukuda Denshi to any I/O connector. If done so by mistake, the DS-7200 system cannot deliver its maximum performance and the connected units may be damaged, resulting in a safety hazard. If the DS-7200 system is used under an environment not fulfilling the specified condition, not only that the equipment cannot deliver its maximum performance, the equipment may be damaged and safety cannot be ensured. If using the equipment under condition other than specified, contact our service representative. Use only the supplied 3-way AC power cable. Use of other cables may result in electric shock to the patient and the operator. The power cable must be connected to the hospital grade outlet. When using multiple ME equipment simultaneously, perform equipotential grounding to prevent potential difference between the equipment. Even a small potential difference may result in electric shock to the patient and the operator. The setup for the alarm mode and display mode remains stored even when the power is turned off or when discharging procedure is performed. Before monitoring, make sure the current monitoring mode is suitable for the patient's condition. The patient classification selection influences the precision of the QRS detection and NIBP measurement. Make sure the correct selection is made. If the QRS pace mask function is set to <u>DFFF</u>. [10ms], or <u>20ms</u>], the pace pulse may not generate due to incorrect HR (counting pace pulse as <u>QRS complex</u>. Select <u>DFFF</u>. [10ms] or <u>20ms</u>] only if you are sure that pacing failure will not occur, or when the patient can be constantly monitored. Be cautious when setting the "SpO₂ alveraging" duration as the SpO₂ alarm is based on th

	 When selecting Silence, Time Disp. Only or OFF (Alarm Pole) for the night mode, pay attention not to miss any important alarm by simultaneously monitoring the bed on other monitors such as central monitor. For the alarm mode, it is recommended to program the alarm mode in rough classification such as patient's age, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient. The RR/APNEA alarm will not be generated unless the parameter key corresponded to the selected RR/APNEA source is displayed. Be sure to display the parameter key for the RR/APNEA source. When lifting this device, hold the handle of the main unit. The "QRS SYNC" signal (No. 1) of the Status II connector is a delay output. (delay: 30 to 75msec, signal width: 100msec). Do not use it as a synchronizing signal for the defibrillator. Make sure the delay time of QRS SYNC signal fulfills the specifications of the connected device. Analog signal is a delay output. (about 35ms for ECG, BP) When connecting to a device using vital signs as trigger signals (ex. IABP), make sure the delay time fulfills the specifications of the connected device. The delay time fulfills the specifications of the connected device. The delay time may differ depending on the waveform shape or artifact interference.
▲ WARNING	 surgery), and if necessary, perform unique setup for each patient. The RR/APNEA alarm will not be generated unless the parameter key corresponded to the selected RR/APNEA source is displayed. Be sure to display the parameter key for the RR/APNEA source. When lifting this device, hold the handle of the main unit. The "QRS SYNC" signal (No. 1) of the Status II connector is a delay output. (delay: 30 to 75msec, signal width: 100msec). Do not use it as a synchronizing signal for the defibrillator. Make sure the delay time of QRS SYNC signal fulfills the specifications of the connected device. Analog signal is a delay output. (about 35ms for ECG, BP) When connecting to a device using vital signs as trigger signals (ex. IABP), make sure the delay time fulfills the specifications of the connected device. The delay time may differ depending on the waveform shape or artifact interference.

CAUTION	 Systems This equipment is intended to be used for only one patient. The installation of this equipment and its option unit should be performed by our service representative or a person who is well acquainted with this equipment. The internal switch setting will be performed by our service representative. Users should not open the maintenance cover. PC Card Slot will be used by our service representative for maintenance purpose. Users should not use it. The software upgrading will be performed by our service representative. The users should not tattempt it. Use only the accessories specified for this device. Otherwise, proper function cannot be executed. Do not reuse a disposable product. For quality improvement, specifications are subject to change without prior notice. When the product is used in regions whose voltage is other than 110-120V, a cable appropriate to the regulations and voltage of the country in which the product is being used shall be used. The display panel utilizes exclusive fluorescent light for the backlight. Since this fluorescent light or other hard-edged instruments. It may cause maffunction or damage the touch panel. In addition, do not apply pressure to any pant of the panel for a prolonged time. Do not use the touch panel with fingers or a touch panel pen. Do not touch with a pen-point or damage the young impact may cause damage. Pay attention not to hit or dorp the touch panel. Do not use the touch panel with the film or adhesive tape attached. Malfunction of the touch panel with the prolenged time. Do not press the touch panel with strength or twis your finger on the panel. It may cause damage. Pay attention not to hit or dorp the touch panel. Do not use the touch panel with the film or adhesive tape attached. Malfunction of the touch panel with strength or twis your finger on the panel. It may cause malfunction or damage the touch panel. Do

	ا ا ا ا ا ا ا ا ا ا ا ا ا ا ا ا ا ا ا
	etc.
▲ CAUTION	 Replace the electrode if the skin contact gets loose due to perspiration, etc. When an electrode is attached to the same location for a long period, some patients may develop skin irritation. Check the patient's skin condition periodically and change the electrode site as required. For stable arrhythmia detection and ECG monitoring, verify proper electrode placement, lead, waveform size, and filter mode selection. If not properly selected, it may cause erroneous detection. Always use the same type of electrodes. If different types of electrodes are used at the same time, the difference between the polarization potential from each electrode may interfere monitoring. The threshold level for arrhythmia detection and QRS detection changes with ECG waveform size. Set a proper waveform size for monitoring. When the waveform size is ×1/4, ×1/2, or ×1, the detection threshold is 250µV. When the waveform size is ×2 or ×4, the detection threshold is 150µV. When arrhythmia is present, HR measurement accuracy may be degraded. Select the appropriate lead for ECG1, 2 to be used for arrhythmia detection, telemeter, central monitor transmission, and recording. The selected lead for ECG1, 2 will be used for recall waveform and recording waveform as well as for arrhythmia analysis. The QRS detection leads, arrhythmia detection leads, monitoring leads on the central monitor, recording leads are fixed as ECG1 and ECG2. Especially for arrhythmia detection, set the most appropriate leads with high QRS for ECG1 and ECG2. Automatic size/position of the ECG is effective only at the time the <u>AUTO</u> key is pressed. This does not continually adjust size and position. There are some cases when the pacemaker pulse can not be detected depending on the pacemaker type, pulse voltage, pulse width, electrode lead type (unipolar, bipolar), or electrode placement which causes the pacemaker pulse amplitude to decrease and disables the pacemak
	AUTO key is pressed. This does not continually adjust size and
ZICAUTION	
	noise and EMG, but it may also reduce the QRS amplitude. The ESIS
	causes the pacemaker pulse amplitude to decrease and disables the
	If signals similar to a pacemaker pulse are present, such as electric
	blanket noise or excessive AC frequency noise, these may be erroneously detected and displayed as a pacemaker pulse.
	 When a spontaneous QRS and pacemaker pulse overlap (as in a fusion beat), QRS detection will be suspended and the heart rate will be
	 reduced. If a pacemaker pulse is continuously detected due to AC frequency interference, QRS detection will be suspended and the heart rate will be reduced. Also arrhythmia detection will not be possible.
	 Respiration Monitoring When the following relay cables are used, respiration cannot be measured.
	 measured. Relay Cable CI-700E-3 (FA) (defibrillation and electrosurgery-proof, 3-electrode)
	•Relay Cable CI-700E-4 (FA) (defibrillation and electrosurgery-proof,
	4-electrode) • Relay Cable CI-700E-5 (FA) (defibrillation and electrosurgery-proof,
	 5-electrode) When a defibrillator is used during respiration monitoring, a large offset voltage will be placed on the ECG electrodes, which may cause
	interruption of monitoring for a few seconds.

CAUTION	 For Masimo[®] sensor, change the sensor attachment site every 4 hours for the reusable sensor, and every 8 hours for the disposable sensor. Exercise extreme caution with poorly perfused patients; skin erosion and pressure necrosis can be caused when the sensor is not frequently moved. Assess site at least every 2 hours with poorly perfused patients. The SpO₂ patient cables (PC04, PC08, and PC12) are intended for Masimo SET sensors only. Connect them only to DS-7210M. If connected to other device, it will not function properly. Measuring on a limb with NIBP cuff, arterial catheter, or intracatheter may result in incorrect measurement. For additional warnings, cautions or contraindications when using sensors with DS-7210 Nellcor[®] model or DS-7210M Masimo[®] model, refer to each SpO₂ sensor instruction manual. If SpO₂measurement failure occurs due to the reason such as sensor detachment from the patient, SpO₂measurement data will be displayed as "". Be cautious as numeric data alarm will not generate in such case. Precautions for DS-7210M Masimo[®] Model The measurable pulse rate range is 25 to 240bpm. "xxx" will be displayed if 25bpm and below or 240bpm and above is measured. If <u>OFF</u> is selected for "PI Display" under the SpO₂ configuration setup, "SpO₂ Low Perfusion" alarm will be indicated by message display only. The alarm sound will not be generated. NIBP Monitoring Select the appropriate cuff size which best fits the arm circumference. If the cuff size is inappropriate, it may cause measurement error. Do not use a cuff which is worn out. The cuff may burst during inflation. Do not use at the the connection is secure. Correct NIBP measurement cannot be performed. Make sure that the connection is secure. Correct NIBP measurement cannot be performed if artificial heart lung machine is used or if the pulse is difficult to detect. Pay attention when measuring the
	 disorders or hyper coagulation. The cuff inflation may cause petechia or circulatory failure by the blood clot. Do not apply the cuff to the arm or thigh where vein is secured. The blood may backflow causing the chemical injection to cease. Pay attention not to bend the cuff hose. Check the condition of cuff-applied part on the patient during
	 measurement so that the blood circulation will not be blocked over long period of time by the squashed or bent cuff hose. Check the patient's condition constantly while measuring over long period of time with interval of 2.5 minutes or less. Also, periodically check the blood circulation while performing periodic measurement over long period of time. Congestion may occur at the measuring site. The following factors may affect the NIBP value. Body motion, arrhythmia, convulsion Continuous noise such as cardiac massage Periodic electromagnetic noise

r	
CAUTION	 When a PTG (SpO₂) sensor is applied to the toe or forehead, the Dyna Alert may not function depending on the patient's condition. When using the Dyna Alert function, be aware of these risks and do not increase the NIBP interval time by relying only on the Dyna Alert function. After the Dyna Alert NIBP measurement, the next Dyna Alert NIBP measurement cannot be performed for 2.5 minutes. The Dyna Alert will not properly function for the following cases. If peripheral circulatory insufficiency or very low BP is developed. If highly-frequent arrhythmia is generated. If a large noise from body movement or electric surgery equipment is interfering. If autonomic nerve or circulatory dynamics is largely affected by medication. For the following situation, measurements will be terminated. When the inflation value has exceeded 300mmHg for adult, 210mmHg for child, and 150mmHg for neonate. If used with the incorrect patient classification, it will not only cause erroneous measurement, but the inflating level for the adult may be applied to child or neonate causing dangerous situation to the patient. The 1-minute interval measurement will automatically stop after 12 minutes. If the mean BP display is set to OFF. BP Monitoring Do not reuse disposable product for BP measurement. When the main power is turned OFF, the previous zero balance information will be maintained, and the BP value will not be displayed or the tabular trend or the NIBP later the power is turned OFF, the previous zero balance is performed. However, if the status is not displayed will be maintained, and the BP value will be displayed or the tabular trend or the NIBP is function if the display is set to OFF. BP Monitoring Do not reuse disposable product for BP measurement. When the main power is turned OFF, the previous zero balance information will be maintained, and the BP value will
	 When the position of the heart has changed due to body movement.
	 When the position of the transducer has changed. When measuring for a long period of time and there is a possibility of
	 measurement error due to change in ambient temperature, etc. When the connector is connected / disconnected, or transducer is
	replaced.
	 When the power has been turned OFF for more than 5 minutes.

CAUTION	 Note that the Systolic Pressure (SYS) = Peak Systolic Pressure (PSP) for the graphic trend, data base, and alarm setup. When ECG is not measured, PDP cannot be calculated. The undisplayed BP data (SYS/DIA/Mean) will not generate a BP alarm or be displayed in the tabular trend. Select the appropriate display type according to the monitoring purpose. CO₂ Monitoring (MGU-722) All FilterLine[®] sampling products are for single patient use only. Perform calibration after Initialization Time (max. 180 seconds has elapsed since the power is turned ON. Do not disconnect the sampling tube during calibration. If disconnected, calibration will cease when the sampling tube is disconnected. Conduct CO₂ calibration for the following case. If the CO₂ gas calibration is not performed at a specified interval, CO₂ measurement accuracy may be affected and also subsequent gas calibration may not be possible. For the following case, a message, "Calibrate the CO₂ unit (MGU-722) is approaching" will be displayed at power ON. Conduct CO₂ calibration. When the accumulated measurement time exceeds 1200 hours from first use. When tBCO₂ measurement is not stable or accuracy is degraded compared with other measuring device conduct CO₂ calibration. CO₂ Monitoring (MGU-721 with CAPNOSTAT[®] 5 CO₂ sensor) The disposable airway adapter should be opened just before use. Do not sterilize the airway adapter should be appended, make sure to set the alarm condition for that unit. The alarm setup is necessary for each measurement unit. Alarm The alarm priority is high for level 1 (life threatening alarm), medium for level 2 (acutionary alarm), and low for level 3 (treatment meeded alarm). Alarm messages for the arrhythmia alarm will continue to be displayed according to the priority. (Level 1 → Level 2 → Level 3 → Level 4)<!--</th-->

CAUTION	 Regardless of ON/OFF setting of "Suspend Arrhy. Analysis during Interference" under Hospital Setup (Preset Menu), the "Cannot analyze" alarm will generate when analysis is suspended for more than 30 seconds. The measurement range and alarm range differs for the following parameters. Be cautious not to set the alarm limit outside the measurement Range: 25 to 240bpm (If 25bpm and below or 240bpm and above is measured, "xxx" will be displayed.) Alarm Range: 20 to 300bpm NIBP Measurement Range: 10 to 280mmHg Alarm Range: 10 to 280mmHg CO₂ for MGU-722 (Microstream[®] CO₂ Unit) Measurement Range: 10 to 15.0kPa CO₂ for MGU-722 (Microstream[®] CO₂ Unit) Measurement Range: 10 to 15.0kPa For the SpO₂ measurement, whether to use the SEC alarm function and its threshold selection should be based on the patient's clinical indication portent and medical evaluation. (For Nellcor[®] SpO₂ unit) If the SpO₂ alarm and SEC alarm setup is set to OFF, the SEC alarm integral value will be set to 0. (For Nellcor[®] SpO₂ unit) The alarm mute ON/OFF setup will remain effective even when the power is turned OFF. Be cautious not to miss any important alarm by leaving the alarm silenced. Pay attention not to set the alarm volume too low to avoid missing any important alarms. System Configuration When the waveform and numeric data display for each parameter is set to OFF, the pulse rate derived from SPO₄ will not be displayed either. When the waveform and numeric data labeled as BP1 or ART is set to OFF, the pulse rate derived from SPO₄ will not be displayed either. When the waveform and numeric data display for SpO₂ is set to OFF, the pulse rate derived from SPO₄ will not be displayed either. When the waveform and numeric data display for SpO₂ is set to OFF, the nulse rate derived from SPO₄ will not be displayed either.
	the same floor. Otherwise, it may cause to remote control more than one monitors at the same time.After the remote control setup, check that the remote control unit is properly operating.
	If you start monitoring a new patient without performing a discharge

▲ CAUTION	 ST Measurement For the lead which the electrode is detached, the reference waveform setup cannot be performed. Check if the electrode is correctly attached, and perform the setup again. CF Card Use only the specified CF card.
	 Use only the CF card formatted with this device. Restart the system after reading the setup data from the CF card. The setup data will become effective after the system is restarted. Reading the patient data from the CF card will erase all previous patient data stored in the patient monitor.
	 TCP/IP Network Connection After setting the IP address, etc. for the laser printer, make sure to turn OFF and back ON the power of the printer. Maintenance
	 The maintenance procedure will be performed by our service representative. Users should not attempt this procedure as malfunction may result to the device.
	 If stains cannot be removed from the touch panel surface, wipe softly with a dry or ethanol dampened cleaning cloth. Never use strong-acidic cleaning solution. (Neither is it recommended that mild acidic or alkaline cleaning solution to be used.)
	• A special coating is applied to the surface of the touch panel. Do not wipe the surface with a cloth or gauze with coarse texture. Wipe the surface with a soft cleaning cloth provided as optional accessory or with an eyeglass cleaning cloth.
	 Clean the equipment frequently so stains can be removed easily. To prevent injury, it is recommended to wear gloves when cleaning the equipment.
	 Do not allow liquids such as alcohol or cleaning solution enter the equipment or connectors.
	 Do not use organic solvents, thinner, toluene and benzene to avoid damaging the resin case.
	 Do not polish the housing with abrasive or chemical cleaner. When sterilizing the entire room using a spray solution, pay close attention not to have liquids get into the equipment or connectors. Use only neutral detergent to clean the housing. Do not use chemical
	cloth, scrub brush, abrasive, polishing powder, hot water, volatile solvent and chemicals (cleanser, thinner, toluene, benzine, benzol, and synthetic detergent for house and furniture), or sharp-edged tools. The surface resin coating may be damaged, resulting in discoloration, scratches, and
	other problems.
	Do not open the housing.If you accidentally wet the device, dry it completely and verify it operates
	safely before usage.
	 If the patient monitor was stored for some while, leave the monitor at the operating environment (10 to 40°C, 30 to 85%) before usage.
	 Replace the components periodically as specified.

Precautions about the Wired Network System (DS-LAN II/DS-LANIII)

A warning	 Do not connect unspecified device to a wired network. Do not mix devices with DS-LANII and DS-LANIII setting in the same wired network. The network may cease and proper monitoring may not be possible. Before setting the bed ID, make sure that the DS-LAN
	(DS-LANII/DS-LANIII) is correctly set on the Monitor Setup menu. If not correctly set, the network may cease which may lead to accidents such as not transmitting life threatening alarms to the central monitor.
CAUTION	 When connecting to the DS-LAN network, perform "DS-LAN Setup" in the Monitor Setup menu and restart the system before connecting the LAN cable. If performing wired network transmission, configure the displays so that the numeric data corresponded to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted. The Bed ID is factory set to 000. If connected to the wired network with the ID unchanged, monitoring on the central monitor will not be possible. When connected to the wired network, make sure that there are no other bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID. If there are more than one bedside monitors with the same Bed ID to 0048 For DS-LANIII network: 001 to 048 For DS-LANIII network: 001 to 100 As the DS-7200 does not have the arrhythmia template display and 12-lead ST display function, these displays on the central monitor will not be corresponded. If connected to a wired network, time/date will be the same with the central monitor. Even if the time/date is changed on the DS-7200 system, it will be corrected to the time/date of the central monitor. On some central monitor depending on the model type or software version, the setups for "HR Low Limit for VT" and "HR Low Limit for Run" cannot be performed. On a wired network, the alarm generated on the DS-7200 will be output to the network with a maximum delay of 1 second, and to the central monitor with a total delay of 2.5 seconds. In case of DS-LANII network, if the HR/PR source is EP (Or, if Auto selects BP for HR/PR source), the EC
	Redrawing the ST display will return the display to normal.

▲ CAUTION	 On the central monitor, the respiration waveform and RR value based on the RR/APNEA alarm source selected on the DS-7200 will be displayed. The RR and APNEA monitored on the central monitor and the DS-7200 will be the same. If the measurement unit of CO₂ concentration is "mmHg", and <u>99mmHg</u> is selected for "CO₂ (mmHg) Upper Limit for LAN, Telemetry" on the monitor setup menu, the CO₂ value of 100mmHg or above will be transmitted as 99mmHg. There are following restrictions when connecting the DS-7200 system to the DS-LANII network. Make sure that the "DS-LAN Setup" on all the bedside monitors and central monitors are set to <u>DS-LANII</u> before connecting the monitors to the network. When DS-5800N/NX/NX^{MB} is used as a central monitor, recall, graphic trend, and tabular trend will not be displayed. Also, Σ recording cannot be performed. For the ST display, overlap waveform will not be displayed on the CS-5800N/X/NX/M^{MB} until 15 minutes have passed since the reference waveform is set on the DS-7200. If the measurement unit for BP (mmHg/kPa) is different between the bedside monitor and the central monitor. If a central monitor which does not support the "kPa" measurement unit is used, and the measurement unit on the bedside monitor is set to "KPa", BP waveform/numeric data, NIBP dat, NIBP dist, etc. in "kPa" unit will be treated as not measured data and will not be displayed on the central monitor. If a central monitor be performed. Arrhythmia alarm of TACHY, BRADY, COUPLET, PAUSE, TRIGEMINY will not be transmitted. Arrhythmia alarm of "SLOW_VT" will be transmitted as "VT". On a wired network, waveform, numeric data, alarm of TEMP3 will not be displayed as "x×x", maximum or minimum value of measurable range will be transmitted. Arrhythmia alarm of "SLOW_VT" will be treated as not measured data. Arrhythmia alarm of "BLOW_VT" will be treated a
	generated on the bedside monitor.
	 If using a HUB for the DS-LANII network construction, make sure to
	use a repeater HUB recommended by Fukuda Denshi.
	use a repeater non recommended by Fukuda Densm.

There are following restrictions when connecting the DS 7200 system to
 There are following restrictions when connecting the DS-7200 system to the DS-LANIII network.
• In order to connect to the DS-LANIII network, the software version
needs to be the version which supports the DS-LANIII. For details,
 refer to our service representative. Make sure that the "DS-LAN Setup" on all the bedside monitors and central monitors are set to DS-LANIII before connecting the monitors to the network. If the measurement unit for BP (mmHg/kPa) and temperature (°C/°F) is different between the bedside monitor and the central monitor, the corresponding waveform and numeric data will not be displayed on the central monitor. If using a HUB for the DS-LANIII network construction, make sure to use a switching HUB recommended by Fukuda Denshi. The displayable waveform, numeric data, alarm will differ depending on the central monitor. There are following restrictions when recording the DS-7200 data on the central monitor recorder. The DS-7200 can not perform the recording with the AU-5500N recorder. If the measurement unit of BP is kPa, the BP waveform, BP numeric data, and NIBP numeric data will be treated as not measured data. If the measurement unit of temperature is °F, the temperature data will be treated as not measured data. When a parameter is not measured, the waveform for that parameter will not be recorded, and measurement data will be recorded as "" or blank. The waveform recording and graphic trend recording, some parameters may not be able to be recorded depending on the central recorder. For the waveform recording and graphic trend recording, some parameters may not be able to be recorded depending on the cantral recorder. For the Mayeform source is BP (Or, if Auto selects BP for HR/PR source), ECG will not be recorded on the central recorder. If the RR/APNEA alarm source is other than [Impedance] (Or, if Auto selects other than impedance for RR/APNEA alarm source), respiration waveform will not be custure to the corder.
selects other than CO_2 for RR/APNEA alarm source), CO_2 waveform
will not be output on the central recorder.
When graphic trend recording, tabular trend recording, or NIBP list
recording is output on the central monitor recorder from the DS-7200, HR measurement value from ECG will be recorded for the HR value
and ST trend.

Precautions about the Wireless Network System

▲ DANGER	When monitoring a patient using wireless telemetry, make sure the patient data is properly received at the central monitor. Pay special attention when the channel ID at the bedside monitor is changed.
▲WARNING	 A password can be set to access the channel ID setup menu to allow only the telemetry channel administrator to change the channel ID. Some wireless combinations of telemetry transmitters may generate interference with other devices. Before selecting a channel, verify that it will not interfere with other channels. Make sure the telemetry manager of your system is aware of any changes to the telemetry channels. If transmitters are used in a neighboring medical facility, your facility and the neighboring facility must make agreements on the setting of the telemetry channels to prevent telemetry interference.
▲CAUTION	 On a wireless network, the alarm generated on the DS-7200 will be transmitted to the central monitor with 15 seconds delay. If the BP unit is kPa and temperature unit is °F, the measurement value will be converted to mmHg and °C respectively when transmitting to the central monitor. If kPa/°F is used as the unit on the central monitor, the measurement value will be reconverted to kPa/°F. If performing telemetry transmission, configure the display so that the numeric data corresponded to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted. The setup of channel ID and group ID should be performed only by our service representative. Users should not perform this procedure as malfunction of the equipment may occur. BP waveform with a scale above the programmed scale can not be properly transmitted. When transmitting the BP waveform, check the displayed BP waveform scale. If the measurement unit of CO₂ concentration is "mmHg", and <u>99mmHg</u> is selected for "CO₂ (mmHg) Upper Limit for LAN, Telemetry" on the monitor setup menu, the CO₂ value of 100mmHg or above will be transmitted as 99mmHg.

Precautions for Use of the Bidirectional Wireless Communications (TCON)

▲ CAUTION	 When using the TCON system, pay attention to the following. The medical institution (hereinafter referred to as "Institution") must execute investigation required to prevent interference including types of radio waves, frequencies, and antenna power if wireless equipment is already installed and being used in the facility. Even if this device is installed within the range of radio communication, the communication may not be possible due to noise or multi-path phasing etc. Always consider this thoroughly before use. Do not install this device in an area where it will be subject to splashing water. Water entering the equipment may cause the equipment to malfunction or be damaged. In managing the TCON system, make sure to follow the precautions below. The Institution should appoint a person (hereinafter referred as the "Overall Manager") to manage the wireless devices for the whole facility. When installing the TCON, the Overall Manager has to receive an explanation of the precautions for use of the TCON from the manufacturer or sales representative. The Overall Manager is responsible for the maintenance and storage of the equipment. The Overall Manager should create a management log (hereinafter referred to as the "log"), which contains a list of the management status of the wireless channels, register it in the log, and give proper instructions to the TCON user. The user needs to verify the transmitting/receiving operation before use. If interference or breakdown occurs in the communication, the TCON user is required to stop using the TCON and to inform the Overall Manager of the problem. The Overall Manager is to deal with the problem properly and/or contacts the nearest Fukuda Denshi representative for service.
	 The Bidirectional Wireless Communications Module (TCON) uses radio waves to transmit data. Therefore, necessary precautions need to be taken for the characteristics and difficulties of using the device that emits radio waves. The TCON user should fully understand these precautions beforehand, and use the TCON device safely. Furthermore, situations in which interference may occur are outlined below. In such cases, pay special attention to the condition of the patient connected to the bedside monitor, and eliminate the cause of interference. 1. The patient's data may become mixed with a different patient's data due to interference. When there are multiple TCON communication devices set to the same TCON ID and channel (group). 2. When symptoms such as being unable to communicate, unstable communication, or poor reception may occur. When the radio communication is bad because there are metal, concrete, or other such obstacles between the Bidirectional Wireless Communications Modules (TCON). When a different wireless device is using the same frequency (channel). When there are other TCON devices nearby using different channels (groups). When a cell telephone or other wireless device is being used nearby. When citizens broadcast bands such as amateur radio or truck radios are used in the vicinity of the TCON operating area.

CAUTION	 When a computer or word processor, or electrical device that has an internal computer, is used near the TCON device antenna. When the TCON device is installed or moved to a location that is outside the radio communication range. If a nearby different TCON group is set with a TCON channel frequency that is too close to the channel frequency set for the current TCON group. Follow the instructions of the Overall Manager for the wireless channel when setting the TCON ID or channel (group) to prevent interference within the same institution. For the TCON ON/OFF setup, if the "OFF" is selected, the message such as "Check TCON Comm." will not be displayed. Check that the TCON radio wave strength between the central monitor and bedside monitor is sufficient. Make sure that "Tu" mark is displayed. Check that the TCON channel (Group) is the same for the bedside monitor and the central monitor in the same TCON group. Do not move the TCON device during operation. Otherwise, symptoms such as being unable to communicate, unstable communication, or poor reception may occur. There are following restrictions when connecting the DS-7200 system to the TCON Network. If the measurement unit for temperature is "°F", the central monitor can not receive the measurement data for temperature. In addition, the alarm settings for temperature can not be operated from the central monitor. If the measurement unit for BP is "kPa", the central monitor can not receive the measurement data for NIBP, BP1, and BP2. In addition, the alarm settings for NIBP, BP1, and BP2. In addition, the alarm settings for NIBP, BP1, and BP2 can not be operated from the central monitor. The NIBP measurement cannot be started from the central monitor via TCON system if the NIBP measurement interval is set to <u>2 min</u> / <u>2.5 min</u> / <u>3 min</u> / <u>5 min</u> or during the 1-minute measurement. However, it can be stopped.
	100mmHg or above will be transmitted as 99mmHg even within measurement range.

Precautions about the Ventilator Monitoring

	• The ventilator alarm on this monitor should be used as supplementary
	function. Check the patient's condition, ventilator alarm sound and
	message occasionally.
	• The ventilator alarm sound is set to OFF at factory default setting.
	The alarm sound can be turned ON on the volume setup menu.
	 If the DS-7200 system does not generate an alarm even though the ventilator is generating an alarm, or if any other malfunction occurs,
	immediately check the ventilator, DS-7200 system, cable, and replace the
	cable if necessary. If the malfunction persists, stop using the device.
	• After connecting the ventilator and the DS-7200, ensure that "Vent.
	Online" message is displayed for the connection status. Otherwise, the
	DS-7200 will not detect the ventilator alarm.
	 The alarm generation on the DS-7200 system is not assured if the alarm other than specified generates at the ventilator.
	See For details of the specified alarms, refer to Δ WARNING on "2. Basic Operation
	Ventilator Alarm Input".
	• The Evita2dura / Evita4 / EvitaXL / Savina acquires alarm information from the serial port. The ventilator alarm that cannot be acquired from the
	serial port is not guaranteed.
	For corresponding alarm, refer to the service representative of the
≜ WARNING	ventilator manufacturer.
	• The DS-7200 system will not correspond to the following alarms
	generated on the Evita 4 / Evita XL / Evita 2 dura.
	• O ₂ monitoring disabled alarm, CO ₂ alarm disabled alarm, Oximeter
	alarm disabled alarm, Neo. volume measurement inoperable alarm, Minute volume alarm disabled alarm, Minute volume alarm low off
	alarm, Tidal volume alarm high off alarm, Apnea alarm off alarm,
	Nebulizer active alarm
	 There is a communication delay of 3 seconds between the DS-7200
	system and the Evita ventilator. Therefore, if the alarm generated at the
	ventilator is resolved within 3 seconds, the ventilator alarm may not be generated at the DS-7200 system.
	 The DS-7200 system will not correspond to the following alarms
	generated on the Savina.
	• O ₂ monitoring disabled alarm, Minute volume alarm disabled alarm,
	Minute volume alarm low off alarm, Tidal volume alarm high off alarm,
	Apnea alarm off alarm, Nebulizer active alarm
	• There is a communication delay of 3 seconds between the DS-7200 system and the Savina ventilator. Therefore, if the alarm generated at the
	ventilator is resolved within 3 seconds, the ventilator alarm may not be
	generated at the DS-7200 system.

	 The ventilator operation should be performed by well-trained and outborized personnel.
	 authorized personnel. For connecting the DS-7200 system and ventilator, use only the specified
	connection cable.
	• Verify that the DS-7200 system and the ventilator are properly connected.
	• When connecting the cable, verify that the main power of the DS-7200
	system and the ventilator is OFF.
	• For the SV-900, PB, Evita, and Savina ventilator alarm factor cannot be
	 transmitted to the central monitor. Depending on the central monitor type and software version, ventilator
	alarm factor may not be displayed. For details, refer to our service
	representative.
	 Check occasionally the communication status of the DS-7200 and the
	ventilator.
	• Verify that the ventilator alarm is not generated, and the "Vent. Online"
	message is displayed.
	• The "Check external alarm" will be displayed until the proper
	communication with the ventilator is resumed. When the communication is
	resumed, the screen will automatically return to the home display.
	• When disconnecting the ventilator and the DS-7200, make sure to select
	OFF on the "Check external alarm" display which appears when the power of the ventilator is turned OFF, or when the cable is disconnected.
	 When connecting the PURITAN-BENNETT ventilator, follow the
	precautions below.
	• The serial port (RS-232C) of the ventilator should be set as follows.
	Refer to the service representative of the ventilator manufacturer.
	Baud Rate : 9600bps
	Data Bit : 8bit
A CAUTION	Parity Bit : None (Stop Bit) : (1bit)
	 The DS-7200 system detects the "ventilator alarm" when the nurse
	call port on the ventilator outputs the alarm signal. For details of
	ventilator setup and alarm signal output condition from the nurse call
	port, refer to the service representative of the ventilator manufacturer.
	• When connecting the Evita2dura / Evita4 / Evita XL / Savina ventilator,
	the serial port (RS-232C) setup of the ventilator should be as follows.
	 Refer to the service representative of the ventilator manufacturer. For Evita 2 dura / Evita 4 / Evita XL
	Protocol : Medibus
	Baud Rate : 19200bps
	Data Bit : 8bit
	Parity Bit : Even
	Stop Bit : 1bit
	For Savina
	Protocol : Medibus
	Baud Rate : 9600bps Data Bit : 8bit
	Parity Bit : None
	Stop Bit : 1bit
	 For PURITAN-BENNETT ventilator, AWP and AWF waveform cannot be
	displayed or recorded. Only the numeric data will be displayed.
	• For SV-300 and Servo-i/s, P-V loop and F-V loop cannot be displayed or
	printed. In addition, Insp Resistance, Exp Resistance, Compliance value
	cannot be displayed or printed on the ventilator numeric data display.
	 For SV-900, P-V loop, F-V loop and numeric data cannot be displayed or printed. Only the alarms will be generated.
	 For PURITAN-BENNETT ventilator, P-V loop and F-V loop cannot be
	displayed or recorded. Only the numeric data will be displayed.
Precautions for Use of SpO₂ Sensor

≜ DANGER	Burn Risk in Using SpO ₂ Sensor In SpO ₂ monitoring, always use the sensor/relay cable specified by Fukuda Denshi. If any other sensor/relay cable is used, a high temperature rise of the sensor may place the patient in danger of burns. If there are any questions regarding the sensor/relay cable use for SpO ₂ measurements of this device, please contact Fukuda Denshi service representative.
-----------------	--

Precautions for Masimo[®] Model: DS-7210M

	No Implied License
A CAUTION	Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized sensors or cables which would, alone, or in combination with this device, fall within the scope of one or more of the patents relating to this device.

Precautions for Use of NIBP Cuff

A CAUTION	This product contains natural rubber latex which may cause allergic reactions. (FDA: Medical Alert on Latex Products, "Allergic Reactions to Latex-Containing Medical Devices", Food & Drug Administration, 9200 Corporate Blvd., Rockville, MD 20850, 1991.)
------------------	--

Disposing of Equipment, Accessories, or Components

MCAUTION When disposing of the equipment, accessories, or components, use an industrial waste distributor. Do not dispose of as ordinary waste.

Precautions about Transportation

For transporting the DS-7200 system, pack with specified packing materials.



Refer to "11. Technical Information Specification / Performance" for environmental condition during transportation.

Precautions about RTC or Data Backup

▲ CAUTION

Precautions for Use of Lithium-Ion Battery Pack

▲ DANGER	 This battery pack is intended for exclusive use with the DS-7200 system (or other specified equipment). Do not use with other equipment. Otherwise, the performance and life cycle of the battery pack deteriorates, and may cause leakage, heating, fuming, ignition, and explosion of the battery. Do not disassemble or remodel the battery pack. If the security apparatus or protector inside the battery pack gets damaged, it may cause leakage, heating, fuming, ignition, and explosion of the battery. Do not use the battery pack if leaked or transformed. If the security apparatus inside the battery pack is damaged, it may cause leakage, heating, fuming, ignition, and explosion of the battery. When installing the battery to the device, ensure the connector direction is correct. If installed in opposite direction, it may cause leakage, heating, fuming, ignition, and explosion. If the leaked solution of the battery gets into the eyes, do not rub the eyes. Wash thoroughly with clean water and immediately receive medical treatment from the doctor. If not treated soon, it may cause serious injury.
₩ARNING	 If the leaked solution of the battery gets on to the skin or clothes, immediately wash down with rinse water. If not treated soon, it may cause serious injury. If the charging operation does not complete within specified time, immediately remove the battery pack from the monitor and unplug the power cable. Otherwise, it may cause leakage or heating of the battery. Do not throw the battery pack into fire or apply heat. The insulator may melt, gas exhaust vent or security apparatus may get damaged, or electrolyte may ignite causing leakage, heating, fuming, ignition, and explosion of the battery. Do not connect the (+) and (-) terminals of the battery with a wire or any other metal. Also, do not carry or store the battery with any metal such as necklace, hairpins, etc. The battery may short causing excessive current flow which may result in leakage, heating, fuming, ignition, and explosion of the battery, or heating of the metal (wire, necklace, hairpin, etc.) Do not directly solder on to the battery pack. The heat may melt the insulator or damage the security apparatus which may result in leakage, heating, fuming, ignition, and explosion of the battery. Do not put the battery pack in microwave oven or a pressure cooker. If heated suddenly or if sealed condition breaks, it may result in leakage, heating, fuming, ignition, and explosion of the battery. Do not drive a nail in, hit with a hammer, step on the battery pack, or peel off or scratch the exterior tube. The battery may explode and transform causing short-circuit which may result in leakage, heating, fuming, ignition, and explosion of the battery.

₩ARNING	 Do not apply strong impact or throw the battery pack. This may result i leakage, heating, fuming, breakage, ignition, and explosion of the batter Also, if the security apparatus incorporated in the battery gets damaged battery charges with abnormal current and voltage, which results in leakage, heating, fuming, ignition, and explosion. Do not get the battery pack wet with water, sea water or chemicals. If the security apparatus incorporated in the battery gets damaged, it may result in leakage, heating, fuming, ignition, and explosion of the battery pack. Do not get the battery pack wet with water, sea water or chemicals. If the security apparatus incorporated in the battery gets damaged, it may result in leakage, heating, fuming, ignition, and explosion of the battery pack. Do not connect the battery pack directly to power outlet or cigarette her socket in a car. A high voltage application will cause excessive current and abnormal chemical reaction inside the battery. This may result in leakage, heating, fuming, ignition, and explosion of the battery. Do not use or leave the battery in a high temperature (80°C or over) seas near the fire or heater. If the resin separator gets damaged by heat, battery pack may become unusable, or may short causing leakage, heatfuming, ignition, and explosion. If the battery is leaking or generating an abnormal odor, immediately remove the battery away from the fire. The leaked electrolyte may cause heating, fuming, ignition, and explosion. 	
	 Do not peel off or scratch the exterior tube. 	
▲ CAUTION	 Do not use or leave the battery in high temperature. It may result in leakage or deterioration of the performance / life cycle of the battery. Immediately stop using the battery if any abnormality is found during use. Do not use / store the battery in reach of infants. If not using the device for a long period of time, turn OFF the power of the monitor and unplug the power cable. Otherwise, it may result in leakage of the battery pack. When disposing of the Lithium-Ion Battery Pack, use an industrial waste distributor. Do not dispose of as ordinary waste. Users should not attempt to install or replace the battery pack. For installation and replacement of the battery pack, contact our service representative. 	

To Prepare for Emergency Use

- 1. Accessories / Optional Accessories
 - (1) The ECG electrodes are consumable products. Always prepare extra supplies of electrodes.(2) Verify that there is no wire break on the patient cable. Check the operation once a week.
- 2. Battery Pack
 - (1) The battery self-discharges even when not in use. If there is any possibility to use the battery in emergency, the power cable should be always connected to the power receptacle. To fully charge the empty battery, it takes approximately 3 hours when the monitor is not operating, and approximately 10 hours when the monitor is operating.



Refer to "2. Basic Operation To Use with the Battery Pack"

(2) The performance of the battery deteriorates with repeated use. To maintain the initial performance, replace the battery at least once a year. It is recommended to indicate the start usage date on the battery so that the replacing date can be easily recognized.

Electromagnetic Compatibility

The performance of this device under electromagnetic environment complies with IEC 60601-1-2 (2007).

Precautions for Safe Operation under Electromagnetic Influence

AUTION	 If any sorts of electromagnetic wave, magnetic field, or static electricity exist around the device, noise interference or malfunction of the device may occur. If any unintended malfunction or noise occurs during monitoring, check the magnetic influence and take appropriate countermeasures. The following are examples of the common cause and countermeasures. Cellular Phone The radio wave may cause malfunction to the device. Cellular phones and radio sets should be turned off in the room (building) where medical device is located. Static Electricity In a dry environment (room), static electricity is likely to occur. Take the following countermeasures. Both operator and patient should remove any static electricity before entering the room. Lightning A lightning nearby may induce excessive voltage to the equipment. If any danger is suspected, use the uninterruptible power supply system. High frequency noise interference from other device through the power outlet Check where the noise is originated and remove it using filtering device, etc. Stop using the device that is originating the noise. Use other power outlet.

EMC Guidance

This equipment complies with IEC60601-1-2 (2007). However, if portable transmitter or wireless LAN equipment is used extremely nearby, the electromagnetic influence may largely exceed the compliance level and may cause unexpected phenomenon such as noise interference on the waveform, etc. Therefore, this equipment should be used in a location specified by each medical institution. If any unexpected noise interference on the waveform or failure to the peripheral device occurs, stop using the equipment and follow the instruction of the technician.

The following is the information relating to EMC (Electromagnetic Compatibility). (When using this equipment, verify that it is used within the environment specified below.)

Compliance to the Electromagnetic Emissions

The DS-7200 system is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-7200 system should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance	
RF Emissions CISPR 11	Group 1	The DS-7200 system uses RF energy only for its internal functioning of the equipment itself. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class A	This DS-7200 system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic Emissions IEC 61000-3-2	Not applicable		
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Not applicable		

•Compliance to the Electromagnetic Immunity (1)

The DS-7200 system is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-7200 system should assure that it is used in such an environment.

Immunity Test	IEC60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV: differential mode ±2kV: common mode	±1kV: differential mode ±2kV: common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_{T} (>95\% dip in U_{T}) for 0.5 cycle 40% U_{T} (60% dip in U_{T}) for 5 cycles 70% U_{T} (30% dip in U_{T}) for 25 cycles <5% U_{T} (>95% dip in U_{T}) for 5 sec.$		Mains power quality should be that of a typical commercial or hospital environment. If the user of the DS-7200 system requires continued operation during power mains interruptions, it is recommended that the DS-7200 system is equipped with an internal battery (option) or is powered from an uninterruptible power supply.
Power Frequency (50/60Hz) Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note U_T is the AC mains voltage prior to application of the test level.

Compliance to the Electromagnetic Immunity (2)

The DS-7200 system is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-7200 system should assure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the DS-7200 system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended Separation Distance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	$d = 1.2\sqrt{P}$
Radiated RF	3V/m	0)///	d = $1.2\sqrt{P}$ 80MHz to 800MHz
IEC 61000-4-3	80MHz to 2.5GHz	3V/m	d = $2.3\sqrt{P}$ 800MHz to 2.5GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^{a)} , should be less than the compliance level in each frequency range ^{b)} .
			Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1: At 80Mhz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DS-7200 system is used exceeds the applicable RF compliance level above, the DS-7200 system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DS-7200 system.
- ^b Over the frequency range 150kHz to 80MHz, field strength should be less than 3V/m.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the DS-7200 System

The DS-7200 system is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the DS-7200 system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DS-7200 system as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum	Separation Distance according to Frequency of Transmitter (m)				
Output Power of Transmitter (W)	$150 \text{ kHz to } 80 \text{ MHz}$ $d = 1.2 \sqrt{P}$	80MHz to 800MHz d = $1.2\sqrt{P}$	800MHz to 2.5GHz d = $2.3\sqrt{P}$		
0.01	0.12	0.12	0.23		
0.1	0.38	0.38	0.73		
1	1.2	1.2	2.3		
10	3.8	3.8	7.3		
100	12	12	23		

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1 : At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note 2 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Blank Page

Contents

Preface		
1. General Description	Describes the outline of this equipment.	1
2. Basic Operation	Describes the basic operation for monitoring.	2
3. Vital Application	Describes the procedure for vital application, etc.	3
4. Monitoring Setup	Describes the procedures to set the monitor according to the monitoring purpose.	4
5. Admit / Discharge of a Patient	Describes the procedure to admit or discharge a patient.	5
6. Parameter Setup	Describes the procedure to set the measurement condition, size, scale, etc. for each parameter.	6
7. Function	Describes about the functions such as arrhythmia analysis, trend, recall, etc.	7
8. System Configuration	Describes about the system configuration such as night mode, alarm mode, display mode, etc.	8
9. Installation	Describes about the environment for use, wireless system, etc.	9
10. Maintenance	Describes about the maintenance, troubleshooting of this equipment.	10
11. Technical Information	Lists the specification, default settings, pin assignments of external connector, etc.	11
12. Accessories	Lists the accessories and optional accessories for this equipment.	12
	 General Description Basic Operation Vital Application Vital Application Monitoring Setup Admit / Discharge of a Patient Parameter Setup Function System Configuration Installation Maintenance Technical Information 	1. General Description Describes the outline of this equipment. 2. Basic Operation Describes the basic operation for monitoring. 3. Vital Application Describes the procedure for vital application, etc. 4. Monitoring Setup Describes the procedure to set the monitor according to the monitoring purpose. 5. Admit / Discharge of a Patient Describes the procedure to admit or discharge a patient. 6. Parameter Setup Describes the procedure to set the measurement condition, size, scale, etc. for each parameter. 7. Function Describes about the functions such as arrhythmia analysis, trend, recall, etc. 8. System Configuration Describes about the system configuration such as night mode, alarm mode, display mode, etc. 9. Installation Describes about the environment for use, wireless system, etc. 10. Maintenance Describes about the maintenance, troubleshooting of this equipment. 11. Technical Information Lists the specification, default settings, pin assignments of external connector, etc. 12. Accessories Lists the accessories and optional

Preface

Composition of This Operation Manual ii		
Safety Precautions iii		
Labels Attached to the Unitiii		
Measurement Unit for Each Parametervi		
Graphic Symbols viii		
Precautions for Safe Operation of		
Medical Electrical Equipmentxi		
Precautions for Safe Operation of		
Medical Telemetry xii		
Precautions about the Maintenance xiii		
Precautions about the Pacemakerxiii		
Non-Explosion Proofxiv		
Defibrillation Safetyxiv		
Electrosurgery Safety xiv		
Precautions about Magnetic Resonance Imaging xiv		
Precautions about Connections to Peripheral Devicesxv		
Precautions about the Fusexv		
Accessories and Optional Accessories		
Precautions about the DS-7200 System xvi Precautions about the Wired Network System		
(DS-LAN II/DS-LAN III)xvii		
Precautions about the Wireless Network System xxx		
Precautions for Use of the Bidirectional Wireless		
Communications (TCON) xxxi		
Precautions about the Ventilator Monitoringxxxiii		
Precautions for Use of SpO2 Sensor xxxv		
Precautions for Masimo [®] Model: DS-7210M ········ xxxv		
Precautions for Use of NIBP Cuff		
Disposing of Equipment, Accessories,		
or Components xxxv		
Precautions about Transportation		
Precautions about RTC or Data Backup xxxv		
Precautions for Use of Lithium-Ion Battery Packxxxvi		
To Prepare for Emergency Usexxxvii		
Electromagnetic Compatibilityxxxviii		
Precautions for Safe Operation under		
Electromagnetic Influencexxviii EMC Guidancexxviii		
Compliance to the Electromagnetic		
Emissionsxxxix		
Compliance to the Electromagnetic		
Immunity (1)xxxix		
Compliance to the Electromagnetic		
Immunity (2)xl		
Recommended Separation Distances between		
Portable and Mobile RF Communications Equipment and the DS-7200 System		

1. General Description

General Description1-2
Combination of Option Unit and Measurement
Parameter ······1-3
Features ······1-4
Names of Parts and Their Functions1-5
DS-7200 Main Unit1-5
[Front Side]1-5
【Rear Side】1-5
[Right Side]1-6
[Left Side] ······1-6
[Bottom]1-7
Option Unit: MGU-721/MGU-722
[MGU-721 RGM Interface Unit]1-8
[MGU-722 Microstream [®] CO ₂ Unit]1-8
Option Unit: HU-71/HU-72/HU-731-9
[Front Side]1-9
【Rear Side】1-9
External Appearance1-10
●DS-72001-10
●Option Unit MGU-721/MGU-7221-10
●Option Unit HU-71/HU-72/HU-731-11

2. Basic Operation Basic Operation for Monitoring ------ 2-2

	_		
Before Use Basic Operation for Monitoring 2-			
Touch Keys2-			
General Key Control 2-			
•Key Control for Each Parameter 2-			
Before Use Home Display 2-			
About the Home Display 2-			
The Description of the Display2-	6		
Numeric Data, Waveform,	~		
Patient Name, etc	6		
Alarm Message for Numeric Data / Arrhythmia2-	6		
 Equipment Status Alarm Message 2- 			
● Ventilator Alarm Message ·······2-			
● Other Bed Alarm Message ······ 2-			
To Return the Display2-			
Preparation for Monitoring	0		
To Turn On the Power 2-	9		
Connecting the Power Cable 2-			
●Main Unit ····································	9		
To Turn ON/OFF the Power Switch 2-1	0		
To Use with the Battery Pack 2-1			
To Start Monitoring 2-1	2		
Discharge Confirmation at Power ON2-1			
To Admit a Patient 2-1			
Basic Operation2-1	5		
Adjusting the Waveform Size, Baseline Position			
(Parameter Key Operation)	5		
Selecting the ECG Lead			
(Parameter Key Operation) 2-1	6		
Scale, Lead, Baseline Position Setup			
(User Key Operation)2-1			
Alarm Setup for Each Parameter 2-1	8		
ON/OFF of Parameter Display	~		
Waveform/Numeric Data Display2-1	9		
Recording			
Start / Stop of Waveform Recording			
To Install the Paper			
To Suspend Monitoring2-2			
Discharging Procedure 2-2	4		
Ventilator Alarm Input			
Connecting the Ventilator 2-2 ●Multiport Relay Cable Connection 2-2	6		
Ventilator Connection			
● Ventilator Selection ······2-2	-		
● Ventilator Alarm Message ······ 2-2			
Check External Alarm			

3. Vital Application

10	Acquire ECG Waveform	3-2
	Before Attaching the Electrodes	3-2
	Electrode Placement	3-3
	Connection to the Patient Monitor	
	About the Arrhythmia Analysis	3-6
	QRS Classification	3-6
	●Arrhythmia Type	
	Filter Selection	
	Filter Mode Setup	
	Procedure for Filter Mode Selection	
	●AC Filter	
	spiration (Impedance Measurement)	3-9
То	Measure the SpO ₂	
	(Nellcor [®] Model: DS-7210)	3-10
То	Measure the SpO ₂ (Masimo [®] Model: DS-7210M)	2.44
	Measure the NIBP	
10		
	Out Organization and Datiant Application	0.40
	Cuff Connection and Patient Application	
	Procedure for Periodic Measurement	3-22
	Procedure for Periodic Measurement Dyna Alert Function Status	····3-22 ····3-23
	Procedure for Periodic Measurement Dyna Alert Function Status	····3-22 ····3-23 ····3-24
То	Procedure for Periodic Measurement Dyna Alert Function Status	3-22 3-23 3-24 3-25
То	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721)	3-22 3-23 3-24 3-25 3-29
To To	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display	3-22 3-23 3-24 3-25 3-29 3-29
To To	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Measure the CO ₂ (Option Unit: MGU-722)	3-22 3-23 3-24 3-25 3-29 3-29 3-2 9
To To	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Patient Application and Display	3-22 3-23 3-24 3-25 3-29 3-29 3-32 3-32
То То То	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Patient Application and Display Procedure for Calibration	3-22 3-23 3-24 3-25 3-29 3-29 3-32 3-32 3-33
То То То То	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Patient Application and Display Procedure for Calibration Measure the Temperature	3-22 3-23 3-24 3-25 3-29 3-29 3-32 3-32 3-33 3-36
To To To To To	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Measure the CO ₂ (Option Unit: MGU-722) Patient Application and Display Procedure for Calibration Measure the Temperature Measure the Cardiac Output (HU-73)	3-22 3-23 3-24 3-25 3-29 3-29 3-32 3-32 3-33 3-36 3-38
To To To To To	Procedure for Periodic Measurement Dyna Alert Function Status Oscillation Graph Display Measure the BP Measure the CO ₂ (Option Unit: MGU-721) Patient Application and Display Patient Application and Display Procedure for Calibration Measure the Temperature	3-22 3-23 3-24 3-29 3-29 3-29 3-32 3-32 3-33 3-36 3-38 3-38

4. Monitoring Setup

Display Configuration For	
To Configure the Display	Standard Mode 34-
To Configure the Display	12-lead Mode 4-7
To Configure the Display	Extended 1 Mode 4-9
To Configure the Display	Extended 2 Mode 4-11
To Configure the Display	Enlarge Mode 4-12
To Configure the Display	Block Cascade 4-14
The Corresponding Key for	
Each Numeric Data Box ····	
To Display the Short Trend	
To Configure the Display	
	timization4-22
●Set the Display Priori	ty 4-22
Optimize the Display	
●Waveform Display Ar	ea 4-24
Numeric Data Box Dis	
(for all parameters)	
Numeric Data Box Dis	splay
(for each parameter) ····	
Alarm Limit Display	4-31
	4-31
	4-32
	Display 4-33
Lead Selection for ST	Data Box 4-34
	age and Alarm Sound 4-35
Equipment Status Ala	arm Message 4-37
Numeric Data Box Me	essage 4-40
	sage 4-44
	or 4-45
Laser Printer Status N	
	sed)
Description of Bidirectional	
	Display 4-46
	 4-47 4-48
To Configure the Menu Dis	ріау 4-50
Erasing the Unnecessary K	
Key Mask Setup	
Recording Setup Waveform	
Manual Recording	4-52
	ecording 4-52
	etup 4-52
	······ 4-53
	4-55
-	ing
	4-56
	Waveform 4-59
	4-60
	c/Tabular Trend, etc.)···· 4-61
	e, Alarm Sound, etc 4-65
Color/Brightness Setup	
	, Waveform) 4-67
Brightness Sotup	, vvaveionn) 4-67
Digitiless Detup	4-07

Sweep Speed
Waveform Display Speed/Time4-68
Remote Control Setup4-69
About the Remote Control Unit4-69
Remote Control Setup4-70
Section Setup
Alarm Pole Setup Notifying the Alarm by Light 4-73
Flashing the Alarm Pole at Alarm Generation4-73
Flashing the Alarm Pole Synchronized with HR4-76
Alarm Setup To Set the Alarm Condition4-77
To Set the System Alarm (ON or Suspend)4-77
To Silence the System Alarm
Precautions about Silencing the Alarm4-78
●To Cancel "Alarm Silence"
Alarm Setup for Each Parameter4-79
Arrhythmia Alarm Setup4-82
To Set ON/OFF of Arrhythmia Alarm4-83
●To Set the Arrhythmia Detection Level ·······4-84
●To Set the HR Low Limit for VT ······4-84
●To Set the HR Low Limit for RUN ······4-84
Alarm Limit for TACHY, BRADY4-85
SpO ₂ SEC Alarm Setup
(Nellcor [®] Model; DS-7210)4-86
Alarm Suspend / Alarm Silence Time4-87
ON / OFF of Alarm Limit Display4-88
Alarm Silence Time for Equipment Status Alarm4-89
ON/OFF of Alarm Occurrence at NIBP Failure4-90

5. Admit / Discharge of a Patient

Adm	it / Discharge of a Patient	5-2
Adm	itting a Patient Name, Sex, and Age	5-3
Р	Patient Name	5-3
Р	atient ID	5-4
Р	Patient Classification	5-5
Р	atient Sex	5-6
Р	acemaker Use	5-6
	Pacemaker Use Selection	5-6
	Impedance Respiration Measurement	5-7
Р	Patient Age	5-7
Α	cquiring Patient Information from the Patient Data	
S	erver (TCON system) ·····	5-8
	harging a Patient	
	rasing Name, Data, etc	
D	Discharging Procedure	5-9
D	Discharge Procedure by User Key 5	-10
	itoring Mode Selection	
Α	larm / Display Mode5	-11
N	Node Selection 5	-11
	pend Monitoring	
S	Suspend / Resume Monitoring5	-12
Т	o Suspend Monitoring5	-12
Т	o Resume Monitoring5	-13

6. Parameter Setup

Par	rameter Setup	
	Setting the Monitoring Condition	
	To Display the Parameter Setup Menu	
	Zero Balance of All Pressure Lines	6-3
	Zero Balance of All Pressure Lines (User Key)	
	G	
	ECG Waveform Size and Lead	6-5
	HR/PR Alarm	6-7
	Arrhythmia Alarm	6-7
	Filter Mode Selection	6-8
	HR Average Selection	6-9
	HR Synchronization Mark and Pulse Tone	
	HR/PR Alarm Source	
	Automatic Lead Switching	
	Pacemaker Pulse	
	QRS Pace Pulse Mask	
	AC Filter	
	ECG Drift Filter	-
	3-lead Override	
	Pace Pulse Detection	
	ON/OFF of Parameter Display	
	spiration	
	Respiration Waveform Size	
	RR Alarm	
	APNEA Alarm	
	CVA Detection	
	RR/APNEA Alarm Source	
	Impedance Respiration Measurement RR Synchronized Indicator	
	ON/OFF of Parameter Display	
	ON/OFF of Parameter Display	0-22
1		
	asive Blood Pressure (BP1 to BP5)	6-23
	asive Blood Pressure (BP1 to BP5) BP Label	•• 6-23 ••6-23
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement	•• 6-23 ••6-23 ••6-25
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale	•• 6-23 ••6-23 ••6-25 ••6-26
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm	•• 6-23 ••6-23 ••6-25 ••6-26 ••6-26
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines	6-23 6-25 6-26 6-26 6-27
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection	6-23 6-25 6-26 6-26 6-27 6-28
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART)	6-23 6-25 6-26 6-26 6-27 6-28 6-28
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display	6-23 6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-29 6-29 6-30
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29 6-30 6-30 6-30
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29 6-30 6-30 6-31
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display 0 2 (Nellcor [®] Model: DS-7210)	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29 6-30 6-30 6-31 6-32
	asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display 0 2 (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size	6-23 6-25 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29 6-30 6-30 6-31 6-33
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display O 2 (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm	6-23 6-25 6-26 6-26 6-28 6-28 6-28 6-28 6-30 6-30 6-30 6-31 6-32 6-33 6-33 6-33
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display C2 (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source	6-23 6-25 6-26 6-26 6-28 6-28 6-28 6-28 6-29 6-30 6-30 6-31 6-32 6-33 6-33 6-33 6-34
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display C2 (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm Muring NIBP Measurement	6-23 6-25 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-29 6-30 6-30 6-31 6-33 6-33 6-33 6-33 6-34 6-35
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display	6-23 6-25 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-28 6-30 6-30 6-31 6-33 6-33 6-33 6-34 6-35 6-36
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display Model: DS-7210M)	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-28 6-28 6-30 6-30 6-31 6-33 6-33 6-33 6-34 6-35 6-36 6-37
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display Convorument of Parameter Display BP Convorument of Parameter Display Convorument of Parameter Display Convorument of Parameter Display Data during NIBP Measurement DN/OFF of Parameter Display Pulse Wave Size	6-23 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-28 6-28 6-30 6-30 6-31 6-33 6-33 6-33 6-34 6-35 6-36 6-37 6-38
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display O ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display O ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm	
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display ON/OFF of Parameter Display HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display O12 (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Monitoring Condition Setup	
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm ON/OFF of Parameter Display ON/OFF of Parameter Display	6-23 6-25 6-25 6-26 6-26 6-27 6-28 6-28 6-28 6-28 6-28 6-28 6-30 6-30 6-31 6-33 6-33 6-34 6-38 6-36 6-36 6-38 6-36 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-38 6-43 6-43 6-43 6-43 6-43 6-43 6-43 6-43
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Monitoring Condition Setup ON/OFF of Parameter Display N-Invasive Blood Pressure	
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display N/OFF of Parameter Display Pulse Wave Size SpO ₂ Alarm SpO ₂ Monitoring Condition Setup ON/OFF of Parameter Display NIBP Automatic Measurement	6-23 6-23 6-25 6-26 6-26 6-28 6-28 6-28 6-28 6-28 6-28 6-28 6-28 6-30 6-30 6-31 6-33 6-33 6-33 6-34 6-35 6-38 6-38 6-38 6-39 6-43 6-45
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display C2 (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display O2 (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm NIBP Automatic Measurement NIBP 1-Minute Interval Measurement	
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display O ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display O ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm NIBP Automatic Measurement NIBP Alarm	
Sp(asive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display C2 (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display O2 (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm SpO ₂ Alarm NIBP Automatic Measurement NIBP 1-Minute Interval Measurement	

Mean BP Display 6-48
Dyna Alert Function 6-49
Pump Setup
Sight Inflation
Pulse Rate Display6-51
Oscillograph Display6-51
Cancelling the NIBP System Error Message
Temperature (T1 to T3)
Temperature Label
Temperature Alarm
∆T Display
ON/OFF of Parameter Display
CO ₂ Concentration (Option Unit: MGU-721) 6-56
CO ₂ Scale6-56
EtCO2 (End-Tidal CO2) Alarm
InspCO ₂ (Inspiratory CO ₂) Alarm 6-58
EtCO ₂ Peak Picking Duration 6-58
Measurement Unit
O ₂ Compensation6-59
N ₂ O Compensation
Anesthetic Gas Compensation 6-60
Atmospheric Pressure Compensation
ON/OFF of Parameter Display6-62 CO ₂ Concentration (Option Unit: MGU-722) 6-63
ON/OFF of Parameter Display6-62
ON/OFF of Parameter Display
ON/OFF of Parameter Display6-62 CO ₂ Concentration (Option Unit: MGU-722) 6-63
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 6-64
ON/OFF of Parameter Display
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 6-65 InspCO2 (Inspiratory CO2) Alarm 6-65 6-65 EtCO2 Peak Picking Duration 6-65 6-65
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 InspCO2 (Inspiratory CO2) Alarm 6-65 EtCO2 Peak Picking Duration 6-66 CO2 Calibration 6-66 Restarting the CO2 Unit 6-67
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 InspCO2 (Inspiratory CO2) Alarm 6-65 EtCO2 Peak Picking Duration 6-66 CO2 Calibration 6-66 Restarting the CO2 Unit 6-67 ON/OFF of Parameter Display 6-68
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 InspCO2 (Inspiratory CO2) Alarm 6-65 EtCO2 Peak Picking Duration 6-66 CO2 Calibration 6-66 Restarting the CO2 Unit 6-67
ON/OFF of Parameter Display
$\begin{array}{llllllllllllllllllllllllllllllllllll$
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 6-65 InspCO2 (Inspiratory CO2) Alarm 6-65 EtCO2 Peak Picking Duration 6-66 CO2 Calibration 6-66 Restarting the CO2 Unit 6-67 ON/OFF of Parameter Display 6-68 Suspending CO2 Measurement 6-69 Ventilator 6-70 AWP Scale 6-70
ON/OFF of Parameter Display
ON/OFF of Parameter Display 6-62 CO2 Concentration (Option Unit: MGU-722) 6-63 CO2 Scale 6-63 EtCO2 (End-Tidal CO2) Alarm 6-64 InspCO2 (Inspiratory CO2) Alarm 6-65 EtCO2 Peak Picking Duration 6-66 CO2 Calibration 6-66 CO2 Calibration 6-66 Restarting the CO2 Unit 6-67 ON/OFF of Parameter Display 6-68 Suspending CO2 Measurement 6-69 Ventilator 6-70 AWP Scale 6-70 AWF Scale 6-70 Tb (Blood Temperature) 6-71 Tb Alarm 6-71
ON/OFF of Parameter Display

7. Function

Arrhythmia Analysis Definition, etc7-2
Arrhythmia Definition7-2
QRS Classification7-2
●Arrhythmia Type
To Set the Arrhythmia Alarm7-3
To Perform Arrhythmia Learning7-5
Graphic Trend Display/Print7-6
To Display the Graphic Trend7-6
Graphic Trend Group Setup7-10
The Description of the Display7-11
Tabular Trend Display/Print7-12
To Display the Tabular Trend7-12
The Description of the Display7-13
Parameter Setup for Tabular Trend
Recall Display/Print7-15
To Display the Recall Menu7-15
• Recall List Display7-15
To Display and Record the Enlarged Recall Waveform
To Set the Recall Condition7-17
NIBP List Display/Print7-19
To Display the NIBP List7-19
The Description of the Display7-20
ST Display ST Display, Alarm Setup, etc7-21
To Display the ST Measurement Menu
To Set the Reference Waveform7-22
●ST Alarm Setup ·······7-23
OCRG Display7-24
CO Measurement Measurement/Editing
To Display the Cardiac Output Menu
Message List
Cardiac Output Setup7-26
To Measure the Cardiac Output7-28
To Edit the Cardiac Output Data7-28
Hemodynamics Calculation/Print7-30
●Calculation Data ······7-30
To Display the Hemodynamic Data7-31
To Calculate the Newly Input Hemodynamic Data 7-31
To Edit the Hemodynamic Data7-33
Ventilator P-V, F-V Loop Display7-34
P-V Loop7-34
F-V Loop7-35
Displaying the Ventilator Measurement7-36
Respiration List Display/Print7-37
To Display the Respiration List7-37
The Description of the Display7-38
Respiration List Setup
ST Graphic Trend Display/Print7-39
To Display the ST Graphic Trend
ST Graphic Trend Group Setup
The Description of the Display
ST Tabular Trend Display/Print
To Display the ST Tabular Trend
The Description of the Display
Other Bed Display/Alarm
Other Bed Alarm Setup
Full Disclosure Waveform Recording
To Record the Full Disclosure Waveform Data

To Display the Full Disclosure Waveform ●About the Waveform Display Duration	7-47
To Shift the Displayed Waveform	7-47
To Change the ECG Waveform Size	7-50
ECG Waveform Amplitude	7-50
To Search by Time	7-50
To Search by Alarm	7-51
To Search by Admit/Discharge Date	7-52
To Print the Waveform	7-52
To Enlarge the Waveform	7-54
To Print the Waveform (Enlarged Display)	7-55
Vigilance/Vigileo List Display/Print	7-56
To Display the Vigilance/Vigileo List	7-56
The Description of the Display	7-57
Vigilance/Vigileo List Setup	

8. System Configuration

Night Mode	Night Mode Monitoring	8-2
About the	Night Mode	8-2
	Night Mode	
●To S	tart/Stop the Night Mode	8-4
●Nigh	t Mode Display Setup	8-4
Alarm Mode		
	ning the Alarm Mode	
	Alarm Mode	
	m the Alarm Mode	
To Initializ	e the Alarm Mode	
Monitor Setu		8-10
About the	Monitor Setup	8-10
	sage Icon	
	ck Discharge at Power ON	
	e/Date Setup	
	word Setup	
	ram Version	
	Setup	
	Mask	
	Key Setup	
	n Pole Setup	
	u Setup	
	lay Optimization Setup	
	kup at Discharge	
	Limit for the Alarm Volume	
	sword	
	itoring Condition after Discharge ht Key	0 40
	Filter display / Exp. clock display	
	Priler display / Exp. clock display PR Source	
	ze Mode Cursor ·····	
	meter Key Operation	
	Alarm Increment	
	(mmHg) Upper Limit for LAN,	0 17
Tele	metry	8-17
	ery Operation	
	e all alarms to "Recall"	
●Buzz	er Tone (Speaker) Failure	8-18
●Built	-in Rec. Status Display	8-18
●Vigila	ance/Vigileo SVR, SVRI, Calc	8-18
	n Level·····	
	n System	
	el 3 Alarm Sound	
	AN Setup	
	Alarm Increment	
	ent Name on Home Display	
	e Setup	
To Progra	m the Display Mode	8-23

Label Setup User Label	8-24
To Set the BP User Label	8-24
To Set the Temperature User Label	8-24
Hospital Setup Setup for Each Hospital	8-25
About the Hospital Setup	8-25
Date Format	
●Alarm Mute	
Arrhythmia Analysis Filter	
External Device Connection Setup	
●NIBP Data Erase Time ·····	
Status Output Setup ······	8-28
Measurement Unit Setup	8-29
Telemeter Setup	
●TCON Setup·····	
Trend Clip	
●BP Recording Scale ·····	8-30
Suspend Arrhy. Analysis during	
MEAN Calculation (BP1, ART, NIBP) Made Cancel	8-31
●Night Mode Cancel ●Asystole, VF, VT	8-31
	8-31
Patient ID Transmission Starting Digit for DS-LAN II	0 00
Admit/Discharge Key Setup ······	
HR/PR Low Limit during Alarm Auto Setting	
Password for Alarm Setup	
Recording Paper	
CF Card Data Transfer	
CF Card Data Transfer	
Full Disclosure Waveform Recording Setup	
Data Transfer (DS-7200 \rightarrow CF Card)	
Data Transfer (CF Card \rightarrow DS-7200)	
Data for Transfer	
Error Message	8-37

9. Installation

Precautions for Installing the Equipment	9-2
Precautions about the Operating Environment	9-2
Power Source and Ground Connection	9-3
Connecting the Slave Monitor	9-4
Attaching the Option Unit	
Attaching the CO ₂ Unit	
Attaching the OAO-12B Battery Pack	9-7
Wireless Network Construction	
To Attach the Telemetry Module	
Channel ID Setup	
Wired Network Construction	
DS-LANII Connection	
DS-LANIII Connection	-
Room/Bed ID Setup	9-18
Bidirectional Wireless Communications (TCON)	
System	9-19
TCON System Configuration Example	
HTC-702 Connection	
TCON ID / TCON Channel Setup	
Setup Item Synchronizing within the Same Netw	ork /
TCON System	
Ventilator Connection	
Oximeter	
Oximeter Connection	
Oximeter Connection	
BIS Monitor	
BIS Monitor BIS Monitor Connection	
External Device Connection Setup TCP/IP Network Connection	9-46
	-
Network Configuration (This Unit)	
To Connect the DS-7200 to the TCP/IP Network	
Network Configuration (Printer)	9-50

10. Maintenance

Handling After Use / Display Panel	10-2
Handling After Use	
Handling the Display Panel	10-2
Storage Device / Recording Paper	10-3
Storing the Device	
Storing the Recording Paper	
Cleaning Touch Panel and Housing	
Cleaning the Display Panel	
Cleaning the Housing	10-5
Cleaning the ECG Lead Cable, Relay Cabl	
Disinfecting the Blood Pressure Transduce	
Cleaning/Disinfecting the Temperature Pro	
Cleaning the Cardiac Output Relay Cable.	
Cleaning and Disinfecting the SpO ₂ Sensor	
Cleaning and Sterilizing the Airway Adapte	
(for Capnostat5)	
Battery Handling the Battery Pack	
Handling the Battery Storing the Battery	
č	
Maintenance Check Daily and Periodic Check	
About the Maintenance Check	
Daily Check	
Periodic Check	
Periodic Replacement Parts	
Software Version	
Troubleshooting	10-12
ECG······	10-12
Respiration	10-14
Invasive Blood Pressure	10-15
SpO ₂ (Nellcor [®] Model: DS-7210)	
SpO ₂ (Masimo [®] Model: DS-7210M) ········	
Non-Invasive Blood Pressure	
Temperature	
Cardiac Output	
CO ₂ Concentration (Option Unit: MGU-722)	
CO ₂ Concentration (Option Unit: MGU-721	
Wired Network (DS-LAN II/ DS-LAN III) ····	
Telemetry	
TCON	
Remote Control Unit	
General	
Battery	
Ventilator	
Oximeter	
BIS Monitor	
Recorder	
Laser Printer	

11. Technical Information

Specification / Performance	
Specification	11-2
Performance	11-3
External Connection Pin Assignments	11-8
RS-232C Connector Output Signal	11-8
COM1 Connector ·····	11-8
COM2 Connector ·····	
COM3 Connector ·····	
Status I/O Signal (Status II Connector).	
ECG Analog Output Connector	
IBP Analog Output Connector	11-10
Setup Item Default and Backup	11-11
Patient Admit / Discharge	
Alarm Setup	
Parameter Setup	11-14
Review Function Setup	11-16
System Configuration Setup	
Tone/Volume Setup	11-19
Display Configuration	11-19
System Configuration Menu	11-22
Hospital Setup	11-24
Monitor Setup	11-26
Alarm Mode Setup	11-29

12. Accessories

Accessories	12-2
Accessories	12-2
Optional Accessories	12-3
ECG, Impedance Respiration Measurement	12-3
Invasive Blood Pressure Measurement	12-3
Non-Invasive Blood Pressure Measurement	12-4
Temperature Measurement	12-4
SpO ₂ Measurement (Nellcor [®] Type: DS-7210)·······	12-4
SpO ₂ Measurement (Masimo [®] Type: DS-7210M) ···	12-5
CO Measurement (HU-73) ·····	12-5
CO2 Concentration Measurement (MGU-722)	12-5
CO ₂ Concentration Measurement	
(MGU-721 with Respironics Novametrix, LLC.	
CAPNOSTAT [®] 5 CO ₂ sensor)	12-7
Others	12-7

Chapter 4

Monitoring Setup

This chapter describes the setup procedure according to the monitoring purpose.

Display Configuration For	Easier View	4-2
To Configure the Display	Standard Mode	
To Configure the Display	12-lead Mode ······	
To Configure the Display	Extended 1 Mode ···	
To Configure the Display	Extended 2 Mode	
To Configure the Display	Enlarge Mode	4-12
To Configure the Display	Block Cascade	4-14
The Corresponding Key fo	r Each Numeric Data	a Box
To Display the Short Trend]	4-19
Waveform Grid Display		
Wave Line Thickness Waveform Clipping		
To Configure the Display		4-21
Auto Mode and Display Op	atimization	1-22
Set the Display Prior		
Optimize the Display		
Description of the Display.		
Waveform Display A		
Numeric Data Box D	isplay (for all parame	eters)
Numeric Data Box D	isplay (for each para	meter)
Alarm Limit Display		
Short Trend Display.		
Cursor Display		4-32
Freeze Mode Cursor	Display	4-33
Lead Selection for S ²		
Description of Alarm Mess		
Vital Alarm Message		4-35
Equipment Status Al		
 Numeric Data Box M 		
●Lead-Off Message…	lessage	/ /2
Ventilator Alarm Mes		
Ventilator Alarm Fac		4-44
Laser Printer Status	Message	
(When laser printer is u	used)	4-44
Description of Bidirectional	I Wireless Communic	cations
(TCON) Display		
Key Setup For Easier Use)	4-46
To Set the User Keys		4-47
To Configure the Menu Dis		
Erasing the Unnecessary F		
	m / Numeric Data …	
Manual Recording		
To Start / Stop the R		
Manual Recording S	etup	4-51
Alarm Recording		4-52
Periodic Recording		4-54
12-Lead Waveform Record	ding·····	4-55
Recording Setup		4-55
To Print the 12-Lead	Waveform	4-58

Recorder Setup ······4-59
Freeze Recording4-60
Graphic Recording (Graphic/Tabular Trend, etc.)4-60
Volume Setup Pulse Tone, Alarm Sound, etc4-64
Color/Brightness Setup4-66
Color Setup (Numeric Data, Waveform)4-66
Brightness Setup4-66
Sweep Speed Waveform Display Speed/Time 4-67
Remote Control Setup4-68
About the Remote Control Unit4-68
Remote Control Setup
Section Setup4-71
Alarm Pole Setup Notifying the Alarm by Light 4-72
Flashing the Alarm Pole at Alarm Generation4-72 Flashing the Alarm Pole Synchronized with HR4-75
· ·
Alarm Setup To Set the Alarm Condition4-76 To Set the System Alarm (ON or Suspend)4-76
To Silence the System Alarm (ON of Suspend)
 Precautions about Silencing the Alarm4-77
To Cancel "Alarm Silence"
Alarm Setup for Each Parameter4-78
Arrhythmia Alarm Setup ······4-81
•To Set ON/OFF of Arrhythmia Alarm
• To Set the Arrhythmia Detection Level4-83
To Set the HR Low Limit for VT
To Set the HR Low Limit for RUN
 Alarm Limit for TACHY, BRADY
SpO_2 SEC Alarm Setup
(Nellcor [®] Model; DS-7210)······4-85
Alarm Suspend / Alarm Silence Time4-86
ON / OFF of Alarm Limit Display4-87
Alarm Silence Time for Equipment Status Alarm 4-88
ON/OFF of Alarm Occurrence at NIBP Failure4-89

Display Configuration

For Easier View

The monitoring display can be configured according to the monitoring purpose. There are 5 types of basic display mode, which are Standard, 12-lead, Extended 1, Extended 2, and Enlarged. For the Standard mode, graphic trend, ventilator, tabular trend, NIBP list, OCRG, block cascade can be simultaneously displayed.



Also, 5 patterns of configured display can be programmed using the display mode setup function. By preprogramming the configuration to each display mode, the display configuration setups at admittance of patient can be simplified by just selecting one of the display modes.



For display mode setup procedure, refer to "8. System Configuration Display Mode To Program the Display Mode"

Standard Mode

To Configure the Display

For the standard mode, maximum of 10 waveforms and 8 numeric data can be displayed. On the waveform display area, graphic trend, ventilator, tabular trend, NIBP list, OCRG can be also displayed. If block cascade is selected, long duration waveform can be displayed. The duration of waveform display is approximately 7.9 seconds.



<Standard>



<Standard (Ventilator)>



<Standard (NIBP List)>



<Standard (Graphic Trend)>



<Standard (Tabular Trend)>



<Standard (OCRG)>

1 Press the Menu \rightarrow Display Config. keys.



Press the Setup key.

2 Select the waveform to display.



Pressing one of the waveform display location will display the waveform selection window. Select the parameter.

Selecting the parameter on the waveform selection window will sequentially assign the parameter from the top. To change the selection, pressing the waveform display location key will allow reselecting the parameter for that location.

<Waveform Display Location Key>

<Waveform Selection Window>



If the same parameter is repeatedly selected, the waveform display area for that parameter will be enlarged.

ĺ	ECG1
[ECG1
[ECG1
	ECG1
[ECG1
[BP Overlap1

3 Select the numeric data to display.



Pressing one of the numeric data display location will display the numeric data selection window. Select the parameter.

Selecting the parameter on the numeric data selection window will sequentially assign the parameter from the top. To change the selection, pressing the numeric data display location key will allow reselecting the parameter for that location.

<Numeric Data Display Location Key>

He de de la de la

<Numeric Data Selection Window>

The numeric data display layout can be changed. By repeatedly selecting the same parameter, the numeric data display area for that parameter will be enlarged. The parameter can be repeatedly selected for up to 3 times.

	HR
ļ	HR
ĺ	HR
į	BP1

Reference

For details of the data which will be displayed by selecting each key, refer to "The Corresponding Key for Each Numeric Data Box" in this chapter.

4 Select the function display.



Press the Func. Disp key and display the function display selection window.

Function Display Trend (3rous) Trend (5rous) List (3rous) List (5rous) NIBP List NIBP List NIBP List	Select the function display from Trend, List, NIBP List, OCRG, VENT 3 rows / 6 rows indicates the size to display the graphic trend and OCRG.
UCRG (3rows) UENT. (5rows) *The display area is fixed. (lowermost part of waveform area) Close	<function display="" tool=""></function>

5 Select "Standard" for the display mode.

Display Config. 1/2 Mode select Page Down Prov. Display Mode Standard Setup 12-lead Setup Ext1 Setup Ext2 Setup Enlarge Setup Auto	Press the <u>Prev. Disp.</u> key to return to the display configuration menu. Then, press the <u>Standard</u> key for the display mode
--	--

If performing telemetry or wired network transmission, configure the display so that the numeric data corresponding to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted.

NOTE	 After configuring the display, press the Home key to verify the programmed display configuration. To maintain the configured display even after the power is turned OFF or after a discharge procedure, save the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode".
------	---

To Configure the Display

12-lead Mode

The 12-lead mode displays 12-lead ECG waveforms and other selected waveforms (respiration, SpO_2 , etc.). ECG waveforms are displayed in segments at the center of the waveform display area. The number of waveforms and numeric data that can be displayed are as follows.

	ECG 12-lead (3.9 sec.)			
	+			
Waveforms (Display Duration)	Other waveforms: max. 2 (approx. 7.9 sec.)			
	Total: Max. 14 waveforms			
Numeric Data	Max. 8 numeric data			
Numeric Data	Max. o humenc uala			
BED-001 FUKUDA DENSHI	02/01 18:42 M			
Kindrahandraha King yang yang yang	60			
Un de de de de 1/2 de 1/2 de 1/2	1			
	e#u % ## % V %			
	V # V = V = V =			
hi had been the the second the				
and a second sec	116/ 77 (92) www.			
All hand have been been been been been been been be	92			
8. II	30			
	129/° 82			
March Record Size/ Adiret/ Essay Key Lock of the NIBP	THE REPORT OF TH			
norm Slerce Start/Starp Scale Decharge Process bittl 2 set inum Auto Hode Start/Starp none				
Press the Menu \rightarrow Display	Config. keys and display the display configuration menu.			

1 Press the Menu \rightarrow Display Config. keys and display the display configuration menu. Then press the Setup key for the 12-lead mode.





<Display Configuration Menu>

2 Select the waveform and numeric data to display.

Prev. Disp.
HR
ST-A
ST-B
ST-C
BP1
Sp02
RR_IMP
NIBP Over lap

Pressing the waveform display location key will display the waveform parameter selection window. Select the parameter.

<Waveform Display Location Key>



Pressing one of the numeric data display location will display the numeric data selection window. Select the parameter.

The numeric data display layout can be changed. By repeatedly selecting the same parameter, the numeric data display area for that parameter will be enlarged.

The parameter can be repeatedly selected for up to 3 times.



<Numeric Data Display Location Key>



For details of the data which will be displayed by selecting each key, refer to "The Corresponding Key for Each Numeric Data Box" in this chapter.

3 Select 12-lead for the display mode.



A CAUTION	If performing telemetry or wired network transmission, display the numeric data corresponding to the waveform. If not, the displayed waveform or numeric data may not be transmitted.
------------------	---

 After configuring the display, press the Home key to verify the programmed display configuration. To maintain the configured display even after the power is turned OFF or after a discharge procedure, save the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode".



The 12-lead waveform can be recorded on the built-in recorder. For setup procedure of 12-lead waveform recording format, refer to "4. Monitoring Setup Recording Setup Recorder Setup"

To Configure the Display

Extended 1 Mode

For The Extended 1 mode, Maximum of 10 waveforms and 16 numeric data can be displayed. The waveform display duration is about 5.5 seconds.



1 Press the Menu \rightarrow Display Config. keys and press the Setup key for the Ext1 mode.



The display configuration menu will be displayed.

<Display Configuration Menu>

2 Select the waveform to display.



Pressing the waveform display location key will display the waveform parameter selection window. Select the parameter.

<Waveform Display Location Key>

3 Select the numeric data to display.

splay (Ext. 1)				Prev Disp
ECG1	HR	- l	JPC+PACE	130
ECG1		T	ST-A	
ECG1				
BP Overlap1	BP		BP2	
BP Overlap1	BP		BP2	L.
BP Overlap1	Sp0	2	TEMP1,2	L.
BP Overlap1	SpO	, [RR_IMP	
BP Overlap1		_		loc Cas
Sp02	NIB	P N	IBP LIST	etu
RESP	NIB	P N	IBP LIST	Over Tap

Pressing one of the numeric data display location will display the numeric data selection window. Select the parameter.

<Numeric Data Display Location Key>

By repeatedly selecting the same parameter, the numeric data display area for that parameter will be enlarged. The same parameter can be repeatedly assigned for up to 3 rows and 2 columns.



Reference

For details of the data which will be displayed by selecting each key, refer to "The Corresponding Key for Each Numeric Data Box" in this chapter.

4 Select Extended 1 mode for the display mode.

Display Config. 1/2 Mode select Page Down Disp.	
Display Mode Standard Setup	
12-lead Setup Ext1 Setup	Press the Prev. Disp. key to return to the display configuration menu.
Enlarge Setup	Then, press the Ext1 key for the display mode.
Auto 📰 / 📰	

▲ CAUTION	If performing wired network transmission, display the numeric data corresponding to the waveform. If not, the displayed waveform or numeric data may not be transmitted.			
NOTE	 After configuring the display, press the Home key to verify the programmed display configuration. To maintain the configured display even after the power is turned OFF or after the discharge procedure, save the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode". 			

To Configure the Display

Extended 2 Mode

For the Extended 2 mode, Maximum of 10 waveforms and 9 numeric data can be displayed. The waveform display duration is about 5.5 seconds.



1 Press the Menu \rightarrow Display Config. keys, and press the Setup key for the Ext2 mode.



The display configuration menu will be displayed.

<Display Configuration Menu>

2 Select the waveform and numeric data to display.

Prev. Disp.

> Block Cas. Setu

Display (Ext.2)				Prev. Disp.
	ECG1				
	ECG1		н	R	
	BP Overlap1				
	BP Overlap1				
	BP Overlap1		NIBP	Sp02	
	BP Overlap1		MIDP	Spoz	
	BP Overlap1		BP1	BP2	
	Sp02	٦	BF1	DFZ	Block Cas.
	Sp02		TEMP 1	TEMP2	Setup
	C02	٦	C02	RR_C02	Over lap

NIBF

BP1 BP2

TEMP 1

CO2 RR_CO2

TEMP2

Pressing the waveform display location key will display the waveform parameter selection window. Select the parameter.

<Waveform Display Location Key>

Pressing one of the numeric data display location will display the numeric data selection window. Select the parameter. The numeric data display layout is fixed.

<Numeric Data Display Location Key>



Display (Ext. 2)

ECG1

ECG1 BP Overlap1 BP Overlap1 BP Overlap1

BP Overlap1 BP Overlap1

Sp02

SpO2

For details of the data which will be displayed by selecting each key, refer to "The Corresponding Key for Each Numeric Data Box" in this chapter.

3 Select Extended 2 mode for the display mode.



Press the Prev. Disp. key to display the display configuration menu.

Then, press the Ext2 key for the display mode.

A CAUTION	If performing telemetry or wired network transmission, display the numeric data corresponding to the waveform. If not, the displayed waveform or numeric data may not be transmitted.		
NOTE	 After configuring the display, press the Home key and verify the programmed display configuration. To maintain the configured display even after the power is turned OFF or after the discharge procedure, store the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode". 		

To Configure the Display

Enlarge Mode

The "Enlarge Mode" displays the numeric data in enlarged format. Maximum of 4 waveforms and 4 numeric data can be displayed. The waveform display duration is about 10.2 seconds.

BED-001 FUKUDA DENSHI	Askult 02/01 13:42 M	٦
the de de de de de	als de de de de de	~
	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	-
		ſ
60	92,	
^{NIBP} <b>129</b> / ^{13:42} <b>82</b>	RR_IMP X	
mmHg		_
Menu Alarm Record Size/ Admit/ Stience START/STOP Scale Discharge Free	20220 Key Lock Alarn NIBP NIBP Hone Start/STOP Hone	

**1** Press the Menu → Display Config. keys and press the Setup key for the Enlarge mode.



The display configuration menu will be displayed.

<Display Configuration Menu>

#### **2** Select the waveform to display.



Pressing the waveform display location key will display the waveform parameter selection window. Select the parameter.

<Waveform Display Location Key>

#### **3** Select the numeric data to display.



Pressing one of the numeric data display location will display the numeric data selection window. Select the parameter.

<Numeric Data Display Location Key>

<b>▲</b> CAUTION	If performing telemetry or wired network transmission, display the numeric data corresponded to the waveform. If not, the displayed waveform or numeric data may not be transmitted.			
NOTE	<ul> <li>After configuring the display, press the Home key and verify the programmed display configuration.</li> <li>To maintain the configured display even after the power is turned OFF or after the discharge procedure, store the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode".</li> </ul>			

## To Configure the Display

For the standard, extended 1, extended 2, and enlarge mode, block cascade waveform can be displayed.

When the display configuration is standard mode with 2 waveforms block cascade, the maximum duration of waveform display is about 39.5 seconds (5 blocks  $\times$  7.9 sec.).

When the display configuration is extended mode with 2 waveforms block cascade, the maximum duration of waveform display is about 27.5 seconds (5 blocks  $\times$  5.5 sec.).

When the display configuration is enlarge mode with 2 waveforms block cascade, the maximum duration of waveform display is about 20.4 seconds (2 blocks  $\times$  10.2 sec.).





**1** Press the Menu → Display Config. keys and press the Setup key for the display mode to set the block cascade.



The display configuration menu will be displayed.

<Display Configuration Menu>

2 Set the block cascade.



Pressing the Block Cas. Setup key will display the block cascade setup menu. Select the waveform quantity and parameter for the block cascade display.

Wave Qty 2 3 4 5 Te "Unue Qty" is channed wheek assende display rear will be set to "UP".				
ECG1	ECG1	ECG2	ECG3	ECG4
ECG2	ECG5	ECG6	ECG7	ECG8
	ECG9	ECG10	ECG11	ECG12
	BP1	BP2		
SpO ₂	SpO ₂	RESP		
RESP	AWF	AWP	OFF	

Select the waveform quantity for the block cascade from 2, 3, 4, 5, 6.

<Block Cascade Setup Menu>



#### **3** Select block cascade for the displaying waveform.



Pressing the waveform location key will display the waveform parameter selection window. Select Block Cascade. Pressing the block cascade key will set the block cascade for the set quantity.

<Waveform Display Location Key>

<Waveform Parameter Selection Window>

#### **4** Select the display mode which the block cascade was set.



Press the Prev. Disp. key to access the display configuration menu.

Then, select the display mode which the block cascade was set.

<b>A</b> CAUTION	If performing wired network transmission, display the numeric data corresponding to the waveform. If not, the displayed waveform or numeric data may not be transmitted.	
NOTE	<ul> <li>After configuring the display, press the Home key and verify the programmed display configuration.</li> <li>The block cascade setup is common for all display mode.</li> </ul>	

## The Corresponding Key for Each Numeric Data Box

The numeric data to be displayed can be selected on the numeric data selection window on the display configuration setup menu. Refer to the following for the corresponding key for each numeric data box.

1	Numeric Dat	a 1/2	Page Down	HR	Prev. Disp.
	HR	UPC+PACE	0FF		
	ST-A	ST-B	st-c		
	Sp02	PR_Sp02	Sp02/PR	BP1	
	NIBP	NIBP LIST		BP2	
	BP1	BP2	PR_IBP	Sp02	Func. Disp
	TEMP1	TEMP2	TEMP1,2	TEMP1,2	Block Cas.
	RR_IMP	C02	RR_C02	RR_IMP	Setup
U		Close		NIBP	Over lap

<Display Configuration Setup Menu Numeric Data (1st Page)>

Reference

For details of the displayed data for each numeric data box, refer to "Display Configuration Description of the Display ●Numeric Data Box Display (for each parameter)".

#### [Numeric Data / 1st Page]

Numeric Data 1/2 Page Down			
HR	UPC+PACE	0FF	
ST-A	ST-B	ST-C	
Sp02	PR_Sp02	Sp02/PR	
NIBP	NIBP LIST		
BP1	BP2	PR_IBP	
TEMP1	Temp2	TEMP1,2	
RR_IMP	C02	RR_CO2	
	Close		

ÍHR Heart Rate HR UPC UPC+PACE VPC, Pace Beat PACE ST (mm) I ≬ II ≬ II ≬ aVR ≬ ST Level ST-C ST-A ST-B Sp02 SpO₂ Value Sp02 PR_Sp02 Pulse Rate (SpO₂) PR_Sp02 Sp02 SpO₂ Value and Pulse Rate Sp02+PR NTRP **NIBP** Value NIBP 13:15 13:10 13:05 13:00 **NIBP** List NIBP LIST **BP** Value BP1 BP2 PR_IBP PR_IBP Pulse Rate (BP) **TEMP** Value TEMP1 TEMP2 TEMP1,2 36.1 RR_IMP RR_IMP **Respiration Rate (Impedance)** CO2 Et Insp EtCO₂ Value / InspCO₂ Value C02 RR_CO2 RR_CO2 Respiration Rate (CO₂)



## [Numeric Data / 2nd Page]

Numeric Data 2/2	Page Up	UENT Ventilator Data	TV Insp 400 Exp 416 TV 6.2
UENT RR_UENT Sv02+C0 BIS		<b>RR_UENT</b> Respiration Rate (Ventilator)	RR_UENT
HEMO SU HEMO-) Sur Ruw Rusw	I SUI Suri Ruwi Ruswi	Sv02+C0 Oximeter Data * Displayed data will differ depending on the used oximeter.	Sv02         83 x           C0 AUG         5.3 L/min           C1 AUG         2.8 L/min/m²           BSA         1.98
Cursor Ref. Cursor	STOP Watch	BIS BIS Value	BIS 58 SQI 87% Emg Odb SR 0%
BP3 TEMP3	Ть	HEno SU HEno-I SUI SUR SURI RUW RUWI RUSW RUSWI Hemodynamic (Based on Vigilance data)	SU         6         5           SUR         1         3         6         3           RUW         0.5         4         RUSU         8.1           SUI         3         8         5         5         5         5         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7         7
Close		Cursor Ref. Reference Cursor	RUUI         0 - 3 2           RUSUI         4 - 2           Cursor (Ref.:BP1)           Omnets
		Cursor Measurement Cursor	Cursor (BP1) Ref. 0
		STOP WATCH Stop Watch	TIMER1 00:00:00 TIMER2 00:00:00
		BP3 BP4 BP5 BP Value	BP3 X 34/21 (26) mHg
		TEMP3 TEMP Value	× 38.3
		Ть Blood Temperature (When CO is measured)	<b>44.9</b>
## To Display the Short Trend

The short trend can be selected to be displayed on the home display.

# **1** Press the Menu $\rightarrow$ Display Config. $\rightarrow$ Page Down keys.



The display configuration menu will be displayed.

Short Trend Selection
ON will display the short trend on the home display.
OFF will not display the short trend on the home display.
Overlap will display the waveform and short trend overlapped.

<Display Configuration Menu>

#### **2** Select the display type for the short trend.

BED-001 FUKUDA DENSHI	02/01 15:58 M
ly_str_str_str_str_str_str_str_s	60
	116/77 ( <u>92)</u> 23/_10
LAMMEL.	
	36.1 37.2 <b>3</b> 8
Norman         Reserve         Reserve <td< td=""><td>NIEP 129/⁰ 82 NIEP NIEP 100 NIEP NIEP NIEP NIEP NIEP NIEP NIEP NIEP</td></td<>	NIEP 129/ ⁰ 82 NIEP NIEP 100 NIEP NIEP NIEP NIEP NIEP NIEP NIEP NIEP
General Source Large (a) much as	PALO TODA STIRITOTOP
BED-001 FUKUDA DENSHI	02/01 15:58 M
" the day of a star	<b>60</b>
100 · PR	116/ 77 ( <u>92)</u>
	<b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b> <b>1</b>
	⇒ 36.1 37.2 ■ 36.1 37.2
	129/ 82

The short trend can be displayed in 5 minutes increments from 0 minute to 30 minutes.

Pressing the waveform display area will change the short trend display duration according to the pressed location.

NOTE	The short trend can be displayed only for the standard display mode. However, if graphic trend, ventilator, OCRG, tabular trend, or NIBP list is displayed, short trend cannot be displayed.
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# Waveform Grid Display



Page Up

ON OFF

ON OFF

The ECG waveform can be selected to be displayed on the grid.

**1** Press the Menu  $\rightarrow$  Display Config.  $\rightarrow$  Page Down keys.

Prev. Disp.

ie is selected

3

The display configuration menu will be displayed.

Grid SelectionON will display the grids on the home display.OFF will not display the grids on the home display.

<Display Configuration Menu>

### **Wave Line Thickness**

Display Config. 2/2

Short Trend

Grid

Wave Line

Thickness

Wave Clip

Invalid if function disp. is selected. Overlap invalid if ECG cascade or block casc

The thickness of the displayed waveforms can be selected from 3 levels.

**1** Press the Menu  $\rightarrow$  Display Config.  $\rightarrow$  Page Down keys.



# **Waveform Clipping**

If the waveform amplitude exceeds the waveform display area, whether or not to clip the exceeded part can be selected.

Press the Menu  $\rightarrow$  Display Config.  $\rightarrow$  Page Down keys. 1

Display Config. 2/2	Page Up Prev. Disp.
Short Trend	
Invalid if function Overlap invalid if E	disp. is selected. CG cascade or block cascade is selected.
Grid	
Wave Line Thickness	
Wave Clip	

ON will display the exceeded part of the waveform in straight line.

OFF will display the whole part of waveform even if it exceeds the display area. However, the exceeded part may not be displayed depending on the sweep speed of the waveform displayed above or below the gas waveform.

### To Configure the Display

# Auto Mode and Display Optimization

The home display layout can be automatically configured by setting the display mode to Auto, or pressing the Optimize Display key preprogrammed as user key.

The display will be automatically configured to either "Standard" mode or "Extended 1" mode depending on the quantity of the measured parameters.

NOTE	<ul> <li>The parameter that is not measured will not be displayed.</li> <li>The low priority parameter may not be displayed.</li> <li>Some parameters (ST, etc.) are not included in the optimization setup selection. To display these parameters on the home display, set the display configuration manually.</li> </ul>
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Set the Display Priority

Set the order of priority to display the parameters for optimized display configuration.

The display configuration optimization setup menu will be displayed.

#### **2** Set the order of priority for displaying the parameters.



**3** Select whether to overlap or separate the BP waveform display.

**BP** Format

Overlap Separate

### Optimize the Display

By selecting Auto for the display mode or by pressing the Optimize Display key preprogrammed as user key, the home display will be automatically configured by arranging the currently measured parameters according to the display priority set on the display optimization setup menu. The display mode will be either "Standard" mode or "Ext. 1" mode depending on the number of parameters.

**1** To optimize the display using the auto mode function, press the Menu  $\rightarrow$  Display Config. keys and display the display configuration menu.



Select Auto for the display mode, and return to the home display.

The home display will be automatically configured by arranging the currently measured parameters according to the display priority set on the display optimization setup menu. Each time the probe or sensor is plugged/unplugged, the home display will be reconfigured.

<Display Configuration Menu>

**2** Another way to optimize the display is to press the Optimize Display key preprogrammed as user key.



The confirmation message will be displayed. Pressing the OK key will optimize the display.

The optimized display configuration will be overwritten to the current display configuration (Standard / Ext. 1).
 To maintain the configured display even after the power is turned OFF or after the discharge procedure, store the configuration to one of the display modes, or select Backup for "Display Config." on the "Backup at Discharge" menu (Monitor Setup). For display mode setup procedure, refer to "8. System Configuration Display Mode".

# **Description of the Display**

This section explains the displayed item on the home display.



#### Bed / Room ID No.

Displays the 4-digit Bed ID and 3-digit (000 to 999) Room ID.

#### **Telemetry Channel ID**

Displays the channel ID of the built-in telemeter.

#### TCON Status (For Bidirectional Wireless Communications Module, HTC-702)

Displays the TCON communication condition, TCON ID, and TCON group (channel).



For more details, refer to "Description of Bidirectional Wireless Communications (TCON) Display" in this section.

#### **Battery Mark**

Displays the battery condition in 3 levels when the monitor is operated by battery.

Battery	Battery Condition	Indication of Operation Time		
Pack	Battery Condition	Normal Mode	Power Saving Mode	
	Full	About 3 to 2 hours	About 4 to 2 hours	
	The remaining battery is less than half.	About 2 hours to 20 minutes	About 2 hours to 20 minutes	
	The battery is almost empty. Connect to the AC power source immediately.	About 20 minutes or less	About 20 minutes or less	



For power saving mode, refer to "8. System Configuration Monitor Setup 

Battery Operation".

<b>A</b> CAUTION	The above operation time indicates the time with a new battery pack performing ECG measurement, NIBP periodic measurement (5-minute interval). Note that the battery pack degrades with continuous use and shortens the usable time.
------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### Date/Time

Displays the current date (month, day) and time (hour, minute).

#### Patient Name / Patient Classification

Displays the patient name and patient classification (adult / child / neonate) selected on the admit menu.

#### Pacemaker Use

When "Used" is selected for "Pacemaker" on the admit menu, "Pacemaker used" will be displayed.

#### **Drift Filter**

When drift filter is set to ON, "Drift-F ON" will be displayed. Depending on the setup, enlarged clock can be displayed instead of drift filter message.

O3/2711:34 ← Enlarged clock display



For procedure to select the drift filter or enlarged clock display, refer to "8. System Configuration Monitor Setup"

#### ECG Filter Mode

Displays the selected ECG filter mode. (M: Monitor mode, E: ESIS mode, D: Diagnosis mode)

#### **Respiratory Sweep Speed**

Displays the sweep speed for the impedance respiration waveform, CO₂ waveform, AWP, AWF waveform.

#### **Event Key**

This touch key will be displayed at alarm occurrence. Even when the alarm is resolved, this key will remained to be displayed until it is pressed. Pressing this key will silence the alarm and display the recall display. The event key display can be selected ON or OFF.

Reference

For ON/OFF of Event Key, refer to "8. System Configuration Monitor Setup".

#### Print Cancel Key a Cancel (For laser printer)

Pressing this Cancel key will cancel printing on the laser printer.



To use the laser printer, TCP/IP network setting is necessary. For procedure, refer to "9. Installation TCP/IP Network Connection".

#### 12LEAD REC key



This key will be displayed when the display mode is "12-lead". Pressing the key will record the 12-lead waveform on the built-in recorder or laser printer.

If laser printer is set for output recorder, laser printer icon will be displayed inside the 12LEAD REC key.



For setup procedure of 12-lead waveform output recorder and recording format, refer to "4. Monitoring Setup Recording Setup 12-lead Waveform Recording"

### Numeric Data Box Display (for all parameters)



Message Icon

When the parameter box size is too small to display the message inside, a message icon will be displayed instead to indicate that message is present. Alarm OFF Mark
Displayed when the alarm is set OFF.



For procedure to select ON/OFF of message icon display, refer to "8. System Configuration Monitor Setup".



### Out of Measurement Range (XXX)

The measurement is out of range.

#### Measurement Error (---)

- When the NIBP measurement ended erroneously. The time of measurement will be displayed.
- When the measurement is ceased manually, the time of measurement will be displayed as "- - -"

### Numeric Data Box Display (for each parameter)



For the corresponding numeric data selection key (display configuration setup menu) for each of the following numeric data box, refer to "Corresponding Key for Each Numeric Data Box" in this chapter.



#### HR Average (Instant / Average)

Displays the averaging method of HR. ("HR Average" on ECG configuration menu)

#### HR / PR Synchronization Mark

Synchronizing to the set HR/PR alarm source, a mark will be displayed. If  $\underline{SpO_2}$  is selected for "Pulse Tone" (ECG config.), the mark will be displayed inside the PR_SpO₂ numeric data box regardless of the HR/PR alarm source setup.

### HR / PR Parameter

Displays HR / PR value.

When the value exceeds the measurable range, "xxx" will be displayed.





#### SpO₂ Value

Displays the arterial oxygen saturation measurement value. - SEC Alarm Indicator (For DS-7210) Displayed when the SEC alarm is set.



#### PI Value (For DS-7210M)

If PI (Perfusion Index) display is set to ON, PI value will be displayed.



#### VPC Value (1 minute, 1 hour)

Displays the VPC rate for the last 1 minute and last 1 hour. "- - -" will be displayed during arrhythmia learning.

#### Pace Beats (1 minute) / Total Beats (1 minute)

Pace beats and total beats for the last 1 minute will be displayed. "- - -" will be displayed during arrhythmia learning.



#### ST Level

ST levels will be displayed. "– – –" will be displayed for the following case.

- $\boldsymbol{\cdot}$  during arrhythmia learning.
- during lead-off condition.
- $\cdot$  when "N" or "S" is not detected for QRS within 30 seconds.
- $\cdot$  when reference waveform is not set for ST measurement.

The leads displayed inside the ST level box can be changed. For procedure, refer to "●Set the Leads for ST Data Box" of this section.



Reference

RR Sync. Mark

#### **RR Value**

Displays the impedance RR /  $CO_2$  RR measurement value corresponded to the respiration synchronization source. When the measurable range is exceeded, "xxx" will be displayed. When the ECG relay cable for electrosurgical knife is used, or when impedance measurement is set to OFF, impedance RR will not be displayed.



#### **RR Synchronization Mark**

Synchronizing to the set RR/APNEA alarm source, a mark will be displayed inside the parameter box.

NOTE The RR synchronization mark will not be displayed when ventilator is the RR/APNEA alarm source, and when the measurement is performed by PURITAN-BENNETT ventilator.

Meas. Interval Meas. Time Cuff Press.



NIBP Value (SYS / DIA)

#### NIBP Measurement Interval

The NIBP measurement interval will be displayed. If the automatic mode is set to OFF, it will be left blank.

If Backup (Resume auto mode by manual measurement.) is selected for "NIBP Auto Mode" on the "Backup at Discharge" menu (Monitor Setup), the NIBP measurement will become in standby condition when a patient is discharged. (When measurement interval is set.) In this case, periodic measurement will resume when the measurement is performed manually, or when the measurement interval is changed.

#### **NIBP Measurement Time**

Displays the start time of NIBP measurement.

#### **NIBP Cuff Pressure**

Displays the cuff pressure during NIBP measurement. **NIBP Value** 

Displays the NIBP measurement value (SYS / DIA / Mean). The mean NIBP display can be set to ON or OFF on the NIBP configuration menu. The value will be displayed as "- - -" when the preprogrammed NIBP erase time has elapsed.



129/

129/

90

90

13:05

13:00

**BP** Label

116/

RP1

#### Oscillograph

The oscillograph can be displayed inside the NIBP numeric data box if the size is 2 boxes or larger and "Oscillograph" is set to ON in the NIBP setup menu.

#### **Dyna Alert Status**

If Dyna Alert function is set to ON, the Dyna Alert function status will be displayed.

#### NIBP List

The latest 4/8/12 data of NIBP list will be displayed. The number of displaying data depends on the size of parameter box.

#### **BP Value**

SYS/DIA/Mean BP value will be displayed.

On the BP configuration menu, ON/OFF of mean BP display can be selected. If the measurable range is exceeded, " $\times \times \times$ " will be displayed. If transducer is not connected or BP zero balance is not performed, "- - -" will be displayed.

### BP Label

The BP label setup for the blood pressure will be displayed.



#### -PCWP Value, PCWP Measured Time

When the BP label is PAP, PCWP (Pulmonary Capillary Wedge Pressure) and measured time can be displayed.

#### PDP Value

When the BP label is IAP, PDP (Peak Diastolic Pressure) of IABP can be displayed.

Systolic Pressure (SYS) = Peak Systolic Pressure (PSP). **CPP Value** 

When the BP label is ICP, and artery pressure is labeled as ART, CPP (Cerebral Perfusion Pressure) can be measured. CPP = Mean Value of Arterial Pressure – Mean Value of Intracranial Pressure If the CPP value is less than 0, ICP or ART is not measured, or zero balance has not been performed for ICP or ART, the value will be displayed as "--". Also, alarm cannot be set for CPP.



### MEAN_WAVE

Displayed when mean waveform is set ON (BP configuration menu.).



### TEMP Label

The label set for the temperature will be displayed. **TEMP Value** 

Displays the temperature measurement value. The 400 series temperature probe can be used. When the measurable range is exceeded, " $\times \times \times$ " will be displayed. When the 700 series is used, "- - -" will be displayed for the measurement value.





### Blood Temperature

By using the thermodilution catheter for the CO measurement, blood temperature can be displayed. When the measurable range is exceeded, " $\times \times \times$ " will be displayed.

#### EtCO₂ Value / InspCO₂ Value

Displays the end-tidal  $CO_2$  concentration and inspiratory  $CO_2$  concentration measurement value.

The measurement unit can be selected from mmHg / kPa / % on the  $CO_2$  configuration menu.



Ventilator Measurement
When ventilator is connect

When ventilator is connected, the ventilator measurement data will be displayed.

Sv02	83	7
CO AVG	5.3	≁ L∕min
CI AVG	2.8	L/min/m ²
Sv02 CO AUG CI AUG BSA	1.98	L/ III 1 11/ III-
L		

#### Oximeter Data

When oximeter (Vigilance / Vigilance CEDV / Vigilance II / Vigileo / OXIMETRIX3 / Q-vue / Q2 Computer) is connected, the oximeter data ( $SvO_2$ , CO, etc.) will be displayed. The displayed data will differ depending on the used oximeter.

Oximeter	Displayed Data				
Vigilance (CCO mode / STAT OFF / Index OFF)	SvO ₂ (ScvO ₂ )	ссо	EDV	ВТ	_
Vigilance (CCO mode / STAT ON / Index OFF)	SvO ₂ (ScvO ₂ )	CCO STAT	EDV STAT	BT	_
Vigilance (CCO mode / STAT OFF / Index ON)	SvO ₂ (ScvO ₂ )	CCI	EDVI	BT	_
Vigilance (CCO mode / STAT ON / Index ON)	SvO ₂ (ScvO ₂ )	CCI STAT	EDVI STAT	BT	_
Vigilance (ICO mode)	SvO ₂ (ScvO ₂ )	CO AVG	CI AVG		_
Oximetrix3	SvO ₂	CO AVG	CI AVG		_
Q-vue (CCO mode)	—	CCO	CCI	BT	_
Q-vue (CCO not measured)	—	CO AVG	—	_	—
Oximetrix3 + Q-vue (CCO mode)	SvO ₂	CCO	CCI	BT	—
Oximetrix3 + Q-vue (CCO not measured)	SvO ₂	CO AVG	—	_	_
Q2 Computer (CCO mode)	SvO ₂	CCO	CCI	BT	_
Q2 Computer (CCO not measured)	SvO ₂	CO AVG	CI AVG	_	BSA

#### Hemodynamic Data (Vigilance)

SU	65
sur	1363
RV₩	0.54
RVS₩	8.1
SVI	3 8
SURI	2304
RŲ₩I	0.32
RUSWI	4.2

Based on the CCO data measured by the Vigilance (or Vigilance CEDV / VigilanceII / Vigileo), the following hemodynamic data are calculated and displayed every second based on the following condition.

 $\cdot$  Measured on CCO mode of Vigilance. (not displayed for ICO mode)

 $\cdot$  SvO2 parameter key (oximeter numeric data box) is displayed.

• BP label is set as ART, PAP, CVP.

(If the unit is "kPa", the data is converted to "mmHg" for calculation.)

Parameter	Description	Equation
SV	Stroke Volume (mL/beat)	$\frac{\text{CCO} \times 1000}{\text{HR}}$
SVR	Systemic Vascular Resistance (dynes-sec-cm ⁻⁵ )	(MAP - CVP) × 79.90 CCO
RVW	Right Ventricular Work (kg⋅m)	$CCO \times (MPAP-CVP) \times 0.0136$
RVSW	Right Ventricular Stroke Work (g·m)	$SV \times (MPAP-CVP) \times 0.0136$
SVI	Stroke Volume Index (mL/beat/m ² )	SV BSA
SVRI	Systemic Vascular Resistance Index (dynes-sec-cm ⁻⁵ ·m ² )	SVR×BSA
RVWI	Right Ventricular Work Index (kg·m/m ² )	RVW BSA
RVSWI	Right Ventricular Stroke Work Index (g·m/m ² )	RVSW BSA

NOTEThe hemodynamic data based on Vigilance data will not be displayed on the<br/>Vigilance / Vigileo list. Only the data directly acquired from the Vigilance will be<br/>displayed on the Vigilance / Vigileo list. For SVR, SVRI, the displaying value can be<br/>selected from Vigilance or DS-7200 on "Vigilance/Vigileo SVR, SVRI<br/>calculation" (Page 4/4 of Monitor Setup)



#### **Stopwatch Key**

Functions as stopwatch.



#### **BIS Data**

When the A-2000 BIS Monitor is connected, BIS data (BIS, SQI, EMG, SR) will be displayed.

If SQI value is below 50%, the background color will turn gray. If SQI value is below 15%, BIS value and SR value will disappear.

### Alarm Limit Display



#### Alarm Limit

The alarm limit can be displayed beside each measurement value. If ON is selected for the individual alarm, the alarm limit will be displayed. The upper and lower limit will be displayed at upper and lower row respectively.

For BP and NIBP, each alarm limit of SYS, DIA, mean BP will be displayed from the top. ON/OFF of alarm limit display can be selected.

For SpO₂, if lower alarm limit is set to 85% or below, the alarm limit will be always displayed regardless of the ON/OFF setting of "Alarm Limit Display".



Reference

Refer to "4. Monitoring Setup Alarm Setup" for ON/OFF of alarm limit display.

### Short Trend Display



#### **Short Trend Display**

Short trend can be displayed beside the measurement data. Pressing the waveform display area will change the displayed trend time to the pressed position.

The trend display is in 5-minute increment from 0 minute to 30 minutes.

#### **Trend Scale**

The short trend scale will be displayed between the short trend and measurement data.

The displayed scale will be in accordance with the scale set on the graphic trend menu.

### Cursor Display

By configuring the Cursor and Cursor Ref to the numeric data display area, measurement cursor (white) and reference cursor (yellow) can be displayed on the selected waveform (ECG or BP) to verify the waveform amplitude.

The cursor position can be moved up or down using the AVV keys.



### •Freeze Mode Cursor Display

FUKUDA DENSHI

By selecting ON for "Freeze Mode Cursor Display" of the monitor setup menu, a vertical cursor (BP data cursor and interval measurement cursor) can be displayed on the home display when the Freeze (user key) is pressed.

The cursor can be moved freely and display the BP value at any cursor position. It can be also used to measure the time interval between each cursor.

#### **1** Press the Freeze key preprogrammed as user key.



Arrow key for BP data cursor

**BP** Data Cursor

Time Interval Cursor



 In the HR data box, time interval (unit: ms) between the BP data cursor and time interval cursor will be displayed.

In the BP data box, BP value at cursor position will be displayed.

#### **2** The cursors can be moved freely using the $\blacksquare \blacksquare \blacksquare$ keys.

The cursors can be also moved by directly pressing the waveform area. The time interval and BP value will be updated each time the cursor is moved.



Press this area to move the time interval cursor.

Press this area to move the BP data cursor.

#### **3** To erase the cursors, press the Freeze (user key) key.

NOTE	<ul> <li>The freeze mode cursor cannot pass over the erase bar (black bar erasing the old waveform).</li> <li>The freeze mode will not be automatically cancelled with the freeze mode cursor displayed.</li> <li>The freeze mode will be automatically cancelled if freeze mode cursor is not displayed for 30 seconds.</li> <li>The display width of short trend cannot be changed with the freeze mode cursor displayed.</li> <li>An error exists between the displayed waveform and cursor position. (For the sweep speed of 25mm/s; an error of 1dot[max.12msec] to the right of the cursor position)</li> </ul>
------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Lead Selection for ST Data Box

The ST value for 4 leads can be displayed in the ST data box. 3 groups (A, B, C) of lead combination can be programmed.

ST (mm)		×	0.5 0.2	
		Â-		
	aVR	₩-	·	

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  ST Disp. Lead Setup keys.



The ST Display Lead Setup menu will be displayed. Set the displaying lead for each group.

**2** Select the group to perform the setup and select the lead for that group.



Select the group from ST-A, ST-B, ST-C.

Select the lead to be displayed. Pressing the lead key will sequentially set the lead to the selected group from the top.

This section explains about the message displayed on the home display.

There are vital alarm message and equipment status alarm message which will be displayed at the top of the home display.

The alarms are classified in level 1, level 2, level 3, level 4, and the alarm message will be displayed according to the priority of level 1>level 2>level 3>level 4. The color of the displayed messages are red for level 1, yellow for level 2, blue for level 3, and white for level 4. The alarm tone will be different depending on the alarm system setting (IEC/FUKUDA DENSHI).

Reference n

The alarm system setting (IEC/FUKUDA DENSHI) can be performed on the "Monitor Setup" menu. For procedure, refer to "8. System Configuration Monitor Setup ●Alarm System" (Default: IEC)

#### [FUKUDA DENSHI]

Alarm Level	Description	Tone	Displayed Color
Level 1	Life Threatening Alarm	<ul><li>(1) Continuous tone with alternate high and low pitch sound</li><li>(2) Continuous rapid tone</li></ul>	Red
Level 2	Cautionary Alarm	<ul> <li>(1) 5 seconds interval alternate high and low pitch sound</li> <li>(2) 5 seconds interval rapid tone</li> </ul>	Yellow
Level 3	Treatment Needed Alarm	<ul> <li>(1) Single beep tone or 15 seconds interval alternate high and low pitch sound</li> <li>(2) Single rapid tone or 15 seconds interval rapid tone</li> </ul>	Blue
Level 4	Notification Alarm	Display Only	White

(1) When the tone setting is set between the 1st and 4th level from the lowest level(2) When the tone setting is set to the 5th level or above from the lowest level

#### [IEC]

Alarm Level	Description	Tone	Displayed Color
Level 1	Life Threatening Alarm	Continuous tone	Red
Level 2	Cautionary Alarm	5 seconds interval beep tone	Yellow
Level 3	Treatment Needed Alarm	Single beep tone (different tone from FUKUDA DENSHI mode) or 15 seconds interval beep tone (*)	Blue
Level 4	Notification Alarm	Display Only	White

▲ CAUTION	<ul> <li>The alarm priority is high for level 1 (life threatening alarm), medium for level 2 (cautionary alarm), and low for level 3 (treatment needed alarm).</li> <li>Alarm messages will be displayed according to the priority. (Level 1 → Level 2 → Level 3 → Level 4)</li> <li>For the same alarm level, the alarm message of the newer alarm will be displayed.</li> </ul>
-----------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



* The time interval for Level 3 alarm sound can be set. The setting can be performed under the "Monitor Setup" menu. For procedure, refer to "8. System Configuration Monitor Setup • Level 3 Alarm System Sound" (Default: One time)

#### Vital Alarm Message

The vital alarm message is generated when a measurement exceeds the alarm limit, or when arrhythmia is detected.

BED-001 FUKUDA DENSHI Check Electrodes	Adult	02/01 14:14 M Lower HR alarm	Numeric Alarm Message
		HR ♥Av.	

BED-001 FUKUDA DENSHI Check Electrodes	Adult	02/01 M	14:14	Arrhythmia Alarm Message
		HR	¥Av.	

There are 2 types of alarm messages, numeric alarm message and arrhythmia alarm message. If the 2 types of alarm generate at the same time, the numeric alarm message and arrhythmia alarm message will be alternately displayed for 2 seconds each. The message will be displayed according to the priority of the alarm level. If the alarms of the same level generate, the message for the newer alarm will be displayed.

selecting User for "Alarm Level" on the 4th page of Monitor Setup, the alarm level set by each user can be applied.	ACAUTION	selecting User for "Alarm Level" on the 4th page of Monitor Setup, the
---------------------------------------------------------------------------------------------------------------------	----------	------------------------------------------------------------------------

#### Life Threatening Alarm (Alarm Level 1)

Parameter	Message
HR	"Lower HR alarm"
пк	"Upper HR alarm"
	"Lower PR alarm"
PR (SpO ₂ , BP)	"Upper PR alarm"
BP (BP1/ART)	"Lower BP1 alarm" or "Lower ART alarm"
BF (BF I/ART)	"Upper BP1 alarm" or "Upper ART alarm"
SpO ₂	"Lower SpO ₂ alarm"
5pO ₂	"Upper SpO ₂ alarm"
Despiration	"Apnea alarm"
Respiration (Impedance, CO ₂ , Ventilator)	"Lower RR alarm"
	"Upper RR alarm"
NIBP	"Lower NIBP alarm"
	"Upper NIBP alarm"
CO ₂	"Upper EtCO₂ alarm"
	"Lower EtCO ₂ alarm"
	"ASYSTOLE"
	"VF"
Arrhythmia	"VT" "SLOW VT"
Arrhythmia	"TACHY"
	"BRADY"
	"RUN"

#### Cautionary Alarm (Alarm Level 2)

Parameter	Message
BP (BP2 to 5)	"Lower BP* alarm" or "Lower (label) alarm"
BF (BF2 10 5)	"Upper BP* alarm" or "Upper (label) alarm"
ST1 to 12	"Lower ST* alarm"
3111012	"Upper ST* alarm"
Temperature (TEMP1 to 3)	"Lower TEMP $*$ alarm" or "Lower (label) alarm"
	"Upper TEMP* alarm" or "Upper (label) alarm"
CO ₂	"Upper InspCO ₂ alarm"
	"PAUSE"
	"COUPLET"
Arrhythmia	"BIGEMINY"
	"TRIGEMINY"
	"FREQUENT"

* indicates the channel no. of BP, ST and TEMP.

#### **Treatment Needed Alarm (Alarm Level 3)**

Parameter	Message
None	

#### Notification Alarm (Alarm Level 4)

Parameter	Message	
All Alarm	"Alarm Suspend (***sec.)"	
Arrhythmic	"LEARN"	
Arrhythmia	"ARRHY. OFF"	

<b>A</b> CAUTION	Even during "LEARN" status, alarm for HR, ASYSTOLE, VF, TACHY, BRADY will be generated.
------------------	-----------------------------------------------------------------------------------------

NOTE	<ul> <li>(***sec) of the "Alarm Suspend (***sec)" message indicates the remaining time to suspend the alarm.</li> </ul>		
	<ul> <li>The "ARRHY OFF" message will be displayed when the ASYSTOLE, VF, VT, SLOW_VT, and HR alarm is OFF.</li> </ul>		

#### • Equipment Status Alarm Message

The equipment status alarm message will be displayed when proper monitoring cannot be performed. The alarm message will be displayed according to the priority of the alarm level. If more than one alarm with the same level is generated, the newest alarm message will be displayed.



#### **Cautionary Alarm (Alarm Level 2)**

Item	Message	
ECG Impedance Respiration	"Check Electrodes"	
Arrhythmia	"Cannot analyze" *1	
SpO ₂	"Check SpO ₂ Sensor"	
(DS-7210, DS-7210M)	"Replace SpO ₂ Sensor"	
SpO ₂ (DS-7210)	"No pulse detect"	
SpO	"Unknown SpO₂ Sensor"	
SpO ₂ (DS-7210M)	"SpO ₂ Low Perfusion"	
, , , , , , , , , , , , , , , , , , ,	"SpO ₂ Pulse Search"	
NIBP	"NIBP measurement failed." *2	
	"Check Sample Line"	
CO ₂	"Check CO ₂ Exhaust Port"	
(MGU-722)	"Check CO ₂ unit"	
	"CO ₂ Cal. Required"	
	"ECG Disconnected"	
	"BP* Disconnected" * ³	
	"BP* Transducer OFF" *3	
	"SpO ₂ Disconnected"	
	"T* Disconnected" * ³	
Connector	"Sample Line Disconnected" (MGU-722)	
	"CO ₂ Disconnected" (MGU-721)	
	"CO Disconnected"	
	"Multiport * Disconnected" * ³	
	"DS-LANII Disconnected" or "DS-LANIII Disconnected"	
	"Printer Cable Disconnected"	
	"Check Backup Battery"	
Others	"Check Equip. Config. (CO ₂ )"	
	"Charge the battery."	
	"Check Memory Card"	

- *1 This alarm will generate when analysis is suspended for more than 30 seconds, regardless of ON/OFF setting of "Suspend Arrhy. Analysis during Noise Interference" under Hospital Setup (Preset Menu).
- Refer to "10. Maintenance Troubleshooting ECG".
   *² "NIBP measurement failed." will be displayed only if "Alarm Occurrence at NIBP Failure" is set to ON on the alarm setup menu. (Default: OFF) * indicates the channel no. of BP and TEMP.
- *3

▲WARNING	When a parameter is in a connector-off condition, the alarm will be generated only on the bedside monitor and not on the central monitor. Make sure that the connector is securely connected. If the waveform/numeric data is not displayed for a monitored parameter, check the patient's condition and pay attention not to miss the connector-off condition.	
ACAUTION	<ul> <li>Even during "Cannot analyze" alarm generation, alarm for HR, ASYSTOLE, VF, TACHY, BRADY will be generated.</li> <li>If OFF is selected for "PI Display" under the SpO₂ configuration setup, "SpO₂ Low Perfusion" alarm will be indicated by message display only. The alarm sound will not be generated.</li> </ul>	
ΝΟΤΕ	The "NIBP measurement failed.", "Connector Disconnected" alarm can be cancelled by pressing the <u>Alarm Silence</u> key. When silencing the alarm, make sure that important alarm is not missed.	

#### **Treatment Needed Alarm (Alarm Level 3)**

Item	Message
NIBP	"Check NIBP Cuff, Air Hose"
Impedance Respiration	"CVA detected"
SpO ₂ (DS-7210)	"No pulse detect"
	"SpO ₂ Low Perfusion"
SpO ₂ (DS-7210M)	"SpO ₂ Pulse Search"
	"SpO ₂ Interference Detected"
	"Check SpO ₂ Sensor"
	"SpO ₂ Low Signal IQ"
ECG	"Pacing detection error"

#### Notification Alarm (Alarm Level 4)

Item	Message
	"Waveform is frozen. (**sec.)"
Operation	"Key Lock"
	"Night Mode Active"
	"ECG Low Amplitude"
ECG	"ECG Noise Interference"
	"ECG Unit Failure"
ECG, Impedance Respiration	
BP	"BP* Transducer OFF" *1
DF	"BP* Zeroing Required" *1
	"Unknown Temp. Sensor"
Temperature	"TEMP Unit Auto Check"
	"TEMP Unit Failure"
SpO ₂ (DS-7210)	"Motion artifact"
SpO ₂ (DS-7210, DS-7210M)	"SpO ₂ unit error"
Microstream [®] CO ₂	"CO ₂ Unit Failure (O)"
(MGU-722)	"CO ₂ Suspended"
	"CO ₂ Unit Failure (C)"
	"CO ₂ Sensor Failure"
Mainstream CO ₂	"CO ₂ Warming Up"
(MGU-721)	"Zero CO ₂ Adapter"
	"Check CO ₂ Adapter"
	"Unknown CO ₂ Sensor"

	"NIBP Unit Failure"
NIBP	"Initializing NIBP"
	"Check Built-in Recorder"
	"Built-in Recorder Unit Failure"
Built-in Recorder	"Check Rec. Paper (Built-in)" *2
	"Check Cassette (Built-in)" *2
	"Recorder busy"
All Alarm	"Alarm is silenced."
	"Check TCON Comm."
TCON	"Chk TCON Receive"
	"TCON Interference"
	"Main Unit Failure"
	"Display Unit Failure"
	"Sub Unit Failure"
Main Unit	"Analog Unit Failure"
Main Onit	"High Internal Temperature"
	"Display Unit Backlight Failure"
	"Check Rotary SW"
	"Check DIP-SW"
External Equipment	"HU Module Failure"
	"Telemetry Unit Failure"
	"Check Central Recorder"
Network Recorder	"Recorder busy (Cent.)"
	"Check Paper (Central)"
	"Check Cassette (Central)"

*1 "*" indicates the channel no. of BP and <u>TEMP</u>.

*² These messages will not be displayed if OFF is selected for "Built-in Rec. Status Display" on the Monitor Setup menu.

#### Numeric Data Box Message

The measurement status of each parameter will be displayed inside the corresponding numeric data box.



#### **HR Data**

Alarm Level	Message
1	"Upper HR alarm"
1	"Lower HR alarm"
2	"Check Electrodes"
2	"Lower ST alarm"
2	"Upper ST alarm"
3	"Pacing detection error"
4	"Out of range"
4	"ECG Low Amplitude"
4	"Artifact"
4	"Check Electrodes"
4	"ECG Unit Error"

ST

Alarm Level	Message	
2	"Lower ST alarm"	
2	"Upper ST alarm"	

### BP1 to 5

Alarm Level	Message
1(*)	"Lower BP alarm"
1(*)	"Upper BP alarm"
4	"Zeroing Required"
4	"Out of range"
te laval 4 fan DD4 and ADT, laval 0 fan ath an lak al	

* : Level 1 for BP1 and ART, Level 2 for other label

### Pulse Rate (BP Source)

Alarm Level	Message
1	"Upper PR alarm" (BP)
1	"Lower PR alarm" (BP)
4	"Out of range"

#### NIBP

Alarm Level	Message	Description
1	"Cannot Measure (C02)"	<ul> <li>Could not measure although the pressure dropped to minimum deflating pressure. (When not quick measurement)</li> <li>Could not measure although the pressure dropped to minimum deflating pressure. (During quick measurement)</li> </ul>
1(*)	"Exhaust Error (C03)"	<ul> <li>Exhaust was suspended for 15 seconds due to body motion.</li> <li>Deflation speed of 1mmHg/sec. continued for 5 seconds.</li> <li>20 seconds elapsed since the exhaust started, but deflation of 30mmHg was not achieved.</li> <li>10 seconds elapsed since the exhaust started, but the target deflation speed was not achieved.</li> </ul>
1(*)	"Insuff. Inflation (C04)"	<ul> <li>The first 5 pulse amplitude was on a decreasing trend.</li> <li>SYS cannot be measured.</li> <li>Pressure difference between the pulses is too large.</li> </ul>
1(*)	"Osc. Pattern Err. (C05)"	<ul> <li>Pressure difference between the pulses is too large.</li> <li>Too much pulse compensation.</li> </ul>
1(*)	"Meas. Error (C06)"	<ul> <li>Set inflation value&gt;systolic&gt;mean&gt;diastolic was not achieved.</li> <li>Systolic pulse amplitude is too small.</li> <li>Pulse pressure is too small.</li> <li>Diastolic pulse amplitude is too small.</li> </ul>
1(*)	"Meas. Timeout (C07)"	Measurement Timeout
1(*)	"PR Over (C08)"	Measured PR value is too large.
1(*)	"Overinflation (C09)"	Maximum pressure is exceeded.
1(*)	"Pulse Amp. Error (C10)"	<ul> <li>Pulse amplitude is too large.</li> <li>Pulse amplitude is too small.</li> </ul>
1(*)	"Check Cuff Size (C11)"	<ul> <li>Neonate cuff is detected with adult mode.</li> <li>Infant cuff is detected with adult mode.</li> <li>Neonate cuff is detected with child mode.</li> </ul>
2	"Lower NIBP Alarm"	
2	"Upper NIBP Alarm"	
2	"Inflation Timeout (C01)"	
2	"System Error (E08-1)"	Communication Error with Sub CPU
2	"System Error (E08-2)"	Watchdog Timeout
2	"System Error (E08-3)"	Pressure Offset Error
2	"System Error (E08-4)"	Pressure Comparison Failure
2	"System Error (E08-5)"	Sub CPU Power Supply Failure
2	"System Error (E08-6)"	Pressure Sensor 2 Power Supply Failure
2	"System Error (E08-7)"	Pressure Sensor 1 A/D Reference Power Voltage Failure
2	"System Error (E09-A)"	Maximum cuff pressure is exceeded.
2	"System Error (E09-B)"	Inflation time is exceeded.

Alarm Level	Message	Description
2	"System Error (E09-C)"	Quick mode timeout
2	"System Error (E09-D)"	Measurement started during the long term mode pause period
2	"System Error (E09-E)"	Measurement Time Over
2	"System Error (E09-F)"	Timeout of pressure data transmission from main CPU
2	"System Error (E09-G)"	Pressure Sensor 1 +5V Power Supply Failure
2	"System Error (E09-H)"	Zero Calibration Timeout
2	"System Error (E09-I)"	ROM Test Error
2	"System Error (E09-J)"	RAM Test Error
2	"System Error (E09-L)"	Clock transmission is ceased.
2	"System Error (E09-M)"	Communication Failure at Power ON
2	"System Error (E09-N)"	Pressure Comparison Failure
2	"System Error (E09-O)"	Maximum inflation time is exceeded.
2	"System Error (E09-Q)"	Measurement was started before zero calibration
2	"System Error (E09-R)"	Zero Calibration Error
2	"System Error (E09-S)"	Watchdog Timeout
2	"System Error (E09-T)"	+5V Digital Power Supply Failure
2	"System Error (E09-U)"	Main CPU Power Supply Failure
2	"System Error (E09-V)"	Pump Control Signal Failure
2	"System Error (E09-W)"	Quick Exhaust Valve Control Signal Failure
2	"System Error (E09-X)"	Sub CPU Constant Exhaust Valve Control Signal Failure
2	"System Error (E09-Y)"	Main CPU Constant Exhaust Valve Control Signal Failure
3	"Check NIBP Cuff, Air Hose"	
4	"NIBP Unit Error"	

4 Display Configuration

*: Level 4 if NIBP measurement is retried, level 1 if measurement is not retried.

"System Error" message can be cleared by pressing the Cancel NIBP System Error key on the second page of the NIBP configuration menu. If it cannot be cleared, failure can be considered.

## SpO₂ (Nellcor[®] Model: DS-7210)

Alarm Level	Message
1	"Lower SpO ₂ alarm"
1	"Upper SpO ₂ alarm"
2	"Replace SpO ₂ Sensor"
2	"Check SpO ₂ Sensor"
2	"No pulse detect"
4	"SpO ₂ Unit Error"
4	"Motion artifact"
4	"SpO ₂ Pulse search"

### SpO₂ (Masimo[®] Model: DS-7210M)

Alarm Level	Message
1	"Lower SpO ₂ alarm"
1	"Upper SpO ₂ alarm"
2	"Replace SpO ₂ Sensor"
2/3*	"Check SpO ₂ Sensor"
2	"SpO ₂ Low Perfusion"
2	"SpO ₂ Pulse search"
3	"SpO ₂ Interference Detected"
3	"Unknown SpO ₂ Sensor"
3	"SpO ₂ Low Signal IQ"
4	"SpO ₂ Unit Error"

*The alarm level differs depending on the cause.

### PR-SpO₂

Alarm Level	Message
1	"Upper PR alarm" (SpO ₂ )
1	"Lower PR alarm" (SpO ₂ )
4	"Out of range"

#### TEMP1 to 3

Alarm Level	Message
2	"Upper TEMP alarm"
2	"Lower TEMP alarm"
4	"Unknown Temp. Sensor"
4	"TEMP Unit Auto Check"
4	"Out of range"

#### Tb

Alarm Level	Message
2	"Upper Tb alarm"
2	"Lower Tb alarm"
4	"Out of range"

### RR (Impedance)

Alarm Level	Message
1	"Apnea alarm"
1	"Upper RR alarm"
1	"Lower RR alarm"
2	"Check Electrodes"
3	"CVA detected"
4	"Check ECG cable"
4	"Out of range"
4	"Suspended"

## RR (Ventilator)

Alarm Level	Message
1	"Apnea alarm"
1	"Upper RR alarm"
1	"Lower RR alarm"

### RR (CO₂)

Alarm Level	Message
1	"Apnea alarm"
1	"Upper RR alarm"
1	"Lower RR alarm"

### CO₂ (When MGU-721 is used)

Alarm Level	Message
1	"Zeroing CO ₂ "
1	"Upper EtCO ₂ alarm"
1	"Lower EtCO ₂ alarm"
2	"Upper InspCO ₂ alarm"
4	"CO ₂ Unit Error"
4	"CO ₂ Sensor Failure"
4	"CO ₂ Warming Up"
4	"Zero CO ₂ Adapter

4	"Check CO ₂ Adapter"
4	"Unknown CO ₂ Sensor"
4	"Out of range"

#### CO₂ (When MGU-722 is used)

Alarm Level	Message
1	"Upper EtCO ₂ alarm"
1	"Lower EtCO ₂ alarm"
2	"Check Sample Line"
2	"Check CO ₂ Exhaust Port"
2	"Check CO ₂ unit"
2	"CO ₂ Cal. Required"
2	"Upper InspCO ₂ alarm"
4	"Initializing CO ₂ "
4	"CO ₂ Suspended"

### Lead-Off Message

If the ECG electrodes are detached, HR alarm and arrhythmia alarm will not be generated. If this condition is left unresolved, a sudden change of the patient may not be noticed. Take prompt action when the lead-off condition is detected.

BED-001 FUKUDA DENSHI Check Electrodes	Adult	02/01 16:43 M
LEAD OFF		HR Av.
	Lead-Off Message	9

<b>A</b> CAUTION	While the "LEAD OFF" message is displayed, HR alarm and arrhythmia alarm will not function. Leaving this condition unresolved may result in missing a sudden change of the patient. Promptly check the electrodes when this message is displayed.
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### Ventilator Alarm Message

When ventilator is connected to the DS-7200, ventilator alarm and the connection status alarm will be generated. The alarm message with the higher alarm level will be displayed.

#### [Ventilator Alarm Message]



#### [Connection Status Alarm Message]

BED-001	FUKUDA DENSH		02/01 14:14 M	Connection Status
Life Th	reatening Alarm (Ala	rm Level 1)		
Life Th	reatening Alarm (Ala Parameter	rm Level 1) Mess	age	

#### **Cautionary Alarm (Alarm Level 2)**

Parameter	Message
None	

#### Treatment Needed Alarm (Alarm Level 3)

Parameter	Message				
None					
Notification Alarm (Alarm Lovel 4)					

#### Notification Alarm (Alarm Level 4)

Parameter	Message		
Ventilator	"Vent. Disable 🔆 "		
	"Vent. Online"		

### Ventilator Alarm Factor

For the SV-300, Servo-i, Servo-s, ventilator alarm factor if specified will be notified and displayed on the central monitor.

Displayed Alarm Message	Description
VENT AWP	Airway Pressure Alarm
VENT MV	Minute Ventilation Alarm
VENT APNEA	Apnea Alarm
VENT CONT. HP	Continuous High Pressure Alarm
VENT Upper FiO ₂	FiO ₂ Upper Limit Alarm
VENT Lower FiO ₂	FiO ₂ Lower Limit Alarm
VENT Upper CO ₂	EtCO ₂ Upper Limit Alarm
VENT Lower CO ₂	EtCO ₂ Lower Limit Alarm
VENT Upper RR	RR Upper Limit Alarm
VENT Lower RR	RR Lower Limit Alarm
VENT PEEP	PEEP Low Alarm
VENT COMM	Power OFF, Cable disconnected, Standby condition, etc.
VENT URGENT	Other high level alarm
VENT	Other ventilator alarm

▲ CAUTION	<ul> <li>For the SV-900, PB, Evita, and Savina ventilator alarm factor cannot be transmitted to the central monitor.</li> <li>Depending on the central monitor type and software version, ventilator alarm factor may not be displayed. For details, refer to our service representative.</li> <li>The ventilator alarm factor listed above are displayed only on the central monitor. These will not be displayed on the bedside monitor.</li> </ul>
-----------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### •Laser Printer Status Message (When laser printer is used)

The laser printer status will be displayed.

LP Com Error

: Printer is not connected to the TCP/IP network, paper-out condition, etc.

Cancel

: In process of printer output. indicates the current output progress.



: Waiting condition for printer output.

## **Description of Bidirectional Wireless Communications (TCON) Display**

This section explains about the message displayed on the home display when performing the bidirectional wireless communications (TCON).



Indications	Ĭ	Y,	¥,	Ť×
Communication Condition	Good	Moderately Good	Bad	Cannot Communicate

# **Key Setup**

# For Easier Use

The DS-7200 system operation is performed through the touch keys displayed on the screen. The touch keys consist of fixed keys (Menu, Home / Enlarge, NIBP Start/Stop) and 6 or 8 user keys which can be programmed according to the monitoring purpose.



<Fixed Keys>



<User Keys>

Menu	Alarm Silence	Rec.	Size/ Scale	Admit/ Discharge	 Freeze	Key Lock (HOLD 2 SEC)	NIBP START/STOP	Home
	<pre></pre>							

Menu	Alarm Silence	Record START/STOP	Size/ Scale	Admit <b>/</b> Discharge	_ Freeze	Key Lock (2 SEC)	Alarm	NIBP Auto Mode	NIBP Start/stop	Home
Small Kov Size: 2 fixed kove 2 user kove										

<Small Key Size: 3 fixed keys, 8 user keys>

By preprogramming the Other Key as user key, 2 pages of user keys can be programmed and pressing the Other Key key allows to switch the user key selection.



Other than the user key setup, the menu key setup allows to change the menu display key configuration, and key mask setup allows to erase the unnecessary keys on the function menu display and system configuration menu display according to user's preference.

### To Set the User Keys

The user keys can be programmed for quick access to the frequently used menu.

1 Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Pre-Set  $\rightarrow$  Monitor Setup  $\rightarrow$ User Key keys. Prev. Disp. Initial DFT Other Key" will be set in same locution for User Key (2/2). The user key setup menu will be displayed. Finne Facerd Facel Reput 
 Orapite Trend
 Tatular Trend
 MBP List
 Recall
 OCR0
 S1
 S1 Transl
 S1 Transl
 Resp. List
 Hean Domain
 Cardeac Ventiletar Night Other Beal Usafance Kave Rec. Phran. Sole HE/PR Br/ann PCHP MIRP Sole See Co-Dee Suspend Alarm HR Alarm Alarm Alarm Auto Suseend Record Record Recorder Setup Config. Optimize University Note Select. Drapky 1 Drapky 2 Depky 8 Key Size ŧ LARGE SMALL 
 I/2
 Surrey Literation
 Surrey Surrey Literation
 Surrey Surrey Literation
 Surrey Literation

#### **2** Select a position to set the user key.



#### **3** Select the function for the user key.



 
 1/2
 Recrit
 Size/ Scale
 Math Recrit
 Freeze
 Key Lock
 Riam
 Output Key

 2/2
 NIBP Rucholse
 Freede Freede
 ST
 Key Lock
 Desky
 Driskey
 Driskey
 Driskey

First, select the location to set the user key. Then, select the function for the user key.

Selecting Other Key key will allow to program 2 pages of user keys.

The Other Key key will be located at the same position for both first and second page.

User Key	Function
Alarm Silence	Silences the alarm for fixed amount of time.
Record START/STOP	Starts/stops the manual recording.
Admit / Disch.	Displays the admit/discharge menu.
Monitor Suspend	Displays the confirmation display whether to suspend monitoring or not.
Freeze	Temporarily stops the waveform trace. By pressing the <u>Record START/STOP</u> key during freeze mode, the waveform in freeze mode can be recorded. When "Freeze Mode Cursor" (monitor setup) is set ON, a cursor will be displayed.
Key Lock	Turns ON/OFF the touch key operation. This function can be used when cleaning the touch screen.
Rapid Discharge	Displays the confirmation screen (latest 12 NIBP list data), and erases patient data, patient information, monitoring condition, etc. (same function as discharge procedure.)
Graphic Trend	Displays the graphic trend.

User Key	Function					
Tabular Trend	Displays the tabular trend.					
NIBP List	Displays the NIBP list.					
Recall	Displays the recall screen.					
OCRG	Displays the OCRG screen.					
ST	Displays the ST measurement menu.					
ST Grap. Trend	Displays the ST graphic trend.					
ST Tab. Trend	Displays the ECG 12-lead ST value in tabular format.					
Resp. List	Displays the respiration measurement value in tabular format.					
Hemodynamic	Hemodynamic calculation menu will be displayed.					
Cardiac Output	CO measurement menu will be displayed.					
Ventilator	P-V, F-V menu will be displayed.					
Night Mode	Turns ON / OFF the night mode.					
Other Bed	Displays the other bed display menu.					
Vigilance List	Displays the Vigilance data list when the Vigilance (or Vigilance CEDV/VigilanceII/Vigileo) oximeter is used.					
Full Disc Wave Rec.	When the CF card for full disclosure waveform recording is inserted in the CF card slot, full disclosure waveform screen will be displayed.					
Parame.	Displays the parameter setup menu.					
Size/Scale	Displays the keys to adjust the size, scale, and the baseline position of the displayed waveform.					
HR/PR Source	Sequentially selects the HR source in the order of $ECG \rightarrow SpO_2 \rightarrow BP1/ART^* \rightarrow Auto \rightarrow ECG.$ * BP1/ART can be selected when ECG/SpO ₂ /BP is set for HR/PR source on the monitor setup menu.					
BP Zero	Performs the zero balance of BP1 to BP5.					
PCWP	Displays the PCWP measurement menu.					
NIBP Auto Mode	Displays the NIBP measurement interval setup menu.					
SpO ₂ Disp ON/OFF	Turns ON/OFF the SpO ₂ display.					
CO ₂ Disp ON/OFF	Turns ON/OFF the CO ₂ display.					
Suspend CO ₂	Suspends the CO ₂ measurement.					
Alarm	Displays the alarm setup menu.					
HR Alarm	Displays the HR/PR alarm setup menu.					
Alarm Auto	Automatically determines the alarm range from the current measurement value.					
Alarm Suspend	Suspends the alarm.					
Record	Displays the recording setup menu.					
Manual Record	Displays the manual record setup menu.					
Recorder Setup	Displays the recorder setup menu.					
Display Config.	Displays the display configuration menu.					
Optimize Display	The confirmation message for optimizing the display will appear. The display will be automatically configured according to the priority set on the "Display Configuration Optimization Setup" (monitor setup).					
Tone/Volume	Displays the tone/volume setup menu.					
Mode Select	Displays the alarm mode/display mode selection menu to set during surgery.					
Display 1 (2, 3)	Switches the home display to Mode 1 (2, 3) preprogrammed on the display mode setup (preset menu).					
Other Key	Switches the first page and second page of the user key.					
OFF	User key will not be displayed.					

### **4** Initialize the user key setup.

Ini	Initial							
CHOLD	2	SEC)						

The factory setting is as follows. User Key 1: Alarm Silence User Key 2: Rec. START/STOP User Key 3: Size / Scale User Key 4: Admit / Discharge Pressing the Initial key for more than 2 seconds will initialize the user key setup to factory default.

User Key 5: Freeze User Key 6: Key Lock User Key 7: Alarm User Key 8: NIBP Auto Mode

# To Configure the Menu Display

# Menu Key Setup

The menu display can be configured for easier use.

Menu			Prev. Disp.
Admit/ Discharge	Alarm	Parameter	
Function	Graphic Trend	Tabular Trend	Recall
	NIBP List	ST	OCRG
	Cardiac Output	Hemodynamics	Ventilator
	Resp. List	ST Graphic Trend	ST Tabular Trend
	Other Bed Display	Night Mode	Vigilance/ Vigileo List
	Full Disc. Waveform Rec.		
System Configuration	Tone/ Volume	Display Config.	Record
	Sweep Speed	Color	Night Mode Setup
	CF Card	Brightness Setup	ST Disp. Lead Setup

Menu Admit/ Discharge	Alarm	Parameter	Prev. Disp.
Function	Graphic Trend	Tabular Trend	Recall
	NIBP List	ST	OCRG
	Cardiac Output		
	Resp.List		
	Other Bed Display	Night Mode	
	Full Disc. Naveform Rec.		
System Configuration	Tone/ Volume	Display Config.	Record
	Sweep Speed	Color	
	CF Card	Brightness Setup	ST Disp. Lead Setup

1 Press the Menu → System Configuration → Pre-Set → Monitor Setup → Menu Setup keys.



First, select the key location.

Select the key for the location.

Select function key for the function key location, and configuration key for the configuration key location.

# **Erasing the Unnecessary Keys**

The unnecessary keys on the function menu and configuration menu can be blanked out.

Function Menu	Prev. Disp.		Function Menu	rev. lisp
Graphic Trend     Tabular Trend       Recall     0006       ST Graphic Trend     ST Tabular ST Tabular       Insh.t Mode     Other Bed Deplay       Full Dec. Haveform Rec.	NBP List ST Resp. List Uentlator Uealance/ Usaleo List		Graphic     Tabular     NIBP List       Trend     Trend     NEP       Recall     Resp. List	
		<function menu=""></function>		
System Configuration	Prev. Disp.		System Configuration	Prev.
				Disp
Display Sweep Config. Speed	_Disp Tone∕ Volume		Display Sweep Tone/ Config. Speed Volume	Disp
Display Config.         Sweep Speed           Record         Color	Tone/		Display Sweep Tone/	Disp
Config. Speed	Tone/ Volume Brightness		Display Sweep Tone/ Config. Speed Dolume Brightness	Disp
Config. Speed Record Color	Tone/ Volume Brightness		Display Sweep Tone/ Config. Speed Dolume Brightness	Disp
Config. Speed Record Color Night Mode Setup Graphic Tabular	Tone/ Uolume Brightness Setup Resp. List		Display Sweep Tone/ Config. Speed Dolume Brightness	Disp.
Confis. Speed Record Color Noth Mode Setup Graphic Trend Setup Recurd Setup Recurd Setup	Tone/ Uolume Brightness Setup Resp. List Setup SI Dep.	>	Display Sweep Tone/ Config. Speed Dolume Brightness	Disp.
Config.         Speed           Record         Color           Night Mode         Statup           Graphic         Trabular           Trend Setup         Trend Setup           Recal Setup         St Graphic           Recal Setup         St Graphic           Set Other         Rev. In	Tone/ <u>Volume</u> <u>Brightness</u> <u>Setup</u> <u>Setup</u> <u>Stopp</u> <u>Lead Setup</u> <u>Ugalance</u>	>	Display Sweep Tone/ Config. Speed Dolume Brightness	Disp.
Config.         Speed           Record         Color           Neiht Rode         Setup           Graphic         Tatulør           Irend Setup         Tatulør           Recall Setup         Sf Graphic           Recall Setup         Sf Graphic           Set Other         Bed ID           BP Über         TBPP Über	Tone/ <u>Volume</u> <u>Brightness</u> <u>Setup</u> <u>Setup</u> <u>Stopp</u> <u>Lead Setup</u> <u>Ugalance</u>	>	Display Sweep Tone/ Config. Speed Dolume Brightness	Disp.
Config.         Speed           Record         Color           Naht Mode         Setup           Graphic         Tratular           Irend Setup         Tratular           Recall Setup         ST Graphic           Set Other         Bed ID           BP User         Izber           Label         Izber	Tone/ <u>Volume</u> <u>Brightness</u> <u>Setup</u> <u>Setup</u> <u>Stopp</u> <u>Lead Setup</u> <u>Ugalance</u>	System Configuration Menu>	Display Sweep Tone/ Config. Speed Dolume Brightness	Disp.

1 Press the Menu → System Configuration → Pre-Set → Monitor Setup → Key Mask keys.

Key Mask (Menu)			Prev. Disp.
Admit/ Discharge	Alarm	Parameter	
Function			
System Configuration			

Select the keys to be displayed on the menu display.

Parameter The key with the LED lit in green will be displayed.

Press the Function key to select the keys to be displayed on the function menu.

Press the System Configuration key to select the keys to be displayed on the function menu.

Key Mask (Function Menu)		Prev. Disp.
Graphic Trend	Tabular Trend	NIBP List
Recall	OCRG	ST ST
ST Graphic Trend	ST Tabular Trend	Resp. List
Hemodynamics	Cardiac Output	Ventilator
Night Mode	Other Bed Display	■ Vigilance/ Vigileo List
Full Disc. Waveform Rec.		

Key Mask (System Config	Menu)	Prev. Disp.
Display Config.	Sweep Speed	Tone/ Volume
Record	Color	Brightness Setup
Night Mode Setup		
Graphic Trend Setup	Tabular Trend Setup	Resp. List Setup
Recall Setup	ST Graphic Trend Setup	<ul> <li>ST Disp.</li> <li>Lead Setup</li> </ul>
Set Other Alarm	Bed ID	Vigilance List Setup
BP User Label	TEMP User Label	
CF Card	<ul> <li>Telemetry Wave Setup</li> </ul>	
Pre-Set		

The green LED will turn off by pressing the key. The key with the LED turned off will not be displayed.

# **Recording Setup**

# Waveform / Numeric Data

On the DS-7200 system, the waveform recording by manual recording, periodic recording, alarm recording, freeze recording, and graphic recording such as graphic trend, tabular trend can be performed. This section describes the procedure for the following recording.

- Manual Recording
- Periodic Recording
- Alarm Recording
- Freeze Recording
- Graphic Recording (Graphic Trend, Tabular Trend, Recall, NIBP List, etc.)

### **Manual Recording**

#### To Start / Stop the Recording

For manual recording (standard recording), pressing the Record START/STOP key will start / stop the recording. Pressing this key during periodic recording, alarm recording, graphic recording, or recall recording will cease the recording in process.



Also, the output recorder status for manual recording will be displayed inside the Record START/STOP key.

Message	Description
None	Normal Operation
PAPER OUT	No recording paper
CASSETTE	Check the cassette.
CHECK?	Other abnormality.



#### Manual Recording Setup

The manual recording can be started from the time the key is pressed, or 8 sec. / 16 sec. prior to the time the key is pressed.

The recording can be set to automatically stop after 24 seconds or continue to record until the Record START/STOP key is pressed again.

1	Press th	ie	Menu	$\rightarrow$	Systen	n Configu	uration	$\rightarrow$	Record	$\rightarrow$	Manual	Record	keys.
	Manual Record				Prev. Disp.	] Th	ne manu	al re	cording se	tup n	nenu will	be display	/ed.
	Rec.Select	Buil	t-in Cent	t.									
	Wave Sel.	ECG 1			 								
		BP1	BP2 BP	3									
	Rec.Dura.	24	sec Cont	t.									
	Delay Time	Nc	one 8se	c 16	Sec								

#### **2** Select the output recorder.

Rec.Select

Built—in Cent.

16sec

Built-in will record on the built-in recorder.

Cent. will record on the 3ch recorder connected to the central monitor.

#### **3** Select the waveform for recording.



Up to 3 waveforms can be selected.

The waveforms will be automatically located at recording.

#### **4** Select the duration for recording.



Select the duration from 24sec or Cont. . 24sec will automatically stop the recording after 24 seconds.

#### **5** Select the delay time for recording.

	·	
Delay Time	 None	8sec

None will start the recording from the time the Rec. START/STOP key is pressed.

8sec, 16sec will start the recording 8 sec. / 16 sec. prior to the time the key is pressed.

	If None is selected for the manual recording dela	ay time,	QR	S classifi	cation will not			
NOTE	NOTE be printed. To record the QRS classification, select 8sec or 16sec for the de							
	time.							

### Alarm Recording

The recording will automatically start at occurrence of a numeric alarm or arrhythmia alarm.

NOTE	<ul> <li>The alarm detection is performed each second, and if more than one alarm occurs at the same time, one data will be stored according to the priority of the alarm factor.</li> <li>Maximum of 3 alarm data can be stored, but if the same or higher priority alarm is newly generated exceeding 3 data, the older recording data will be replaced with the newly generated alarm recording data. The stored data will be erased when recorded.</li> <li>Priority of alarm recording factor ; ASYSTOLE &gt; VF &gt; VT &gt; SLOW VT &gt; TACHY &gt; BRADY &gt; RUN &gt; HR( HR / PR_SpO₂ / PR_IBP ) &gt; APNEA &gt; BP1(or ART) &gt; SpO₂ &gt; NIBP &gt; RR (RR_IMP / RR_CO₂ / RR_VENT)&gt; EtCO₂ &gt; PAUSE &gt; COUPLET &gt; BIGEMINY &gt; TRIGEMINY &gt; FREQUENT &gt; BP2 &gt; BP3 &gt; BP4 &gt; BP5 &gt; ST &gt; TEMP &gt; Tb &gt; InspCO₂</li> <li>If recording on the central monitor recorder, alarm recording and recall recording cannot be performed for the following alarm factor; T3, TACHY, BRADY, SLOW_VT, COUPLET, PAUSE, TRIGEMINY</li> </ul>

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Record  $\rightarrow$  Alarm Record keys.



The alarm recording setup menu will be displayed.

#### 2 Select the output recorder



Built-in will record on the built-in recorder.

Cent. will record on the 3ch recorder connected to the central monitor.

The central monitor recorder will be the one with the smallest central ID.

#### **3** Select the waveform for recording.



Up to 3 waveforms can be selected.

The waveforms will be automatically located at recording. Alarm will record the waveform which generated the alarm.

#### **4** Select the recording factor.



Select the recording factor for alarm recording.

HR will start the alarm recording when a HR or PR alarm is generated.

Other will start the alarm recording when a numeric alarm other than HR and PR alarm is generated.

Arrhy. will start the alarm recording when an arrhythmia alarm is generated.

#### **5** Select the recording duration.

Rec.Dura.

12sec 24sec

Select the recording duration from <u>12sec</u>, <u>24sec</u>. The recording will automatically stop after the selected time.

NOTE	The delay time differs depending on the recording duration.					
	Describer		Dela	Delay Time		
	Recording Duration	Adult	Child	Neonate		
				Numeric Alarm	Arrhythmia Alarm	
	12 sec.	12 sec.	12 sec.	8 sec.	12 sec.	
	24 sec.	16 sec.	16 sec.	16 sec.	16 sec.	
	Ex) If 24sec is selected for adult, the recording will start from the data 16					
	seconds before the alarm and ends 8 seconds after the alarm.					

#### **6** Select the arrhythmia type.

If arrhythmia is selected for the recording factor, select the arrhythmia type.



Pressing the Arrhythmia Factor key will display the arrhythmia selection window for alarm recording.

Selected as alarm recording factor.

Not selected as alarm recording factor.

#### 7 Start the alarm recording.

Alarm Rec.	ON	OFF	ON If al
			in ai

ON will automatically start the recording at alarm occurrence. f alarm recording is not required, select OFF.

NOTE The data at the time of alarm occurrence will be recorded.

### **Periodic Recording**

The recording will be automatically performed with the selected interval. Periodic recording can be performed on the central monitor connected on the wired network system.

NOTE	<ul> <li>If the periodic recording was interrupted due to paper out, etc., only the latest data will be printed when it becomes available again.</li> <li>QRS judgment will not be printed for periodic recording.</li> </ul>
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Record  $\rightarrow$  Periodic Record keys.



The periodic recording setup menu will be displayed.

#### **2** Select the output recorder.



Built-in will record on the built-in recorder.

Cent. will record on the 3ch recorder connected to the central monitor.

The central monitor recorder will be the one with the smallest central ID.

#### **3** Select the waveform for recording.

Wave Sel.	ECG 1	ECG2	 Sp <b>0</b> ₂	RESP	
	 BP1	BP2	BP3	BP4	BP5
	-J AWF	-J AWP			

Up to 3 waveforms can be selected.

The waveforms will be automatically located at recording.

#### **4** Select the periodic interval.



er.	Timer	ſ





Interval recording or timer recording can be selected for periodic recording.

- Timer : Recording will automatically start at the programmed time. Select the time to start recording.
- Inter. : The recording will automatically start with the selected interval. 5min. will start the recording at 10:00, 10:05, ... 10:25. 60 min. will start the recording at 10:00, 11:00, ... 12:00.
#### **5** Select the recording duration.

Rec.Dura.

a. 6sec 12sec 24sec

Select the duration from 6sec, 12sec, 24sec keys. The recording will automatically stop after the selected time.

### 6 Start the periodic recording.

Periodic ON OFF

ON will activate the periodic recording with the selected interval. If periodic recording is not necessary, select OFF.

## 12-Lead Waveform Recording

The monitored 12-lead waveform can be output on the built-in recorder or laser printer.

$\sim$
Reference

To use the laser printer, TCP/IP network setup is necessary. For procedure, refer to "9. Installation TCP/IP Network Connection".

### Recording Setup

The procedure to set the 12-lead waveform recording condition is explained below.

1	Press the	Menu	→	System Configuration	→	Record	→	12-Lead	keys.

The 12-Lead Waveform Recording Setup menu will be displayed. The setup screen differs depending on the laser printer usage.

12-Lead Record Setup Rec. Select Built-in Laser							
Rec. Format							
3 Waves x 4 x 4 + Rh	s 6 Waves x 2	12 Waves					
Position							
Center	Proportional	OFF					
Wa∨e Format	Regular	Reverse					
Recorder Auto Scale	ON	OFF					
Print Calibration	ON	OFF					
Lead Boundary		OFF					

<When laser printer is used>

12-Lead Record Setup Rec. Select	Built-in	_aser						
Rec. Format								
3 Waves × 4 × 6								
Position	Position							
Center	Proportional	OFF						
Wa∨e Format	Regular	Reverse						
Recorder Auto Scale	ON	OFF						
Print Calibration	ON	OFF						

<When laser printer is not used>

### **2** When laser printer is used, select the output printer for 12-lead waveform.

Rec. Select

Built—in Laser

### **3** Select the output format for the laser printer.

Rec. Format

-			-
3 Waves	3 Waves	6 Waves	12 Waves
× 4	$\times$ 4 + Rhy.	×2	1∠ waves
~ 1	<u>× • • • • • • • • • • • • • • • • • • •</u>		

3 Waves $\times$ 4 will print 3 waveforms  $\times$  4 columns. The length of each waveform is 2.5 seconds.



For the example shown on left; the waveforms are as follows; 1st column: I, II, III 2nd column: aVR, aVL, aVF 3rd column: V1, V2, V3 4th column: V4, V5, V6

<u>3 Waves×4+Rhythm</u> will print 3 waveforms  $\times$  4 columns along with10 seconds of rhythm waveform (ECG1 lead on the home display).



For the example shown on left; the waveforms are as follows; 1st column: I, II, III 2nd column: aVR, aVL, aVF 3rd column: V1, V2, V3 4th column: V4, V5, V6

6 Waves $\times 2$  will print 6 waveforms  $\times 2$  column. The length of each waveform is 5 seconds.

	缺形記錄			007/05/02 10	:11 ID:12841			a 709%		
R 60	Sp0+	. 12	20 m		V, 0.03	N 0.10	V. 0.11	A 0.00 V. 0.10	.F. 0.00 # Y. 0.05 #	No 0.01 No 0.01
Lis	l	~	en			1-	-Y~	~	~~	~
Ids	l	<u> </u>	h	h	h	r	7	-	-	-1
In		~	e	A	-1	n	r	r	-r-	-r
lq	-r	~	m	-h-	-y~	m	m	h	h	+
Iq_		-	•••			en	h	h	h	h
I.	h	~	h	h	h	de	h	h	h	h
Sem/s MA	SEAL REC		111-		711.9:1	= * ^4:	on 新住用:0	**		

For the example shown on left, the waveforms are as follows; 1st column: I, II, III, aVR, aVL, aVF 2nd column: V1, V2, V3, V4, V5, V6

12 Waves will print 12 waveforms. The length of each waveform is 10 seconds.

1	2 誘導波	15EM	1	ED-013	2007/	05/02 10:11	ID:12841		종, 가운 로스	フクダ電子		
		Sp01	• 88	20 10		÷.	0.01 - 1	0.01 - 1	0.01 = - W	0.00Vi 0.10 - Vi	0.00V. 0.05 - V	0.01 -
	11	~	d		La	da	da	da	da	dr	ha	-l-
		~	1		In	.la	.lr	In	h	.lr	.hr	l
		~	1		h		h		h	_l_		
ŀ	ar	-	n	-	r	m	m	m	m	m	m	m
Ļ	a.,	~	-									
L	a.t.	~	l		h	h	h	h	h	h	h	l
	e-y	-	r		r	m	r	r	m	~~	m	r
	ny	-	-fr		m	m	mp	m	m	m	m	-y-
L	any	~	÷	<u> </u>	h	in	in	in	in	-j-	-in	-ir
1	e j	1	i		in	in	in	in	in	in	in	-i
ļ	.1	~	L	-	In	th	h	h	In	h	h	V
1	.1	~	1	4	In	th	In	h	In	h	In	1
_								7 MA:05				

For the example shown on left, the waveforms from the top are as follows; I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6

### **4** Select the output format for the built-in recorder.

Rec. Format	3 Waves×4 will print 3 waveforms × 4 columns.
3 Waves 2 Waves	1st column: I, II, III
× 4 × 6	2nd column: aVR, aVL, aVF
	3rd column: V1, V2, V3
	4th column: V4, V5, V6
	2 Waves $\times$ 6 will print 2waveforms $\times$ 6 columns.
	1st column: I, II
	2nd column: III, aVR
	3rd column: aVL, aVF
	4th column: V1, V2,
	5th column: V3, V4
	6th column: V5, V6

#### [Example of "3 Waves×4"]

BED-013 2007/05/02 09:55 10 フクダ電子	0:12841 SEX: AGE:0	ADULT	HR EObpa	RR 208pm	ST (I) ST (avr)	0.01mV 0.00mV 5	ST (III) ST (aVL)	0.01mV ST(00) 0.00mV ST(aVF)	0, 01mV 0, 01mV	ST (V1) 0.03 ST (V4) 0.10	InV ST (V2) InV ST (V5)	0. 10m¥ 0. 05m¥	ST (¥3) ST (¥6)	0. 11mV 0. 01mV
	1n	-l-	-1	-l-				$\neg \gamma$			$\sim$		$\sim$	
	A	h	h	h							~		~	
	1	mhn				L		h	h	h	~		~	L

If the output recorder is built-in recorder, the length of each waveform is fixed as 6 seconds for both output formats.

#### **5** Set the waveform position.

Position

Center Proportional OFF

Center will equalize the printing width of each lead so that the waveform baseline will be at the center. The printing scale of the waveform will be also automatically adjusted.

Proportional will equalize the blank space between each lead to avoid overlapping of the waveforms. The recording scale of the waveform will be also automatically adjusted.

OFF will not adjust the waveform position when printing.

### **6** Set the printing order of the waveforms.

Wa∨e Format	Regular       Regular       Regular       will start printing from the limb leads.         (In the order of I, II, III, aVR, aVL, aVF, V1, V2, V3, V4,
	V5, V6) Reverse will start printing from the chest leads. (In the order of V1, V2, V3, V4, V5,V6, I, II, III, aVR, aVL, aVF)
	and the sector of the Breat the sector

**7** Select whether or not to automatically adjust the scale.

Recorder Auto Scale	Γ ON	OFF	When position adjustment is OFF, select whether or
			not to automatically adjust the scale.
			ON will automatically adjust the printing scale.
			OFF will not automatically adjust the printing scale
			and prints with the displayed scale.

NOTE	The printing scale will be adjusted in the range of $\times 1$ , $\times 1/2$ , $\times 1/4$ . It will not be adjusted to $\times 2$ or $\times 4$ even if the amplitude is small.
------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### 8 Select whether or not to print the calibration waveform.

Print	Calibration

ON will print the calibration waveform. OFF will not print the calibration waveform.

### **9** Select whether or not to print the lead boundary when recording on the laser printer.

Lead Boundary	OFF	ON will print the lead boundary between the leads.
	 	OFF will not print the lead boundary.

### To Print the 12-Lead Waveform

When the display mode is "12-Lead", 12LEAD REC key will be displayed on the home display. Press this key to start printing.

### **1** Set the display mode to "12-Lead".

12LEAD REC key will be displayed on the home display.



When laser printer is set as the output recorder, laser printer icon will be displayed inside the 12LEAD REC key.



### **2** Print the 12-lead waveform.

Press the 12LEAD REC key to start recording. For the laser printer, printing can be cancelled by pressing the Cancel key.

- · The waveform for each lead will be printed in the same phase.
- · The printing duration of the waveforms for each format are as follows.

Output Recorder	Output Format	Printing Duration	
Built-in Recorder	3 Waves×4	6 sec.	
Built-III Recorder	2 Waves×6	0 360.	
	3 Waves×4	2.5 sec.	
Laser Printer	3 Waves×4+Rhythm	2.0 560.	
	6 Waves ×2	5 sec.	
	12 Waves	10 sec.	

### **Recorder Setup**

The recording condition common to manual, periodic, alarm recording, and output recorder for graphic recording can be set.

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Record  $\rightarrow$  Setup keys.

Recorder Setup		Prev. Disp.	
QRS Classification	ON	OFF	
Graphic Recording	Built-in	Cent.	
Recall Recording	Graphic Recording	Manual Recording	
Recorder			Th
Feed	Both Top	End OFF	
Speed	50mm/s 25	mm/s	
Print Calibration	Tob	 Each Page	
	OFF		
-			

The recorder operation setup menu will be displayed.

### 2 Select ON/OFF for QRS classification symbol recording.

QRS Classification

OFF ON will waveform.

ON will record the QRS classification symbol on to the ECG waveform.

Symbol	Description
N (Normal)	Normal QRS beat
V (VPC)	Ventricular Extrasystole
S (SVPC)	Supraventricular Extrasystole
P (Pacing Beat)	Pacing beat
F (Fusion Beat)	Fusion beat of pacing and spontaneous beat.
? (Undetermined beat)	Learning arrhythmia, or beat not matching the pattern

The QRS classification symbol cannot be recorded for the manual recording without delay time and periodic recording. To record the QRS classification symbol, set the delay time to 8 seconds or 16 seconds for manual recording.
 The "S" (QRS symbol) will be printed as "N" on the central recorder.

### **3** Select the output recorder for the graphic recording.

Cent.

Graphic Recording

Feed

Built-in will record on the built-in recorder.

Cent. will record on the central monitor recorder. The central monitor recorder of the smallest ID will be used.

### **4** Select the output recorder for the recall recording.

Built—in

Recall Recording Graphic Recording Manual Recording will record on the recorder selected at procedure 3. Manual Recording will record on the recorder selected for manual recording.

5 Set the paper feed operation for the recorder.

Both Top End OFF

Both will start the recording from the perforation, and feeds the paper to the next perforation, so that the paper can be easily cut off. Top will start the recording from the perforation, and will not feed the paper after recording.

End will start the recording from the position where the previous recording ended, and feeds the paper to the next perforation after recording so that the paper can be easily cut off.

OFF will start the recording from the position where the previous recording ended, and will not feed the paper after recording.



### **Freeze Recording**

The waveform display can be frozen and recorded from 12 seconds prior to the frozen point. The recording duration is fixed as 12 seconds.

To freeze the waveform display, the Freeze key needs to be assigned as a user key.

### **1** Freeze the waveform display.

Press the Freeze key on the user key.



### **2** Start freeze recording.

Press the Record START/STOP key to record the displayed waveform. The freeze recording will be printed on the built-in recorder. The waveform set for the manual recording will be recorded.

# Graphic Recording (Graphic/Tabular Trend, etc.)

Graphic trend recording, tabular trend recording, NIBP list recording, recall recording, and Vigilance list recording can be performed.

Reference

Refer to respective section of "7. Function " for recording procedure.

### **Graphic Trend**

3 (three) parameters from the following can be selected for graphic recording.

Parameter	Description
HR	HR, PR (PR_SpO ₂ ), PR (PR_IBP)
ST	ST(I), ST(II), ST(III), ST(aVR), ST(aVL), ST(aVF), ST(V), ST(V2), ST(V3), ST(V4), ST(V5), ST(V6)
VPC	VPC beats
BP1	BP1 (SYS / Mean / DIA)
BP2	BP2 (SYS / Mean / DIA)
BP3	BP3 (SYS / Mean / DIA)
BP4	BP4 (SYS / Mean / DIA)
BP5	BP5 (SYS / Mean / DIA)
PDP	Peak Diastolic Pressure

Parameter	Description
CPP	Cerebral Perfusion Pressure
NIBP	NIBP (SYS / Mean / DIA)
SpO ₂	SpO ₂ value
TEMP1, TEMP2	TEMP1, TEMP2
TEMP3	TEMP3
Tb	Blood Temperature
RR	Impedance Resp. (RR), CO ₂ Resp. (RR_CO ₂ ), Ventilator Resp. (RR_VENT)
APNEA	Apnea Time (Impedance, CO ₂ , Ventilator)
CO ₂	EtCO ₂ /InspCO ₂
SvO ₂	Mixed Venous Oxygen Saturation
ScvO ₂	Central Venous Oxygen Saturation
CCO	Continuous Cardiac Output
CCI	Continuous Cardiac Index
BT	Blood Temperature (Vigilance Data)
BIS	BIS Monitor Data
EVENT1	ASYSTOLE, VF, VT, SLOW_VT, RUN, BIGEMINY
EVENT2	TRIGEMINY, PAUSE, COUPLET, TACHY, BRADY, FREQUENT

### **Tabular Trend**

17 parameters from the following can be displayed / recorded as tabular trend.

Parameter	Description
HR	HR, PR (PR_SpO ₂ ), PR (PR_IBP)
ST	ST(I), ST(II), ST(III), ST (aVR), ST (aVL), ST (aVF), ST (V), ST (V2), ST (V3), ST (V4), ST (V5), ST (V6)
VPC	VPC beats
BP	BP (SYS): BP1-S, BP2-S, BP3-S, BP4-S, BP5-S, BP (Mean): BP1-M, BP2-M, BP3-M, BP4-M, BP5-M BP (DIA): BP1-D, BP2-D, BP3-D, BP4-D, BP5-D
PDP	Peak Diastolic Pressure
CPP	Cerebral Perfusion Pressure
PCWP	Pulmonary Capillary Wedge Pressure
NIBP	NIBP (SYS / Mean / DIA) NIBP-S, NIBP-M, NIBP-D
SpO ₂	SpO ₂ value
TEMP	T1, T2, T3
Tb	Blood Temperature
RR	Impedance Resp. (RR), CO ₂ Resp. (RR_CO ₂ ), Ventilator Resp. (RR_VENT)
APNEA	Apnea Time (Impedance, CO ₂ , Ventilator)
EtCO ₂	EtCO ₂ value
InspCO ₂	InspCO ₂ value
SvO ₂	Mixed Venous Oxygen Saturation
ScvO ₂	Central Venous Oxygen Saturation
CCO	Continuous Cardiac Output
CCI	Continuous Cardiac Index
BT	Blood Temperature (Vigilance Data)
BIS	BIS, BIS_SQI, BIS_EMG, BIS_SR (BIS Monitor Data)

#### **NIBP List**

The numeric data for the following parameters can be recorded.

Parameter	Description
HR	HR
PR_SpO ₂	PR (PR_SpO ₂ )
SpO ₂	SpO ₂ value
NIBP	NIBP (SYS / Mean / DIA)

#### Recall

The waveform and numeric data at alarm occurrence will be recorded with the following setup.

Parameter	Description
Waveform (Max. 2)	ECG1, ECG2, BP1 to BP5, SpO ₂ , CO ₂ , RESP
Numeric	HR, ST, NIBP, BP1 to BP5, RR, APNEA, SpO ₂ , T1 to T3, Tb, CO ₂ ,
Arrhythmia Factor	ASYSTOLE, VF, VT, SLOW VT, RUN, BIGEMINY, TRIGEMINY, PAUSE,
	COUPLET, TACHY, BRADY, FREQUENT
<b>Recording Duration</b>	12 sec.

### **Respiration List**

17 parameters from the following can be selected for list recording.

Parameter	Description
RR_IMP	Impedance RR
RR_CO ₂	CO ₂ RR
RR_VENT	Ventilator RR
SpO ₂	SpO ₂ value
APNEA	Apnea time (Impedance, CO ₂ , Ventilator)
E-TV	Expiratory Tidal Volume
I-TV	Inspiratory Tidal Volume
MV	Minute Ventilation
SMV	Spontaneous Minute Ventilation
P_PEAK	Maximum Airway Pressure
P_PAUSE	Pause Airway Pressure
PEEP	Peak End Expiratory Pressure
P_MEAN	Mean Airway Pressure
D-RES	Dynamic Resistance
S-RES	Static Resistance
D-COMP	Dynamic Compliance
S-COMP	Static Compliance
P_Min	Minimum Airway Pressure
S_RR	Spontaneous Respiration
FIO ₂	Fractional Concentration of Inspiratory Gas
EtCO ₂	EtCO ₂ value
SvO ₂	Mixed Venous Oxygen Saturation
ScvO ₂	Central Venous Oxygen Saturation
CCO	Continuous Cardiac Output
CCI	Continuous Cardiac Index

NOTE	The displayed name of the parameter differs depending on the used ventilator. Also, some parameters may not be displayed depending on the used ventilator. For example, Resistance (Insp/Exp) and Compliance will not be displayed for the SV-300 and Servo-i.
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Vigilance List

17 parameters from the following can be selected for list recording.

Parameter	Description
SvO ₂	Mixed Venous Oxygen Saturation
ScvO ₂	Central Venous Oxygen Saturation
SaO ₂	Arterial Oxygen Saturation
O ₂ EI	Oxygen Uptake Index
DO ₂	Oxygen Transport
VO ₂	Oxygen Consumption
SV	Stroke Volume
SV_STAT	Stroke Volume (STAT Mode)
SVI	Stroke Volume Index
SVI_STAT	Stroke Volume Index (STAT Mode)
HR	Heart Rate
MAP	Mean Arterial Pressure
CVP	Central Venous Pressure
CCO	Continuous Cardiac Output
CCO_STAT	Continuous Cardiac Output (STAT Mode)
CCI	Continuous Cardiac Index
CCI_STAT	Continuous Cardiac Index (STAT Mode)
SVR	Systemic Vascular Resistance
SVRI	Systemic Vascular Resistance Index
B_Temp	Blood Temperature
EF	Ejection Fraction
EF_STAT	Ejection Fraction (STAT Mode)
EDV	End-Diastolic Volume
EDV_STAT	End-Diastolic Volume (STAT Mode)
EDVI	End-Diastolic Volume Index
EDVI_STAT	End-Diastolic Volume Index (STAT Mode)
ESV	End-Systolic Volume
ESVI	End-Systolic Volume Index
SVV	Stroke Volume Variance

# **Volume Setup**

# Pulse Tone, Alarm Sound, etc.

This menu allows volume setup of the pulse tone, alarm sound, key sound and other bed alarm sound. The ON/OFF of ventilator alarm sound can be also selected.

#### **1** Press the Menu $\rightarrow$ Tone / Volume keys.

The tone/volume setup screen will be displayed.

#### **2** On the first page, set the tone/volume for the alarm sound.



The tone and volume can be set for each alarm level. However, if <u>IEC</u> is selected for "Alarm System" under the "Monitor Setup" menu, tone for level 2 and level 3 cannot be set. The tone setting for level 1 will be applied.



Reference

For the night mode setup, refer to "8. System Configuration Night Mode"

# **3** Press the Page Down key and adjust the tone/volume of the other bed alarm sound and other sound.



Adjust the tone/volume of the pulse, key, other bed, and other sound.



# **Color/Brightness Setup**

This menu allows the setup of the colors of numeric data / waveform and brightness of the display.

# Color Setup (Numeric Data, Waveform)

The displayed color for each parameter can be selected from the 32-color palette.

### **1** Press the Menu $\rightarrow$ System Configuration $\rightarrow$ Color keys.



<Color Setup Menu>

### **2** Select the color from the 32-color palette.



Press the key of the desired color.

### **3** Assign the selected color to the parameter.



Press the parameter key to assign the selected color. The selected color for the parameter will be applied to the waveform, numeric data, graphic trend, and tabular trend.

## **Brightness Setup**

The brightness of the display can be adjusted.



# **Sweep Speed**

# Waveform Display Speed/Time

This section describes the procedure to set up the sweep speed of the waveform display.

The sweep speed can be set separately for ECG/BP waveform and respiration waveform.



Display Unit	Display Time		
Display Config. Selection	Standard	Ext.1, Ext.2	Enlarged
25mm/s	7.9sec.	5.5sec.	10.2sec.
12.5mm/s	15.8sec.	11.0sec.	20.4sec.
6.25mm/s	31.6sec.	22.0sec.	40.8sec.

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Sweep Speed keys.



**2** Set the sweep speed for ECG, BP,  $SpO_2$  waveform.

ECG, BP, SpO2 (mm/s) 12.5 6.25 25

**3** Select the sweep speed for respiration waveform.

RESP (mm/s)

6.25	<b>12.5</b>	<b>2</b> 5
------	-------------	------------

# **Remote Control Setup**

This section explains the setup procedure on how to use the optional remote control unit, CF-700.

### About the Remote Control Unit



#### **Monitor Selection Dial**

The remote control bed ID to control can be selected. The remote control bed ID should be preprogrammed on each bedside monitor. Maximum of 8 monitors can be controlled with one remote control unit.

#### **Check Key**

Displays the remote control bed ID on the bedside monitor. Also, pressing one of the function keys (ex. NIBP START/STOP key) will display a message on the monitor which will ask you to press this Check key.

### Alarm Silence Key

This key has the same function as the Alarm Silence key on the bedside monitor.

#### Home Key

This key functions the same as the Home key on the bedside monitor.

#### User Key (F1 to F8)

The function for each key can be assigned on the bedside monitor.

Key	Default	Function
F1	ECG1 Size	Switches the ECG1 size each time the key is pressed.
		$\times 1/4 \rightarrow \times 1/2 \rightarrow \times 1 \rightarrow \times 2 \rightarrow \times 4 \rightarrow \times 1/4$
F2	ECG1 Lead	Switches the ECG1 lead each time the key is pressed.
		3-electrode: I→II→III→I
		4-electrode: I→II→III→aVR→aVL→aVF→I
		5-electrode: I→II→III→aVR→aVL→aVF→V→I
		10-electrode: $I \rightarrow III \rightarrow III \rightarrow aVR \rightarrow aVL \rightarrow aVF \rightarrow V1 \rightarrow V2 \rightarrow V3 \rightarrow$
		V4→V5→V6→I
F3	NIBP START/STOP	Starts/stops the NIBP measurement.
		Pressing this key will display a message on the monitor to press
		the "Check" key.
		When the "Check" key is pressed, the measurement will start.
		To cancel the process, press the "Home" key.
		Pressing this key during the measurement will stop the
		measurement.
F4	Record START/STOP	Starts/stops the manual recording.
F5	Night Mode	Sets ON/OFF of the Night Mode.
F6	Tabular Trend	Displays the tabular trend data.
F7	NIBP List	Displays the NIBP list.
F8	BP Zero	Starts the BP zeroing.

# **Remote Control Setup**

Functions can be assigned to 8 user keys on the remote control unit.

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  Pre-Set  $\rightarrow$  Monitor Setup  $\rightarrow$  R.C. Setup keys.

The remote control unit setup menu will be displayed.



### **2** Set the remote control bed ID.

Set the bed ID that will respond to the monitor selection dial on the remote control unit.



### **3** Select the key location.



Key Location on the Remote Control Unit

#### **4** Select the function.



Select a function to be assigned for the selected key location.

Pressing the function key will assign the function to the selected key location.

Use the 1/2, 2/2 keys to switch the page for function selection.

Switch page



Select R.C. OFF if not using the remote control function.

Select the key location of the remote control unit to change or assign a function.

The "F1" key on the remote control unit is the same key as the "1" key on the remote control setup menu.

### Functions that can be assigned to the User Keys

Function	Description
ECG1 Size	Switches the ECG1 (ECG2) size each time the key is pressed.
ECG2 Size	$\times 1/4 \rightarrow \times 1/2 \rightarrow \times 1 \rightarrow \times 2 \rightarrow \times 4 \rightarrow \times 1/4$
ECG1 Lead	Switches the ECG1 (ECG2) lead each time the key is pressed.
ECG2 Lead	3-electrode: I→II→III→I
	4-electrode: I→II→III→aVR→aVL→aVF→I
	5-electrode: I→II→III→aVR→aVL→aVF→V→I
	10-electrode: $I \rightarrow II \rightarrow III \rightarrow aVR \rightarrow aVL \rightarrow aVF \rightarrow V1 \rightarrow V2 \rightarrow V3 \rightarrow V4 \rightarrow V5 \rightarrow V6 \rightarrow I$
ECG Auto Size	Automatically adjusts the ECG size to 10mm. This function is effective only
	when the key is pressed.
BP1 to BP5 Scale	Switches the BP1 to BP5 scale each time the key is pressed.
	20→50→75→100→150→200→250→300→20 (mmHg)
	4→8→12→16→20→24→32→40 (kPa)
PCWP	If the BP label is PAP, PCWP input menu will be displayed.
BP Zero	Starts zeroing for all BP.
	It will not function unless the transducers for all BP is opened to air.
NIBP START/STOP	Starts/stops the NIBP measurement.
	Pressing this key will display a message on the monitor to press the "Check"
	key. When the "Check" key is pressed, the measurement will start.
	To cancel the process, press the "Home" key.
<u> </u>	Pressing this key during the measurement will stop the measurement.
Record	Starts/stops the manual recording. The recording duration set on the manual
START/STOP	recording setup menu will be applied.
Monitor Resume	Resumes monitoring when the monitoring is suspended.
Alarm Suspend	Suspends the alarm for fixed amount of time.
	When pressed during the alarm is suspended, the alarm function will
Freeze	resume. Temporarily stops the waveform trace. Pressing the key again will resume
116626	the waveform trace.
Graphic Trend	Displays the graphic trend.
Tabular Trend	Displays the tabular trend.
NIBP List	Displays the NIBP list.
OCRG	Displays the OCRG.
Recall	Displays the recall data.
ST	Displays the ST measurement display.
ST Graphic Trend	Displays the ST graphic trend.
ST Tabular Trend	Displays the ST tabular trend.
Cardiac Output	CO measurement menu will be displayed. This key will not start the CO
ourdido output	measurement.
Hemodynamic	Hemodynamic calculation menu will be displayed.
VENT (P-V)	Displays P-V loop on the ventilator display.
VENT (F-V)	Displays F-V loop on the ventilator display.
VENT (Numeric)	Displays numeric data on the ventilator display.
Resp. List	Displays respiration list.
Vigilance List	Displays Vigilance list/
Night Mode	Turns ON/OFF the Night Mode.
Config. Enlarge	Switches the display configuration to "Enlarge" mode. Pressing the key
3	again will return the display to previous configuration.
Config. 12LEAD	Switches the display configuration to "12-lead" mode. Pressing the key again
Config. 12LEAD	Switches the display configuration to "12-lead" mode. Pressing the key again will return the display to previous configuration. Turns OFF the key operation.

### **5** Check the setting.

▲ CAUTION	<ul> <li>Do not set the same remote control bed ID to more than one monitors on the same floor. Otherwise, it may cause to remote control more than one monitors at the same time.</li> <li>After the remote control setup, check that the remote control unit is properly operating.</li> </ul>
-----------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **Section Setup**

One remote control unit can control a maximum of 8 monitors. When controlling more than 8 monitors, it is necessary to set a section for the remote control. If the section is not properly set, the remote control may unintentionally control 2 monitors at the same time. The set section number must be the same between the remote control unit and the monitor.



# **Alarm Pole Setup**

# Notifying the Alarm by Light

This monitor is equipped with an alarm pole.

The alarm generation can be indicated with the alarm pole flashing and can be seen from distance. The alarm pole can be also synchronized with HR.

The alarm pole function setup can be performed for (1) flashing at alarm generation, and for (2) flashing synchronized with HR.

When <u>IEC</u> is selected for "Alarm System" on the Monitor Setup menu, not all items can be set on the Alarm Pole Setup menu.

When FUKUDA DENSHI is set	When IEC is set
- Setting is possible.	Setting is not possible. (Fixed to "ON")
	Setting is not possible. (Fixed to "Level 1, 2 and 3")
	Setting is not possible. (Fixed to "ON")
	Setting is not possible.
	Setting is possible.
	_

**CAUTION** The alarm priority is high for level 1 (life threatening alarm), medium for level 2 (cautionary alarm), and low for level 3 (treatment needed alarm).

### Flashing the Alarm Pole at Alarm Generation

The alarm level to flash the alarm pole and the flash pattern for each alarm level can be set.

Press theMenu $\rightarrow$ SystemAlarm Polekeys.	$\begin{array}{c c c c c c c c c c c c c c c c c c c $
Alam Pole Setue	
Sync. with Alarm ON OFF	
Uentilator Alarn         ON         OFF           Pattern Setup         >When "Sync. with HR" is 0N, pattern is fixed.	The alarm pole setup menu will be displayed.
Sync. with HR ON OFF *"Sync. with Alara" will have priority if set to ON.	
Select whether the alarm pole	e should flash at alarm generation.
Sync. with Alarm	To flash at alarm generation, select ON. To not flash at alarm generation, select OFF.
Select the alarm type to flash	the alarm pole.
Alarm Type Level1 Level1 and 2 Level1, 2 a Ventilator Alarm ON OFF	Select the alarm level to flash the alarm pole from Level 1, Level 1 and 2, and Level 1, 2, and 3. Whether to flash the alarm pole at ventilator alarm generation can be also selected. ON will flash the alarm pole ventilator alarm generation.
CAUTION By selectin	level shown in next page is the standard level set by Fukuda Denshi. Ig User for "Alarm Level" in the 4th page of Monitor Setup menu, rogrammed alarm level will be applied.

### Alarm Level 1 (Life Threatening Alarm)

Parameter	Message
HR	"Lower HR alarm"
	"Upper HR alarm"
	"Lower PR alarm"
PR (SpO ₂ , BP)	"Upper PR alarm"
BP (BP1/ART)	"Lower BP1 alarm" or "Lower ART alarm"
	"Upper BP1 alarm" or "Upper ART alarm"
SpO ₂	"Lower SpO ₂ alarm"
	"Upper SpO₂ alarm"
Respiration	"Apnea alarm"
(Impedance, CO ₂ , ventilator)	"Lower RR alarm"
	"Upper RR alarm"
NIBP	"Lower NIBP alarm"
	"Upper NIBP alarm"
CO ₂	"Upper EtCO ₂ alarm"
	"Lower EtCO ₂ alarm"
	"ASYSTOLE"
	"VF"
Arrhythmia	"VT" "SLOW VT"
Arrhythmia	"TACHY"
	"BRADY"
	"RUN"

### Alarm Level 2 (Cautionary Alarm)

Parameter	Message
BP (BP2 to 5)	"Lower BP* alarm" or "Lower (label) alarm"
	"Upper BP* alarm" or "Upper (label) alarm"
OT1 to 10	"Lower ST* alarm"
ST1 to 12	"Upper ST* alarm"
Temperature	"Upper TEMP $*$ alarm" or "Lower (label) alarm"
(TEMP1 to 3)	"Lower TEMP* alarm" or "Upper (label) alarm"
CO ₂	"Upper InspCO ₂ alarm"
	"PAUSE"
	"COUPLET"
	"BIGEMINY"
Arrhythmia	"TRIGEMINY"
	"FREQUENT"
	"Cannot analyze"
ECG, Impedance RESP	"Check Electrodes"
SpO ₂	"Check SpO ₂ Sensor"
(DS-7210, DS-7210M)	"Replace SpO ₂ Sensor"
SpO ₂ (DS-7210)	"No pulse detect"
6-0	"Unknown SpO₂ Sensor"
SpO ₂ (DS-7210M)	"SpO ₂ Low Perfusion"
	"SpO ₂ Pulse search"
	"Check Sample Line"
CO ₂	"Check CO ₂ Exhaust Port"
(MGU-722)	"Check CO ₂ unit"
	"CO ₂ Cal. Required"

Parameter	Message
	"ECG Disconnected"
	"BP* Disconnected"
	"SpO ₂ Disconnected"
Connector	"T* Disconnected"
Connector	"CO ₂ Disconnected"
	"CO Disconnected"
	"Multiport* Disconnected"
	"DS-LANII Disconnected" or "DS-LANIII Disconnected"
Others	"Check Backup Battery"
	"Check Equip. Config. (CO ₂ )"
Others	"Charge the battery."
	"Check Memory Card"

 $\star$  indicates the channel no. of BP and TEMP.

#### Alarm Level 3 (Treatment Needed Alarm)

Parameter	Message
NIBP	"Check NIBP Cuff, Air Hose"
Impedance RESP	"CVA detected"
SpO ₂ (DS-7210)	"No pulse detect"
	"SpO ₂ Low perfusion"
	"SpO ₂ Pulse Search"
SpO ₂ (DS-7210M)	"SpO ₂ Interference Detected"
	"Check SpO ₂ Sensor"
	"SpO ₂ Low Signal IQ"
ECG	"Pacing detection error"

#### Ventilator Alarm

Parameter	Message	
Ventileter	"Vent. Alarm"	
Ventilator	"Vent. Invalid"	

### **4** Set the alarm pole flash pattern for each alarm level.



Press the Pattern Setup key to display the alarm pole flash pattern setup menu.

Assign the flash pattern for each alarm level.

1) First, select the flash pattern from Pattern 1 to Pattern 10.

The alarm pole consists of 3 blocks.

(Left: Red, Middle: Red/Orange/Green, Right: Red) The ways these 3 blocks flash are different for each pattern.

By pressing the Pattern Test key, the flashing of the selected pattern can be verified.

2) Next, press the alarm level key to assign the selected flash pattern.

#### **Flash Pattern**

("*" indicates that light is OFF.)

Pattern	Flash
Pattern 1	(Red, Red, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Red, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Red, Red)
Pattern 2	(Red, Orange, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Orange, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Orange, Red)
Pattern 3	(Red, Green, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Green, Red) $\rightarrow$ (***) $\rightarrow$ (Red, Green, Red)
Pattern 4	$(*, Red, *) \rightarrow (* * *) \rightarrow (*, Red, *) \rightarrow (* * *) \rightarrow (*, Red, *)$
Pattern 5	$(*, Orange, *) \rightarrow (* * *) \rightarrow (*, Orange, *) \rightarrow (* * *) \rightarrow (*, Orange, *)$
Pattern 6	$(*, Green, *) \rightarrow (* * *) \rightarrow (*, Green, *) \rightarrow (* * *) \rightarrow (*, Green, *)$
Pattern 7	$(Red, Red, *) \rightarrow (* * *) \rightarrow (*, Red, Red) \rightarrow (* * *) \rightarrow (Red, Red, *)$
Pattern 8	(Red, Orange, $*$ ) $\rightarrow$ ( $*$ $*$ ) $\rightarrow$ ( $*$ , Orange, Red) $\rightarrow$ ( $*$ $*$ ) $\rightarrow$ (Red, Orange, $*$ )
Pattern 9	(Red, Green, $*$ ) $\rightarrow$ ( $*$ $*$ ) $\rightarrow$ ( $*$ , Green, Red) $\rightarrow$ ( $*$ $*$ ) $\rightarrow$ (Red, Green, $*$ )
Pattern 10	$(Red, *, *) \rightarrow (* * *) \rightarrow (*, Red, *) \rightarrow (* * *) \rightarrow (*, *, Red)$

### **Default Setting**

Alarm Level	Pattern	Flash
Level 1	Pattern 1	$\begin{array}{ccc} (Red,  Red,  Red) \rightarrow (\ast \ast \ast) \rightarrow (Red,  Red,  Red) \rightarrow (\ast \ast \ast) \rightarrow (Red,  Red,  Red) \\ 480 \mathrm{ms} & 480 \mathrm{ms} & 480 \mathrm{ms} & 480 \mathrm{ms} \\ \end{array}$
Level 2	Pattern 10	(Red**)→(***)→(*Red*)→(***)→(**Red) 480ms 240ms 480ms 240ms 480ms
Level 3	Pattern 4	(*Red*)→(***)→(*Red*)→(***)→(*Red*) 480ms 480ms 480ms 480ms 480ms

### Flash Pattern when IEC is selected for "Alarm System" (Fixed)

Alarm Level	Flash			
Level 1	$(\text{Red, Red, Red}) \rightarrow (* * *) \rightarrow (\text{Red, Red, Red}) \rightarrow (* * *) \rightarrow (\text{Red, Red, Red})$ 320ms 320ms 320ms 320ms 320ms 320ms			
Level 2	$(*Orange*) \rightarrow (***) \rightarrow (*Orange*) \rightarrow (***) \rightarrow (*Orange*)$ 800ms 800ms 800ms 800ms 800ms 800ms			
Level 3	(*Orange*) Lights constantly.			

# Flashing the Alarm Pole Synchronized with HR

The alarm pole can be flashed synchronizing with HR.

### **1** Select whether to flash the alarm pole synchronizing to the heartbeat.

Sync. with HR

ON	OFF

To synchronize with HR, select ON. The middle LED will flash in green synchronized with HR. To not synchronize with HR, select OFF.

ΝΟΤΕ	<ul> <li>When asystole is generated, the green LED at middle part of the alarm pole will remain lighted. However, if "Sync. with Alarm" is set to ON, the fixed flash pattern for the asystole alarm (Level 1: Pattern 1) will have priority.</li> <li>If the parameter other than ECG is selected as "HR/PR Alarm Source" on ECG (SpO₂, BP) setup menu, the alarm pole will not flash synchronizing with the pulse.</li> </ul>
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# **Alarm Setup**

# To Set the Alarm Condition

This section explains the setup of the alarm suspension and upper / lower alarm limit. On the alarm setup menu, ON/suspend of system alarm, ON/OFF and upper / lower alarm limit of each parameter can be set.

Also, 5 patterns of alarm setting can be programmed using the alarm mode setup function. By preprogramming the alarm setting to each alarm mode, the alarm setups at admittance of patient can be simplified by just selecting one of the alarm modes.



For alarm mode setup procedure, refer to "8. System Configuration Alarm Mode Setup To Program the Alarm Mode"

# To Set the System Alarm (ON or Suspend)

The system alarm can be set to ON or suspend, but it cannot be turned OFF.

<b>A</b> WARNING	• When the system alarm is suspended, all the alarm will be suspended even if the parameter alarm is set to ON. Also, the alarm event will not be stored as recall.
	<ul> <li>If the upper/lower alarm limit of the parameter is set to OFF, or arrhythmia alarm is set to OFF, alarm will not function even if the system alarm is set to ON. Pay attention when setting them OFF.</li> </ul>

### **1** Press the Menu $\rightarrow$ Alarm keys.



The alarm setup menu will be displayed.

**2** When the **Suspend** key LED is extinguished, the system alarm is ON.

The alarm ON/OFF function and alarm limit for each parameter is effective.



### **3** When the Suspend key LED is lighted, the system alarm is suspended.

Pressing the Suspend key when the LED is extinguished will light the LED and temporarily suspends the system alarm.

While the alarm is suspended, "Alarm Suspend (***sec)" message will be displayed. (***sec.) indicates the remaining time.

The alarm will turn ON when the suspended time completes.

Supend		Alarm Suspend Message
FUKUDA DENSHI	Adult	02/01 16:51 M Alarm Suspend(115sec.) HR ♥Av.

### To Silence the System Alarm

The alarm sound can be silenced for fixed amount of time. This setting will not affect the alarm message. If the alarm cause still remains at completion of the silence duration, the alarm sound will generate again. Also, if another alarm with the same or higher priority occurs during the alarm silence duration, the alarm sound for the new alarm will generate.



### Precautions about Silencing the Alarm

- Alarm silence function is effective for each parameter. If an alarm condition for the selected parameter is resolved for a moment but is generated again during the alarm silence time, the alarm will remain silenced. The recall and alarm recording will not function at this time.
- If another alarm with the lower priority occurs during the alarm silence time, alarm sound will not generate. The recall and alarm recording will function.
- If the <u>Alarm Silence</u> key is pressed for the alarm of another parameter which occurred during the alarm silence time, the alarm silence time for the first alarm will not be extended.
- The alarm silence state for all parameters will cease in the event of any of the following.
  - $\boldsymbol{\cdot}$  When the main power is turned ON.
  - · When the system alarm status (ON / suspend) is changed.
  - When Resume All Al. Sound key is pressed on the alarm setup menu.
  - When monitoring is suspended on the patient admit / discharge menu.
  - · When the alarm mode is changed on the patient admit / discharge menu.
  - When the patient has been discharged.
- The alarm silence state for each parameter will cease in the event of any of the following.
  - · When the alarm silence duration for the parameter is completed.
    - · When automatic alarm is selected for the parameter.
    - · When the alarm is turned OFF for the parameter.
- If <u>Linked to each new occurrence</u> is selected for "Status Alarm Control" in the alarm setup menu, the status alarm sound will not resume after the alarm silence time unless a new status alarm generates.

### To Cancel "Alarm Silence"

The alarm silence state can be cancelled by pressing the Resume All Al. Sound on the Alarm Setup menu. The alarm silence state for all parameters and equipment status will be cancelled and alarm sound will resume if alarm factor exists.

Alarm	■ Suspend	Mode Select	Alar Reco		
HR/PR _{Recall}	0	40-120	ST 🕅 Recall	STIL	F
Arrhy. _{Recall}	Asystole Run	UF UT : Tachy Bra	Slow UT ady		
NIBP Recall	S D	80-180 0FF-0FF			
SpO ₂ _{Recall}	Sec	90-0FF 0FF	RESP _{Recall}	rr Apnea	5- 30 15
BP <b>1</b> _{Recall}	S D M	80-180 0FF-0FF 0FF-0FF	BP <b>2</b> _{Recall}	XAS ⊓	OFF OFF OFF
T <b>1 )</b> Recall	X	OFF	T <b>2</b> _{Recall}	X	OFF
BP Alar	m Tem	p Alarm	S	etup	

### Alarm Setup for Each Parameter

The alarm for each parameter can be turned ON or OFF, and upper and lower alarm limit can be set. The numeric data alarm generates when a value exceeds the upper and lower alarm limit. But, for adult and child, HR/PR, BP, RR, SpO₂, TEMP, EtCO₂/InspCO₂, TACHY, BRADY alarm generates after 5 seconds from the point when the limit is exceeded. For neonate, alarm will generate at the point when the limit is exceeded.

### **1** Press the Menu $\rightarrow$ Alarm keys.



The alarm setup menu will be displayed.

**2** Select the parameter to display the alarm setup menu.



The alarm setup menu for each parameter will be displayed.

<HR/PR Alarm Setup Menu>

Display	Description
0 50 100 150 200 250 300 ▲	<ul> <li>Displays the upper and lower alarm limit and the current measurement (▲).</li> <li>The limits can be adjusted by directly pressing the bar display or using the arrow keys.</li> </ul>
Lower Upper 40 ← 60 → 120	Displays lower limit←current value→upper limit.

Key	Item	Description	
ON OFF	Individual Alarm	Selecting ON will generate the alarm.	
		Selecting OFF will not generate the alarm.	
Lower Lower Alarm Limit		Sets the lower alarm limit. The lower limit will be turned OFF when a value	
		below the range is selected.	
Upper	Upper Alarm LimitSets the upper alarm limit.Upper Alarm LimitThe upper limit will be turned OFF when a v above the range is selected.		
Auto	Automatic Setup	Automatically sets the limits corresponding to the current value. If the limit is turned OFF, it will be remained OFF. The system alarm and parameter alarm will be in a ON condition.	

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu. (Monitor Setup).



For alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

### **3** Select ON/OFF and set upper and lower alarm limit for the parameter alarm.

Numeric Data Key	Item	Description
HR/PR 40-120 Recall	HR / PR / BPR	ON, OFF 20 to 300bpm*
STXXST Alarm OFF Recall	ST1 to ST12	ST All Alarm ON, OFF ST1 to ST12 ±2.0mV Individual Alarm ON, OFF
BP1 \$ 80-180 D 0FF-0FF Recall H 0FF-0FF	BP1	ON, OFF 0 to 300mmHg 0 to 40.0kPa
BP2 X 0FF 0FF Recall 11 0FF	BP2	ON, OFF 0 to 300mmHg 0 to 40.0kPa
SpO _{2 Sec} 90-0FF Recall	SpO ₂	ON, OFF 50 to 100%
RESP RR 5- 30 APNEA 15 Recall	RR	ON, OFF 5 to 150bpm (Adult) 2 to 150bpm (Child, Neonate)
	APNEA (Upper Limit)	ON, OFF 5 to 20 sec.
T 1 🕅 OFF Recall	TEMP1	ON, OFF 30 to 50°C 86 to 122°F
T2 X 0FF Recall	TEMP2	ON, OFF 30 to 50°C 86 to 122°F
NIBP S 80-180 D 0FF-0FF Recall	NIBP	ON, OFF 10 to 300mmHg* 1.5 to 40.0kPa
CO ₂ Et 30-45 Insp 3	EtCO ₂	ON, OFF 1 to 115mmHg* 0.1 to 15.0kPa 0.1 to 15.0%
Recall	InspCO ₂ (Upper Limit)	ON, OFF 1 to 24mmHg 0.1 to 3.0kPa 0.1 to 3.0%

▲ CAUTION	<ul> <li>NIBP Measurement Range: 10 to 280mmHg Alarm Range: 10 to 300mmHg</li> <li>CO₂ for MGU-722 (Microstream[®]CO₂ Unit)</li> </ul>
	<ul> <li>CO₂ for MGU-722 (Microstream[®] CO₂ Unit) Measurement Range: 0 to 99mmHg/0 to 13.3kPa Alarm Range: 1 to 115mmHg/0.1 to 15.0kPa</li> </ul>

### **4** Set ON/OFF and upper and lower limit for each BP alarm (BP1 to 5).



Pressing the BP Alarm key will display the alarm setup menu for BP1 to BP5.

Set ON/OFF and upper and lower limit (0 to 300mmHg / 0 to 40.0kPa) for each BP alarm by pressing the BP1 to BP5 keys.

### **5** Set ON/OFF and upper and lower limit for each temperature alarm (T1 to T3).



Pressing the Temp Alarm key will display the alarm setup menu for T1 to T3.

Set ON/OFF and upper and lower limit (30 to  $50^{\circ}C / 86$  to  $122^{\circ}F$ ) for each temperature alarm by pressing the T1 to T3 keys.

### **6** Set ON/OFF and upper and lower limit for blood temperature. (For HU-73)



Press the Tb key to set ON/OFF and upper and lower limit (30 to  $45^{\circ}$ C / 86 to  $113^{\circ}$ F) for the blood temperature.

# Arrhythmia Alarm Setup

The arrhythmia alarm can be turned ON or OFF, and arrhythmia detection level can be set.

**1** Press the Menu  $\rightarrow$  Alarm keys.



The alarm setup menu will be displayed. Press the Arrhy. key to display the arrhythmia alarm setup menu.

**2** Set the alarm condition for each arrhythmia. Pressing the Page Down key will display the next page.

Arrhythmia ala	rm 1/3	Page Do	wn Prev. Disp.
Asystole	5 ON	OFF	Ree
VF	ON	OFF	Recall
VT (HR > 120 bea	ats/min) ON	OFF	Recall
Slow VT		OFF	Recall
Tachy	ON	OFF	Recall
Brady	ON	OFF	Recall
			HR Alarm
Run (HR> 40 bea	3 ON	OFF	Recal
Run			
(HR > 40 loea Bigeminy	ats/min)		
- /			Recall
Trigeminy		OFF	Recall
Pause	3.0 ⁻ ON	OFF	Recall
Couplet	- ON	OFF	Recall
Frequent	10 ON	OFF	Recall
			HR Alarm
Arrhythmia ala	rm 3/3 Page	Цр	Prev. Disp.
HR Low Lin	nit for VT 120	140	

HR Alarm

HR Low Limit for RUN 40

The alarm setup menu for Asystole, VF, VT, Slow_VT, Tachy, and Brady will be displayed.

On the second page, the alarm setup menu for Run, Couplet, Bigeminy, Trigeminy, Pause, and Frequent will be displayed.

On the third page, the analysis condition (HR Low Limit) setup menu for VT and RUN will be displayed.

Page	Item
Page 1/3	ASYSTOLE, VF, VT, SLOW_VT, TACHY, BRADY
Page 2/3	RUN, BIGEMINY, TRIGEMINY, PAUSE, COUPLET, FREQUENT
Page 3/3	HR Low Limit for VT, HR Low Limit for RUN

### ●To Set ON/OFF of Arrhythmia Alarm

Selecting ON will generate the arrhythmia alarm. Selecting OFF will not generate the arrhythmia alarm.

Page	Item	Selection
	Asystole	ON (ON, OFF)
	VF	ON (ON, OFF)
Dogo 1/2	VT	ON (ON, OFF)
Page 1/3	Slow_VT	ON (ON, OFF)
	Tachy	ON, OFF
	Brady	ON, OFF
	Run	ON, OFF
	Couplet	ON, OFF
Dogo 2/2	Bigeminy	ON, OFF
Page 2/3	Trigeminy	ON, OFF
	Pause	ON, OFF
	Frequent	ON, OFF

<b>▲</b> CAUTION	<ul> <li>In case of the wired (DS-LANII) network, there are the following restrictions.</li> <li>Arrhythmia alarm of TACHY, BRADY, COUPLET, PAUSE, TRIGEMINY will not be transmitted.</li> <li>"SLOW_VT" will be transmitted as "VT".</li> <li>The setups for "HR Low Limit for VT" and "HR Low Limit for Run" cannot be performed on some central monitors and on some monitors with old software version.</li> </ul>		
	The "Arrhythmia alarm OFF" message will be displayed when the ASYSTOLE,		
NOTE	VF, VT, SLOW VT, and HR alarm is OFF.		



If ON is selected for "Asystole, VF, VT" on the hospital setup menu, Asystole, VF, VT, Slow_VT alarm can not be set to OFF.  $\rightarrow$  "8. System Configuration Hospital Setup Asystole, VF, VT"

### To Set the Arrhythmia Detection Level

- Select the level to detect each arrhythmia.
- 1 Pressing the detection level key (ex. 5 for Asystole) for each arrhythmia will display the window to adjust the detection level.



Set the detection level.

ltem	Range
Asystole	3 to 10 sec.
Run	2 to 8 beats
Pause	1.5 to 5 sec.
Frequent	1 to 50 beats/min.

#### **2** Close the window to adjust the detection level.

Close

Press the Close key.



Set the analysis condition to detect VT.

Arrhythmia alarm 3/3 Page Up	Pre	Υ. p.
HR Low Limit for VT 120 HR Low Limit for RUN 0	<b>1</b> 40	
	HR Alar	m
HR Low Limit for	VT 🗖	120

Select the HR low limit to detect VT from 120 or 140bpm. If HR is below the set value, it will be detected as Slow_VT.

### To Set the HR Low Limit for RUN

Set the analysis condition to detect RUN.

120

140

#### 1 Pressing the detection level key (ex. 40 for the following screen) will display the window to adjust the detection level.



### **2** Set the detection level.



Use the arrow keys to set the detection threshold. It can be set in the range from 0 to 100bpm.

### **3** Close the window to adjust the detection level.

Close

Press the Close key.



N The settings for the "HR Low Limit for VT" and "HR Low Limit for RUN" will be compared with the average HR of continuous VPC. Therefore, the displayed HR value at alarm generation may be lower than the settings if it is just after the VT detection, or if RUN with few continuous VPC is detected.

### Alarm Limit for TACHY, BRADY

The arrhythmia detection level for tachycardia (Tachy) and bradycardia (Brady) alarm links with the upper and lower alarm limit for HR / PR.

The tachycardia (Tachy) alarm generates when the measurement exceeds the HR / PR upper alarm limit. When the upper alarm limit is OFF, alarm will not generate.

The bradycardia (Brady) alarm generates when the measurement is below the HR / PR lower alarm limit. When the lower alarm limit is OFF, alarm will not generate.

# SpO₂ SEC Alarm Setup

# (Nellcor[®] Model; DS-7210)

When the  $SpO_2$  value is unstable around the lower alarm limit, the frequently generated alarm may be annoying. The SEC alarm function controls these frequent alarms.

This function generates the alarm only when the integral value (the accumulation of difference between the alarm limit and  $SpO_2$  value at every second) reaches the preprogrammed SEC alarm threshold value.



The integral value of the SEC alarm is calculated as follows.

On this graph, the SEC alarm threshold value is set as 100.

The  $SpO_2$  value begins to fall below the alarm limit at approximately 10 seconds. At the same time, the integral value begins to increase.

(Alarm limit) – (SpO₂ value) is accumulated each second.

At around 25 seconds, the integral value reaches 100 and the alarm is generated.

At approximately 36 seconds, the SpO₂ value returns to the level within the alarm limit, and at the same time, the integral value begins to decrease.  $\{(Alarm limit) - (SpO_2 value)\} \times 2$  is subtracted each second.

Also, there is a safety net when setting the SEC alarm function. This safety net is for the case when the  $SpO_2$  value frequently falls below the alarm limit but does not last long enough to reach the SEC alarm threshold.

If the SpO₂ value falls below the limit 3 times or more during the last 60 seconds, an alarm will be generated even if the SEC alarm threshold is not reached.

### **1** Press the Menu $\rightarrow$ Alarm keys.



The alarm setup menu will be displayed.

**2** Press the  $[SpO_2]$  key to display the alarm setup menu.



**3** Select the SEC alarm value according to the alarm frequency.

Indicator	2. If 10 / 25 / 50 / 100 is selected, a circular SEC alarm indicator will be displayed inside the parameter key. As the integral value increases, the indicator will begin to fill, and when it is completely filled, an alarm will be generated. If OFF is selected, this SEC alarm indicator will not be displayed.
<b>▲</b> CAUTION	<ul> <li>Whether to use the SEC alarm function and its threshold selection should be based on the patient's clinical indication portent and medical evaluation.</li> <li>If the SpO₂ alarm and SEC alarm setup is set to OFF, the SEC alarm integral value will be set to 0.</li> </ul>

# Alarm Suspend / Alarm Silence Time

The time for suspending the system alarm and suspending the alarm sound can be selected.

1	Press the M	$ enu  \rightarrow  A arm$	n → Se	tup keys.
	Alarm Setup	Prev. Disp.		
	Alarm Suspend Time	1min 3min 5min		
	Alarm Silence Time	1min 3min 5min		
	Alarm Limit Display	OFF	The al	arm setup menu will be displayed.
	Status Alarm Control	Linked to alarm silence time		
	Alarm Occurrence at NIBP Failure			
2	Select the tim	ne for "Alarm Su	spend Tin	ne".
	Alarm Suspend Ti	me 1min 3min	<b>5</b> min	Select the appropriate time for alarm suspend time.
3	Select the tim	ne for "Alarm Sil	ence Time	
	Alarm Silence Tim	e <mark>1</mark> min 3min	<b>5</b> min	Select the appropriate time for alarm silence time.

# **ON / OFF of Alarm Limit Display**

The alarm limit can be selected to display or not display on the home display.



NOTE	<ul> <li>The alarm limit for the parameter with the alarm turned OFF will not be displayed regardless of this setup.</li> <li>If the lower limit setting for SpO₂ alarm is 85% or below, alarm limit will be displayed even if OFF is selected for "Alarm Limit Display".</li> </ul>
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Alarm Silence Time for Equipment Status Alarm

The alarm silence time for the level 2 equipment status alarm ("Check electrodes", "Check SpO₂ sensor", etc.) can be set.

Reference

For details of level 2 equipment status alarm, refer to " –Display Configuration – Description of Alarm Message and Alarm Sound ●Equipment Status Alarm Message" in this chapter.

1	Press the M	lenu → Alarm	$\rightarrow$ Setup keys.
	Alarm Setup	Prev. Disp.	
	Alarm Suspend Time	1min 3min 5min	
	Alarm Silence Time	1min 3min 5min	
	Alarm Limit Display		
	Status Alarm Control	Linked to alarn silence time Linked to each new occurrence	The alarm setup menu will be displayed
	Alarm Occurrence at NIBP Failure		-

2 Set the "Status Alarm Control".

-
inked to each ew occurrence
e₩

Linked to alarm silence time will silence the alarm when the Alarm Silence key is pressed for fixed amount of time set on "Alarm Silence Time".

- If the alarm cause still remains at completion of silence time, the alarm sound will generate again.
- If the same alarm occurs during the alarm silence time, the alarm sound will not generate.
- If the new alarm occurs during the alarm silence time, the alarm sound for the new alarm will generate.

Linked to each new occurrence will silence the alarm when the Alarm Silence key is pressed until the situation changes.

- The alarm will be silenced as long as the alarm cause remains.
- If the alarm cause is resolved during the alarm silence time, the alarm silence will be cancelled.
- If the same alarm generates again during the alarm silence time, the alarm sound will generate.

# **ON/OFF of Alarm Occurrence at NIBP Failure**

The NIBP measurement failure can be notified by alarm.

Press the Me	enu → Alarm	$\rightarrow$ Setup keys.
Alarm Setup	Prev. Disp.	Alarm Setup menu will be displayed.
Alarm Suspend Time	1min 3min 5min	
Alarm Silence Time	1min 3min 5min	
Alarm Limit Display	OFF	
Status Alarm Control	Linked to alarn Linked to each silence time	
Alarm Occurrence at NIBP Failure		

2 Select ON or OFF for "Alarm Occurrence at NIBP Failure".

Alarm Occurrence at NIBP Failure

|--|

ON will display a "NIBP measurement failed." message (equipment status alarm, level 2) and generates an alarm sound when NIBP measurement fails.

This alarm can be cancelled by pressing the Alarm Silence key.

OFF will not generate alarm ("NIBP measurement failed." message and alarm sound) even if NIBP measurement fails.

Blank Page
# Chapter 5 Admit/Discharge of a Patient

This chapter describes the procedure to admit or discharge a patient to the monitor.

Admit / Discharge of a Patient5-2
Admitting a Patient Name, Sex, and Age5-3
Patient Name5-3
Patient ID5-4
Patient Classification5-5
Patient Sex5-6
Pacemaker Use5-6
Pacemaker Use Selection
Impedance Respiration Measurement5-7
Patient Age5-7
Acquiring Patient Information from the Patient Data
Server (TCON system)
Discharging a Patient Erasing Name, Data, etc 5-9
Discharging Procedure5-9
Discharge Procedure by User Key5-10
Monitoring Mode Selection
Alarm / Display Mode5-11
Mode Selection
Suspend Monitoring
Suspend / Resume Monitoring5-12
To Suspend Monitoring

# Admit / Discharge of a Patient

This menu allows admitting, discharging, suspend monitoring of a patient, and selection of the display configuration mode and alarm mode according to the monitoring purpose.

<b>A</b> CAUTION	If you start monitoring a new patient without performing a discharge procedure for the previous patient, new data will be added to the previous data which will result in inaccuracy.
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### **1** Press the Menu $\rightarrow$ Admit / Discharge keys.



The admit / discharge menu will be displayed.

There are two ways to enter the patient information.

- 1. Manually using the alphanumeric keypad displayed on the screen.
- 2. Automatically acquiring patient information from the patient data server using the patient ID via TCON communication with the central monitor.

# **Admitting a Patient**

# Name, Sex, and Age

This menu allows entering patient's name, ID, age, and selection of the patient classification (adult, child, neonate) and pacemaker use (used, not used), which affects the monitoring accuracy.



Whether to display or not display the following keys on the admit menu can be selected.

Impedance Mode key
Filter Mode key (Monitor / ESIS / DIAG.)
Bed ID key

For procedures, refer to "8. System Configuration Hospital Setup"

# **Patient Name**

For patient's name, up to 16 characters can be used.



**1** Press the Name key.



Enter the name using the alphanumeric keypad.

**2** The entered patient's name will be displayed on the home display.



# **Patient ID**

Up to 20 characters of alphabets, numbers, or symbols can be used for the patient ID. Enter the ID according to the monitoring purpose. The entered ID will be printed on the recording paper.



#### **1** Press the ID key.



Enter the ID using the alphanumeric keypad. 20 digits can be input, but only 10 digits can be transmitted through the DS-LAN II network. On the hospital setup under the preset menu, set which 10 digits to send to the central monitor.

If DS-LANIII network is used, all 20 digits can be transmitted to the central monitor.



Refer to "8. System Configuration Patient ID Transmission Starting Digit for DS-LANII" for procedure to set the 10-digit patient ID.

# **Patient Classification**

The selection of the patient classification affects the accuracy of NIBP, HR, and RR measurement. Also the delay time to generate the measurement data alarm will change according to the patient classification.

		Adult	Adult Child	
	SYS	60 to 250mmHg	30 to 180mmHg	40 to 120mmHg
NIBP	MEAN	45 to 235mmHg	15 to 160mmHg	30 to 100mmHg
	DIA	40 to 200mmHg	10 to 150mmHg	20 to 90mmHg
HR		0bpm, 12 to 300bpm		0bpm, 30 to 300bpm
	Monitor	0.5 to 40Hz		1.6 to 40Hz
Filter Mode	ESIS	1.6 to 15Hz		1.6 to 15Hz
	Diagnosis	3-lead: 0.05 to 100Hz		
Diagnosis		4, 5, 10-lead: 0.05 to 150Hz		
Impedance Respiration		1.5Hz		2.5Hz
Alarm Delay Time		5 sec.		0 sec.

The alarm delay time is a function to prevent frequent generation of the measurement data alarm by holding the alarm generation for the corresponding alarm delay duration.

The alarm delay time applies to the measurement data alarm for the following parameters; HR / PR, BP, RR, SpO₂, TEMP, EtCO₂ / InspCO₂, TACHY, BRADY.

**MARNING** The monitor determines the detection algorithms for QRS and NIBP according to the selected patient classification. Make sure the proper selection is made.

#### **1** Select the patient classification.



#### **2** The selected patient classification will be displayed on the home display.

	Patient Classification
FUKUDA DENSHI	oult 02/01 13:42 M
	HR Av.

# **Patient Sex**

Select the patient's sex from male or female. The default is set as undetermined. The selected sex will be printed on the recording paper.

#### **1** Select Male or Female.



This selection will not affect the measurement accuracy of the monitoring.

# Pacemaker Use

#### Pacemaker Use Selection

If the patient is wearing a pacemaker, the monitor will identify the pacemaker pulse and insert an artificial pulse onto the ECG waveform for easy identification. By detecting the pacemaker pulse, it prevents to erroneously detect QRS as pacemaker pulse when pacing waveform does not appear (pacing failure). The arrhythmia analysis analyzes pacing beat as P (pacemaker beat) or F (fusion beat) to prevent erroneous judgment of VPC.

MARNINGThe pacemaker use selection influences the precision of the QRS detection<br/>and arrhythmia analysis. Make sure the correct selection is made.

**1** Select Used or Not used for pacemaker use.



**2** The pacemaker use will be displayed on the home display.



#### Impedance Respiration Measurement

The respiration measurement using the impedance method conducts high-frequency and weak current between the ECG electrodes attached to the patient, and measures the potential difference between the electrodes caused by thoracic movement using the synchronous rectification system. For the patient using the adaptive (minute ventilation) pacemaker, the pacemaker measurement signal and the high-frequency current of this equipment interferes with each other which causes incorrect respiration measurement.

If the patient is using an adaptive (minute ventilation) pacemaker, set the impedance respiration measurement OFF.

#### **1** Press the Impedance Mode key.



# **Patient Age**

There are two ways to set patient's age. One is to enter the birth date which will automatically calculate the age, and the other is to directly enter the age using the numeric keypad.

If <u>Neonate</u> is selected as patient classification, the age in days will be displayed.



**1** Press the Age key, and enter patient's age.



To directly enter the age, use the numeric keypad to enter the age and press the Yrs key. The entered age will be displayed inside the key.

**2** Enter the patient's birth date using the numeric keypad. The age will be automatically calculated.



Enter the year, month, day using the numeric keypad, and press the Yr, Mo, Dy keys respectively. The entered year, month, day will be displayed inside the Yr, Mo, Dy keys respectively.

# Acquiring Patient Information from the Patient Data Server (TCON system)

If the TCON system is used and the central monitor is connected to the Patient Data Server, the patient information can be automatically acquired from the Patient Data Server via the central monitor.

- **1** Press the Menu  $\rightarrow$  Admit / Discharge  $\rightarrow$  ID keys.
- **2** Enter the patient ID.
- **3** Press the Search Patient key to search patient information on the patient data server.

ID	Prev. Disp.		
<u>111111112</u>	Erase		
123456 QWERTYUIOP			
	<b>×</b> ← →		
ABC<-> QVERTY			
	Search Patient		
	3		
Based on the entered pa server via TCON comm The acquired patient inf	unication with the	e central monitor.	
Admit Patient Information	]		
Current Information ID 11111111111			
Name FUKUDA			
New Information			
1111111111 Sex : No Input			
Birth Date: Yr Name	Mo Dy( OYrs)		
Change only patient info, (Keep current meas. data/settings.)	rching		
Admit as new patient. (Initialize curr. data/settings.)	Cancel		

Change only patient info. : Only the patient information will be changed to the new information. This function is mainly used to correct the patient information.

Admit as new patient. : Initializes the current patient data/monitoring condition and performs the admit process with the newly acquired information.

Cancel : Cancels the acquired data.

The item not acquired from the patient data server will be left blank. For the blank item, manually input the information.

Reference	For procedure on how to manually enter data using the alphanumeric keypad, refer to the above "Patient Name", "Patient Classification", "Patient Sex", "Pacemaker Use", and "Patient Age".

After the information for a new patient is acquired by searching the patient data server, make sure to perform the admit process by pressing the Admit as new patient. key.

# **Discharging a Patient**

# Erasing Name, Data, etc.

This menu allows to clear the patient name, ID, age, and past measurement data such as tabular trend, graphic trend, and recall data.

By pressing the Rapid Discharge key preprogrammed as user key, a quick discharge process can be performed.

# **Discharging Procedure**



#### **1** Press the Discharge key.

The confirmation display will appear. To cancel the discharge process, press the **Cancel** key to return to the previous display.

	•		
Discharge			
All	Data for this pa	itient will be erased	
	ок	Cancel	

#### **2** Erase the patient data.

To initialize the patient data and patient information, press the OK key. The data will be initialized and the screen will return to the home display selected for the display mode. The alarm setup will be initialized to the selected alarm mode.

Data	Description
Patient Data	Erases the data of graphic trend, tabular trend, NIBP list, recall, ST Display, OCRG, cardiac output, hemodynamic, P-V / F-V control data, vigilance list. The setup condition of recall setup, tabular trend setup, graphic trend setup, vigilance list will remain.
Patient Erases the data of patient name, ID, sex, age.	
Information The patient classification will not be initialized.	
Measuremen Condition	Pacemaker use will be set to unused, and impedance respiration measurement will be set to ON. The BP zero-balance condition will be cleared.
NOTE	<ul> <li>Depending on the setup of "Backup at Discharge" on the monitor setup menu, some data may not be initialized.</li> <li>If discharge procedure is performed during stopwatch operation, the counting will stop and the stopwatch time will be reset to "00:00:00".</li> </ul>

Monitor Suspend Patient has been discharged. Admit	When Suspend is selected for "Discharge Mode" (Monitor Setup), the window as shown in the left will be displayed. Pressing the Admit key will start monitoring.
Reference Refer to "To Resume Mo	nitoring" (p.5-13).

# **Discharge Procedure by User Key**

The patient's monitoring data (trend data), patient information, monitoring condition can be initialized by pressing the Rapid Discharge key preprogrammed as user key. This key will function the same as the discharge procedure.

**1** Press the Rapid Discharge key preprogrammed as user key.



The confirmation message for erasing the data will be displayed.

**2** Pressing the OK key will initialize the patient data.

# **Monitoring Mode Selection**

# Alarm / Display Mode

This menu allows to select the alarm mode and display mode.

The alarm setting and display configuration can be each selected from 5 modes depending on the monitoring purpose. Select the appropriate mode when admitting a patient.

<b>A</b> warning	The setup for the alarm mode and display mode remains stored even when the power is turned off or when discharging procedure is performed. Before monitoring, make sure the current monitoring mode is suitable for the patient's condition.		
<b>CAUTION</b> If the built-in backup battery is depleted when the power is turned ON, the alarm mode setting will be initialized to default setting.			
Refer to "11. Technical Information Setup Item Default and Backup" for default settings of the display modes and alarm modes.			

# **Mode Selection**

Admit	Prev. Disp.
Name	FUKUDA DENSHI
Sex	Male Female Class Adult Child Neonate
Height Weight	CM Blood type Rh kg ( BSA m ² )
Age ID	OYrs yr Mo Day
Pacemaker	Used Not used Impedance Mode
Filter Mode	Monitor ESIS DIAG. Admit date
Discharg	Monitor Bed ID Bed ID



Mode select	Prev. Disp.
Alarm Mode	Display Mode
ALM MODE1	DISP MODE1
ALM MODE2	DISP MODE2
ALM MODES	DISP MODE3
ALM MODE4	DISP MODE4
ALM MODE5	DISP MODE5

The mode selection menu for alarm mode and display mode will be displayed.

#### **2** Select an alarm mode from the selection.

Alarm Mode	Display Mode
ALM MODE1	DISP MODE1
ALM MODE2	DISP MODE2
ALM MODE3	DISP MODE3
ALM MODE4	DISP MODE4
ALM MODE5	DISP MODE5

Select from the	ALM MODE1	to	ALM MODE5	keys or	
DISP MODE1	to DISP MC	DE5	keys to set th	e mode which	
meets the monitoring purpose.					

The alarm mode and display mode setting remains stored even when the power is turned off or when discharging procedure is performed. The previously selected mode will be effective if the selection is not made.

# **Suspend Monitoring**

# **Suspend / Resume Monitoring**

This menu allows to suspend and resume monitoring when a patient temporarily leaves the bed. Turning the power OFF will erase the recall and ST measurement data, but with this suspend monitoring function, data measurement, alarm generation, automatic measurement, and automatic recording can be suspended without erasing any data and setup condition.

# **To Suspend Monitoring**

Admit	Prev. Disp.
Name	FUKUDA DENSHI ( )
Sex	Male Fenale Class Adult Child Neonate
Height Weight	cm Blood type Rh kg ( BSA m²)
Age ID	OYrs yr Mo Day
Pacemaker	Used Not used Impedance Mode
Filter Mode	Monitor ESIS DIAG. Admit date
Discharge	Mode select Monitor Bed ID

#### **1** Press the Monitor Suspend key.

The confirmation menu will be displayed. If the Monitor Suspend key was pressed by mistake, press the Cancel key to return to the previous display.

Monitor Susp	end		
	Monitorina wi	l be suspended.	
	normorning wi	i be suspended.	
	ок	Cancel	
		Garicer	

#### **2** Suspend monitoring.

Pressing the OK key will return to the home display with the Resume key displayed. The numeric data display and waveform display on the home display will be suspended and only the Resume key will be effective.

Monitor S	Suspen	al							
Μ	0	n	i	t	0 S	r	i	ng	
s	u	s	р	ė	-	d	е	d.	
			[	Res	ume	]			



**CAUTION** Resuming monitoring will also resume the suspended alarm.

#### **1** Press the Resume key.

The monitor suspend display will be cleared and monitoring will resume.



2 When the discharge procedure is performed, and <u>Suspend</u> is selected for "Discharge Mode" (Monitor Setup), the following window will be displayed. Pressing the Admit key will clear the window and start monitoring.



Blank Page

# **Chapter 6**

# **Parameter Setup**

#### **Parameter Setup**

Setting the Monitoring Condition	6-2
To Display the Parameter Setup Menu	6-2
Zero Balance of All Pressure Lines	6-3
Zero Balance of All Pressure Lines (User Key) -	6-3
ECG	
ECG Waveform Size and Lead	6-5
HR/PR Alarm	6-7
Arrhythmia Alarm	6-7
Filter Mode Selection	
HR Average Selection	6-9
HR Synchronization Mark and Pulse Tone	6-9
HR/PR Alarm Source	6-10
Automatic Lead Switching	
Pacemaker Pulse	6-12
QRS Pace Pulse Mask	6-13
AC Filter ·····	
ECG Drift Filter	
3-lead Override	6-15
Pace Pulse Detection	6-15
ON/OFF of Parameter Display	6-16
Respiration	6-17
Respiration Waveform Size	
RR Alarm	6-18
APNEA Alarm	6-18
CVA Detection	6-19
RR/APNEA Alarm Source	6-20
Impedance Respiration Measurement	6-21
RR Synchronized Indicator	6-21
ON/OFF of Parameter Display	6-22
Invasive Blood Pressure (BP1 to BP5) BP Label	••••• 6-23
Invasive Blood Pressure (BP1 to BP5)	••••• 6-23
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale	••••• 6-23 •••• 6-23 •••• 6-25 •••• 6-26
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm	••••• 6-23 ••••• 6-23 ••••• 6-25 ••••• 6-26 ••••• 6-26
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines	••••• 6-23 •••• 6-25 •••• 6-26 •••• 6-26 •••• 6-26 •••• 6-27
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection	<b>6-23</b> 6-23 6-25 6-26 6-26 6-27 6-28
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART)	<b>6-23</b> <b>6-23</b> <b>6-25</b> <b>6-26</b> <b>6-26</b> <b>6-26</b> <b>6-27</b> <b>6-28</b> <b>6-28</b>
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data	<b> 6-23</b> 6-25 6-26 6-26 6-26 6-27 6-28 6-28 6-29
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display	6-23           6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter	6-23           6-25           6-26           6-27           6-28           6-28           6-29           6-30
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-30           6-31
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210)	6-23           6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Waveform Size	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-32           6-33
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm	6-23           6-23           6-25           6-26           6-26           6-27           6-28           6-29           6-30           6-31           6-33           6-33
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-34
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm Source SpO ₂ Alarm during NIBP Measurement	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-34
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm Muring NIBP Measurement ON/OFF of Parameter Display	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-34           6-35           6-36
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M)	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-32           6-33           6-33           6-34           6-35           6-36
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-32           6-33           6-33           6-34           6-35           6-36           6-37           6-38
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-38
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-38           6-39
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm	6-23           6-25           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-38           6-39           6-38           6-39           6-39
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Qaveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm Non-Invasive Blood Pressure	6-23           6-25           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-39           6-38           6-39           6-43
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Qaveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm Non-Invasive Blood Pressure NIBP Automatic Measurement	6-23           6-25           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-33           6-34           6-35           6-36           6-37           6-38           6-38           6-39           6-34           6-35           6-36           6-37           6-38           6-39           6-44           6-45
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Qaveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm NO/OFF of Parameter Display Non-Invasive Blood Pressure NIBP Automatic Measurement NIBP 1-Minute Interval Measurement	6-23           6-25           6-26           6-26           6-27           6-28           6-28           6-29           6-30           6-31           6-32           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-38           6-39           6-43           6-44           6-45
Invasive Blood Pressure (BP1 to BP5) BP Label PCWP Measurement BP Scale BP Alarm Zero Balance of Pressure Lines Filter Selection HR/PR Alarm Source (BP1 or ART) Display Selection of Numeric Data Mean BP Waveform Display Respiration Rejection Filter ON/OFF of Parameter Display SpO ₂ (Nellcor [®] Model: DS-7210) SpO ₂ Qaveform Size SpO ₂ Alarm HR/PR Alarm Source SpO ₂ Alarm during NIBP Measurement ON/OFF of Parameter Display SpO ₂ (Masimo [®] Model: DS-7210M) Pulse Wave Size SpO ₂ Alarm SpO ₂ Alarm Non-Invasive Blood Pressure NIBP Automatic Measurement	6-23           6-25           6-26           6-26           6-27           6-28           6-29           6-30           6-31           6-32           6-33           6-33           6-33           6-34           6-35           6-36           6-37           6-38           6-38           6-39           6-43           6-44           6-45           6-46

End of Measurement Tone Mean BP Display Dyna Alert Function Pump Setup Sight Inflation Pulse Rate Display Oscillograph Display Cancelling the NIBP System Error Message <b>Temperature (T1 to T3)</b> Temperature Label Temperature Alarm	6-48 6-50 6-50 6-51 6-51 6-52 6-53 6-53
∆T Display	6-54
ON/OFF of Parameter Display CO ₂ Concentration (Option Unit: MGU-721)	
CO ₂ Scale	6-56
EtCO ₂ (End-Tidal CO ₂ ) Alarm	
InspCO ₂ (Inspiratory CO ₂ ) Alarm	
EtCO ₂ Peak Picking Duration	6-58
Measurement Unit	6-59
O ₂ Compensation	
N ₂ O Compensation	6-60
Anesthetic Gas Compensation	6-60
Atmospheric Pressure Compensation	6-61
ON/OFF of Parameter Display	
CO ₂ Concentration (Option Unit: MGU-722)	6-63
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale	<b>6-63</b>
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm	•••••6-63 •••••6-63 •••••6-64
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm	•••••6-63 ••••6-63 ••••6-64 ••••6-65
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration	<b>6-63</b> 6-63 6-64 6-65 6-65
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit	<b>6-63</b> 6-64 6-65 6-65 6-66
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration	<b>6-63</b> 6-63 6-64 6-65 6-65 6-66
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale Tb (Blood Temperature) Tb Alarm	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale Tb (Blood Temperature) Tb Alarm	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale Tb (Blood Temperature) Tb Alarm Vigilance Data STAT Mode / Index Display	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale Tb (Blood Temperature) Tb Alarm Vigilance Data STAT Mode / Index Display Stopwatch	<b>6-63</b> 
CO ₂ Concentration (Option Unit: MGU-722) CO ₂ Scale EtCO ₂ (End-Tidal CO ₂ ) Alarm InspCO ₂ (Inspiratory CO ₂ ) Alarm EtCO ₂ Peak Picking Duration Measurement Unit CO ₂ Calibration Restarting the CO ₂ Unit ON/OFF of Parameter Display Suspending CO ₂ Measurement Ventilator AWP Scale AWF Scale Tb (Blood Temperature) Tb Alarm Vigilance Data STAT Mode / Index Display	<b>6-63</b> 

# 6 Parameter Setup

# **Parameter Setup**

# **Setting the Monitoring Condition**

This menu allows setting the measurement condition, waveform size, scale, etc. of ECG, BP, NIBP, SpO₂, RESP, TEMP, CO₂, and Vigilance data.

# To Display the Parameter Setup Menu

Press the Menu  $\rightarrow$  Parameter keys to display the parameter setup menu, and select the parameter. On the parameter setup menu, BP zero balance can be performed.

The parameter setup menu for each parameter can be also accessed by pressing the parameter key where the numeric data is displayed.



By selecting Store for "Parameter Key Operation" in the monitor setup menu, the previously displayed screen can be directly accessed when pressing the parameter key.

Reference

For setup procedure, refer to "8. System Configuration Monitor Setup".

# Zero Balance of All Pressure Lines

If all the displayed BP's are opened to air, the zero balance procedure for all BP can be performed. If any of the BP is in progress of measurement, perform the zero balance on each BP parameter setup menu.

- **1** Open the three-way valve of all the pressure transducers to air.
- **2** Press the Zero All BP key when "Perform zero balance" message is displayed.



Verify the BP waveform is positioned at zero, and "0" is displayed for the BP value. A message, "BP zero complete" will be displayed when the procedure is complete. A message, "BP zero failed" will be displayed when the process fails. The three-way valve may not be opened to air, artifact may be present, or the transducer may be defective.

Check the cause and try the zero balance procedure again. A message, "BP zero drift" will be displayed when the interface cable is not connected. Check if the cable is connected correctly.

**3** Close the three-way valve when the zero balance is complete.

▲ CAUTION	<ul> <li>Each time the blood pressure transducer or tubing is replaced, the zero balance procedure is required to ensure accurate measurements.</li> <li>"Perform zero balance" message will not be displayed unless the three-way valve of all pressure transducers are opened to air. If the status is not</li> </ul>
	displayed, or if "Open stop cock to air" message is displayed, check if the three-way valve of pressure transducers are opened to air.

# Zero Balance of All Pressure Lines (User Key)

The zero balance procedure of all pressure lines can be performed using the user key. If any of the BP is in progress of measurement, perform the zero balance on each BP parameter setup menu.



**1** Open the three-way valve of all pressure transducers to air.

#### **2** Press the **BP Zero** key when "READY" message is displayed inside the user key.

Verify the BP waveform is positioned at zero, and "0" is displayed for the BP value. A message, "COMPLETE" will be displayed when the procedure is complete. A message, "FAILED" will be displayed when the process fails. The three-way valve may not be opened to air, artifact may be present, or the transducer may be defective.

Check the cause and try the zero balance procedure again.

A message, "DRIFT" will be displayed when the interface cable is not connected. Check if the cable is firmly connected.

# **3** Close the three-way valve when the zero balance is complete.

# [BP zero status displayed inside the user key]

BP Zero DRIFT BP zero status BP zero status BP zero status BP zero status	<ul> <li>Open transducer to air</li> <li>Open transducer to air</li> <li>Ready to perform zero balance.</li> <li>BP zero in progress</li> <li>BP zero failed</li> <li>BP zero complete</li> <li>BP zero drift</li> </ul>
------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

▲CAUTION
----------

# **ECG**

EC6 Lead·Size	ECG1:	×1,		Prev. Disp.
HR/PR Alarm	ON	40-1	20	Graphic Trend
Arrhy.Alarm	Arrhy, Learn			
Config.		Pace Pulse	0FF	
Filter	Monitor	Pace Pulse Mask Time	Auto	
HR Average Pulse Tone Alarm Source	Average ON Auto	ECG Drift Filter AC Filter	off On	Display ON/OFF ON
3lead Override Auto Lead	off off	Pace Pulse Detection	Med.	

The ECG measurement condition can be set on this menu.

Lead, Size : Sets the waveform size and lead for ECG display and recording. HR/PR Alarm : Sets ON/OFF of the HR/PR alarm, and sets the upper and lower alarm limit. Arrhy. Alarm : Sets ON/OFF and the detection threshold for each arrhythmia alarm. Configuration : Sets the condition for measuring ECG and HR. Arrhy. Learn : The monitor learns the normal QRS at ECG electrode replacement or at misdetection of the arrhythmia analysis.

# ECG Waveform Size and Lead

#### **1** Press the Lead, Size key to display the size / lead setup menu.

Select the lead and size for the waveform to display on the home display.



<6 Waveforms Display>



#### **2** Select the ECG channel to set the lead, size, and the baseline position.

ECG1 ×1 ,	ÉCG2	×1 , aVR
-----------	------	----------

Select the ECG channel by pressing the  $\boxed{ECG*}$  key. The selected channel key LED will light.

#### **3** Select the waveform size.

Auto Size ×Ķ ×2 ×4

Select the waveform size for displaying and recording. Pressing the Auto key will automatically adjust the ECG amplitude to 10mm. The automatic adjustment is effective only when the key is pressed.

Size	×1/4	×1/2	×1	×2	×4
Voltage (10mm)	4mV	2mV	1mV	500μV	250μV

▲ CAUTION	<ul> <li>The threshold level for arrhythmia detection and QRS detection changes with the ECG waveform size. Set a proper waveform size for monitoring. When the ECG waveform size is ×1/4, ×1/2, or ×1, the detection threshold is 250µV.</li> <li>When the ECG waveform size is ×2 or ×4, the detection threshold is 150µV.</li> <li>Automatic size/position of the ECG is effective only at the time the Auto key is pressed. This does not continuously adjust the size and position.</li> </ul>
-----------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### **4** Select the lead for ECG1, ECG2.

The leads can be selected from 3 leads, 6 leads, 7 leads, or 12 leads depending on the connected ECG relay cable.

ECG Relay Cable	Lead
3-electrode	
4-electrode	
5-electrode	
10-electrode	$\begin{bmatrix} I \\ I \\ V_1 \end{bmatrix} \begin{bmatrix} I \\ V_3 \end{bmatrix} \begin{bmatrix} I \\ V_4 \end{bmatrix} \begin{bmatrix} I \\ V_5 \end{bmatrix} \begin{bmatrix} I \\ V_6 \end{bmatrix}$

When the 12-lead display configuration is used, press the ECG1, 2 key to display the lead selection menu.



#### **5** Set the baseline position.



If the waveform is difficult to see due to ECG amplitude, set the 0 mV baseline position.

The baseline position for the waveform display and recording will be adjusted.

# **HR/PR Alarm**

#### **1** Press the HR/PR Alarm key to display the alarm setup menu.



Select ON/OFF for HR/PR alarm, and set the upper and lower alarm limit.

The common alarm value for HR measured from ECG, PR measured from  $SpO_2$ , PR measured from BP can be set. The upper and lower limit can be set in 5 bpm increments.

Key	ltem	Description	
ON OFF	Individual Alarm	SelectingONwill generate the HR/PR alarm.SelectingOFFwill not generate the HR/PR alarm.	
← Lower →	Lower Alarm Limit	Sets the lower alarm limit (20 to 295bpm). Setting a value 20bpm or below will turn OFF the alarm.	
🗲 Upper 🔿	Upper Alarm Limit	Sets the upper alarm limit (25 to 300bpm). Setting a value 300bpm or above will turn OFF the alarm.	
Auto	Automatic Setup	Automatically sets the upper limit to +40bpm, and the lower limit to -40bpm from the current value. The lower limit will be restricted by "HR/PR Low Limit during Alarm Auto Setting (OFF/30bpm/40bpm)" under the Hospital Setup.	

To maintain the alarm setting even after the power is turned OFF or after a discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# Arrhythmia Alarm

**1** Press the Arrhy. key to display the arrhythmia alarm setup menu.

 $\ensuremath{\mathsf{ON/OFF}}$  for each arrhythmia alarm and analysis threshold level can be set.

Arrhythmia alarm 1/3	Page Down Prev. Disp.			
Asystole 5	ON OFF Recall			
VF [		Reference	Refer to "4. Monitoring Setup	Alarm Setup" for details.
VT (HR > 120 beats/min)		0		
Slow VT				
Tachy				
Brady				
	HR Alarm			

# **Filter Mode Selection**

The waveform frequency characteristic can be selected from Monitor mode, ESIS mode, or Diagnosis mode according to the monitoring purpose. Each mode has different frequency characteristic. The selected filter mode will be displayed on the home display. (Monitor: M, ESIS: E, Diagnosis: D)



The selected filter mode will be printed when recording.

#### **1** Press the Config. key.

ECG Configuration 1/3	Page Down Prev. Disp.
Filter	Monitor ESIS
HR Average	Instant Average
Pulse Tone	
HR/PR Alarm Source	
	Auto

The configuration menu for selecting a filter will be displayed.

2 Select the filter mode from Monitor, ESIS, or Diag.

#### **Monitor Mode**

#### (Frequency Characteristic: Adult/Child 0.5 to 40Hz, Neonate 1.6 to 40Hz)

This is the standard mode for ECG monitoring. The upper frequency is set to 40Hz to reduce artifact caused by EMG, etc.

# ESIS Mode (When a defibrillation and electrosurgery-proof ECG relay cable is used) (Frequency Characteristic: Adult/Child/Neonate 1.6 to 15Hz)

The upper frequency is set to 15Hz, so that it can largely reduce the high-frequency artifact.

#### **Diagnosis Mode**

#### (Frequency Characteristic for 3-electrode Lead: Adult/Child/ Neonate 0.05 to 100Hz Frequency Characteristic for 4, 5, 10-electrode Lead: Adult/Child/ Neonate 0.05 to 150Hz)

Select this mode if ST measurement or high frequency ECG monitoring is performed. As the lower frequency is set to 0.05Hz, ST level can be accurately measured.

<b>▲</b> CAUTION	The ESIS mode can largely reduce the artifact such as electrosurgery noise and EMG, but it may also reduce the QRS amplitude. The ESIS mode should be selected only during electrosurgery.
------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



# **HR Average Selection**

The averaging method of HR measured from ECG can be selected.

1 Press the Config. key.



The configuration menu for selecting HR average method will be displayed.

2 Select Instant or Average.

Selecting Instant will display the HR measured from RR interval of each heartbeat. Selecting Average will display the HR measured from 6 seconds of heartbeat for adult and child, and 3 seconds of heartbeat for neonate.

# HR Synchronization Mark and Pulse Tone

The HR mark synchronized to ECG or PR can be displayed inside the parameter key. ON/OFF of HR pulse tone can be also set.



**1** Press the Config. key.



The configuration menu for setting the pulse tone will be displayed.

# 2 Select from ON / OFF / SpO₂.

OFF will not display the synchronized mark. The pulse tone will not be generated.
 ON will display the synchronized mark. The pulse tone will be generated.
 SpO₂ will always synchronize to SpO₂ for synchronized mark/tone regardless of the "HR/PR Alarm Source" setup.

# **HR/PR Alarm Source**

The parameter to display the HR synchronized mark and to generate the HR/PR alarm can be selected from ECG, SpO₂, and BP (BP1 or ART).

# **1** Press the Config. key.

ECG Configuration 1/3	Page Down Disp.
Filter	Monitor ESIS
	DIAG.
HR Average	Instant Average
Pulse Tone	
HR/PR Alarm Source	
	Auto
-	

The configuration menu for selecting the HR/PR alarm source will be displayed.

#### 2 Select a parameter.

ECG	SpO2	BP
Auto		

Selecting ECG will generate the alarm based on HR measured				
from ECG. HR synchronized mark will be displayed.				
Selecting SpO ₂ will generate the alarm based on PR measured				
from $SpO_2$ . $SpO_2$ synchronized mark will be displayed.				
Selecting BP will generate the alarm based on PR measured from				
BP (BP1 or ART). BP synchronized mark will be displayed.				
However BP can be selected only when ECG/SpO ₂ /BP is				
selected for "HR/PR Source" of the monitor setup menu.				
Selecting Auto will automatically select the measurable HR				
source in the priority of ECG>SpO ₂ >BP. If the corresponding				
parameter key is displayed, alarm generation will be also effective.				

▲WARNING	<ul> <li>The HR/PR alarm will not be generated unless the parameter key corresponding to the selected HR/PR source is displayed. Be sure to display the parameter key for the HR/PR source.</li> <li>The alarm for the parameter not selected for the "HR/PR Alarm Source" (ECG/SpO₂/BP) will be set to OFF on the DS-7600 Central Monitor.</li> <li>The "HR/PR Alarm Source" setting will synchronize between the bedside monitor and the central monitor.</li> <li>For example, if PR is set as the HR/PR alarm source on the DS-7200, HR alarm will be set to OFF on the central monitor.</li> </ul>
<b>▲</b> CAUTION	In case of DS-LANII network, if the HR/PR source is BP (Or, if Auto selects BP for HR/PR source), the ECG waveform will not be transmitted on the network. On the central monitor, PR_IBP value will be displayed instead of HR. However, on some central monitor depending on the model type, the HR value from ECG will be displayed on the NIBP list and ST measurement list. Refer to the operation manual for the respective central monitor. In case of DS-LANIII network, refer to the operation manual for the central monitor.
NOTE	If $[SpO_2]$ is selected for "Pulse Tone" (ECG config.), the synchronized mark will be displayed inside the PR_SpO ₂ numeric data box regardless of the HR/PR alarm source setup.

# **Automatic Lead Switching**

When the lead off condition occurs, the "Check Electrodes" message will be displayed and a new ECG lead will be automatically set if the Automatic Lead Switching is set to ON. The automatic lead switching will be performed for ECG 1 and ECG 2.

#### Lead Switching

Lead Type	Electrode Off	Auto Lead Selected	
Leau Type	Electrode On	ECG1	ECG2
4-electrode	RA	III	III
Lead	LA	II	II
5-electrode Lead	RA / RA+V	III	III
	LA / LA+V	II	II
	V	II	aVR
10-electrode	RA / RA+V	III	III
	LA / LA+V	II	II
Loud	V, V2 to V6	II	aVR

#### **1** Press the Config. $\rightarrow$ Page Down keys.

ECG Configuration 2/3	Page Up Page Down Prev. Disp.
Auto Lead Switch	OFF OFF
Pacemaker Pulse	
Pace Pulse Mask Time	Auto 10ms 20ms 40ms 0FF
AC Filter	
ECG Drift Filter	
3lead Override	OFF OFF

The configuration menu for selecting the automatic lead switching will be displayed.

# 2 Select ON or OFF.

OFF will not switch the lead when an electrode comes off. ON will automatically switch to another lead when an electrode comes off.

# **Pacemaker Pulse**

The pacemaker pulse can be displayed in yellow superimposing on the ECG waveform.

#### Pacemaker Pulse Detection Algorithm



#### **1** Press the Config. $\rightarrow$ Page Down key.

ECG Configuration 2/3	Page Up	Page Dow	n Disp.
Auto Lead Switch	ON	OFF	
Pacemaker Pulse	ON	OFF	
Pace Pulse Mask Time	Auto 40ms	10ms OFF	20ms
AC Filter	ON	OFF	
ECG Drift Filter	ON	OFF	
3lead Override	ON	OFF	

The configuration menu to select pacemaker pulse will be displayed.

# 2 Select ON or OFF.

OFF will not display the pacemaker pulse.

ON will display the pacemaker pulse in a different color from the ECG waveform. This will automatically set to ON when "Used" is selected for pacemaker use on the patient admit/discharge menu.

# **QRS Pace Pulse Mask**

For patients using pacemakers, there are cases when the paced QRS may not occur in spite of the pacing stimulus. This condition is called "pacing failure", or "failure to capture". To avoid detecting pacemaker pulses as a QRS complex, the monitor has a function to suspend QRS detection for a fixed amount of time starting from the detection of the pacing stimulus. This function is called "pace pulse mask". But if the pacemaker does not detect the patient's spontaneous heartbeat (sensing failure), and the pacing stimulus is applied at the same time as the QRS, this "pace pulse mask" function may erroneously mask the QRS and cause the heart rate measurement to decrease. To avoid this, QRS pace pulse mask function can be set to OFF, 10ms, or 20ms for correct measurement of the heart rate. (Default setting : Auto)

<b>∕∆</b> warning	If the QRS pace mask function is set to OFF, 10ms, or 20ms, the pace pulse may be erroneously be detected as a QRS complex and HR/Asystole Alarms may not generate due to incorrect HR (counting pace pulse as QRS complex). Select OFF, 10ms, or 20ms only if you are sure that pacing failure will not occur, or when the patient can be constantly monitored.
-------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



#### **1** Press the Config. $\rightarrow$ Page Down keys.

ECG Configuration 2/3	Page Up Page Down Prev. Disp.
Auto Lead Switch	OFF OFF
Pacemaker Pulse	
Pace Pulse Mask Time	Auto 10ms 20ms 40ms 0FF
AC Filter	ONOFF
ECG Drift Filter	
3lead Override	OFF

The configuration menu for selecting pace pulse mask time will be displayed.

#### **2** Select the pace pulse mask time.

Select from 10ms, 20ms, 40ms depending on the pace pulse amplitude or presence of fusion beat.

Selecting Auto will automatically select the pace pulse mask time from 20ms, 30ms, or 40ms depending on the pace pulse amplitude.

Selecting OFF will set the mask time to 0ms.

# **AC Filter**

#### **1** Press the Config. $\rightarrow$ Page Down keys.

ECG Configuration 2/3	Page Up Page Down Disp.
Auto Lead Switch	
Pacemaker Pulse	
Pace Pulse Mask Time	Auto 10ms 20ms
	40ms OFF
AC Filter	
ECG Drift Filter	
3lead Override	

The configuration menu for selecting the AC filter will be displayed.

#### **2** Select ON or OFF for the AC filter.

Selecting ON will set the AC filter which attenuates the AC noise (50 to 60Hz). Selecting OFF will not set the AC filter.

# **ECG Drift Filter**

By setting the ECG drift filter ON, only the amplitude with frequency component under 1Hz will be attenuated to prevent the ECG baseline to drift.

#### **1** Press the Config. $\rightarrow$ Page Down keys.

EC6 Configuration 2∕3	Page Up Page Down Disp.
Auto Lead Switch	OFF
Pacemaker Pulse	
Pace Pulse Mask Time	Auto 10ms 20ms 40ms 0FF
AC Filter	
ECG Drift Filter	
3lead Override	ON OFF

The configuration menu to select the ECG drift filter will be displayed.

#### 2 Select ON or OFF for the ECG drift filter.

Selecting ON will set the ECG drift filter and controls the baseline drift. When the ECG drift filter is set, the waveform will be delayed about 0.5 seconds. The message, "Drift-F ON" will be displayed on the home display. Selecting OFF will not set the ECG drift filter.

FL	KUDA DENSHI Adult 12/11 14:51 M Drift-F ON Message
	stead of "Drift-F ON" message, the enlarged time can be displayed depending on the setup. or procedures, refer to "8. System Configuration Monitor Setup".
NOTE	<ul> <li>When the "ECG Drift Filter" is set to ON, the waveform display will be delayed about 0.5 seconds.</li> <li>When the defibrillation and electrosurgery-proof ECG relay cable (CI-700E-3 (FA), CI-700E-4 (FA), CI-700E-5 (FA)) is used, ECG Drift Filter will be always set to OFF.</li> </ul>

# **3-lead Override**

When a 5-lead or 10-lead relay cable is used with a 3-lead cable, the device will judge as lead-off condition and display the "LEAD OFF" message. This can be avoided by selecting ON for "3-lead Override".

NOTE If 4-lead, 5-lead, or 10-lead relay cable is used with all the lead cables and electrodes connected, but ON is selected for "3-lead Override", the device will aknowledge that only 3 electrodes are used and displays only one waveform. Also, artifact may interfere on the waveform, and lead-off condition cannot be correctly displayed. When setting ON for "3-lead Override", use only 3 electrodes of LA, RA, and LL.

- Press the Config.  $\rightarrow$  Page Down keys. 1 ECG Configuration 2/3 Page Down Disp. Page Up Auto Lead Switch ON OFF OFF Pacemaker Pulse ON Pace Pulse 20ms The configuration menu to set the "3-lead Override" will be . Auto Mask Time displayed. OFF 40ms AC Filter ON OFF ECG Drift Filter ON OFF 3lead Override ON OFF
- 2 If using the 3-lead cable, select ON.

# **Pace Pulse Detection**

The pacemaker pulse detection sensitivity level can be selected from Low, Medium, and High.



# **ON/OFF of Parameter Display**

1	Press the Display ON/OFF	] key.
	ECG Prev. Lead-Size ECG1:×1,    HR/PR Alarm ON 4 0 - 1 2 0 Graphic Arrhy.Alarm Arrhy.Learn Trend Config. Pace Pulse OFF Filter Nontor Pace Pulse Pulse Tone ON Drift Filter Rem Source Auto ACF Filter Stead Overnako OFF Auto Lead OFF Detection Hed.	The confirmation screen for ON/OFF of ECG display will appear.
2	Select Display ON or D	isplay OFF.         Pressing the Display ON key will display the waveform and numeric data.         Pressing the Display OFF key will not display the waveform and numeric data.
	BED-001 FUKUDA DENSH	Result: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Prodult: Product Product Pro

When ECG electrodes are attached to the patient with the ECG display set to OFF, the ECG waveform and numeric data will be automatically displayed after 10 seconds.

<b>▲</b> CAUTION	When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend input will also be suspended.
------------------	----------------------------------------------------------------------------------------------------------------------------------------

# Respiration

This menu allows setting the impedance respiration measurement,  $CO_2$  respiration measurement or ventilator respiration measurement.



Size: Selects the waveform size to display the impedance respiration.RR Alarm: Selects ON/OFF for respiration rate alarm, and sets the upper and lower alarm limits.APNEA Alarm: Selects ON/OFF for apnea alarm and sets the upper alarm limit.Configuration: Sets the respiration monitoring configuration.				
<b>▲</b> CAUTION	<ul> <li>When the following relay cables are used, respiration cannot be measured.</li> <li>Relay Cable CI-700E-3 (FA) (Defibrillation and electrosurgery-proof, 3-electrode)</li> <li>Relay Cable CI-700E-4 (FA) (Defibrillation and electrosurgery-proof, 4-electrode)</li> <li>Relay Cable CI-700E-5 (FA) (Defibrillation and electrosurgery-proof, 5-electrode)</li> <li>When a defibrillator is used during respiration monitoring, a large offset voltage will be placed on the ECG electrodes, which may cause interruption of monitoring for a few seconds.</li> </ul>			

# **Respiration Waveform Size**

**1** Press the Size key to display the size setup menu.



Select the waveform size and baseline position to display and record the impedance respiration waveform.

#### **2** Select the waveform size.

Select the size from	×1/4 ,	×1/2 ,	×1 ,	×2 ,	×4 .
----------------------	--------	--------	------	------	------

# 3 Set the baseline position using the $\frown$ , $\checkmark$ keys.

Adjust the baseline position for  $0\Omega$  if the waveform is hard to see due to the waveform amplitude.

#### **1** Press the RR Alarm key to display the alarm setup menu.



Set ON/OFF of the RR alarm and upper and lower alarm limit. The alarm will be set common to RR measured from the impedance respiration waveform or RR measured from the  $CO_2$  waveform.

The increment will be according to the "RR Alarm Increment" setting. (Normal / Small).

Normal		Small
Adult	5Bpm increment	1Bpm increment
Child/Neonate	2Bpm increment	1Bpm increment

Reference

For the RR alarm increment setup, refer to "8. System Configuration Monitor Setup".

If the alarm is based on the RR measured from CO₂ waveform, RR alarm will not generate unless 2 or more respiration is detected within 30 seconds after power ON or after discharge.

Key	Item	Description
ONOFF	Individual Alarm	Selecting ON will generate the RR alarm. Selecting OFF will not generate the RR alarm.
Lower	Lower Alarm Limit	Sets the lower alarm limit (5 to 145Bpm / 2 to 148Bpm). Setting the value to 5Bpm/2Bpm or below will turn OFF the alarm.
🗲 Upper 🔿	Upper AlarmSets the upper alarm limit (10 to 150Bpm / 4 to 150Bpm).LimitSetting the value to 150Bpm or above will turn OFF the alar	
Auto	Automatic Setup	Automatically sets the upper limit to +20Bpm, and the lower limit to -20Bpm from the current value.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# **APNEA** Alarm

#### **1** Press the **APNEA** Alarm key to display the alarm setup menu.



Set ON/OFF of the apnea alarm and upper limit of apnea time. Apnea will be set common to apnea time measured from the impedance respiration waveform or apnea time measured from the  $CO_2$  waveform.

The upper limit can be set in 1-second increment. There is no lower limit.

<ul> <li>The purpose of this apnea alarm is to alert the user to evaluate for the possible occurrence of apnea events by identifying the absence of respiration. It is not intended to be classified as an "Apnea Monitor" and will not identify the condition creating the possible event. (Central, Obstructive or Mixed.)</li> <li>When PURITAN-BENNETT Ventilator is used, APNEA alarm will not generate if ventilator is selected for RR/APNEA Alarm Source.</li> </ul>

Key	Item	Description
ON OFF	Individual Setup	Selecting ON will generate the apnea alarm.
		Selecting OFF will not generate the apnea alarm.
Upper	Upper Alarm Limit	Sets the upper alarm limit (5 to 20sec.). Setting the value equal to or above 20sec. will turn OFF the alarm.
Auto	Automatic Setup	Sets the apnea alarm value set for the currently selected alarm mode.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).

For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# **CVA Detection**

When the amplitude of the respiration waveform decreases due to causes such as respiratory pause, the ECG waveform may be superimposed on to the respiration waveform, making the RR equal to the HR. This condition is called CVA (Cardio-Vascular Artifact), and is detected using the CVA detection function.

If the ECG waveform is superimposed on to the respiration waveform, with HR (RR) 30bpm, for 20 seconds or over (10 seconds or over for neonates) and the CVA detection function set to ON, the "CVA detected" message will be displayed, and an alarm sound will be generated.

This function will be effective only when Impedance is set as the RR/APNEA alarm source.

#### **1** Press the Config. key on the RESP setup menu.

_	RESP	Prev.
C	CVA detect	
	RR/APNEA Alarm Source	Impedance     CO2       Ventilator     Auto
	Impedance measure	
	RR sync, Indicator	

The configuration menu to set the CVA detection will be displayed.

2 Select ON or OFF.

ON will generate an alarm and display a message when CVA is detected. OFF will not perform CVA detection.

# **RR/APNEA Alarm Source**

The parameter to display the RR synchronized mark and to generate the RR/APNEA alarm can be selected from impedance RR,  $CO_2$  RR, and ventilator RR.

#### **1** Press the Config. key on the RESP setup menu.



The configuration menu to select RR source, APNEA source will be displayed.

#### **2** Select the parameter.

 Impedance
 CO2

 Ventilator
 Auto

Selecting Impedance will generate RR alarm based on the impedance respiration curve. Impedance synchronized mark will be displayed.

Selecting  $\boxed{CO_2}$  will generate RR alarm based on the  $CO_2$ 

waveform.  $CO_2$  synchronized mark will be displayed.

Selecting Ventilator will generate RR alarm based on RR measured by the ventilator. Ventilator synchronized mark will be displayed.

Selecting Auto will automatically select the measurable parameter in the priority of  $CO_2$ >ventilator>impedance, and generates alarm if the corresponding parameter key is displayed on the home display.

<b>M</b> WARNING	The RR/APNEA alarm will not be generated unless the parameter key corresponding to the selected RR/APNEA source is displayed. Be sure to display the parameter key for the RR/APNEA source.	
▲ CAUTION	In case of DS-LANII network, if the RR/APNEA alarm source is other than Impedance (Or, if Auto selects a setting other than impedance for RR/APNEA alarm source), the respiration waveform will not be transmitted on the network. In addition, if the RR/APNEA alarm source is other than $CO_2$ (Or, if Auto selects a setting other than $CO_2$ for RR/APNEA alarm source), the $CO_2$ waveform will not be transmitted on the network. In case of DS-LANIII network, refer to the operation manual for the central monitor.	

# **Impedance Respiration Measurement**

The respiration measurement using the impedance method conducts high-frequency and weak current between the ECG electrodes attached to the patient, and measures the potential difference between the electrodes caused by thoracic movement using the synchronous rectification system. For a patient using the adaptive (minute ventilation) pacemaker, the pacemaker measurement signal and the high-frequency current of this equipment interferes with each other which causes incorrect respiration measurement.

If the patient is using an adaptive (minute ventilation) pacemaker, set the impedance respiration measurement OFF.

#### **1** Press the Config. key on the RESP setup menu.



The configuration menu to set the impedance respiration measurement will be displayed.

#### 2 Select ON or OFF.

ON will perform standard impedance respiration measurement.

OFF will stop the impedance respiration measurement and will not display the impedance respiration waveform and RR. A high frequency electric discharge which is a measurement signal will be also ceased.

# **RR Synchronized Indicator**

The RR mark synchronized to impedance respiration or CO₂ waveform will be displayed inside the parameter key.

NO	тс	If PURITAN-BEN	NETT Vent	ilator is used, RR synchronized indicator will not be
	OTE	displayed when	Ventilator	is selected for "RR/APNEA Alarm Source".

#### **1** Press the Config. key on the RESP setup menu.

RESP Configuration CVA detect	ON OFF
RR/APNEA Alarm Source	Impedance     Impedance       Ventilator     Auto
Impedance measure	ON
RR sync, Indicator	

The configuration menu to set the RR synchronized indicator will be displayed.

#### 2 Select ON or OFF.



RR. JIP RR Sync. Indicator

# **ON/OFF of Parameter Display**

**1** Press the Display ON/OFF key.



# 2 Select Display ON or Display OFF.

RESP Prov. Disp. RESP display can be turned ON or OFF.	Pressing the Display ON key will display the waveform and numeric data. Pressing the Display OFF key will not display the waveform and numeric data.
Display ON Display OFF	Etc02       33         InspC02       1         RR_ITP       The Display OFF message will be displayed inside the parameter key.         Disp.OFF       9:19         \$129/0       82         metric       metric

When ECG electrodes are attached to the patient with the respiration display set to OFF, the respiration waveform and numeric data will be automatically displayed after 10 seconds.

When the waveform and numeric data display is set to OFF, the alarm ▲CAUTION generation and tabular/graphic trend will also be suspended.
# Invasive Blood Pressure (BP1 to BP5)

This menu allows setting the measurement condition for BP1 to BP5.

BP1	BP1	Prev. Disp.
Scale	0 - 200 mmHg	Graphic Trend Tabular Trend
BP Alarm	ONS 80 - 180 DOFF - OFF MOFF - OFF	
BP Zero	BP zero drift	Display ON/OFF
Config.	Filter 12Hz	ON
	Display Type S/D/M Mean Wave OFF Resp.Reject OFF	

Scale BP Alarm	<ul> <li>Selects the scale for BP waveform display.</li> <li>Sets the upper and lower alarm limit of systolic, diastolic, mean blood pressure and ON/OFF of the alarm.</li> </ul>
BP Zero	: Performs zero balance.
Configuration	: Sets the BP monitoring condition.

▲ CAUTION	The BP value will not be displayed until the zero balance is performed after discharging or when the main power is turned ON. If BP interface cable or transducer is not connected for 5 minutes or more, the zero balance information will be cleared. Make sure to perform the zero balance. However, if the power is turned ON within 5 minutes after the power is turned OFF, the previous zero balance information will be maintained, and the BP value will be displayed.
-----------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# **BP Label**

#### 1 Press the Label key.



The BP label setup menu will be displayed. Select the BP label for display and recording.

#### 2 Select the label.

Select from	BP*,	ART,	RAP ,	RVP,	PAP ,	CVP,	ICP,	UAP ,	IAP ,	LAP ,	LVP,
US1, US	2.										

#### [The Description of Each Label]

- ART (Arterial Pressure)
- RAP (Right Atrial Pressure)
- RVP (Right Ventricular Pressure)
- PAP (Pulmonary Artery Pressure)
- CVP (Central Venous Pressure)
- ICP (Intra-cranial Pressure)
- UAP (Umbilical Artery Pressure)
- IAP (Intra-aortic Balloon Pumping Pressure)
- LAP (Left Atrial Pressure)
- LVP (Left Ventricular Pressure)
- US1 (BP User Label 1)
- US2 (BP User Label 2)

Reference

Refer to "8. System Configuration Label Setup" for procedure to set the BP user label.

#### When the BP Label is IAP

When the BP label is IAP, PDP (Peak Diastolic Pressure) will be displayed in addition to systolic, diastolic, and mean pressure.

Please be aware that the Systolic Pressure (SYS) = Peak Systolic Pressure (PSP)



<ul> <li>Note that the Systolic Pressure (SYS) = Peak Systolic Pressure (PSP) for the graphic trend, data base, and alarm setup.</li> </ul>
<ul> <li>When ECG is not measured, PDP cannot be calculated.</li> </ul>

#### When the BP Label is CVP

When the BP label is CVP, the measurement unit can be selected from "mmHg", "kPa" or "cmH₂O". The measurement unit can be selected on the preset menu. The selected unit will be displayed on the BP parameter key.



Refer to "8. System Configuration Hospital Setup" for CVP measurement unit setup.

#### When the BP Label is ICP

Reference

When the BP label is ICP, labeling the artery pressure as ART will allow measuring the CPP (Cerebral Perfusion Pressure). CPP = Mean Arterial Pressure - Intracranial Pressure. If CPP is a negative value, the value will not be displayed. Also, the alarm cannot be set for CPP.



#### **PCWP Measurement**

When PAP is set as BP label, the mean value can be displayed as PCWP (Pulmonary Capillary Wedge Pressure).

#### 1 Set the BP label to PAP, and display the BP menu.



If the BP label is set to PAP, the PCWP key will be displayed on the BP menu.

# 2 Press the PCWP key.



The PCWP measurement display will appear.

#### **3** Press the WAVE READ key.

PCWP waveform and respiration waveform will be displayed. The cursor point indicates the current mean pressure.

#### **4** Use the **↑**, **↓** keys to manually set the PCWP value.

Move the cursor to manually set the PCWP value. Select the waveform scale from 20 or 50mmHg.

#### 5 Press the INPUT key.

Press the INPUT key when the PCWP value is correctly set. The PCWP value will be displayed inside the PAP (BP label) parameter box with the measurement time. Also, the value will be input to the graphic trend.



### **BP Scale**

#### 1 Press the Scale key



The BP scale setup menu will be displayed. Select the full scale for displaying and recording. The scale selection will differ depending on the label.

BP Label	Scale
BP1 to BP5, UAP,	20, 50, 75, 100, 150, 200, 250, 300 (mmHg)
User Label	4, 8, 12, 16, 20, 24, 32, 40 (kPa)
ART, LVP, IAP	50, 75, 100, 150, 200, 250, 300 (mmHg)
	8, 12, 16, 20, 24, 32, 40 (kPa)
PAP, RVP	20, 40, 50, 75, 100, 150 (mmHg)
FAF, NVF	4, 6, 8, 12, 16, 20 (kPa)
	10, 20, 30, 40, 50 (mmHg)
CVP, RAP, LAP	2, 4, 5, 6, 8 (kPa)
	20, 40 (cmH ₂ O) * only for CVP
ICP	5, 10, 15, 20, 50, 100 (mmHg)
	1, 2, 3, 4, 8, 16 (kPa)

#### **2** Select one scale from the displayed selection.

<b>A</b> CAUTION	The BP waveform with a scale above the programmed scale will not be properly transmitted on a wireless network. Select a proper scale for the waveform.
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------

#### **BP** Alarm

#### **1** Press the BP Alarm key.



The BP alarm setup menu will be displayed. Select ON/OFF for BP alarm and set the upper and lower alarm limit for systolic (S), diastolic (D), and mean (M) BP.

The alarm value should be set for each unit. (mmHg / kPa) The adjustable increment will be according to the "BP Alarm Increment" setting. (Normal / Small). The adjustable increment for upper and lower limit changes from

The adjustable increment for upper and lower limit changes from 50mmHg / 7kPa.

	"BP Alarm Increment" Setup					
	When Normal When Small					
0 to 50mmHg	2mmHg increment	1mmHg increment				
50 to 300mmHg	5mmHg increment					
0 to 7kPa	0.2kPa increment	0.1kPa increment				
7 to 40.0kPa	0.5kPa increment	0. TKPa Increment				



For the BP alarm increment setup, refer to "8. System Configuration Monitor Setup".

Key	ltem	Description		
ONOFF	Individual Alarm	SelectingONwill generate BP alarm.SelectingOFFwill not generate BP alarm.		
S 		Select from S (systolic BP), D (diastolic BP), M (mean BP).		
	Lower Alarm Limit	Sets the lower alarm limit (0 to 295mmHg / 0 to 39.5kPa). Setting the value equal to or below 0mmHg/0kPa will turn OFF the alarm.		
Upper	Upper Alarm Limit	Set the upper limit (2 to 300 mmHg / 0.2 to 40.0kPa). Setting the value equal to or above 300 mmHg / 40.0kPa will turn OFF the alarm.		
Auto	Automatic Setup	When the BP label is BP1/ART, the upper and lower lim will be automatically set to +40mmHg / +5kPa ar -20mmHg / -3kPa respectively from the current value. When the BP label is other than BP1/ART, the upper ar lower limit will be automatically set to the current value +20%, -20% respectively from the current value.		

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).

For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# Zero Balance of Pressure Lines

**1** Open the three-way valve of the pressure transducer to air.

#### **2** Press the BP zero key.



Verify that the BP waveform is positioned at zero, and "0" is displayed for the BP value. A message, "BP zero complete" will be displayed when the procedure is completed. A message, "BP zero failed" will be displayed when the process fails. The three-way valve may not be opened to air, artifact is present, or the transducer may be defective. Check the cause and try the zero balance procedure again.

A message, "BP zero drift" will be displayed when the interface cable is not connected. Check if the cable is firmly connected.

**3** Close the three-way valve when the zero balance is completed.

<b>∧</b> CAUTION	Each time the blood pressure transducer or tubing is replaced, the zero balance procedure is required to ensure accurate measurements.
	balance procedure is required to ensure accurate measurements.

Reference

### **Filter Selection**

Noise may interfere on the BP waveform depending on the combination of BP measurement circuit. Select the appropriate filter from the low-pass filter of 6Hz, 8Hz, 12Hz, or 40Hz.

#### 1 Press the Config. key.

BP1 Configuration	Prev. Disp.
Filter	6Hz 8Hz 12Hz 40Hz
HR/PR Alarm Source	ECG SpO ₂ BP
Display Type	S/D/M S/D M
Mean Wave	
Resp.Reject	OFF

The configuration menu to select the filter will be displayed.

#### 2 Select the filter.

Select the appropriate filter from 6Hz, 8Hz, 12Hz, 40Hz.

# HR/PR Alarm Source (BP1 or ART)

The parameter to display the HR synchronized mark and to generate the HR/PR alarm can be selected from ECG, SpO₂, or BP (BP1 or ART). If BP1 and ART are measured simultaneously, ART will be prioritized.

NOTE	This setting will be displayed only when	ECG/SpO ₂ /BP	is selected for "HR/PR
NOTE	Alarm Source" on the Monitor Setup (2/4	) menu.	

# 1 Press the Config. key.



The configuration menu to select the HR/PR alarm source will be displayed.

#### **2** Select a parameter.

ECG SpO₂ BP BF AR Auto Hoy

ECG will generate the alarm based on HR measured from ECG. HR synchronized mark will be displayed.

 $SpO_2$  will generate the alarm based on PR measured from SpO₂. SpO₂ synchronized mark will be displayed.

BP will generate the alarm based on PR measured from BP (BP1 or ART). BP synchronized mark will be displayed.

However BP can be selected only when ECG/SpO₂/BP is selected for "HR/PR Alarm Source" of the monitor setup menu.

Auto will automatically select the measurable HR source in the priority of ECG>SpO₂>BP. If the corresponding parameter key is displayed, alarm generation will be also effective.

▲WARNING	<ul> <li>The HR/PR alarm will not be generated unless the parameter key corresponding to the selected HR/PR source is displayed. Be sure to display the parameter key for the HR/PR source.</li> <li>The alarm for the parameter not selected for the "HR/PR Alarm Source" (ECG/SpO₂/BP) will be set to OFF on the DS-7600 Central Monitor.</li> <li>The "HR/PR Alarm Source" setting will synchronize between the bedside monitor and the central monitor.</li> <li>For example, if PR is set as the HR/PR alarm source on the DS-7200, HR</li> </ul>
	alarm will be set to OFF on the central monitor.
▲ CAUTION	In case of DS-LANII network, if the HR/PR source is BP (Or, if Auto selects BP for HR/PR source), the ECG waveform will not be transmitted on the network. On the central monitor, PR_IBP value will be displayed instead of HR. However, on some central monitor depending on the model type, the HR value from ECG will be displayed on the NIBP list and ST measurement list. Refer to the operation manual for the respective central monitor. In case of DS-LANIII network, refer to the operation manual for the central monitor.
ΝΟΤΕ	<ul> <li>If the HR/PR alarm source is BP, the PR will be displayed as "" if the corresponding BP (ART or BP1) is not measured.</li> <li>If <u>SpO</u>₂ is selected for "Pulse Tone" (ECG config.), the synchronized mark will be displayed inside the PR_SpO₂ numeric data box regardless of the HR/PR alarm source setup.</li> </ul>

# **Display Selection of Numeric Data**

The display type of BP numeric data can be selected from systolic / diastolic / mean, systolic / diastolic, or mean BP. The BP alarm will not be generated unless the data is displayed.

If the BP label is CVP, IAP, ICP or PAP, the display type is fixed. The selection is possible only for the BP labels other than above.

#### 1 Press the Config. key.

BP 1 Configuration	Prev. Disp.
Filter	6Hz 8Hz 12Hz 40Hz
HR/PR Alarm Source	ECG SpO ₂ BP
Display Type	Б/D/М
Mean Wave	OFF

The display type selection will be displayed.

#### **2** Select the display type.



S/D/M will display the systolic / diastolic / mean pressure.



S/D will display the systolic / diastolic pressure.



M will display only the mean pressure.

<b>▲</b> CAUTION	The undisplayed BP data (SYS/DIA/Mean) will not generate a BP alarm or be displayed in the tabular trend. Select the appropriate display type according to the monitoring purpose.
------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# Mean BP Waveform Display

The mean BP waveform can be selected to be continuously displayed on the home display.

**1** Press the Config. key.



The mean waveform display selection will be displayed.

2 Select ON/OFF of mean BP waveform display.



Selecting ON will display the mean BP waveform and "MEAN_WAVE" will be displayed inside the parameter key.

# **Respiration Rejection Filter**

The BP waveform baseline drift caused by the respiration influence can be prevented by setting the respiration rejection filter ON.

**1** Press the Config. key.



The respiration rejection filter selection will be displayed.

**2** Select ON or OFF for the respiration rejection filter.

# **ON/OFF of Parameter Display**

CAUTION

•

**1** Press the Display ON/OFF key. BP 1 Prev. Disp. Label BP 1 Graphic Trend 0 - 200 mmHg Scale Tabular Trend 80 - 180 OFF - OFF OFF - OFF BP Alarm ON S D The confirmation screen for ON/OFF of BP display will appear. M OFF BP Zero BP zero drift Display ON/OFF ON Filter Config. 12Hz S/D/M Off Off Display Type Mean Wave Resp.Reject A 2 Select Display ON or Display OFF. Prev. Disp. BP1 display can be turned ON or OFF. Pressing the Display ON key will display the waveform and numeric data. Display OFF Display ON Pressing the Display OFF key will not display the waveform and numeric data. 60 The Display OFF message will be displayed inside the parameter key. Disp, ÓFF nHa 23/ 10 When waveform and numeric data display is set to OFF, the alarm 

generation and tabular/graphic trend will be suspended.

OFF, the pulse rate derived from BP will not be displayed either.

If the display of waveform / numeric data labeled as BP1 or ART is set to

# SpO₂

This menu allows setting the SpO₂ monitoring condition.



Size : Sets the SpO₂ waveform size.

SpO₂ Alarm : Sets ON/OFF of the alarm, upper and lower alarm limit, and SEC alarm. Configuration : Sets the SpO₂ monitoring configuration.

▲ CAUTION	<ul> <li>Take the following precautions when monitoring over a long period of time.</li> <li>To avoid skin rash or low-temperature burn, it is recommended to change the measurement position several times a day, which is specified for each SpO₂ sensor. Be especially careful when continuously using on neonates, infants, or patients with peripheral circulatory disturbance.</li> <li>Excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.</li> <li>When not performing the measurement, unplug the relay cable and sensor from the SpO₂ connector. Otherwise, the measurement data may be erroneously displayed by the ambient light.</li> <li>The pulse wave is normalized for SpO₂ measurement. It does not indicate perfused blood volume. Check proper probe attachment by observing the pulse wave.</li> <li>Precautions for Reusable Type Sensor (DS-100A)</li> <li>The DS-100A is intended for use on finger of adults weighing over 40 kg (approximate). Do not use them on children or neonates. Also do not apply them on the thumb or toe.</li> <li>The DS-100A must be moved to a new site at least every 4 hours. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients. If skin integrity changes, move the sensor to another site.</li> <li>Precautions for Single-Patient-Use Type Sensors</li> <li>Do not wind the tape too strong. It may obstruct the blood flow.</li> <li>The sensor can be reused on the same patient as long as the adhesive tape attaches without slippage. But do not reuse on other patients. It is intended for single patient use only.</li> <li>For the sensor to another site.</li> <li>Do not wind the tape too strong. It may obstruct the blood flow.</li> <li>The sensor can be reused on the same patient as long as the adhesive tape attaches without slippage. But do not reuse on other patients. It is intended for single patient use only.</li> <li>For the single pati</li></ul>

▲ CAUTION	<ul> <li>Measuring on a limb with NIBP cuff, arterial catheter, or intracatheter may result in incorrect measurement.</li> <li>For additional warnings, cautions, or contraindications when using sensors with the DS-7210 Nellcor[®] model, refer to each SpO₂ sensor instruction manual.</li> <li>If measurement failure occurs due to the reason such as sensor detachment from the patient, measurement data will be displayed as "". Be cautious as numeric data alarm will not generate in such case.</li> </ul>
-----------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# SpO₂ Waveform Size

1 Press the Size key.



#### **2** Select the waveform size.

Select the size from	×1/4 ,	×1/2 ,	×1 ,	×2 ,	×4 .
----------------------	--------	--------	------	------	------

# SpO₂ Alarm

#### **1** Press the $SpO_2$ Alarm key to display the alarm setup menu.

Select ON/OFF of SpO₂ alarm, and set the upper and lower alarm limit. Also, when the SpO₂ value is unstable around the lower alarm limit, the frequently generated alarm can be corrected by setting the SEC (second) alarm function.



Reference	Refer to "4. Monitoring Se details of SEC alarm setu	etup SpO ₂ SEC Alarm Setup" for procedure.

The upper and lower limits can be set in 1% increment.

Key	ltem	Description		
ONOFF	Individual Alarm	SelectingONwill generate the SpO2 alarm.SelectingOFFwill not generate the SpO2 alarm.		
← Lower →	Lower Alarm Limit	Sets the lower alarm limit (50 to 99%). Setting the value below 50% will turn OFF the alarm.		
🗲 Upper 🗭	Upper Alarm Limit	Sets the upper alarm limit (51 to 100%). Setting the value above 100% will turn OFF the alarm.		
Auto	Automatic Setup	Automatically sets the upper limit to OFF, and the lower limit to 90%.		

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

### **HR/PR Alarm Source**

The parameter to display the HR synchronized mark and to generate the HR/PR alarm can be selected from ECG, SpO₂, and BP (BP1 or ART).

ΝΟΤΕ	The BP key will be displayed only when ECG/SpO ₂ /BP is selected for "HR/PR	
NOTE	Alarm Source" on the Monitor Setup (2/4) menu.	

#### **1** Press the Config. key.

Sp02 Configuration	Prev. Disp.
HR/PR Alarm Source	ECG     SpO2     BP       Auto
Ignore NIBP	ON OFF

The configuration menu to select the HR/PR source will be displayed.

#### 2 Select a parameter.

ECG SpO2 BP

Selecting ECG will generate the alarm based on HR measured from ECG. HR synchronized mark will be displayed.
Selecting SpO₂ will generate the alarm based on PR measured from SpO₂. SpO₂ synchronized mark will be displayed.
Selecting BP will generate the alarm based on PR measured from BP (BP1 or ART). BP synchronized mark will be displayed.
BP can be selected only when ECG/SpO₂/BP is selected for "HR/PR Alarm Source" of the monitor setup menu.
Selecting Auto will automatically select the measurable HR source in

the priority of ECG>SpO₂>BP. If the corresponding parameter key is displayed, alarm generation will be also effective.

▲WARNING	<ul> <li>The HR/PR alarm will not be generated unless the parameter key corresponding to the selected HR/PR source is displayed. Be sure to display the parameter key for the HR/PR source.</li> <li>The alarm for the parameter not selected for the "HR/PR Alarm Source" (ECG/SpO₂/BP) will be set to OFF on the DS-7600 Central Monitor.</li> <li>The "HR/PR Alarm Source" setting will synchronize between the bedside monitor and the central monitor.</li> <li>For example, if PR is set as the HR/PR alarm source on the DS-7200, HR alarm will be set to OFF on the central monitor.</li> </ul>
----------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<b>▲</b> CAUTION	In case of DS-LANII network, if the HR/PR source is BP (Or, if Auto selects BP for HR/PR source), the ECG waveform will not be transmitted on the network. On the central monitor, PR_IBP value will be displayed instead of HR. However, on some central monitor depending on the model type, the HR value from ECG will be displayed on the NIBP list and ST measurement list. Refer to the operation manual for the respective central monitor. In case of DS-LANIII network, refer to the operation manual for the central monitor.
NOTE	If $[SpO_2]$ is selected for "Pulse Tone" (ECG config.), the synchronized mark will be displayed inside the PR_SpO ₂ numeric data box regardless of the HR/PR alarm source setup.

# SpO₂ Alarm during NIBP Measurement

This setup is to be made when the  $\mbox{SpO}_2$  sensor and NIBP cuff is placed on the same limb for measurement.

During the NIBP measurement, the cuff inflation restricts the blood flow which disables the correct detection of the  $SpO_2$  value and PR, and may generate an improper alarm.

Selecting OFF will not generate the alarm until the NIBP measurement is complete. Similarly, when the HR source is set as  $SpO_2$ , the PR alarm will not be generated during NIBP measurement.

#### **1** Press the Config. key.



The configuration menu to set "Ignore NIBP" will be displayed.

#### 2 Select ON or OFF.

ON will generate the alarm during NIBP measurement.

OFF will not generate the SpO₂/PR alarm during NIBP measurement.

# **ON/OFF of Parameter Display**

**1** Press the Display ON/OFF key.



### 2 Select Display ON or Display OFF.



When  $SpO_2$  sensor is attached to the patient with the  $SpO_2$  display set to OFF, and  $SpO_2$  is measured for 10 seconds, the  $SpO_2$  waveform and numeric data will be automatically displayed.

<b>▲</b> CAUTION	<ul> <li>When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend input will also be suspended.</li> <li>When the waveform and numeric data display is set to OFF, the pulse rate derived from SpO₂ will not be displayed either.</li> </ul>
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# SpO₂

# (Masimo[®] Model: DS-7210M)

This menu allows the setup of SpO₂ monitoring condition for DS-7210M (Masimo^{$\degree$} Model).

This menu allows t	he setup of SpO ₂ monitoring condition for DS-7210M (Masimo ^{$\circ$} Model).
\$p02	Prev. Disp.
Size x 1	
Size x 1	Graphic Trend
SpO ₂ alarm ON	Tabular Trend
SpO ₂ alarm ON 9 0 -	OFF
Config.	Display
Alarm Source Auto Puls	se Sens. Low
-	T SAT OFF
	Xaplay ON Nal IQ Have OFF
0	
	Sets the pulse wave size
	Sets ON/OFF of SpO ₂ alarm, upper and lower alarm limit. Sets the SpO ₂ monitoring condition.
Configuration . C	
	• Take the following processions when menitoring ever long periods of time
	• Take the following precautions when monitoring over long periods of time.
	<ul> <li>To avoid skin rash or low-temperature burn, it is recommended to shange the measurement position several times a day, which is</li> </ul>
	change the measurement position several times a day, which is specified for each SpO ₂ sensor.
	Be especially careful when continuously using on neonates, infants, or
	patients with peripheral circulatory disturbance.
	Excessive light may cause inaccurate measurements. In such cases,
	cover the sensor with opaque material.
	<ul> <li>When not performing the measurement, unplug the relay cable and</li> </ul>
	sensor from the SpO ₂ connector. Otherwise, the measurement data
	may be erroneously displayed by the ambient light.
	• The pulse wave is normalized for SpO ₂ measurement. It does not
	indicate perfused blood volume. Check proper probe attachment by
	observing the pulse wave.
	<ul> <li>Precautions for Reusable Type Sensor (LNOP[®] DCI)</li> <li>The light amitting part of the sensor abound be over the root of the</li> </ul>
	<ul> <li>The light-emitting part of the sensor should be over the root of the fingernail. Do not insert the finger too far into the sensor as it may hurt</li> </ul>
	the patient.
	<ul> <li>Precautions for Single-Patient-Use Type Sensors</li> </ul>
	<ul> <li>Do not wind the tape too strong. It may obstruct the blood flow.</li> </ul>
	The sensor is contraindicated for use on patients who exhibit allergic
<b>A</b> CAUTION	<ul> <li>The Masimo[®] LNOP sensor can be reused on the same patient as long</li> </ul>
	as the light emitting and receiving part is clean, and if it is still adhesive
	to the skin. But do not reuse it on other patients. It is intended for single
	patient use only.
	<ul> <li>Do not reuse the sensor by resterilizing it.</li> <li>Dispass the sensor of the uses in the superior of demose to the starile</li> </ul>
	<ul> <li>Dispose the sensor after use. In the event of damage to the sterile packaging, do not use it.</li> </ul>
	<ul> <li>Change the sensor attachment site every 4 hours for the reusable sensor,</li> </ul>
	and every 8 hours for the single patient use type sensor. Exercise extreme
	caution with poorly perfused patients; skin erosion and pressure necrosis
	can be caused when the sensor is not frequently moved. Assess site at
	least every 2 hours with poorly perfused patients.
	• Measuring on a limb with NIBP cuff, arterial catheter, or intracatheter may
	result in incorrect measurement.
	• For additional warnings, cautions, or contraindications when using
	sensors with the DS-7210M Masimo [®] model, refer to each SpO ₂ sensor
	instruction manual.
	<ul> <li>The measurable pulse rate range is 25 to 240bpm. However, "xxx" will be displayed if 25bpm and below or 240bpm and above is measured.</li> </ul>
	<ul> <li>displayed if 25bpm and below or 240bpm and above is measured.</li> <li>If measurement failure occurs due to the reason such as sensor</li> </ul>
	detachment from the patient, measurement data will be displayed as
	" $ -$ ". Be cautious as numeric data alarm will not generate in such case.
	. De datione de hamene data alarm wil het generate in such case.

NOTE	The pulse wave for the Masimo [®] model (DS-7210M) will be displayed with approximately 630msec delay from the actual pulse.
	approximately boomsec delay nom the actual pulse.

#### **Pulse Wave Size**

**1** Press the Size key.



Select the waveform size for displaying and recording from  $[\times 1/4]$ ,  $[\times 1/2]$ ,  $[\times 1]$ ,  $[\times 2]$ ,  $[\times 4]$ .

# SpO₂ Alarm

#### **1** Press the $SpO_2$ Alarm key to display the alarm setup menu.

Select ON/OFF of SpO₂ alarm, and set the upper and lower alarm limit.

Sp02 Alarm SpO ₂ Alarm		ON	]	)FF	Prev. Disp.
Lower		Auto	] [	🗲 Upp	oer 🔿
70 75	80	85	90	95	
Lower		%		Upper	
90	÷	98	÷	OFF	
					Recall

The upper and lower limits can be set in 1% increment.

Key	Item	Description
ONOFF	Individual Alarm	Selecting ON will generate the SpO ₂ alarm. Selecting OFF will not generate the SpO ₂ alarm.
Lower	Lower Alarm Limit	Sets the lower alarm limit (50 to 99%). Setting a value 50% or below will turn OFF the alarm.
Upper 🔿	Upper Alarm Limit	Sets the upper alarm limit (51 to 100%). Setting a value above 100% will turn OFF the alarm.
Auto	Automatic Setup	Automatically sets the upper limit to OFF, and the lower limit to 90%.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# SpO₂ Monitoring Condition Setup

Pressing the Config. key on the  $SpO_2$  setup menu will allow to set the  $SpO_2$  measurement condition.

#### HR/PR Alarm Source

The parameter to display the HR synchronized mark and to generate the HR/PR alarm can be selected from ECG, SpO₂, and BP (BP1 or ART).

NOTE The BP key will be displayed only when ECG/SpO₂/BP is selected for "HR/PR Alarm Source" on the Monitor Setup (2/4) menu.

#### **1** Press the Config. key.



Selecting ECG will generate the alarm based on HR measured from ECG. HR synchronized mark will be displayed. Selecting  $SpO_2$  will generate the alarm based on PR measured from  $SpO_2$  SpO₂ synchronized mark will be

measured from  $SpO_2$ .  $SpO_2$  synchronized mark will be displayed. Selecting BP will generate the alarm based on PR

measured from BP (BP1 or ART). BP synchronized mark will be displayed.

BP can be selected only when ECG/SpO₂/BP is selected for "HR/PR Alarm Source" of the monitor setup menu.

Selecting Auto will automatically select the measurable HR source in the priority of  $ECG > SpO_2 > BP$ . If the corresponding parameter key is displayed, alarm generation will be also effective.

AWRNING	<ul> <li>The HR/PR alarm will not be generated unless the numeric data box corresponding to the selected HR/PR source is displayed. When HR/PR source selection is changed, be sure to display the numeric data box corresponding to the selected HR/PR source.</li> <li>The alarm for the parameter not selected for the "HR/PR Alarm Source" (ECG/SpO₂/BP) will be set to OFF on the DS-7600 Central Monitor.</li> <li>The "HR/PR Alarm Source" setting will synchronize between the bedside monitor and the central monitor.</li> <li>For example, if PR is set as the HR/PR alarm source on the DS-7200, HR alarm will be set to OFF on the central monitor.</li> </ul>
ACAUTION	In case of DS-LANII network, if the HR/PR source is BP (Or, if Auto selects BP for HR/PR source), the ECG waveform will not be transmitted on the network. On the central monitor, PR_IBP value will be displayed instead of HR. However, on some central monitor depending on the model type, the HR value from ECG will be displayed on the NIBP list and ST measurement list. Refer to the operation manual for the respective central monitor. In case of DS-LANIII network, refer to the operation manual for the central monitor.
NOTE	<ul> <li>If PR source is set to BP, make sure to measure the corresponding BP (BP1 or ART). Otherwise PR will be displayed as "".</li> <li>If <u>SpO</u>₂ is selected for "Pulse Tone" (ECG config.), the synchronized mark will be displayed inside the PR_SpO₂ numeric data box regardless of the HR/PR alarm source setup.</li> </ul>

#### ●SpO₂ Alarm during NIBP Measurement (Ignore NIBP)

This setup can be used when the  $SpO_2$  sensor and the NIBP cuff is placed on the same limb for measurement. During the NIBP measurement, the cuff inflation restricts the blood flow which disables the correct detection of  $SpO_2$  and may generate an improper alarm.

Selecting OFF for "Ignore NIBP" will not generate a SpO₂ alarm until the NIBP measurement is complete. Similarly, when the PR source is SpO₂, the PR alarm will not be generated during NIBP measurement.

#### 1 Press the Config. key.



ON will generate an alarm during NIBP measurement. OFF will not generate a SpO₂/PR alarm during NIBP measurement.

### ●SpO₂ Averaging

The averaging time for SpO₂ value can be selected.

<b>A</b> warning	E t
	c

Be cautious when setting the "SpO₂ Averaging" duration as the SpO₂ alarm is based on the displayed SpO₂ value which is averaged from the duration set in SpO₂ Averaging". The alarm occurrence time will be affected or may not occur for the transient value of SpO₂ depending on the set duration.

#### **1** Set the SpO₂ averaging time.

HR/PR Alarm Source ECG Sp02 BP Auto Ignore NIBP ON OFF Sp02 Averaging Z-4sec 4-6sec 8sec 10sec 12sec 14sec 16sec	Sp02 Configuration 1/2	Page Down Disp.
SpO ₂ Averaging 10sec 12sec 14sec 14sec	HR/PR Alarm Source	
10sec 12sec 14sec	Ignore NIBP	
	SpO ₂ Averaging	

								4-6sec	
8sec	/	10sec	/	12sec	/	14sec	/	16 sec	

# Pulse Wave Detection Sensitivity

The sensitivity to detect the pulse wave can be selected from high or low.

Sp <b>0</b> 2 Configuration 2/2	Page Up Prev. Disp.	]
^D ulse Sens,	High Low	
FAST SAT		For standard use, select Low.
^o I Display		To improve the low perfusion condition, or to perform fas tracking when the SpO ₂ value changes abruptly, select
Signal IQ Wave		High.

CAUTION If High is selected for pulse wave sensitivity, sensor-detached detection will become somewhat inaccurate.	วท
--------------------------------------------------------------------------------------------------------------------	----

# ●FAST SAT Setup

By selecting ON for "FAST SAT", abrupt change of the  $\ensuremath{\text{SpO}}_2$  value can be monitored.

	To pick up the abrupt change of the value sooner, and to take advantage of the
	qualities of FAST SAT mode, it is recommended to set 2-4sec for SpO ₂
	averaging time when FAST SAT is set ON.

#### **1** Select ON/OFF for "FAST SAT".



#### PI Display

Whether or not to display the PI (Perfusion Index) data can be selected.

The perfusion index is calculated by pulsatile signal  $\div$  apulsatile signal  $\times 100$ , and indicates patient's circulation condition. This can be used to find a good perfusion site to attach the sensor. Also, it can be used as diagnosis index to predict the patient's critical condition when at low perfusion.

#### **1** Select ON/OFF of PI Display.

(	Sp0z Configuration 2/2 Pulse Sens. FAST SAT PI Display Signal IQ Wave	Pase UP High Low ON OFF ON OFF ON OFF	Prev. Disp.	The setup menu to set the PI display will be displayed.
2	Select from Sp02 PI 7.00 Sp02	ON or OF 92 ^x 92 ^x	F]. [ON	
	<b>≜</b> CAUT			ected, "SpO ₂ Low Perfusion" alarm will be indicated by y only. The alarm sound will not be generated.

#### Signal IQ Wave Display

Whether or not to display the signal IQ wave can be selected.



The signal IQ wave indicates the signal force and pulse wave timing. The vertical length indicates the signal quality. A low vertical line

NOTE

The signal IQ wave cannot be recorded.

#### **1** Set ON/OFF of Signal IQ Wave display.

Sp02 Configuration 2/2	Page Up	Prev. Disp.	
Pulse Sens.	High Low		
FAST SAT	OFF		
PI Display			Select ON (display) or OFF (not display).
Signal IQ Wave	OFF		

indicates a bad signal quality.

# **ON/OFF of Parameter Display**



When  $SpO_2$  sensor is attached to the patient with the  $SpO_2$  display set to OFF, and  $SpO_2$  is measured for 10 seconds, the pulse wave and numeric data will be automatically displayed.

▲ CAUTION	<ul> <li>When the waveform and numeric data display is set to OFF, the alarm generation and list/trend input will also be suspended.</li> <li>When the waveform and numeric data display is set to OFF, the pulse rate derived from SpO₂ will not be displayed either.</li> </ul>
-----------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# **Non-Invasive Blood Pressure**

This menu allows setting the NIBP monitoring condition.

NIBP				Prev. Disp.
Auto Mod	de OFF			Graphic Trend
NIBP alar	m ON S D C	80 FF	— 180 — OFF	NIBP list
Config.				
	Quick Measurement	OFF	Sight Inflation	ON
	End Tone	ON	PR	OFF
	Mean	OFF	Oscill. Display	OFF
	Dyna Alert	OFF		
	Pump Setup	Norma	I	

Auto Mode	: Sets the automatic interval measurement and starts the 1-minute interval
	measurement and quick measurement.
NIBP Alarm	: Sets ON/OFF of the NIBP alarm and upper / lower limit of systolic, diastolic, and
	mean value of NIRP

mean value of NIBP. Configuration : Sets the NIBP monitoring configuration.

<b>▲</b> CAUTION	<ul> <li>Pay attention when measuring the NIBP of patient with bleeding disorders or hypercoagulation. The cuff inflation constricting the arm may cause petechia or circulatory failure with blood clot.</li> <li>Check the patient's condition constantly while measuring over a long period of time with interval of 2.5 minutes or less. Also, periodically check the blood circulation while performing periodic measurement over a long period of time. Congestion or rash may occur at the measuring site.</li> <li>For the following situation, measurements will be terminated.</li> <li>When the measurement time has exceeded 160 seconds for adult and child, 80 seconds for neonate.</li> <li>When the inflation value has exceeded 300mmHg for adult, 210mmHg for child, and 150mmHg for neonate.</li> <li>If used with the incorrect patient type, it will not only cause erroneous measurement, but the inflating level for the adult may be applied to child causing a dangerous situation to the patient.</li> </ul>

Non-invasive blood pressure can be measured automatically at selected time intervals. If Quick measurement is performed during the NIBP automatic measurement, the automatic measurement will automatically resume when Quick measurement completes.



When NIBP automatic measurement is set, the set interval time will be displayed inside the parameter key.

### 1 Press the Auto Mode key.



The measurement interval setup menu for the automatic measurement will be displayed.

#### 2 Select the interval time.

Select from 2min /	2.5min /	3min /	5min /	10min	/ 15min	/ 20min /
30min / 60min /	120min .					

Select OFF if not performing interval measurement.

The measurement starting time will be the integral multiple of the selected interval time beginning with 0 minute.

Ex.) If the present time is 13:14, the measurement starting time will be as follows for each interval time.

2 min. : 13:16, 13:18, 13:20, . . . 2.5 min. : 13:15, 13:17:30, 13:20, . . . 3 min. : 13:15, 13:18, 13:21, . . . 5 min. : 13:15, 13:20, 13:25, . . . 120 min. : 14:00, 16:00, 18:00, . . . (The measurement will start at every even hours.)

**Mathematical CAUTION** The NIBP measurement cannot be started from the central monitor via TCON system if the NIBP measurement interval is set to 2 min / 2.5 min / 3 min / 5 min. However, it can be stopped.

#### **3** Press Start early to start the measurement 5 minutes early.

NIBP Auto Mode				Prev. Disp.
	1min start			
Periodic	OFF			
2min		 3min	<b>5</b> min	]
<b>1</b> Omin	<b>15</b> min	20min	- 30min	
60min	120min		art early:	

If 60min or 120min is selected for interval time, the measurement will start 5 minutes early.

If outputting the data to PC or other external equipment using the PC communication function of this device, an error may be generated to the NIBP measurement time depending on the input interval of the external equipment. As this device outputs the data at completion of NIBP measurement, if the external equipment inputs the data at 60 minutes interval, 60 minutes time lag will occur. By starting the measurement 5 minutes early, this time lag between the external equipment can be minimized.

NOTE	When using the DS-LANIII network or TCON system and if "Timer" is selected for NIBP measurement on the central monitor, NIBP auto mode will be set OFF on the DS-7200, but the measurement will start at fixed time according to the central monitor setting.
------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

The NIBP measurement interval condition at discharge and at power ON can be set on the monitor setup menu.



- "NIBP Auto Mode" on "Backup at Discharge" menu
  - OFF / Backup / Backup (Resume auto mode by manual measurement)
- "NIBP Measurement Interval at Power ON" : According to Setup / 2.5min. when OFF

# **NIBP 1-Minute Interval Measurement**

The 1-minute interval measurement will automatically stop after 12 minutes and returns to the previous interval mode setup.

#### **1** Press the Auto Mode key.

NIBP Auto Mode		
1min start       Periodic       OFF       2min       2.5min       3min       5min	The NIBP auto mode menu will be displa	ayed.
Image: 2.5 min     Image: 3 min       Image: 10 min     15 min       Image: 2.5 min     Image: 3 min       Image: 10 min     120 min       Image: 2.5 min     Image: 3 min       Image: 3 min		

**2** Press the <u>1min Start</u> key to start the 1-minute interval measurement.

Pressing the NIBP START/STOP key will not stop the 1-minute interval measurement. To cancel the measurement, press the 1min Start key again.

ACAUTION	<ul> <li>The 1-minute interval measurement will always start from 00 second. Pressing the <u>1min Start</u> key will start the measurement from the next 00 second.</li> <li>The 1-minute interval measurement will automatically stop after 12 minutes and returns to the previous interval mode setup.</li> <li>The NIBP measurement cannot be started from the central monitor via TCON system during the 1-minute interval measurement. However, it can be stopped.</li> </ul>
----------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# **NIBP Alarm**

**1** Press the NIBP Alarm key.



The alarm setup menu will be displayed. Set ON/OFF of the NIBP alarm, upper and lower alarm limits of systolic (S), diastolic (D), mean (M) NIBP.

Set the alarm value for each measurement unit (mmHg / kPa). The upper and lower limit can be set in 5mmHg / 0.5kPa increment.

Key	Item	Description	
ONOFF	Individual Alarm	SelectingONwill generate a NIBP alarm.SelectingOFFwill not generate a NIBP alarm.	
S D M		Select from S (systolic), D (diastolic), or M (mean).	
Lower	Lower Alarm Limit	Sets the lower alarm limit (10 to 295mmHg/1.5 to 39.5kPa). Setting the value 10mmHg/1.5kPa or below will turn OFF the alarm.	
🗲 Upper 🔿	Upper Alarm Limit	Sets the upper limit (15 to 300mmHg/2.0 to 40.0kPa). Setting a value 300bpm/40.0kPa or above will turn OFF the alarm.	
Auto	Automatic Setup	Automatically sets the upper limit to +40mmHg/+5kPa from the current value, and the lower limit to -20mmHg/3kPa from the current value.	
<b>A</b> CALLER NIBP alarm will not generate if the alarm limit outside the NIBP measurement			

CAUTION NIBP alarm will not generate if the alarm limit outside the NIBP measurement range (10 to 280mmHg) is set.

To maintain the alarm setting even after the power is turned OFF or after a discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode". Whether or not to generate an alarm at NIBP measurement failure can be set. (Default: OFF) For details, refer to "4. Monitoring Setup Alarm Setup ON/OFF of Alarm Occurrence at NIBP Failure".

# **Quick Measurement**

The duration of quick measurement can be selected from 3 min., 5 min., or 10 min. The long duration of continuous measurement may congest the blood stream of the measured location. Set the duration according to the patient condition.

#### **1** Press the Config. key.



The NIBP configuration menu to set the Quick Measurement will be displayed.

#### **2** Select the measurement duration.

Select ON, OFF for quick measurement. ON will perform the measurement in duration of about 20 to 25 seconds in case of adult patient.

NOTE	The quick measurement setting will be effective only when the patient classification is adult or child. If the patient classification is neonate, standard NIBP measurement will be performed regardless of the quick measurement setting.
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# **End of Measurement Tone**

By selecting ON for the "End Tone", a tone will be generated when the NIBP measurement completes.

#### 1 Press the Config. key.



The NIBP configuration menu to set the "End Tone" will be displayed.

#### 2 Select ON or OFF.

ON will generate a tone when the measurement completes. OFF will not generate a tone when the measurement completes.

# Mean BP Display

The ON/OFF of mean BP display can be selected.

#### **1** Press the Config. key.



The NIBP configuration menu for setting the mean BP display will be displayed.

#### 2 Select ON or OFF.



ON will display the mean BP.

OFF will not display the mean BP.

**A**CAUTION If the mean BP display is set to OFF, the mean BP alarm will not be generated. Also, the mean BP will not be displayed for the tabular trend or the NIBP list function if the display is set to OFF.

# **Dyna Alert Function**

Using a cuff allows to measure the blood pressure noninvasively, but on the other hand, there is a demerit of not being able to perform the measurement continuously. Therefore, there is always a risk of sudden blood pressure change in between the periodic measurements.

The Dyna Alert function is a technology to prevent accidents which may occur by this sudden BP change during the non-measured duration by estimating the variation of circulatory dynamics using the parameters obtained from ECG and PTG (photoplethysmograph), and initiating a new NIBP measurement if a change in the circulatory dynamics is detected.

#### Parameters used for Dyna Alert Function



The Dyna Alert function will operate on the following measurement condition.

Patient Classification: Adult (Cuff Applied Site: UpperSpO2 Sensor Applied Site: FingerNIBP Measurement Interval: 5 to 60

: Adult (20kg or above) : Upper arm : Fingertip : 5 to 60 minutes

CAUTION When a PTG (SpO₂) sensor is applied to the toe or forehead, the Dyna Alert may not function depending on the patient's condition.

#### **1** Select ON/OFF of Dyna Alert function.



ON will set the Dyna Alert function ON. OFF will set the Dyna Alert function OFF.

# **Pump Setup**

The pump inflation operation can be selected from normal or silent operation.

#### **1** Display the second page of the NIBP setup menu, and set the "Pump Setup".



Silent will inflate with decreased speed to reduce the pump inflating sound.

# Sight Inflation

When "Sight Inflation" is set to ON, the maximum blood pressure level will be estimated during the inflation process, then 35mmHg will be added and the measurement will start. The inflation speed is slow but allows to detect any sudden increase of blood pressure to prevent re-inflation and discomfort for the patient.

When "Sight Inflation" is set to OFF, it will inflate to the target level set according to the previous measurement result.

NOTE	• The "Sight Inflation" will not function if the patient classification is "Neonate".
NOTE	<ul> <li>The "Sight Inflation" will not function during 1-minute measurement.</li> </ul>

#### **1** Select ON/OFF for "Sight Inflation".

NIBP Configuration 1/2		Page Down	Prev. Disp.
Quick Meas.	ON	OFF	
End Tone	ON	OFF	
MEAN	ON	OFF	
Dyna Alert	ON	OFF	
Pump Setup	Normal	Silent	
Sight Inflation	ON	OFF	
PR	ON	OFF	

ON will set the Sight Inflation ON.

OFF will set the Sight Inflation OFF.

# **Pulse Rate Display**

The pulse rate measured during NIBP measurement can be displayed. This pulse rate is only for display. It will not generate alarm, or be displayed for the list function.

#### **1** Press the Config. keys.



The configuration menu to set the pulse rate display will appear.



# **Oscillograph Display**

When the NIBP numeric data box size is 2-box size or larger, and "Oscillograph" is set to ON on the NIBP setup menu, the oscillation graph will be displayed inside the NIBP numeric data box.



The horizontal axis shows the cuff pressure, and vertical axis shows the pulse amplitude with reference to maximum pulse amplitude.

The bar graph shown at left indicates the size of maximum pulse amplitude compared with the reference value. For example, if the maximum pulse amplitude at measurement is 1/2 of the reference value, the bar graph will be half filled in.



#### **1** Press the Config. $\rightarrow$ Page Down keys.



# **Cancelling the NIBP System Error Message**

The displayed NIBP system error message can be cleared.





Refer to "10. Maintenance Troubleshooting" for the cause of air hose check.

# Temperature (T1 to T3)

This menu allows setting the monitoring condition for TEMP 1 to TEMP 3.

TEMP 1	Prev. Disp.
TEMP1 Label T1	Graphic Trend
TEMP1 Alarm OFF	Tabular Trend
TEMP2	
TEMP2 Label T2	
TEMP2 Alarm OFF	TEMP 1 Display ON/OFF ON
⊿T Display ON OFF	TEMP2 Display ON/OFF ON

TEMP* Alarm : Sets ON/OFF of the temperature alarm, and upper and lower alarm limits. TEMP* Label : Set the temperature measuring location. (* indicates T1 to T3)

# **Temperature Label**

**1** Press the **TEMP*** Label key.

TEMP 1 Label Prev. Disp.	
TEMP1 TEMP2 T1 T2 T1 Tsk Tre Tes Tco US1 US2	The temperature label setup menu will be displayed.

#### 2 Select a label.

Select from T*, Tsk, Tre, Tes, Tco, US1, US2.

#### [Description of Each Label]

T1-T3 (Default)

- Tsk (Skin Temperature)
- Tre (Rectal Temperature)
- Tes (Esophageal Temperature)
- Tco (Core Temperature))
- US1 (Temperature User Label 1)
- US2 (Temperature User Label 2)

Reference

Refer to "8. System Configuration Label Setup" for temperature user label setup.

# **Temperature Alarm**

**1** Press the **TEMP*** Alarm key.



The alarm setup menu for temperature will be displayed. Select ON/OFF of the temperature alarm, and set the upper and lower alarm limit.

The alarm limit can be set for each measurement unit (°C / °F). The upper and lower limit can be set in increments of 0.5°C / 1.0°F.

Кеу	Item	Description
ON OFF	Individual Alarm	Selecting ON will generate a TEMP alarm. Selecting OFF will not generate a TEMP alarm.
Lower	Lower Alarm Limit	Sets the lower alarm limit (30.0 to 49.0°C / 86.0 to 120.0°F). Setting the value 30.0°C / 86.0°F or below will turn the alarm OFF.
Upper	Upper Alarm Limit	Sets the upper alarm limit (31.0 to 50.0°C / 88.0 to 122.0°F). Setting the value 50.0°C / 122.0°F or above will turn the alarm OFF.
Auto	Automatic Setup	Automatically sets the upper limit to $+2.0^{\circ}C / +4.0^{\circ}F$ from the current value, and lower limit to $-2^{\circ}C / -4.0^{\circ}F$ from the current value.

To maintain the alarm setting even after the power is turned OFF or after a discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).

Reference

For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# $\Delta T$ Display



NOTE

 The temperature difference of TEMP1 and TEMP3, TEMP2 and TEMP3 cannot be displayed for ∠T.

• The alarm cannot be set for  $\angle T$ .

# **ON/OFF of Parameter Display**

**1** Press the Display ON/OFF key.



# CO₂ Concentration

# (Option Unit: MGU-721)

This menu allows to set the monitoring condition of  $CO_2$  concentration when mainstream  $CO_2$  unit (MGU-721) is used.

Reference

• For parameter setup of the microstream CO2 unit (MGU-722), refer to P6-63 "CO₂ Concentration (MGU-722)".

002			Prev. Disp.
Scale	50 mm	⊣g	
EtCO ₂ Alarm	ON 30 —	45	Graphic Trend Tabular Trend
InspCO ₂ Alarm	ON		
	з		
Config,	EtCO2 Peak Picking Duration Unit O2 Comp. N2O Comp. Anesthetic Comp. Atmos. Pressure	10sec mmHg 21% 0FF 0.0 % 760 mmH	Display ON/OFF ON
Cal. Airway Adpt			

Scale	: Sets the CO ₂ waveform scale.
EtCO ₂ Alarm	: Sets ON/OFF of the EtCO ₂ alarm, and upper and lower alarm limits.
InspCO ₂ Alarm	: Sets ON/OFF of the InspCO ₂ alarm and upper alarm limit.
Configuration	: Sets the CO ₂ monitoring conditions.

# CO₂ Scale

**1** Press the Scale key.

 Description
 Scale
 Image: Scale key.

 Scale
 Image: Scale key.
 The scale setup menu will be displayed.

 50
 100
 mmHg

 50
 100
 mmHg

 Scale setup menu for the unit in mmHg>

#### **2** Select the CO₂ waveform scale for displaying and recording.

```
For the measurement unit in mmHg, select the scale from 50, 100.
For the measurement unit in kPa and %, select the scale from 4, 8, 10.
```

# EtCO₂ (End-Tidal CO₂) Alarm

#### **1** Press the $EtCO_2$ Alarm key.



The alarm setup menu will be displayed.

Select ON/OFF of  $\ensuremath{\mathsf{EtCO}}_2$  alarm, and set the upper and lower alarm limits.

Alarm condition should be set for each measurement unit (mmHg / kPa / %).

Upper and lower alarm limits can be set in increments of 1mmHg, 0.1kPa, and 0.1%.

NOTE EtCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after power ON or after discharge.

Кеу	ltem	Description
	Individual Alarm	Selecting ON will generate the $EtCO_2$ alarm.
		Selecting OFF will not generate the EtCO ₂ alarm.
Lower	Lower Alarm Limit	Sets the lower alarm limit (1 to 98mmHg, 0.1 to 13.1kPa, 0.1 to 13.1%). Setting a value equal to or below 1mmHg, 0.1kPa, 0.1% will turn the alarm OFF.
🗲 Upper 🍑	Upper Alarm Limit	Sets the upper alarm limit (3 to 115mmHg, 0.3 to 15.0kPa, 0.3 to 15.0%). Setting a value equal to or above 115mmHg, 15.0kPa, 15.0% will turn the alarm OFF.
Auto	Automatic Setup	Automatically sets the upper alarm limit to +10mmHg, +1.3kPa, +1.3% to the current value, and the lower alarm limit to -10mmHg, -1.3kPa, -1.3% to the current value.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# InspCO₂ (Inspiratory CO₂) Alarm

**1** Press the InspCO₂ Alarm key.



The alarm setup menu will be displayed. Select ON/OFF of  $InspCO_2$  alarm, and set the upper alarm limit. Set the alarm condition for each measurement unit (mmHg / kPa / %). Upper alarm limit can be set in increments of 1mmHg, 0.1kPa, 0.1%. Lower alarm limit cannot be set.

InspCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after power ON or after discharge.

Key	ltem	Description
ON OFF	Individual Alarm	Selecting ON will generate the InspCO ₂ alarm.
		Selecting OFF will not generate the InspCO ₂ alarm.
🗲 Upper 🗭	Upper Alarm Limit	Sets the upper alarm limit (1 to 24mmHg, 0.1 to 3.0kPa, 0.1 to 3.0%). Setting a value equal to or above 24mmHg, 3.0kPa, 3.0% will turn the alarm OFF.
Auto	Automatic Setup	Automatically sets the upper alarm limit to 3mmHg, 0.4kPa, 0.4% to the current measurement.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).

Reference

For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

# **EtCO₂ Peak Picking Duration**

The duration to pick the EtCO₂ maximum value can be selected from 10 sec., 20 sec., or OFF.

#### **1** Press the Config. key.

CO2 Configuration	Prev.
EtCO ₂ Peak Picking Duration	10sec 20sec OFF
Unit	mmHg kPa %
O ₂ Comp.	21%
N ₂ O Comp.	OFF OFF
Anesthetic Comp.	0.0 %
Atmos, Pressure	760 mmHg

The  $CO_2$  configuration menu to select  $EtCO_2$  peak picking duration will be displayed.

#### 2 Select the duration.

10sec, 20sec will display the maximum  $EtCO_2$  value for the selected duration. If OFF is selected,  $EtCO_2$  value for each respiration will be displayed. As the  $EtCO_2$  value display is updated each second,  $EtCO_2$  value for each respiration cannot be displayed if respiration rate is above 60 Bpm.
### **Measurement Unit**

The measurement unit can be selected from mmHg, kPa, or %.

NOTE	<ul> <li>If the measurement unit is changed frequently, the continuity of the graphic trend and tabular trend may be lost.</li> <li>The alarm setup is necessary for each measurement unit. When a measurement unit is changed, make sure to set the alarm condition for the changed unit.</li> </ul>
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### Press the Config. key. 1

C <b>O</b> 2 Configuration	Prev. Disp.	
EtCO ₂ Peak		
Unit	mmHg kPa %	)
O ₂ Comp.	21 %	
N ₂ O Comp.	ON OFF	
Anesthetic Comp.	0.0 %	
Atmos, Pressure	760 mmHg	

The CO₂ configuration menu to select the measurement unit will be displayed.

2 Select the measurement unit from mmHg, kPa, %.

The graphic trend and tabular trend will be displayed with the selected measurement unit.



### O₂ Compensation

By inputting the O₂ concentration value, compensation can be made to display more accurate value.

#### Press the Config. key. 1

C <b>O</b> 2 Configuration	Prev. Disp.
EtCO ₂ Peak Picking Duration	10sec 20sec OFF
Unit	mmHg kPa %
O ₂ Comp.	21%
N ₂ O Comp.	
Anesthetic Comp.	0.0 %

The  $CO_2$  configuration menu to set the  $O_2$  compensation will be displayed.

#### **2** Set the O₂ compensation.

Pressing the key will display the tool box to set the value.

O ₂ Comp. 2 1 % ↓ (0 - 100%) Close	Use the $\bigwedge$ $\checkmark$ keys to adjust the O ₂ compensation (O ₂ concentration). The O ₂ compensation can be set in 1% increment for the value up to 30%, and 5% increment for the value above 30%. Press the Close key after setting the O ₂ compensation.
(0 - 100%)	The O ₂ compensation can be set in 1% increment for the value up to $30\%$ , and $5\%$ increment for the value above $30\%$ .

NOTE	The value cannot be changed if the total value of $O_2$ compensation and anesthetic gas compensation exceeds 100%. In such case, change the $O_2$ compensation value after changing the anesthetic gas compensation value.
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### N₂O Compensation

If  $N_2O$  is present in the respiration circuit, the  $CO_2$  value tends to be displayed higher than the actual value. By setting the  $N_2O$  compensation ON, this can be adjusted.

#### 1 Press the Config. key.



The  $\text{CO}_2$  configuration menu to set the  $N_2\text{O}$  compensation will be displayed.

**2** Select ON/OFF for " $N_2O$  Comp".

### **Anesthetic Gas Compensation**

By inputting the anesthetic gas concentration value, compensation can be made to display more accurate value.

#### **1** Press the Config. key.

C <b>O</b> 2 Configuration	Prev. Disp.
EtCO ₂ Peak Picking Duration	10sec 20sec OFF
Unit	mmHg kPa %
O ₂ Comp.	21%
N ₂ O Comp.	OFF
Anesthetic Comp.	0.0 %
Atmos, Pressure	760 mmHg

The  $CO_2$  configuration menu to set the anesthetic gas compensation will be displayed.

#### **2** Set the anesthetic gas compensation.

Pressing the key will display the tool box to set the value.



Use the **(k**eys to adjust the anesthetic gas compensation (anesthetic gas concentration), and press the **Close** key.

NOTE	The value cannot be changed if the total value of $O_2$ compensation and anesthetic gas compensation exceeds 100%. In such case, change the $O_2$
	compensation value after changing the anesthetic gas compensation value.

### **Atmospheric Pressure Compensation**

The atmospheric pressure can be adjusted to compensate for pressure differences.

#### **1** Press the Config. key.

C <b>O</b> 2 Configuration	Prev. Disp.
EtCO ₂ Peak Picking Duration	10sec 20sec OFF
Unit	mmHg kPa %
O ₂ Comp.	21%
N ₂ O Comp.	OFF
Anesthetic Comp.	0.0 %
Atmos. Pressure	760 mmHg

The  $\text{CO}_2$  configuration menu to set the atmospheric pressure will be displayed.

Pressing the key will display the tool box to set the atmospheric pressure.

#### **2** Set the atmospheric pressure.

Pressing the key will display the tool box to set the value.



Use the **1**, **1** keys to adjust the atmospheric pressure (mmHg), and press the **C**lose key.

### **ON/OFF of Parameter Display**

1 Press the Display ON/OFF key. C**O**2 Prev. Disp. Scale 50 mmHg Graphic Trend EtCO₂ Alarm ON Tabular Trend зо – 45 InspCO₂ Alarm ON The confirmation display for ON/OFF of CO₂ display will appear. З EtC**O**2 Peak Picking Duration Config. 10sec Unit Oz Comp. N20 Comp. Anesthetic Comp. Atmos. Pressure mmHg 21% OFF 0.0 % 760 ∎ Display ON/OFF ON H Cal. Airway Adıpt **2** Select Display ON or Display OFF. C**O**2 Prev. Disp. CO₂ display can be turned ON or OFF. Pressing the Display ON key will display the waveform and Display ON numeric data. Display OFF Pressing the Display OFF key will not display the waveform and numeric data. The patient is equipped with a filter line in the state "Display OFF", and it changes to the state of "Display ON" automatically that two breathing or more are detected in 30 seconds. **36.1** 37.2 The Display OFF message will be displayed inside the EtC02 parameter key. InspCODisp.OFF mmH≋ RR_IMP

If the filter line is applied to the patient and more than 2 respirations are detected in 30 seconds during "Display OFF" condition, the display will automatically return.

▲ CAUTION	<ul> <li>When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend will also be suspended.</li> <li>Whe the waveform and numeric data display is set to OFF, the respiration rate measured by CO₂ will not be displayed either.</li> </ul>
-----------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## CO₂ Concentration

# (Option Unit: MGU-722)

This menu allows to set the monitoring condition of  $CO_2$  concentration when Microstream[®]  $CO_2$  unit (MGU-722) is used.

Reference

• For parameter setup of the mainstream CO₂ unit (MGU-721), refer to P6-56 "CO₂ Concentration (MGU-721)".

C02			Prev. Disp.
Scale	50	mmHg	Graphic Trend
EtCO ₂ Alarm	ON 3 O	- 45	Tabular Trend
InspCO ₂ Alarm	ON 3		Display ON/OFF ON
Config.	EtCO2 Peak Picking Duratik Unit	on ^{10sec} mmHg	Suspend CO ₂

### CO₂ Scale

**1** Press the Scale key.



#### **2** Select the CO₂ waveform scale for displaying and recording.

For the measurement unit in mmHg, select the scale from	50	, 1	00.	
For the measurement unit in kPa and %, select the scale fr	om	4,	8,	10.

### EtCO₂ (End-Tidal CO₂) Alarm

**1** Press the EtCO₂ Alarm key.



The alarm setup menu will be displayed. Select ON/OFF of  $EtCO_2$  alarm, and set the upper and lower alarm limits. Alarm condition should be set for each unit (mmHg / kPa / %). Upper and lower alarm limits can be set in increments of 1mmHg, 0.1kPa, and 0.1%.

NOTE EtCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after power ON or after discharge.

Кеу	Item	Description		
ON OFF	Individual Alarm	Selecting ON will generate the EtCO ₂ alarm.		
	Lower Alarm Limit	SelectingOFFwill not generate the EtCO2 alarm.Sets the lower alarm limit (1 to 98mmHg, 0.1 to 13.1kPa, 0.1 to 13.1%). Setting a value equal to or below 1mmHg,		
	Upper Alarm Limit	0.1kPa, 0.1% will turn the alarm OFF. Sets the upper alarm limit (3 to 115mmHg, 0.3 to 15.0kPa, 0.3 to 15.0%). Setting a value equal to or above 115mmHg, 15.0kPa, 15.0% will turn the alarm OFF.		
Auto	Automatic Setup	Automatically sets the upper alarm limit to +10mmHg, +1.3kPa, +1.3% to the current value, and the lower alarm limit to -10mmHg, -1.3kPa, -1.3%.		

CAUTION For MGU-722, the upper EtCO₂ alarm will not generate if the upper limit is set to 100mmHg/13.4kPa and above as the measurement range is 0 to 99mmHg / 0 to 13.3kPa.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).

Reference

For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

### InspCO₂ (Inspiratory CO₂) Alarm

**1** Press the InspCO₂ Alarm key.



The alarm setup menu will be displayed. Select ON/OFF of  $InspCO_2$  alarm, and set the upper alarm limit. Set the alarm condition for each measurement unit (mmHg / kPa / %). Upper alarm limit can be set in increments of 1mmHg, 0.1kPa, 0.1%. Lower alarm limit cannot be set.

**NOTE** InspCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after power ON or after discharge.

Key	ltem	Description
ON OFF	Individual Alarm	Selecting ON will generate the InspCO ₂ alarm.
ONOFF	Inulviuual Alann	Selecting OFF will not generate the InspCO ₂ alarm.
Upper	Upper Alarm Limit	Sets the upper alarm limit (1 to 24mmHg, 0.1 to 3.0kPa, 0.1 to 3.0%). Setting a value equal to or above 24mmHg, 3.0kPa, 3.0% will turn the alarm OFF.
Auto	Automatic Setup	Automatically sets the upper alarm limit to 3mmHg, 0.4kPa, 0.4%.

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

### **EtCO₂ Peak Picking Duration**

The duration to pick the EtCO₂ maximum value can be selected from 10 sec., 20 sec., 30 sec., or OFF.

#### 1 Press the Config. key.

C <b>O</b> 2 Configuration	Prev. Disp.
EtCO ₂ Peak Picking Duration	10sec 20sec 30sec
Unit	mmHg kPa %
CO ₂ Cal.	Restart CO2 (NOLB 2 DEC)

The  $CO_2$  configuration menu to select  $EtCO_2$  peak picking duration will be displayed.

#### **2** Select the duration.

Select the duration to display the maximum  $EtCO_2$  value from 10sec, 20sec, 30sec. If OFF is selected,  $EtCO_2$  value for each respiration will be displayed. As the  $EtCO_2$  value display is updated each second,  $EtCO_2$  value for each respiration cannot be displayed if respiration rate is above 60 Bpm.

### **Measurement Unit**

The measurement unit can be selected from mmHg, kPa, or %.

NOTE	<ul> <li>If the measurement unit is changed frequently, the continuity of the graphic trend and tabular trend may be lost.</li> <li>The alarm setup is necessary for each measurement unit. When a measurement unit is changed, make sure to set the alarm condition for the changed unit.</li> </ul>
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **1** Press the Config. key.

C <b>O</b> 2 Configuration	Prev. Disp.
EtCO ₂ Peak Picking Duration	10sec 20sec 30sec
Unit	mmHg kPa %
CO ₂ Cal.	Restart CO2 (HOLD 2 SEC)

The CO₂ configuration menu to select the measurement unit will be displayed.

**2** Select the measurement unit from mmHg, kPa, %.

The graphic trend and tabular trend will be displayed with the selected measurement unit.



### **CO₂ Calibration**

CO₂ calibration can be performed using calibration gas. Calibration should be conducted every specified period or when any measurement error is found.

Perform calibration when 1 year has elapsed from the last calibration, or accumulated  $EtCO_2$  measurement time exceeds 4000 hours, or any measurement error is found.

	If the $CO_2$ gas calibration is not performed at a specified interval, $CO_2$ measurement accuracy may be affected and also subsequent gas calibration may not be possible.
--	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

#### **1** Press the Config. key.

CO ₂ Configuration	
EtCO2 Peak Picking Duration	
Unit mmHg kPa %	The $CO_2$ configuration menu with the $CO_2$ Cal. key will be displayed.
CO ₂ Cal.	

**2** Press the  $CO_2$  Cal. key to display the calibration menu.

- **3** Press the **Start Cal** key and conduct calibration according to the displayed messages.
- **4** The message, "Feed CAL. GAS" will be displayed. Press the injection button to inject the calibration gas.
- **5** The message, "Cal. Gas can be removed" will be displayed. Stop pressing the injection button to cease the injection.
- **6** The message, "CAL. OK" will be displayed. "Last Cal. Date" will be updated to the current date.

If any of the following messages is displayed, start the procedure again from step 2. "CAL. error", "CAL GAS error", "Auto Zero fail", "No stable gas flow", "CAL. failure"

o Cal Complete	CAL. OK	
	Press "Cal Comp Last Cal Date:	olete"key 2003 06/12 07:00

7 Press the Cal Complete key to end the calibration.

▲ CAUTION	<ul> <li>Perform calibration after Initialization Time (max. 180 seconds) has elapsed since the power is turned ON.</li> <li>Do not disconnect the sampling tube during calibration. Calibration will cease when the sampling tube is disconnected.</li> <li>For the following case, a message, "Calibrate the CO₂ unit (MGU-722)" or "The periodic calibration of the CO₂ unit (MGU-722) is approaching" will be displayed at power ON. Conduct CO₂ calibration.</li> <li>When the accumulated measurement time exceeds 1200 hours from first use.</li> <li>When 1 year has elapsed from the last calibration date.</li> <li>When the accumulated measurement time exceeds 4000 hours from the last calibration date.</li> <li>When EtCO₂ measurement is not stable or accuracy is degraded compared with other measuring device, conduct CO₂ calibration.</li> </ul>
-----------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Restarting the CO₂ Unit

The sampling tube will cease functioning when erroneous condition such as blocking of exhaust tube, sampling tube or nasal prong is detected. When the pump ceases functioning, "Check  $CO_2$  unit" message will be displayed. After resolving the problem, press the Restart  $CO_2$  key and restart the measurement.

#### **1** Press the Config. key.

_Disp
10sec 20sec 30sec
mmHg KPa %
Restart CO2 (HOLB 2 SEC)

The  $CO_2$  configuration menu with the Restart  $CO_2$  key will be displayed.

**2** Press the Restart  $CO_2$  key for 2 seconds.

If the "Check  $CO_2$  unit" message is not displayed, the Restart  $CO_2$  key will not function.

#### **3** Check that the unit is restarted.

The sampling pump will start to function, and the "Check  $CO_2$  unit" message will disappear. Check that the message has disappeared and the measurement data is displayed.

NOTE	If the "Check Sample Line", "Check $CO_2$ Exhaust Port", "Check $CO_2$ unit", $CO_2$ Unit Error" message does not disappear after restarting the unit, equipment failure can be considered, or parts replacement of $CO_2$ unit may be necessary. Contact our service representative.
------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **ON/OFF of Parameter Display**

**1** Press the Display ON/OFF key.



**2** Select Display ON or Display OFF.





The Display OFF message will be displayed inside the parameter key.

If the filter line is applied to the patient and more than 2 respirations are detected in 30 seconds during "Display OFF" condition, the display will automatically return.

<b>▲</b> CAUTION	<ul> <li>When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend will also be suspended.</li> <li>Whe the waveform and numeric data display is set to OFF, the respiration rate measured by CO₂ will not be displayed either.</li> </ul>
------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Suspending CO₂ Measurement

**1** Press the Suspend  $CO_2$  key.



**2** To resume the  $CO_2$  monitoring, press the Resume  $CO_2$  key.



Pressing the Resume CO₂ key will display the waveform and numeric data.

### Ventilator

This menu allows to set the ventilator monitoring condition.



AWP Scale : Sets the scale for AWP (airway pressure) waveform. AWF Scale : Sets the scale for AWF (airway flow) waveform.

CAUTION For PURITAN-BENNETT Ventilator, AWP and AWF waveform cannot be displayed or recorded. Only the numeric data will be displayed.

### **AWP Scale**

**1** Select the scale.



### **AWF Scale**

**1** Select the scale.



Select the appropriate scale from 5, 10, 20, 50, or 180.

### **Tb (Blood Temperature)**

When thermodilution catheter is used to measure the cardiac output, the blood temperature can be monitored.



Tb Alarm : Sets ON/OFF of blood temperature alarm and upper and lower alarm limits.

#### **Tb Alarm**

**1** Select ON/OFF of blood temperature alarm and set upper and lower alarm limits.



The alarm setup should be performed for each measurement unit (°C / °F). The upper and lower limit can be set in increments of 0.5°C / 1.0°F.

Key	Item	Description	
ON OFF	Individual Alarm	Selecting ON will generate the TEMP alarm.	
		Selecting OFF will not generate the TEMP alarm.	
	Lower Alarm Limit	Sets the lower alarm limit (30.0 to 44.0°C / 86.0 to	
🗲 Lower 🔿		111.0°F). Setting a value 30.0°C / 86.0°F or below will	
		turn the alarm OFF.	
	Upper Alarm Limit	Sets the upper alarm limit (31.0 to 45.0°C / 88.0 to	
🗲 Upper 🔿		113.0°F). Setting a value 45.0°C / 113.0°F or above will	
		turn the alarm OFF.	
		Automatically sets the upper limit to +2.0°C / +4.0°F to	
Auto	Automatic Setup	the current value, and lower limit to $-2^{\circ}C$ / $-4.0^{\circ}F$ to the	
<u> </u>		current value.	

To maintain the alarm setting even after the power is turned OFF or after the discharge procedure, store the setting to one of the alarm modes, or select Backup for "Alarm" on the "Backup at Discharge" menu (Monitor Setup).



For the alarm mode setup procedure, refer to "8. System Configuration Alarm Mode".

### **Vigilance Data**

When Vigilance, Vigilance CEDV, VigilanceII, or Vigileo (oximeter/CCO measurement device manufactured by Baxter) is used, Vigilance data display can be selected from several modes.

#### [ICO Mode]

Vigilance				Prev. Disp.
STAT Mo	ode [	ON	OFF	
Index Di	sp	ON	OFF	
CO_AVG		CI_AVG		
ICO[1]		ICI[1]		Graphic Trend
IC0[2]		IC1[2]		Vigilance List
IC0[3]		IC1[3]		List
IC0[4]		IC1[4]		
IC0[5]		IC1[5]		
IC0[6]		IC1[6]		

STAT Mode : When Vigilance is in CCO mode, STAT mode display can be set ON or OFF. Index Disp. : When Vigilance is in CCO mode, Index display can be set ON or OFF.

When the Vigilance is in ICO mode, the 6 latest data of ICO (Intermittent Cardiac Output) and ICI (Intermittent Cardiac Index) will be displayed.

### **STAT Mode / Index Display**

**1** ON / OFF of STAT Mode and ON / OFF of Index display can be selected on the Vigilance display.

#### [STAT Mode OFF, Index Display OFF]

\$v02	83	
000	5.0	L/min
EDV	160	mL
RI	38.5	°C

 $SvO_2$  (or  $ScvO_2$ ), CCO, EDV, BT data will be displayed inside the  $SvO_2$ +CO numeric data box.

#### [STAT Mode OFF, Index Display ON]

Sv02	83	
CCI		L/min/m2
EDAI		mL/m²
BI	38.5	°C

By setting the Index display ON, CCI, EDVI data will be displayed instead of CCO and EDV.

#### [STAT Mode ON, Index Display OFF]

\$v02	83	×	
Sv02 CCO_STAT		L/min	
EDV_STAT	_1_6 0	mL	
BT	38.5	°C	

By setting the STAT mode ON, CCO_STAT, EDV_STAT data will be displayed instead of CCO and EDV.

#### [STAT Mode ON, Index Display ON]

n/m2
2

By setting the STAT mode and Index display ON, CCI_STAT, EDVI_STAT will be displayed instead of CCO and EDV.

### Stopwatch

By setting a stopwatch key on the home display, a stopwatch function can be used.

Stopwatch	Prev. Disp.
TIMER1	00:00:00 START STOP RESET
TIMER2	00:00:00 START STOP RESET
LA	BEL1 LABEL2

TIMER1, 2 : Starts/stops the stopwatch function. LABEL1, 2 : A label can be set for each timer.

### Label Setup

**1** Press the LABEL key.



The stopwatch label setup menu will be displayed. Enter 8 characters using alphanumeric keypad.

### **START/STOP** of Stopwatch



	• If discharge procedure is performed during stopwatch operation, the counting will
ΝΟΤΕ	stop and the stopwatch time will be reset to "00:00:00".
	• The stopwatch will continue counting even when the monitoring is suspended.

Blank Page

# Chapter 7

# **Function**

This chapter describes the function such as arrhythmia analysis, trend, and recall.

Arrhythmia Analysis Definition, etc.	• 7-2
Arrhythmia Definition	· 7-2
QRS Classification	
●Arrhythmia Type ·····	
To Set the Arrhythmia Alarm	• 7-3
To Perform Arrhythmia Learning	
Graphic Trend Display/Print	• 7-6
To Display the Graphic Trend	· 7-6
Graphic Trend Group Setup	7-10
The Description of the Display	
Tabular Trend         Display/Print	7-12
To Display the Tabular Trend	7-12
The Description of the Display	7-13
Parameter Setup for Tabular Trend	7-13
Recall Display/Print	7-15
To Display the Recall Menu	7-15
Recall List Display	
To Display and Record the Enlarged Recall Wavefor	orm
To Set the Recall Condition	
NIBP List Display/Print	7-19
To Display the NIBP List	7-19
The Description of the Display	
ST Display ST Display, Alarm Setup, etc	7-21
To Display the ST Measurement Menu	7-21
To Set the Reference Waveform	
●ST Alarm Setup ······	
OCRG Display	7-24
CO Measurement Measurement/Editing	7-25
To Display the Cardiac Output Menu	7-25
Message List	7-26
Cardiac Output Setup	7-26
To Measure the Cardiac Output	7-28
To Edit the Cardiac Output Data	7-28
Hemodynamics Calculation/Print	7-30
Calculation Data	
To Display the Hemodynamic Data	7-31
To Calculate the Newly Input Hemodynamic Data .	7-31
To Edit the Hemodynamic Data	7-33
Ventilator P-V, F-V Loop Display	7-34
P-V Loop	7-34
F-V Loop	7-35
Displaying the Ventilator Measurement	7-36

Respiration List Display/Print	7_37
To Display the Respiration List	
The Description of the Display	
Respiration List Setup	7-38
ST Graphic Trend Display/Print To Display the ST Graphic Trend	
ST Graphic Trend Group Setup	
The Description of the Display	7 41
ST Tabular Trend Display/Print	
To Display the ST Tabular Trend	7-42
The Description of the Display	
Other Bed Display/Alarm	7-44
Other Bed Display	7-44
Other Bed Alarm Setup	
Full Disclosure Waveform Recording	7-46
To Record the Full Disclosure Waveform Data	
To Display the Full Disclosure Waveform	
About the Waveform Display Duration	
To Shift the Displayed Waveform	
To Change the ECG Waveform Size	7-50
ECG Waveform Amplitude	7-50
To Search by Time	7-50
To Search by Alarm	7-51
To Search by Admit/Discharge Date	7-52
To Print the Waveform	
To Enlarge the Waveform	7-54
To Print the Waveform (Enlarged Display)	7-55
Vigilance/Vigileo List Display/Print	
To Display the Vigilance/Vigileo List	
The Description of the Display	7-57
Vigilance/Vigileo List Setup	7-57

### **Arrhythmia Analysis**

This section explains the arrhythmia analysis, alarm setup procedure, etc.

### **Arrhythmia Definition**

The arrhythmia detection is performed by learning the normal waveform of the patient and by determining VPC by comparing the waveform (QRS pattern) and R-R interval for each heartbeat. A pattern matching is performed with each VPC detected from R-R interval, QRS amplitude, QRS area, QRS polarity, etc., and determines as VPC after discriminating the noise from VPC.



#### QRS Classification

The QRS analysis is performed by comparing with the learned waveform and QRS pattern matching.

N (Normal)	Normal QRS beat	
V (VPC)	Ventricular extrasystole	
? (Undetermined beat)	Learning arrhythmia, or unmatched beat	
P (Pacing beat) Pacing beat		
F (Fusion beat)	Fusion beat of pacing and spontaneous beat	
S (SVPC) Supraventricular extrasystole		

### Arrhythmia Type

The alarm is generated according to the arrhythmia classification of a normal QRS and VPC.

Туре	Meaning	Detection Criteria
ASYSTOLE	Cardiac Arrest	Cardiac arrest is detected for more than the preprogrammed time.
VF	Ventricular Fibrillation	A random, rapid electrical activity of the heart is detected.
VT	Ventricular Tachycardia	HR is same or above the preprogrammed value (140bpm or 120bpm), and 9 or more continuous ventricular beats are detected.
SLOW_VT		9 or more continuous ventricular beats are detected. (HR: below 140bpm / 120bpm)
TACHY	Tachycardia	HR is over the upper alarm limit.
BRADY	Bradycardia	HR is below the lower alarm limit.
RUN	Consecutive VPC	HR is same or above the preprogrammed value, and continuous VPC exceeding the preprogrammed value (2 to 8beats) is detected.
COUPLET	Couplet Ventricular Extrasystole	2 continuous VPC beats are detected.
PAUSE		Cardiac arrest is detected for more than the preprogrammed time.
BIGEMINY	Ventricular Bigeminy	3 or more continuous QRS pattern of V-N is detected.
TRIGEMINY	Ventricular Trigeminy	3 or more continuous QRS pattern of V-N-N is detected.
FREQUENT	Frequent VPC	VPC exceeding the preprogrammed value is detected within 1 minute.



For details of Arrhythmia Alarm Setup, refer to " 4. Monitoring Setup Alarm Setup ● To Set the Arrhythmia Detection Level".

▲WARNING	Objective and constant arrhythmia detection is possible through the fixed algorithm incorporated in this monitor. However, excessive waveform morphology change, motion artifact, or the inability to determine the waveform pattern may cause an error, or fail to madequate detection. Therefore, physicians should make final decisions using manual recording, alarm recording and recall waveform for evaluation.	
▲ CAUTION	For proper arrhythmia detection and ECG monitoring, verify proper electrode placement, lead selection, and ECG waveform size. If necessary, turn ON the AC filter. Improper electrode placement, lead selection, and ECG waveform size can cause errors in detection.	

# To Set the Arrhythmia Alarm

ON/OFF of arrhythmia alarm and reference of arrhythmia analysis can be set.

<b>1</b> Press the Menu $\rightarrow$	larm $\rightarrow$ Arrhy. keys.
Arrtvrthnia alarn 1/3     Page Down       Asystole     5     ON     OFF        VF     ON     OFF        VT     ON     OFF        Slow VT     ON     OFF        Tachy     ON     OFF        Brady     ON     OFF	Cal Cal Cal Cal Cal Cal Cal Cal
Arrtvithila alarin 2/3       Page lip       Page Down         Run       3       ON       OFF       Image Down         Bigeminy       ON       OFF       Image Down       Image Down         Trigeminy       ON       OFF       Image Down       Image Down         Pause       3.0       ON       Image DoFF       Image Down         Couplet       ON       Image DoFF       Image Down         Frequent       10       ON       Image DoFF       Image Down	Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call Call
Arrhvthnia alarn 3/3 Page llp HR Low Limit for VT 120 140 HR Low Limit for RUN 40	On page 3/3, the analysis condition (HR Low Limit) setup menu for VT and RUN will be displayed.

Page	Arrhythmia
Page 1/3	ASYSTOLE, VF, VT, SLOW_VT, TACHY, BRADY
Page 2/3	RUN, COUPLET, BIGEMINY, TRIGEMINY, PAUSE, FREQUENT
Page 3/3	HR Low Limit for VT, HR Low Limit for RUN

#### **2** Set the reference range.



#### <Arrhythmia Reference Range>

Arrhythmia	Reference Range	Default
ASYSTOLE	3 to 10 sec.	5 sec.
RUN	2 to 8 beats	3 beats
PAUSE	1.5 to 5 sec.	3 sec.
FREQUENT	1 to 50 beats/min.	10 beats/min.

On the 3rd page, HR low limit to perform the arrhythmia analysis for VT and RUN can be set.

Arrhythmia	HR Low Limit	Default
VT	120bpm, 140bpm	120bpm
RUN	0 to 100bpm	40bpm

Reference

For details of Arrhythmia Alarm Setup, refer to " 4.Monitoring Setup Alarm Setup ●To Set the Arrhythmia Detection Level".

#### **3** Select ON or OFF for the alarm.

Asystole [ Aları	5 ON OFF Recall m will generate. Alarm will not generate.	
Reference	For details of Arrhythmia Alarm Setup, refer to " 4. Monitoring Setup ON/OFF of Arrhythmia Alarm".	Alarm Setup ●To Set

#### **4** Select ON or OFF for recall factor.

ON/OFF of recall factor can be set on the alarm setup menu.



Pressing the Recall key will switch the ON/OFF selection.

### **To Perform Arrhythmia Learning**

Learning of normal ECG largely affects the accuracy of arrhythmia analysis.

If any error occurs in arrhythmia detection and QRS judgment, performing arrhythmia learning will recover the original analyzing accuracy. Arrhythmia learning will be performed for about 20 beats for the normal ECG, but it may take longer if the heartbeat is unstable. During arrhythmia learning, arrhythmia alarm other than ASYSTOLE, TACHY, BRADY will not be

During arrhythmia learning, arrhythmia alarm other than ASYSTOLE, TACHY, BRADY will not be generated.



**2** Start arrhythmia learning.

Arrhy, Learn Learning arrhythmia

Pressing the key while learning arrhythmia will not stop the learning.

#### **3** During arrhythmia learning, a message will be displayed.



NOTE If Used is selected for "Pacemaker" on the Admit menu, LED on the Arrhy. Learn key will not light and the "LEARN" message will not be displayed although the arrhythmia learn procedure is performed.

This section explains the graphic trend function and recording procedure.

### To Display the Graphic Trend

The graphic trend menu can be accessed from the menu, or from the preprogrammed user key. If the data is displayed on the home display, the 24 hours of graphic trend data in 1-minute interval will be automatically stored and displayed.

- 1 Press the Menu  $\rightarrow$  Function  $\rightarrow$  Graphic Trend keys. Graphic Trend Prev. Disp. 4H 🗲 Time Cursor 200/19 10:15 HR 60 bpm Scale 0 150 The graphic trend menu will be displayed. BP1 Scale 116/77 (92) mmHg The display will switch by pressing the A, B, C, D keys. 0 50 BP2 Scale 23/ 10 (15) nnHg 9:45 Group A B C D Group Setup Print
- **2** The Group D display will allow selecting the parameter for numeric data on the trend menu.



Parameter	Description
HR	Heart Rate
VPC	VPC beats
EVENT1	ASYSTOLE, VF, VT, SLOW_VT, RUN, BIGEMINY
EVENT2	TRIGEMINY, PAUSE, COUPLET, TACHY, BRADY, FREQUENT
SpO ₂	SpO ₂ value
PR_SpO ₂	SpO ₂ pulse rate
ST(I), ST(II), ST(III), ST(aVR), ST(aVL),           ST(aVF), ST(V1), ST(V2), ST(V3),           ST(V4), ST(V5), ST(V6)	ST level
NIBP	Noninvasive Blood Pressure (Systolic/Mean/Diastolic)*
TEMP1,2, TEMP3	Temperature
Tb	Blood Temperature (Cardiac Output Measurement)
BP1, BP2, BP3, BP4, BP5	Blood Pressure (Systolic/Mean/Diastolic)
PR_IBP	Blood Pressure Pulse Rate (BP1 or ART)
PDP	Peak Diastolic Pressure of IABP
CPP	Cerebral Perfusion Pressure
SvO ₂	Mixed Venous Oxygen Saturation
ScvO ₂	Central Venous Oxygen Saturation
CCO	Continuous Cardiac Output
CCI	Continuous Cardiac Index
BT	Blood Temperature
RR_IMP	Impedance Respiration Rate
APNEA	Apnea Time (Impedance, CO ₂ , ventilator)
CO ₂	EtCO ₂ / InspCO ₂
RR_CO ₂	CO ₂ Respiration Rate
RR_VENT	Ventilator Respiration Rate
BIS	BIS Monitor Data

CAUTION If the mean BP display is set to OFF on the NIBP configuration setup, the mean BP will not be displayed for the tabular trend or the NIBP list.

ΝΟΤΕ

The apnea time will be stored when it exceeds the alarm threshold level. If lower than the alarm threshold level, it will be stored as "0 (zero)".

#### **3** Select the scale for display.



Pressing the Scale key will switch the scale according to the displayed parameter as shown below.

Parameter	Scale	Unit
HR	100, 200, 300	bpm
VPC	20, 50, 100	beat
EVENT	none	
SpO ₂	0–100, 50–100, 80–100	%
PR_SpO ₂	100, 200, 300	bpm
ST	±0.2, ±0.5, ±1.0, ±2.0	mV
	±2, ±5, ±10, ±20	mm
NIBP	100, 150, 200, 300	mmHg
	16, 20, 24, 40	kPa
	20–45, 30–40	°C
TEMP1,2, TEMP3	68–113, 86–104	°F
TL	20–45, 30–40	°C
ТЬ	68–113, 86–104	°F
	20, 50, 100, 150, 200, 300	mmHg
BP1, BP2, BP3, BP4, BP5	4, 8, 16, 20, 24, 40	kPa
PR_IBP	100, 200, 300	bpm
PDP	20, 50, 100, 150, 200, 300	mmHg
FDF	4, 8, 16, 20, 24, 40	kPa
СРР	20, 50, 100, 150, 200, 300	mmHg
GFF	4, 8, 16, 20, 24, 40	kPa
SvO ₂	0–100, 50–100, 80–100	%
ScvO ₂	0–100, 50–100, 80–100	%
CCO	6.0, 12.0, 20.0	L/min
CCI	6.0, 12.0, 20.0	L/min/m ²
BT	20-45, 30-40	°C
RR_IMP	50, 100, 150	bpm
APNEA	15, 30	
CO ₂	50, 100	mmHg
	4.0, 8.0, 10.0	kPa
RR_CO ₂	50, 100, 150	bpm
RR_VENT	50, 100, 150	bpm
BIS	0–100 fixed	(no unit)

**4** Select the display time range.



Time Range 1H 2H 4H 8H 12H 24H Close

Pressing the display time range key will display the time range selection tool.

Select the time range for the graphic trend.

Time Range	Resolution
1 hour	1 min.
2 hour	1 min.
4 hour	1 min.
8 hour	2 min.
12 hour	3 min.
24 hour	6 min.

#### 5 Select the time span.



Cursor 06/04 10:06

Scrolls to present

Scrolls to past Scrolls to present

Scrolls the graphic trend display to past or present data with the selected time range.

Pressing the  $\leftarrow$  key will scroll the display to the past data with the selected time range.

Pressing the  $\rightarrow$  key will scroll to the present data with the selected time range.

The data of selected time can be displayed by moving the cursor.

Pressing the  $\leftarrow$  key will scroll to the past data Pressing the  $\rightarrow$  key will scroll to the present data.

### 7 Enlarge the display.

Scrolls to past

6 Move the cursor.



#### 8 Print the graphic trend data.

	Print The displayed graphic trend data will be printed.
NOTE	<ul> <li>The following graphic trend data cannot be printed on the central monitor recorder.</li> <li>BP trend, NIBP trend, PDP trend, CPP trend when the BP measurement unit is kPa.</li> <li>TEMP trend, Tb trend when the temperature measurement unit is °F.</li> </ul>

### **Graphic Trend Group Setup**

Each trend group displays combination of 3 parameters simultaneously. 3 types of trend group can be programmed.



**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Graphic Trend  $\rightarrow$  Group Setup keys.



The graphic trend group setup menu will be displayed. On this menu, parameters to display for Group A, Group B, and Group C can be selected.

#### **2** Select the trend group to set the parameters.

Graphic T	oond Grou	o Sotup	[	Prev. Disp.
	Group A HR BP1 BP2		Group B HR RR_IMP APNEA	_DTSP
		-Tr	end Items	
HR	UPC	EVENT1	EUENT2 Sp02 PR_Sp02	Ł
ST(I)	ST (II)	ST (III)	ST (aUR) ST (aUL) ST (aUF)	
ST (U)	ST (U2)	ST (U3)	ST (U4) ST (U5) ST (U6)	C
NIBP	TEMP1,2	TEMP 3	ТЬ	
BP1	BP2	BP3		
PR_IBP	PDP	CPP		
Sv02	Scv02	CCO	CCI BT BIS	
RR_IMP	APNEA	CO2	RR_CO2 RR_UENT	

Select the trend group by pressing the Group A, Group B, or Group C key.



Select the parameter by pressing the displayed parameter keys. Pressing the parameter key will sequentially set the 3 (three) parameters from the top.

## The Description of the Display



The measured data will be compressed for the 8-hour / 12-hour / 24-hour display.

Parameter	Compressed Form
HR	Mean Value
VPC	Maximum Value
EVENT	Logical Sum
SpO ₂	Mean Value
PR_SpO ₂	Mean Value
ST	Mean Value
NIBP	Current Value
TEMP	Mean Value
BP	Mean Value
PR_IBP	Mean Value
PDP	Mean Value
CPP	Mean Value
SvO ₂	Mean Value
ScvO ₂	Mean Value
CCO	Mean Value
CCI	Mean Value
BT	Mean Value
RR_IMP	Mean Value
APNEA	Maximum Value
CO ₂	Mean Value
RR_CO ₂	Mean Value
RR_VENT	Mean Value
BIS	Mean Value

### **Tabular Trend**

This section explains the tabular trend function and recording procedure.

### To Display the Tabular Trend

The tabular trend menu can be accessed from the menu, or from the preprogrammed user key. The 24 hours of data in 1-minute interval will be automatically stored and displayed if the data is displayed on the home display.

Tekular Tekular 09/19 HR UPC/n in \$1(10) mn 8(P1.5 milg 	60 0.5 0.2 116 77 92 23 10 15 92 30 33 30 30 10 36.1	9:10 60 0.5 0.2 23 116 77 92 23 10 15 92 30 10 35.1 37.2	Inter∨ 9:20 60 0.5 0.2 116 77 92 23 10 15 92 30 33 30 10 36.1 37.2	9:30 60 0.5 116 777 92 23 10 15 92 30 30 10 36.1 37.2	9:40 60 0.5 0.2 116 77 92 23 10 15 92 30 33 30 10 36.1 37.2	0 0.5 0.2 0 116 1 77 92 92 92 92 92 93 10 15 92 92 93 30 33 30 3 30	0         60           0         0           5         0.5           2         0.2           6         116           7         77           2         92           3         23           0         10           5         15           5         15           2         92           0         30           3         33           0         30           1         36.1	The tab	ular trer	nd of 17 j	paramet	ers wi	ll be dis	play
--------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------	-----------	------------	---------	--------	-----------	------

		_
Interval		
<b>1</b> min	<b>5</b> min	10min
15min	30min	60min
	Close	

Interv.60M

Pressing the time interval key will display the time interval selection tool.

Select the time interval for the tabular trend display.

5M will display the data in real time such as 10:00, 10:05, 10:25. Selecting Selecting 60M will display the data in real time such as 10:00, 11:00, 12:00. If the list is displayed at 10:35, the data from 10:00 will be displayed.

#### 3 Shift the page.



The page can be shifted past or present by page with the displayed time interval.

Shift to past Shift to present

Pressing the **F** key will shift one page to the past data. Pressing the  $\rightarrow$  key will shift one page to the present data.

The data will be listed in 8 columns.

If 5-minute time range is selected and the starting time on the list is 10:00, 35 minutes from 10:00 to 9:25 will be listed in 1 page.

Pressing the **F** key will display the list from 9:20 to 8:45.

#### Shift the displayed column. 4

The tabular trend data can be shifted. Pressing the *key* will shift the display to past. Pressing the  $\rightarrow$  key will shift the display to present. Shift to past Shift to present

#### 5 Print the list data.

	Print The displayed tabular trend data will be printed.
NOTE	<ul> <li>The following tabular trend data cannot be recorded on the central monitor recorder.</li> <li>BP trend, NIBP trend, PDP trend, CPP trend, PCWP trend when the BP measurement unit is kPa.</li> <li>TEMP trend, Tb trend when the temperature measurement unit is °F.</li> </ul>

### The Description of the Display

Latest Date

Tabular Trend	Shift	⇒	Interv	~10M	•	Page	➡	Prev. Disp.	Latest Time
09/19	9:00	9:10	9:20	9:30	9:40	9:50	10:00	10:1 <del>0</del>	
HR	60	60	60	60	60	60	60	60	
UPC/min	0	0	0	0	0	0	0	0	
ST(I) mm	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5	
ST(II) mm	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2	
BP1_S_mmHg	116	116	116	116	116	116	116	116	
_D nnHg	77	77	77	77	77	77	77	77	
_M mmHg	92	92	92	92	92	92	92	92	
BP2_S mnHg	23	23	23	23	23	23	23	23	
_D nnHg	10	10	10	10	10	10	10	10	
_M nnHg	15	15	15	15	15	15	15	15	
Sp02	92	92	92	92	92	92	92	92	
RR_IMP	30	30	30	30	30	30	30	30	
EtCO ₂ mmHg	33	33	33	33	33	33	33	33	
RR_C02	30	30	30	30	30	30	30	30	
APNEA	10	10	10	10	10	10	10	10	
T1 ℃	36.1	36.1	36.1	36.1	36.1	36.1	36.1	36.1	
T2 C	37.2	37.2	37.2	37.2	37.2	37.2	37.2	37.2	
					[	List setup	] [	^{&gt;} rint	

For the data when the measurement was not performed (before admittance) or when the monitoring was suspended, the time will be displayed as "--:--". Also, if the data is not displayed on the home display, or the BP is not zero balanced, the data will be displayed as "--".

### Parameter Setup for Tabular Trend

#### **1** Press the List Setup key on the tabular trend menu to display the tabular trend setup menu.



List setup	_		Page [	own Prev.	
HR HR	1	OFF		Disp.	
⊣ UPC/min		HR	UPC/min	UPC/hour	
⊐ stan	1	ST(I)	ST (II)	STOLD	
⊐ BP1_S	l	ST (aUR)	ST (aUL)	ST (aUF)	
<u> </u>					
<u>۳</u>	]	ST (U)	ST (U2)	ST (Ų3)	
⊐ BP2_S		ST (U4)	ST (U6)	ST (U6)	
)		NIBP_S	NIBP_D	NIBP_M	
<u>н</u>	1	HIDP_3	MIDP_D	HIBP_II	
→ Sp02	i	BP1_S	BP1_D	BP1_M	
RR_IMP	1	BP2_S	BP2_D	BP2_M	
- EtCO2	1				
- RR_C02	i	BP3_S	BP3_D	BP3_M	
- APNEA	j –				
- 11	]				
- T2	]				

**2** Select the position on the list.

List setup HR		OFF	Page (	own Prev. Disp.
UPC/min		HR	UPC/min	UPC/hour
ST (II)		ST(I)	ST (II)	ST (III)
BP1_S		ST (aUR)	ST (aUL)	ST (aUF)
		ST (Ų)	ST (U2)	ST (Ų3)
BP2_S	4	ST (U4)	ST (U6)	ST (U6)
	<b>•</b>	NIBP_S	NIBP_D	NIBP_M
- Sp02		BP1_S	BP1_D	BP1_M
RR_IMP		BP2_S	BP2_D	BP2_M
- EtC02 - RR_C02		BP3_S	BP3_D	BP3_M
APNEA				
	ļ			

Select the position.

There are 17 positions on the list to set the parameter.

#### **3** Select the parameter for display.

List setup		Page	Jown Pro Disp
<u>HR</u>	OFF		
UPC/min	HR	UPC/min	UPC/hour
	ST(I)	ST (II)	ST CIID
BP1_S	ST (aUR)	ST (aUL)	ST (aUF)
	ST (U)	ST (U2)	ST (Ų3)
BP2_S	ST (U4)	ST (Us)	ST (U6)
	NIBP_S	NIBP_D	NIBP_M
	BP1_S	BP1_D	BP1_M
- RR_IMP	BP2_S	BP2_D	BP2_M
- EtC02 RR_C02	BP3_S	BP3_D	BP3_M
- APNEA			
<u> </u>			
		<b>-</b> 1	

Select the parameter to display for the previously selected position. The position will automatically shift downward to allow consecutive parameter selection.

**NOTE** The apnea time will be stored when it exceeds the alarm threshold level. If lower than the alarm threshold level, it will be stored as "0 (zero)".

This section explains the recall menu function and recording procedure.

### To Display the Recall Menu

The recall menu can be accessed from the menu, or from the preprogrammed user key.

Recall	1/20	OEvent	s 🔳	▼	Prev. Disp.
	Asystole				
16:21	UT	Slow UT			
L					
Disp Selec	lay tion_				Recall Setup

When the assigned alarm factor occurs, the waveform (12 seconds) and numeric data at alarm occurence will be stored for up to 200 data.

The recall data to be displayed can be selected on the display selection menu.

On the recall list display, 5 compressed recall waveform will be displayed. Pressing one of the compressed recall waveform will enlarge the waveform.



Time at Alarm Occurrence Recall Factor Recall Waveform (Compressed: 7se

If the recall data exceeds 200, the data will be erased from the oldest one.

#### Recall List Display

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Recall keys.



The recall list will be displayed. The compressed waveform of about 7 seconds will be displayed.

The alarm occurrence time, the recall factor occurred at the same time, and the compressed waveform of recall waveform 1 will be displayed.

#### **2** Select the recall factor to display on the recall list.



Select the numeric data, arrhythmia to display as recall factor.

HR --If the key LED is lighted, recall data will be displayed. If the key LED is extinguished, recall data will not be displayed.



### **3** Shift the recall list display.



## To Display and Record the Enlarged Recall Waveform

On the recall list display, pressing one of the recall factor will display the enlarged recall waveform. On the enlarged recall waveform display, the recall waveform will be displayed in 25mm/s and by using the cursor, the data before and after the alarm occurrence can be checked.



Press one of the recall factors on the recall list.

1 Pressing one of the recall factors will display the enlarged recall waveform.



**2** Shift the waveform left or right.



The recall waveform display can be shifted to left or right. ← key will shift to the older data. ➔ key will shift to the newer data.

**3** The alarm factor occurred at the same time will be displayed.



### To Set the Recall Condition

On the recall menu, the storing condition at alarm occurrence can be set. The recall waveform and recall factor (numeric data, arrhythmia) can be selected. The displayed parameters will differ depending on the used option unit.

RR_CO:

Recall Setup	Prev Disp	
Wave	12sec.,Wave1 OFF ,Wave2 EC62	
Numeric	HR , ST , NIBP , BP1 , BP2 , RR APNEA , Sp0 $_2$ , T1 , T2 , C0 $_2$	
Arrhy.	Asystole , VF , VT , Slow VT , Run Bigeminy , Trigeminy , Pause Couplet , Tachy , Brady VPC Frequent	

NOTE	By setting "Store all alarms to Recall" (Monitor Setup) to ON, all alarms (numeric and arrhythmia alarm) will be stored regardless of this Recall Setup. (Default: OFF) At this time, Numeric, Arrhy. keys will not be displayed on the "Recall Setup" menu.
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



For details, refer to "8. System Configuration Monitor Setup Store all alarms to "Recall"".

#### 1 Select the recall waveform.

	Wave
Recall Wave	Numeric Arrhy. Prev. Disp.
Wave Tin	ne 12second
Wave1	ECG1     ECG2     BP1     BP2       Sp02     RESP     C02
Wave2	ECG1         ECG2         BP1         BP2           Sp0z         RESP         CO2           orF

-

Pressing the Wave key will display the menu to select the recall waveform.

Up to 2 waveforms can be selected for recall waveform. Select the recall waveform from Wave 1 and Wave 2. The key with the LED lighted is the selected waveform.

#### **2** Select the recall factor (numeric data).

Numeric	Pressing the Numeric key will display the menu to select the numeric data recall factor.
Recall Factor (numeric) Wave Arrhy. Prov. HR ST M NIBP BP1 BP2 M RR APNEA Sp02 T1 M T2 M CO2	Select the recall factor by pressing the keys. The key with the LED lighted will be the recall factor. The alarm OFF mark will be displayed inside the parameter key if the alarm is set to OFF for that parameter.

#### **3** Select the recall factor (arrhythmia).

Arrhy.

Pressing the Arrhy. key will display the menu to select the arrhythmia alarm factor.

Recall Factor Carrhy	rthmia)	Wave	Numeric	Prev. Disp.
Asystole	VF		VT	
Slow VT	<b>R</b> un		Bigeminy	X
Trigeminy 💥	Pause	X	Couplet	X
Tachy	Brady		Frequent	X
			Arrhy. A	Alarm

Select an arrhythmia for recall factor. The key with LED lighted will be the recall factor.

The alarm OFF mark will be displayed inside the arrhythmia key if the alarm is set to OFF for that arrhythmia.

	The recall waveform will start with the following delay time tracing back from the alarm occurrence.						
NOTE		Adult	Child	Neonate Meas. Data Alarm Arrhy. Alarm			
	Delay Time	12 sec.	12 sec.	8 sec.	12sec.		

### **NIBP List**

This section explains the NIBP list function and recording procedure.

### To Display the NIBP List

The NIBP list display can be accessed from the menu, or from the preprogrammed user key.



On the NIBP list, NIBP data and HR,  $SpO_2$  pulse rate,  $SpO_2$  value at the commencement of NIBP measurement will be stored and displayed for 120 NIBP measurements. If the data exceeds 120, the data will be erased from the oldest one.

1 Press the | Menu |  $\rightarrow$  | Function | NIBP List keys. NIBP List  $\mathbf{\overline{}}$ ▼ Prev. Disp. The NIBP list will be displayed. Date Time 09/19 10:15 09/19 10:10 
 HR
 PR_Sp02
 Sp02

 76
 76
 96

 76
 76
 96
 No. 1 NIBP mnHg 128/ 89 120/ 85 2 09/19 10:05 129/ 90 76 76 96 T 76 76 78 96 96 95 09/19 10:00 09/19 9:55 129/ 90 129/ 90 76 76 Displays the newest 12 data of the NIBP list. 6 09/19 9:50 128/ 91 78 96 98 96 7 09/19 9:45 129/ 90 76 76 Shifts the display to newer data by 1 page (12 data). 09/19 9:40 09/19 9:35 129/ 90 129/ 90 76 76 76 10 09/19 9:30 11 132/ 93 76 76 96 Shifts the display to older data by 1 page (12 data). ▼ Print All Print

#### 2 Print the NIBP list.



All the data stored on NIBP list will be printed on the built-in recorder.





Currently displayed NIBP list will be printed.

# The Description of the Display

BP Lis	it			▼	Pre Dis
No.	Date Time	NIBP mmHg	HR	PR_Sp02	Sp02
1	09/19 10:15	128/ 89	76	76	96
2	09/19 10:10	120/ 85	76	76	96
3	09/19 10:05	129/ 90	76	76	96
4	09/19 10:00	129/ 90	76	76	96
5	09/19 9:55	129/ 90	76	76	96
6	09/19 9:50	128/ 91	78	78	95
7	09/19 9:45	129/ 90	76	76	96
8	09/19 9:40	129/ 90	76	76	98
9	09/19 9:35	129/ 90	76	76	96
10	09/19 9:30	132/ 93	76	76	96
11					
12					
			P	rint All	Prin

The mean BP will be displayed on the NIBP list only if it is displayed on the home display. If HR or SpO₂ is not measured, or not correctly measured at the commencement of NIBP measurement, the measured data will be displayed as "--".

For Quick SYS measurement, only the SYS (highest BP value) will be displayed.

NOTE	If the NIBP measurement was not completed, the data will not be displayed on the NIBP list. At some telemetry central monitor (ex. DS-5700), the time and measurement will be displayed as "00:00" and "" respectively.
------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------
# **ST Display**

# ST Display, Alarm Setup, etc.

This section describes the operation procedure for the ST display and alarm setup.

### To Display the ST Measurement Menu

The ST display can be accessed from the menu, or from the preprogrammed user key.



On the ST display, the averaged ECG waveform of 16 beats will be superimposed for 5 minutes. 3 frames of superimposed waveform will be displayed. Also, HR and ST level will be simultaneously displayed as graphic trend. ST1 will be measured for each lead. On the ST display, ST alarm limit and ST reference point / measurement point can be set.

NOTE
------

### **1** Press the Menu $\rightarrow$ Function $\rightarrow$ ST Display keys.



The ST display will appear.

### 2 Select the superimposed waveform.



#### **3** Select the waveform size for the superimposed waveform.



Pressing the key will sequentially change the key as follows;  $\times 1/4 \rightarrow \times 1/2 \rightarrow \times 1 \rightarrow \times 2 \rightarrow \times 4 \rightarrow \times 1/4$ .

NOTE

The selection of displayed waveform size for the superimposed waveform synchronizes with the ECG waveform size on the home display.

### **4** Select the trend scale.



Select the displaying scale for the trend.

Trend	Scale	Unit
HR	100, 200, 300	bpm
ST	±0.2, ±0.5, ±1.0, ±2.0	mV
	±2, ±5, ±10, ±20	mm

### **5** Print the ST display.



Prints the currently displayed graphic trend, waveform on the ST display.

### To Set the Reference Waveform

The reference waveform and reference / measurement point for measuring the ST level can be set on this menu.

1 Press the Menu → Function → ST Display → Reference keys to display the reference waveform setup menu.



### **2** Read the waveform by pressing the Wave Set key.



REF.PT.

-80ms

-

16 beats average of the ECG judged as normal QRS by arrhythmia analysis will be read. If during arrhythmia learning, or if VPC is present, the reference waveform setup will take for more than 16 beats.

During the reference waveform setup, the key LED will light.

### **3** Set the reference point on the ST display.

►

The reference point can be set in the range of -240 to 0ms in increments of 10ms from the peak of QRS to the P wave direction.

### **4** Set the measurement point on the ST display.



The measurement point can be set in the range of 0 to 560ms in increments of 10ms from the peak of QRS to the T wave direction. Moving the cursor will display the currently measured ST value.

**CAUTION** For the lead which the electrode is detached, the reference waveform setup cannot be performed. Check if the electrode is correctly attached, and perform the setup again.

### ST Alarm Setup

The ST upper value and lower value compared with the reference waveform will be set. The alarm value is to be set for each measurement unit (mm / mV). The upper and lower limit can be set in 1 mm / 0.1 mV increment.

**1** Press the Menu → Function → ST Display → Alarm keys to display the alarm setup menu.





2 Select ON/OFF for "ST All Alarm".

ST All Alarm	OFF

ON will generate ST alarm. However, the alarm will not generate for the lead which individual alarm is set to OFF. OFF will not generate ST alarm.

### **3** Select the lead to set the alarm limit.



Press one of the lead keys to set the alarm limit.

4. Set the upper and lower alarm limit.



Use the  $\uparrow$ ,  $\checkmark$  keys to adjust the alarm limit.

ltem	Description
Lower Alarm Limit	Select the lower alarm limit ( $\pm$ 20mm / $\pm$ 2.0mV). Alarm will be set to OFF if the value –20mm / –2.0mV or lower is selected.
Upper Alarm Limit	Select the upper alarm limit ( $\pm$ 20mm / $\pm$ 2.0mV). Alarm will be set to OFF if the value +20mm / +2.0mV or above is selected.

### 5 Select ON/OFF of "Indiv. Alarm".



ON will generate the ST alarm for the selected lead. OFF will not generate the ST alarm for the selected lead.

### **6** Select "Auto" for automatically setting the alarm limit.



Indiv. Alarm

Pressing the Auto key will automatically set the upper alarm limit to current ST value +0.2mV (+2mm), and lower alarm limit to current ST value -0.2mV (-2mm). Selecting "Auto" will automatically turn ON the ST alarm. If the upper or lower limit is OFF, the limits will remain to be OFF.

# OCRG

Display

This section describes the procedure for OCRG display. The OCRG display can be accessed from the menu, or from the preprogrammed user key.

ocre 5min 10min		Prev. Disp.
	HR	Sp02
	300	100
	250	95
	200	90
	150	30 85
	100	80
	50	75
10min 8min 6min 4min 2min	0	70
10min 8min 6min 4min 2min Hali Balan (Hali Balan) Hali Balan (Hali Balan) RESP CO2	RESF	×1 Print

On the OCRG display, compressed respiration waveform, HR trend and  $\text{SpO}_2$  trend are displayed simultaneously.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  OCRG keys to display the OCRG menu.



- **2** Select the respiration waveform.
  - RESP CO2

Select RESP or  $CO_2$  to display the compressed respiration waveform from impedance respiration (RESP) or  $CO_2$  waveform.

**3** Select the displaying duration.

Select a displaying duration from 5min or 10min.

#### **4** Select the waveform size for compressed respiration waveform.



Pressing the size key will sequentially change the waveform size.

Respiration Waveform	Size, Scale
Impedance, RESP	$\times 1/4 \to \times 1/2 \to \times 1 \to \times 2 \to \times 4 \to \times 1/4$
CO ₂	$100 \rightarrow 50 \rightarrow 100$ (unit : mmHg)
	$4 \rightarrow 8 \rightarrow 10 \rightarrow 4$ (unit : % or kPa)

### **5** Print the OCRG display.

Print

The currently displayed graphic trend and compressed waveform on the OCRG display will be printed.

**NOTE** The OCRG cannot be printed on the central monitor recorder.

# **CO Measurement**

**Measurement/Editing** 

This section explains about the cardiac output measurement using the thermodilution method, setup procedure for catheter type, and procedure for editing the measurement result.

# To Display the Cardiac Output Menu

The cardiac output menu can be accessed from the menu, or from the preprogrammed user key. Status Message Cardiac Output (CO) Cardiac Output Edit Prev. Disp. Config. Cardiac Constant (CC) CO Start WAIT CO (L/min) cc 0.542 Cardiac Index (CI) Auto Start: ON 5.32 Scale  $(L/min/m^2)$ СΙ CO OK (°C) 43.9 2.80 Thermodilution Curve **Blood Temperature** 44 . 9°C Tb Ti 0.0°C Injection Temperature 44.9↓ 0 Print 20 30 sec 10 **Press the** Menu  $\rightarrow$  Function Cardiac Output keys. 1  $\rightarrow$ Cardiac Output Config. Prev. Disp. Edit CO Start WAIT CO (L/min) CC 0.542 Auto Start: 0 32 5 Scale The cardiac output menu will be displayed. СΙ  $(L/min/m^2)$ CO OK (°C) 43.9 The message will be displayed depending on the status. 2.80 The measurement can be started when "READY" is displayed. 44 . 9°C Тb Ti 0°0.0°C Print



Message	Description
Status Message	
WAIT	Preparing for measurement. Also displayed when the catheter relay cable is not connected to the CO module, or when the thermodilution catheter is not connected.
READY	Ready to begin the measurement.
BUSY	In process of measurement.
END	End of measurement.
Result Message	
CO_OK	CO is correctly measured.
	Measurement error
	<ul> <li>The blood temperature is out of the measurable range after the injection.</li> </ul>
UPPER_FAULT	<ul> <li>The thermistor connector and relay cable is not properly connected.</li> </ul>
	<ul> <li>The line is cut on the sensor or relay cable.</li> </ul>
	Measurement error
	<ul> <li>The peak of the thermodilution curve can not be detected.</li> </ul>
PEAK_FAULT	<ul> <li>The thermistor connector and relay cable is not properly connected.</li> </ul>
	<ul> <li>The line is cut on the sensor or relay cable.</li> </ul>
	Measurement error
	<ul> <li>The blood temperature has not returned to stable condition after</li> </ul>
LOWER_FAULT	measurement.
	<ul> <li>The thermistor connector and relay cable is not properly connected.</li> </ul>
	<ul> <li>The line is cut on the sensor or relay cable.</li> </ul>
	Measurement error
SENSOR_ERROR	<ul> <li>The thermistor connector and relay cable is not properly connected.</li> </ul>
	<ul> <li>The line is cut on the sensor or relay cable.</li> </ul>
	Measurement error
OVER RANGE	<ul> <li>The CO value is out of measurable range.</li> </ul>

The result status will be displayed for 30 seconds after completion of measurement.

### **Cardiac Output Setup**

Before measuring the cardiac output, set the measurement condition such as ON/OFF of auto start, injection condition, etc.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Cardiac Output  $\rightarrow$  Config. keys.

Cardiac Output		Prev. Disp.	
Auto Start	ON	_OFE_	
Time Scale	30sec	60sec	The cardi Set the m
CC 0.54	42 Auto Input	Manual Input	displaying

The cardiac output configuration menu will be displayed. Set the measurement condition such as ON/OFF of auto start, displaying scale of thermodilution curve, CC value for injection, etc.

### 2 Set ON/OFF of "Auto Start".

		ON will automatically start the measurement without pressing the
		CO Start key.
Auto Start	ON	OFF Will start the measurement when CO Start key is pressed.
		Even when ON is selected, the measurement can be manually
		started by pressing the CO Start key.

### **3** Set the time scale.



When using the CJ-382 catheter relay cable, make sure to set the "Injectate Temperature"

Ice Room

Cancel

Inj. Temp.

CC

INPUT

After setting the CC value, press the Input key to finalize the value.

### To Measure the Cardiac Output

The measurement can be started when "READY" is displayed.

If "WAIT", "BUSY", "CO_OK" message is displayed, the measurement cannot be started. Wait until it is ready for measurement. Particularly when "WAIT" message is continuously displayed, verify that catheter relay cable is properly connected to cardiac output module, and thermodilution catheter is securely connected.

### **1** Display the cardiac output menu.

#### **2** Start the measurement.



Press the <u>CO Start</u> key, and inject as soon as the beep sound generates

If "Auto Start" is set to ON, injecting without pressing the

CO Start key will automatically start the measurement by detecting the blood temperature change.

When the measurement is complete, CO and CI value will be displayed.

#### **3** Print the measurement result.

Print
-------

Pressing the Print key will print the displayed thermodilution curve, cardiac output, and cardiac index on the recorder.

NOTE	<ul> <li>Before injecting, check that the Ti (injectate temperature) setting is correct.</li> <li>When repeatedly performing the measurement, inject at intervals of 30 to 60 seconds</li> <li>In the following cases, measurements may be inaccurate. <ul> <li>Shunt disease, tricuspid regurgitation or pulmonic regurgitation.</li> <li>During exercise stress</li> <li>As body temperature differs sequentially by exercise, constant CO value cannot be measured.</li> <li>Excessive Arrhythmia</li> <li>Body Temperature varies non-continuously as a result of arrhythmia. Accurate CO value cannot be measured.</li> </ul> </li> <li>The CI value will not be displayed unless height/weight or BSA value is input on the admit menu.</li> </ul>
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### To Edit the Cardiac Output Data

By performing the CO measurement continuously, mean CO and mean CI can be calculated by editing the measurement result.



**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Cardiac Output  $\rightarrow$  Edit keys.



The cardiac output edit menu will be displayed. The mean CO and mean CI value obtained from the measurement result will be displayed. The data can be omitted from the averaging by turning OFF the LED of the corresponded data.

### **2** Delete the measurement result.



Pressing the Delete key will delete all the measurement data of the thermodilution curve with the key LED turned OFF.

### **3** Input the data to the list.



Pressing the Mean CO Input key will input the displayed mean CO data to the list.

NOTE	If the height, weight, and BSA are changed on the patient admit menu, the average CI will be recalculated. As the CI will not be recalculated after the hemodynamic calculation, store the average CI by hemodynamic calculation before changing the height, weight, and BSA.
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

# Hemodynamics

# **Calculation/Print**

This section explains the procedure for hemodynamic calculation and printing.



### Calculation Data

Data	Description	Formula
BSA	Body Surface Area (m ² )	$h^{0.725} \times w^{0.425} \times 71.84 \times 10^{-4}$ (Dubois Formula)
СІ	Cardiac Index (L/min/m ² )	CO BSA
SV	Stroke Volume (mL/beat)	$\frac{\text{CO} \times 1000}{\text{HR}}$
SVI	Stroke Volume Index (mL/beat/m ² )	SV BSA
SVR	Systemic Vascular Resistance (dynes-sec·cm ⁻⁵ )	(MAP - CVP)× 79.90 CO
SVRI	$\begin{array}{c c} Systemic Vascular Resistance Index \\ (dynes \cdot sec \cdot cm^{-5} \cdot m^2) \end{array} SVR \times BSA \end{array}$	
PVR	Pulmonary Vascular Resistance (dynes·sec·cm ⁻⁵ )	(MPAP – PCWP)× 79.90 CO
PVRI	Pulmonary Vascular Resistance Index (dynes-sec·cm ⁻⁵ ·m ² )	$PVR \times BSA$
LVW	Left Ventricular Work (kg·m)	$CO \times (MAP - PCWP) \times 0.0136$
LVWI	Left Ventricular Work Index (kg·m/m ² )	LVW BSA
LVSW	Left Ventricular Stroke Work (g·m)	$SV \times (MAP - PCWP) \times 0.0136$
LVSWI	Left Ventricular Stroke Work Index (g·m/m ² )	LVSW BSA
RVW	Right Ventricular Work (kg·m)	$CO \times (MPAP - CVP) \times 0.0136$
RVWI	Right Ventricular Stroke Work Index     RVW       (kg·m/m ² )     BSA	
RVSW	Right Ventricular Stroke Work (g·m)	$SV \times (MPAP-CVP) \times 0.0136$
RVSWI	Right Ventricular Stroke Work Index (g·m/m ² )	RVSW BSA

<b>NOTE</b> The blood pressure unit for hemodynamics is mmHg or kPa. The unit, cmH ₂ O cannot be used.	
-------------------------------------------------------------------------------------------------------------------	--

### To Display the Hemodynamic Data

The latest 5 hemodynamics data will be displayed.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Hemodynamic keys.



The hemodynamic menu will be displayed.

### **2** Select the data to display.



Select the data to display the calculation result by pressing the data selection key.

On the data selection key, calculated date and time will be displayed.

**3** Print the calculation data.



The currently displayed hemodynamic calculation data will be printed.

### To Calculate the Newly Input Hemodynamic Data

The hemodynamic calculation can be performed using the newly input data. The data can be manually entered using the numeric keys, or the current measurement data can be automatically entered.



### 2 Automatically enter the current measurement data.



### **3** Enter the data using the numeric keys.



Input the data using the numeric keys, and press the corresponded	
key from the Height, Weight, BSA, CO, HR, MAP,	
CVP, MPAP, PCWP.	

BSA will be automatically calculated when height and weight is entered, but it can be also manually entered using the numeric keys.

### [Input Data]

Data	Description (Unit)	
Height	(cm)	
Weight	(kg)	
BSA	Body Surface Area (m ² )	
CO	Cardiac Output (L/min)	
HR	Heart Rate (bpm)	
MAP	Mean Artery Pressure (mmHg)	
MPAP	Mean Pulmonary Artery Pressure (mmHg)	
CVP	Central Venous Pressure (mmHg)	
PCWP	Pulmonary Capillary Wedge Pressure (mmHg)	

NOTE	If the height, weight, BSA is changed on the patient admit/discharge menu, mean CI will be recalculated. However, the hemodynamic will not be recalculated with the new CI data.
------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **4** Execute the hemodynamic calculation.



## To Edit the Hemodynamic Data

The hemodynamic data can be edited.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Hemodynamic keys



Select the hemodynamic data to perform editing.

**2** Press the Edit key to edit the data.



Enter the value using the numeric keys, and press the corresponded
key from Height, Weight, BSA, CO, HR, MAP, CVP,
MPAP, PCWP keys.

### **3** Recalculate the hemodynamic data.



data.

After entering the data, press the Calc. key. The calculation result will be displayed.

To cancel the calculation, press the Prev. Disp. key. The date/time will not change after recalculation.

# Ventilator

By connecting the ventilator, P-V loop (airway pressure / ventilation) and F-V loop (airway flow / ventilation) can be monitored on the ventilator display.

<b>▲</b> CAUTION	<ul> <li>cannot be displayed or printed on the ventilator numeric data display.</li> <li>For SV-900, P-V loop, F-V loop and numeric data cannot be displayed or</li> </ul>
	printed. Only the alarms will be generated.

### P-V Loop

The P-V loop is sampled each 60ms and displayed for each respiration. The beginning of the loop trace is displayed in white, and the rest of the loop is displayed in cyan. The horizontal axis shows AWP (Unit: cmH₂O), and vertical axis shows Volume (Unit: mL).



**2** Select the scale for P-V loop.



The scale will change as the Scale key is pressed.

- Vertical Axis (Volume) Select from 250 / 500 / 750 / 1000 (mL). Horizontal Axis (Pressure)
  - Select from 10 / 20 / 30 / 50 / 120 (cmH₂O).

**3** Freeze the loop drawing.



Pressing the Freeze key will stop the P-V loop drawing. Pressing the key again will resume the waveform trace.

### **4** Program the reference loop.

The control loop can be stored to see the change in P-V loop.



### **F-V Loop**

The F-V loop is sampled each 60ms and displayed for each respiration. The beginning of the loop trace is displayed in white, and the rest of the loop is displayed in cyan. The horizontal axis shows AWF (Unit: L/min), and vertical axis shows Volume (Unit: mL).

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Ventilator  $\rightarrow$  F-V keys.



The F-V loop will be displayed.

### 2 Select the scale for F-V loop.



- The scale will change as the Scale key is pressed.
  - Vertical Axis (Flow)
  - Select from  $\pm 5 / \pm 10 / \pm 20 / \pm 50 / \pm 180$  (L/min). Horizontal Axis (Volume)
  - Select from 250 / 500 / 750 / 1000 (mL).

**3** Freeze the loop drawing.



Pressing the Freeze key will stop the F-V loop drawing. Pressing the key again will resume the waveform trace.

### **4** Program the reference loop.

The control loop can be programmed to see the change in F-V loop.



Pressing the Store key will store the displayed F-V loop as the control loop. Pressing the Review key will display the stored control loop. The control loop 1 will be displayed in yellow, and control loop 2 will be displayed in green.

5 Print the F-V loop.

The currently displayed F-V loop will be printed.

NOTE The F-V loop cannot be printed on the central monitor recorder.

### **Displaying the Ventilator Measurement**

Print

The numeric data measured by the ventilator can be displayed.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  Ventilator  $\rightarrow$  Numeric keys.

entilator	P-U		V	Prev. Disp.	
TV Insp	400	mL.	RESISTANCE		
Exp	416	mL	Insp	- cnH20/L/sec	
MU	6.2	L/nin	Ехр	- cnH20/L/sec	
FI <b>0</b> 2		%	COMPLIANCE		
PRESSURE			]	- mL/cnH20	
PEAK	2	cnH20	RR 20	)	The ventilator display will appear.
PAUSE		cnH ₂ 0			· · · · · · · · · · · · · · · · · · ·
PEEP	0	cnH20			
MEAN	1	cmH ₂ 0			

**2** Print the measurement data.



The currently displayed ventilator measurement data will be printed.

## **Respiration List**

This section explains about the respiration list display and recording procedure.

### To Display the Respiration List

The respiration list display can be accessed from the menu, or from the preprogrammed user key. If the respiration data is displayed on the home display, 24 hours of data will be automatically stored and displayed in 1-minute interval.

1 Press the Menu  $\rightarrow$  Function  $\rightarrow$  Resp. List keys.



The respiration list of 17 parameters will be displayed.

2 Select the displaying interval.



Pressing the time interval key will display the time interval selection tool.

Select the time interval for the respiration list display.

- If 5M is selected, the time will be displayed in real time as follows. 10:00, 10:05, ...10:25.
- If 60M is selected, it will be displayed as 10:00, 11:00, 12:00.

If the respiration list is displayed at 10:35, the data will be displayed from 10:00.

#### 3 Switch the page.



The page will be switched by one page with the displayed time interval.

key will display the previous page listing the older data.

older data

← →

key will display the next page listing the newer data.

The respiration list will be displayed in 8 columns.

If 5-minute interval is selected and if the list starts from 10:00, 35 minutes data from 10:00 to 9:25 will be displayed in 1 page.

Pressing the **←** key will display the list of 9:20 to 8:45.





The displayed list will be shifted by one column.

← key will shift the display to older data by one column.

 $\rightarrow$  key will shift the display to newer data by one column.

Print the respiration list.

5



The currently displayed respiration list will be printed.

### The Description of the Display



If the time is before the admitted time or if monitoring is suspended, the time will be displayed as "---:---". Also, if the measured data is not displayed on the home display, or BP zero balance is not performed, the data will be displayed as "---:---".

### **Respiration List Setup**

### **1** Press the List Setup key on the respiration list display.



2 Select the display position on the list.



Maximum of 17 parameters can be displayed on the list.

### **3** Select the parameter to display.



Select the parameter by pressing the corresponded key. The display will automatically shift downward to allow continuous parameter selection. This section explains about the ST graphic trend display and recording procedure.

### To Display the ST Graphic Trend

The ST graphic trend display can be accessed from the menu, or from the preprogrammed user key. If the ST data is displayed on the home display, 24 hours of data will be automatically stored and displayed in 1 minute interval.

- 1 Press the Menu  $\rightarrow$  Function  $\rightarrow$  ST Graphic Trend keys. ST Graphic Trend Prev. Disp. 4H 🗲 Time Cursor 09/19 10:15 ➡ Zoom ST(I) Scale 0.5 nm -5 5 The ST graphic trend will be displayed. ST (II) Pressing one of the Group A, B, C, D key will switch the Scale 0.2 nn display. ST (Ⅲ) Scale 0.0 m Group Setup Print
- **2** Displaying Group D will allow parameter selection on the ST graphic trend menu.



Select the position to set the parameter.

Select the parameter and press the Close key.

### **3** Select the displaying scale.



Pressing the Scale key will sequentially switch the scale depending on the displayed parameter as follows.

Parameter	Scale	Unit
HR	100, 200, 300	bpm
ет	±0.2, ±0.5, ±1.0, ±2.0	mV
51	±2, ±5, ±10, ±20	mm

#### **4** Select the displaying time span.



Time Range

 1H
 2H
 4H

 8H
 12H
 24H

 Close

Pressing the display time range key will display the time range selection tool.

Select the time range for the ST graphic trend.

Time Span	Sample Rate
1 hour	1 min.
2 hours	1 min.
4 hours	1 min.
8 hours	2 min.
12 hours	3 min.
24 hours	6 min.

### 5 Shift the display.



The ST trend display can be shifted to the older or newer data with the displayed time span.

← key will display the older data by the selected time span.

 $\rightarrow$  key will display the newer data by the selected time span.

### 6 Move the cursor.



The cursor can be moved to older or newer data. The time and data at cursor point will be displayed.

key will display the older data.

key will display the newer data.

**7** Enlarge the display.



Pressing the Zoom key will display 1 hour of data with the cursor point at center.

Directly pressing the graph area will move the cursor position.

8 Print the ST graphic trend.

Print

The currently displayed ST graphic trend will be printed.

### ST Graphic Trend Group Setup

Each trend group displays combination of 3 parameters simultaneously. 3 types of trend group can be programmed.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  ST Graphic Trend  $\rightarrow$  Group Setup keys.



The ST graphic trend group setup menu will be displayed. On this menu, parameters to display for Group A, Group B, and Group C can be selected.

**2** Select the trend group to set the parameters.



Select the trend group by pressing the Group A, Group B, or Group C key.

Select the parameter by pressing the displayed parameter keys. Pressing the parameter key will sequentially set the 3 (three) parameters from the top.

## The Description of the Display



The measured data will be compressed for the 8-hour / 12-hour / 24-hour display.

Parameter	Compressed Form
HR	Mean Value
ST	Mean Value

# **ST Tabular Trend**

This section explains about the ST tabular trend display and printing procedure.

### To Display the ST Tabular Trend

The ST tabular trend display can be accessed from the menu, or from the preprogrammed user key. If the ST data is displayed on the home display, 24 hours of data will be automatically stored and displayed in 1 minute interval.

**1** Press the Menu  $\rightarrow$  Function  $\rightarrow$  ST Tabular Trend keys.



Each ST level of I, II, III, aVR, aVL, aVF, V₁, V₂, V₃, V₄, V₅, V₆, and HR will be displayed in list format.

### **2** Select the display interval.



Pressing the time interval key will display the time interval selection tool.

Select the time interval for the tabular trend display.

If 5M is selected, the time will be displayed in real time as follows. 10:00, 10:05, ...10:25.

If 60M is selected, it will be displayed as 10:00, 11:00, 12:00.

If the ST tabular trend is displayed at 10:35, the data will be displayed from 10:00.

### **3** Switch the page.



The page will be switched by one page with the displayed time interval.

key will display the previous page listing the older data.
 key will display the next page listing the newer data.

The ST tabular trend will be displayed in 8 columns.

If 5-minute interval is selected and if the ST trend starts from 10:00, 35 minutes of data from 10:00 to 9:25 will be displayed in 1 page.

Pressing the key will display the ST trend from 9:20 to 8:45.

### 4 Shift the displayed columns.



The displayed tabular trend can be shifted by one column.
★ key will shift the display to older data by one column.
★ key will shift the display to newer data by one column.

Print the tabular trend.



The currently displayed ST tabular trend will be printed.

Latest Measurement Date 🔍



For the data when the measurement was not performed (before admittance) or when the monitoring was suspended, the time will be displayed as "--:--".

If the lead can not be monitored depending on the ECG lead cable, the data display will be left blank, and if the ST reference point is not set, the data will be displayed "- - -".

### **Other Bed**

This section explains about the function to display the waveform and numeric data and to set alarms for other bedside monitors. To use this function, DS-LAN II or DS-LAN III wired network connection is required.

NOTE The DS-7200 system cannot connect to a wired network of AU-5500N 8ch Re set as administrator (1:N network). Even if connected, other bed display, record and other network function cannot be used.
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **Other Bed Display**

The other bed display can be accessed from the menu or from the preprogrammed user key. Also, by setting the other bed alarm ON, Other Alarm key will be displayed when other bedside monitor generates an alarm. By pressing this Other Alarm key, the display for the other bed can be accessed.



**1** Press the Menu → Function → Other Bed Display keys to display the other bed selection menu.

Other B Display	Bal			Set Ot	ner Aları	n	Prev. Disp.
BED-001	BED-002	」 BED-003	교 BED-004	山 BED-005	∟ BED-006	교 BED-007	6ED-008
LW -009 ch5001	LW -010 ch5002	LW -011 [ch5003]	LW -012 ch5004_	L₩ -013 ch5005_	L₩ -014 ch5006_	LW -015 ch5007_	LW -016 ch5008_
	11	1	1	1	1	1	
			Other	Alarm	^{-'} ΟΝ		OFF

On the other bed selection menu, select the Room / Bed ID to display. For the DS-LANII network, there are 48 beds selection, and for the DS-LANIII network, there are 100 beds selection.

The Room / Bed ID for the alarm generating bed will be displayed in red.

The bed displaying this menu will be displayed in gray.

The key LED for the bed selected as the other bed alarm generating bed will be lighted.

### **2** Press the Room / Bed ID key and access the display for the other bed.

HR 75 BP1 140/ 95(110) ¹¹ 36.8 PR 75 CUP 0 cmto 12 36.3 So2 96 MBP 138/ 93(108)
RR 33 ELCO2 3.5 * ECG2 BP1 CVP Sp02 RESP C02 Silence

ECG waveform and numeric data for the selected bed will be displayed.

If an alarm is generated for this bed, the physiological alarm / arrhythmia alarm message will be displayed.

By pressing the Alarm Silence key on the other bed display, the alarm sound for the displayed bed can be silenced.

NOTE	<ul> <li>In case of DS-LANII network, when the temperature unit is °F, the temperature data will not be displayed.</li> <li>In case of DS-LANIII network, if the measurement unit for BP (mmHg/kPa) and temperature (°C/°F) is different between the bedside monitor and the central monitor, the numeric data will not be displayed.</li> </ul>
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### **Other Bed Alarm Setup**

From the bedside monitors connected to the wired network, the bed to generate the other bed alarm and ON/OFF of other bed alarm display can be performed.



**1** Select the bed to generate the other bed alarm.



Press the Set Other Alarm key to display the other bed alarm setup menu.

Select the bed to generate the other bed alarm.

The key LED for the bed selected as the other bed alarm generating bed will be lighted.

### **2** Select ON for the other bed alarm.



Press the Prev. Disp. key and display the other bed selection menu.

Selecting ON will generate the other bed alarm when an alarm generates at the other bed.

Selecting OFF will not generate the other bed alarm.

# **Full Disclosure Waveform Recording**

This section explains about the full disclosure waveform recording function.

Maximum of 48 hours of waveforms and numeric data can be recorded on the CF card. The data recorded on the CF card can be displayed on the monitor screen and printed on the built-in recorder.

### To Record the Full Disclosure Waveform Data

To record the full disclosure waveform data, a specified CF card (FCF-16GA) is required. When using the CF card on this equipment for the first time, make sure to format the CF card in advance.

**1** Press the Menu  $\rightarrow$  System Configuration  $\rightarrow$  CF Card keys.

	Prev. Disp.	The "CF Card" screen will be displayed.
CF Card Format. (HOLD 2 SEC)		
32/64/128MB/16GB		
DS-7200→CF Card		
CF Card→DS-7200		
	(H0LD 2 SEC) 32/64/128HB/166B DS-7200→CF Card CF Card→DS-7200	LPisp. CF Card Format (HOLD 2 SEC) 32/64/128MB/1668 DS-7200→CF Card CF Card→DS-7200 Full Disc. Have Rec.

**2** Insert the CF card (FCF-16GA) to the CF card slot.

### **3** Format the CF card.



Press the CF Card Format key for more than 2 seconds. The system will automatically detect the CF card type and starts the format process.

When the format process completes, full disclosure waveform data of 48 hours/6 waveforms will be recorded according to the setup.



For details of full disclosure waveform recording setup, refer to "8. System Configuration CF Card".

### To Display the Full Disclosure Waveform





### About the Waveform Display Duration

The waveform display duration will differ depending on the waveform quantity.

Waveform Quantity 1:60 seconds per waveform

Waveform Quantity 2:30 seconds per waveform

Waveform Quantity 3: 20 seconds per waveform

Waveform Quantity 6: 10 seconds per waveform

### To Shift the Displayed Waveform

◀

◀

The displayed waveform data can be shifted to older data.



(60 seconds if 1 waveform is displayed)



The waveform will shift to older data in interval of half the displayed duration.

The waveform will shift to older data in interval of one displayed duration.

(30 seconds if 1 waveform is displayed)

The displayed waveform data can be shifted to newer data. ₩

The waveform will shift to newer data in interval of one displayed duration.



(60 seconds if 1 waveform is displayed)

The waveform will shift to newer data in interval of half the displayed duration. (30 seconds if 1 waveform is displayed)

### 2 Select the quantity of numeric data to be displayed.



Pressing the Meas. Qty key will sequentially change the quantity in the order of  $1 \rightarrow 2 \rightarrow 3 \rightarrow 4 \rightarrow 5 \rightarrow 6 \rightarrow 1$ .

#### 3 Select the numeric data box.

Press the numeric data box area to open the parameter selection window.

Full Disc. Waveform Rec.	fukuda dens	Slow UT	Upper HR a	larm HR	Print Prev. Disp.	
12/20 12:23	48 Alarm Search	24	₀ Adm./Dis. ➡ Search			
ala pla	da da da		landanda landanda		30	
ula da	de de de		la da da			
als_als_	als als als		lslsls_	ECG (II)		D
ECG ×1	Time Search Enlar	ge Meas.	1 Have 1	]	-	

Select the parameter on the selection window.

Full Disc. Waveform	Rec. FU	(UDA DENS	Slow U	T Upper HR ala	rm Print Prev. Disp.
Meas.Sele	ct				HR
HR	UPC+PACE			Adm. /Dis.	
ST-A	ST-B	ST-C		Search	
Sp02	PR_Sp02	Sp02/PR	NIBP		
RR_IMP	C02	RR_C02		lllll	
BP1	BP2	PR_IBP		plslsls	
TEMP1	E.	TEMP1,2	TEMP3	r_dr_dr_dr_	
BP3	BPA	BP5	ТЬ	r_dr_dr_dr_	
				n dradradr	
	Clo	ose		leas. 1 Wave 1 Qty 1 Qty 1	

The selected numeric data will be displayed.



### **4** Select the waveform quantity.

Wave	1
Цtу	

Pressing the Wave Qty key will sequentially change the waveform quantity in the order of  $1 \rightarrow 2 \rightarrow 3 \rightarrow 6 \rightarrow 1$ .

### **5** Select the waveform to display.

Press the waveform parameter key to open the parameter selection window.

Full Disc. Waveform Rec.	fukuda dens	Slow UT	Upper HR a	arm	Print Prev. Disp.
• •	48	24	⊒₽₩	BP1	
12/20 12:23	Alarm Search		Adm./Dis.		160
eco(")	als_als_als_		ls_ls_ls_		
als_als_	als_als_als_		lalala_	-	83
ala_ala_	als als als		la da da	-	
ul-ul-	als_als_als_		lılılı_	-	( <b>95</b> )
ala_da_	als_als_als_		h_ala_da_		200
ala_ala_	als_als_als_		lslsls_		
ECG ×1 Size ×1	Time Search Enlai	∿9e Meas. Qty	1 Wave 1 Qty	]	

Select the parameter on the selection window.

Full Disc. Waveform Rec Waveform	c. FUKUDA D	ENS <mark>Sto</mark>	W UT	Upper HR alarm	Print Prev. Disp.
ECG(1)	ECG (II)	ECG (III)	ECG (aUR)		
ECG (aVL)	ECG (aUF)	ECG (V/V1)	ECG (U2)		
ECG (U3)	ECG (U4)	ECG (US)	ECG (U6)	ll_	
Sp02	RESP	C02		- la	83
BP1	Con la	BP3	BP4	landa	
BP5	</td <td></td> <td></td> <td></td> <td>(95)</td>				(95)
	Ű			[	nnHg
					ECG (II)
				hin -	
Close				Wave 1 Qty 1	

The selected waveform will be displayed.

Full Disc. Waveform Rec. FUKUDA DENS <b>Stow UT Upper HR alarm</b>	Print Prev. Disp.
	1
48 24 0 12/20 12:23	
Search Search Search	05
	00
	<b>0</b> 3
	(95)
	\$p02
	- Shoz
ECG 1 Time February Meas. 1 Wave 1	
ECG     ×1     Time Search     Enlarge     Meas. Qty     Have Qty	

### To Change the ECG Waveform Size

#### 1 Select the ECG waveform size.



Pressing the ECG Size key will sequentially change the displayed ECG waveform size in the order of  $\times 1 \rightarrow \times 2 \rightarrow \times 4 \rightarrow \times 1/4 \rightarrow \times 1/2 \rightarrow \times 1$ . This setting is not synchronized with the ECG waveform size on the home display.

### ECG Waveform Amplitude

The ECG waveform amplitude will differ depending on the ECG waveform size.

- $\times 1/4$  : ± 6.0mV
- ×1/2 : ± 3.0mV
- : ± 1.5mV ×1
- ×2 : ± 0.75mV
- ×4 : ± 0.33mV

### To Search by Time

The full disclosure waveform of a specified time can be displayed.

	If the time information on the CF card is not correct due to the following cause, the time search operation may fail. •When the time setting is changed during the full disclosure waveform recording.
	<ul> <li>recording.</li> <li>When the time setting on the DS-7200 and the central monitor do not match when connected to the wired network.</li> </ul>

### **1** Press the Time Search key.

Time

Press the Time Search key to set the date/time for searching. Search



**2** Enter the date/time, and press the Search key to start searching.



The waveform of the specified date/time will be displayed.

Full Disc. Waveform Rec.	FUKUDA DENS	\$100 UT		BP1	Print Prev. Disp.
12/20 12:23	Alarm Search		Adm./Dis.	<b>)</b>	<b>163</b>
	ala_ala_ala			~ ~ ~	(95)
dr_dr_	ala da da				nnHg ECG (II)
ECG ×1	Time Search Enla	rge Mea	us. 1 Wave	1	

## To Search by Alarm

The full disclosure waveform data of an alarm generated point can be searched.

#### Press the arrow keys for "Alarm Search". 1

Alarm Search

keys will start the searching process. Pressing the |⇒|

The generated arrhythmia alarm or numeric data alarm will be displayed at the upper part of the screen. The alarm generated numeric data will be highlighted.



### To Search by Admit/Discharge Date

The full disclosure waveform data of an admitted/discharged date/time can be searched.

**1** Press the arrow keys for "Adm./Dis. Search".



The full disclosure waveform data at the admitted/discharged date/time will be displayed.

The patient name will be also displayed if entered at admittance.

Patient Name	
Full Disc. Haveform Rec. FUKUDA DENS STORUT Discer HR at 48 24 12/20 12:23  Rarm Search  FCG(13)  FC	RI Print Prev. Disp. HR
ECG ×1 Time Enlarge Meas. 1 Wave 1 Oty 1	ECG (II)

### To Print the Waveform

The enlarged full disclosure waveform can be printed on the built-in recorder.

#### **1** Select the quantity of waveform to print.

Press the Wave Qty key to select the quantity of waveform to display and print. When 6 waveforms are selected, only the 3 waveforms from the top will be printed.

Full Disc. Waveform Rec.	fukuda dens	Slow UT	Upper	• HR alarm	Print Prev. Disp.
◀ ◀	]⊢				
12/20 12:23	48 Alarm Search	24	Adm./Dis. Search	€	
	مام مام مام		_ courcin		
	als als als	lslslslslslsls_			
	·····				
~l~_~l~_			~l^l^		
l-l-			.ll	-l-	ECG (II)
ala_da_		ll.	.ll	-l-	
ECG ×1	Time Search Enla	arge Meas Qty	3. 1 Wa	ty 1	
•				- Jes	C h
				R	20

#### **2** Select the waveform to print.

On the parameter selection window, select the waveform to display and print.

Full Disc. Waveform Rec.	FUKUDA DENS	\$100 UT	Jøper HR a	larm BP1	Print Prev. Disp.
12/20 12:23	Alarm Search	➡ €	Adm./Dis. Search		63
l			المصام الم		83
als_ds_			la da da	. (	( <b>95</b> )
ntr_dr_ ntr_dr_	ala ala ala		ىلى_ملى_ملى باي_ملى_ملى_	ECG (II	
ECG ×1	Time Search Enla	irge Meas Qty	· 1 Wave 1 Qty 1		E

Full Disc. Waveform Rec Waveform	:. FUKUDA D	ENS <mark>St</mark> o	UT wo	Upper HR alarm	Print Prev. Disp.
Advertorin				▶ ▶ BP1	
ECG(I)	ECG (II)	ECG (III)	ECG (aUR)	Dis.	
ECG (aVL)	ECG (aVF)	ECG (Ų/Ų1)	ECG (U2)		<b>[[63</b> ]
ECG (U3)	ECG (Ų4)	ECG (US)	ECG (U6)	LL	
Sp02	RESP	C02	]	lsls	83
BP1		BP3	BP4	l-1-	
BP5	S			llr_	(95)
					nnHg
				hand [	ECG (II)
				l-i-	
	Clo	ose	Wave 1 Qty 1		

### **3** The displayed waveform can be printed.

Print

Pressing the Print key will print the displayed waveform in enlarged format.

Full Disc. Waveform Rec.	Fukuda dens	Slow UT	Upper HR alar	Print Prev. Disp.
•	48	24		BP1
12/20 12:23	Alarm Search	→ ←	Adm./Dis. 🗲 Search	
\$P02		~~~~	$\sim \sim \sim$	
~~~	~~~~~	~~~~	~~~~	<b>83</b>
~~~	~~~~	~~~~	~~~~	
~~~	~~~~	~~~~	~~~~	(95)
~~~	~~~~	~~~~		nnHg
~~~~	~~~~~	~~~~		Sp02
ECG Size ×1	Time Search Enla	rge Meas. Qty	1 Wave 1 Qty 1	

[Output Example of Full Disclosure Waveform]



To Enlarge the Waveform

The full disclosure waveform is displayed in compressed format. On the enlarged display, the selected waveform on the full disclosure waveform screen will be displayed enlarged.

QRS Classification	Full Disc. <u>Waveforr</u> Enlarge 12/20 N EC6(III)	Slow UT	KUDA DENS		Upper HR alar	Print Prev. Disp.	
Scale in Seconds	<u>*</u> µЪ-Л		^&_^	ላቤሊ		180	
Measurement Data	0 HR Sp02 PR_Sp02 ST(1) ST(1) ST(1) ST(1) ST(aUR) ST(aUF)	0.00 ST(U4	P 23 PR_IBP NT 0 0 0.00 BP1 0 0.00 BP2 0 0.00 BP3 0 0.00 BP4 0 0.00 BP5	23/10(15) / ()	o o EtC02 49 InspC02 1 RR_C02 20 T1 36.1 T2 37.2 T3 Tb	ECG (II)	
Reference 2 The disp	Defi	nition".	QRS classific m can be sl			n Arrhythmia Analysis	Arrhythmia

1 Press the Enlarge key to display the enlarged waveform.

🗲 Wave 🔿

Press the arrow keys to shift the waveform.

 \leftarrow key will display older data, and \rightarrow key will display newer data.

NOTE The enlarged waveform can be shifted in 1-minute interval.

To Print the Waveform (Enlarged Display)

The enlarged waveform (displayed) can be printed on the built-in recorder.

1 Print the displayed waveform.



Press the Print key to print the displayed waveform.



[Output Example of Full Disclosure Waveform]

BED-001 Oct. 25 2011 19:10 11 Fukuda densi) millig (++ +++)			1 1 1 1 1 1		
U x1								II 1				
n_nln_	Mr.	A h h	 21	<u> </u>	l	$^{-}$	- Al-	1	_^			
N 25mm/s FULL DISC R		- M	N				N 1			N		H

Vigilance/Vigileo List

By connecting the Vigilance, Vigilance CEDV, VigilanceII, Vigileo (oximeter /CCO measurement device manufactured by Baxter), the Vigilance data such as SvO_2 (mixed venous oxygen saturation), CO (cardiac output) can be displayed in list format.

NOTE The Vigilance list data will be erased if the power has been turned OFF for more than 5 minutes.

To Display the Vigilance/Vigileo List

The Vigilance/Vigileo list can be accessed from the menu, or from the preprogrammed user key. If the data is displayed on the home display, the 24 hours of data in 1-minute interval will be automatically stored and displayed.

lgilance Vigileo ist	Shift	⇒	Interv	/-10M	+	Page	⇒	Prev. Disp.
09/22	8:20	8:30	8:40	8:50	9:00	9:10	9:20	9:30
SV02	85	86	84	85	86	84	85	86
CC0	4.3	4.4	4.2	4.3	4.4	4.2	4.3	4.4
EDU	121	122	120	121	122	120	121	122
B_Temp C	37.1	37.2	37.0	37.1	37.2	37.0	37.1	37.2
HR	61	65	63	61	65	63	61	65
EF	61	62	60	61	62	60	61	62
SU	61	62	60	61	62	60	61	62
CC I	2.4	2.5	2.3	2.4	2.5	2.3	2.4	2.5
EDUI	66	67	65	66	67	65	66	6
ESU	46	47	45	46	47	45	46	4
SUR	1617	1617	1617	1617	1617	1617	1617	161
Sa02	96	97	95	96	97	95	96	9
SUI	41	42	40	41	42	40	41	4
ESUI	28	29	27	28	29	27	28	2
SURI	2247	2248	2246	2247	2248	2246	2247	224
CCO_STAT	4.6	4.7	4.5	4.6	4.7	4.5	4.6	4.
EDU_STAT	121	122	120	121	122	120	121	12
		rint						

1 Press the Menu \rightarrow Function \rightarrow Vigilance/Vigileo List keys.

2 Select the display interval.



Pressing the time interval key will display the time interval selection tool.

Select the time interval for the tabular trend display.

If <u>5M</u> is selected, the time will be displayed in real time as follows. 10:00, 10:05, ...10:25.

If 60M is selected, it will be displayed as 10:00, 11:00, 12:00.

If the Vigilance/Vigileo list is displayed at 10:35, the data will be displayed from 10:00.

3 Switch the page.



→ key will display the next page listing the newer data.

The Vigilance/Vigileo list will be displayed in 8 columns.

If 5-minute interval is selected and if the list starts from 10:00, 35 minutes of data from 10:00 to 9:25 will be displayed in 1 page.

Pressing the key will display the Vigilance/Vigileo list from 9:20 to 8:45.

4 Shift the displayed columns.



The displayed list can be shifted by one column.

key will shift the display to older data by one column.

→ key will shift the display to newer data by one column.

Print the Vigilance/Vigileo list.



The currently displayed Vigilance/Vigileo list will be printed.
The Description of the Display

Latest Measurement Date	Vigilance /Vigileo List	Shift	⇒	Interv	~10M	•	Page	→ [Prev. Disp.	
	- 09/22	8:20	8:30	8:40	8:50	9:00	9:10	9:20	9: 30	Latest Measurement Time
	\$v02	85	86	84	85	86	84	85	86	
	000	4.3	4.4	4.2	4.3	4.4	4.2	4.3	4.4	
	EDŲ	121	122	120	121	122	120	121	122	
	B_Temp C	37.1	37.2	37.0	37.1	37.2	37.0	37.1	37.2	
	HR	61	65	63	61	65	63	61	65	
	EF	61	62	60	61	62	60	61	62	
	SU	61	62	60	61	62	60	61	62	
	CCI	2.4	2.5	2.3	2.4	2.5	2.3	2.4	2.5	
	EDŲI	66	67	65	66	67	65	66	67	
	ESŲ	46	47	45	46	47	45	46	47	
	SUR	1617	1617	1617	1617	1617	1617	1617	1617	
	Sa02	96	97	95	96	97	95	96	97	
	SUI	41	42	40	41	42	40	41	42	
	ESUI	28	29 2248	27 2246	28 2247	29 2248	27	28 2247	29 2248	
	SURI CCO_STAT	4.6	4.7	4.5	4.6	4.7	2246 4.5	4.6	4.7	
	EDU_STAT	4.6	4.1	4.5	4.6	4.1	4.5	4.6	4.1	
	CDO-21H1	121	122	120	121	122	120	121	122	<u>-</u>
							List setup	F	Print	

For the data when the measurement was not performed (before admittance) or when the monitoring was suspended, the time will be displayed as "---:---".

Also, if the data is not displayed on the home display, it will be displayed as "- - -" on the list.

Vigilance/Vigileo List Setup

The parameter to display on the Vigilance/Vigileo list can be selected.

1 Press the List Setup key on the Vigilance/Vigileo list display.



ligilance/Vigile	o List S	etup		Prev.
\$v02				_Disp
- CC0	1	OFF		
		Sv02	Scv02	Sa02
B_Tenp	Į			
⊐ HR		02E1	D02	Ų02
J EE		SU	SU_STAT	SUI
⊐ su				
- ccı		SUI_STAT	HR	MAP
⊐ EDVI		CUP	CCO	CCO_STAT
ESU		CCI	CCL_STAT	SUR
⊐' SUR				
J Sa0₂	1	SURI	B_Temp	EF
⊐ sui	l	EF_STAT	EDU	EDU_STAT
- ESŬI				
⊐' SURI		EDVI	EDULSTAT	ESU
→ CCO_STAT		ESVI	SUU	
- EDV_STAT	1			
	-			

2 Select the display position on the list.



Maximum of 17 parameters can be displayed on the list.

3 Select the parameter to display.

Vigilance/Vigileo Lis	t Setup		Prev. Disp.
Sv02	OFF		Lorsp.
- EDU			
- B_Temp	Sv02	Scv02	Sa02
HR	02EI	D02	U02
	SU	SU_STAT	SUI
	SUI_STAT	HR	MAP
EDVI 1	CUP	CCO	CCO_STAT
ESU SUR	CC I	CC I_STAT	SUR
- Sa02	SURI	B_Tenp	EF
- SUI	EF_STAT	EDV	EDU_STAT
	EDŲI	EDUI_STAT	ESU
CCO_STAT	ESUI	SUU	
EDIL STOT	-		

Select the parameter by pressing the corresponded key. The display will automatically shift downward to allow continuous parameter selection. Blank Page

Chapter 8 System Configuration

Night Mode Night Mode Monitoring	8-2
About the Night Mode	8-2
To Set the Night Mode	
●To Start/Stop the Night Mode	
Night Mode Display Setup	
Alarm Mode Setup	
Programming the Alarm Mode About the Alarm Mode	
To Program the Alarm Mode	
To Initialize the Alarm Mode	
Monitor Setup Setup for Each Monitor	
About the Monitor Setup	
Message Icon	
Check Discharge at Power ON	
●Time/Date Setup ······	8-12
Password Setup	
Program Version	
●R.C. Setup	
●Key Mask ······	
●User Key Setup	8-13
●Alarm Pole Setup ······	8-13
Menu Setup	
Display Optimization Setup	8-14
Backup at Discharge	8-14
●Low Limit for the Alarm Volume	
Password	
Monitoring Condition after Discharge ···	
●Event Key ·····	
Drift Filter display / Exp. clock display ··	
●HR/PR Source·····	
Freeze Mode Cursor	
Parameter Key Operation	
●BP Alarm Increment ·····	8-17
OCO_2 (mmHg) Upper Limit for LAN,	o (7
Telemetry	
 Battery Operation Store all alarms to "Recall" 	
Buzzer Tone (Speaker) Failure	
Built-in Rec. Status Display	0 10
●Vigilance/Vigileo SVR, SVRI, Calc. ······	
Alarm Level	
Alarm System	
Alarm System Sound	
●DS-LAN Setup······	
RR Alarm Increment	
 Patient Name on Home Display 	
Display Mode Setup	
To Program the Display Mode	

Label Setup User Label	8-24
To Set the BP User Label	
To Set the Temperature User Label	8-24
Hospital Setup Setup for Each Hospital	
About the Hospital Setup	
Date Format	
Alarm Mute	
Arrhythmia Analysis Filter	
External Device Connection Setup	8-27
NIBP Data Erase Time	
Status Output Setup	
Measurement Unit Setup	
Telemeter Setup	
●TCON Setup ······	
●Trend Clip	8-30
BP Recording Scale	8-30
Suspend Arrhy. Analysis during	
Noise Interference	
MEAN Calculation (BP1, ART, NIBP)	
Night Mode Cancel	
●Asystole, VF, VT ······	8-31
Patient ID Transmission Starting Digit	
for DS-LAN II ·····	
Admit/Discharge Key Setup	8-32
HR/PR Low Limit during Alarm Auto Setting	
Password for Alarm Setup	
Recording Paper	8-32
CF Card Data Transfer	8-33
CF Card Format	
Full Disclosure Waveform Recording Setup	
Data Transfer (DS-7200 \rightarrow CF Card)	
Data Transfer (CF Card \rightarrow DS-7200)	
Data for Transfer	
Error Message	8-37

Night Mode

This section explains the procedure to set the night mode.

About the Night Mode

The display brightness and alarm volume can be set to night mode when turning off the light of the ward or when the patient is asleep.

The night mode can be manually set to ON, or automatically set to ON by preprogramming the time to turn ON/OFF the night mode.

Operation flow when the night mode is set to "Time Display Only"



Operation flow when the night mode is set to "Slightly Dark" or "Dark"

BED-001 FUKUDA DENSHI AAK 02/04 15:16 M 150000000 PI 116/277 Provide Trende Recall Function Prevale Trende Recall NUMP Let 01 0000 Rece-Let 01 freede Trende Denket Net Denket Net	(1) The night mode can be turned ON manually by pressing the <u>Night Mode</u> key on the menu display or the <u>Night Mode</u> key preprogrammed as user key. It can be also automatically turned ON at the preprogrammed time.
(2) (1)	 (2) During the night mode, a message, "Night Mode Active" will be displayed. To cancel the night mode, select <u>Any Key</u> for "Night Mode Cancel" on the hospital setup of the preset menu. Touching anywhere on the screen will cancel the night mode. Selecting <u>Night Mode Key</u> will cancel the night mode by pressing again the <u>Night Mode</u> key on the menu display or <u>Night Mode</u> key preprogrammed as user key.
	Hospital Setup" for procedure to cancel the night mode.

NOTE	 Even when the automatic night mode is set, the night mode can be manually set to ON using the user key. In this case, the night mode will automatically set to OFF at the preprogrammed time. The night mode cannot be set when the ventilator alarm is generated.
------	---

To Set the Night Mode

The time to start and stop the night mode, and the night mode display can be set.



2 Select Manual or Auto to start the night mode.

Auto

Manual

Niaht	Mode

Manual key will start the night mode manually by pressing the Night Mode | key set as user key. Auto will start the night mode automatically at the preprogrammed time. If automatic night mode is set to ON, the night mode can be manually turned ON using the user key or remote control.

3 Set the "Auto Start Time" and "Auto End Time" of the night mode. (Only for the automatic night mode)



Enter the hour and minute using the numeric keypad and press the Hr key, Min key for the start time and end time.

Night Mode Display Setup

1 Press the Page Down key on the Night Mode Setup (1/2).



The second page of the night mode setup menu will be displayed.

2 Set the volume for the night mode.

This volume setup will be effective for all sounds such as key sound and alarm sound.

/ol,No Change	Quiet Very Quiet Silence	
Selection	Actual Volume	
No Change	Standard volume	
Quiet	Third level from the minimum	
Very Quiet	Minimum volume	
Silence	No sound	

∧ warning	When selecting Silence, pay attention not to miss any important alarm by simultaneously monitoring the bed on other monitors such as central monitor.
	simultaneously monitoring the bed on other monitors such as central monitor.

3 Select the display brightness of the Night Mode.

Disp. No Change Slightly Dark	D ark	Time Only
-------------------------------	--------------	-----------

Selection	Actual Brightness	
No Change	Standard display	
Slightly Dark	80% of maximum brightness	
Dark	50% of maximum brightness	
Time Disp. Only	Only the time will be displayed. The message will disappear after 1 minute from starting the night mode.	

M warning	When selecting Time Disp. Only, pay attention not to miss any important alarm by simultaneously monitoring the bed on other monitors such as central monitor.
------------------	---

4 Select ON/OFF of alarm pole for the night mode.

Alarm Pole

ON OFF

Selection	Alarm Pole
ON	The alarm pole will light during the night mode.
OFF	The alarm pole will not light during the night mode.

≜ WARNING	When selecting OFF for "Alarm Pole" for the night mode, pay attention not to miss any important alarm by simultaneously monitoring the bed on other monitors such as central monitor.
------------------	---

5 Set the slave monitor display for the night mode.

		Selection	Operation			
Slave Monitor		ON	The waveforms/numeric data from the home			
	OFF(Time Only)		display will be displayed on the slave monitor			
			during night mode.			
		OFF	The slave monitor display will be OFF during night			
			mode.			
		OFF (Time	When Time Only is selected for night mode			
		Only)	display setup, only the time will be displayed on			
			the slave monitor.			
			If No Change, Slightly Dark or Dark is			
			selected, then the slave monitor display will be			
			OFF.			

		Slave monitor Setup			Operation	
		ON	OFF	Time Only	Home Display	Slave Monitor
Night Mode Setup	Time Only	0	_	_	Display will depend on the night mode setup.	The display brightness will not be dimmed and the waveforms/numeric data will be displayed.
	-	_	0	—		Display will be OFF.
		_		0		Only the time will be displayed.
	No Change/		_	_		Same as the home display.
		_	0	—		Display will be OFF.
	Dalk	_	_	0		Display will be OFF.

This section explains the procedure to program the alarm mode.

About the Alarm Mode

On the DS-7200 system, 5 patterns of alarm mode can be programmed according to the monitoring purpose.

Setting all the alarm condition for each time the patient is admitted may be troublesome. To simplify this procedure, 5 patterns of alarm mode other than default setting can be programmed according to the monitoring purpose.

By preprogramming the setups to each alarm mode, alarm setups at admitting procedure can be simplified by just selecting the alarm mode.

MARNING It is recommended to program the alarm mode in rough classification such as patient's age, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient.



To Program the Alarm Mode

Programming the alarm condition for each alarm mode can be performed on the standard alarm setup menu. The default setting can be changed for each alarm mode.



2 Select a mode for programming.



3 Program the alarm condition.





The programmed alarm condition will be displayed inside the numeric data key.

Select the mode to program the alarm condition.

The programmed value and alarm OFF mark will be displayed simultaneously.

Pressing the numeric data key will display the alarm setup menu which allows to change the alarm condition.

Pressing the Register to Mode key will register the current alarm condition to the alarm mode (1 to 5) selected at procedure 2.

4 Enter a comment.



Pressing the Comment key will display the keyboard display. Enter the comment using the keyboard.

To Initialize the Alarm Mode

The alarm mode setting can be initialized to factory default setting.



More Mode Setup Prov. Comment ALARM 1 Mode1 Mode2 Mode3 Mode4	Initialize Current Mode will initialize the currently selected alarm mode to factory default setting.
HZ/PR 40-120 ST XX CT Marm or Arrhy, épotelie 0F 3 low 0T Indo Irado irado irado NIEP 3 irado irado SpG2 isc. 00-17 irado BP1 5 0 00-180 in 00-191 BP2 X5 00-11 1 X 00-11 12 X 00-1	Initialize All Mode will initialize all alarm modes (Mode1 to Mode5) to factory default setting.
Initialize Current Mode (HOLD 2 SEC) IP /kern Tere /fere	

Reference

For factory default setting of the alarm mode, refer to "11. Technical Information".

Monitor Setup

Setup for Each Monitor

This section describes the setup items that can be set individually for each monitor.

)	
Monitor Setup 1/5 Page Down Di	ev. Monitor Setup 2/5	Page Up Page Down Disp.	Monitor Setup 3/5	Page Up Page Down Prev. Disp.
Message Icon	Password	ON OFF	Parameter Key Operation	Store Not Store
Check discharge at power ON	Discharge Mode	Suspend Admit	BP Alarm Increment	Normal
Time/Date Set Program Password Version	Event Key		CO ₂ (mmHg) upper limit for LAN, telemetry	No limit 99mmHg
R.C.Setup Key Mask	Drift Filter display/ Exp clock display	Drift Filter Disp.	Battery Operation	Normal Power Save
User Key Alarm Pole Menu Setup	HR/PR Alarm Source	ECG/ SpO ₂ SpO ₂ /	Store all alarms to "Recall".	
Display Optim, Backup at Low Limit			Buzzer Tone (speaker) Failure	Enable Disable
Setup discharge Alarm Vol.	Freeze Mode Cursor		Built—in Rec. Status Display	
Monitor Setup 4/5 Page Up Page Down Pr	ev. sp. Monitor Setup 5/5	Page Up Page Down Prev. Disp.		
Vigilance/Vigileo SVR, SVRI Calc.	RR Alarm Increment	Normal		
Alarm Level Standard User	Patient Name on Home Display			
Level				
"User" cannot be selected if connected to DS-LAN II/III If "User" is already selected when the DS-LAN II/III is connected, it will automatically change to "Standard".				
Alarm System				
DS-LAN Setup				
*To validate the setup, you need to restart the system.				
Level 3 Alarm Sound One time				

About the Monitor Setup

The monitoring condition can be set for each monitor on the monitor setup menu.

First Page (1/5)

- Message Icon
- Time/Date
- Program Version
- Key Mask
- Alarm Pole Setup
- Display Optim. Setup
- Low Limit Alarm Vol.

Second Page (2/5)

- · Password
- Event Key
- HR/PR Alarm Source

Third Page (3/5)

- Parameter Key Operation
- \cdot CO₂ (mmHg) upper limit for LAN
- Store all alarms to "Recall"
- Built-in Rec. Status Display

Fourth Page (4/5)

- · Vigilance/Vigileo SVR, SVRI Calc. · Alarm Level
- Alarm System
- · Level 3 Alarm Sound

Fifth Page (5/5)

RR Alarm Increment

- Check Discharge at Power ON
- Set Password
- · R.C. Setup
- User Key
- Menu Setup
- Backup at Discharge
- · Discharge Mode
- Drift Filter Display / Exp Clock Display
- Freeze Mode Cursor
- BP Alarm Increment
- Battery Operation
- · Buzzer Tone (speaker) Failure
- · DS-LAN Setup

Patient Name on Home Display

Message Icon



Check Discharge at Power ON

Check discharge at power ON	ON	OFF

Select ON/OFF to display message icon.

When there are many numeric data display, the parameter key size will be reduced which may disable the message to be displayed inside the parameter key.

In such case, an icon will be displayed inside the parameter key to indicate that there is a message.

The message can be checked by pressing the parameter key with the icon, and displaying the parameter setup menu.

If the parameter key size is large, a message will be displayed inside the parameter key.

The trend data and NIBP list data will remain stored even when the power is turned OFF. To start monitoring a new patient, it is necessary to perform discharge procedure on patient admit/discharge menu, and clear the data of previous patient.

This function allows to select ON/OFF of discharge confirmation display when previous data remains at power ON.

To immediately display the waveform and numeric data at power ON, select OFF. The discharge confirmation display will not be displayed and monitoring will be immediately started. Selecting ON will display the discharge confirmation display when the previous data remains at power ON.

Previous patient data still in memory.
Press "Discharge" to clear.
Discharge Initialize patient data/info, monitoring parameters, etc.
Continue Continue monitoring.

< Discharge Confirmation at Power ON>

Time/Date Setup

Time/Date	Time / Date key will display the time/date setup menu.
Time/Date Prev. Disp. 2004 09/17 09:03:00 09 Hr 03 2004 Yr 45 9 Mo 1 23 Setup 0	Enter the time/date using the numeric keypad, and press the corresponded key. For example, to change the time from 2min to 5min, enter 5 on the numeric keypad. Next, press the 02 Min key. Then, press the Setup key to finalize the setup.
malfuncti NIBP list	ne/date is not correctly set, or changed during monitoring, on may occur to NIBP measurement, periodic recording, trend, data, and age calculation from the birth date.

	The list data, and ago baloalation normano birth dato.	
▲ CAUTION	• When connected to a wired network, the same time/date with the central	
	monitor will be set. Even if the time/date is changed on the DS-7200 system,	
	it will be corrected to the time/date of the central monitor.	
	 If the time/date is changed, the time/date for the trend, NIBP list, recall data 	
	will also change.	
	• If the time/date is changed during monitoring, the patient's age will not be	
	recalculated.	

Password Setup



A 4-digit password to be entered on the preset menu can be set.

If the password is set to ON in the monitor setup menu, a password will be required to access the preset menu.

Enter the numbers, and press the Setup key.

The programmed password will be displayed when the setup is complete.

Program Version

Program Version	

ran Versen					
	Version	Date	Contrent		
DS-7200	V01-01(#0134)	2007/07/20	DYNASCOPE	DS-7200	
Bout Bersan	QQ1-Q1				
Display	V01-01(#0001)				
NIBP MAIN-CPU		0000/00/00			
NIBP SUB-CPU		0000/00/00			
HU Module	HU-73				
CO ₂ Module	Capnostat5				
Notale versi					

The software information will be displayed.

Pressing the Program Version key will display software version of the monitor, produced date, and comment.

The software version required for the DS-7200 system will be displayed.

- DS-7200 Software
- Display Unit Software
- NIBP MAIN-CPU Software
- NIBP SUB-CPU Software
- HU Module Type
- · CO₂ Module Type

The boot version will be also displayed.

Pressing the Module Version key will display the equipment information of the equipment connected to serial connector of the main unit.



Set the user key, ID and zone for remote control unit (CF-700) of option. Pressing the R.C. Setup key will display R.C. Setup display to select the ID function.

Reference For details, refer to "4. Monitoring Setup R.C. Setup"

Key Mask



Unnecessary keys on the Menu display can be erased. Pressing the Key Mask key will display Key Mask display to select the key to erase from the Menu display.



For details, refer to "4. Monitoring Setup Key Setup Erasing the Unnecessary Keys"

User Key Setup



User Key

6 or 8 user keys can be programmed to be displayed on the home display. Pressing the User Key key will display the user key setup menu.

Refer to "4. Monitoring Setup Key Setup" for details.

Alarm Pole Setup



Sets the function for the alarm pole located at the top of the monitor. Press the Alarm Pole Setup key to display the alarm pole setup menu.



Refer to "4. Monitoring Setup Alarm Pole Setup" for details.

Menu Setup

Menu Setup	Sets the key and location to display on the menu display. Press the Menu Setup key to display the Menu Setup display.
View Torus Constraint State	Reference Reference Refer to "4. Monitoring Setup Key Setup To Configure the Menu Display" for details.
Display Optimization Setur Sets the display priority when t) the display configuration is optimized.
Display Optim. Setup	Press the Display Optim. Setup to display the Display Configuration Optimization Setup Screen.
Instruction denotion forme Image: Set of the providy for degring configuration optimization. Image: Set of the providy for degring configuration optimization. Image: Set of the providy for degring configuration optimization. Image: Set of the providy for degring configuration optimization. Image: Set of the providy for degring configuration optimization. Image: Set of the provide degring configuration. Image: Set of the provide degring configuration.	Reference For details, refer to "4. Monitoring Setup Display Configuration Optimizing the Display Configuration".
Backup at Discharge Backup at discharge	Press the Backup at discharge key to display the "Backup at Discharge" screen.
Backup at discharge 1/2 Page Doum Prev. Display Config. Backup Initial Alarm Backup Initial ECG1.ECG2 Lead Backup Off Impedance Resp. Backup Initial DN/OFF Backup Initial BP Scale Backup Initial NIBP Auto Mode Backup Off Umpedance Resp. Off Initial BP Scale Backup Off Impedance Resp. Off Initial Bischup Off Initial	Switch the page using the Page Down / Page Up keys. Select whether to backup or to initialize after discharge for each item. Select Backup if you do not want to initialize the item after the discharge procedure.
Backup at discharge 2/2 Page Up Prev. EtCO2 Peak Backup 10scc Picking Duration Backup Initial	 Selecting Initial will initialize the item to factory default setting after the discharge procedure. If Initial is selected for "Display Config." and "Alarm", the currently selected display mode and alarm mode will be initialized. Select Backup or OFF for "CVA Set". Select Backup or 10sec for "EtCO₂ Peak Picking Duration".

[NIBP Auto Mode]

Whether or not to back up the NIBP Auto Mode after discharge can be selected.

Selecting OFF will turn off the NIBP auto mode after the discharge procedure.

Selecting Backup will back up the NIBP auto mode even after the discharge procedure. It will function even if the patient is not admitted.

NIBP Auto Mode

⊐ Backup	OFF	
- Backup		
(Resume auto mode		
by manual measurement.)		

Backup (Resume auto mode by manual measurement) will resume the NIBP auto mode when the next admitted patient has first started the manual measurement. NIBP will be periodically measured at the same interval with the previous patient. Until the NIBP auto mode is resumed or the interval is changed, "Standby" will be displayed inside the NIBP numeric data box.

•Low Limit for the Alarm Volume

Low Limit for the Alarm Volume The low limit for the alarm volume range on the Prev. Disp. Alarm Priority "Tone/Vol." screen can be set. Hi The alarm volume range can be changed for each alarm Level 1 Sound Test Inc. Dec. level. The adjustable alarm volume range will be indicated by a yellow underline in the "Tone/Vol." screen. Level 2 Sound Test Inc. Dec. Level 3 Sound Test Inc. Dec. Lo Decreases the low limit. Current low limit. Increases the low limit. Tests the set Low Limit Alarm Vol. Level 1 Dec. Inc. Sound Test

Reference

For the night mode setup, refer to "8. System Configuration Night Mode

Night Mode Display Setup"

Password

Sets the password requirement to access the preset menu.

Password

ON	OFF

Selecting ON will require to enter password to access the preset menu.



A 4-digit password can be set on the monitor setup menu. Also, "7200" can be used as maintenance password.

Monitoring Condition after Discharge Sets the monitoring condition after the patient has discharged.

Sets the monitoring		Admit will cont Suspend will s Numeric data dis	tinue monitoring after of suspend monitoring aft play will be erased and periodic measuremen be performed.	er discharge. d alarm
Discharge Mode	Suspend Admit	м	Monitoring is suspended. Resume	
Event Key		_		_
•	ON OFF ala	N will display th rm occurrence.	e event key on the hor ay the event key on th	
Pressing the event	key will suspend the alarm	E	ay the recall menu. vent Key	
^{BED-001} FU	KUDA DENSHI	Aault	1	2/11 14:51 M
●Drift Filter display Drift Filter display/ Exp clock display	y / Exp. clock display	· · · · · · · · · · · · · · · · · · ·] will display the drift f .] will not display the enlarged time instead.	drift filter status
02/04 16:1 M Drift-		^{02/04} 1 MR	6:16 Enlarge Display	ed Time
HR/PR Source This setup will allow	w HR/PR source selection o	of ECG/SpO ₂	or ECG/SpO ₂ /BP.	
HR/PR Source	ECG/ ECG/ SpO2. SpO2 BP ECG/Sp SpO2, or B SpO2, or B	O ₂ /BP will allow	R/PR source selection	

The HR/PR source selection can be performed on each configuration menu of ECG/SpO₂/BP1 (or ART). For the setup procedure, refer to "6. Parameter Setup HR/PR Alarm Source".

A CAUTION	If HR/PR alarm source is BP (Or, if Auto selects BP for HR/PR source), ECG waveform will not be transmitted on a wired network.
NOTE	 If HR/PR alarm source is BP, and ART is not selected as the first BP label, BP1 will be the HR/PR alarm source. If BP1 is not measured at this time, PR_IBP value will be blank. If HR/PR alarm source is BP, PR_IBP value will be displayed as HR on the wired network central monitor.

Freeze Mode Cursor

When the Freeze key preprogrammed as user key is pressed, Freeze Mode Cursor will be displayed on the home display. By moving this cursor, BP value at cursor position can be displayed, and interval time between the cursors can be measured.

Freeze Mode Cursor



For details of Freeze Mode Cursor, refer to "4. Monitoring Setup Display Configuration Freeze Mode Cursor Display".

Parameter Key Operation

Whether or not to store the screen accessed from the parameter key (numeric data box) can be selected. By storing the parameter key operation, the previously accessed screen can be directly displayed.

Parameter Key Operation

■ Not Stor

ON

Store will store the screen accessed from the parameter key. The next time the parameter key is pressed, the previously accessed screen will be directly displayed.

Not Store will not store the parameter key operation. The next time the parameter key is pressed, the initial screen (ex. ECG menu if HR parameter key is pressed) will be displayed.

BP Alarm Increment

The BP alarm increment can be selected from Normal or Small.

	When Normal is selected	When Small is selected		
0 to 50mmHg	2mmHg increment	1mmHg increment		
55 to 300mmHg	5mmHg increment	mining increment		
0 to 7kPa	0.2kPa increment	0.1kPa increment		
7.5 to 40.0kPa	0.5kPa increment	0. IKPa increment		

BP Alarm Increment

Normal		Select	Normal	or	Small	according	to	the
Normai	Small	monitori	ng purpose.					

CO₂ (mmHg) Upper Limit for LAN, Telemetry

When the measurement unit of CO_2 is mmHg, whether or not to limit the CO_2 value to 99mmHg when transmitting to the central monitor can be selected.

CO₂ (mmHg) upper limit for LAN, telemetry No limit will transmit the actual CO₂ value to the central monitor even if the value is 100mmHg or above. 99mmHg will transmit the CO₂ value as 99mmHg if the value is 100mmHg or above.

Battery Operation

The battery operation mode can be selected from "Normal" or "Power Save".

The battery operation mode can be selected from "Normal" or "Power Save".
Battery Operation Image: Power Save Save Save Save Save Save Save Save
●Store all alarms to "Recall"
Store all alarms to "Recall". ON OFF OFF Will store all alarms to "Recall". OFF Will store only the selected alarm factors to "Recall".
Buzzer Tone (Speaker) Failure
Whether or not to generate a buzzer tone during speaker failure can be selected. Buzzer Tone (speaker) Failure
 Enable will generate a buzzer tone instead of an alarm sound under the following condition. Speaker failure Alarm sound level is not set to the lowest level, or "Alarm Mute" (Hospital Setup) is set OFF. If the monitor is in a night mode, the night mode volume is not set to Silence. Alarm (level 1, 2, 3 or ventilator alarm) is generating.
If the buzzer tone is generated at alarm generation, it can be silenced by pressing the Alarm Silence key. Disable will not generate a buzzer tone even during speaker failure.
Built-in Rec. Status Display
Built—in Rec. Status Display ON OFF OF OFF ON will display built-in recorder status message on the home display. OFF will not display built-in recorder status message on the home display.
●Vigilance/Vigileo SVR, SVRI, Calc.
The source of SVR, SVRI, MAP, CVP value for the Vigilance/Vigileo list can be selected from Vigilance or DS-7200.
Vigilance/Vigileo SVR, SVRI Calc. Vigilance/Vigileo ist. Vigilance/Vigileo Vigileo
●Alarm Level
The alarm level for numeric data alarm (HR/PR, SpO ₂ , BP, NIBP, RR, EtCO ₂) and arrhythmia alarm can be changed.
Alarm Level Standard Level "User" cannot be selected if connected to DS-LAN II. If "User" is already selected when the DS-LAN II.
is connected, it will automatically change to "Standard".

Reve Level texture Level 1 Level 2 SpO2 Level 1 Level 2 BP1 Level 1 Level 2 BP2 Level 1 Level 2 BP3 Level 1 Level 2 BP4 Level 1 Level 2 BP5 Level 1 Level 2 RR Level 1 Level 2 EtCO2 Level 1 Level 2	Tachy. Sync. To HR/PR Level Brady Even 1 Level 2 Run Level 1 Level 2 Run	If User is selected, press the Level key to display the Alarm Level Setup menu. Select Level 1 or Level 2 for each alarm.
A CAUTION	If the monitor is cor	nnected to wired network, User cannot be selected.

Alarm System

The alarm system such as alarm sound and alarm indicator will differ depending on this selection.

Alarm System

FUKUDA Denshi iec

WARNING When "Alarm System" setting (IEC/FUKUDA DENSHI) is changed on the Monitor Setup menu, make sure to check the alarm sound and alarm indicator.

[Tone/Volume Setup]

	•			
	When FUKUDA DENSHI is set	When IEC is set		
Sound (1) When the tone setting is set between the 1st and 4th level from the lowest level (2) When the tone setting is set to the 5th level or above from the lowest level				
Level 1	 (1) Continuous tone with alternate high and low pitch sound (2) Continuous rapid tone 			
Level 2	(1) 5 seconds interval alternate high and low pitch sound(2) 5 seconds interval rapid tone	5 seconds interval beep tone		
Level 3	 Single beep tone or 15 seconds interval alternate high and low pitch sound Single rapid tone or 15 seconds interval rapid tone (*) 	15 seconds interval beep tone (*)		
Volume Setup				
Level 1 Level 2 Level 3 Setting is possible.		The volume for low level alarm cannot be set higher than the high level alarm.		
Tone Setup				
Level 1		Setting is possible.		
Level 2 Setting is possible. Level 3		Setting is not possible. (Setting for Level 1 will be applied.)		
		Setting is not possible. (Setting for Level 1 will be applied.)		
Other Setup				
Other Bed Alarm	Setting is possible.	Setting is not possible. (Setting for Level 1 will be applied.)		
Ventilator Alarm	Setting is possible.	Only ON/OFF setting is possible. (For tone and volume, setting for Level 1 will be applied.)		
Alarm Mute (Hospital Setup)	Setting is possible.	Setting is not possible. (Fixed to "OFF")		



* The time interval for Level 3 alarm sound can be set. The setting can be performed under the "Monitor Setup" menu. For procedure, refer to "8. System Configuration Monitor Setup ●Level 3 Alarm System Sound" (Default: One time)

[Alarm Pole Setup]

Setup	When FUKUDA DENSHI is set	When IEC is set
Sync. with Alarm		Setting is not possible. (Fixed to "ON")
Alarm Type	Sotting in possible	Setting is not possible. (Fixed to "Level 1, 2 and 3")
Ventilator Alarm	Setting is possible.	Setting is not possible. (Fixed to "ON")
Pattern Setup		Setting is not possible.
Sync. with HR		Setting is possible.

A CAUTION	The alarm priority is high for level 1 (life threatening alarm), medium for level 2 (cautionary alarm), and low for level 3 (treatment needed alarm).
------------------	---

Level 3 Alarm Sound

The alarm generating time interval for Level 3 alarm can be set.

Level 3 Alarm Sound

```
One time 15s inter∨.
```

If One time is selected, the alarm sound will generate only one time.

If <u>15s interv.</u> is selected, the alarm sound will generate in 15 seconds interval.

DS-LAN Setup

The DS-LAN network type can be selected.

DS-LAN Setup

*To validate the setup, you need to restart the system.

NOTE	 When the "DS-LAN Setup" is changed, make sure that the same setting is made on the central monitor. If the setting is different, proper communication cannot be performed. The following central monitors can connect to DS-LANII network only. When connecting these central monitors, make sure all monitors in the same wired network is set to DS-LANII. DS-5700, DS-5800N/NX/NX^{MB}, DS-7600/7600W (software version of V05 and prior) To validate the "DS-LAN Setup", it is necessary to restart the system. Make sure to restart the system when the setting is changed for "DS-LAN Setup".
------	--

ORR Alarm Increment

The RR alarm increment can be setup as Normal or Small.

	Normal	Small
Adult	5Bpm increment	1Bpm increment
Child/Neonate	2Bpm increment	1Bpm increment

RR Alarm Increment	Normal	Small	Select	Normal	or	Small.
--------------------	--------	-------	--------	--------	----	--------

Patient Name on Home Display

ON/OFF of the patient name on home display can be selected.

```
Patient Name
on Home Display
```

ON OFF ON will

ON will display the patient name. OFF will not display the patient name.

Display Mode Setup

On the DS-7200 system, 5 patterns of display configuration can be programmed according to the monitoring purpose.

Setting the display configuration for each time the patient is admitted or each time the parameter is added or deleted may be troublesome.

To simplify this procedure, 5 patterns of display configuration other than default setting can be programmed according to the monitoring purpose.

By preprogramming the configuration to each display mode, display configuration setups at admitting procedure can be simplified by just selecting the display mode.

It is recommended to program the display mode in rough classification such as patient's condition, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient.



To Program the Display Mode

Programming the display configuration for each display mode can be performed on the display configuration setup menu. The default setting can be changed for each display mode.

1 Press the Menu \rightarrow System Configuration Pre-Set Display Mode Setup keys. \rightarrow \rightarrow Display Hode Sets Frer. Enter a comment Comment for each mode. Mode1 Mode2 Mode3 Mode4 Mode5 Change Setup Register to Mode Program the display configuration. Standard, 12-lead, Extension 1, Veve CCGLOP Overlap1.5p02.855P Manamic H8, ST-11, ST-8, S1-0, 591, Sp02, 59, 109 Select a mode for Extension 2, and Enlarge can be NAME BELIEVES programmed. programming. Nonemo URUBPI Spozi NUPI UPO POGLISTI NUPI. ILIPTI 2188, INP. NUPI LIS ve FCGL.8P Overlap1.5p0+.855P Manaent: HR.NIBP.Selly.BP1.6P2.TETP1.TENP2.00/.BR.03 Det 0 kave ECG1.8P Overlap1.5p0g.03g Numeric IR: Sp02.N3IP.3R_IRP Nove FEGL, Sp07, RESP Short Irensi OH Grai OFF Navo Lino Tridenoso _____ 2 ECG1.LCG2 BP1.UP2 Have Uty Veveforn BP Overlap I BP Durchus?

2 Select a mode for programming.

<u>Mode1</u> <u>Mode2</u> <u>Mode3</u> <u>Mode4</u> <u>Mode5</u> Select a mode to program the display configuration.

NOTE On the display mode setup menu, the setup of currently selected mode will be displayed. Changing the mode and returning to the home display will set the display configuration to the setup of the last selected mode.

3 Program the display configuration for the mode.

Change Setup	Press the Change Setup key.
Display Kode Standard Setup	Set the display configuration.
Ext1 Ext2 Enlarge Auto	For display configuration setup procedure, refer to "4. Monitoring Setup Display Configuration"
<pre>view of the second second</pre>	Pressing the Register to Mode key will register the current display configuration setup to the display mode (1 to 5) selected at procedure 2.
Comment $\underline{D \mid S P \mid A Y 1}$ Comment $\underline{D \mid S P \mid A Y 1}$ Comment $\underline{D \mid S P \mid A Y 1}$ Comment $\underline{D \mid S P \mid A Y 1}$ Frase $1 2 3 4 5 6 7 8 9 0$ $0 W \in \mathbb{R} T Y \cup 1 0 P$ $A S D F G H J K \sqcup X \leftarrow \rightarrow$ $Z \times C \vee B N M, . \downarrow$ \underbrace{R}_{R}^{R}	Pressing the <u>Comment</u> key will display the keyboard to enter the comment.

Label Setup

This section describes the procedure to program the user label for BP and TEMP.

To Set the BP User Label

Any 3 letters can be programmed as BP user label. Up to 2 BP user labels can be set.



2 Set the user label for Label 1 or Label 2.

Enter 3 letters using alphabet, number, and symbol keys.

To Set the Temperature User Label

Any 3 letters can be programmed as TEMP user label. Up to 2TEMP user labels can be set.



2 Set the user label for Label 1 or Label 2.

Enter 3 letters using alphabet, number, and symbol keys.

Hospital Setup

Setup for Each Hospital

This section explains about the different setup for each hospital.

Hospital Setup 1/3		Page Down Prev. Disp.	
Date	05/09 M	ay.09 09 May.	
Alarm Mute		OFF	
Arrhy, Analysis Filter	Disp Waveform	Fixed	
Ext. Device Connection	NIBP Data Erase Time	Status Output Setup	
Unit	Telemeter Setup	TCON Setup	
Hospital Setup 3/3	Page Up	Prev. Disp.	
DS-LAN Pat. ID	Tx 🔶	1 char. 븆	

+

Full

OFF

_____ 30bpm

ON

۵4

1 char. 🔿

Light

40bpm

OFF

Hospital Setup 2/3	age Up P	age Down Disp
Trend Clip	ON	OFF
BP Record Scale	40mm	
Suspend Arrhy, Analysis during Noise Interference	ON	OFF
MEAN Calculation	- Wave	Calc.
Night Mode Cancel	Any Key	⊐ Night Mode Key
Asystole , VF, VT (Neonate, only Asystole)	ON	ON/OFF

About the Hospital Setup

Admit/Discharge

HR/PR Low Limit during Alarm Auto Setting

Key Setup

Password for

Alarm Setup

Rec. Paper

A different monitoring condition can be set for each hospital.

First Page (1/3)

- · Date
 - Arrhythmia Analysis Filter
 - NIBP Data Erase Time
 - Unit
 - TCON Setup

Second Page (2/3)

Trend Clip

· BP Record Scale

· Admit/Discharge Key Setup

- · Suspend Arrhy. Analysis during Noise Interference
- MEAN Calculation Night Mode Cancel
- · Asystole, VF, VT

Third Page (3/3)

- DS-LAN Pat. ID Tx
- · HR/PR Low Limit during Alarm Auto Setting
- Password for Alarm Setup · Rec. Paper

- · Alarm Mute
- Ext. Device Connection
- Status Output Setup
- · Telemeter Setup

Date Format

Date

05/09 May.09 09 May.

The date format for display and recording can be selected from 3 formats.

Alarm Mute

Alarm Mute

ON	OFF

The alarm sound can be silenced at time of alarm occurrence. To avoid shocking the patient from the alarm sound, or if alarm sound is not necessary, this function can be used. This setting will not affect recall and alarm recording.

When <u>IEC</u> is selected for "Alarm System" on the Monitor Setup menu, ON/OFF setting for "Alarm Mute" is not possible.



•Arrhythmia Analysis Filter

Arrhy, Analysis Filter Disp Waveform Fixed The ECG filter to perform arrhythmia analysis can be set. Disp. Waveform will set the filter mode selected on admit menu or ECG configuration menut.

Fixed will set the filter to 1.0 to 30Hz regardless of the filter mode selection.

External Device Connection Setup



Status II 1
Same function cannot be set to multiple ports. If function is duplicated, one port turns OFF. (10ximetrix33 and IO-vue] can be selected at the same time.)
OFF
Vent- SU-300 Servo-1/s SU-900 PB-740/ Alanm 780/840
-Sv02/CC0
Usilance Oximetrix3 0-vue 02 Computer
Others-
BIS
Close
[Setun Screen]

(Setup Screen)

Port A Same Function cannot be set to multiple ports. If function is duplicated, one port turns OFF. (DXinetrx3) and ID=vuel can be selected at the same time.) OFF -Vent SU=300 Servo=1/s PB=740/ 760/840 Evita4/XL Savina PB=7200 -SVO2/CCO -Usinetrix3 0=vue D2 Computer Close	
IF function is duplicated, one port turns OFF. (Usinetrix 3) and (D-vue) can be selected at the same time.) OFF Vent SU-300 Servo-v/s PB-740/ 760/840 Evita 4/XL Savna PB-7200 SvO2/CCO Usalance Dximetrix 3 O-vue D2 Computer	Port A
Vent SU-300 Servo-v/s PB-740/ 760/840 Evite 4/XL Savina PB-7200 SvO2/CCO Usalance Dximetrix3 0-vue D2 Computer	If function is duplicated, one port turns OFF.
SU-300 Servo-r/s PB-740/760/840 Evita 4/XL Savna PB-7200 SVO2/CCO Valence Dximetrix3	
Sub-300 Server/s 760/840 Exit ad /XL Savna PB-7200 CSVO2/CCO Usalance Dxinetrix 3 Usalance Dxinetrix 3 0-vue 02 Computer	-Vent
Sub-300 Server/s 760/840 Exit ad /XL Savna PB-7200 CSVO2/CCO Usalance Dxinetrix 3 Usalance Dxinetrix 3 0-vue 02 Computer	
Zaura Savina PD-7200 SvO2/CCO Usakance Diametrix3 Usakance Diametrix3 Diametrix3	SU-300 Servo-i/s PB-740/ 760/840
Zaura Savina PD-7200 SvO2/CCO Usakance Diametrix3 Usakance Diametrix3 Diametrix3	
Usakance Oxinetrix 3 O-vue 02 Computer	
Usakance Oxinetrix 3 O-vue 02 Computer	-Sv0 ₂ /CC0
	Vigilance Ovimetrix 3 Onvine 0.2 Computer
Close	/Vigileo
Close	
Liose	Olara I
	Close

For StatusII 5 port and Multiport B, only one of either can be used. Press the Setup key to select the using port from StatusII 5 or Multiport B.

	COM1	COM2	COM3	atus]]	Port A	Port B	
	COMI	CONIZ	CONIS	1	2 to 5	FUILA	FULL
PC Comm.	Yes	Yes	Yes	No	No	No	No
SV-300	No	No	Yes	Yes	Yes	Yes	Yes
Servo-i/s	No	No	Yes	Yes	Yes	Yes	Yes
SV-900	No	No	Yes	Yes	Yes	No	No
PB-740/760/840	No	No	Yes	Yes	Yes	Yes	No
Evita4/XL/2dura	No	No	Yes	No	Yes	Yes	Yes
Savina	No	No	Yes	No	Yes	Yes	Yes
PB-7200	No	No	No	No	No	Yes	No
Vigilance/Vigileo	No	No	Yes	Yes	Yes	Yes	Yes
Oximetrix3	No	No	Yes	Yes	Yes	Yes	Yes
Q-vue	No	No	Yes	Yes	Yes	Yes	Yes
Q2 Computer	No	No	Yes	Yes	Yes	Yes	Yes
BIS	No	No	No	Yes	Yes	No	No
HLX-561/HLX-801	Yes	No	No	No	No	No	No
TCON (HTC-702)	No	Yes	No	No	No	No	No

[Equipments that can be connected to each port]

•NIBP Data Erase Time



Set the time to erase the NIBP data. Press the NIBP Data Erase Time key to display the setup menu for the erase time.

Select the time from	10min ,	30min ,	60min ,	24hour.
When the selected tir	ne passes	, the NIBF	data will	be erased.

Status Output Setup



atus Output Se	etup		Prev Disp
Super Module(H	S-700) Status Po	rt	
Sync Sienal	Signal output	OFF	
Output	Output logic	Negative logic	
DSC-7300 Sta	tus Port		
DSC-7300 Sta	tus Port		
Alarm	Alarm output	OFF	
		OFF Negative logic	
Alarm	Alarm output		
Alarm	Alarm output		

This device is capable to output the synchronized signal (HR, RR) and alarm.

Also, this device can output the generating alarm to external device. Pressing the <u>Status Output Setup</u> key allows to set the details of synchronized signal output and alarm output.

Pressing the Sync Signal Output key will display the menu to select the synchronized signal (HR, RR) and output logic (positive logic, negative logic).

Pressing the Alarm Output key will display the menu to select the alarm to output and output logic (positive logic, negative logic).

Synchronized Signal Output Setup

Status Output Setup(Sync Signal) Prev. Disp.	Select the output signal from HR, RR.
Genal output HR RR	HR will output the synchronized signal according to the selected HR source (ECG, SpO ₂ , BP1).
hutput losic Positive losic Negative losic	RR will output the synchronized signal according to the selected
	RR source (impedance, CO ₂). Select the output logic from Positive Logic , Negative Logic .
	Positive logic outputs the synchronized signal in plus, and negative logic outputs the synchronized signal in minus.

Alarm Output Setup

Select the alarm to output.

	Level 1 will output the signal when the level 1 alarm generates.
Status Output Setup(Alarm)	Level 1 and 2 will output the signal when level 1 or level 2 alarm
filarm output	generates. Level 1, 2 and 3 will output the signal when level 1, level 2, or level 3 alarm generates.
Level1 and 2.3	APNEA will output the signal when apnea alarm generates. Select OFF if not necessary to output the alarm.
Pulse Result ogc	Select the output logic from Positive Logic, Negative Logic, Pulse.
	Positive logic outputs the signal in plus, and negative logic outputs the signal in minus.
	A square wave of 440ms cycle is output for Pulse.
Refer to "11. Technical output signal.	nformation External Connection" for connector pin assignments of the

The equipment status alarm will be output as level 3. Select Level 1, 2, and 3 when outputting the equipment status alarm.

Measurement Unit Setup

Select the unit for the measurement.

Changing the unit will erase the tabular trend data and graphic trend data.

Also, perform the alarm setup again as alarm condition should be set for each measurement unit.



Unit

BP : Changes the unit for BP and NIBP.

CVP : Changes the unit for CVP. (When BP label is CVP.)

- TEMP : Changes the unit for temperature.
- ST : Changes the unit for ST measurement.

Telemeter Setup



Suspend Arrhy. Analysis during Noise Interference

	Suspend Arrhy, Anal during Noise Interfere	will v ysis ON OFF OFF OFF OFF OFF OFF Wh OFF Wh	ring the arrhythmia analysis, the following conditions be detected (every second) as noise. When poor signals are frequently detected from the bicked up QRS candidates. When the baseline drift duration is too long. When spike noises similar to QRS and baseline noises are frequently detected. DFF will suspend arrhythmia analysis for 1 second en a noise is interfering. N will suspend arrhythmia analysis for 5 seconds en a noise is interfering. In addition, more noises will detected when ON is selected.
	▲ CAUTION	under Heepitel Setur (Dress)	
●ME	AN Calculation		e mean blood pressure value of BP and NIBP can be

Wave Calc.

The mean blood pressure value of BP and NIBP can be selected to be measured from the waveform from calculation.



- Mean BP = (Systolic BP + Diastolic BP \times 2) \div 3
- Wave : Measures the mean BP as follows.



Night Mode Cancel

Night Mode Cancel

(ART,NIBP)

Select the procedure to cancel the night mode when "Slightly Dark" or "Dark" is set.

Any Key will cancel the night mode by pressing anywhere on the screen.

Night Mode Keywill cancel the night mode by pressing theNight Modekey preprogrammed as user key orNight Modekey on the menu display.

●Asystole, VF, VT

Asystole , VF, VT	turned off in the arrhythmia
(only asystole for neonate)	F ON/OFF will allow the al
	be turned ON or OFF

J Night Mode Key

Any Key

ON will not allow the alarms for asystole, VF, VT, Slow_VT to be turned off in the arrhythmia alarm setup menu.

ON/OFF will allow the alarms for asystole, VF, VT, Slow_VT to be turned ON or OFF.

Reference

Refer to "4. Monitoring Setup To Set ON/OFF of Arrhythmia Alarm" for arrhythmia alarm setup.

Patient ID Transmission Starting Digit for DS-LAN II

On the DS-7200 system, patient ID of up to 20 digits can be set, but only 10 digits can be transmitted on the DS-LANII network. This setup allows to set the starting digit of the 10 digits to be transmitted. On the DS-LANIII network, if Center is selected for the recorder and recording is started on the DS-7200, the central monitor recorder can print only up to 10 digits. This setup will set the starting digit to be printed. However, all 20 digits can be displayed on the central monitor.

	DS-LAN Pat. ID Tx	(♠ 1 char. ♣ 🛄	e starting digit will shift to left. e starting digit will shift to right.			
	NOTE	On the DS-7200 system, pa can be transmitted on the D	atient ID of up to 20 digits can be set, but only 10 digits oS-LANII network.			
		<u>Mode</u> key Monitor / ESIS / DIAG	on the admit menu can be selected.			
HR	HR/PR Low Limit during Alarm Auto Setting					
			OFF will not restrict the HR/PR low limit			





Password for Alarm Setup

Whether or not to require password for alarm setup menu can be selected.

Password for Alarm Setup

ON	OFF

ON will require password for alarm setup menu. OFF will not require password for alarm setup menu.

Recording Paper

Recording paper size can be selected.

Rec. Paper



A4 will record on the A4 size paper.

LETTER | will record on the letter size paper.

CF Card

This section explains about transferring the setup data and recording the full disclosure waveform using the optional CF card. To set all the monitors in the same ward to the same alarm setup and display configuration may take large amount of time. However this process can be simplified by performing the setup on one monitor, and copying the data to all the other monitors using the CF card.

Also, by using the CF card for full disclosure waveform recording, full disclosure waveform (48 hours, 6 waveforms) can be recorded and displayed on the monitor.

▲ CAUTION	 Use only the specified CF card. For data transfer: FCF-128 For full disclosure waveform recording: FCF-16GA When using the CF card for the first time, make sure to format the CF card on the using equipment. When the CF card is formatted, all the data stored in the CF card will be erased. The CF card intended for full disclosure waveform recording cannot be used for data transfer. The CF card intended for data transfer cannot be used for full disclosure waveform recording.
	When using a CE card with write protect function, make sure to cancel the write

- **NOTE** When using a CF card with write protect function, make sure to cancel the write protect function before data transfer.
- **1** Press the Menu \rightarrow System Configuration \rightarrow CF Card keys.



CF Card Format

1 Insert the CF card to CF card slot 1.

For data transferFCF-128For full disclosure waveformFCF-16GA

2 Format the CF card.

CF Card Format
(HOLD 2 SEC)

Pressing the CF Card Format key for more than 2 seconds will automatically determine the CF card type (32/64/128Mbyte, 16Gbyte), generates a beep sound and starts formatting.

3 A beep sound will generate again to notify that the format process is complete. </pr

The CF card can be used for data transfer for setup data and patient data.

Refer to "Data Transfer (DS-7200 \rightarrow CF Card) or "Data Transfer (CF Card \rightarrow DS-7200) in this section.

<In case of CF card for full disclosure waveform recording> The CF card can be used for full disclosure waveform recording. As the full disclosure waveform will be automatically recorded, leave the CF card inserted in the CF card slot.

To select the waveform for recording, refer to the next section.

Full Disclosure Waveform Recording Setup

Select the waveform for full disclosure waveform recording.

Full Disc. Have Rec.	Press the \uparrow Up key to increase the priority of the
Set the waveform priority. *6 waveforms of high priority will be recorded. Priority [11] $ECG(II)$ [11] $ECG(II)$ [21] $Sp02$ [12] $ECG(III)$ [31] $RESP$ [13] $ECG(II2)$ [41] $BP1$ [14] $ECG(II2)$ [51] $BP2$ [15] $ECG(II4)$ [61] $CO2$ [16] $ECG(II5)$ [7] $ECG(II/21)$ [17] $ECG(II6)$ Set the $CO2$ measurement unit.	Press the <u>→Down</u> key to increase the priority of the selected waveform. Press the <u>→Down</u> key to decrease the priority of the selected waveform. The 6 waveforms of high priority will be recorded on the CF card.
[8] ECG (aUR) [18] BP3 mnHg [9] ECG (aUF) [19] BP4 kPa [10] ECG (aUL) [20] BP5 2%	



For the CO_2 waveform, select the measurement unit from mmHg/kPa/%.

Data Transfer (DS-7200 \rightarrow CF Card)

The data can be transferred from the monitor to the CF card.

- **1** Insert the CF card to CF card slot 1.
- **2** Write data to the CF card

	DS-7200→CF Card	The data will be transferred from the monitor to the CF card. Press the DS-7200 \rightarrow CF Card key to select the data type to transfer.
3	Select the data type.	Setup data key will transfer the setup data such as alarm setup, display configuration, parameter setup to the CF card. Patient data key will transfer the patient data such as admitting data (name, age, ID, etc.), graphic trend and tabular trend. Recall data, OCRG, ST measurement and short trend data cannot be transferred. If the CF card is not inserted to the card slot, a message will be displayed.
	CF card slot:0UT Insert card into CF card slot.	Message Display

4 Confirm if OK to write the data to the CF card.

Press the YES key if you are sure to overwrite the CF card data with the monitor data.

Write setup data to CF card,	Write patient data
OK?	OK?
YES NO	YES [
<setup data=""></setup>	<patient< td=""></patient<>

Data Transfer (CF Card \rightarrow DS-7200)

The data can be transferred from the CF card to the monitor.

- **1** Insert the IC card to IC card slot.
- **2** Read the data from the CF card.



The data will be transferred from the CF card to the monitor. Press the CF Card \rightarrow DS-7200 key to select the data type to transfer.

NO

Data>

3 Select the data type.



4 Confirm if OK to read the data from the CF card.

Press the YES key if you are sure to overwrite the monitor data with the CF card data.

					_
Read setup data from CF card.			Read patient da	ita from CF card.	
OK?			c	K?	
YES	NO		YES	NO	
*If BP or TEMP unit of measurement is different from current setup, all the corresponded trend data and recall data will be initialized. *After reading setup data, power the unit again.			*If BP or TEMP unit of measur current setup, all the corre recall data will be initialized.		
<setup data=""></setup>		1	<patier< td=""><td>nt Data></td><td>-</td></patier<>	nt Data>	-
NOTE	When the data retrieving process completes, the display will return to the home display.				
▲ CAUTION	 Restart the system after reading the setup data from the CF card. The setup data will become effective after the system is restarted. Reading the patient data from the CF card will erase all previous patient data stored in the patient monitor. 				

Data for Transfer

The setup data such as monitoring condition, alarm setup, and patient data such as graphic trend data and tabular trend data can be transferred.

Setup Data

	Data	Description	
Parameter Setup		Stores the monitoring condition (size, lead, etc.)	
		for all the monitoring parameters.	
Alarm Setup	_	Stores the alarm threshold level.	
	ST Meas. Condition		
	Record		
	Sweep Speed		
	Tone / Volume		
	Color / Brightness		
	Display Configuration		
	Night Mode Setup		
System	Graphic Trend Setup		
Configuration	Tabular Trend Setup	Stores the current setup.	
5	Resp. List Setup	7	
	Vigilance List Setup		
	Recall Setup		
	ST Graphic Trend Setup		
	ST Display Lead Setup		
	Set Other Alarm		
	BP User Label		
	TEMP User Label		
	Alarm Mode Setup		
	Display Mode Setup		
Pre-Set	Hospital Setup	Stores the current setup.	
	(Not all data)		
	Monitor Setup		
	(Not all data)		

Patient Data

Data	Description
Patient Information	Stores patient information such as name, ID, age, sex, pacemaker use, patient type.
Graphic Trend Data	Stores 24 hours of graphic trend data.
Tabular Trend Data	Stores 24 hours of tabular trend data.

Error Message

CF card slot : OUT

Cause : CF card is not inserted or not correctly set in the CF card slot. Solution : Set the CF card into the CF card slot.

Invalid CF card.

Cause : Unspecified CF card is used.

Solution : Set the specified CF card into the CF card slot.

No data on the CF card.

Cause : There is no data on the CF card to be read out.

Solution : Check if the correct CF card is being used, or rewrite the data on the CF card.

CF card error.

Cause 1 : An error has been detected when writing/reading data on the CF card.

- Solution : If the error has been detected during writing, try again.
 - If the error has been detected during reading, data might not be correctly written on the CF card. Rewrite the data after formatting and try the procedure again.
- Cause 2 : The software version of the DS-7200 is older than that of the data stored in the CF card.
- Solution : The data of newer software version cannot be read. Update the software version of the DS-7200.
- Cause 3 : There is no more space on the CF card to write the data.
- Solution : Format the CF card.

Blank Page

FUKUDA DENSHI CO., LTD.

3-39-4 Hongo, Bunkyo-ku, Tokyo, Japan Tel: +81-3-5684-1455 Fax: +81-3-3814-1222 http://www.fukuda.com

Printed in Japan 4L010029I 201305