DYNASCOPE 8000 Series Patient Monitor

DS-8400 system

Ver. 05

Operation Manual



- Before using the product, please read this manual thoroughly.
- * Store this manual where it can be always referred to.



This manual is for the DS-8400 System Version 05.



MONITORING EQUIPMENT

AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY ANSI/AAMI ES 60601-1 (2005)+AMD1 (2012), IEC 60601-2-25 (2011), IEC 60601-2-27 (2011), IEC 80601-2-30 (2009)+AMD1 (2013), IEC 60601-2-34 (2011), IEC 60601-2-49 (2011), ISO 80601-2-55 (2011), ISO 80601-2-56 (2009) AND ISO 80601-2-61 (2011).

Federal Law restricts this device to sale by or on the order of a physician.

CAUTION

- Only physician or persons instructed by physicians are allowed to use the equipment.
- The information contained in this document is subject to change without notice due to improvement in the equipment.

© 2017 Fukuda Denshi Co., Ltd.

No part of this document may be reproduced or transmitted in any form without the prior written permission of Fukuda Denshi Co., Ltd.

If this manual has pages missing or out of order, contact Fukuda Denshi for replacement.

Contents

Preface

Introduction	i
Important Notice	
For Safe Operation of the Equipment	
Intended Use of this Equipment	
Copyright	i
Maintenance, Repair, Replacement	
Contact	i
About This Manual	iii
Expression Used in This Manual	iii
Composition of This Manual	
Composition of This Manual	iv

Safety

The Meaning of Each Safety Precaution	
Warping Labels Attached to the Unit	
Warning Labels Attached to the Unit	
Graphic Symbols	
Precautions for Safe Operation of Medical Electrical Equipment	
Precautions about the Maintenance	V
Precautions about the Network System	vi
Medical Telemetry	
Precautions when Using with Other Equipment	
Pacemaker	
Non-Explosion Proof	
Defibrillator	
Electrosurgical Instrument MRI (Magnetic Resonance Imaging)	
Precautions about Connections to Peripheral Devices	
Precautions for Using the Equipment	
This System	
Wired Network (DS-LANII/ DS-LANIII)	
Wireless Network System	
RTC and Data Backup	
Precautions about the Ventilator Monitoringx	
Precautions about the SpO2 Sensorx	
Precautions about the Masimo Modelx	
Precautions about the NIBP Cuffx	
Precautions about Disposing of the Equipment, Accessories, or Components	xiii
Precautions about Transportationxx	xiii
Monitoring after Power Failurexx	xiii
To Prepare for Emergency Usexx	
Electromagnetic Compatibilityxx	
Precautions for Safe Operation under Electromagnetic Influence	

Chapter 1 General Description

Composition of the System	1-1
Features	1-3
Menu Configurations	1-5

Chapter 2 Name of Parts and Their Functions

Name of Parts and Their Functions	
Main Unit: DSC-8400 series	2-1
Display Unit: LC-8018TC/LC-8016TC	
Super Unit: HS-8000 Series	
HS Adapter: HSA-80	
HS Adapter for DS-8400: HSA-81	
Transport Monitor: DS-8007 Series	
Adapter for DS-8007: DSA-82	2-8
Multi Module: HM-800	
HM-801 Multi Module	2-10
SpO2 Module (HG-810/HG-820)	
Multiport Module: HP-800	2-11
Recorder Unit: HR-800	2-12
Input Box: IB-8004	2-13
CO2 Gas Unit: HCP-800	2-13
CO2 Gas Unit: HCP-810	
CO2 Gas Unit: HCP-820	2-15
Gas Unit I/F: HPD-800	
Gas Unit I/F: HPD-810	
Gas Unit I/F: HPD-820	
BISx I/F Unit: HBX-800	2-17
Multigas Unit: MGU-800/810 Series	2-18

Chapter 3 Operation Procedure and Screen Examples

Operation Procedure	
Fixed Keys	
Touch Key	
Mouse	
Home Display	
About the Home Display	
Displayed Items	
Description of the Display	
Messages and Sound	
Window Display	
About the Window Display	
Display	

Floating Window Screen Display	
Minimize Window	
Transfer Window	
Operation Restriction	
Procedure to Return the Display	
To Enter Characters	
For Easier Use	
User Key	
Menu Screen	
To Delete the Unnecessary Keys (Key Mask)	
Display on the External Monitor and Extended Display Unit	
External Monitor Display	
Extended Display 1	
Extended Display 2	

Chapter 4 Preparation

Daily Check	4-1
To Start Monitoring	4-1
Check Discharge When Start Monitoring a New Patient	4-2
Data Transfer Function Using the Super Unit	
Data Transfer Function Using the Transport Monitor (DS-8007)	4-4
To Stop Monitoring	4-10
Clock Setup	4-11
Installing the Recording Paper	4-12

Chapter 5 Admit/Discharge

To Display the "Admit/Discharge" Screen	5-1
Admit	5-1
Entering the Patient Information	5-1
Entering Patient Information from the Magnetic Card	
Entering Patient Information from the Patient Data Server (When DS-LANIII is used)	
Discharge	
Discharging Procedure	5-6
User Mode	5-8
To Select the User Mode	5-8
Suspend Monitoring	5-9
To Suspend Monitoring	5-9
To Resume Monitoring	

Chapter 6 Alarm Function

Alarm	6-1
To Set the Arrhythmia Alarm	6-1
SpO2 Second Álarm Setup	
ST Alarm Setup	6-4
List of Alarm Settings	6-5
Detail Setup	6-6
Alarm Limit Setup	6-7
To Set the System Alarm (ON or Suspend)	6-7
To Silence or Suspend the System Alarm Sound	6-8
Alarm Limit Setup for Each Parameter	
About the Alarm Threshold Limit	6-11
Alarm Assist Screen	6-13

Chapter 7 Monitoring

To Display the Parameter Setup Screen	7-1
ECG	7-1
Before Attaching the Electrodes	
Electrode Placement	
Type of Electrodes and Lead Cable	
Connection to the Patient Monitor	
ECG Parameter Setup	
Respiration	
Respiration Monitoring (Impedance Method)	
RESP Parameter Setup	
BP	7-19
BP Monitoring	
Zero Balance of All Pressure Lines (User Key)	
Zero Balance of All Pressure Lines ([BP Zero] Key)	
Zero Balance for Each Pressure Line	
BP Parameter Setup	
Pulse Oximetry	
SpO2 Monitoring	
SpCO, SpMet, SpHb, SpOC Measurement (Masimo)	
Precautions about the Masimo Sensors and Cables	
SpO2 Parameter Setup (Nellcor)	
RR_SpO2 Parameter Setup (Nellcor)	
SpO2 Parameter Setup (Masimo) Non-Invasive Blood Pressure	
Lineup of Cuffs.	
NIBP Monitoring	
Inflation Mode Setup NIBP Auto Mode Setup	
Oscillation Graph Display	
Dyna Alert Function Status	
NIBP Parameter Setup	
Temperature	
, en per acare	

TEMP Monitoring	
TEMP Parameter Setup	
Cardiac Output and Blood Temperature	
Connecting the Super Unit	
Cardiac Output Measurement Algorithm	
Blood Temperature Alarm Setup	
CO2 Concentration (Mainstream Method)	
Patient Application and Display	
CO2 Parameter Setup	
CO2 Concentration (Sidestream Method)	
Patient Application and Display	
CO2 Parameter Setup	
Multigas Unit/SPIRO	
Connecting to the Respiration Circuit	
CO2 Measurement Unit Setup GAS Display during Undetected Breath	
Multigas Unit Data Setup (Multigas Concentration/Spirometry)	
BIS Data (HBX-800 with BISx)	
Preparation for Monitoring	
BIS Setup	
BIS Data (A-2000/A-3000)	
Ventilator	
AWP/AWF/AWV Scale Setup	
P-V/F-V Loop Display	
FLOW-i Data	
CO2 Measurement Unit Setup	
FLOW-i Setup	
SvO2/CCO Data	
STAT Mode / Index Display	
INVOS Data	
Channel Number Setup for INVOS Data	
Stopwatch	
Label Setup	
Start/Stop	
Multiparameter Connector Setup for BP, TEMP, CO Measurement	
Multiparameter Connector Setup	

Chapter 8 Review Function

Common Operation	
Arrhythmia Analysis	
Arrhythmia Definition	
Arrhythmia Alarm Setup	
Arrhythmia Learn	
Graphic Trend	
Graphic Trend Setup	
Description for Each Parameter	
Short Trend	8-10

Tabular Trend	
To Display/Print the Tabular Trend	8-13
The Description of the Display	
Parameter Setup for Tabular Trend	
Recall	
To Display the Recall Waveform	
To Display/Print the Enlarged Recall Waveform	8-18
Recall Setup	
OCRG	
Alarm History	
Alarm History Setup	
Description for Each Item	
Zoom Wave	
ST Measurement	
To Display/Print the ST Measurement	
Reference Waveform Setup	
ST Alarm Setup	
12-Lead Analysis	
12-Lead ECG Display	
12-Lead Analysis Setup	
12-Lead ECG Analysis	
12-Lead Analyzed Result Output Example	
Full Disclosure Waveform (Optional Function)	
To Format the CF Card	
Waveform Setup	
Description of the Full Disclosure Waveform Display	
To Search by Time	
Hemodynamics	
Calculation Data	
To Display/Print the Hemodynamics Data	
New Input of Hemodynamics Calculation	
To Edit the Hemodynamics Input Data	
Lung Function	
Calculation Data	
To Display/Print the Lung Function Data	
New Input of Lung Function Calculation	
To Edit the Lung Function Input Data	
Cardiac Output (CO)	
To Display the CO Measurement Screen	
Cardiac Output Setup	
CO Measurement	
To Edit the CO Measurement Result	
Drug Calculation	
Other Bed Display	
Other Bed Display/Alarm	8-54

Chapter 9 Printing

inting Setup9-1

Manual Printing (Basic)	. 9-1
Manual Printing (12-Lead)	
Manual Printing (Other Setup)	
Automatic Printing (Alarm Printing)	
Automatic Printing (Periodic Printing)	
Common Setup for Printing	
Freeze Printing	9-10
12-lead Waveform Printing	

Chapter 10 System Configuration

Display Configuration	
Numeric Data Selection	
To Configure the Display	
Waveform Selection	
User Key Selection	
Changing the Display Layout from the Home Display	
Tone/Volume	10-20
Color	10-21
Brightness	10-23
Night Mode	
Night Mode	10-25

Chapter 11 Troubleshooting

Message List	11-1
Vital Alarm Message	11-1
Vital Alarm Message (DS-LAN Standard Setup)	11-4
Equipment Status Alarm Message	11-8
Numeric Data Box Message	11-14
Ventilator Alarm Message	11-20
Ventilator Alarm Factor	11-20
Cardiac Output Message	11-21
Troubleshooting	11-22
ECG	
Respiration	
Invasive Blood Pressure	
SpO2 Measurement (HS-8312N, DS-8007N, HG-820)	11-31
SpO2 Measurement (HS-8312M, DS-8007M, HG-810)	11-33
Non-Invasive Blood Pressure	11-37
Temperature	
Cardiac Output (CO)	11-42
CO2 Measurement (HCP-800/HCP-810/HCP-820)	11-45
CO2 Measurement (HPD-800/HPD-810/HPD-820)	11-46
Recorder Unit (HR-800)	11-48
Network Printer	11-49
Wired Network (DS-LANII/ DS-LANIII)	

Telemeter (HLX-801)	
Remote Control	
General	
Super Unit	
Data Transfer Function	
IB-8004 Input Box	
Expansion Module	
Extended Display Unit	
Ventilator	
Multigas Unit	
SvO2/CCO Monitor	
BIS Monitor (A-2000/A-3000)	
BIS (When HBX-800 is used)	
INVOS	
FLOW-i	
PC Communication	11-73
TCM4/TCM5 FLEX	11-74
Magnetic Card Reader/Barcode Reader	
External Media	11-75

Chapter 12 Setup Item/Default Value

Chapter 13 Accessories

Accessories	
Optional Accessories	13-1
ECG, Impedance Respiration Measurement (Manufactured by Fukuda Denshi)	13-1
Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)	
Non-Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)	13-2
Temperature Measurement (Manufactured by Fukuda Denshi)	13-3
Pulse Oximetry Measurement (Manufactured by Covidien)	13-3
Pulse Oximetry Measurement (Manufactured by Masimo)	
CO Measurement (Manufactured by Fukuda Denshi)	
CO2 Concentration Measurement (Manufactured by Philips)	
CO2 Concentration Measurement (Manufactured by Covidien)	
Anesthetic Gas Concentration Measurement (Manufactured by Mindray Medical Sweden	AB)13-7
BIS Measurement (Manufactured by Covidien)	
Others (Manufactured by Fukuda Denshi)	
· · ·	

Chapter 14 Specification

Specification	14-1
Main Unit: DSC-8400 Series	14-1
Transport Monitor: DS-8007 and DSA-82	14-2
Display Unit: LC-8018TC/LC-8016TC	14-3
Super Unit: HS-8312N/8312M and HSA-80/HSA-81	14-4
Expansion Unit: MGU-800/810 Series and HR-800	14-5
Expansion Module: HM-800/801, HP-800, HG-810/820	14-6
Gas Unit I/F: HPD-800/HPD-810 and CO2 Gas Unit: HCP-800/HCP-810	14-7
Gas Unit I/F: HPD-820 and CO2 Gas Unit: HCP-820	
Input Box: IB-8004	14-9
BISx I/F Unit: HBX-800	4-10
Performance	4-11
Measurement Unit for Each Parameter	4-21
Alarm Limit Range for Each Parameter14	4-26
About the SpO2 Clinical Test	4-28

Contents

Preface

Introductioni
Important Noticei
For Safe Operation of the Equipmenti
Intended Use of this Equipmenti
Copyright ii
Maintenance, Repair, Replacement ii
Contact ii
About This Manualiii
Expression Used in This Manualiii
Composition of This Manual iv

Preface

Introduction

Thank you for purchasing this product. Read the "Safety Precautions" thoroughly before use to ensure correct and safe use of the product.

Before using or installing this product, read this manual thoroughly.

Important Notice

For Safe Operation of the Equipment

- (1) Before using this equipment, read this operation manual.
- (2) Fukuda Denshi cannot predict all the dangers which may be caused by misusage of this product or environmental condition.
- (3) For using this equipment, there are many items that "should be performed", "should not be performed", and "cannot be performed". It is not possible to cover all these items in this manual or warning labels. Therefore, it is necessary to also follow the general safety precaution other than the items described in this manual.
- (4) To prevent accidents, usage other than intended, or usage, cleaning, and maintenance not described in this manual should not be performed.
- (5) When using this equipment, follow the respective regulation to minimize the probability of accidents.

Intended Use of this Equipment

This equipment is designed for the following <Intended Use>.

<Intended Use>

This equipment is intended for measuring parameters such as ECG, respiration, NIBP, pulse rate, SpO₂, SpCO, SpMet, pulse wave, temperature, BP, CO, respiration gas (concentration of CO_2 , N₂O, volatile anesthetic agent, O₂), spirometry, BIS, brain wave, and monitors patient condition by displaying/printing the measurement data on this equipment or central monitor and generates alarm as required.

This equipment is intended for monitoring one patient. It is not intended for monitoring multiple patients. The 12-lead ECG analysis function is intended for adult and pediatric patients.

For specification of this equipment, refer to "Chapter 14 Specification" of this manual.

This equipment is intended to be used by healthcare professionals. Users should have a thorough knowledge of the function and operation before using this equipment. The maintenance of this equipment should be performed by skilled personnel who received a training of possible hazards and measures to avoid those hazards. Also, your local regulation must be followed. If this equipment is used for the purpose other than intended, or if the user does not follow the safety instructions, the following hazard may result.

- Hazard to the Life and Health of the Patient or the User
- A Problem Related to Medical Practice

• Damage to the Equipment

Copyright

- (1) The copyright of this manual is owned by Fukuda Denshi. No part of this document may be copied or transmitted in any form without the prior written permission of Fukuda Denshi Co., Ltd.
- (2) This manual includes the description for the optional equipments that can be connected.
- (3) The illustration in this manual may differ with the actual equipment.
- (4) If you lose or damage this manual, contact your nearest sales representative. Using the equipment without this manual may cause accidents.
- (5) When handing over this equipment, make sure to also pass this manual to the next owner.

Maintenance, Repair, Replacement

Fukuda Denshi is liable for the safety, reliability, and performance of its equipment only if;

- Maintenance, modifications, and repairs are carried out by authorized personnel or organization.
- Components are used in accordance with Fukuda Denshi operating instructions.

A full technical description of the DS-8400 System is available from your local Fukuda Denshi sales representative.

Contact

If you need more detailed information, please contact following.

(1) Fukuda Denshi Co., Ltd., Head Office

3-39-4 Hongo, Bunkyo-ku, Tokyo, Japan Phone:+81-3-5684-1455 Fax:+81-3-3814-1222 E-mail: info@fukuda.co.jp Home Page: http://www.fukuda.com

(2) Fukuda Denshi USA, Inc.

17725-C NE 65th Street Redmond, WA 98052 USA Toll Free: +1-800-365-6668 Local: +1-425-881-7737 Fax: +1-425-869-2018 Home Page: http://www.fukuda.com/index_usa.html

About This Manual

Expression Used in This Manual

Meaning of the Symbols

Type of Precaution	Description
A DANGER	Failure to follow this message may cause immediate threat of death or serious injury.
	Failure to follow this message may result in death or serious injury.
▲ CAUTION	Failure to follow this message may cause injury or failure to the equipment.
NOTE	"Note" is used to emphasize important information.
REFERENCE	"Reference" is used to provide useful information.
Ē	Indicates the reference page for the procedure and precaution.
*	Used in a table which indicates that there is detailed explanation outside the table.

□Indications for the Screens and Keys

The keys displayed on the monitor screen are indicated by []. (Ex.: [Display Config.], [Manual Printing], etc.)

The expressions displayed on the monitor screen are indicated by " ". (Ex.: "Volume", "Admit/Discharge", etc.)

The messages displayed on the screen are indicated by < >. (Ex: <Searching>, <Alarm Suspend>, etc.)

Composition of This Manual

Chapter Title	Description
Preface	Outline and purpose of this manual (Important Notice, About This Manual)
Safety	Warning, Precautions for Safety
1.General Description	Composition, features, menu configuration of this equipment
2.Name of Parts and Their Functions	Name and function of each part, external appearance
3.Operation Procedure and Screen Examples	Operation procedure, home display, window, procedure to return to the previous display, user key setup
4.Preparation	Installing the recording paper, power ON/OFF, time/date, daily checks
5.Admit/Discharge	Entering patient information (name, age, etc.) at admittance, discharging the patient, user mode selection, suspend monitoring
6.Alarm Function	General description of alarm function, alarm-related setups
7.Monitoring	Measurement condition setup of the monitoring parameters, size/scale setup, etc. Setup of the stop watch, connector
8.Review Function	Arrhythmia analysis, 12-lead analysis, trend, recall, NIBP list, ST measurement, hemodynamics, lung function, alarm history, other bed display, full disclosure waveform
9. Printing	Recorder output function
10. System Configuration	Setup of the display configuration, tone/volume, color, etc.
11. Troubleshooting	Message list, maintenance and troubleshooting of this equipment
12. Setup Item/Default Value	Setup details and default value
13. Accessories	List of accessories and optional accessories of this equipment
14. Specification	Specification and performance of this equipment

The operation manual is composed of the following chapters.

The maintenance manual is composed of the following chapters.

Chapter Title	Description
Preface	Outline and purpose of this manual (Important Notice, About This Manual)
Safety	Warning, Precautions for Safety
1.Installation of the Unit	Precautions about the operating environment, system construction, mouse connection
2.Network System Construction	Network connection and setup
3.Using the External Media	Procedure to use the external media
4.Connection to the External Devices	External equipment connection/setup, magnetic card reader usage
5.Initial Settings	Initial setup, administrator setup, alarm/measurement setup, user I/F, user mode registration
6.Setup Item/Default Value	Default and backup of setup items
7.Replacement Parts	Precautions about the periodic replacement parts, consumable parts
8.Cleaning/Disinfecting/Storing	Procedure to handle, clean, store this equipment
9. Maintenance Check	Daily and periodic checks, self-diagnosis function, software version software install

Safety

Safety

About the Safety Precautions

The Meaning of Each Safety Precaution

Read this manual thoroughly before use to ensure correct and safe use of the product.

Be sure to follow the precautions indicated below, as these are important messages related to safety.

Type of Precaution	Description
	Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.
	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
	Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Warning Labels Attached to the Unit

Make sure to read the warning label attached to the equipment and comply with the requirements while operating the equipment.

 Do not damage or erase the warning label attached to the equipment. This warning label contains important descriptions for handling and operating the equipment properly and safely. A damaged label may compromise safe operation.

DS-8400 System Main Unit



Warning Label Attached to the Equipment



Contents of the Label





Warning Label Attached to the Equipment (HS-8312)



Contents of the Label

Graphic Symbols

Symbol	Description
	Follow operating instructions (Warning); indicated in blue. Failure to follow operating instructions could place the patient or operator at risk.
[]ii	Follow operating instructions (Information); Indicates the need to refer to the related accompanying documents before operation.
	General precaution
	Caution, refer to accompanying documents Indicates the need to refer to the related accompanying documents before operation.
Ą	Potential Equalization Terminal Indicates the terminal to equalize the potential difference when interconnecting the devices.
÷	Protective Earth Indicates the protective earth inside the equipment.
~	Alternating Current (Main Power Input Indicator)
\odot	Indicates that the equipment is in normal operation.
Ó	Indicates that the equipment is in standby mode.
A LA	Electrostatic Sensitive Part Directly touching this connector part with hands should be avoided.
i 🗨 i	Type CF Applied Part with Defibrillation-Proof Indicates the degree of protection against electric shock is Type CF Applied Part with defibrillation- proof.
۱ <u>۴</u> ۱	Type BF Applied Part with Defibrillation-Proof Indicates the degree of protection against electric shock is Type BF Applied Part with defibrillation- proof.
⇔	Signal Output
G	GAS Input Part
⇒	GAS Output Part
÷	Signal Input Part
공문공	TCP/IP Network Connector Connects to TCP/IP network.
	RS-232C Connector Connects the related device.
	Eject Indicates the switch to pull out the paper tray.
\otimes	Indicates prohibited actions. Refer to the instruction.
0	Indicates mandatory or instructed actions. Refer to the instruction.
	Battery
	Date of Manufacture Indicates the date of manufacture.
	Name and Address of Manufacturer Indicates the name and address of manufacturer.
	WEEE (Waste Electrical and Electronics Equipment) Indicates a separate collection for electrical and electronic equipment.

Refer to the following for the meaning of the symbol indicated on the equipment.

Symbol	Description
IP32	Dustproof (IP3X): Protection against tips of tools. Waterproof (IPX2): Protection against water drops falling vertically over 15 degrees range.
IPX1	Waterproof (IPX1): Protection against water drops.
×	Alarm Silence
\$ \$ \$	NIBP Start/Stop
Se	NIBP Periodic Measurement
Ĩ	Lock
1	Unlock

Precautions for Safe Operation of Medical Electrical Equipment

- Users should have a thorough knowledge of the operation before using this equipment.
- Do not use the equipment in an environment where protective earth and wiring is questionable.

Precautions about the Location of Installation and Storage of the Equipment

- Set the monitor to the user's intended position where the user can easily recognize the visual and audible monitoring conditions. Normally it is recommended to set at a distance of one (1) m from the user.
- Install or store in a place where the equipment will not be exposed to splashing water.
- Install or store in an area where environmental conditions such as atmospheric pressure, temperature, humidity, ventilation, sunlight, dust, sodium, and sulfur will not adversely affect the system.
- Place the equipment on a stable surface where there is no inclination, vibration, or shock (including during transportation).
- Do not install or store in an area where chemicals are stored or gases are evolved.
- Verify the power frequency, voltage and allowable current (or power consumption).
- Ensure the grounding is proper by connecting the accompanying power cable to the hospital grade outlet.
- Make sure to secure the equipment during usage.
- Do not place the equipment or accessories in any position that might cause it to fall on the patient.

Precautions Before Using the Equipment

- Verify the power voltage. Charge the battery pack fully before operating the system with the battery pack.
- Check the cable connection and polarity to ensure proper operation of the equipment.
- Make sure the power system has adequate earth ground.
- Ensure that all cables are firmly and safely connected.
- Pay special attention when the equipment is used in conjunction with other equipment as it may cause erroneous diagnosis and danger.

Precautions During Using the Equipment

- Always observe the equipment and patient to ensure safe operation of the equipment.
- If any abnormality is found on the equipment or with the patient, take appropriate measures under the safe conditions, such as ceasing operation of the equipment.
- Do not allow the patient to come in contact with the equipment.
- On start-up of the system, verify that the start-up tone generates and alarm indicator lights.
- For the connectors which are not Type BF, CF applied part, do not touch them and the patient at the same time.

Precautions After Using the Equipment

- Unplug all the cables from the patient before turning off the power.
- When unplugging the cables, do not apply excessive force by pulling on the cord. Pull by the connector part of the cable.
- Clean the accessories and cables, and keep them together in one place.
- Keep the equipment clean to ensure proper operation for the next use.

Precaution when Equipment Failure Occurs

• If the equipment is damaged and in need of repair, the user should not attempt service. Label the unit "OUT OF ORDER" and contact your nearest service representative.

Precaution about Disassembling/Remodeling the Equipment

- Do not disassemble or remodel the equipment.
- If water or other liquids enter the equipment, cease using the equipment and contact your nearest service representative.

Precautions about Maintenance Check

- Make sure to periodically check the equipment, accessories, and cables.
- Before reusing the equipment that has been left unused for a while, make sure that the equipment operates normally and safely.

Precautions when Using with Other Equipment

• To prevent patient from burn injury, verify proper attachment of patient ground plate, ECG electrode type when using the electrosurgical knife, and verify paste volume, output energy when using the defibrillator. Also, verify that each equipment is properly grounded.

Precautions about the Maintenance

WARNING

• Never open the housing while the equipment is in operation or connected to hospital grade outlet as it may result in electric shock.

CAUTION Precautions about Safety Check

- For safe operation of the equipment, regular inspection and maintenance are required. Once a year, check all cables, devices, and accessories for damage, earth impedance, earth and leakage currents, and all alarm functions. Also, ensure that all safety labels are legible. Maintain a record of these safety inspections.
- Immediate maintenance has to be carried out for the following case.
 - When the equipment was subjected to extreme mechanical stress, e.g. after a heavy fall.
 - When the equipment was subjected to liquid spill.

- When the monitoring function is interrupted or disturbed.
- When parts of the equipment enclosure are cracked, removed, or lost.
- When any connector or cable shows signs of deterioration.

Precautions about the Network System

Medical Telemetry

CAUTION Precautions about the Installation

- The medical institution (hereinafter referred as "Institution") must decide the telemetry installation plan for the medical institution in order to prevent interference between transmitters (telemetry based on destination country's radio law). When telemetry has already been installed and been used, radio format, frequency, and antenna power are required to be examined to prevent interference.
- When using telemetry which requires zone location, the institution is to set up the zones as an operation unit for each transmitter to prevent electronic interference between telemetry throughout the Institution.
- When using telemetry which requires zone location, display and identify each prepared zone in the equipment.
- When laying receiver antenna for each transmitter, the Institution has to examine the installation so that electronic interference does not occur.
- Based on the above examination result, the Institution should place each receiver antenna as required.

CAUTION Precautions about the Management

- The institution appoints a person to manage the wireless channels for the whole medical institution. And when using telemetry which requires zone location, the Institution should nominate a person to manage the wireless channels in each zone (a "Coordinator"). However, when using such telemetry in a local medical institution, one person can perform both functions.
- Select a telemetry coordinator who understands the characteristics and functionality of telemetry systems, and is skilled in operating telemetry.
- When installing telemetry, the Coordinators have to understand the precautions for use of the telemetry in advance.
- The Coordinator takes responsibility of wireless channel management and transmitter storage for the whole Institution by giving proper instruction.
- The Coordinator should create a management log (hereinafter referred to as the "log"), which contains a list of the management status of the wireless channels for the whole Institution. When changing a wireless channel, register it in the log and give proper instructions to the user.
- The Coordinator assumes responsibility for managing the wireless channels, storing, and managing telemetry.
- The Coordinator assigns the transmitter to the user, and provides enough education for use inside the zone.
- The telemetry user verifies operation of the transmitter/receiver before use.
- The telemetry user, if using the telemetry in a zone location, follows the instructions of the Coordinator for the zone and gives instructions to the patient if required.
- When interference or breakdown occurs in telemetry communication, the user is required to inform the Coordinators of the problems. The Coordinators are to deal with the problem properly and/or contact their nearest Fukuda Denshi representative for service.

Precautions when Using with Other Equipment

Pacemaker

WARNING

- Minute ventilation rate-adaptive implantable pacemakers can occasionally interact with certain cardiac
 monitoring and diagnostic equipment, causing the pacemakers to pace at their maximum programmed rate. The
 cardiac monitoring and diagnostic equipment may possibly send wrong information. If such event occurs,
 please disconnect the cardiac monitoring and diagnostic equipment, or follow the procedures described in the
 operation manual of the pacemaker. For more details, contact FUKUDA DENSHI personnel, your institution's
 professionals, or your pacemaker distributors.
- Rate meters may continue to count the pacemaker rate during occurrences of cardiac arrest or some arrhythmias. Do not rely entirely upon rate meter alarms. Keep pacemaker patients under close surveillance.

Reference

"Minute Ventilation Rate-Adaptive Pacemakers"

FDA alerts health professionals that minute ventilation rate-adaptive implantable pacemakers can occasionally interact with certain cardiac monitoring and diagnostic equipment, causing pacemakers to pace at their maximum programmed rate.

[Based on a safety bulletin issued by FDA Center for Devices and Radiological Health on October 14, 1998]

Non-Explosion Proof

DANGER

• Never operate the equipment in the presence of flammable anesthetics, high concentration of oxygen, or inside hyperbaric chamber. Also, do not operate the equipment in an environment in which there is a risk of explosion. Explosion or fire may result.

Defibrillator

WARNING

When defibrillating, keep away from the electrodes or medicament applied to the patient chest. If this is not
possible, remove the electrodes or medicament before defibrillating.
If the defibrillator paddles are directly in contact with the electrodes or medicament, an electrical shock may

If the defibrillator paddles are directly in contact with the electrodes or medicament, an electrical shock may result by the discharged energy.

• When defibrillating, make sure that the electrodes, sensor cables, or relay cables are firmly connected to the device.

Contacting the metal part of the disconnected cable may result in electrical shock from the discharged energy.

- When defibrillating, do not touch the patient and the metal part of the device or cables. Electric shock may result from the discharged energy.
- This equipment will return to standard operating mode within 10 seconds after defibrillating. However, when in diagnosis mode, it may require 10 seconds or more after defibrillation to display the normal ECG waveform as the time constant setting is large.

The stored data will not be affected. The measurement accuracy will temporarily decrease during defibrillation, but it will not compromise the safety of patient and the equipment.

• The QRS synchronized signal is not intended to be used as synchronized signal for defibrillator.

Electrosurgical Instrument

WARNING

• The monitoring system contains protection against interference generated by electrosurgical instruments. However, depending on the operating conditions, surgery site with respect to the location of ECG electrodes, ground plate attachment condition, or the type of instrument used, it may cause burn injury at the electrode site or noise on the ECG. The noise is generated at the tip of the electrosurgical knife and is difficult to completely eliminate because of the frequency components of the ECG. To reduce electrosurgical interference, take the following precautions:

Location:

Locate the electrosurgical unit as far as possible from this equipment and the patient cable. This will help reduce interference on the ECG through the monitor or cables.

Power Supply:

Connect the electrosurgical unit to a power supply that is different from that of this equipment. This will help prevent interference through the power cable.

Electrode Placement

The amount of noise interference is considerably different depending on the electrode position and surgery site. Place the ECG electrodes as far away as possible from the surgery site and the ground plate. Do not place electrodes in the path between the surgery site and the ground plate. If the electrodes are placed in this path, the amount of interference will be quite large. Position (+) and (-) electrodes as close as possible to each other.

Ground Plate

When using electrosurgical instruments, make sure the contact between the patient and the ground plate is secure. If the connection is incomplete, the patient may suffer from burn at the electrode site.

- The stored data will not be affected. The measurement accuracy will temporarily decrease during electrosurgery, but it will not compromise the safety of patient and the equipment.
- When using the electrosurgery-proof type ECG relay cable, the impedance respiration cannot be measured, and its numeric data and waveform will not be displayed. When measuring in an environment where electrosurgery is not performed, make sure to use the standard ECG relay cable.
- As this equipment utilizes capacitive touch panel, the energy from the electrosurgical knife may pass through the cable to the touch panel causing unintentional touch panel control. Locate the cables as far away as possible from the touch panel.

MRI (Magnetic Resonance Imaging)

WARNING

MR Unsafe-Keep away from magnetic resonance imaging (MRI) equipment.

- Do not use this equipment in magnetic resonance imaging (MRI) environments.
- When conducting MRI test, remove the electrodes and sensors connected to the patient (test subject). This equipment may be pulled towards the MRI device. Also, the local heating caused by the induced electromotive force may cause burn injury to the patient or performance degradation, failure, damage of this equipment.

For details, refer to the operation manual for the MRI testing device.

Precautions about Connections to Peripheral Devices

To use the equipment safely and to ensure maximum performance of the equipment, connection of other manufacturer's equipment to this equipment is not authorized, unless the connection is explicitly approved by Fukuda Denshi. It is the user's responsibility to contact Fukuda Denshi to determine the compatibility and warranty status of any connection made to another manufacturer's equipment.

WARNING

- When multiple equipments are connected to the patient, it may be necessary to take measures for connection (use of separation device), power supply (use of isolation power), grounding (additional protective earth). If these measures are not properly taken, a leakage current may flow between the equipments, or the total amount of leakage current may exceed the limit specified on IEC 60601-1-1.
- Only the peripheral devices specified by Fukuda Denshi should be connected with the given procedure. Use of an unspecified device may cause electric shock to the patient and/or operator due to excessive leakage current.

- Although the peripheral device connectors on the DS-8400 System are, with some exceptions, isolated from the power supply, the connecting peripheral devices should comply with IEC 60601-1. It is the user's responsibility to verify that the overall system complies with IEC 60601-1-1.
- To prevent danger of electric shock, always position the peripheral devices away from the patient.
- Network equipment including printer and hub should be located outside the "Patient Environment". If located inside the "Patient Environment", it may result in electric shock to the patient or the operator.
- Combinations of medical equipment with non-medical equipment must comply with IEC 60601-1-1 or IEC 60601-1. Never use a multiple portable socket-outlet or extension cable when connecting the equipments unless it is supplied specifically for use with that equipment.

Precautions for Using the Equipment

This System

A DANGER

• When connecting to other equipments, contact your nearest representative. Danger such as electric shock may result to the patient and operator.

WARNING Warnings about the System

- Do not connect any equipment or cable not authorized by Fukuda Denshi to any I/O connector. Also, do not connect any damaged equipment or cable. If done so by mistake, not only that the equipment cannot deliver its maximum performance, the equipment may be damaged and safety cannot be ensured.
- If the equipment is used under an environment not fulfilling the specified condition, not only that the equipment cannot deliver its maximum performance, the equipment may be damaged and safety cannot be ensured. If using the equipment under condition other than specified, contact your nearest representative.
- Use only the supplied 3-way AC power cable. Use of other cables may result in electric shock to the patient and the operator.
- The power cable must be connected to a hospital grade outlet.
- When using multiple ME equipment simultaneously, perform equipotential grounding to prevent potential difference between the equipments. Even a small potential difference may result in electric shock to the patient and the operator.

- Carefully route the cables to avoid patient entanglement and strangulation.
- When lifting this equipment, hold the bottom part of the main unit and not the display unit.
- Damage to the LCD may cause leakage of liquid crystal. In such case, do not touch the leaked liquid crystal with bare hands, or put it into your mouth as it may cause intoxication. If the liquid crystal accidentally enters the eyes or mouth, wash off immediately with water and consult a physician.
- When attaching the display unit to the main unit, slowly insert the rail on the rear side of the display unit to the guides on both sides of the main unit. Then, secure it with the specified screws.

WARNING Warnings about the monitoring

- The patient classification selection influences the precision of the QRS detection and NIBP measurement. Make sure the proper selection is made.
- The pacemaker usage setting influences the precision of the QRS detection and arrhythmia analysis. Make sure the correct selection is made.
- If the QRS pace mask function is set to [OFF], [10ms]/[20ms], the pace pulse may be erroneously detected as a QRS complex and HR alarm or asystole alarm may not generate due to incorrect HR (counting pace pulse as QRS complex). Select [OFF], [10ms]/[20ms] only if you are sure that pacing failure will not occur, or when the patient can be constantly monitored.
- When measuring the SpO₂ of patient with high fever or peripheral circulatory insufficiency, check the sensor attachment periodically and change the attachment site. The temperature of the attachment site will rise due to the sensor heat which may result in burn injury.
- For the following case, accurate measurement of SpO2 may not be possible.
 - Patient with excessive abnormal hemoglobin (COHb, MetHb)
 - Patient with the pigment injected to the blood
 - Patient receiving CPR treatment
 - When a sensor is applied to a limb with NIBP cuff, arterial catheter, or intracatheter
 - When measuring at site with venous pulse
 - Patient with body motion
 - · Patient with small pulse
- When a patient is receiving a photodynamic therapy, measuring SpO₂ on a same site for a long duration may cause blisters from the irradiation light of the SpO₂ sensor. Make sure to periodically change the sensor attachment site.
- Before the measurement, make sure the patient classification (Adult/Child/Neonate) is properly selected. Otherwise, correct measurement cannot be performed, and congestion or other injury may result.
- When the system alarm is suspended, all the alarm will be suspended even if the parameter alarm is set to [ON]. Also, the alarms will not be stored as recall events.
- If the upper/lower alarm limit of the parameter is set to [OFF], or arrhythmia alarm is set to [OFF], alarm will not function even if the system alarm is set to [ON]. Pay attention when setting them [OFF].
- Objective and constant arrhythmia detection is possible through the fixed algorithm incorporated in this monitor. However, excessive waveform morphology change, motion artifact, or the inability to determine the waveform pattern may cause an error, or fail to make adequate detection. Therefore, physicians should make final decisions using manual printing, alarm printing and recall waveform for evaluation.
- The RR/APNEA alarm will not be generated unless the numeric data box corresponded to the selected RR/ APNEA alarm source is displayed. Make sure to display the numeric data box for the RR/APNEA alarm source.
- The SpO₂ respiration measurement function is not intended for use as an APNEA monitor.
- When selecting [0] for "Volume" or [Timer] for "Display" for the Night Mode, pay attention not to miss any important alarm by simultaneously monitoring the patient on central monitor or other monitors.
- When the alarm sound is suspended, the alarm sound will not generate for the fixed amount of time. Pay

Safety

attention not to miss any important alarm by simultaneously monitoring the patient on central monitor or other monitors.

- If the safety of the patient cannot be ensured, do not suspend the alarm or decrease the alarm volume.
- The oxygenator mode is intended to prevent alarms during cardiopulmonary bypass surgery. Pay special attention when using this mode as the alarm generation will not be the same as to the standard monitoring mode.
- If the "Alarm Setting" under the Oxygenator Mode Setup is set to [All OFF], all vital alarm will not generate regardless of the alarm setting of each parameter. Also, if [Sel. Parameter] is set, vital alarm for unselected parameter will not generate. Pay attention to not miss any significant change of the patient's vital sign as the alarms will not be generated during the Oxygenator Mode.
- Once the cardiopulmonary bypass is finished, make sure to cancel the Oxygenator Mode and return to the standard monitoring mode.

WARNING Warnings about the CO₂ Monitoring

(HCP-800/HCP-810/HCP-820, HPD-800/HPD-810/HPD-820)

- Only one of either HCP-800/HCP-810/HCP-820/HPD-800/HPD-810/HPD-820 can be connected.
- When using a sampling line for intubated patients with a closed suction system, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.
- To prevent cross-infection, do not allow the sampling gas to return to the breathing system.
- To protect the hospital staffs from unnecessary anesthetic agent when using the HCP-800/HCP-810/HCP-820, it is strongly recommended to connect the exhaust hole to the gas exhaust system in the hospital.
- Loose or damaged connections of the sampling line may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.
- Do not cut or remove any part of the sampling line. It could lead to erroneous readings.
- If too much moisture enters the sampling line (i.e., from ambient humidity or breathing of unusually humid air) when using the HCP-800/HCP-810/HCP-820, <Check Sample Line> will be displayed in the message area. Replace the sampling line once this message is displayed.
- Carefully route the sampling line to reduce the possibility of patient entanglement or strangulation.
- Do not lift the HCP-800/HCP-810/HCP-820 by the sampling line, as the sampling line could disconnect from the equipment, causing the equipment to fall on the patient.
- CO₂ readings and respiratory rate can be affected by sensor application errors, certain ambient environmental conditions, and certain patient conditions.

WARNING Warnings about the Gas Monitoring (MGU-800/810)

- Make sure to use only the specified Mindray Medical Sweden AB product.
 (Planesthetic Gas Concentration Measurement (Manufactured by Mindray Medical Sweden AB)" P13-7)
- Be careful not to damage the water trap during operation as bacteria and/or mucus may contaminate the MGU-800/810 series.
- The airway adapter, sampling line, flow sensor are disposable products that are intended for single patient use only. Do not reuse them on other patients as it may cause cross-infection.
- To prevent cross-infection, do not allow the sampling gas to return to the breathing system.
- When using a sampling line for intubated patients with a closed suction system, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.
- Do not use the MGU-800/810 series with the flammable anesthetic agents.
- To protect the hospital staffs from unnecessary anesthetic agent, it is strongly recommended to connect the exhaust hole to the gas exhaust system in the hospital.
- The sampling line may get clogged by internal condensation.

- - Do not use adult/pediatric type water traps and/or sampling lines with neonates to avoid high sampling flow.
 - Connect only DRYLINE gas sampling lines to the water trap. Note that there may be other compatible tubing present, e.g. IV-lines.
 - Do not use DRYLINE neonatal sampling lines (blue luer lock nuts) with DRYLINE adult water traps as this could result in incorrect measurement data.
 - Do not use DRYLINE adult sampling lines (colorless luer lock nuts) with DRYLINE neonatal water traps as this could result in incorrect measurement data.
 - Only combine the SPIRIT Flow Sensors and DRYLINE Water Traps as specified. Other combinations might lead to incorrect measurements.

(Connecting to the Respiration Circuit" P7-85)

- Use the adult flow sensor for a patient whose tidal volume is above 150 mL.
- Use the pediatric flow sensor for a patient whose tidal volume is below 300 mL.

• The contents of the water trap should be handled as a potential infection hazard.

- Make sure to use the correct flow sensor depending on the patient conditions, adult or pediatric and the tidal volume.
- Do not confuse the gas sampling line with other compatible tubing, e.g. IV-lines.

WARNING Warnings about the 12-Lead ECG Analysis Function

- The 12-Lead ECG analysis function is designed to acquire and interpret ECG data from a resting, supine patient. If ECG signals from moving or shaking patients are acquired, erroneous 12-lead interpretation result may occur. Always ensure that the patient is kept motionless during 12-lead ECG signal acquisition and analysis.
- The 12-lead ECG analysis function is intended for use with adult and pediatric patients.
- All computerized ECG analysis results should be reviewed by a physician before making decision for the patient treatment.

WARNING Warnings about the BIS Monitoring (HBX-800)

- Clinical judgment should always be used when interpreting BIS in conjunction with other available clinical signs. Reliance on BIS alone for intra-operative anesthetic management is not recommended.
- BIS values should also be interpreted cautiously with certain anesthetic combinations, such as those relying primarily on either ketamine or nitrous oxide/narcotics to produce unconsciousness.
- As with any monitored parameter, artifacts and poor signal quality may lead to inappropriate BIS values. Potential artifacts may be caused by poor skin contact (high impedance), muscle activity or rigidity, head and body motion, sustained eye movements, improper sensor placement and unusual or excessive electrical interference.

CAUTION Precautions about the System

- Do not assess the patient's condition only with the information from this equipment. A clinical judgment based on the information from the equipment should be made by a doctor who fully understands functions of the equipment, in a comprehensive manner combined with clinical findings and other test results.
- Do not assess the patient's condition only with the alarm from this equipment. When the alarm is set to OFF or if the alarm priority is low, a sudden change of the patient may not be noticed.
- If an alarm generates, check the patient's condition first and ensure the safety. Depending on the alarm, take appropriate measures to remove the problem. If the problem lies with the alarm setting, set the alarm properly.
- When measuring for a long period of time, make sure not to compress the patient with the lead cables and the electrodes. Compressing the same site for a long duration may inhibit the blood flow and generate compression necrosis and burn injury.
- Use only the spare parts specified for this equipment. Otherwise, proper function cannot be executed.
- Do not use the touch panel with the film attached. It may cause malfunction or damage the touch panel.
- For quality improvement, specifications are subject to change without prior notice.

- This equipment utilizes LED for the backlight. Since this LED deteriorates by the life cycle, the display may become dark, scintillate, or may not light by the long term use. In such case, contact your nearest service representative.
- This equipment is intended to be used for only one patient.
- The installation of this equipment should be performed by our service representative or a person who is well acquainted with this equipment.
- If not using the equipment for a long period, disconnect the power cable and lithium-ion battery.
- Although the LCD utilizes highly accurate picture elements, occasionally, there may be a few pixels which do not light or constantly light. Please note that this is not an equipment failure, and will not affect monitoring operation.
- Exposing LCD panel to intense light may deteriorate display property. Do not expose the panel to direct sunlight or strong ultra violet (UV) light.
- The lithium-ion battery can only be charged in the specified operational temperatures of the equipment. Refer to the operation manual of the lithium-ion battery (BTO-005) for details.

CAUTION Precautions about the ECG Monitoring

- If any electrodes get detached from the patient after being connected to the lead cable and patient monitor, pay attention that the metal part of the electrode does not get in touch with any metal parts of the bed or any conductive parts. Also, the operator should not touch any conductive parts with bare hands. Otherwise, it may cause electric shock to the patient and/or operator due to excessive leakage current.
- The indication for continuous use of the electrode is about one day.
- Replace the electrode if the skin contact gets loosen due to perspiration, etc.
- When an electrode is attached to the same location for a long period, some patients may develop skin irritation. Check the patient's skin condition periodically and change the electrode site as required.
- For stable arrhythmia detection and ECG monitoring, verify proper electrode placement, lead, waveform size, and filter mode selection. If not properly selected, it may cause erroneous detection.
- The threshold level for arrhythmia detection changes with ECG waveform size. Set a proper waveform size for monitoring.
 - When the ECG waveform size is x1/4, x1/2, or x1, the arrhythmia detection level is 250 μ V.
 - + When the ECG waveform size is x2 or x4, the arrhythmia detection level is 150 μ V.
- The leads for arrhythmia detection, central monitor display, printing are fixed as ECG1 and ECG2. Set the most appropriate leads with high QRS for ECG1 and ECG2, especially for arrhythmia detection. If the QRS amplitude for the set lead is low, it may cause erroneous arrhythmia detection.
- In ESIS Mode, artifacts such as electrosurgical noise or EMG can be largely reduced, but QRS amplitude attenuation, waveform distortion, or ST segment change may occur compared with other filter modes.
- The ESIS mode cannot completely reduce the electrical noise, and may erroneously detect the pacemaker spike. This mode should be selected only when a high frequency noise largely affects the HR measurement.
- There are some cases when the pacemaker pulse cannot be detected depending on the pacemaker type, pulse voltage, pulse width, electrode lead type (unipolar, bipolar), or electrode placement which causes the pacemaker pulse amplitude to decrease, and disables the pacemaker pulse detection.
- If signals similar to a pacemaker pulse are present, such as electric blanket noise or excessive AC frequency noise, these may be erroneously detected and displayed as a pacemaker pulse.
- When a spontaneous QRS and pacemaker pulse overlap (ex. fusion beat, etc.), QRS detection cannot be performed properly. In this case, the heart rate is degraded.
- If a pacemaker pulse is continuously detected due to AC frequency interference, QRS detection will be suspended and the heart rate will be reduced. Arrhythmia will not be detected either.

CAUTION Precautions about the ST Measurement

• The ST algorithm has been tested for accuracy of the ST segment data. The significance of the ST segment

changes need to be determined by a clinician.

• For the lead which the electrode is detached, the reference waveform setup cannot be performed. Check if the electrode is appropriately attached, and perform the setup again.

CAUTION Precautions about the 12-Lead Analysis	
---	--

• Interpretation and Minnesota codes given by this equipment do not instruct the physician as to the kind and degree of cardiac disease. Accordingly, four value judgments are given for ECG waveform and although "abnormal" indicates a large possibility of organic cardiac disease, there are cases where no cardiac disease exists despite an abnormal ECG (that is, an abnormal ECG may be caused by something other than heart). On the other hand, care should be taken in the event any preclinical coronary arteriosclerosis could be present despite a normal ECG interpretation.

Therefore, for a proper diagnosis, the ECG should be integrated with other interpretations.

• ECG Recording by the Mason-Likar System The 12-lead ECG recorded with the torso placement of the limb leads (Mason-Likar 12-lead system) may differ from the standard 12-lead ECG. Moreover, waveforms may differ somewhat also in a supine position and a standing position (sitting position).

Fukuda Denshi recommends to carry out the recording of the ECG by taking into consideration the waveform differences according to electrode positions or postures.

- For the model installed with ECG analysis program The ECG analysis program is intended to analyze the standard 12-lead ECG waveforms. Therefore, the analyzed result for waveforms recorded with the torso placement of the limb leads (Mason-Likar 12-lead system) may differ from that of the standard 12-lead ECG waveforms.
- Select "Used" for the pacemaker setting on the patient admit/discharge menu if a patient has a pacemaker.
- The threshold values for classification of 12-lead ECG interpretation and Minnesota code are set by age and sex as follows:
 - 1. Male and Female of ages 19 years old and above
 - 2. Male of age 12 through 18 years old
 - 3. Female of age 12 through 18 years old
 - 4. Male and Female of ages 3 through 11 years old
 - 5. Male and Female of ages below 2 years old
- If no patient information (i.e. Default : "Class." [Adult], "Sex": undetermined, and "Age" [0]) has been entered, the system algorithm will handle the patient as a "35 years old male".
- Before the analysis, make sure the patient classification ([Adult] / [Child]) is properly selected. If [Neonate] is selected, the 12-lead ECG analysis will not function.
- Enter the age of patient if known. If no age information (i.e. Default: [0]) has been entered, the system algorithm will handle the patient as "35 years old".
- Enter the sex of patient if known. If no sex information (i.e. Default: undetermined) has been entered, the system algorithm will handle the patient as "Male".
- If the patient classification is set as [Child] and no age (i.e. Default: [0]) has been entered, the system algorithm will handle the patient as "less than 2 years old""

CAUTION Precautions about the SpO₂ Monitoring

- Use only the sensor/relay cable specified by Fukuda Denshi. Otherwise, it may cause measurement error. If the sensor is damaged, stop using it.
- If the nail is rough, dirty, or manicured, accurate measurement will not be possible. Change the finger or clean the nail before attaching the sensor.
- If irritation such as skin reddening appears with the sensor use, change the attachment site or stop using the sensor.
- Do not apply the sensor too tight. At the same time, check the blood flow constantly so that congestion is not generated at the peripheral site.
- Do not use tape to attach the sensor.

- Even attachment for a short duration may inhibit the blood flow and generate compression necrosis or burn injury. Also, blood flow inhibition may prevent correct measurements.
- Check the sensor attachment site constantly in every 4 hours when probes or reusable sensor are used, and at least every 8 hours when single patient use sensors are used. Be especially careful of a patient with bad perfusion. If the sensor attachment position is not changed constantly, skin irritation or skin necrosis due to compression may be developed. For the patient with bad perfusion, check the sensor attachment position at least every 2 hours.
- As skin for neonate, premature infant is immature, change the sensor attachment site more frequently depending on the condition.
- Direct sunlight to the sensor area can cause a measurement error.Place a black or dark cloth over the sensor if using in direct sunlight.
- When not measuring, unplug the relay cable and sensor from the SpO₂ connector. Otherwise, the outside light may affect to falsely display measurements.
- The pulse wave is normalized for SpO₂ measurement, and does not indicate perfused blood volume. Check proper probe attachment by observing the pulse wave.
- Precautions for Reusable Sensors

The light-emitting part of the sensor should be over the root of the fingernail or as instructed per the related sensor instruction manual. Do not insert the finger too far into the sensor as it may hurt the patient. For details, refer to the SpO_2 sensor instruction manual.

• Precautions for Single-Patient-Use Type Sensors

The sensor can be reused on the same patient as long as the adhesive tape attaches without slippage. But do not reuse on other patients to avoid cross contamination. It is intended for single patient use only. For details, refer to the SpO_2 sensor instruction manual.

- If "---" is displayed for the SpO₂ numeric data, make sure that the sensor is properly attached.
- Measuring on a limb with NIBP cuff, arterial catheter, or intracatheter may result in incorrect measurement.
- Venous congestion may cause under reading of actual oxygen saturation. Therefore, assure proper venous outflow from monitored site. Sensor should not be below heart level (e.g. sensor on hand of a patient in a bed with arm dangling to the floor).

CAUTION Precautions about the NIBP Monitoring

- Do not apply the NIBP cuff to site of injury. An injury may be worsened by the measurement.
- Do not apply the NIBP cuff to the arm on side treated axillary lymph nodes dissection. It may lead to lymphatic edema by the cuff pressure.
- Measuring on a limb with SpO₂ sensor, arterial catheter, or intracatheter may result in incorrect measurement.
- An operator must not get away from a patient during the NIBP measurement. However, when getting away from the patient is necessary, do not activate the Alarm Suspend and Silence functions in order not to miss any sudden changes in the patient's condition.
- Pay attention when measuring the NIBP of patient with bleeding disorders or hyper coagulation. The cuff inflation may cause petechia or circulatory failure by the blood clot.
- For the following situation, measurements will be terminated.
 - When the measurement time has exceeded 160 seconds for adult and child, 80 seconds for neonate.
 - When the inflation value has exceeded 300mmHg for adult, 210mmHg for child, and 150mmHg for neonate.
- If used with the incorrect patient classification, it will not only cause erroneous measurement, but the inflating level for the adult may be applied to child or neonate causing dangerous situation to the patient.
- The continuous measurement and 1-minute interval measurement will automatically stop after 12 minutes (maximum 15 minutes).
- If the mean MAP display is set to OFF, the MAP alarm will not be generated. Also the MAP data will not be displayed for the tabular trend or the NIBP list.

CAUTION Precautions about the BP Monitoring

- Do not reuse / re-sterilize the disposable type transducers.
- If using a reusable blood pressure transducer, disinfect it according to the manufacturer's guidelines.
- The long-term use of the blood pressure transducer, tube and catheter may increase the risk of infection. Perform periodic replacement with new one. The guidelines of the CDC (Disease Control and Prevention) recommend replacing within 96 hours.
- If the ambient temperature of the blood pressure transducer has changed greatly, the zero balance may cause the drift. Perform the zero balance again.
- An operator must not get away from a patient during the BP measurement. However, when getting away from the patient is necessary, do not activate the Alarm Suspend and Silence functions in order not to miss any sudden changes in the patient's condition.
- Be sure to perform Daily Check. Use of faulty equipment might harm the patient or operator.
- If the Equipment Status Alarm occurs or if you feel the unusual operation of the equipment, perform the inspections to confirm the safety or contact our service representative.
- If the transducer get disconnected, pay attention that the metal part of the transducer does not get in touch with any metal parts of the bed or any conductive parts. Also, the operator should not touch the conductive parts with bare hands. Otherwise, it may result in electric shock to the patient and/or operator.
- When the power is turned ON, the BP value will not be displayed until zero balance is performed. Make sure to perform the zero balance. Once the zero balance is performed, the zero balance information will be maintained, and the BP value will be displayed.
- Each time the blood pressure transducer or tubing is replaced, the zero balance procedure is required to ensure accurate measurements.
- The zero balance procedure is required for the following case.
 - When starting the measurement.
 - When the position of the heart has changed due to body movement.
 - When the position of the transducer has changed.
 - When measuring for a long period of time and there is a possibility of measurement error due to change in ambient temperature, etc.
 - When a connector is connected/disconnected, or a transducer is replaced.
- Note that Systolic Pressure (SYS) = Peak Systolic Pressure (PSP) for the graphic trend, data base, and alarm setup.
- When ECG is not measured, Peak Diastolic Pressure (PDP) cannot be calculated.
- The undisplayed BP data (SYS/DIA/Mean) will not generate a BP alarm or be displayed in the tabular trend. Select the appropriate display type according to the monitoring purpose.

CAUTION Precautions about the CO₂ Monitoring (HCP-800/HCP-810/HCP-820)

• Conduct CO₂ calibration for the following case. If the CO₂ gas calibration is not performed at a specified interval, CO₂ measurement accuracy may be affected and also subsequent gas calibration may not be possible.

- When the accumulated measurement time exceeds 1,200 hours from the first use. However, if the first calibration was performed before the accumulated measurement time reaches 720 hours, another calibration is required when the accumulated measurement time exceeds 1,200 hours from the first calibration.
- When 12 months has elapsed or the accumulated measurement time has exceeded 4,000 hours from the previous calibration.
- When EtCO₂ measurement is not stable or accuracy is degraded compared with other measuring device.
- When the patient monitor was not used for a while, or when EtCO₂ was not measured for a while.
- Perform the calibration 5 minutes after turning ON the power on the HCP-800/HCP-810/HCP-820.
- Do not disconnect the sampling tube during calibration. If disconnected, calibration will cease.
- Dispose of calibration gas according to the regulation of each medical institution.
- Microstream[®] EtCO₂ sampling tubes are designed for single patient use, and are not to be reused. Do not attempt to clean, disinfect, sterilize or flush any part of the sampling tube as this can cause damage to the monitor or lead to cross-infection.
- Dispose of sampling tubes according to standard operating procedures or local regulations for the disposal of contaminated medical waste.
- Before use, carefully read the Directions for Use for the Microstream[®] EtCO₂ sampling tube.
- Only use Microstream[®] EtCO₂ sampling lines to ensure the monitor functions properly.

CAUTION Precautions about the CO₂Monitoring (HPD-800/HPD-810/HPD-820 Gas Unit I/F)

- The disposable airway adapter should be opened just before use.
- Do not reuse the disposable airway adapter. If sterilized, it will become unusable.
- The airway adapter should be attached with the thicker side facing to the patient. If attached oppositely, it may damage the CO₂ sensor or airway adapter.
- Dispose of calibration gas according to the regulation of each medical institution.

CAUTION Precautions about the BIS Monitoring (HBX-800 BISx I/F Unit)

- The conductive parts of sensors and connectors should not contact other conductive parts, including earth.
- To reduce the hazard of burns in the high-frequency surgical neutral electrode connection, the BIS sensor should not be located between the surgical site and the electrosurgical unit return electrode.
- The BIS sensor must not be located between defibrillator pads when a defibrillator is used on a patient connected to the patient monitor.
- Continuous impedance checking may need to be disabled if the 1 nA, 128 Hz impedance check signal interferes with other equipment (e.g., evoked potential monitors).
- Considerations when using Electro-Convulsive Therapy (ECT) equipment during BIS monitoring: Place ECT electrodes as far as possible from the BIS sensor to minimize the effect of interference. Certain ECT equipment may interfere with the proper function of the BIS monitoring system.

CAUTION Precautions about the Alarm

- Alarm messages will be displayed according to the priority. (Level S > Level H > Level M > Level L > Level N)
- For the same alarm level, the alarm message for the newer alarm will be displayed.
- The arrhythmia alarm message other than Tachy, Brady, Ext Tachy, Ext Brady will continue to be displayed for 30 seconds after the alarm is resolved.
- When "LEAD OFF", "Check Electrodes" is displayed, HR alarm or arrhythmia alarm will not function. If this condition is left unresolved, a sudden change of the patient may not be noticed. Take prompt action when the lead-off condition is detected.
- For the HPD-800/HPD-810/HPD-820 Gas Unit I/F and HCP-800/HCP-810/HCP-820 CO₂ Gas Unit, the upper EtCO₂ alarm will not generate if the upper limit is set to 100 mmHg/13.4 kPa and above as the measurement range is 0 to 99 mmHg / 0 to 13.3 kPa.
- Whether to use the SpO₂ second alarm function and its threshold selection should be based on the patient's clinical indication/portent and medical evaluation.
- If the SpO₂ alarm and second alarm setup is set to [OFF], the second alarm integral value will be set to 0.
- Pay attention not to set the alarm volume too low to avoid missing any important alarms.
- On a wired network, the alarm generated on the bedside monitor will be output to the network with a maximum delay of 1 second, and to the central monitor with a total delay of 2.5 seconds.
- If the NIBP alarm is turned OFF under the Oxygenator Mode, NIBP auto mode measurement and NIBP measurement at alarm occurrence will not be performed.

• If the same or similar equipments with different alarm settings are used in the same facility or same department, pay attention not to misjudge the alarms.

CAUTION Precautions about the System Setup

- When the waveform and numeric data display for each parameter is set to OFF, the alarm and trend input will be also suspended.
- If the HR/PR source is set to [BP], and if BP waveform/numeric data is set to [Disp. OFF], the PR value will not be displayed.
- If the HR/PR source is set to [SpO₂], and if SpO₂ waveform/numeric data is set to [Disp. OFF], the PR value will not be displayed.
- If the RR source is set to [CO₂/GAS], and if CO₂ waveform/numeric data is set to [Disp. OFF], the RR value will not be displayed.
- If the RR source is set to [CO₂/GAS], and if GAS waveform/numeric data is set to [Disp. OFF], the RR value will not be displayed.
- Do not set the same remote control ID to multiple monitors in the same floor. Otherwise, the remote control operation may control multiple monitors at the same time.
- After the remote control setup, check that the remote control unit is properly operating.
- If the time/date is not correctly set, or changed during monitoring, malfunction may occur with NIBP measurement, periodic printing, trend, NIBP list data, and age calculation from the birth date.
- If the time/date is changed, the time/date for all the saved patient data (trend, list, recall, etc.) will also change. The printed time/date before changing and the displayed time/date after changing will differ. Also, the data transmitted to the central monitor before the time/date is changed will be displayed on the central monitor with the previous time/date.

Precautions about the Multigas Unit Data Monitoring (MGU-800/MGU-810 Series)

- When the multigas unit (MGU-800/MGU-810 series) and HPD-800/HPD-810/HPD-820, HCP-800/HCP-810/ HCP-820 are simultaneously used, the CO₂ concentration measurement will be performed by the equipment selected for the "CO₂ Source Priority" under ([Menu] > "Parameter" [CO₂]).
- The MGU-800/MGU-810 series require warm up of about 10 minutes to correctly measure the data.
- If the power supply is interrupted due to power failure, etc., MGU-800/810 series multigas unit will initialize and enter into warm-up mode even if the power interruption is within 30 seconds.
- About the Gas Calibration

The gas calibration will automatically start when the MGU-800/810 series multigas unit is connected. After the warm-up completes, gas calibration will be performed every 4 hours during stable operation. During warm-up, gas calibration interval will become shorter than during normal operation. During gas calibration, measurement data will not be updated. Calibration gas is not required during gas calibration.

- Dispose of calibration gas according to the regulation of each medical institution.
- Make sure the sampling line and flow sensor is securely connected to prevent any leakage.
- The gas leakage inside or outside the equipment will cause measurement error. Make sure to check for leakage before usage.
- An environment with alcoholic vapor may adversely affect the measurement readings.
- CO₂, N₂O or anesthetic agent in the atmosphere around the MGU-800/810 series may adversely affect the measurement readings.
- SPIRO and ventilator cannot be used simultaneously.
- During the warming up process, the date of the last measurement accuracy check cannot be updated. Perform the measurement accuracy check after the warming up process is completed.
- If the accuracy check is performed using a low pressure gas, the accuracy of gas measurement will be reduced.

Make sure to perform the accuracy check using the specified calibration gas before its expiration date.

- If the error persists, refer to your nearest service representative.
- About the MGU-810 Series
 - The adult flow sensor dead space is 6.9 mL and the flow resistance is 1.8 cmH₂O at 60 L/min. The pediatric flow sensor dead space is 0.75 mL and the flow resistance is 0.9 cmH₂O at 10 L/min. Adjust ventilation accordingly.
 - To prevent condensation, the patient breathing circuit, flow sensor and pressure tubing should not be directly exposed to cooling equipment such as fans or cooling blankets.
 - Leakage of gas from the patient breathing system may occur if the pressure or gas sampling lines are not connected to the MGU-810.
 - The pressure tube and gas sampling lines of the flow sensor should always be routed from the patient circuit to the MGU-810 such a way as to avoid kinking.
 - Flow sensors that have suffered damage to sensor head, tubing or tubing connector must not be used.
 - If liquid has entered the pressure tubes, it can be removed by gently tapping or shaking the flow sensor.

CAUTION Precautions about the Patient Admit/Discharge

- If monitoring of a new patient is started without discharging the previous patient, data of the new patient will be added to the data of the previous patient which will result in inaccuracy.
- The user mode setting (alarm/display configuration) will remain effective even when the power is turned OFF or when the patient is discharged. Before monitoring, make sure the current user mode is suitable for the patient's condition.
- Resuming monitoring will also resume the alarm in suspension.
- **CAUTION** Precautions about the External Media

• Use only the specified external media.

- Use only the external media formatted on this equipment.
- Make sure to power cycle the system after the setup data is read from the CF card. By power cycling the system, the read data will become effective.
- Reading the patient data from the CF card will erase all previous patient data stored in the patient monitor.

CAUTION Precautions about the Maintenance

- When cleaning the touch panel, never use strong-acidic cleaning solution.
- To clean the touch panel, use an optional cleaning cloth, eyeglass cleaning cloth, soft cotton cloth, or nonwoven cloth (pulp, rayon, polyethylene, etc.).
- Clean the equipment frequently so stains can be removed easily.
- To prevent injury, it is recommended to wear gloves when cleaning the equipment.
- Pay attention not to allow chemical solution to enter the equipment or connectors.
- Do not use organic solvents, thinner, toluene or benzene to avoid damaging the resin case.
- Do not polish the equipment with abrasive or chemical cleaner.
- When disinfecting the entire room using a spray solution, pay close attention not to get any solution into the equipment or connectors.
- Use only neutral detergent to clean the equipment. The surface resin coating may damage, resulting in discoloration, scratches, and malfunction. Example:

chemical cloth, scrub brush, abrasive, polishing powder, hot water, volatile solvent and chemicals (cleanser, thinner, benzine, benzol, and synthetic detergent for house and furniture), or sharp-edged tools

- Do not open the housing.
- Do not allow alcohol or other liquids to enter the equipment.

• Replace the periodic replacement parts periodically as specified.

Wired Network (DS-LANII/ DS-LANIII)

WARNING

- Do not connect unspecified device to the wired network.
- Do not mix devices with DS-LANII and DS-LANIII setting in the same wired network. The network may cease and proper monitoring may not be possible.

- If performing wired network transmission, configure the display so that the numeric data corresponded to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted.
- The default setting of bed ID is "000". If connected to a wired network with the bed ID unchanged, monitoring on the central monitor will not be possible.
- When connecting to a wired network, make sure that there are no other bedside monitors with the same ID. If there is more than one bedside monitor with the same bed ID, the duplicated bedside monitors cannot be monitored on the central monitor.
- When connected to the DS-LAN II network, set the bed ID in the range from "001" to "048".
- When connected to the DS-LAN III network, set the bed ID in the range from "001" to "100".
- The alarms that can be notified to the central monitor depend on the model type and software version of the central monitor. For details, refer to the operation manual of the central monitor.
- There are following restrictions when connecting the DS-8400 System to the wired network.
 - The BP measurement unit setting should be the same for all central monitors and bedside monitors. If the setting is different among the monitors, data such as BP waveform, BP numeric data, NIBP numeric data, NIBP list will not be transmitted. It will be treated as not measured data, and will not be displayed on the central monitor. The alarm limit setup from the central monitor cannot be performed either.
 - On the DS-LAN II network, the following arrhythmia alarms will not be transmitted. TACHY, BRADY, COUPLET, PAUSE, TRIGEMINY, TRIPLET, EXT TACHY, EXT BRADY, R on T, MULTIFORM, VENT RHYTHM, SVT, IRREGULAR RR, PROLONGED RR, S FREQUENT, S COUPLET, VPC, SVPC, PACER NOT CAPTURE, PACE NOT PACING
 - On the DS-LAN II network, arrhythmia alarm of "SLOW VT" will be transmitted as "VT" .
 - On the DS-LAN II network, waveform, numeric data, and alarm of BP7, BP8, TEMP3 to 8 will not be transmitted. (These can be transmitted on the DS-LAN III network). Also, the displayable waveform, numeric data, and alarm will differ depending on the central monitor model type. Refer also to the operation manual for the respective central monitor.
 - The PR_IBP alarm will not be transmitted to the central monitor.
 - If the "RR/APNEA alarm source" is other than [Impedance] (or, if [Auto] selects a setting other than [Impedance]), the RESP waveform will not be transmitted on a wired network.
 - If the "RR/APNEA alarm source" is other than [CO₂/GAS] (Or, if [Auto] selects a setting other than [CO₂/GAS]), the CO₂ waveform will not be transmitted on a wired network.
 - For the numeric data displayed as "xxx", maximum or minimum value of measurable range will be transmitted.
 - The numeric data displayed as "--- " will be treated as not measured data.
 - If the measurement unit of CO₂ concentration is mmHg and [99mmHg] is selected for "Upper Limit of CO₂ (mmHg) Transmission" ([Initial Settings]>[DS-LAN]), the CO₂ value of 100 mmHg or above will be transmitted as 99 mmHg.
- As the DS-8400 System do not have the arrhythmia template display and 12-lead ST display function, waveforms and other data will not be displayed for these displays on the central monitor connected to the DS-

LAN network.

- When connected to the wired network, the time/date will synchronize with the central monitor. Even if the time/date is changed on the DS-8400 System, it will be corrected to the time/date of the central monitor.
- The ST display will be distorted on the central monitor if the ECG lead (ECG1 or ECG 2) is changed on the DS-8400 System. Redrawing the ST display will return the display to normal.
- On the central monitor, the respiration waveform and RR value based on the RR/APNEA alarm source selected on the DS-8400 System will be displayed. The monitored RR and APNEA will be the same for the central monitor and the DS-8400 System.

Wireless Network System

DANGER

• When monitoring a patient using wireless telemetry, make sure the patient data is properly received at the central monitor. Pay special attention when the channel ID at the bedside monitor is changed.

WARNING

- A password can be set to access the channel ID setup menu to allow only the telemetry channel administrator to change the channel ID.
- Some type of wireless combinations may generate interference with other telemetry.
- Before selecting a channel, verify it will not interfere with other channels.
- Inform the supervisor of the use of telemetry channels to avoid interference with other telemetry.
- If transmitters are used in a neighboring medical facility, your facility and the neighboring facility must make agreements on the setting of the telemetry channels to prevent telemetry interference.

CAUTION Precautions about the Telemetry

- When performing telemetry transmission, configure the display so that the numeric data corresponded to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted.
- The setup of channel ID and group ID should be performed only by the telemetry channel administrator or our service representative. Users should not perform this procedure as malfunction may occur.
- When the measurement unit of BP is "kPa", corresponding waveform and numeric data will not be transmitted. When using a wireless network, use "mmHg" for the BP measurement unit.
- BP waveform with a scale above the set scale can not be properly transmitted. When transmitting the BP waveform, pay attention to the displayed BP scale.
- If the measurement unit of CO₂ concentration is mmHg and [99mmHg] is selected for "CO₂ (mmHg) Upper Limit of Transmission" ([Initial Settings]>[System]>[Telemeter]), the CO₂ value of 100mmHg or above will be transmitted as 99mmHg.

RTC and Data Backup

- This equipment is equipped with a built-in clock. When the power of this equipment is turned OFF, this clock is backed up by a lithium primary battery. If incorrect time is displayed when turning ON the power, a low battery may be the cause. In such case, contact Fukuda Denshi service representative for replacing the battery.
- To protect the data during voltage dip, short interruptions and voltage variations on power supply input lines or during short duration of power turned OFF, this equipment performs 5-minute (approx.) data backup using the secondary battery. The data may not be protected if the power is turned OFF within 30 minutes from power ON.

Precautions about the Ventilator Monitoring

WARNING

- The ventilator alarm sound is set to OFF (factory default). The alarm sound can be turned ON on the Tone/Volume setup screen.
- If the DS-8400 System does not generate an alarm even though the ventilator is generating an alarm, or if any other malfunction occurs, immediately check the ventilator, this equipment, cable, and replace the cable if necessary. If the malfunction persists, stop using the equipment.
- The alarm generation on the DS-8400 System is not guaranteed if the alarm other than the specified one generates at the ventilator.

(@Maintenance Manual "Ventilator Measurement and Alarm Input" P4-1)

- The ventilator operation should be performed by well-trained and authorized personnel.
- When connecting this equipment and a ventilator, use only the specified connection cable.
- Verify that this equipment and the ventilator are properly connected.
- When connecting the cable, verify that the main power of this equipment and the ventilator are OFF.

Precautions about the SpO₂ Sensor

A DANGER Danger of Burn Injury Caused by the SpO₂ Sensor

• When monitoring SpO₂, make sure to use only the specified sensor/relay cable. If any other sensor/relay cable is used, a high temperature rise of the sensor may place the patient in danger of burns. If there are any questions regarding the sensor/relay cable use for SpO₂ measurements of this equipment, please contact Fukuda Denshi service representative.

Precautions about the Masimo Model

• Possession or purchase of this device does not convey any express or implied license to use the device with unauthorized sensors or cables which would, alone or in combination with this device, fall within the scope of one or more of the patents relating to this device.

Precautions about the NIBP Cuff

• Some of the NIBP cuffs used for this equipment contain natural rubber latex which may cause allergic reactions.

(FDA: Medical Alert on Latex Products, "Allergic Reactions to Latex-Containing Medical Devices", Food & Drug Administration, 10903 New Hampshire Avenue, Silver Spring, MD 20993, 1991.)

Precautions about Disposing of the Equipment, Accessories, or Components

- When disposing of this equipment, accessories, or components, use an industrial waste distributor. Do not dispose of as ordinary waste.
- When disposing of the battery, separate it from other wastes and contact your nearest service representative.

Precautions about Transportation

 When transporting this equipment, pack it with specified packing materials. Also, transport it under appropriate environment condition.
 (Operation Manual "Specification" P14-1)

Monitoring after Power Failure

When the power failure is less than 30 seconds, monitoring will resume with the display mode and patient information unchanged. When the power failure is 30 seconds or more, monitoring will resume with the default display mode set by the user, or the display mode which was last set.

However, this is only if the equipment was operated for 30 minutes or more before the power failure.

MGU-800/810, HCP-800/810/820, HPD-800/810/820 will start up from the warm-up mode. The warm-up time differs for each unit.

To Prepare for Emergency Use

Accessories/Optional Accessories

- The ECG electrodes are consumable products. Always prepare extra supplies of electrodes.
- Verify that there is no wire break on the patient cable once a week.

Battery Pack

- Even if the battery pack is not in use, the remaining capacity decreases due to self-discharge.Make sure to verify once a week that the battery pack is fully charged.
- To fully charge the empty battery pack, it takes 5.0 hours during operation, and 2.5 hours when the power is OFF and AC cable is connected.
- The performance of the battery deteriorates with repeated use. To ensure performance of the battery, it is recommended to replace it once a year.

Electromagnetic Compatibility

The performance of this equipment under electromagnetic environment complies with IEC 60601-1-2: 2007.

Do not use any unauthorized equipment or cables as they may not comply with the EMC standard.

Precautions for Safe Operation under Electromagnetic Influence

If any sorts of electromagnetic wave, magnetic field, or static electricity exist around the device, noise interference or malfunction of the device may occur. If any unintended malfunction or noise occurs during monitoring, check the electromagnetic influence and take appropriate countermeasures.

The following are examples of the common cause and countermeasures.

A DANGER Static Electricity

In a dry environment (room), static electricity is likely to occur. Take the following countermeasures.

- Both operator and patient should remove any static electricity before entering the room.
- Humidify the room.

WARNING Cellular Phone

• The radio wave may cause malfunction to the device. Cellular phones and radio sets should be turned off in the room (building) where medical device is located.

CAUTION Lightning

A lightning nearby may induce excessive voltage to the equipment. If any danger is suspected;

- Use the uninterruptible power supply system.
- Use the battery.
- **CAUTION** High frequency noise interference from other device through the power outlet
- Check where the noise is originated and remove it using filtering device, etc.
- Stop using the device that is originating the noise.
- Use other power outlet.

EMC Guidance

This device complies with the Safety Standard IEC 60601-1-2: 2007. However, if portable transmitter or wireless LAN equipment is used extremely nearby, the electromagnetic influence may largely exceed the compliance level and may cause unexpected phenomenon such as noise interference on the waveform, etc.

Also, if this equipment is installed close to, or stacked with other equipment, malfunction may occur. Make sure to verify that the equipment operates properly in a used location.

This equipment should be used in a location specified by each medical institution.

If any unexpected noise interference on the waveform or failure to the peripheral device occurs, stop using the equipment and follow the instruction of the technical engineer.

The following is the information relating to EMC (Electromagnetic Compatibility). (When using this equipment, verify that it is used within the environment specified below.) This equipment complies with IEC 60601-1-2: 2007 for the following system configuration.

- Main Unit: DSC-8410
- Display Unit: LC-8018TC
- Patient Monitor: DS-8007
- Multi Module: HM-801
- Multi Module: HM-800 x 2
- Recorder Unit: HR-800
- Telemetry Transmitter Module: HLX-801
- Gas Unit I/F: HPD-820
- Input Box: IB-8004
- Multiport Module HP-800
- Multigas Unit: MGU-810
- BISx I/F Unit: HBX-800
- BISx Module: BISx
- Lithium-Ion Battery: BTO-005

Compliance to the Electromagnetic Emissions

The DS-8400 System is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-8400 System should assure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions				
Emissions Test	Compliance	Electromagnetic Environment - Guidance		
RF Emissions CISPR 11	Group 1	The DS-8400 System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.		
RF Emissions CISPR 11	Class A			
Harmonic Emissions IEC 61000-3-2	Class A	The DS-8400 System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network which supplies buildings		
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	used for domestic purposes.		

Compliance to the Electromagnetic Immunity (1)

The DS-8400 System is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-8400 System should assure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity				
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance	
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV: contact ±8kV: air	±6kV: contact ±8kV: air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Electrical fast transient/burst IEC 61000-4-4	±2kV: power supply lines ±1kV: input/output lines	±2kV: power supply lines ±1kV: input/output lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surge IEC 61000-4-5	±1kV: differential mode ±2kV:common mode	±1kV: differential mode ±2kV:common mode	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-11			Mains power quality should be that of a typical commercial or hospital environment. If it is required to continuously operate the DS-8400 System during power failure, it is recommended to operate on an uninterrupted power supply.	
Power Frequency (50/60Hz) Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	

*: U_T is the AC mains voltage prior to application of the test level.

Compliance to the Electromagnetic Immunity (2)

The DS-8400 System is intended for use in the electromagnetic environment specified below. The customer or the user of the DS-8400 System should assure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity				
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance	
			Portable and mobile RF communications equipment should be used no closer to any part of the DS-8400 System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance	
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	d = 1.2 √₽	
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	d = 1.2√戸 80MHz to 800MHz d = 2.3 √戸 800MHz to 2.5GHz	
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^{*1} , should be less than the compliance level in each frequency range ^{*2} . Interference may occur in the vicinity of equipment marked with the following symbol:	
Note 1	At 80MHz and 800MHz_th	e senaration dist	-	
	At 80MHz and 800MHz, the separation distance for the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
*1:	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DS-8400 System is used exceeds the applicable RF compliance level above, the DS-8400 System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DS-8400 System.			
*2:	Over the frequency range	150kHz to 80MH	z, field strength should be less than 3V/m.	

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the DS-8400 System

The customer or the user of the DS-8400 System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DS-8400 System as recommended below, according to the maximum output power of the communications equipment.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the DS-8400 System				
Rated Maximum Output	Separation Distance according to Frequency of Transmitter (m)			
Power of Transmitter (W)	150kHz to 80MHz d = 1.2 √P	80MHz to 800MHz d = 1.2 √₽	800MHz to 2.5GHz d = 2.3 √p	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Chapter 1 General Description

Composition of the System	1-1
Features	1-3
Menu Configurations	1-5

Chapter 1 General Description

Composition of the System

The DS-8400 system is composed of Display Unit (LC-8016TC/LC-8018TC), Main Unit (DSC-8410), Super Unit (HS-8000 series/DS-8007 series), Multigas Unit (MGU-800 series/MGU-810 series), Recorder Unit (HR-800), expansion modules and Input Box (IB-8004).





Lineup of Main Unit

Model Type	External Monitor Output	Extended Display Unit Output	LAN (TCP/IP)
DSC-8410	Yes	No	No
DSC-8410 (with CC-84)	Yes	2ch	1ch

Lineup of Display Unit

Model Type	Display Size
LC-8016TC	15.6 inch
LC-8018TC	18.5 inch

Lineup of Super Unit	
----------------------	--

Model Type	Fixed Parameter	SpO ₂ Unit	Multiparameter Measuring Items	CO ₂ Measurement (Optional)
HS-8312N	ECG (Max.12 leads), RESPx1, NIBPx1 SpO ₂ x1	Nellcor	3 ports Temperature x6 (maximum) BP x6 (maximum) CO Measurement x1 (maximum)	Yes
HS-8312M	ECG (Max.12 leads), RESPx1, NIBPx1, SpO ₂ x1, SpCO x1 [*] , SpMet x1 [*] , SpHb x1 [*] , PVI x1 [*]	Masimo		Yes
DS-8007N	ECG (Max.12 leads), RESPx1, NIBPx1, SpO ₂ x1, RR_SpO ₂ x1 [*] , TEMPx2	Nellcor	2 ports TEMP x4+2 phone plugs (maximum) BP x4 (maximum) CO Measurement x1 (maximum)	Yes
DS-8007M	ECG (Max.12 leads), RESPx1, NIBPx1, SpO ₂ x1, TEMPx2, PI, SpMet x1 [*] , SpCO x1 [*] , SpHb x1 [*] , SpOC x1 [*] , PVI x1 [*]	Masimo		Yes

*: RR_SpO₂, SpMet, SpCO, SpHb, SpOC, and PVI measurements are optional functions.

Lineup of Adapter

Model Type	Remarks
HSA-80	Adapter for HS-8000 series (Connects to the patient monitor with module connection cable.)
HSA-81	Adapter for HS-8000 series (Connects to the patient monitor by sliding on the rail at the rear side of the patient monitor.)
DSA-82	Adapter for DS-8007 series (Connects to the patient monitor with module connection cable.)

Lineup of Multigas Unit

Model Type	CO ₂ /N ₂ O Measurement	Anesthetic Agent Measurement	O ₂ Measurement	Spirometry Function Assessment
MGU-801P	Yes	Yes	Yes Paramagnetic	No
MGU-811P	Yes	Yes	Yes Paramagnetic	Yes

Lineup of Recorder

	Model Type	Remarks	
Ī	HR-800	50 mm Roll Paper, 3 waveforms printing	

Lineup of Input Box

Model Type	Number of Slots
IB-8004	4

Model	Module	Parameter
HM-800	Multi Module	(IBP, TEMP, CO) x 2
HM-801	Multi Module	(IBP, TEMP, CO) x 1 AUX Connector x1
HP-800	Multiport Module	Serial Communication 2ch, Analog Communication 1ch
HG-810	SpO ₂ Module M (Masimo)	SpO ₂ , SpCO [*] , SpMet [*] , SpHb [*] , PVI [*] , PR
HG-820	SpO ₂ Module N (Nellcor)	SpO ₂ , PR

Lineup of Expansion Module

*: SpCO, SpMet, SpHb, PVI measurements are optional functions.

Lineup of Gas Module

Model	Module	Parameter
HPD-800	Gas Unit I/F	
HPD-810		CO ₂ measurement by Mainstream method with a connection to Philips Capnostat 5
HPD-820		
HCP-800	CO ₂ Unit	Incorporates Microstream technology developed by
HCP-810		Covidien
HCP-820		Sidestream Method

Lineup of EEG Module

Model	Module	Parameter
HBX-800	BISx I/F Unit	BIS, SQI, EMG, SR, EEG, SEF, TOTPOW

Features

• Various displays such as enlarged numeric data, trend, or ventilator can be selected according to monitoring conditions.

The operation can be performed with the touch panel. Also, frequently used keys can be programmed as user key.

- It is possible to connect two types of display units in addition to the main display unit to extend the display.
- An optional mouse can be connected allowing mouse operation. (However, some operations are not possible.)
- The alarm indicator notifies the alarm with different flashing patterns corresponding to the alarm level so that the users can easily identify the alarm level of the generating alarm.
- Remote control is possible using the optional remote control unit. Using the CF-820 IR Remote Control Unit allows to remotely control the patient monitor.
- By using the multiparameter amplifier, the HS-8000 series Super Unit is capable of monitoring parameters in combination of BP (maximum 6 channels), temperature (maximum 6 channels), and CO (maximum 1 channel). In addition to ECG, respiration, SpO₂ (pulse wave), BP, NIBP, temperature and CO, it is also possible to measure CO₂ and BIS as optional function.
- By using the multiparameter amplifier, the DS-8007 is capable of monitoring basic parameters of BP (max. 4 ch.), temperature (max. 6ch.), CO (max. 1ch.), and optional parameters of CO₂ and BIS.
- 27 types of arrhythmia can be analyzed.

- By using the optional SpO₂ Module (HG-810/HG-820), arterial oxygen saturation can be also measured. By using the system with the Super Unit, arterial oxygen saturation measured at 2 different sites can be monitored as additional parameter.
- This system uses pulse oximetry to measure and display functional oxygen saturation in the blood. There are two model types with different built-in SpO₂ modules, which are Covidien/Nellcor and Masimo.
- SpCO, SpMet, SpHb, PVI measurement are optional function available when using the HS-8312M/DS-8007M and HG-810 with built-in Masimo SpO₂ module.
 - SpOC measurement is optional function available when using the DS-8007M.
- RR_SpO₂ is optional parameter which can be measured on the DS-8007N with the built-in Covidien SpO₂ module.
- By using the optional HP-800 Multiport Module, or connecting the ventilator to Status II port on the main unit, airway flow, airway pressure waveform, minute ventilation, airway resistance, etc. can be monitored. Also, ventilator alarm can be notified to the central monitor via wireless and wired network. The following ventilators can be connected.
 - SV-300/300A
 - SERVO-i, SERVO-s, SERVO-U/n/air
 - PURITAN-BENNETT Ventilator 740/760, 840
 - Evita 4/Evita XL/Evita 2 dura
 - Velia, Ultra, Astral
- Wired network (DS-LANII/DS-LANIII) is possible via the Ethernet LAN cable. DS-LAN II is a network based on 10BASE-T with transmission speed of 10 Mbps and maximum transmission distance of 100 m. DS-LANIII is a network based on 100BASE-TX with transmission speed of 100 Mbps and maximum transmission distance of 100 m.
- Wireless network construction is possible using the optional telemetry transmitter module (HLX-801 (FA) / HLX-801 (G)).
- By attaching the module to the built-in slot, monitoring parameters can be added without the Input Box.
- By using the optional Multi Module (HM-800/HM-801), the monitoring parameters can be added. To use more than one expansion module, optional Input Box (IB-8004) is required.
- By using the optional Recorder Unit (HR-800), the measurement data can be printed.
- By connecting the optional Multigas Unit (MGU-800/810 series), CO₂ concentration, anesthetic gas concentration, O₂ measurement, N₂O concentration can be measured. The following anesthetic agents can be measured.
 - Halothane
 - Isoflurane
 - Sevoflurane
 - Enflurane
 - Desflurane
- By connecting the Gas Unit I/F (HPD-800/HPD-810/HPD-820) or CO₂ Gas Unit (HCP-800/HCP-810/HCP-820), CO₂ concentration can be measured.
- By using the HP-800 Multiport Module, or by connecting the FLOW-i Anesthesia Delivery System to Status II port or to COM1 to COM4 port on the main unit, CO₂ concentration, anesthetic gas concentration (ISO, SEV, DES), O₂ concentration, N₂O concentration, airway flow, airway pressure waveform, minute ventilation, airway resistance, etc. can be monitored.
- By using the HP-800 Multiport Module, or connecting the Oximeter to Status II port or to COM 1 to 4 port on the main unit, SvO₂, CO, etc. can be monitored. The following Oximeter/CCO measurement device can be connected.
 - Vigilance
 - Vigilance CEDV

- Vigilance II
- Vigileo
- PiCCO2
- EV-1000
- PulsioFlex (connects to COM1 to COM4 port)
- By using the HP-800 Multiport Module, or by connecting the A-2000 BIS Monitor/A-3000 BIS Vista (Covidien) to Status II port or to COM 1 to 4 port on the main unit, the patient's wakeful state can be monitored.
- By using the HP-800 Multiport Module, or by connecting the INVOS 5100C Cerebral Oximeter (Covidien) to Status II port or to COM1 to 4 port on the main unit, regional cerebral oxygen saturation data can be monitored.
- By connecting the following transcutaneous blood gas monitors to COM1 to COM4 ports, transcutaneous blood gas partial pressure can be monitored.

• TCM4

• TCM5 FLEX

Menu Configurations

The menu configuration of this equipment is as follows.

Menu Screen

The menu screen is a group of shortcut keys to jump to each menu. The menu is composed of the following 10 groups and can be accessed from the menu screen.

Function Groups	Displayed Menu
Admit/Discharge	Admit/Discharge
Basic Setup	Maximum of 9 functions are displayed.
Alarm	Maximum of 9 functions are displayed.
Parameter	Maximum of 18 functions are displayed.
Data Review	Maximum of 9 functions are displayed.
Waveform Review	Maximum of 9 functions are displayed.
Calculation	Maximum of 5 functions are displayed.
Other Bed	Other bedside monitors connected to the DS-LAN will be displayed.
Initial Settings	Initial settings menu will be displayed.
Maintenance	Maintenance menu will be displayed.

REFERENCE

• Other than the "Initial Settings", the items to be displayed on the menu screen can be customized by groups.

(@Maintenance Manual "Menu Setup" P5-21)

Admit/Discharge

Admit/Discharge	Mode Selection
	ID, Name, Classification, Sex, Team, Birth Date, Age, Height, Weight, BSA, Blood Type (ABO, Rh), Pacemaker, Impedance Meas., Admit Date/Time
	Monitor Suspend
	Discharge

Basic Setup

Display Configuration	Layout, Background, Palette, Detail Setup, Meas., Wave (Sweep Speed), Short Trend, User Key
Manual Printing	Basic (Printer, Waveform, Print Duration, Delay Time), 12-Lead (12-Lead Waveform/ Analysis Format, Position, Wave Format, Print Calibration, Printer Auto Scale, Lead Boundary), Other Setup (Graphic Printing, Recall Printing), Common (QRS Classific., Speed, Calibration: Print Calibration, Print NIBP Data)
Auto Printing	Alarm Printing (Print, Printer, Waveform, Print Duration), Periodic Printing (Print, Printer, Periodic Interval, Waveform, Print Duration), Common (QRS Classific, Speed, Print Calibration, Print NIBP Data)
Tone/Volume	Vital Alarm Sound, Ventilator Alarm Sound, Status Alarm Sound, Tone Source, Key Sound, Other Bed Alarm Sound, Boot/Shutdown, Other
Time/Date	Time, Date
Color	Waveform/Numeric Data, Background, Palette, User Key
Brightness	Brightness
Night Mode	Night Mode, Detail Setup (Volume, Display, Alarm Indicator)

Alarm

Basic	The parameters to be displayed are selectable.
Dasie	
	Alarm Suspend, Mode Select, Print, All Auto
Circulatory	Alarms for HR, Ext Tachy, Ext Brady, PR, SpO ₂ , Ext SpO ₂ , NIBP, BP, TEMP can be set.
	Alarm Suspend, Mode Select, Print, All Auto
Respiratory/Gas	Alarm for RR, APNEA and gas can be set.
	Alarm Suspend, Mode Select, Print, All Auto
Arrhythmia	Arrhythmia Alarm, Detail Setup
ST	ST Alarm, Waveform Review (ST), Update Ref. Wave
List	List of alarm ON/OFF setting and lower/upper limits, Meas. List/All List, Print Setup, Recall Setup
Detail Setup	Alarm Suspend Time, Alarm Silence Time, Alarm Silence, Alarm Sound Suspend, Status Alarm Control, Alarm Limit Display

Parameter

ECG		Arrhy., Arrhy. Alarm Setup, ST Setup, HR
		Size, Lead, Optimize Size, Alarm Assist, Disp. ON/OFF
		Detail Setup Filter, Synchronized Mark/Tone, Pacemaker, Pacemaker Pulse, Pace Pulse Mask Time, HR Average, HR Delay, ECG Drift Filter, AC Filter, Auto Lead, 3Lead Override, ST/VPC/ Arrhy. Alarm Display, ECG Analog Output, ECG waveform display during Lead-OFF, Noise Detection
RESP		Size, Common Setup (RR Synchronized Mark, RR/Apnea Alarm Source), RR, Apnea, Alarm Assist, Display ON/OFF
		Impedance Setup CVA Detection, Impedance Measurement, Impedance Detection Lead, Impedance Detection Level
NIBP		NIBP Auto Mode, NIBP S, M, D, Alarm Assist, Cancel Error
		Detail Setup Patient Classification, Dyna Alert, Oscillograph, PR Display, NIBP Erase Time, Measure at Alarm, Quick Measurement, Sight Inflation, MAP, End Tone, User Interval, Auto Mode with Start/Stop key, Time Display, Periodic Measurement Starting Time, Oscill. Print, Target Inflation Value
BP		BP Zero (BP1 to BP8), Scale Selection, Label, Alarm Assist, Display ON/OFF
		Detail Setup Synchronized Mark/Tone, Display Type, Wave Filter, Mean Wave, Respiration Filter, IBP Analog Output, Alarm during NIBP, ART Catheter Check Message
SpO ₂		Size, Label, Alarm Assist, Disp. ON/OFF
	HS-83xxN/HG-820 /DS-8007N	Detail Setup Alarm during NIBP, Synchronized Mark/Tone, Second Alarm
	HS-83xxM/HG-810 /DS-8007M	Detail Setup Alarm during NIBP, Synchronized Mark/Tone, SpO ₂ Averaging, Pulse Sensitivity, FAST SAT, PI Display, Signal IQ Wave
Sp*		SpCO, SpMet, SpHb Setup
Sp*	RR_SpO ₂	Common Setup RR Synchronized Mark, RR/Apnea Alarm Source), Display ON/OFF, Alarm Assist, RR, Apnea
TEMP		Label, ∆T Setting, Alarm Assist, T1 to T8 Disp. ON/OFF
GAS		Scale, Gas Calibration, Detail Setup (Flow Rate, Wave Clip, CO ₂ Source Priority)
CO ₂		Scale, Airway Adapter Calibration, Alarm Assist, Disp. ON/OFF
		Detail Setup EtCO ₂ Peak Duration, N ₂ O Compensation, Atmospheric Pressure, O ₂ Compensation, Anesthetic Compensation, CO ₂ Source Priority
BIS		Scale, Alarm Assist, TREND-E
		Common Setup Short Trend 2nd Parameter, Smoothing Rate, Continuous Impedance Check, EEG Filter
External Device		SvO ₂ /CCO, VENT, INVOS

Data Review

Graphic Trend	Latest Data, Alarm Review, Trend Group, Alarm Disp. Sel., Print
Tabular Trend	Latest Data, Alarm Review, List Group, Setup, Print, Print (All)
Recall	Latest Data, Display Selection, Setup, Delete Sel.
OCRG	Latest Data, Resp. Wave, Impedance, Resp. Wave Size, Print
Alarm History	Latest Data, Display Selection, Print

Uwaveform Review

Zoom Wave	Latest Data, Alarm Review, Meas., Print, Delete
ST	ST Waveform, Reference Waveform, Setup, Slide Show, Size, Latest Data, Print
12-lead	Latest Data, Review, Start Analyze, Setup, Print
Full Disclosure Waveform	Latest Data, Alarm Review, Slide Show, Time Search, Size/Scale, Setup, Alarm Display, Print

Hemodynamics	Input Data, Edit, calculation results list, New Regist., Index Display, Print
Lung Function	Input Data, Edit, calculation results list, New Regist., Index Display, Print
СО	Meas., Edit, Setup, Hemodynamics, Average CO Input, Delete Sel.
Drug Calculation	Drug Name, Drug Amount, Diluent Amount, Dosing Rate, Flow Rate, Weight

Other Bed

Other Bed Display	Area Selection (Area 1 to 5), Other Bed Alarm Sound, Alarm Display, Area Setup (Area 1 to 5), Bed List
Other Bed Display	Area Selection (Area 1 to 5), Other Bed Alarm Silence, ON/OFF of menu title display, Waveform Selection

Initial Settings

Alarm	-	Alarm System, Basic Alarm Parameter, Oxygenator Mode Setup, Buzzer Tone at Speaker Failure, Suspend Arrhy. Analysis during Noise Interference, Lower Limit for Alarm Volume, Alarm Indicator, Alarm Level, HR/PR Lower Limit during Alarm Auto Setting, Alarm Threshold Limit, Auto Alarm Setup
Measurement	User Label	BP User Label, TEMP User Label
	Unit	CO ₂ , BP, CVP, TEMP, ST, Height/Weight
	Other	NIBP Start 5min. early, MAP Calc.(ART, NIBP), Arrhythmia Analysis Filter, Synchronized Mark/Tone Priority, HR/PR Source Priority, Gas Display during Undetected Breath, Catheter Manufacturer for CC Input
User I/F	Display/Print	Date Format, BP Alarm Increment, Trend Clip, BP Printing Scale, Night Mode Cancel, ST Display Lead Setup, VENT Display Parameters, Hemo, etc Display Parameters, Auto Display Configuration, Dim All Data Other than Numeric, All Window Opaque, Printer Message Display, Message Icon, Operation Guide Display, Notification when Changing Equipment Configuration, Sync wave size/scale of extended display with main unit, 12-lead Analysis Filter Display, Waveform Size Display, Shift Time (Day Shift, Twilight Shift, Night Shift), Key Group Setup, Event Label Setup, RR Alarm Increment, Patient Name on the Information Display Area, External Device Numeric Data Box Operation, Drug Calculation, Monitor Suspend Label, Monitor Suspend Timer
	Power ON/ Discharge	Check Discharge at Power ON, Discharge Mode, NIBP Resume Auto Mode with Manual Measurement., Backup Setting at Power ON/Discharge
	Menu	Items to be displayed on the menu screen can be selected.
	Key Mask	Items not to be displayed on the menu screen can be selected.
	Remote Control	Remote Control Key Function, Room ID, Bed ID
	Operation	Mouse, Auto Hide Window, Auto Minimize
External Device	Main Unit Port HP-800	COM, Status II, U-LINK, Numbering of HP-800 Ventilator (SV-300, SERVO-i/s, SERVO-U/n/air, PB, Evita), SvO ₂ /CCO (Vigilance, PiCCO, PulsioFlex), GAS/SPIRO (MGU-800, MGU-810), Other (PC Comm., HLX, Barcode, Magnetic Card, BIS, INVOS, FLOW-i, TCM4, TCM5), PC Comm. (DS-5000), Trend Data Setup
	Magnetic Card Reader	Data digits for each patient information, Auto Reference to Central Monitor when Reading Patient ID
	Network	Main Unit (IP Address, Sub-Network Mask, Default Gateway), Printer (network printer, IP address, MAC address, printer specification, paper size)
	Status Output	Alarm Output
	Analog Output	Analog Output Setup, Sync. Signal Output
System	DS-LAN	DS-LAN Setup, Room ID, Bed ID, DS-LAN Pat. ID Tx Transmission Start Position, Synchronize Hemodynamic Data with the Central Monitor, CO ₂ (mmHg) Upper Limit of Transmission
	Telemeter	Tele. ON/OFF, Channel/Group ID, Telemetry Wave, CO ₂ (mmHg) Upper Limit of Transmission
	Unit Module	Multiparameter connector setup of HS-8000/DS-8007, HM-800/HM-801, channel setup of ${\rm SpO}_2$
	Other	AC filter, Search Patient ID, Data Transfer, Data Selection for Transfer, Numeric Data External Output, HLX, HR-800, Extended Display Unit
User Mode Registration	-	Registering/Changing/Initializing the User Mode, Link with Patient Classification
Administrator Setup	Key Lock	Key lock for each function can be set.
	Password Setup	Password for each administrator level can be registered/changed.

Maintenance

Maintenance	Program Version, External Media, GAS Accuracy Check, Parts Usage Time, Install,
	Module Install, Test Menu

Chapter 2 Name of Parts and Their Functions

Name of Parts and Their Functions	2-1
Main Unit: DSC-8400 series	2-1
Display Unit: LC-8018TC/LC-8016TC	
Super Unit: HS-8000 Series	2-4
HS Adapter: HSA-80	
HS Adapter for DS-8400: HSA-81	2-6
Transport Monitor: DS-8007 Series	2-6
Adapter for DS-8007: DSA-82	
Multi Module: HM-800	2-9
HM-801 Multi Module	2-10
SpO2 Module (HG-810/HG-820)	2-11
Multiport Module: HP-800	
Recorder Unit: HR-800	
Input Box: IB-8004	2-13
CO2 Gas Unit: HCP-800	2-13
CO2 Gas Unit: HCP-810	2-14
CO2 Gas Unit: HCP-820	2-15
Gas Unit I/F: HPD-800	2-15
Gas Unit I/F: HPD-810	2-16
Gas Unit I/F: HPD-820	2-17
BISx I/F Unit: HBX-800	
Multigas Unit: MGU-800/810 Series	2-18

Chapter 2 Name of Parts and Their Functions

Name of Parts and Their Functions

WARNING

 Do not connect a unit or cable not authorized by Fukuda Denshi to any I/O connector. If done so by mistake, the equipment cannot deliver its maximum performance and the connected units may be damaged, resulting in a safety hazard.

Main Unit: DSC-8400 series

Generation Front Side

- 1 Display Unit Attaching Position Attaches the display unit (LC-8018TC/LC-8016TC).
- 2 Display Unit Connector Connects the display unit.



Rear Side

- 1 HLX Fixing Position Fixes the Telemetry Transmitter Module (HLX-801 (FA) / HLX-801 (G)).
- 2 Potential Equalization Terminal Used for equipotential connection.
- 3 Power Supply Connector Connects the power cable.
- 4 Battery Attaches the battery (BTO-005).



Left Side

- 1 Super Unit Connector Connects the specified equipment.
- 2 Module Connector Connects the specified equipment.
- Super Unit Release Lever
 Releases the Super Unit from the main unit.



Right Side

The illustration shown on right is when the optional CC-84 and HR-800 are installed. No. 11 to No. 14 described below can be used only when the CC-84 is installed.

- 1 Card Slot (CF, CFast) A slot to insert the specified external media
- 2 Card Access Indicator (CF, CFast) Indicates the external media access status.
- 3 DS-LAN Connector Connects to the wired network using the Branch Cable (CJ-522).
- 4 I/O Connector Connects the specified equipment.
- 5 External Monitor Connector Connects the external monitor.
- 6 AUX Connector For future function enhancement.
- 7 Serial Connector (COM1 to 4) Connects the specified equipment.
- 8 Status II Input/Output Connector Connects the specified equipment.
- 9 U-LINK Connector Connects the Recorder Unit (HR-800) and Multigas Unit (MGU-800/MGU-810 series).
- 10 module-LAN Connector Connects the HS Adapter (HSA-80), DS-8007 Adapter (DSA-82), or the Input Box (IB-8004).
- 11 LAN (TCP/IP) Connector Connects the specified equipment.
- 12 I/O Connector (Alarm) Connects the specified equipment.
- 13 Serial Connector (COM A, B) Connects the specified equipment.



- 14 Extended Display Unit Connector Connects the specified equipment.
- 15 Recorder Unit SlotA slot to insert the HR-800.

Display Unit: LC-8018TC/LC-8016TC

This section explains about the LC-8018TC (18.5 inch), LC-8016TC (15.6 inch).

Generation Front Side

1 Remote Control Sensor

Receives the infrared remote control signal.

2 Battery Charging LED

Indicates the battery-charging status. During battery operation, the LED will not light. Orange: Charging is in process Green: Charging is complete Light Off: During battery operation, or when battery is not installed, or when battery charging is ceased (due to temperature, etc.) Flash: Battery charging error



3 Power Supply LED

Indicates the power supply status.

Lights when the AC power is supplied to the main unit and links with the standby switch.

- Orange: Standby Mode
- Green: In normal operation

The light will be OFF when the AC power is not supplied to the main unit or during battery operation.

4 Standby Switch

Sets ON/OFF the standby condition.

5 Alarm Indicator

Lights/blinks when the alarm generates. Red: Level H (Urgent Alarm, Alarm Priority/High) Yellow: Level M (Cautionary Alarm, Alarm Priority/Medium) Blue: Level L (Status Alarm, Alarm Priority/Low)

6 Fixed Keys

(@"Fixed Keys" P3-1)

NOTE

• If the battery charging LED flashes, battery charging error is occurring. Remove the battery and install it again. If the error persists, contact your nearest service representative.

Rear Side

- 1 Display Unit Attaching Position Attaches to the main unit.
- 2 Mouse/Keyboard Connector Connects the optional mouse.
- 3 Main Unit Connector Connects the main unit.



• The display unit utilizes LED for the backlight. Since this LED deteriorates by the life cycle, the display may become dark, scintillate or may not light by the long term use. In such case, contact your nearest service representative.

Super Unit: HS-8000 Series

Generation Front Side

1 NIBP Start/Stop Key

Starts/stops the NIBP measurement. The indicator lights during the NIBP measurement.

2 BP Zero Balance Indicator

Performs BP zero balance.The indicator lights during the BP zero balancing.

- 3 Alarm Silence Key Silences the Alarm. The indicator lights during the alarm silence condition.
- 4 Power Supply LED Indicates the power supply status.
- 5 ECG Connector Connects the ECG cable.
- 6 AUX Connector

Connects the Gas Unit I/F (HPD-800/HPD-810), CO₂ Gas Unit (HCP-800/HCP-810), or BISx I/F Unit (HBX-800).

- 7 Multiparameter Connector Connects the input cables for BP, TEMP or CO.
- 8 NIBP Connector Connects the NIBP air hose.
- 9 SpO₂ Connector

Connects the $\ensuremath{\text{SpO}}_2$ sensor, or relay cable (patient cable).



Rear Side

- 1 HS Adapter Connector Connects the HSA-80 HS Adapter.
- 2 Analog Output Connector Outputs the ECG and BP waveforms.



HS Adapter: HSA-80

Generation Front Side

- 1 Super Unit Connector Connects the HS-8000 series.
- 2 Release Lever Releases the HS-8000 series from the HS Adapter.



Rear Side

1 module-LAN Connector Connects the Main Unit or IB-8004.



HS Adapter for DS-8400: HSA-81

Generation Front Side

- 1 Super Unit Connector Connects the HS-8000 series.
- 2 Release Lever Releases the HS-8000 series from the HSA-81.



Rear Side

1 DS I/F Connector Connects the HSA-81 to the DSC-8410.

Transport Monitor: DS-8007 Series

NOTE

- When the DS-8007 is connected to this equipment via DSA-82, the DS-8007 fixed keys (MENU, PRINT START/STOP, NIBP START/STOP) and alarm silence key are enabled on the host monitor.
- When this equipment is connected, the standby switch, HOME key, alarm indicator and alarm sound suspend key on the DS-8007 will not function.
- For the operation procedure when DS-8007 is used without connection to other system, refer to the DS-8007 system Operation Manual.

Generation Front Side

- 1 Standby Switch Sets ON/OFF the standby condition.
- 2 Power Supply LED

Indicates the power supply status. Light will be off when the AC power is not supplied to the monitor. Orange: Standby Mode Green: In normal operation Light Off: During battery operation (AC power cable is not connected.)



3 Battery Charging LED

Indicates the battery-charging status. During battery operation, the LED will not light. Orange: Charging is in process Green: Charging is complete Light Off: During battery operation, or when battery is not installed, or when battery charging is ceased (due to temperature, etc.) Flash: Battery charging error



- If the battery charging LED flashes, battery charging error is occurring. Remove the battery and install it again. If the error persists, contact your nearest service representative.
- 4 Fixed Key HOME Key MENU Key PRINT START/STOP Key NIBP START/STOP Key
- 5 Alarm Indicator, Alarm Silence Key/Alarm Sound Suspend Key Lights/blinks when the alarm generates.
 Red: Level H (Urgent Alarm, Alarm Priority/High) Yellow: Level M (Cautionary Alarm, Alarm Priority/Medium)
 Blue: Level L (Status Alarm, Alarm Priority/Low)
 By pressing this key during alarm generation, the alarm will temporarily silence.
 When "Alarm Sound Suspend" setting is ON, the alarm sound can be suspended by holding down this key.
 (Per "Detail Setup" P6-6)

Rear Side

1 CO₂ I/F Connector Attaches the Gas Unit I/F (HPD-820) or CO₂ Gas Unit (HCP-820).



Right Side

- 1 USB Memory Slot Insert a USB memory.
- 2 Battery Cover Stores the specified lithium-ion battery.
- 3 DS I/F Connector Connects to the DSA-81 AC Unit, DSA-82 DS-8007 Adapter, or DSC-8410.



Left Side

- 1 Analog Output Connector Outputs the ECG, BP, synchronized signal.
- 2 ECG Connector Connects the specified ECG relay cable.
- 3 NIBP Connector Connects the NIBP air hose.
- 4 AUX Connector Connects the Gas Unit I/F (HPD-810), CO₂ Gas Unit (HCP-810), BISx I/F Unit (HBX-800).
- 5 Temperature Connector Connects the TEMP sensor cable.
- 6 Multiparameter Connector Connects the input cables for BP, TEMP or CO.
- 7 SpO $_2$ Connector Connects the SpO $_2$ sensor, or relay cable (patient cable).



NOTE

 When operating as an input module by connecting to the rear side of the host monitor (DS-8400 system), standby switch, fixed key, alarm indicator will not function..

Bottom Side

The illustration shown on right is when the cover is removed.

- 1 SD Card Indicator Lights when the SD card is accessed.
- 2 SD Card Slot Card slot for the specified SD card



Adapter for DS-8007: DSA-82

1 DS I/F Connector

Connects the DSA-82 to the DS-8007.

2 DS-8007/HS-8000 Switch

This is a switch to change between DS-8007 and HS-8000.The selected unit and the host monitor will communicate.The LED for the selected unit will light in green.Switching to the left will communicate with the DS-8007.Switching to the right will communicate with the HS-8000.

3 Lock Lever

Lever to release the DS-8007 series from the DSA-82



- 4 Host Monitor Connector Connects the host monitor.
- 5 HS-8000 Connector Connects the HS-8000 via HSA-80.



NOTE

 When the DS-8007 is connected via DSA-82, the NIBP auto mode setting of the DS-8007 will be displayed on the DS-8007, but the actual measurement will be performed according to the NIBP auto mode setting of the DS-8400.

Multi Module: HM-800

Generation Front Side

- 1 Power Supply Indicator Indicates the power status.
- 2 BP Zero Balance Indicator Lights during BP zero balancing.
- 3 BP Zero Balance Key Starts BP zero balance.
- 4 Multiparameter Connector Connects the relay cables for BP, TEMP or CO.
- 5 Release Lock Button Press to lock the release lever.
- 6 Release Lever Press here to remove the expansion modules from the Input Box.

Rear Side

- 1 Power Input Connector Supplies power while connecting to the Input Box.
- 2 Infrared Communication Port Communicates with the Input Box via IrDA.





HM-801 Multi Module

Generation Front Side

- 1 Power Supply LED Indicates the power ON/OFF status.
- 2 BP Zero Balance Indicator Lights during BP zero balancing.
- 3 BP Zero Balance Key Starts BP zero balance.
- 4 Multiparameter Connector Connects the input cables for BP, TEMP or CO.
- 5 AUX Connector

Connects the Gas Unit I/F (HPD-810), CO_2 Gas Unit (HCP-810), BISx I/F Unit (HBX-800).

- 6 Release Lock Button Press to lock the release lever.
- 7 Release Lever

Press here to remove the expansion modules from the input box.

Rear Side

- 1 Power Input Connector Supplies power while connecting to the input box.
- 2 Infrared Communication Port Communicates with the input box via IrDA.



1

2 3


SpO₂ Module (HG-810/HG-820)

There are HG-810 (Masimo model) and HG-820 (Nellcor model) for the SpO_2 modules. The following shows the example of HG-810.

Generation Front Side

- 1 Power Supply LED Indicates the power ON/OFF status.
- 2 SpO₂ Connector
 Connects the SpO₂ sensor, or relay cable (patient cable).
- 3 Release Lock Button Press to lock the release lever.
- 4 Release Lever Press here to remove the expansion modules from the input box.



Rear Side

- 1 Power Input Connector Supplies power while connecting to the input box.
- 2 Infrared Communication Port Communicates with the input box via IrDA.



Multiport Module: HP-800

Generation Front Side

- 1 Power Supply Indicator Indicates the power status.
- 2 Status Input/Output Connector

Performs serial communication with the external device, and inputs the alarm status of the external device.

- 3 Analog Input Connector Inputs analog signal of the external device.
- 4 Release Lock Button Press to lock the release lever.
- 5 Release Lever

Press here to remove the expansion modules from the Input Box.



Rear Side

- 1 Power Input Connector Supplies power while connecting to the Input Box.
- 2 Infrared Communication Port Communicates with the Input Box via IrDA.



Recorder Unit: HR-800

Generation Front Side

- 1 Power Supply Indicator Indicates the power status.
- 2 Printing Indicator Lights during printing.
- 3 Print Key Starts/stops the printing.
- 4 Paper Feed Indicator Lights during paper feeding.
- 5 Paper Feed Key Feeds the paper.
- 6 Open/Close Lever Press to open the paper holder.

Rear Side

1 U-LINK Connector

Connects to the MGU-800/810 series Multigas Unit or Main Unit.





Input Box: IB-8004

Generation Front Side

- 1 LAN-ID Setting Indicator Indicates the assigned LAN-ID.
- 2 Power Output Connector Supplies power to the expansion module.
- 3 Infrared Communication Port Communicates with the expansion module via IrDA.
- 4 Expansion Module Connection Slot Connects maximum of 4 expansion modules to slot 1 to 4.

Rear Side

- 1 LAN ID Setting Dial Sets the LAN ID to 1.
- 2 module-LAN Connector x3 Connects the Main Unit, IB-8004, or HSA-80 HS Adapter.





REFERENCE

· For the connection procedure, refer to the operation manual of the IB-8004.

CO₂ Gas Unit: HCP-800

Generation Front Side

1 Power Supply LED

Indicates the power ON/OFF status. It will light in green while the power is ON.

- 2 Sampling Tube Connector Connects the sampling tube manufactured by Covidien.
- 3 Clip

Attaches to the bedside rail or headboard for bedside use.



Rear Side

- 1 AUX connection cable Connects to the AUX connector of the HS-8000.
- 2 Exhaust Hole

Connects the gas exhaust system and exhausts sampling gas.



• Do not block the exhaust hole as it may cause damage to the equipment.

CO₂ Gas Unit: HCP-810

Generation Front Side

- 1 Power Supply LED Indicates the power ON/OFF status. It will light in green while the power is ON.
- 2 Sampling Tube Connector Connects the sampling tube manufactured by Covidien.
- 3 Clip

Attaches to the bedside rail or headboard for bedside use.

Rear Side

1 AUX Connector

Connects to the AUX connector of HS-8000/DS-8007/HM-801 with AUX connection cable.

2 Exhaust Hole

Connects the gas exhaust system and exhausts sampling gas.



NOTE

 The usable AUX connection cable differs depending on the connecting equipment. For the combination of the AUX connection cable and the connecting equipment, refer to the section on "Optional Accessories".

Do not block the exhaust hole as it may cause damage to the equipment.

CO₂ Gas Unit: HCP-820

Generation Front Side

1 Sampling Tube Connector Connects the sampling tube manufactured by Covidien.



Rear Side

- 1 Gas Unit Connector Connects to the CO₂ I/F connector of the DS-8007.
- 2 Exhaust Hole Connects the gas exhaust system and exhausts sampling gas.



· Do not block the exhaust hole as it may cause damage to the equipment.

Gas Unit I/F: HPD-800

Generation Front Side

1 Power Supply LED

Indicates the power ON/OFF status. It will light in green while the power is ON.

- 2 CO₂ Connector Connects to the Capnostat 5 (Philips).
- 3 Clip

Attaches to the bedside rail or headboard for bedside use.



Rear Side

1 AUX Connection Cable Connects to the AUX connector of the HS-8000.



Gas Unit I/F: HPD-810

Generation Front Side

1 Power Supply LED

Indicates the power ON/OFF status. It will light in green while the power is ON.

2 CO₂ Connector

Connects to the Capnostat 5 (Philips).

3 Clip

Attaches to the bedside rail or headboard for bedside use.



Rear Side

1 AUX Connector

Connects to the AUX connector of HS-8000/DS-8007/HM-801 with AUX connection cable.



NOTE

 The usable AUX connection cable differs depending on the connecting equipment. For the combination of the AUX connection cable and the connecting equipment, refer to the section on "Optional Accessories".

Gas Unit I/F: HPD-820

Generation Front Side

1 CO₂ Connector Connects to the Capnostat 5 (Philips).

Rear Side

1 Gas Unit Connector Connects to the CO_2 I/F connector on the DS-8007.



BISx I/F Unit: HBX-800

Generation Front Side

- 1 BISx Module Connector Connects the BISx (Covidien).
- 2 Power Supply LED

Indicates the power ON/OFF status. Lights in orange when the power is ON, and BISx is not connected.Lights in green when the BISx is connected.





Rear Side

1 AUX Connector

Connects to the AUX connector of HS-8000/DS-8007 with AUX connection cable.

NOTE

• The usable AUX connection cable differs depending on the connecting equipment. For the

combination of the AUX connection cable and the connecting equipment, refer to the section on "Optional Accessories".

Multigas Unit: MGU-800/810 Series

Front Side (MGU-801P/MGU-811P)

- 1 Power Supply LED Indicates the power ON/OFF status.
- 2 Inhale Port Connects the sampling tube to inhale sampling gas.
- 3 Water Trap with Reservoir

Removes the water from the sampling tube connected to the patient. When the reservoir is more than half full with water, empty the water. (@Maintenance Manual "Water Trap (Multigas Unit)" P8-5)



4 Exhaust Hole

Connects gas exhaust system and exhausts sampling gas.

WARNING

• To protect the hospital staffs from unnecessary anesthetic agent, it is strongly recommended to connect the exhaust hole to the gas exhaust system in the hospital.

- Do not block the exhaust hole as it may cause damage to the equipment.
- 5 Flow Sensor Connector (For MGU-810 series only) Connects the flow sensor cable.

Rear Side

- 1 External Equipment Connector 1 Connects the main unit.
- 2 External Equipment Connector 2 Connects the Recorder Unit (HR-800).



Chapter 3 Operation Procedure and Screen Examples

Operation Procedure	3-1
Fixed Keys	3-1
Touch Key	
Mouse	3-2
Home Display	3-3
About the Home Display	
Displayed Items	
Description of the Display	3-16
Messages and Sound	
Window Display	
About the Window Display	
Display	3-19
Floating Window Screen Display	3-21
Minimize Window	3-22
Transfer Window	3-23
Operation Restriction	3-24
Procedure to Return the Display	3-24
To Enter Characters	3-25
For Easier Use	3-25
User Key	3-26
Menu Screen	
To Delete the Unnecessary Keys (Key Mask)	3-27
Display on the External Monitor and Extended Display Unit .	
External Monitor Display	
Extended Display 1	
Extended Display 2	

2

3

Chapter 3 Operation Procedure and Screen Examples

Operation Procedure

All operation of this equipment is performed using the fixed keys, touch keys, and mouse (optional). The equipment can also be operated remotely using the optional remote control.

Δ

Fixed Keys

1 Alarm Silence Key

Silences the alarm. The key LED will light during the alarm silence condition.

2 NIBP Start/Stop Key

Starts/stops the NIBP measurement. If pressed during the measurement, the measurement will cease. The key LED will light during the NIBP measurement.

3 NIBP Auto Mode Key

The NIBP Auto Mode setup window will be displayed. The key LED will light when the NIBP Auto Mode is set.

4 Standby Switch

Enters into standby mode or cancels the standby mode.

Touch Key

General Key Control



- 1 Pressing the [Menu] key will switch the display with a pip sound.
- 2 The touch key will respond by pressing any part of the key.
- 3 Pressing the [Home] key at any time will return the display to the home display.

REFERENCE

- The above is an example of the screen. The user keys can be customized and can be placed to any position.
 - (☞ "To Configure the Display" P10-4)

Key Control for Each Parameter



- 1 Press the numeric data box area. The touch key will respond by pressing any part of the numeric data box.
- 2 Pressing the [Home] key at any time will return the display to the home display.

REFERENCE

Frequently used touch keys can be programmed as user key. The user keys can be positioned to the user key display area and also to the numeric data area.
 (@"For Easier Use" P3-25)

Mouse

An optional mouse can be connected allowing touch key control using the mouse.

By moving the pointer on the displayed keys, and left-clicking the mouse, the operation can be performed just the same as by directly touching the displayed keys.



The pointer will be hidden if the mouse is not used for 5 minutes.(default operation) The hidden mouse pointer will be displayed again by moving the mouse.

NOTE

- It is necessary to set the mouse function (ON/OFF, pointer shape, moving speed) in advance. (maintenance Manual "Operation Related Setup" P5-26)
- · Some mouse operation such as dragging operation is not possible on this equipment.

Home **Display**

About the Home Display

The display can be configured according to the monitoring purpose.

The display layout can be selected from the following combinations of display location and size of numeric data box.

Display Location of Numeric Data	Numeric Data Box Size	Features
Right/Left	Right/Left: 1 column/ 2 columns By selecting 2 columns, more numeric data can b	
Right/Left+Bottom	Bottom: 1 row/2 rows	monitored. Also, 12-lead ECG can be displayed on the waveform display area.
Bottom	2 rows to 6 rows	The waveform display area or numeric data can be enlarged.

By connecting maximum of 2 extended display units, maximum of 2 displays independent from the main display can be used. (extended display function)

Also, the user key location can be selected from left, right, or bottom. By selecting the layout of 2 columns on left/ right, more user keys can be displayed.

REFERENCE

• The display layout can be configured/registered according to the monitoring waveform and numeric data as necessary.

(${}_{\bigcirc}$ "To Configure the Display" P10-4)

Display Example







12-Lead (Box Layout: Right)



On this system, 9 main modes and 6 sub modes can be preprogrammed according to the monitoring purpose. By preprogramming the configuration to each mode, the display configuration setups at admittance of patient can be simplified by just selecting one of the modes.

("To Select the User Mode" P5-8)

It is recommended to program the mode in rough classification such as patient's condition, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient.

Oxygenator Mode

WARNING

- The oxygenator mode is intended to prevent alarms during cardiopulmonary bypass surgery. Pay special attention when using this mode as the alarm generation will not be the same as to the standard monitoring mode.
- If the "Alarm Setting" under the Oxygenator Mode Setup is set to [All OFF], all vital alarm will
 not generate regardless of the alarm setting of each parameter. Also, if [Sel. Parameter] is
 set, vital alarm for unselected parameter will not generate. Pay attention to not miss any
 significant change of the patient's vital sign as the alarms will not be generated during the
 Oxygenator Mode.
- Once the cardiopulmonary bypass is finished, make sure to cancel the Oxygenator Mode and return to the standard monitoring mode.

• If the NIBP alarm is turned OFF under the Oxygenator Mode, NIBP auto mode measurement and NIBP measurement at alarm occurrence will not be performed.

Oxygenator mode can be used to prevent frequent alarm generation when oxygenator is used for extracorporeal circulation during cardiac surgery.

During oxygenator mode, "Oxygenator Mode" will be displayed on the screen, alarm generation will be stopped, and low priority parameter will be displayed with decreased brightness.

The main difference	of standard monitor	oring mode and	oxygenator mod	e is as follows.
		0	10	

	Standard Monitoring Mode	Oxygenator Mode	
Vital Alarm	will be generated.	will not be generated, or only the alarm for specified parameter will be generated. [*]	
Equipment Status Alarm	will be generated. will be generated for specified parameter.		
NIBP Periodic Measurement	······································		
Night Mode Night mode can be used. Night mode cannot be used.			
*It is also possible to set the same alarm function with the standard monitoring mode.			

REFERENCE

 The oxygenator mode setup can be performed under [Menu>Initial Settings>Alarm>Oxygenator Mode Setup].
 (@Maintenance Manual "Alarm Related Setup" P5-4)

- **1** Press the [Oxygenator] key on the user key.
 - The confirmation window will be displayed.

Press the [OK] key to change the monitoring mode to oxygenator mode.





Press the [Standard Mode] inside the message window to return to the standard monitoring mode.

REFERENCE

- The message window can be dragged to any position within the waveform area.
- The message window will not be displayed by selecting [OFF] for "Oxygenator Mode Message" (Menu>Initial Settings>Oxygenator Mode Setup).

Displayed Items

Other than waveforms and numeric data, patient name, alarm message, status message, etc. will be displayed on the screen.

Numeric Data, Waveform, Patient Name, etc.

- Information Display Area Room/Bed ID, Patient Name, Patient Class., current time, messages, etc., will be displayed.
- 2 Numeric Data Area
- 3 User Key Area
- 4 Waveform Area



Information Display Area



1 Telemetry Channel (When HLX-801 (FA) / HLX-801 (G) is connected) Displays the telemetry channel ID.

2 Room/Bed ID

Displays the 4-digit Room ID and 3-digit (000-999) Bed ID.

3 Nurse Team Color

Displays the color of the nurse team set on the "Admit/Discharge" menu.

4 Patient Name

The patient name set on the "Admit/Discharge" menu will be displayed. The patient name can be hidden from the display area by selecting [OFF] for "Patient Name on the Information Display Area" (Initial Settings>User I/F>Display/Print).

5 Pacemaker Usage

When [Used] is set for "Pacemaker" on the "Admit/Discharge" menu, <Pacemaker> will be displayed.

6 Patient Classification

The patient classification (Adult, Child, Neonate) set on the "Admit/Discharge" menu will be displayed.

7 Date/Time

Displays the current date (month, day) and time (hour, minute).

8 Set Mode

The currently selected user mode will be displayed. Sub mode will be also displayed if selected.

When using the data transfer function with the Super Unit, alarm settings and parameter settings can be also transferred. When the settings are changed by the data transfer function, the mode name will be highlighted to notify that the setting has been changed. Pressing the highlighted mode name will display the confirmation message window (shown on right). Pressing the [OK] key will clear the highlight.



When the alarm settings are changed, the alarm settings list will be displayed.

9 Ventilator Connection Status

Displays the connection status of the ventilator. <Vent. Comm.>: Communication with the ventilator is in progress. <Vent. Offline>: Communication with the ventilator is interrupted. <Vent. Disable.>: Communication with the ventilator is disabled. No display: Ventilator is not set for "External Device" setting.

10 Drift Filter

When drift filter is set to ON, <Drift-F ON> will be displayed.

11 Message Area

When an alarm generates, a message will be displayed.

By pressing the message display area, the alarm message history can be verified.

12 Other Bed Status

Displayed when connected to central monitor. Pressing the [Other Bed] key will display the Other Bed display.

13 Power Supply/Battery

The power supply status and battery charging status (when battery is installed) will be displayed.

Waveform Area

- 1 ECG
- 2 ECG Lead
- 3 ECG Size

The waveform size display of ECG, RESP, SpO₂ can be selected from [Numeric]/[Bar]/[Bar (10 mm)]. [Initial Settings > User I/F > Display/Print > Waveform Size Display]

(Plaintenance Manual "Display/Print Setup" P5-13)

- 4 SpO₂ Waveform
- 5 SpO₂ Size
- 6 BP Label
- 7 BP Scale
- 8 BP Waveform
- 9 Respiration Waveform
- 10 RESP Size
- 11 AWF, AWP, AWV Waveform and Scale
- 12 CO₂ Scale
- 13 CO₂ Waveform
- 14 EEG Waveform
- 15 EEG Scale
- 16 ECG Drift Filter/AC Filter Display
 - AC: AC Filter ON, DF: Drift Filter ON M: Monitor Mode, E: ESIS Mode, D: Diagnosis Mode





Enlarged Waveform

By selecting [ON] for "Zoom" under [Menu>Basic Setup>Display Config.], the waveform display can be enlarged. The waveform display can be enlarged for ECG, pulse wave, and respiration waveforms. <200%> will be displayed for the enlarged waveform. Also, the sweep speed will be doubled.



Graphic/Tabular Trend Display

By selecting [ON] for "Graphic/Tabular Trend" under [Menu>Basic Setup>Display Config.>Detail Setup], graphic/tabular trend can be displayed below the waveform display area. It cannot be displayed if "Bottom" is set for the display layout.

(Review Function" P8-1)



Numeric Data Box Display (for all parameters)

1 Message Icon

When the numeric data box size is too small to display the message inside, a message icon will be displayed instead to indicate that message is present.

(@Maintenance Manual "Display/Print Setup" P5-13)

2 Alarm OFF Mark

Indicates that the alarm is set to OFF.

- 3 Alarm Silence Mark Indicates that the alarm is silenced.
- 4 Out of Measurement Range (XXX) Indicates that the measurement is out of range.
- 5 Invalid (---)

Indicates that the NIBP measurement ended erroneously, or NIBP erase time has elapsed.

UNumeric Data Box Display (for each parameter)

|--|

 The following numeric data box is displayed when the corresponding parameter is selected on the "Numeric Data Selection" window under "Display Config.". (@"Numeric Data Selection" P10-2)

HR, HR/PR

1 HR/PR Synchronization Mark

When HR or PR according to the setting of "Synchronized Mark/ Tone" is detected, HR/PR synchronized mark will be displayed inside the corresponding numeric data box.

2 HR/PR Value

The HR/PR value will be displayed. When the value exceeds the measurable range, "xxx" will be displayed.

3 HR Average (Instant / Average, or no display)

Displays the averaging method of HR. ("HR Average" setting on ECG setup.) When the patient classification is [Adult] or [Child], and "HR Delay" is set to [ON], "Inst." or "Av." will not be displayed.







PR, HR/PR

- 1 Pulse Rate (BP)
- 2 Pulse Rate (SpO₂)
- 3 PR_IBP Source

SpO_2

1 SpO₂ Value

The arterial oxygen saturation will be displayed.

- 3 Second Alarm Indicator (Nellcor Only)

When the second alarm is set, the second alarm indicator is displayed.

- 4 Pulse Rate The pulse rate is displayed. When the value exceeds the measurable range, "xxx" will be displayed.
- 5 PI Value (Masimo only) The perfusion index will be displayed.
- 6 PVI Value (Masimo only, optional) The pleth variability index will be displayed.
- 7 SpCO Value (DS-8007M only, optional) The arterial oxygen content will be displayed.

SpCO/SpMet/SpHb (Masimo only, optional)

- 1 SpCO Value: The carboxyhemoglobin concentration will be displayed.
- 2 SpMet Value: The methemoglobin concentration will be displayed.
- 3 SpHb Value: The total hemoglobin concentration will be displayed.

VPC

1 VPC (1 min)

The VPC rate for the last 1 minute will be displayed. "---" will be displayed during arrhythmia learning.

2 Pace Beats (1 minute) / Total Beats (1 minute)

Pace beats and total beats for the last 1 minute will be displayed. <---> will be displayed during arrhythmia learning.

ST

ST Level

The ST value for 4 leads can be displayed in the ST data box. 3 groups (A, B, C) of lead combination can be programmed. For the following case, "---" will be displayed.

- •During arrhythmia learning
- During lead-off condition







ST (mm)	Ι	×	0.5
(mm)	II	×	0.2
	III	\boxtimes	
	a¥R	\boxtimes	

- •When "N" or "S" is not detected for QRS within 30 seconds.
- •When reference waveform is not set for ST measurement.

	 REFERENCE The leads displayed inside the ST level box can be (A maintenance Manual "Display/Print Setup" P 		
1 Impeda	(Impedance Respiration) ance Detection Lead et detection lead (I/II) will be displayed.	1 RR(i-I) 2 X	30. ³

2 RR Synchronized Mark

When the respiration of the set RR source is detected, a synchronized mark will be displayed inside the corresponding numeric data box.

3 Respiration Rate

The impedance RR, CO₂ RR, and ventilator RR will be displayed. When the value exceeds the measurable range, "xxx" will be displayed.

When the impedance measurement is set to OFF, impedance RR will not be displayed.

NOTE

• To display the respiration rate for each respiration source (impedance/CO₂/ventilator/ SpO₂), set the corresponding numeric data box on the display configuration setup menu.

NIBP

1 NIBP Value/Cuff Pressure

The NIBP measurement value (SYS / DIA / MAP) will be displayed.

The mean NIBP display can be set to ON or OFF on the NIBP setup menu. The value will be displayed as "---" when the preprogrammed NIBP erase time has elapsed.

During measurement, a cuff pressure will be displayed.

2 Dyna Alert Message

This message will be displayed when the Dyna Alert is effective.

3 NIBP Measurement Interval

The NIBP measurement interval will be displayed.

4 Elapsed Time/Measured Time

The elapsed time or measured time will be displayed. The display can be selected under [Menu>Parameter>NIBP>Detail Setup>Time Display].

5 Oscillation Graph

The horizontal axis in the graph shows the cuff pressure, and vertical axis shows the pulse amplitude with reference to maximum pulse amplitude.

6 NIBP List

The NIBP list of the latest 3/6/9/12/18 data and measured date/ time will be displayed. The number of displaying data depends on the size of numeric data box.

Blood Pressure

1 BP Label

The label set for the blood pressure will be displayed.

2 <MEAN_WAVE>



~	- 4/	7	08:15	129/	82 (98)
6			08:10	1207	82 (95)
			08:05	129/	90 (98)



<MEAN_WAVE> is displayed when [ON] is set for "Mean Wave" under [Menu>Parameter>BP>Detail Setup].

3 Blood Pressure

The BP measurement value (SYS/DIA/MEAN) will be displayed. On the BP setup menu, the display type (S/ D/M, S/D, M) can be selected. When the value exceeds the measurable range, "xxx" will be displayed. If BP zero balance is not performed, "---" will be displayed, and if transducer is not connected, nothing will be displayed.

PAP/ IAP/ ICP

1 PCWP Value, PCWP Measured Time

When the BP label is PAP, PCWP (Pulmonary Capillary Wedge Pressure) and measured time can be displayed.

2 PDP Value

When the BP label is IAP, PDP (Peak Diastolic Pressure) of IABP can be measured. Note that Systolic Pressure (SYS) = Peak Systolic Pressure (PSP).



T2 (10) 🖄

37.2

38

T1 (c) 🛛

Tb (C)

CO2 (MGU)

(mmHg)

X

Insp

17

X

36.1

3 CPP Value

When the BP label is ICP, labeling the artery pressure as ART will allow measuring the CPP (Cerebral Perfusion Pressure). CPP = Mean Arterial Pressure – Mean Intracranial Pressure If the CPP value is negative value, or zero balance has not been performed for ICP or ART, "---" will be displayed, and if ICP or ART has not been measured, nothing will be displayed. Also, alarm cannot be set for CPP.

Temperature

1 TEMP Label

The label set for the temperature will be displayed.

2 TEMP Value

The temperature will be displayed.400 series temperature sensor can be used. When the value exceeds the measurable range, "xxx" will be displayed. When 700 is used, "---" will be displayed.

Blood Temperature

By using the thermodilution catheter for the CO measurement, blood temperature can be displayed. When the value exceeds the measurable range, "xxx" will be displayed.

EtCO₂/InspCO₂

InspCO₂ Value/EtCO₂ Value

The end-tidal CO_2 concentration and inspiratory CO_2 concentration measurement value will be displayed.

The measurement unit can be selected from mmHg / kPa / % under the "Initial Settings" menu.

Ventilator

Ventilator Data

When ventilator is connected, the ventilator measurement data will be displayed.

TV i 400 TV e 418 MV e 6.2 PEAK 2 PEEP O (emH>0) MEAN 1



P-V, F-V

P-V, F-V Loop

By connecting the ventilator, multigas unit (MGU-810 with SPIRO unit), or FLOW-i, P-V loop (airway pressure / ventilation) and F-V loop (airway flow / ventilation) can be monitored on the ventilator display.

SvO₂/CCO Measurement Device



SvO₂/CCO Data

When SvO_2/CCO Measurement Device (Vigilance/Vigilance CEDV/VigilanceII/Vigileo/PiCCO2/EV-1000/ PulsioFlex) is connected, the measurement data (SvO_2 , CO, etc.) acquired from these devices will be displayed. The displayed data will differ depending on the used SvO_2/CCO Measurement Device and measurement mode.

Oximeter/CCO Measurement Device	Displayed Data			
Vigilance (CCO mode/STAT OFF/Index OFF)	SvO ₂ (ScvO ₂)	ссо	EDV	BT
Vigilance (CCO mode/STAT OFF/Index OFF)	SvO ₂ (ScvO ₂)	CCO STAT	EDV STAT	BT
Vigilance (CCO mode/STAT OFF/Index OFF)	SvO ₂ (ScvO ₂)	CCI	EDVI	BT
Vigilance (CCO mode/STAT OFF/Index OFF)	SvO ₂ (ScvO ₂)	CCI STAT	EDVI STAT	BT
Vigilance (ICO mode)	SvO ₂ (ScvO ₂)	CO AVG	CI AVG	-
PiCCO2/PulsioFlex	ScvO ₂	CCO	CCI	BT

Hemodynamic Data

Hemodynamic Data (Vigilance)

Based on the CCO data measured by the Vigilance (or Vigilance CEDV/ Vigilance/Vigileo), the following hemodynamic data are calculated and displayed every second based on the following condition. However the following condition should be met.

- It is measured on Vigilance with CCO mode. (It will not be displayed during ICO mode.)
- SvO₂ parameter key (oximeter numeric data box) is displayed.
- BP label is set as ART, PAP, CVP.

SV	65	SVR 1363
RV₩	0.54	^{RVS₩} 8.1
S¥I	38	^{SVRI} 2304
RV₩I	0.32	RVSWI 4.2

Data	Description	Formula
SV	Stroke Volume (mL/beat)	CCO x 1000 HR
SVR	Systemic Vascular Resistance (dynes*sec*cm ⁻⁵)	(MAP-CVP) x 79.90 CCO
RVW	Right Ventricular Work (kg*m)	CCOx(MPAP-CVP)x0.0136
RVSW	Right Ventricular Stroke Work (g*m)	SVx(MPAP-CVP)x0.0136
SVI	Stroke Volume Index (mL/beat/m ²)	SV BSA
SVRI	Systemic Vascular Resistance Index (dynes*sec*cm ⁵ •m ²)	SVRxBSA
RVWI	Right Ventricular Work Index (kgm/m ²)	RVW BSA
RVSWI	Right Ventricular Stroke Work Index (g*m/m ²)	RVSW BSA

NOTE

• The hemodynamic data based on Vigilance data will not be stored in the list.For the Vigilance list, the data directly acquired from the Vigilance will be stored.

Multigas Unit Data

Multigas Unit Data

When multigas unit or mainstream module is connected, the numeric data measured by the connected unit or module (CO₂ /anesthetic gas/O₂ /N₂O concentration) will be displayed.



Stopwatch Key

Functions as stopwatch.

BIS

BIS Value

By connecting the BISx module through the HBX-800, or by connecting the BIS monitor, BIS data (BIS, SQI, EMG, SR) will be displayed.

If SQI value is below 50%, the BIS value will be displayed in gray. If SQI value is below 15%, the BIS value and SR value will disappear.

EMG and SQI will be displayed in bar graph.

Bar Graph	SQI (0 to 100) [%]	EMG (30 to 55) [dB]
1	1 to 20	30 to 38
2	21 to 40	39 to 47
3	41 to 60	48 to 55
4	61 to 80	55 and above
5	81 to 100	-

The alarm bar will be displayed only when measurement is performed on BISx using the HBX-800.

INVOS

INVOS 5100C Measurement Data

When connected to INVOS 5100C, regional cerebral oxygen saturation value will be displayed.

Lt- indicates left brain, and Rt- indicates right brain.



CO2 (mnHg)	1/	38
AGN 図 3 5/ 1 2	128	⊈
(%) 図 3 5/ 1 2	(%) 28	⊈ 30/ 0

TIMER1	00:00:00
TIMER2	00:00:00



User Selectable Numeric Data (Ventilator, Hemodynamics, Transcutaneous Blood Gas)

For the following numeric data, the data to be displayed in the numeric data box are selectable by the users.

- Ventilator Data
- Hemodynamics Data
- Transcutaneous Blood Gas Partial Pressure Data

The number of displaying data depends on the size of numeric data box.Small: 2Medium: 4Large: 6Two types of user selectable numeric data (A, B) can be set

REFERENCE

 The parameters to be displayed can be set under [Initial Settings>User I/F>Display Print] for [VENT Display Parameters] or [Hemo/etc Display Parameters].

Extended Function (Recall List)

The recall data will be displayed in a list format. By selecting a data from the list, an enlarged waveform of the corresponding data will be displayed.



Alarm Limit Display



The alarm limit can be displayed beside each numeric data.The display type can be selected from [Graph]/ [Numeric]/[OFF] for "Alarm Limit Display" under [Menu>Alarm>Detail Setup].

If ON is selected for the individual alarm, the alarm limit will be displayed.

The upper and lower limit will be displayed at upper and lower row respectively.

For BP and NIBP, each alarm limit of systolic BP (SYS), diastolic BP (DIA), mean BP (MEAN/MAP) will be displayed from the top.

ON/OFF of alarm limit display can be selected. (@"List of Alarm Settings" P6-5)

- 1 Upper Alarm Limit
- 2 Lower Alarm Limit
- 3 Current Measurement Value

NOTE

• If the alarm limit display for BP is [Graph], systolic value will be displayed.



- Depending on the numeric data box type, alarm limit may not be displayed.
- If the SpO₂ lower alarm limit is set to 85%SpO₂ or below, the alarm limit value will be displayed regardless of the "Alarm Limit Display" setting.

Short Trend Display

1 Short Trend Display

On the waveform display area, short trend can be displayed.

The parameters to be displayed for short trend can be set under [Menu>Basic Setup>Display Config.].

The short trend width can be enlarged/reduced to the pressed position on the waveform area. The short trend width can be selected from 7 levels.

The graph displayed in red indicates the alarm occurrence point. Pressing the short trend of an alarm generated parameter will display the recall screen.



2 Trend Scale

The short trend scale will be displayed between the short trend and numeric data. The displayed scale will be in accordance with the scale set on the "Trend" screen.

For the following parameters, the short trend scale can be synchronized with the corresponding waveform scale by selecting [Waveform] for "Short Trend Scale" under [Menu>Display Config.>Detail Setup].

▶ BP, PEAK, TV, CO₂, O₂, Agent

For operation procedure on the short trend display, refer to @"Short Trend" P8-10.

Displayed number of waveform and numeric data

Layout	Maximum Displayed Waveforms	Display Duration (25 mm/s)	Maximum Displayed Boxes
Right/Left, 1 column	27 (22)	13 (11)	11
Right/Left, 1 column+bottom	20	13 (11)	15
Right/Left, 2 columns	27 (22)	10 (8)	22
Right/Left, 2 columns+bottom	20	10 (8)	25
Bottom, 2 rows to 6 rows	18 to 10	16 (3)	10 to 30

*: () indicates when LC-8016TC is connected.

Description of the Display

Symbol	Description	
×	Alarm OFF Indicates the alarm is OFF.	
×	Alarm Silence Indicates the alarm is silenced when the alarm system is IEC mode.	
•	HR Synchronized Mark This mark flashes synchronizing to the heartbeat.	
Λ	RR Synchronized Mark This mark flashes synchronizing to the inspiration.	
0	Message Icon Indicates that an alarm message is present for that parameter. Whether or not to display this icon can be selected under "Initial Settings".	
Ê	Key Lock Mark Indicates that the item requires a password to change the setting.	
2 1	Key Unlocked Mark Indicates that the key is unlocked	
AC Power =	Indicates that the equipment is connected to AC power source.	
1771	Indicates the remaining battery level. This icon (full green) indicates that the battery is fully charged. *While charging, the corresponding battery level icon flashes.	
	This icon (2/3 green) indicates that the battery is less than full, but still usable.	
- Z Z	This icon (1/3 yellow) indicates that the battery level is low and needs to be charged.	
	This icon (1/3 red) indicates that the battery level is very low and flashes to alert the low battery status. Immediate battery charge is required. Technical alarm will generate.	
	This icon (red frame) indicates that the battery is almost depleted and it flashes to alert that charging is necessary. Make sure to charge the battery immediately.	
Ċ.	This icon (black frame with a slash) indicates that the battery is not installed. Pay attention as power will not be supplied if AC power cable is disconnected during this state.	

The following symbols are used for this equipment.

Messages and Sound

This section explains about the message displayed on the home display.

There are vital alarm message and equipment status alarm message which will be displayed at the top of the home display.

The alarms are classified to Level S (top priority), Level H (high priority, urgent), Level M (medium priority, caution), Level L (low priority, status), and Notification, and the message will be displayed according to the priority of Level S > Level H > Level M > Level L > Notification.

The displayed messages will flash in red and white for Level S, red for Level H, yellow for Level M, blue for Level L, and white for Notification.

Alarm Priority, Le	vel	Description	Sound	Displayed Color
Top Priority	S	Top Priority Alarm	Continuous	Red/White
High Priority	н	Life Threatening Alarm	Continuous	Red
Medium Priority	М	Cautionary Alarm	5 seconds interval	Yellow
Low Priority	L	Status Alarm	15 seconds interval	Blue
Notification	N	Notification Alarm	Display Only	White

• When more than one alarms of the same priority are generated, the newer alarm message will be displayed.

Vital Alarm Message

The vital alarm message is generated when a measurement exceeds the alarm limit, or when arrhythmia is detected.



1 Numeric Data Alarm Message

2 Arrhythmia Alarm Message

There are 2 types of vital alarm messages; numeric data alarm and arrhythmia alarm. If both alarms occur at the same time, the numeric alarm message and arrhythmia alarm message will be displayed alternately in 2-seconds intervals. The alarm messages will be displayed according to the priority. If the priority level is the same, the newer alarm message will be prioritized.



• The arrhythmia alarm message other than Tachy, Brady, Ext Tachy, Ext Brady will continue to be displayed for 30 seconds after the alarm is resolved.

Equipment Status Alarm Message

The equipment status alarm message will be displayed when proper monitoring cannot be performed. The alarm messages will be displayed according to the priority. If the priority level is the same, the newer alarm message will be prioritized.



Numeric Data Box Message

The measurement status of each parameter will be displayed inside the corresponding numeric data box.



Lead-Off Message

If the ECG electrodes used for HR measurement or arrhythmia analysis are detached, the status will be notified.

WARNING

 When <Lead-Off> is displayed, HR alarm or arrhythmia alarm will not be generated.If this condition is left unresolved, a sudden change of the patient may not be noticed. Take prompt action when the lead-off condition is detected.



□ Ventilator Alarm Message

When a ventilator is connected to this equipment, ventilator alarm and connection status alarm will be displayed on the equipment status alarm message area.

The alarm message with the higher alarm level will be displayed.

WARNING

- The ventilator alarm sound is set to OFF (factory default).
- The alarm sound can be turned ON under [Menu>Tone/Volume].
 (@"Tone/Volume" P10-20)



□Ventilator Alarm Factor Message

For the Servo ventilators (SV-300, SERVO-i/s, SERVO-U/n/air), ventilator alarm factor if specified will be notified and displayed on the central monitor.

- For the ventilators other than Servo ventilators, ventilator alarm factor will not be notified to the central monitor.
- Depending on the central monitor type and software version, ventilator alarm factor may not be displayed. For details of the central monitor type and software version, refer to your nearest service representative.
- The ventilator alarm factors are displayed only on the central monitor. These will not be displayed on the bedside monitor.

Uventilator Disconnected Confirmation Window

A confirmation window will be displayed when a ventilator cable is disconnected from the DS-8400, or when the power of the ventilator is turned OFF.

[ON] will continue communication with the ventilator during ventilator alarm condition. Check the ventilator power and cable connection.

[Alarm Silence] will silence the ventilator alarm for 2 minutes.

If the ventilator alarm condition remains after 2 minutes, the alarm will generate again.

[OFF] will disable the ventilator alarm until the ventilator connection status recovers.

Vent.	0N
Alarm	Alarm Silence 0FF
ci	Ventilator is in power off or standby condition. Check ventilator. Or, check if cable is properly pomected to DS-8400 and ventilator.

- Check occasionally the communication status of this equipment and the ventilator.
- Verify that a ventilator alarm is not generated, and that the <Vent. Comm.> message is displayed.

This confirmation window will be displayed until the displayed key is pressed or proper communication with the ventilator is resumed. When the communication is resumed, the window will automatically close.

When disconnecting the ventilator and this equipment, make sure to select [OFF] on the confirmation window which will be displayed when the power of the ventilator is turned OFF, or when the cable is disconnected.

Window Display

About the Window Display

The screens that are displayed when operating this system are referred to as windows. (The windows that appear by pressing the numeric data area are called floating windows, as they can be moved to any desired position.) The target window can be displayed by using various method, such as selecting the menu items, pressing a parameter key or using a short cut key such as user key.

Display

The common items on the window are explained below.

1 Hierarchical Level Display

The hierarchical level of the current window is displayed. The level is expressed using the ">" symbol.

This area also functions as keys, making it possible to return from the lowermost to topmost window in a one-touch operation.

2 Previous Display

Pressing this key will return the display to the previous window.

3 Upper Level Key

Returns to the upper level display.



4 Minimize Key

Pressing this key will minimize the currently displayed window and will be stored to the user key. To restore the minimized window, press the Restore key in the user key and select the window to restore.

5 Key Lock Icon

Key lock icon will be displayed for the setup item that is locked.

To unlock the setup item, enter the password.

It will return to locked condition after 30 seconds if no key operation is performed.

- 🚹 : Locked
- 🔁: Unlocked

NOTE

• The color of each key lock icon indicates its administrative level, and a higher level password must be entered to unlock it.

6 Setup Item

Most of the setups can be performed by selecting from the dropdown list. The dropdown list will close when a selection has been made.

Pressing the item again or selecting a different item will also close the dropdown list.

Some menu may display a subwindow to perform the setup.

To close the subwindow, press either the \mathbf{X} key, [Home] or [Prev. Disp.] key.

Pressing the key with the "
"
icon will display another window.To return to the original display, press the
key.

• Example of an item which displays another screen



7 Dropdown List

Arrhy. Alarm Setup

Select one from the displayed selection list.

8 Page Switch Key

This key will appear when the setup items or display data are on multiple pages.

The currently displayed page is indicated by "•".

9 Tab Display Area

These are the tabs to display the screens under the same menu level. The screens under the same menu level can be switched by one-touch operation of these tabs without returning to the main menu.

For example, to change the blood pressure scale after changing the ECG waveform size, it is not necessary to return to the main menu.

For the review screens, the date/time of each review data are linked which allows to switch the display of the tabular trend, graphic trend, waveform of the same date/time in one-touch operation.

3

Δ

Floating Window Screen Display

The descriptions of the floating window which is displayed by pressing the numeric data area are as follows. The displayed items on the floating window depends on the parameter, but there are some common items as follows.

1

2

Parameter

Size

Optimiz

×1

rhyth

Learn

Arrhy. Alarm

Lead

0

HR/PR

Auto

T Setu

িচা

HR

0 N

 $\overline{\mathbf{X}}$

Ext Tac Ext Bra

Disp. ON

1 Window Title

The windows can be moved to any desired position by dragging the window title.

2 Alarm Assist Key

The alarm assist screen will be displayed. On the alarm assist screen, maximum of 24 hours of trend data for the corresponding parameter will be displayed, and alarm threshold can be adjusted by checking the trend data.

3 Close Key

Press the \mathbf{X} key to close the window.

4 Detail Key

On the floating window, minimum items are displayed. Press the 🔄 key to display more detailed items.



Minimize Window

To temporarily display the home display during the setup, press the \square (Minimize) key. The current window will be minimized. By pressing the [Restore Window] key, the window will be redisplayed.

NOTE
• To use the restore window function, the [Restore Window] key needs to be preprogrammed as user key.

1 Press ☑.



• The window will be minimized.

 $\mathbf{2}$ Press the minimized window.

▶ The original window will be displayed again.

D-001 FUKUDA DENSHI - MHT R INITIAL	14:49 14:49
we also also also also also also also also	Home
	Nenu
n [200]	Alarm Silence
	Admit/ BP Zero Disch. 19875
	NIBP Start/Stop
	NIBP Alarm Auto Node (ALL)
	NIBP Cont. Alarn History NIBP List Recall
	Tabular Graphic Trend Trend
129/°82 °98). 116/ 77(92). 30	Print Start/Stop

NOTE

- · Maximum of 8 windows can be minimized. If exceeded, the oldest window will be deleted.
- To delete all minimized window, press the [Delete All] key which will be displayed when
 is pressed for more than 1 second.
- The window which has been automatically erased after fixed amount of time can be remained minimized by selecting [ON] for "Auto Minimize" under [Initial Settings>User I/F >Operation].

Transfer Window

The floating window displayed by pressing the numeric data area can be moved by dragging the window title. This operation is possible on the touch panel.

1 Press the window title.



2 Place the finger on the window title and drag to the desired position.



NOTE

- The floating window cannot be overlapped to the numeric data area or information display area.
- The window which is displayed from "Menu" cannot be moved.
- The displayed position of the floating window will be stored until the power is turned OFF.

Operation Restriction

To restrict the operator to change the setup items, key lock function can be used.

(B Maintenance Manual "Key Lock" P5-2)

For the items that are key locked, the settings cannot be changed unless the password is entered.

The unlocked condition will return to locked condition if operation has not been performed for about 30 seconds.

For the key locked item, **1** icon will be displayed. (shown on right)

When the password is entered and key is unlocked, the icon will change to \mathbf{f} .

NOTE

- There are 3 key lock levels.
- The level is distinguished by the color of Administrator)" > "Green (User)", and the upper level password can unlock the lower level key lock.

Procedure to Return the Display

To Return to Home Display

Press the [Home] key to return to the home display.



To Return to the Previous Display

Press the (5) key of each setup window to return to the previous display.

_		
З	-24	
3	-24	

Menu > Alarm Basic Circ. Resp./Gas Explanation Area	Arrhy. ST List Petal Gettip
Asystole ON 5 sec.	Tachy ON
VF DN	Brady DN
(HR > 120bpm) ON	Run (HR > 40 bpm) ON 3 beats
Ext Tachy OFF 150 bpm	Pause OFF 3.0 sec.
Ext Brady & OFF 30 bpm	Triplet OFF
SLOW VT ON	Couplet OFF

To Enter Characters

Alphanumeric characters and symbols can be entered using the displayed keyboard. The procedure to enter characters is explained below using the example of patient admit menu.

Entering Alphanumeric Characters

Enter alphabets, numerics, or symbols.

1 Press [ABC] or [QWERTY] to switch the displayed keyboard. Enter the alphanumeric characters.



Entering Numerics

For age, telemetry channel ID, etc., only numbers can be entered. In such case, only numeric keys will be displayed.

Enter the numbers.



For Easier Use

The user keys and menu can be customized according to the monitoring purpose.

REFERENCE

- From the preprogrammed user mode, the display configuration and alarm settings can be selected according to the monitoring purpose.
 - (@Maintenance Manual "User Mode Registration" P5-30)

User Key

The user keys can be customized according to the monitoring purpose. (@"To Configure the Display" P10-4)



By assigning the $[\diamondsuit]$ to the user key area, 2 pages of user keys can be registered, and pressing the $[\diamondsuit]$ allows to switch the pages. The user key can be enlarged by using 2 display areas.

The user key can be also assigned to the numeric data area. It is useful if the key related to numeric data is assigned near the numeric data.

Menu Screen

The key position can be changed and unnecessary keys can be deleted on the "Menu" screen. (@Maintenance Manual "Display/Print Setup" P5-13)




To Delete the Unnecessary Keys (Key Mask)

Unused keys, items, tabs can be masked. (@Maintenance Manual "Key Mask" P5-22)

Menu 🕽	> Admit/	Discharge		ۍ	Menu 🕽	> Admit/	Discharge		(ح)
				(†					
	= Require	ed Iten]			■ Require	ed Iten		(\mathbf{P})
BED- 001	Hode Select	INITIAL	■Class. Adult		BED- 001			■Class. Adult	
	=ID		Sex			=ID		Sex	
	Name	FUKUDA DENSHI	Team 📗			Name	FUKUDA DENSHI		
			Blood AB0						
			Blood Rh						
			■Pacemaker Not Used					■Pacemaker Not Used	
Bi	rth Date	Age 0	Impedance Meas, ON				Age 0	Impedance Meas. ON	
	Height (cm)	0.0 Weight 0.0 BSA 0.00 (m ²)		Nonitor Suspend					Nonitor Suspend
Admit D	ate/Time								
				Discharge					Discharge



-		
Menu >	Basic Setup	3
	Display Manual Auto Config. Minung Printing Tone/ Volume Time/Date Color	
U	Explanation Area	Ø
Basic	Printer Bedside)	
	Print Duration 24 sec.	
12-Lead	• • Print ON Position Proportion	
	Have Format Regular	
Other Setup	Graphic Printer Printing Sel. Printing Printing	
Common	• • • • • • • • • • • • • • • • • • •	

Ĺ	Sasic Setup Declay Menual Config. Menual Antenarity nee	
Basic	Printer Built-in) Waveform Select ECG1 II Have	
	Print Duration 24 sec.	
12-Lead	• c	
	Wave Format Regular	
Other Setup	Graphic Printer Printing Sel. Recall Braphic Printing Printing	
Common	•• Implication of the set of the	

Example on Tab Display

Display on the External Monitor and Extended Display Unit

For the DS-8400 system, in addition to the main display, another display unit can be used for extended display.

• Use only the specified 19-inch display unit. For details, refer to our service representative..

Model		Displayable Screen	
Model	External Monitor Display	Extended Display 1	Extended Display 2
DSC-8410	Yes	No	No
DSC-8410 (with optional CC-84)	Yes	Yes	Yes

External Monitor Display

The same monitoring display can be displayed on other display unit. However, only the messages related to arrhythmia alarm will be displayed.

Other messages, menu, and setup screens are not displayed. Also, operation is not possible on the external monitor.



Display on the Main Display

Display on the External Monitor

REFERENCE

Extended Display 1

On the extended display 1, the independent display from the main display can be displayed.



Display on the Main Display

On the extended display 1, the following operations are possible on the touch panel.

*Selection of Preprogrammed Display Layout (3 types)

*Graphic Trend Display

*Tabular Trend Display

*ON/OFF of Short Trend Display

:*ON/OFF of Enlarged Numeric Data Display

*Display Configuration Setup for Extended Display 1, 2

*Waveform Size Selection

*Alarm Silence

*Alarm Suspend

*Parameter Selection for HR/PR Numeric Data Box

- *NIBP Start/Stop
- *NIBP Auto Mode Interval Selection

Display on the Extended Display 1

- *Print Start/Stop
- *BP Zero Balance
- *Lead Selection
- *ON/OFF of Oxygenator Mode
- *Alarm Setup
- *Parameter Setup

- As a speaker is not equipped for the extended display unit, key sound and alarm sound will not be generated. Alarm sound will be generated from the main display.
- The same setup window cannot be opened simultaneously on the main display and the extended display. If the same setup window is opened, the previously opened window will close.

REFERENCE

- The waveform size for the main display and the extended display is independent. On the "Initial Settings" menu, whether or not to synchronize the waveform size/scale of extended display with the main unit can be selected.
 (Paintenance Manual "Display/Print Setup" P5-13)
- The waveform size for the extended display 1 and 2 is common.
- By setting the [Scale (Extended Display)] key as user key, the waveform scale on the extended display can be changed on the main unit.
 (@"User Key Selection" P10-15)
- The display configuration for the extended display 1 can be also changed on the main display. (Menu > Initial Settings > User Mode Regist. > Select mode for "Extended Display 1" > Change the setting)

• The key mask function, key lock function, auto hide function cannot be used.

Extended Display 2

On the extended display 2, the independent display from the main display can be displayed.



On the extended display 2, the following operations are possible on the touch panel.

- Selection of Preprogrammed Display Layout (3 types)
- ON/OFF of Short Trend Display
- ON/OFF of Enlarged Numeric Data Display
- Waveform Size Selection
- Alarm Silence
- Alarm Suspend
- NIBP Start/Stop
- Print Start/Stop
- Parameter Selection for HR/PR Numeric Data Box
- BP Zero Balance

• As a speaker is not equipped for the extended display unit, key sound and alarm sound will not be generated. Alarm sound will be generated from the main display.

(REFERENCE

- On the extended display unit 2, operation such as alarm setup is not possible.
- The waveform size for the main display and the extended display is independent. On the "Initial Settings" menu, whether or not to synchronize the waveform size/scale of extended display with the main unit can be selected.
 (maintenance Manual "Display/Print Setup" P5-13)
- The waveform size for the extended display 1 and 2 is common.
- By setting the [Scale (Extended Display)] key as user key, the waveform scale on the extended display can be changed on the main unit.
 (@"User Key Selection" P10-15)
- The display configuration for the extended display 2 can be changed on the main display. (Menu > Initial Settings > User Mode Regist. > Select mode for "Extended Display 2" >

Change the setting)

Chapter 4 Preparation

Daily Check	4-1
To Start Monitoring	4-1
Check Discharge When Start Monitoring a New Patient	4-2
Data Transfer Function Using the Super Unit	4-3
Data Transfer Function Using the Transport Monitor (DS-80	07)4-4
To Stop Monitoring	4-10
Clock Setup	4-11
Installing the Recording Paper	4-12

Chapter 4 Preparation

Daily Check

Before using the equipment, perform the daily check.

Take necessary measures for the items with the "NG" judgment, and use the equipment only if the judgments for all the items are "OK".

To Start Monitoring

This section explains about the procedure to turn the power ON and start monitoring.

- If not using the equipment for a long period, disconnect the power cable and lithium-ion battery.
- During transportation, firmly grasp the handle and make sure that the equipment does not fall. Otherwise, it may cause injury to the operator or damage to the equipment.

1 If operating with AC power supply, verify that the power supply cable is properly connected to the main unit.

If operating with battery, verify that the lithium-ion battery

(BTO-005) is properly installed in the main unit.

(☞ Maintenance Manual "Power Connection of the Main Unit" P1-18)

(@Maintenance Manual "Installing the Battery Pack (BTO-005)" P1-21)

➤ When connected to the AC power source with battery installed, charging will automatically start.

1 Rapid Charge (when the equipment is not in operation): 2.5 hours

2 Normal Charge (when the equipment is operating): 5 hours



WARNING

• Do not connect a battery other than the lithium-ion battery (BTO-005).

 $\mathbf{2}$ Turn ON the standby switch on the main unit.

- The system will turn ON and monitoring will start.
- The power supply LED on the front side of the main unit will light.
 1 Power Supply LED

Green: Power ON Orange: Standby Mode Light Off: During battery operation 2 Battery Charging LED Green: Charging is complete Orange: Charging is in process



Light Off: During battery operation, or when battery is not installed, or when battery charging is ceased (due to temperature, etc.)

Flash: Battery Charging Error

NOTE

- The operation after the power is turned ON will be according to the setting made on [Initial Settings] > [User I/F] > [Power ON/Discharge]. However, if the power was turned OFF for less than 30 seconds, the setting before the power was turned OFF will remain.
- If the battery charging LED flashes, battery charging error is occurring. Remove the battery and install it again. If the error persists, contact your nearest service representative.

REFERENCE

 The power ON/OFF operation of the main unit, Super Unit, expansion unit, expansion module, and input box synchronizes with the standby switch operation (ON/OFF) on the display unit.

Check Discharge When Start Monitoring a New Patient

The trend data, tabular trend data, recall, ST measurement, OCRG data will be stored for 5 minutes even after the standby switch is turned OFF. If the previous data is remained when the standby switch is turned ON again, the discharge confirmation screen will be displayed.

	Dis	scharge		
Patient da	ta/info., nonitor	ing paraneters,	etc. will be i	initialized.
	Cont i nue	monitoring.		7
Monitoring	will continue.			2

Check Discharge

Select from [Discharge] / [Continue].

- [Discharge]: The previous data will be deleted.
- ▶ [Continue]: The monitoring will start with the previous data retained.

NOTE

- If the standby switch was turned OFF for less than 30 seconds, this discharge confirmation screen will not be displayed. To perform the discharge procedure, press the [Discharge] key on the "Admit/Discharge" screen.
 (P5-6)
- To start monitoring a new patient, select [Discharge] and enter the new patient information on the "Admit/Discharge" screen.

REFERENCE

Whether or not to display the discharge confirmation screen can be selected.
 (@ Maintenance Manual "Power ON/Discharge" P5-18)

Periodic Replacement Message

When the periodic replacement period approaches for each part, a message will be displayed to notify the user.

Continue			
Periodic Check	Replace the following		
HCP-800 Calibration	NIBP Unit		

REFERENCE

- The parts which the replacement period will be notified are the NIBP unit in the Super Unit and the CO₂ unit in the HCP-800/HCP-810/HCP-820.
 (Priodic Replacement" P7-1)
- Even if it is set not to display the discharge confirmation screen, the confirmation message for parts replacement will be displayed when the replacement period approaches.

Data Transfer Function Using the Super Unit

The patient data and settings are stored in the HS-8000 Series Super Unit.

When transferring the patient to another bed, the same monitoring condition such as patient data and settings can be used on new bed by transferring the HS-8000 along with the patient.

This function can be used when the HS-8000 series Super Unit is connected. When the DS-8007 is connected, refer to provide the Transfer Function Using the Transport Monitor (DS-8007)" P4-4, as the operation will differ.

1 Turn OFF the power of the DS-8400 system.

 $\mathbf 2$ Connect the HS-8000 to the DS-8400 system of the new bed, and turn ON the power of the DS-8400 system.

- > The "Patient Selection" window will be displayed.
 - [Continue monitoring data of ***]: The patient data will be transferred and monitoring will resume.
 - 2 This will be displayed when the data is damaged and cannot be transferred.
 - 3 The data that can be transferred will be displayed. To change the transferring data, press the [Change Data] and change the setting.

tient Selection	1
▲Depending on the selection, patient data on the bedside monitor and/or HS-8000 will be deleted.	
Patient data was found on the HS-8000.	1
Continue monitoring data	<u> </u>
of .	
*The following data is corrupted and cannot be updated.	<u>+</u> 2
Part of Recall Data,	
The folloving data will be updated with the HS-8000 data.(to 2000/01/01 00:00)	<u>+3</u>
Patient Admit/Discharge Data, Trend Data, Initial Settings,	
Change Data	4
Discharge	5
	-
*Patient data/info., nonitoring parameters, etc. will be initialized.	
	6
Cancel	<u> </u>
*The HS-8000 patient data vill be deleted.	1
	1

- 4 [Change Data] : The data to be transferred can be changed.
- 5 [Discharge]: The data will not be transferred and monitoring of new patient will start.
- 6 [Cancel]: The patient data and settings stored on the main unit will be used.

· After the data transfer process, make sure that the setting and patient data are correct.

NOTE

- During the data update process, the patient name on the home display will flash.
- When [Continue monitoring] is selected, the stored data on the main unit will be overwritten with that of the HS-8000.
 If a central monitor is connected, the data on the central monitor will be also deleted.
 The alarm settings and parameter settings can be also transferred. When the settings are changed by the data transfer function, the mode name will be highlighted to notify that the setting has been changed. Pressing the highlighted mode name will display the confirmation message window, and pressing the [OK] key will clear the highlight. When the alarm settings are changed, the alarm settings list will be displayed.
- When [Discharge] is selected, both data on the main unit and the HS-8000 will be deleted/ initialized.
- When [Cancel] is selected, the stored data on the HS-8000 will be overwritten with that of the main unit.
- The data on the HS-8000 will be updated if any of the [Continue monitoring]/[Discharge]/ [Cancel] is selected. Do not disconnect the HS-8000 during the update process. If disconnected, the data consistency may be lost.
- The BP zero balance value on the HS-8000 will not be cleared. After transferring the data, make sure to verify the BP zero balance value.
- The recall event generated during the data update process will not be stored.
- If the time setting is different between the data transferring monitors, the time of the recall
 data and trend data may not be correctly displayed on the monitor which the data was
 transferred.
- Do not disconnect the HS-8000 while setting up the extended display.

REFERENCE

 ON/OFF of data transfer function and the data selection to be transferred can be performed on the "Initial Settings" menu.

(@Maintenance Manual "System Setup" P5-28)

Data Transfer Function Using the Transport Monitor (DS-8007)

The DS-8007 system can be used as a transport monitor by using the stored patient data and settings. For the monitoring system using the DS-8007, DS-8400, DS-8500, DS-8900, the data during transport (when not connected to the central monitor) can be transferred between the monitors using the transport monitor (DS-8007).

 This function can be used only when the monitoring system is constructed with the DS-8007, DS-8400, DS-8500, DS-8900. To use this function, refer also to the DS-8900 Operation Manual. <General Description of Data Transfer Function>



- While Monitoring on the Bedside Monitor Remove the DS-8007 from the DS-8400/DS-8500.
 When the confirmation window is displayed, select [Monitor Suspend] or [Discharge].
- 2 While Leaving the Bed During transfer or examination, monitoring on the DS-8007 will continue. The monitoring data will be saved on the DS-8007.
 - Resume Monitoring on the Bedside Monitor
 When the patient returns to bed, attach the DS-8007 to the DS-8400/DS-8500.
 By selecting [Monitor Patient of the Transport Monitor] on the patient selection window, the saved data on the DS-8007 while in transfer will be uploaded to the DS-8400/DS-8500.
 - 4 Uploading to the Central Monitor The data will be automatically uploaded to the central monitor.

Condition to Use the Data Transfer Function

To use the data transfer function, the following conditions need to be satisfied.

- On the DS-8007, full disclosure waveform recording on the SD card is required.
- The software version of the DS-8007 should be V03-01 and newer.
- On the DS-8400/DS-8500, the data transfer function needs to be enabled. Select [Transport] for "Data Transfer" under [Menu > Initial Settings > System > Other].
- To transfer the DS-8007 data (full disclosure waveform, trend, recall) to the DS-8400/DS-8500, full disclosure waveform recording on the CFast Card is required.
- To transfer the setup data (alarm, parameter), the data transfer function needs to be enabled and the setup data to be transferred needs to be selected on the DS-8400/DS-8500.
 Select [ON] for "Alarm Setup" and "Parameters" under [Menu > Initial Settings > System > Other > Data Selection for Transfer].

For uploading to the central monitor, the following additional conditions needs to be satisfied.

- The central monitor needs to be compatible with the data transfer function.
- The bed is registered on the central monitor.

When the DS-8007 is disconnected from the DS-8400

When the DS-8007 is disconnected from the DS-8400, a confirmation window will be automatically displayed.

- [Monitor Suspend]: The confirmation window will close, and monitoring will be suspended.
 This is to be selected when returning to the same bed after transferring.
- ► [Discharge]: The patient will be discharged on the DS-8400, DS-8900. The monitoring condition after discharge will be according to the "Discharge Mode" setting (

discharge will be according to the "Discharge Mode" setting. (@Maintenance Manual "Power ON/ Discharge" P5-18)

This is to be selected when returning to the other bed after transferring.

• [Close]: A confirmation window will close.

When the DS-8007 is connected to the DS-8400

The operation will differ depending on the patient information (patient name, patient ID, patient identification code) saved separately on the DS-8400 and DS-8007.

NOTE

- The patient identification code will be updated when the patient is discharged. If the patient information such as patient name, ID is changed without performing the discharge process, the patient identification code will not be updated, and the patient will be treated as the same patient.
- When the EMR link function is used, the patient selection window shown below will not be displayed.

When the patient information matches between the DS-8400 and DS-8007

When the patient information matches between the DS-8400 and DS-8007, the patient selection window will not be displayed, and monitoring will continue.

2 When the patient information does not match between the DS-8400 and DS-8007

- The patient selection window will be automatically displayed. Select the monitoring patient.
- [Monitor Patient of the Transport Monitor] Starts monitoring the patient of the DS-8007 by discharging the patient of the DS-8400.
 All data on the DS-8400 will be deleted. If a central monitor is connected, the data on the central monitor will be also deleted.Part of the patient information, review data will be overwritten with the data of the DS-8007.

L	Wonitor the patient of the transport monitor
L	Patient ID:
	Patient Name:
	Deletes the patient review data of this equipment, and starts monitoring the patient of the transport monitor.
ſ	Monitor the patient of this equipment
L	Patient ID:
	Patient Name: FUKUDA DENSHI
	Deletes the patient review data of the transport monitor, and starts monitoring the patient of this equipment.
ſ	Start monitor the new patient
-	Deletes the patient review data of the transport monitor/this equipment, and starts monitoring the new patient.

 [Monitor Patient of This Equipment] Starts monitoring the patient of the DS-8400 by discharging the patient of the DS-8007.
 All data on the DS-8007 will be deleted. Only the patient information will be overwritten with the data of the DS-8400.

• [Monitor New Patient]: Starts monitoring a new patient by discharging the patients for both DS-8400 and DS-8007.

Both data on the DS-8400 and DS-8007 will be deleted/initialized.

> During data transfer, <Uploading> will be displayed.



CAUTION

· After the data transfer process, make sure that the patient information is correct.

REFERENCE

- To transfer the alarm settings and parameter settings when [Monitor Patient of the Transport Monitor] or [Monitor Patient of This Equipment] is selected, it is necessary to set [Transport] for "Data Transfer" under [Menu > Initial Settings > System > Other], and set [ON] for "Alarm Setup" and "Parameters" under [Menu > Initial Settings > System > Other > Data Selection for Transfer].

When EMR Link Function is Used

When the EMR link function is used on the central monitor, the patient admit/discharge process will be performed through the EMR.

When Starting the DS-8007 Transfer

- As the discharge process is performed through the EMR, monitor suspend confirmation window will be displayed when starting the DS-8007 transfer.
- ▶ To start the DS-8007 transfer, press the [Yes] key.



 $\mathbf{2}$ When Ending the DS-8007 Transfer

- ▶ When the transport monitor (DS-8007) with the patient unadmitted is connected to the host monitor (DS-8400) which is connected to the central monitor with the EMR link function, a confirmation window to upload to the central monitor will be displayed.
- When a patient is admitted to the transport monitor, uploading to the DS-8900 will automatically start.

Transport monitor is disconnected. Do you want to suspend monitoring?
OK Cancel
The data of an unadmitted patient is present. Do you want to upload the data to the central monitor?
2017/01/01 08:00:00 to 2017/01/01 10:00:00
Yes No

Data Transfer of Patient Data, Alarm Settings, Parameter Settings

By selecting [Transport] for "Data Transfer" (Menu > Initial Settings > System > Other), alarm settings, parameter settings, patient review data can be transferred.

Also, the patient data while in transfer can be transmitted to the DS-8900.

While the data is transmitted from the DS-8007, <Uploading > will be displayed on the central monitor and DS-8400.

CAUTION

 Do not disconnect the DS-8007 while <Uploading> is displayed. Otherwise, upload process cannot be completed.

Transferring the Alarm Setting, Parameter Setting

- **1** For the following case, a confirmation window to transfer the settings will be displayed.
- When the patient information matches between the DS-8400 and DS-8007
- When [Monitor Patient of the Transport Monitor] is selected on the confirmation window

Do you want to apply the following settings of the transport monitor to this monitor?
Alarm, Parameter
Yes No
Setting

2 [Yes]: Alarm settings, parameter settings will be transmitted according to the settings made on "Data Selection for Transfer" screen ([Menu > Initial Settings > System > Other]). [No]: Alarm settings, parameter settings will not be transmitted.

When [No] is selected, alarm settings, parameter settings will not be transmitted until the DS-8007 is disconnected and reconnected to the DS-8400.

[Change Setting]: "Data Selection for Transfer" screen will be displayed to change the settings.

Data Selection for Transfer
Alara Setup OFF
Parameters ON ON OFF
Set

- If [OFF] is set for alarm setting and parameter setting on the "Data Selection for Transfer" screen, a confirmation window to transfer the settings will not be displayed, and the settings will not be transferred.
- When the alarm settings and parameter settings are transferred, make sure the settings are properly transferred to the DS-8400.

NOTE

 When the alarm settings of the transport monitor is transferred: Even if the alarm threshold of the transport monitor exceeds the alarm threshold limit of the DS-8400, the exceeded alarm threshold will be applied to the DS-8400. Make sure to check the alarm setting on the DS-8400 as the alarm threshold limit status will be changed to "Limit Deactivating Mode".

2 Transferring the Patient Review Data

- 1 For the following case, the review data of the DS-8007 will be saved to the DS-8400.
- When the patient information matches between the DS-8400 and DS-8007 Trend, recall, full disclosure waveform, numeric data of the latest transferring data
- When [Monitor Patient of the Transport Monitor] is selected on the confirmation window Trend, recall, full disclosure waveform, numeric data for the selected patient from admittance
- When the bed is transferred using the EMR link function Trend, recall, full disclosure waveform, numeric data while the patient has been transferred
- When the EMR link function is used, and the DS-8007 is temporarily disconnected Trend, recall, full disclosure waveform, numeric data of the latest transferring data
- 2 The saved waveform data will be displayed on the full disclosure waveform screen as shown below. On the explanation area on the review screen, the starting and ending time of transfer will be displayed.

Menu	Vaveform Review Zoom Wave ST 12-Lead Full Disc. Data Raview Trend ► Transport 01/01 00:00 - 01/01 00:00 -	
[24h	■	▶ 24h Latest
9.51.00	-lash-da-da-da-da-da-da-da-da-da-da-da-da-	
9:51:30	had the first of t	
9.52.00	-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-	Tine Search Size/Scale
9:52:30	-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-	Setup
9.53.00	-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-dr-	Alarn Display
9:53:30	-dr. dr. dr. dr. dr. dr. dr. dr. dr. dr.	Print

- (Orange): Indicates the starting time (>) and ending time (<) of transfer. For starting point and ending point, one (1) second blank display is shown.
- (Purple): Indicates the starting time (>) and ending time (<) of data while the patient was monitored on another bed.



 To transfer the patient review data, it is necessary to use the CFast Card to save the full disclosure waveform data. If the CFast Card is not used, the patient review data will not be transferred.

Precautions when Starting the Data Transfer

- The data transfer process will not start for approximately 10 seconds after the DS-8007 is connected to the DS-8400. Do not disconnect the DS-8007 or turn OFF the power of the DS-8400 during this time.
- When the DS-8007 in standby mode is connected to the DS-8400, the data transfer process will not start for approximately one (1) minute after the DS-8007 is started. Do not disconnect the DS-8007 or turn OFF the power of the DS-8400 during this time.
- When a patient is discharged while DS-8400 and DS-8007 are connected, it will take approximately 10 seconds to synchronize the discharge information between the DS-8400 and DS-8007. Do not disconnect the DS-8007 or turn OFF the power of the DS-8400 for 10 seconds after the discharge process.

• When a patient is admitted/discharged on the DS-8400 while the DS-8007 is transferred, patient selection window will be displayed when the DS-8007 is connected.

When [Monitor Patient of the Transport Monitor] is selected, the patient on the DS-8400 will be discharged which will delete all data.

When [Monitor Patient of This Equipment] is selected, the patient on the DS-8007 will be discharged which will delete all data.

When [Monitor New Patient] is selected, both patients on the DS-8007 and DS-8400 will be discharged which will delete all data.

Cancellation of Uploading

The uploading will be canceled under the following condition.Once canceled, the uploading will not resume. To upload, disconnect and reconnect the DS-8007 to the DS-8400, and manually upload the data by selecting the corresponding data on the central monitor. For procedure, refer to the operation manual of the central monitor. <Cancellation of Uploading to the Central Monitor>

- The DS-8007 was disconnected from DS-8400 during uploading.
- The patient of the DS-8400 was discharged during uploading.
- On the central monitor, discharge process or bed transfer was performed, or bed registration was canceled for the uploading bed.
- The power of the DS-8400 or central monitor was turned OFF.

- DS-LAN cable was disconnected.
- <Cancellation of Uploading to the Host Monitor (DS-8400)>
- The DS-8007 was disconnected from DS-8400 during uploading.
- The patient of the DS-8400 was discharged during uploading.
- The power of the DS-8400 was turned OFF.
- Uploading to the central monitor was canceled.

About the Uploading Process

- The uploading will be performed in the order of central monitor and DS-8400. During the process, <Uploading> will be displayed on the central monitor and DS-8400.
- The uploading process is performed one at a time within one DS-LAN network. If the DS-8007 is connected while uploading for other bed is in process, <Upload Standby> will be displayed until the uploading can be started. When the uploading completes for the other bed, uploading for the connected bed starts, and the displayed message will change to <Uploading>. If the DS-8007 is disconnected, or a patient is discharged on the DS-8400, or bed transfer/exchange is performed on the central monitor during the standby condition, uploading will not be performed.
- If the MPDR function of the DS-8007 is used, and uploading data to the central monitor is selected, the data will be uploaded in the following order.
 - Selected data on the MPDR data list (Upload 1) (Upload to the central monitor)
 - Selected data on the MPDR data list (Upload 2) (Upload to the central monitor)
 - Transferring data (Upload to the central monitor)
 - Transferring data (Upload to the bedside monitor)
 - Other bed data (Upload to the bedside monitor)

REFERENCE

• For details of the MPDR function of the DS-8007, refer to the DS-8007 System Operation Manual.

To Stop Monitoring

This section explains about the procedure to stop monitoring.

 $m{7}$ Turn OFF the standby switch on the display unit.

• A standby confirmation message will appear.

2 Press [OK] to enter into standby mode.

- A 10-seconds progress bar will be displayed.
- Press the [Cancel] key to stop entering into standby mode. Only the [Cancel] key will be effective while the progress bar is displayed.

When 10 seconds has elapsed without pressing the [Cancel] key, the display will turn OFF and monitoring will stop.

• The operation of the Super Unit and the Input Box will also stop.

Monitor will enter into standby mode.
 Monitor will enter into standby mode.
To continue monitoring, press the [Cancel] key.
Cancel

Using the standby switch to stop monitoring will allow to easily resume monitoring by turning ON the standby switch again.

Æ CAUTION

- · If the remaining battery capacity becomes extremely low during battery operation, monitoring will automatically stop.
- If not using the equipment for a long period, disconnect the power cable and lithium-ion battery.

NOTE

When the power is turned OFF, graphic/tabular trend data (Vigilance, ventilator), recall, ST measurement, OCRG data will be erased after 5 minutes.

Clock Setup

This section explains about the time/date setup procedure.

CAUTION ∕¶∖

- If the time/date is not correctly set, or changed during monitoring, malfunction may occur with the NIBP measurement, periodic printing, graphic/tabular trend data, and age calculation from the birth date.
- The time/date can not be set while connected to a wired network system. The time/date will synchronize with the central monitor.
- If the time/date is changed, the time/date for all the saved patient data (trend, list, recall, etc.) will also change.

The printed time/date before changing and the displayed time/date after changing will differ. Also, the data transmitted to the central monitor before the time/date is changed will be displayed on the central monitor with the previous time/date.

Press the [Menu], [Time/Date] ("Basic Setup") keys.

Or, press the time/date on the information display area at the upper part of the screen.

The Time/Date setup window will be displayed.





 $m{Z}$ Press on the area to perform the setup.

- A blue frame will be displayed on the selected area.
- > When the screen is first displayed, the blue frame will be positioned on "hour".

 $\mathbf{3}$ Use the numeric keys to enter the numerics.

• The blue frame will automatically move to the next item.

4 Enter the current date/time and press the [Set] key.

- > The entered date/time will be set.(The number of seconds will be set to "00".)
- ▶ Press [Cancel] to cancel the time/date setup.

Installing the Recording Paper

- About the Recording Paper
 - Use only "OP050-01TDR" for the recording paper. If the surface treatment and thickness of the recording paper are different, it may result in poor print quality.
- Storing the Recording Paper Since the recording paper is thermal type, inappropriate storage may change the quality of the printed content, and make it illegible.

When storing the recording paper, follow the precautions below.

- Store in a place where light is shut off and avoid direct sunlight.
- Do not leave the paper in a high temperature (50 °C/122 °F and above).
- Do not store the paper in a polyvinyl chloride bag.
- · Do not superpose the papers until the diazo copy is completely dried.
- Do not expose the paper to alcohol, hydrochloric acid, or ester ketone.
- Avoid using adhesive agents other than water based glue.
- Installing the Recording Paper
 - When installing the recording paper, pay attention not to touch the thermal head or sensor. The temperature of those parts rises immediately after printing and may cause burn injury. Also, it may cause failure to the thermal head and sensor.
 - · Do not operate the equipment with wet hand. Doing so may short the thermal head.

Install the recording paper with the following procedure.

Press the Open/Close Lever.

• The paper holder will open.



2 Set the Paper.

The outside surface of the paper is heat-sensitive. Make sure to place the outside surface of the paper facing up.





Chapter 5 Admit/Discharge

To Display the "Admit/Discharge" Screen	5-1
Admit	
Entering the Patient Information	5-1
Entering Patient Information from the Magnetic Card	5-4
Entering Patient Information from the Patient Data Server (W	Vhen DS-
LANIII is used)	5-4
Discharge	5-6
Discharging Procedure	5-6
User Mode	5-8
To Select the User Mode	5-8
Suspend Monitoring	5-9
To Suspend Monitoring	
To Resume Monitoring	

Chapter 5 Admit/Discharge

This menu allows setup of admitting, discharging, suspend monitoring of a patient, and selection of the user mode (display configuration) according to the monitoring purpose.

• If monitoring of new patient is started without performing a discharge procedure of the previous patient, new data will be added to the previous data which will result in inaccuracy.

To Display the "Admit/Discharge" Screen

Press the [Menu], "Admit/Discharge" icon.

• The "Admit/Discharge" screen will be displayed.

Menu > Admit/	Discharge		<u>)</u>
= Require	:d Iten		
BED- Hode 001 Select	INITIAL	■Class. Adult	<u>ן</u>
		Sex 📃	
Name	FUKUDA DENSHI	Team 📕	
		Blood ABD	
		Blood Rh	
		■Pacemaker Not Used	
Birth Date	Age O	Impedance ON Meas.	
Heisht (cm)	0.0 Weight 0.0 BSA 0.00 (kg) (m²)]	Nonitor Suspend
Admit Date/Time]	
			Discharge

Admit

This section explains the admit procedure.

This menu allows entering of patient's name, ID, age, and selection of patient classification (adult, child, neonate) and pacemaker usage (used, not used) which affects the monitoring accuracy.

Entering the Patient Information



Up to 16 characters of alphabets, numbers, or symbols can be used. The entered name will be displayed on the home display. (PTo Enter Characters' P3-25)

$\mathbf{2}_{\mathsf{Enter}}$ the patient ID.

Up to 20 characters of alphabets, numbers, or symbols can be used. After entering the ID, press the [Set] key. If the [Set] key is not pressed, the entered ID will not be finalized. (To Enter Characters' P3-25)

NOTE

On a wired network (DS-LANII/III), up to 10 digits of ID can be transmitted.
 (Plantenance Manual "DS-LAN Setup" P2-2)

 ${f 3}$ Select the patient classification from [Adult] / [Child] / [Neonate].

> The selected patient classification and icon will be displayed on the home display.



- The patient classification affects the accuracy of NIBP measurement, HR measurement, and RR measurement. It also affects the delay time to generate the measurement data alarm.
- The alarm delay time is the function to prevent frequent generation of the measurement data alarm by holding the alarm generation for the duration of each delay time.

The alarm delay functions for HR/PR, BP, RR, SpO₂, TEMP, EtCO₂/InspCO₂, Tachy, Brady, Ext Tachy, Ext Brady.

			Adult	Child	Neonate		
NIBP Measurement Range MAP		SYS	30 mmHg to 280 mmHg	30 mmHg to 180 mmHg	30 mmHg to 130 mmHg		
		MAP	15 mmHg to 235 mmHg	15 mmHg to 160 mmHg	15 mmHg to 100 mmHg		
		DIA	10 mmHg to 200 mmHg	10 mmHg to 150 mmHg	10 mmHg to 90 mmHg		
HR			0 bpm, 12 bj	0 bpm, 30 bpm to 300 bpm			
	Monitor ESIS		0.5 Hz to 40 Hz		1.6 Hz to 40 Hz		
Filter Mode			1.6 Hz	1.6 Hz to 15 Hz			
	Diagnosis		3 electrodes: 0.05 Hz to 100 Hz				
Diagnosis			4, 5, 10 electrodes: 0.05 Hz to 150 Hz				
Impedance Respiration			1.:	2.5 Hz			
Alarm delay time			5 sec. 0 sec.				

WARNING

- The patient classification selection influences the precision of the QRS detection and NIBP measurement. Make sure the correct selection is made.
- The NIBP air hose corresponded to the set patient classification must be used to perform NIBP measurement. (However, if the patient classification is child, NIBP air hose for adult

can be used.)

 When "Link with Patient Class." is set to [ON], and patient classification is changed, the main mode will change to the selected mode on the "Link Settings". (Phaintenance Manual "To Program the User Mode" P5-31)

4 Select the sex from [Male]/[Female].

5 Select the color of the nurse team.

Enter the patient's age.

There are two ways to enter the patient's age. One is to enter the birth date which will automatically calculate the age, and the other is to directly enter the age using the numeric keypad. If [Neonate] is selected for patient classification, age will be displayed in days.

When Pacemaker is Used

WARNING

• The pacemaker usage setting influences the precision of the QRS detection and arrhythmia analysis.Make sure the correct selection is made.

If [Used] is selected for "Pacemaker", the monitor will detect the pacing pulse (pacemaker pulse) to perform the following process.

- The artificial pacemaker pulse will be displayed.
- When pacing waveform does not appear (pacing failure), erroneously detecting the pacemaker pulse as QRS will be prevented.
- •The arrhythmia analysis will detect pacing beat as P (Pacemaker Beat) or F (Fusion Beat) to prevent erroneous judgment of VPC.
- Menu 📏 Admit/Discharg 5 1 F = Required Iten Hode INITIAL BED Adult =ID Name FUKUDA DENSH Not Used Birth Date Age Height 0.0 Weight 0.0 (cm) (kg) 0.0 BSA 0.00 Date/Time

Press the key for "Pacemaker", and select from [Used]/[Not Used].

> When [Used] is selected, <Pacemaker> will be displayed at the upper part of the home display.

BED-001 CH0801	FUKUDA DENSHI P138240	- Adult Pacenaker	INIT Other IAL Ded Vent. Drift-F Enable ON
		l'l'_	
allR	r~~~		

Entering Patient Information from the Magnetic Card

By using the magnetic card reader, patient information can be entered from the magnetic card. The admittance process will speed up compared to manually entering each information.

NOTE
 To automatically enter the patient information from the magnetic card or barcode, it is necessary to perform the setup in advance. (Maintenance Manual "Using the Magnetic Card Reader" P4-22)

1 Read the data from the magnetic card or barcode.

> The acquired data will be displayed.

2 Press the [Change only patient info.]/[Cancel] key.

- ▶ [Change only patient info.] : Replaces the current patient information with the newly acquired information.
- ▶ [Cancel] : Cancels the acquired data.

NOTE

- Make sure the patient is discharged before replacing the patient information.
- The item which the information was not acquired from the magnetic card or barcode will be left blank. For the blank item, manually enter the information.

Entering Patient Information from the Patient Data Server (When DS-LANIII is used)

When the central monitor is connected to the patient data server, and the DS-LANIII network is used, patient information can be automatically entered by searching on the patient data server.

NOTE

 When a DS-LAN II network is used, patient information cannot be entered from the patient data server.

 When "Link with Patient Class." is set to [ON], and patient classification is changed by acquiring patient information from the patient data server, the Main Mode will change to the selected mode on the "Link Settings". (PMaintenance Manual "To Program the User Mode" P5-31)

When Using the Patient Data Server and Magnetic Card Reader (or Barcode Reader)

NOTE

- Select [ON] for "Auto Reference to Central Monitor when Reading Patient ID" under [Initial Settings>Magnetic Card Reader] in advance.
 - (@Maintenance Manual "Magnetic Card Reader Setup" P4-22)

Read the data from the magnetic card or barcode.

 ${f 2}$ The acquired patient information from the patient data server will be displayed in the "New Information" area.



If there is no applicable patient information, current patient information will be displayed in the "New Information" area.

Select whether or not to enter the searched patient information.

Select from [Change only patient info.] / [Discharge and admit as new patient.] / [Cancel].

[Change only patient info.] will replace the current patient information to the newly acquired information. [Discharge and admit as new patient.] will initialize the current patient data/monitoring condition and admit the searched patient as new patient.

[Cancel] will invalidate the acquired data.

\square	NOTE	\square							
	 1 ··			 		 			

- The item not acquired from the patient data server will be left blank.
- · For the blank item, manually enter the information.

When Magnetic Card Reader (or Barcode Reader) is not Used

1 Press the [Menu], "Admit/Discharge" icon, [ID]. "ID" window will be displayed.



 $\mathbf{2}$ Enter the patient ID.

3 Press the [Search ID] key and start searching on the patient data server.

- 1 Use the touch keys to enter the ID.
- 2 Based on the entered patient ID, patient information will be searched on the patient data server through the **DS-LANIII** network.

	ID	X)
ID [101234567		
12345	67890-	_]
QWERT	TYUIOP]
		h
ZXCV	V B N M , . /	J
ABC OWERTY	■ Delete Input C:	ancel

The searched patient information will be displayed under "New Information" .

E	xplanation Area			
Ne# Info. Search ID	ID: P123456 Name: FURUDA DERSHI Class.: Adult Sex: D00: Age: 55 Height(ca): 70.0 Beight(ka): 78.0 BSA(af): 1.08 Pacesaker: Hot Used		Change only patient info. Durront ness. data/cettings will remain. Sicharge, and admit as may pair of the source of the sour	
Current Info.	ID: Name:	ו⊾ן	Cancel	

4 Select whether or not to enter the searched patient information.

Select from [Change only patient info.] / [Discharge and admit as new patient.] / [Cancel]. [Change only patient info.] will replace the current patient information to the newly acquired information. [Discharge and admit as new patient.] will initialize the current patient data/monitoring condition and admit the searched patient as new patient. [Cancel] will invalidate the acquired data.

The item not acquired from the patient data server will be left blank. For the blank item, manually enter the information.

(NOTE

- The item not acquired from the patient data server will be left blank.
- For the blank item, manually enter the information.

Discharge

This section explains about the discharge process.

This procedure will erase the patient name, ID, age, and past measurement data such as tabular / graphic trend, and recall.

By pressing the [Rapid Discharge] key preprogrammed as user key, a discharge process can be performed.

Discharging Procedure

- If monitoring of new patient is started without discharging the previous patient, the measurement data of the previous and new patient will become mixed up on the recall and trend data.
- When the discharge process is performed, patient data such as recall and trend will be initialized. The parameter and alarm settings will be reset according to the settings made under [Menu>Initial Settings>User I/F>Power ON/Discharge).
 When the discharge process is performed on the central monitor, alarm will be reset according to the setting on "Admit Setup" of the central monitor.
 (Paintenance Manual "Power ON/Discharge" P5-18)
- If the power is turned OFF or if the system enters into standby mode soon after the discharge procedure, the patient may not be discharged on the central monitor.
 If it is necessary to turn OFF the power or enter into standby mode after the discharge procedure, select [Standby] for "Discharge Mode" under [Initial Settings>User I/F>Power ON/Discharge].

NOTE

- Depending on the setting made for "At Discharge" under ([Initial Settings>User I/F>Power ON/Discharge], some items may not be initialized.
 (Power ON/Discharge" P5-18)
- If the discharge procedure is performed during stopwatch operation, the counting will stop and will be reset to "00:00:00".
- The monitoring condition after discharge can be set on "Discharge Mode" under [Initial Settings>User I/F>Power ON/Discharge].

Press the [Discharge] key on the "Admit/Discharge" screen.

- > The discharge confirmation window will be displayed. (shown on right)
- ➤ To cancel the discharge process, press the [No] key or close the discharge confirmation window.



Press the [Yes] key.

• The patient data, patient information will be initialized.

Data	Description
Patient Data	Trend, Tabular Trend, Recall, ST, OCRG, CO, Hemodynamics, Lung Function, P-V/F-V control data will be erased. The settings for recall, tabular trend, graphic trend, vigilance list will remain.
Patient Information	Erases the data of patient name, ID, sex, age. The patient classification will not be initialized.
Measurement Condition	The learned arrhythmia waveform data will be deleted. The BP zero-balance condition will be initialized. The NIBP target inflation value will be initialized to the default value of each patient classification.

 If [Standby] is selected for "Discharge Mode" under [Initial Settings>User I/F>Power ON/Discharge], standby progress window will be displayed.
 Pressing the [Cancel] key will cancel the process to enter into standby mode.

After 10 seconds, discharge procedure will be performed and the system will enter into standby mode.



To continue monitoring, press the [Cancel] key.

Cancel

REFERENCE

For details, refer to [Initial Settings>User I/F>Power ON/Discharge].
 (@Maintenance Manual "Power ON/Discharge" P5-18)

User Mode

This section explains about the user mode selection.

From the preprogrammed user mode, an appropriate user mode can be selected according to the monitoring purpose.

 The selected user mode will be stored even after the power is turned OFF or discharge process is performed.

Before monitoring, make sure the current user mode is suitable for the patient's condition. (@Maintenance Manual "User Mode Registration" P5-30)

REFERENCE

For the user mode, up to 9 main modes of display configuration and alarm settings can be registered according to the patient's age and monitoring purpose.
 Also, for temporarily changing the display configuration (ex. when checking the 12-lead ECG), 6 sub modes of display configuration can be registered.
 For the extended display, 3 modes for each extended display (1, 2) can be registered.
 (Patient Content Conten Content Content Conten Content Content Content Content Con

To Select the User Mode

1 Press the [Menu], "Admit/Discharge" icon, "Mode Select" key. Or, press the mode key on the information display area at the upper part of the screen.

The "Mode Select" window will be displayed. (shown on right)

T INITIAL A HENO. FUL Á LOCA Å Ŵ A int Return to Wain Node Sub Mode 12LEAD SUB SUB6 MACHIN STAFF EXT1-3 Extended CIRCU ECG CA EXT2-Extended

Select the main mode or sub mode appropriate for the patient.

When the extended display is used, select the mode for the extended display.

WARNING

After changing the mode, make sure that the monitoring setting is appropriate.
 When the mode is changed, patient classification, alarm settings, etc. will be changed.

REFERENCE

- The selected user mode will be stored even after the power is turned OFF. If a new patient is admitted without changing the user mode, the monitoring will start with the previous user mode.
- The mode setting after the discharge operation can be set under [Initial Settings>User I/ F>Power ON/Discharge].
- · To change from the sub mode to the main mode, press [Return to Main Mode].
- Refer to "Setup Item/Default Value" for the default setting of each mode.
 (@ Maintenance Manual "User Mode Registration" P5-30)

Suspend Monitoring

This section explains about the monitoring suspend/resume function.

Monitoring suspend function can be used when a patient temporarily leaves the bed. If the monitoring is ceased by turning the power OFF, recall and ST data will be erased.

By using the monitoring suspend function, measurement, alarm, printing will be suspended but data and settings will remain, which allows to resume monitoring smoothly.

By using the monitor suspend label function, different labels in different colors according to the patient's destination can be displayed during the monitoring suspended condition.

To remind the user to resume monitoring, alarm will generate after the preprogrammed duration (15 min./30 min./1 hr/1.5 hr/2 hr) for "Monitor Suspend Timer".

REFERENCE

The monitor suspend label can be set on the Initial Settings menu.
 (@ Maintenance Manual "Monitor Suspend Setup" P5-16)

To Suspend Monitoring

When "Monitor Suspend Label" is not set:

1 Press the [Menu], "Admit/Discharge" icon, [Monitor Suspend] keys.

- The monitor suspend confirmation window will be displayed.
- ► If [Cancel] is pressed, monitoring will not be suspended and the confirmation window will close.

2 Press the [OK] key.

- ► The screen will automatically return to the home display with "Monitoring is suspended" message and [Resume] key.
- ➤ On the home display, numeric data and waveform display will be suspended.

Monitoring will be suspended.
Monitoring is suspended. Resume

REFERENCE

- When the monitoring is suspended, telemetry transmission will cease. Note that the square wave will be displayed on the central monitor indicating the too far condition of the telemetry.
- The stopwatch counting will continue even when the monitoring is suspended.
- The setting can be changed even when the monitoring is suspended.

Uhen Both "Monitor Suspend Label" and "Monitor Suspend Timer" are set

Press the [Menu], "Admit/Discharge" icon, [Monitor Suspend] keys.

- ▶ The "Monitor Suspend" screen will be displayed.
- Z Select the label to be displayed during the monitoring suspended condition.
 - The monitoring suspend duration selection will be displayed after selecting the monitoring suspend label.
- 3 Select the monitoring suspend duration from [15Min.]/ [30Min.]/[1Hr.]/[1.5Hr.]/[2Hr.]/[Continuous]. [Continuous] will start to suspend monitoring without setting the duration.
 - Confirmation window to suspend monitoring will be displayed. (shown on right)
 - Pressing the [Suspend] key will suspend the monitoring.

Verify that the monitoring is suspended on the home display.

The selected label with the set color will be displayed on the home display.

- On the home display, the time will start counting for the set duration.
- ➤ When the set duration completes, alarm sound will generate (5 sec. interval), and alarm indicator will light.

REFERENCE

• To extend the monitoring suspended duration, press [Extend] to display the timer selection.

When "Monitor Suspend Label" is set, but "Monitor Suspend Timer" is not set:

Press the [Menu], "Admit/Discharge" icon, [Monitor Suspend] keys.

▶ The "Suspend" screen will be displayed.

- Select the label to be displayed during the monitoring suspended condition.
 - A confirmation message will be displayed.
 - > Pressing the [Suspend] key will suspend the monitoring.
 - The selected monitor suspend label with the set color will be displayed on the home display.

	Menu > Admit/Discharge > Suspend Monitor
2.	Explanation Area
	SISYSACC INCR IZAN IN REMAN BATHTAG OUT SIREENY RESTRONG

	Menu > Admit/Discharge > Suspend Monitor
	(T)
2	SUSPENDED UNDER EXAM IN REHAB BATHING
<u> </u>	OUT SURGERY RESTROOM
3	When the selected time has classed, it will be notified by an alarm.
0	15 min. 30 min. 1 hr. 1.5 hr. 2 hr. Continuous





To Resume Monitoring

• Resuming monitoring will also resume the suspended alarm.

1 Press the [Resume] key.

• The "Monitoring is suspended" message will disappear and monitoring will resume.
Chapter 6 Alarm Function

Alarm	6-1
To Set the Arrhythmia Alarm	6-1
SpO2 Second Alarm Setup	
ST Alarm Setup	6-4
List of Alarm Settings	6-5
Detail Setup	6-6
Alarm Limit Setup	6-7
To Set the System Alarm (ON or Suspend)	6-7
To Silence or Suspend the System Alarm Sound	6-8
Alarm Limit Setup for Each Parameter	6-10
About the Alarm Threshold Limit	6-11
Alarm Assist Screen	6-13

Chapter 6 Alarm Function

Alarm

To Set the Arrhythmia Alarm

The arrhythmia alarm can be turned ON or OFF, and arrhythmia detection level can be set.

WARNING

- Set the appropriate upper and lower alarm limit for each parameter according to the monitoring condition.
- When the system alarm is suspended, all the alarms will be suspended even if the parameter alarm is set to ON. Also, the alarms will not be stored as recall events.
- If the upper/lower alarm limit of the parameter is set to OFF, or if arrhythmia alarm is set to OFF, alarm will not function even if the system alarm is enabled. Pay attention when setting them OFF.

1 Press the [Menu], [Arrhy.] ("Alarm") key.

> The arrhythmia alarm setup screen will be displayed.



 $\mathbf{2}$ Set ON/OFF of each arrhythmia.

- [ON]: Arrhythmia alarm will generate.
- ▶ [OFF]: Alarm will not generate.

NOTE

• The <ARRHY OFF> message will be displayed when the Slow_VT, Tachy, Brady, Ext Tachy, Ext Brady, and HR alarm is OFF.

REFERENCE

- The arrhythmia detection level for tachycardia (Tachy), bradycardia (Brady), extreme tachycardia (Ext Tachy), extreme bradycardia (Brady) alarms link with the upper and lower alarm limit for HR/PR.
 - The tachycardia (Tachy) alarm generates when the value exceeds the HR/PR upper alarm limit. When the upper alarm limit is OFF, alarm will not generate.
 - · For the Ext Tachy alarm, the alarm threshold level cannot be set below that of Tachy

alarm.

- The bradycardia (Brady) alarm generates when the value exceeds the HR/PR lower alarm limit. When the lower alarm limit is OFF, alarm will not generate.
- For the Ext Brady alarm, the alarm threshold level cannot be set above that of Brady alarm.

 $\mathbf{3}$ Select the level to detect each arrhythmia.

Item	Description	ltem	Description
Asystole	3 sec. to 10 sec.	R on T	200 ms to 600 ms
Run	2 beats to 8 beats	SVT	2 beats to 10 beats
Pause	1.5 sec. to 5 sec.	Irregular RR	10, 15, 20%
Frequent	1 bpm to 50 bpm	S Frequent	1 bpm to 50 bpm
Ext Tachy	22 bpm to 300 bpm	Pacer Not Capture	80 ms to 480 ms
Ext Brady	20 bpm to 295 bpm	Pacer Not Pacing	20 bpm to 200 bpm

4 Press the [Detail Setup] key, and set the HR Lower Limit for VT and RUN.

- 1 "HR Lower Limit for VT"
 - Select the lower limit of HR value from 120 bpm / 140 bpm to generate VT.
 - ▶ If the HR is below the selected value, Slow_VT will generate.
- 2 "HR Lower Limit for RUN"
 - ▶ If the HR is same or above the selected value, RUN will generate.
- **3** "HR Lower Limit for SVT"
 - ▶ If the HR is same or above the set value, SVT alarm will generate.



SpO₂ Second Alarm Setup

The SpO₂ second alarm function is available when HS-8312N or DS-8007N, HG-820 is connected.

When the SpO_2 value is unstable around the lower alarm limit, the frequently generated alarm may be bothersome. The second alarm function controls these frequent alarms.

This function generates the alarm only when the integral value (the accumulation of difference between the alarm limit and SpO_2 value at every second) reaches the preprogrammed second alarm threshold value.

NOTE

 The SpO₂ second alarm function utilizes SatSecondsTM technology of Covidien. SatSecondsTM is a trademark of Covidien. The integral value of the second alarm is calculated as follows.



- 1 Integral Value
- 2 Alarm Generation
- 3 SpO₂ Value
- 4 Alarm Limit

On this graph, the second alarm threshold value is set as 100.

The SpO_2 value begins to fall below the alarm limit at approximately 10 seconds. At the same time, the integral value begins to increase. (Alarm limit) – (SpO_2 value) is accumulated each second.

At approximately 25 seconds, the integral value reaches 100 and the alarm is generated.

The SpO_2 value begins to fall below the alarm limit at approximately 36 seconds. At the same time, the integral value begins to decrease. [(Alarm limit) – (SpO_2 value)]x 2 is subtracted each second.

Also, there is a safety net when setting the second alarm function. This safety net is for the case when the SpO_2 value frequently falls below the alarm limit but does not last long enough to reach the second alarm threshold.

If the SpO_2 value falls below the limit 3 times or more during the last 60 seconds, an alarm will be generated even if the second alarm threshold is not reached.

- Whether to use the second alarm function and its threshold selection should be based on the patient's clinical indication portent and medical evaluation.
- If the SpO₂ alarm and second alarm setup is set to [OFF], the second alarm integral value will be set to 0.

Press the [Menu], [SpO₂] ("Parameter") keys to display the "SpO₂" setup screen.



- **2** Set the "Second Alarm".
 - [10]/ [25]/ [50]/ [100]: A circular second alarm indicator will be displayed inside the numeric data box.
 As the integral value increases, the indicator will begin to fill, and when it is completely filled, an alarm will be generated.



• [OFF]: Second alarm indicator will not be displayed.

ST Alarm Setup

Set the ST upper limit and lower limit for the reference waveform.

The alarm limit should be set for each measurement unit (mm/mV). The upper and lower limit can be set in 1 mm / 0.1 mV increment.

1 Press the [Menu], [ST] ("Alarm") key.

> The ST alarm setup screen will be displayed.



2 Select [ON]/[OFF] for "ST All Alarm" .

• [OFF]: Alarms will not generate even if the alarm for each lead is set to ON.

 $\mathbf{3}$ Select the lead to set the alarm limit.

• The selected lead will be displayed large at the right.

4 Select [ON]/[OFF] of ST alarm for each lead.

5 Slide the \boxed{XXX} / \boxed{XXX} and set the upper, lower limit (±20 mm / ±2.0 mV).

- Alarm will be set to OFF if the value -20 mm / +2.0 mV or lower is selected.
- Alarm will be set to OFF if the value +20 mm / +2.0 mV or above is selected.

The alarm settings can be verified in list format. The alarm settings for each parameter can be changed on this list.

- **1** Press the [Menu], [List] ("Alarm") key.
 - The alarm settings list will be displayed.



2 Select from [All List] / [Meas. List].

- [All List]: The settings for all the parameters will be displayed.
- [Meas. List]: The settings for only the measured parameters will be displayed.

 $\mathbf{3}$ Change the alarm threshold.

- 1 Select a parameter.
 - The alarm setup window will be displayed.



2 Press xxx / xxx to set the threshold level.

Detail Setup

The alarm-related setup such as alarm suspend time and alarm silence time can be performed.

1 Press the [Menu], [Detail Setup] ("Alarm") keys.

• The alarm detail setup screen will be displayed.



2 Select [1 min.] / [2 min.] for "Suspend Time".

3 Select [1 min.] / [2 min.] for "Silence Time".

4 Set the "Alarm Sound Suspend" function.

- ▶ [ON]: The alarm sound suspend function will turn ON.
- ▶ [OFF]: The alarm sound suspend function will turn OFF.

5 Select the "Alarm Sound Suspend Time" from 1 min., 2min., 5 min., 10min., 30 min., 60 min., 90 min., 120 min., 240 min., 360 min.

6 Set the "Status Alarm Control".

REFERENCE

- The alarm silence time for the level L equipment status alarm ("Check electrodes", "NIBP Check patient type, air hose", etc.) can be set.
 (Prequipment Status Alarm Message" P11-8)
- [Link to Alarm Silence Time]: When the [Alarm Silence] key is pressed at occurrence of equipment status alarm, alarm will be silenced for fixed amount of time set for "Silence Time". If the alarm factor still remains at completion of silence time, the alarm sound will generate again.
- If the same alarm occurs during the alarm silence time, the alarm sound will not generate.
- If a new alarm occurs during the alarm silence time, the alarm sound for the new alarm will generate.
- [Link to each new occurrence]: When the [Alarm Silence] key is pressed at occurrence of equipment status alarm, the alarm will be silenced as long as the alarm factor remains regardless of the "Silence Time" setting. While the same equipment status alarm is generated, the alarm will remain silenced.
 If the alarm factor is resolved during the alarm silence time, the alarm will be canceled.
 If the same alarm generates again during the alarm silence time, the alarm sound will generate.

If the same alarm generates again during the alarm silence time, the alarm sound will generate.

Zselect [Graph]/[Numeric]/[OFF] for "Alarm Limit Display".

• The upper and lower alarm limit will be displayed on the home display.



NOTE

- The alarm limit for the parameter with the alarm turned OFF will not be displayed regardless of this setup.
- If the alarm limit display for BP is [Graph], systolic value will be displayed.
- Depending on the numeric data box type, alarm limit may not be displayed.
- If the SpO₂ lower alarm limit is set to 85%SpO₂ or below, the alarm limit value will be displayed regardless of the "Alarm Limit Display" setting.

Alarm Limit Setup

This section explains the procedure to enable/suspend the system alarm, and to set the upper/lower alarm limit for each parameter.

On this system, 9 modes can be preprogrammed according to the monitoring purpose. By preprogramming the alarm setting to each mode, the alarm setups at admittance of patient can be simplified by just selecting a mode. It is recommended to program the mode in rough classification such as patient's age, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient.

To Set the System Alarm (ON or Suspend)

The system alarm can be enabled or suspended.

When the system alarm is enabled, the alarm suspended condition is canceled, and alarm limit and alarm ON/OFF setting for each parameter are effective. The system alarm cannot be disabled.

WARNING

- When the system alarm is suspended, all the alarm will be suspended even if the parameter alarm is set to ON. Also, the alarms will not be stored as recall events.
- If the upper/lower alarm limit of the parameter is set to OFF, or if arrhythmia alarm is set to OFF, alarm will not function even if the system alarm is enabled. Pay attention when setting them OFF.

Press the [Menu], [Basic] or [Circ.] or [Resp./Gas] ("Alarm") keys.

• The alarm setup screen will be displayed.



 $\mathbf{2}$ Select whether to enable or suspend the alarm.

To Suspend the Alarm

1 Press the [Alarm Suspend] key.

- The key will change to blue.
- ▶ The alarm will suspend temporarily.
- Alarm Suspend (xxx sec.)> will be displayed.

BED-001 FUKUDA DENSHI - ▲ 4011 ALT

(xxx sec.) indicates the remaining time. The system alarm will be enabled when the suspended time completes.

To Enable the System Alarm

1 Press the [Alarm Suspend] key while in alarm suspended condition.

- The key will change to gray.
- > The alarm limits and ON/OFF settings for each parameter will become effective.
- > The alarm suspended condition will be canceled.

To Silence or Suspend the System Alarm Sound

The alarm sound can be suspended for fixed amount of time. There are two ways to suspend the alarm sound, which are "Alarm Silence" and "Alarm Sound Suspend".

The "Alarm Silence" function suspends the alarm sound for fixed amount of time (1 min. / 2 min.).

The "Alarm Sound Suspend" function suspends the alarm generation in advance such as during surgery when the alarm generation is expected. Alarm monitoring will continue even while the alarm sound is suspended. The alarm sound suspend duration can be selected from 1 min., 2min., 5 min., 10min., 30 min., 60 min., 90 min., 120 min., 240 min., 360 min.

To silence the alarm, press the [Alarm Silence] key (fixed key).

- The alarm sound will be silenced for fixed amount of time.
- > If the alarm factor still remains at completion of silence time, the alarm sound will generate again.
- > The [Alarm Silence] can also be operated on user keys or remote control.

 $\mathbf{2}$ To suspend the alarm sound, press the Alarm Silence key (fixed key) for more than 3 seconds.

> The alarm sound will be suspended for fixed amount of time.

• During the alarm sound suspended duration, the alarm sound will not generate.



- If the [Alarm Silence] key is pressed while the alarm sound is generated, it will bring the system to "Alarm Silence" condition and not the "Alarm Sound Suspend" condition.
- During the "Alarm Sound Suspend" duration, other bed alarm sound will not generate.

Precautions about Silencing the Alarm

The alarm silence function is effective for each parameter. Once the alarm cause is resolved, the alarm silence condition for that parameter will be canceled.

When [Fukuda Tone] is set for "Alarm System" under [Menu>Setup>Initial Settings], and if another alarm with the lower priority occurs during the alarm silence duration, alarm sound will not generate. The recall and alarm printing will function.

When [Fukuda Tone] is set for the "Alarm System" and equipment status alarm is silenced, the alarm sound for the lower priority numeric and arrhythmia alarm will generate.

When [Melodic Tone] or [Standard Tone] is set for the "Alarm System" and if another alarm with lower priority occurs, the alarm sound will generate.

If the [Alarm Silence] key is pressed for the alarm of another parameter which occurred during the alarm silence condition, the alarm silence duration for the first alarm will not be extended.

The alarm silence condition for all parameters will be canceled for the following case.

- When the power is turned ON.
- When the system alarm status (enable/suspend) is changed.
- When the monitoring is suspended on the "Admit/Discharge" screen.
- When the user mode is changed.
- When the patient is discharged.
- When [Resume All Al. Sound] key on the alarm setup screen is pressed.

The alarm silence condition for each parameter will be canceled for the following case.

- When the alarm cause is resolved for that parameter.
- When the alarm silence time for the parameter is completed.
- When automatic alarm is set for the parameter.
- When the alarm is turned OFF for the parameter.

If [Link to each new occurrence] is set for "Status Alarm Control" (Menu>Alarm>Detail Setup), the alarm sound will not generate until the alarm condition changes even the set alarm silence duration completes.

Precautions about Suspending the Alarm Sound

During the alarm sound suspended duration, recall and alarm printing will function.

The alarm sound suspended condition will cease in the event of any of the following.

- Discharge
- When OFF is set for "Alarm Sound Suspend".
- When the ventilator alarm is generated.
- When resumed from monitor suspend condition.
- When the [Alarm Silence] key is pressed.

Alarm Limit Setup for Each Parameter

The alarm for each parameter can be turned ON or OFF, and upper and lower alarm limit can be set.

- Set the appropriate upper and lower alarm limit for each parameter according to the monitoring condition.
- When the system alarm is suspended, all the alarms will be suspended even if the parameter alarm is set to ON. Also, the alarms will not be stored as recall events.
- If the upper/lower alarm limit of the parameter is set to OFF, or if arrhythmia alarm is set to OFF, alarm will not function even if the system alarm is enabled. Pay attention when setting them OFF.
- When the numeric data acquired from FLOW-i is displayed, the following alarms cannot be set. Also, these alarms will not generate.
 InspCO₂/EtCO₂, InspO₂/ExpO₂, InspN₂O/ExpN₂O, InspAgent/ExpAgent, MAC, ExpMV, PEAK, PEEP

1 Press the [Menu], and then the key for "Alarm".

▶ The alarm setup menu will be displayed.



 $\mathbf{2}$ Select the parameter group from the tab.

REFERENCE

 The standard parameters will be displayed on the Menu screen. The parameters to be displayed here are selectable.
 (Control Manual "Alarm Related Setup" P5.4)

(@Maintenance Manual "Alarm Related Setup" P5-4)

 $\mathbf{3}$ Select ON/ OFF for the individual alarm.

- ▶ [ON]: Alarm of the corresponding parameter will generate.
- [OFF]: Alarm of the corresponding parameter will not generate.

4 Set the upper/ lower limit.

- 1 Slide the (xxx) (xxx) keys on the right side of the bar.
 - $\angle xxx$: Adjusts the upper limit.
 - XXX : Adjusts the lower limit.
 - ▶ By releasing the finger from the key, fine-tune keys will appear for a fixed period of time.



About the Alarm Threshold Limit

By setting the alarm threshold limit ("Initial Settings") in advance, the alarm threshold can be limited within the preprogrammed range. When the alarm threshold limit function is enabled, threshold limit will be displayed beside the alarm bar.

(@Maintenance Manual "Alarm Related Setup" P5-4)



Above is an example of alarm threshold limit setting where HR is set to [Enable], and upper and lower limits are set to 180 bpm and 40 bpm respectively.

- NOTE
- The alarm threshold limit can be set for each parameter. When enabling this function, make sure the upper and lower limits are set appropriately.
- When the alarm threshold limit function is enabled, pressing the [Auto] key for alarm settings will set the alarm threshold within the limit range.

Limit Deactivating Mode

Even when the alarm threshold limit function is enabled, the alarm threshold outside the limit can be temporarily set. This is called the "Limit Deactivating Mode."

By pressing the up arrow key \blacktriangle for 2 seconds at the upper threshold limit, the limit can be deactivated. The arrow keys will turn to blue indicating that the upper threshold limit can be exceeded.

In the same way, by pressing the down arrow key \checkmark for 2 seconds at the lower threshold limit, the limit can be deactivated. The arrow keys will turn to blue indicating that the lower threshold limit can be exceeded.

When the alarm threshold is set within the limit range, the limit deactivating mode will end.



Above is an example of HR upper threshold limit being deactivated. The upper limit keys are turned to blue indicating that the upper limit 180 bpm can be exceeded.

NOTE

- If the alarm threshold set on the central monitor exceeds the threshold limit set on the bedside monitor, the alarm threshold set on the central monitor will be applied. In such case, the threshold limit is deactivated.
- If the alarm threshold of the transport monitor exceeds the alarm threshold limit of the DS-8400, the exceeded alarm threshold will be applied to the DS-8400. Make sure to check the alarm setting on the DS-8400 as the alarm threshold limit status will be changed to "Limit Deactivating Mode".
- If the alarm threshold of "Setup at Discharge" exceeds the alarm threshold limit, the
 exceeded alarm threshold will be applied. Make sure to check the alarm setting at
 admittance as the alarm threshold limit status will be changed to "Limit Deactivating Mode".
- If the monitor mode is changed, and the alarm threshold of the current monitor mode exceeds the threshold limit, this alarm threshold will be applied. In such case, the threshold limit is deactivated.

Alarm Assist Screen

On the alarm assist screen, maximum of 24 hours of trend data for the corresponding parameter will be displayed. Alarm limit can be set by using the past trend data as reference.

1 To display the alarm assist screen, press [Menu], select a parameter, and press **a** on the corresponding parameter setup screen.

Or, press the numeric data box on the home display, and press **mathematical** on the corresponding parameter setup screen.

> The alarm assist screen will be displayed.



2 Select the time range on the time bar.

- > Dragging the slider to the right will display newer data, and dragging it to the left will display older data.
- ▶ Pressing [24h] will switch the display by 24 hours.



4 Set the upper and lower alarm limit.

- **1** Press /xxx / xxx on the right of the bar.
 - Alarm zone will be displayed on the trend.



- The displayed alarm zone will slide by sliding the \boxed{XXX} or \boxed{XXX} .
- The displayed alarm zone will also slide by pressing the ▲ / ▼
- 2 Set the alarm limit by using the alarm trend as reference.

Chapter 7 Monitoring

To Display the Parameter Setup Screen	7-1
ECG	7-1
Before Attaching the Electrodes	
Electrode Placement	
Type of Electrodes and Lead Cable	
Connection to the Patient Monitor	
ECG Parameter Setup	
Respiration	
Respiration Monitoring (Impedance Method)	
RESP Parameter Setup	
BP	
BP Monitoring	
Zero Balance of All Pressure Lines (User Key)	
Zero Balance of All Pressure Lines ([BP Zero] Key)	
Zero Balance for Each Pressure Line	
BP Parameter Setup	
Pulse Oximetry	
SpO2 Monitoring	
SpCO, SpMet, SpHb, SpOC Measurement (Masimo)	
Precautions about the Masimo Sensors and Cables	
SpO2 Parameter Setup (Nellcor)	
RR_SpO2 Parameter Setup (Nellcor)	
SpO2 Parameter Setup (Masimo)	
Non-Invasive Blood Pressure	
Lineup of Cuffs	
NIBP Monitoring	
Inflation Mode Setup	
NIBP Auto Mode Setup	
Oscillation Graph Display	
Dyna Alert Function Status NIBP Parameter Setup	
Temperature	
TEMP Monitoring	
TEMP Parameter Setup	
Cardiac Output and Blood Temperature	
Connecting the Super Unit	
Cardiac Output Measurement Algorithm	
Blood Temperature Alarm Setup	
CO2 Concentration (Mainstream Method)	7_71
Patient Application and Display	
CO2 Parameter Setup	
CO2 Concentration (Sidestream Method)	
Patient Application and Display	
CO2 Parameter Setup	
Multigas Unit/SPIRO	
Connecting to the Respiration Circuit	
CO2 Measurement Unit Setup	
GAS Display during Undetected Breath	7-88
Multigas Unit Data Setup (Multigas Concentration/Spirometry)	
BIS Data (HBX-800 with BISx)	
Preparation for Monitoring	
BIS Setup	
BIS Data (A-2000/A-3000)	
	1-30

Ventilator	7-98
AWP/AWF/AWV Scale Setup	7-99
P-V/F-V Loop Display	7-100
FLOW-i Data	7-101
CO2 Measurement Unit Setup	7-101
FLOW-i Setup	
SvO2/CCO Data	7-107
STAT Mode / Index Display	7-107
INVOS Data	7-108
Channel Number Setup for INVOS Data	7-108
Stopwatch	7-109
Label Setup	
Start/Stop	
Multiparameter Connector Setup for BP, TEMP, CO Measure	ement.7-110
Multiparameter Connector Setup	7-111

Chapter 7 Monitoring

To Display the Parameter Setup Screen

This section explains how to display the "Parameter Setup" of monitoring parameters.

Press the [Menu], and then select the parameter to perform the setup.

Or, press the numeric data box on the home display, and press () on the corresponding parameter setup screen.

> The "Parameter Setup" screen will be displayed.



• When the numeric data box on the home display is pressed, a floating window for the basic setup such as size/scale will be displayed.

To display the "Parameter Setup" screen for detailed setup, press $\textcircled{\sc star}$.



ECG

This section explains the procedure for ECG measurement preparation and monitoring condition setup.

Before Attaching the Electrodes

- Make sure to use electrodes of the same type.
 If different types of electrodes are used at the same time, the difference between the polarization potential from each electrode may interfere monitoring.
- ECG measurement part is Type CF applied part, but it is not intended to directly apply on patient's heart.



Electrode Placement

Depending on the lead cable type, 3-electrode/4-electrode/5-electrode/10-electrode placements are available. Using the 4-electrode, 5-electrode or 10-electrode application allows simultaneous monitoring of 2 ECG waveforms, and high accuracy of arrhythmia analysis can be attained. (1 to 12 waveforms can be displayed depending on the number of electrodes.)

Also, the displayed lead type can be changed.

□ For 3-electrode lead cable (1 waveform monitoring)

Lead Type: [I]/[II]/[III]

Symbol	Color	Electrode Site
RA	White	On the right infraclavicular fossa
LA	Black	On the left infraclavicular fossa
LL	Red	On the left midclavicular line, near the supracrestal line.



□ For 4-electrode lead cable (Maximum 6 waveforms monitoring) Lead Type: [I]/[II]/[aVR]/[aVL]/[aVF]

Symbol	Color	Electrode Site
RA	White	On the right infraclavicular fossa
LA	Black	On the left infraclavicular fossa
LL	Red	On the left midclavicular line, near the supracrestal line.
RL	Green	On the right midclavicular line at the same height as LL.



□ For 5-electrode lead cable (Maximum 7 waveforms monitoring) Lead Type: [I]/[II]/[aVR]/[aVL]/[aVF]/[V]

Symbol	Color	Electrode Site
RA	White	On the right infraclavicular fossa
LA	Black	On the left infraclavicular fossa
LL	Red	On the left midclavicular line, near the supracrestal line.
RL	Green	On the right midclavicular line at the same height as LL.
V	Red/Brown	Chest electrodes (V1 to V6)



General For 10-electrode lead cable (Maximum 12 waveforms monitoring)

Lead Type: [I]/[II]/[aVR]/[aVL]/[aVF]/[V1]/[V2]/[V3]/[V4]/[V5]/[V6]

	I	
Symbol	Color	Electrode Site
RA	White	On the right infraclavicular fossa
LA	Black	On the left infraclavicular fossa
LL	Red	On the left midclavicular line, near the supracrestal line.
RL	Green	On the right midclavicular line at the same height as LL.
V	Red/Brown	The fourth intercostal space at the right sternal border.
V2	Yellow/Brown	The fourth intercostal space at the left sternal border.
V3	Green/Brown	On the midway between V2 and V4.
V4	Blue/Brown	The fifth intercostal space on the left midclavicular line.
V5	Orange/ Brown	On the left anterior axillary line at the same horizontal level as V4.
V6	Violet/Brown	On the left midaxillary line at the same horizontal level as V4.



```
NOTE
```

Electrode Placement for 12-Lead ECG Analysis

When acquiring 12-lead ECG signals, Fukuda Denshi recommends placing the limb electrodes anywhere along the arms and legs as shown below.
However if it is difficult, use the Mason-Likar 12-lead system.
To reduce the waveform differences from the standard 12-lead, Fukuda Denshi recommends that the torso placement of the RA and LA electrodes be near as possible to each arm, in the infraclavicular fossae, within the area unaffected by myoelectricity.)

Symbol	Color	Electrode Site
RA	White	On the right arm
LA	Black	On the left arm
LL	Red	On the left leg.
RL	Green	On the right leg.
V	Red/Brown	The fourth intercostal space at the right sternal border.
V2	Yellow/Brown	The fourth intercostal space at the left sternal border.
V3	Green/Brown	On the midway between V2 and V4.
V4	Blue/Brown	The fifth intercostal space on the left midclavicular line.
V5	Orange/ Brown	On the left anterior axillary line at the same horizontal level as V4.
V6	Violet/Brown	On the left midaxillary line at the same horizontal level as V4.



Type of Electrodes and Lead Cable

There are various types of disposable electrodes for ECG measurement depending on the connection method with the lead cable and materials which the electrodes are made of. Make sure to use the appropriate electrodes which will make full use of the characteristics.

Do not reuse/resterilize the disposable electrodes.

For details of usable lead cables, refer to *F*"ECG, Impedance Respiration Measurement (Manufactured by Fukuda Denshi)" P13-1

Connection to the Patient Monitor

- The indication for continuous use of the electrode is about one day.
- · Replace the electrode if the skin contact gets loosen due to perspiration, etc.
- When an electrode is attached to the same location for a long period, some patients may develop skin irritation. Check the patient's skin condition periodically and change the electrode site as required.
- When using the electrosurgery-proof type ECG relay cable, the impedance respiration cannot be measured, and its numeric data and waveform will not be displayed. When measuring in an environment where electrosurgery is not performed, make sure to use the standard ECG relay cable.

(NOTE

- Use only the specified relay cables, lead cables, and electrodes.
- The conductive parts of electrodes and associated connectors for applied parts, including the neutral electrode, should not contact other conductive parts including earth.

1 Clip on the lead cable end to the electrode convex part.

 $\mathbf{2}$ Turn right and left to verify that it is securely connected.









• ECG waveform and HR data will be displayed on the monitor.



5 Adjust the waveform size and position, and change the monitoring lead as necessary. ([] "ECG Parameter Setup" P7-6)

ECG Parameter Setup

Press the [Menu], [ECG] keys to display the "ECG" setup screen.



Adjustment of Waveform Size and Baseline Position

Adjust the waveform size and baseline position.

- The threshold level for arrhythmia detection changes with ECG waveform size. Set a proper waveform size for monitoring. When the ECG waveform size is x1/4, x1/2, or x1, the arrhythmia detection level is 250 μ V. When the ECG waveform size is x2 or x4, the arrhythmia detection level is 150 μ V.
- Automatic size/position of the ECG is effective only at the time the [Auto] key is pressed. This does not continuously adjust the size and position.
- The waveform size and position cannot be set if the waveform is not displayed. Refer to "To Configure the Display" P10-4, and change the display configuration as necessary.
- The threshold level for HR detection changes with ECG waveform size. Set a proper waveform size for monitoring.

REFERENCE

- By setting the [ECG Size (All Leads)] key as user key, ECG size for all leads can be changed at once.
 - (@"User Key Setup" P10-10)

Press the key for "ECG1" to "ECG12", and display the "Size" selection window.

When the display layout is "12-Lead", the waveform size can be set differently for limb leads and chest leads.

 $\mathbf{2}$ Select the waveform size for displaying/printing.

▶ [Auto]: ECG amplitude will be automatically adjusted to 10 mm. The automatic adjustment is effective only when the [Auto] key is pressed.

Waveform Size	x1/4	x1/2	x1	x2	x4
Voltage (10 mm)	4mV	2mV	1mV	500µV	250µV



If the waveform is difficult to see due to ECG amplitude, press ▲/ ▼ and set the baseline position to 0 mV. The baseline position for the waveform display and printing will be adjusted.

Lead Selection

Set the monitoring lead.

- The leads for arrhythmia detection, central monitor display, printing are fixed as ECG1 and ECG2. Set the most appropriate leads with high QRS for ECG1 and ECG2, especially for arrhythmia detection.
- The alarms for HR, Tachy, Brady will not be generated when the electrode for ECG1 or ECG2 lead is detached, and for 30 seconds after the electrode is reattached.

Press the key for "ECG1" to "ECG12", and display the "Lead" selection window.

When the display layout is "12-Lead", select the lead for ECG1 and ECG2 on the lead selection window.

Z Select the ECG monitoring lead.



HR Alarm Setup

Set the HR alarm.

("Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Set the upper limit in the range of 22 bpm to 300 bpm. The upper limit alarm will become OFF if the value exceeds 300 bpm.
- Set the lower limit in the range of 20 bpm to 295 bpm. If a value below 20 bpm is set, the lower alarm will turn OFF.
- Ext Tachy alarm threshold cannot be set below HR upper alarm limit, and Ext Brady alarm threshold cannot be set above HR lower alarm limit.

REFERENCE

- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Baintenance Manual "Alarm Related Setup" P5-4)
- When [Auto] is set, the upper and lower limit will be automatically set to +40 bpm and -40 bpm to the current value respectively. The lower limit will be clipped to the setting made for "HR/PR Lower Limit during Alarm Auto Setting" (Menu>Initial Settings>Alarm).
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit. However if "HR/PR Lower Limit during Alarm Auto Setting" is also set, HR lower alarm limit will be clipped to the larger value.
- Ext Tachy will be set to HR upper limit+10 bpm, Ext Brady will be set to HR lower limit-10 bpm. When the set value exceeds 300 bpm for the upper limit and 20 bpm for the lower limit, the setting will be clipped to 300 bpm and 20 bpm respectively
- When [Auto] is set for Ext Tachy, Ext Brady, the same setting, HR upper limit+10 bpm, HR lower limit-10 bpm, will be set respectively.

Arrhythmia Alarm Setup

Set the arrhythmia alarm.

(To Set the Arrhythmia Alarm" P6-1)

Detail Setup



Set the filter mode.

- Select from [Monitor]/[ESIS]/[Diag.] according to the monitoring purpose.
- Each mode has different frequency characteristic.
- The selected filter mode will be printed along with other data.
- ▶ On the waveform area, "M" (Monitor), "E" (ESIS), or "D" (Diagnosis) will be displayed.

CAUTION Ŕ

- · The ESIS mode cannot completely reduce the electrical noise, and may erroneously detect the pacemaker spike.
- The ESIS mode should be selected only when a high frequency noise largely affects the ٠ HR measurement.
- In ESIS Mode, artifacts such as electrosurgical noise or EMG can be largely reduced, but QRS amplitude attenuation, waveform distortion, or ST segment change may occur compared with other filter modes.

Monitor Mode	This is the standard mode for ECG monitoring.
(Frequency Characteristic: Adult/Child 0.5 Hz to 40 Hz, Neonate	The highest frequency is set to 40 Hz to reduce the
1.6 Hz to 40 Hz)	artifact caused by EMG, etc.
ESIS Mode (Frequency Characteristic: Adult/Child/Neonate 1.6 Hz to 15 Hz)	By selecting this mode during electrosurgery, noise can be largely reduced.
Diagnosis Mode	Select this mode if ST measurement or high
(Frequency Characteristic: 3-electrode Adult/Child/Neonate 0.05	frequency ECG monitoring is performed.
Hz to 100 Hz	As the lowest frequency is set to 0.05 Hz, ST level
4, 5,10-electrode Adult/Child/Neonate 0.05 Hz to 150 Hz)	can be accurately measured.

(NOTE

 When the filter mode is changed, a notch will appear on the ECG waveform due to the change in frequency characteristic as shown below.



2 Select [Used]/[Not Used] for "Pacemaker".

- [Used]: Pacemaker pulse will be detected and pace pulse mask function will be performed for set duration.
- [Not Used]: Pacemaker pulse will not be detected.

3 Set the "Synchronized Mark/Tone".

- [OFF]: Synchronized mark will not be displayed.
- [Auto]: The priority will be according to the setting of "Synchronized Mark/Tone Priority" [Menu>Initial Settings>Meas.>Other]. (Maintenance Manual "Other Setup" P5-11)
 [ECG]: The priority of synchronized mark/tone will be set in the order of ECG>SpO₂-1>SpO₂-2>BP.
 [SpO₂]: The priority of synchronized mark/tone will be set in the order of SpO₂-1>SpO₂-2>ECG>BP.
- ▶ [ECG]: HR synchronized mark will be displayed. The synchronized tone will turn ON.
- ► [SpO₂-1]/[SpO₂-2]: SpO₂ synchronized mark will be displayed. The synchronized tone will turn ON.
- [BP]: BP synchronized mark will be displayed. The synchronized tone will turn ON.



4 Set the "Pacemaker Pulse".

- · Precautions about Pacemaker Pulse Detection
 - There are some cases when the pacemaker pulse cannot be detected depending on the pacemaker type, pulse voltage, pulse width, electrode lead type (unipolar, bipolar), or electrode placement which causes the pacemaker pulse amplitude to decrease, and disables the pacemaker pulse detection.
 - If signals similar to a pacemaker pulse are present, such as electric blanket noise or excessive AC frequency noise, these may be erroneously detected and displayed as a pacemaker pulse.
 - When a spontaneous QRS and pacemaker pulse overlap (ex. fusion beat, etc.), QRS detection cannot be performed properly. In this case, the heart rate is degraded.
 - If a pacemaker pulse is continuously detected due to AC frequency interference, QRS detection will be suspended and the heart rate will be reduced. Arrhythmia will not be detected either.

REFERENCE
Pacemaker Pulse Detection Algorithm
ECG Signal Input ECG signal will be input.
Pacemaker Pulse Detection and Suspension of QRS Detection
Detects the high frequency and large amplitude signal as pacemaker pulse.
When pacemaker pulse is detected, QRS detection will be suspended for fixed amount of time to avoid erroneous detection of pacemaker pulse as QRS.
Canceling of Arrhythmia Detection Arrhythmia detection of the waveform following the pacemaker pulse will be canceled.

Select from [ON] or [OFF].

[ON]: The pacemaker artificial pulse will be displayed on to the ECG waveform with a different color. "Pacemaker Pulse" will be automatically set to [ON] when [Used] is selected for "Pacemaker" on the "Admit/ Discharge" screen.

[OFF]: The pacemaker artificial pulse will not be displayed.

5 Set the "Pace Pulse Mask Time".

• WARNING

 If the QRS pace mask function is set to [OFF]/[10ms]/[20ms]/[40ms], the pace pulse may be erroneously be detected as a QRS complex and HR alarm or asystole alarm may not generate due to incorrect HR (counting pace pulse as QRS complex). Select [OFF]/ [10ms]/[20ms]/[40ms] only if you are sure that pacing failure will not occur, or when the patient can be constantly monitored.

REFERENCE

 For the patients using pacemakers, there are cases when the pacing waveform may not occur in spite of the pacing stimulus. This condition is called "pacing failure". To avoid detecting pacemaker pulses as a QRS complex, this monitor has a function to suspend QRS detection for a fixed amount of time starting from the detection of the pacing stimulus. This function is called "pace pulse mask".

But if the pacemaker does not detect the patient's spontaneous heartbeat (sensing failure), and the pacing stimulus is applied at the same timing as QRS, this pace mask function may erroneously mask the QRS and cause the heart rate measurement to decrease.

To avoid this, QRS pace pulse mask function can be set to [OFF]/[10ms]/[20ms] for correct measurement of the heart rate. (Default: Auto)



- 1 Pacemaker Pulse
- 2 Pacing waveform caused by pacemaker pulse
- 3 No waveform in spite of pacing stimulus
- 4 Pacemaker pulse and spontaneous heartbeat occurring at the same time
- 1 Press the key for "Pace Pulse Mask Time".
- 2 Select the mask time depending on the pace spike amplitude or presence of fusion beat.
 - [Auto]: Pace pulse mask time will be automatically set according to the pace pulse amplitude.
 - [OFF]: Pace pulse mask time will be set to 0 ms.



• [ON]: Only the amplitude with frequency component under 1 Hz will be attenuated to prevent the ECG baseline drift.

The patient signal display will delay about 0.5 seconds.

П

On the home display, "Drift-F ON" will be displayed in the information area, and "DF" will be displayed in the waveform area.

• [OFF]: ECG drift filter will not be set.

7Set the "HR Average".

▶ [Instant]: HR measured from RR interval of each heartbeat will be displayed.

FUKUDA DENSH

• [Average]: HR measured from 6 seconds of heartbeat for adult and child, and 3 seconds of heartbeat for neonate will be displayed.

8 Set the "HR Delay".

- ▶ [OFF]: HR will be calculated based on the "HR Average" setting.
- ▶ [ON]: HR will be calculated based on the arrhythmia analysis. 5 seconds delay will occur compared to when [OFF] is selected. It may improve the HR detection when T wave or noise is interfering. When two ECG waveforms (ECG1 and ECG2) are measured, HR will be calculated by merging ECG1 and



ECG2.

If artifact is present on one of the waveforms, HR will be calculated using only the stable ECG waveform. If artifact is present on both of the waveforms, HR value will be displayed as "---".

When ECG electrodes are detached, arrhythmia analysis cannot be performed, and <Lead OFF> message will be displayed. Alarm sound will also generate.

NOTE

- · When the patient classification is set to [Neonate], "HR Delay" will be set to [OFF].
- When the patient classification is set to [Adult] or [Child], and "HR Delay" is set to [ON], "Inst." or "Av." will not be displayed inside the HR or HR/PR numeric data box.

9Set the "AC Filter".

If the ECG waveform is interfered with AC noise, the AC filter cuts off the frequency component (50 Hz/60 Hz).

- ▶ [ON]: AC filter which attenuates the AC noise of 50 Hz to 60 Hz will be set. "AC" will be displayed in the waveform area.
- ▶ [OFF]: AC filter will not be set.

10 Set the "Auto Lead". The automatic lead switching will be performed for ECG 1 and ECG 2.

During Lead OFF

Lead Cable Type	Detached Electrode	Auto Lead Selected	
		ECG1	ECG2
4-electrode	RA		
	LA	II	II
5-electrode	RA/RA+V		III
	LA/LA+V	II	II
	V	II	aVR
10-electrode	RA/RA+V		III
	LA/LA+V	II	II
	V,V2 to V6	11	aVR

 [ON]: When lead-off condition occurs, the lead will automatically switch. Also, <Check Electrodes> will be displayed.

▶ [OFF]: The lead will not automatically switch even when lead-off condition occurs.

Set the "ST/VPC/Arrhy. Alarm Display".

- [ON]: If 2 or more boxes are used for ECG numeric data display, ST level, VPC, arrhythmia alarm factor will be displayed inside the ECG numeric data box.
- ▶ [OFF]: ST level, VPC, arrhythmia alarm factor will not be displayed inside the ECG numeric data box.

12 Set the "3lead Override".

\square	NOTE	

- When a relay cable for 5-lead or 10-lead is used with a 3-lead cable, it will be judged as lead-off condition and <LEAD OFF> message will be displayed.
 If a 3-lead cable is intentionally used, select [ON] for "3lead Override" to avoid displaying the <LEAD OFF> message.
- If [ON] is selected for "3-lead Override" even though 4-lead, 5-lead, or 10-lead relay
 cable is used with all the lead cables and electrodes connected, it will be acknowledged
 as only 3 electrodes are used and only one waveform will be displayed.
 Also, artifact may interfere to the waveform or lead-off information may become incorrect.

When using the "3lead Override" function, use only 3 electrodes of LA, RA and LL.

Select from [ON] or [OFF].

13 Select the lead for ECG analog output.

- [Disp. Lead]: The lead of the displayed waveform will be output.
- ▶ [Selected Lead]: The lead selected on "Output Lead Sel." window will be output.
- **14** Set the "ECG Waveform Display during Lead-OFF".

When the lead-OFF condition is detected, whether or not to display the waveform for detached lead can be selected.

- > [ON]: The input waveform will be displayed even during lead-off condition.
- [OFF]: Baseline will be displayed during lead-off condition.

15 Set the "Chest Lead-OFF".

Whether or not to detect the chest lead OFF condition can be selected. If set to [Enable], chest lead OFF condition will be notified by an alarm generation.

- [Enable]: Chest lead OFF condition will be notified by an alarm generation.
- ▶ [Disable]: Chest lead OFF condition will not be notified by an alarm generation.



 If chest lead is set for ECG1/ECG2, chest lead OFF condition will be notified by an alarm generation even if [Disable] is set for "Chest Lead-OFF".

16 Set the "Noise Detection".

When a noise generating from electrosurgery, body motion, etc. is detected, whether or not to retain the HR data before the noise detection and to switch the synchronizing source to SpO₂/BP can be selected.

- [ON]: HR data before the noise detection will be retained, and synchronizing source will switch to SpO₂, BP.
- ▶ [OFF]: HR data before the noise detection will not be retained, and synchronizing source will not switch to SpO₂, BP.

NOTE

 Even if the synchronizing source is switched to SpO₂, the ECG tone will remain and not change.

ON/OFF of Parameter Display

Select ON/OFF for parameter display.

• When the waveform and numeric data display is set to OFF, the alarm generation and tabular trend input will also cease.



1 Press the [Disp. ON] key.

> The "Display ON/OFF" confirmation window will be displayed.

2 Select from [Display ON] or [Display OFF].

- [Display ON]: Waveform and numeric data will be displayed.
- [Display OFF]: Waveform and numeric data will not be displayed. A message will be displayed inside the numeric data display area.



Display ON/OFF	\mathbf{X}		
ECG display can be turned ON or OFF.			
Display ON Display OFF If the electrodes are attached to the patient during "Display OFF" condition, the setup will automatically switch to "Display ON" after 10 seconds. Close			

When ECG electrodes are attached to the patient with the ECG display set to OFF, the ECG waveform and numeric data will be automatically displayed after 10 seconds.

Respiration

This section explains about the respiration measurement by the impedance, CO_2 , ventilator, or SpO_2 method and the measurement condition settings.

WARNING

• The SpO₂ respiration measurement is not intended for use as an APNEA monitor.

- When a defibrillator is used during respiration monitoring, a large offset voltage will be placed on the ECG electrodes, which may cause interruption of monitoring for a few seconds.
- When using the electrosurgery-proof type ECG relay cable, the impedance respiration cannot be measured, and its numeric data and waveform will not be displayed.When measuring in an environment where electrosurgery is not performed, make sure to use the standard ECG relay cable.

Respiration Monitoring (Impedance Method)

1 Check that the displayed ECG waveform is stable.

The respiration waveform is detected from ECG II or ECG I lead explained in the previous section. Therefore, a stable ECG waveform is necessary to acquire respiration waveform.

 $\mathbf 2$ Verify that the respiration waveform and respiration rate is displayed on the home display.



NOTE

 Adjust the detection lead, waveform size, baseline position, and sweep speed for optimum waveform display.

(@"To Configure the Display" P10-4)

• To change the lead, press the lead name on the waveform area, and display the lead selection window.

RESP Parameter Setup

Press the [Menu], [RESP] keys to display the "RESP" setup screen. The example when the HS-8000 is used is shown below.



Press the key for "Size" to adjust the waveform size and baseline position.

- Select from [x1/4] / [x1/2] / [x1] / [x2] / [x4].
- If the waveform is difficult to see due to impedance waveform amplitude, set the baseline position to 0Ω. The baseline position for printing will not change.

Use the $\blacktriangle/\checkmark$ keys to adjust the baseline position.



2 Set the RR alarm.

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The same RR alarm setting will be applied for impedance, CO₂, ventilator, gas unit, and SpO₂ measurement.
- For RR measured from CO₂ waveform, alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit within the following range for each patient classification. Adult: 10 Bpm to 150 Bpm Child/Neonate: 4 Bpm to 150 Bpm The upper limit alarm will turn OFF if the value above 150 Bpm is set.
- Set the lower limit within the following range for each patient classification. Adult: 5 Bpm to 145 Bpm Child/Neonate: 2 Bpm to 148 Bpm If a value below 5 Bpm / 2 Bpm is set, the lower alarm will turn OFF.
- For the impedance respiration, RR alarm will not generate when the electrode of detection lead is detached and for 30 seconds after the electrode is reattached.
- When Capnostat 5 is used, RR alarm will not generate unless 2 or more respiration is detected after completion of the airway adapter calibration.

REFERENCE

- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Plantenance Manual "Alarm Related Setup" P5-4)
- When [Auto] is set, the upper and lower limit will be automatically set to +20 Bpm and 20 Bpm to the current value respectively.
- The adjustable increment for upper and lower limit depends on the patient classification and "RR Alarm Increment" setting under "Initial Settings" > "User I/F".

	Alarm Increment (Initial Settings > User I/F)		
	Normal	Small	
Adult	5 Bpm increment	1 Bpm increment	
Child/Neonate	2 Bpm increment	1 Bpm increment	

• When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

3 Set the APNEA alarm.

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The same APNEA alarm setting will be applied for impedance, CO₂, SpO₂ module, and ventilator measurement.
- If the alarm is based on the apnea time measured from CO₂ waveform, apnea alarm will
 not generate unless 2 or more respiration is detected within 30 seconds after the power
 is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit in the range of 10 sec. to 60 sec. If a value above 60 sec. is set, the upper alarm will turn OFF.
- For the impedance respiration, apnea alarm will not generate when the electrode of detection lead is detached and for 30 seconds after the electrode is reattached.
- When Capnostat 5 is used, apnea alarm will not generate unless 2 or more respiration is

detected after completion of the airway adapter calibration.

REFERENCE

- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Plaintenance Manual "Alarm Related Setup" P5-4)
- When [Auto] is set, the apnea alarm setting registered for the currently selected mode will be applied.
- The upper limit can be set in 1 second increment. There is no lower limit.
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

4 Set the "RR Synchronized Mark".

▶ [ON]: The mark synchronized to impedance respiration or CO₂ waveform will be displayed.



• [OFF]: Synchronized mark will not be displayed.

5 Set the "RR/APNEA Alarm Source".

The parameter to display the RR synchronized mark and to generate the RR/APNEA alarm can be selected from impedance, CO₂/multigas unit, ventilator, and SpO₂.

WARNING

- The RR/APNEA alarm will not be generated unless the numeric data box corresponded to the selected RR/APNEA alarm source is displayed.Make sure to display the numeric data box for the parameter set as the RR/APNEA alarm source.
- The SpO₂ respiration measurement function is not intended for use as an APNEA monitor.

- If the "RR/APNEA Alarm Source" setting is other than [Impedance] (or, if [Auto] selects a setting other than [Impedance]), the RESP waveform will not be transmitted on a wired network.
- RR_SpO₂ measurement is optional function available on DS-8007N.
- [Impedance]: RR alarm will be generated based on the impedance respiration curve. The RR synchronized mark based on impedance respiration will be displayed.
- ▶ [CO₂/GAS]: When multigas unit/FLOW-i is used, RR alarm will be generated based on the RR measured by the multigas unit/FLOW-i.

If multigas unit is not used, RR alarm will be generated based on the RR measured by the HPD-800/HPD-810/HPD-820 (Capnostat 5) or HCP-800/HCP-810/HCP-820. The RR synchronized mark based on CO_2 waveform will be displayed.

- [Ventilator]: RR alarm will be generated based on the RR measured by the ventilator. The RR synchronized mark based on ventilator measurement will be displayed.
- [SpO₂]: RR alarm will be generated based on the RR measured by the SpO₂ module. The RR synchronized mark will not be displayed.
- [Auto]: The measurable parameter will be selected in the priority of CO₂/GAS or FLOW-i > ventilator > impedance > SpO₂, and generates the alarm if the corresponded numeric data box is displayed on the home display.

6 Set the "CVA Detect".

REFERENCE

- When the amplitude of the respiration waveform decreases due to causes such as respiratory pause, the ECG waveform may be superimposed on to the respiration waveform, making the RR equal to the HR. This condition is called CVA (Cardio-Vascular Artifact), and is detected using the CVA detection function.
- This function will be effective only when [Impedance] is set as the "RR/APNEA Alarm Source" or, when [Auto] selects impedance respiration.
- If the ECG waveform is superimposed on to the respiration waveform with HR (RR) of 30 Bpm or above for 20 seconds (10 seconds for neonates) or more and if the "CVA Detect" is set to [ON], the <CVA detected> message will be displayed, and an alarm sound will be generated.
- [ON]: When CVA is detected, alarm will generate and message will be displayed.
- [OFF]: CVA detection will not be performed.

Zset the "Impedance Measurement".

WARNING

 If a patient is using an adaptive (minute ventilation) pacemaker, "Impedance Measurement" should be set to OFF.
 The respiration measurement using the impedance method conducts high-frequency and weak current between the ECC electrodes attached to the patient, and measures the

weak current between the ECG electrodes attached to the patient, and measures the potential difference between the electrodes caused by thoracic movement using the synchronous rectification system. For the patient using the adaptive (minute ventilation) pacemaker, the pacemaker measurement signal and the high-frequency current of this equipment interferes with each other which causes incorrect respiration measurement.

- [ON]: Standard impedance respiration measurement will be performed.
- [OFF]: Impedance respiration measurement will not be performed and impedance respiration waveform and RR data will not be displayed. A high-frequency current which is a measurement signal will not be conducted.

Set the "Impedance Detection Lead".

Select the respiration detection lead from [I] or [II].

9 Select ON/OFF for parameter display.

(@"ECG Parameter Setup" P7-6)
BΡ

This section explains about the procedure of BP1 to BP 8 measurement preparation and measurement condition setup.

AUTION

- · Do not reuse / re-sterilize the disposable type transducers.
- If using a reusable blood pressure transducer, disinfect it according to the manufacturer's guidelines.
- The long-term use of the blood pressure transducer, tube and catheter may increase the risk of infection. Perform periodic replacement with new one. The guidelines of the CDC (Disease Control and Prevention) recommend replacing within 96 hours.
- If the ambient temperature of the blood pressure transducer has changed greatly, the zero balance may cause the drift. Perform the zero balance again.
- An operator must not get away from a patient during the BP measurement. However, when getting away from the patient is necessary, do not activate the Alarm Suspend and Silence functions in order not to miss any sudden changes in the patient's condition.
 (P"To Set the System Alarm (ON or Suspend)" P6-7)
 (PTo Silence or Suspend the System Alarm Sound" P6-8)
- Be sure to perform Daily Check. Use of faulty equipment might harm the patient or operator.
 (P "Daily Check" P4-1)
- If the Equipment Status Alarm occurs or if you feel the unusual operation of the equipment, perform the inspections to confirm the safety or contact our service representative.
 (P11-8)
- The BP value will not be displayed until zero balance is performed after the power is turned ON. Make sure to perform the zero balance.
 Once the zero balance is performed, the zero balance information will be maintained, and the BP value will be displayed.

BP Monitoring

The HS-8000 utilizes multiparameter amplifier input method which allows monitoring of 2 channels of BP through the 2ch BP conversion cable, CJO-P01B-DJ0.5. The BP relay cable can be directly connected to the multiparameter connector.

The measurement is also possible using the HM-800/HM-801 Multi Module

(@"Multiparameter Connector Setup for BP, TEMP, CO Measurement" P7-110)

Connect the 2ch BP interface cable to the HS-8000.

For Connection via 2ch BP Conversion Cable (CJO-P01B-DJ0.5):

7 Connect the interface cable to the multiparameter connector via 2ch BP conversion cable (CJO-P01B-DJ0.5).

- 1 Multiparameter Connector
- 2 2ch BP Conversion Cable CJO-P01B-DJ0.5
- 3 1ch BP Relay Cable CJO-P01B-S**



For Direct Connection:

1 Connect the BP relay cable directly to the multiparameter connector.

 1ch BP Relay Cable CJO-P01B-S**

 2ch BP Relay Cable CJO-P01B-D**



2 Assemble the BP measurement device.

REFERENCE

- The following procedure explains the case when a BP transducer (LS575 series) is used. If using other transducers, refer to the operation manual for the corresponding transducer.
- 1 Inspect transducer packaging for damage prior to opening.
- 2 Verify that each connector is securely connected.



3 Connect the BP relay cable to the transducer.



4 Align the bracket to patient's heart position (about 1/2 of the chest depth).



5 Inject 1000 units of heparin into the saline bag, mix thoroughly and puncture the infusion line through the same hole.



6 Set the saline bag to pressure bag, and hang from the infusion device. Fill saline to about 1/3 of the drip.

7 After loosening the zero-port plug, push the flash button to perform priming to remove air bubbles.



 $m{8}$ Verify that all air bubbles are removed, and tighten the zero-port plug. Turn on the zero-port plug side of

the open-air three-way valve.



 ${\bf 9}\,$ Inflate the pressure bag to 300 mmHg.



 $10\,$ Set the BP device and wait for about 5 minutes.

3Perform zero balance.

1 Loosen the zero-port plug on open-air three-way valve one-half turn.



 $2\,$ Press the BP numeric data box (parameter key) on the home display.

NG-001 FURUDA DENSHI - 10011	t Power • 10:08 2017/01/19
T 1 HR brokk. Sp02 96	Home
	Nenu
ST I COVIC SO	Aların Silence
·····································	Admit/ BP Zero Disch. 1987
AP1 control 12 12 12 12 12 12 12 12 12 12 12 12 12	NIBP Start/Stop
Menter 116/ 77 (92)	NIBP Alarm Setup ito Node (ALL)
0 BP2 Include A Control A	BP Cont. Alarn History
	BP List Recall abular Graphic
20/12(18).	Trend Trend
	Start/Stop Mathematics Key Lock
58 58 58 50 36.1 37.2 Marcine 53 37.5	KEY LOCK

> The BP floating window will be displayed.



3 Press the [Zero] key on the BP floating window or BP parameter setup screen.



- > Zero balance will start.
- When the BP zero balance is complete, the completed date/time will be displayed inside the [Zero] key.
- **4** Turn off the zero-port plug side of the open-air three-way valve.



5 Connect the catheter to the end of monitoring line.



> The measurement preparation is completed, and BP measurement will start.

4 Press the [Home] key on the user key.

5 Verify that the BP waveform and numeric data is displayed on the home display.

لمت المتحدة: 10:08 من 10:08 م
60 Rttoni 60 Meru
0.2 VPC 30 Alarm Silence
He 29/° 82 (* 98),
29/ OZ (98/, NIBP Start/Stop
NIBP Alarm
16/ 77 (92)
m ^{Mgl} 6/ 6 (6) BP List Recall
abular Graphic Trend
2/ 38 3102 (*) CCI LI/MR/M*I Strict
30. 0 5.3 37.5 Key Lock

- · The zero balance procedure is required for the following case.
 - · When starting the measurement.
 - · When the position of the heart has changed due to body movement.
 - · When the position of the transducer has changed.
 - When measuring for a long period of time and there is a possibility of measurement error due to change in ambient temperature, etc.
 - · When a connector is connected/disconnected, or a transducer is replaced.

Zero Balance of All Pressure Lines (User Key)

The zero balance for all the displayed BP can be performed using the user key. If any of the BP is in progress of measurement, perform the zero balance on each BP parameter setup screen.



1 Open the three-way valve of all the pressure transducers to air.

READY> will be displayed inside the user key.

Press the [BP Zero] key on the user key.

 ${f 3}$ Verify the BP waveform is positioned at zero, and "0" is displayed for the BP value.

- COMPLETE> will be displayed when the procedure is complete.
- FAILED> will be displayed when the process fails.
- > <DRIFT> will be displayed when the BP relay cable is not connected.

(NOTE

- If <FAILED> is displayed, the three-way valve may not be opened to air, artifact is
 present, or the transducer may be defective. Check the cause and try the zero balance
 procedure again.
- If <DRIFT> is displayed, verify that all the connections are secure.

4 Close the three-way valve when the zero balance is complete.

- When the transducer or tubing is replaced, make sure to perform the zero balance. Otherwise, accurate measurement will not be performed.
- <READY> will not be displayed unless the three-way valves of all pressure transducers are opened to air. If the status is not displayed, or if <MEASURE> is displayed, check if the three-way valve of pressure transducers are opened to air.

BP zero status displayed inside the user key

	No display	: Open transducer to air
DD 7		
BP Zero	MEASURE	: Open transducer to air
READY	READY	: Ready to perform zero balance.
	BP ZERO	: BP zero in progress
	FAILED	: Zero failed
	COMPLETE	: Zero complete
	DRIFT	: Zero drift

Zero Balance of All Pressure Lines ([BP Zero] Key)

By using the [BP Zero] key on the HS-8000 or Multi Module, zero balance can be performed for all the BP even if not displayed.

- When the BP zero balance properly completes, a beep sound will generate for 1 second and LED will light in blue.
- When the BP zero balance fails, a beep sound will generate for 3 seconds and LED will flash in blue.

NOTE Using the [BP Zero] key will allow to perform zero balance for all the BP even if not displayed on the home display. For the BP channel with the transducer in progress of measurement, zero balance will not

be performed.

Zero Balance for Each Pressure Line

1 Open the three-way valve of the pressure transducer to air.
2 Verify that "Zero ready" is displayed on the BP parameter setup screen for BP1 to BP8, and press the [Zero] key.
${f 3}$ Verify the BP waveform is positioned at zero, and "0" is displayed for the BP value.
A message, "Zero complete" will be displayed when the procedure is complete. When the BP zero balance is complete, the completed date/time will be displayed at the lower part of the [Zero] key.
▶ A message, "Zero failed" will be displayed when the process fails.
▶ A message, "Zero drift" will be displayed when the BP relay cable is not connected.
 NOTE If a message, "Zero failed" is displayed, the three-way valve may not be opened to air, artifact is present, or the transducer may be defective. Check the cause and try the zero balance procedure again.
 If a message, "Zero drift" is displayed, verify that all the connections are secure.
4 Close the three-way valve when the zero balance is complete.

BP Parameter Setup

Channel Settings of BP, TEMP

The default settings of BP and TEMP channels are as follows.

For HS-8000 For DS-8007



When the DS-8007 is used, the temperature measured on temperature jack is fixed as T1, T2.

NOTE

 The channel settings can be changed on the "Multi-amplifier" screen ([Initial Settings > System > Unit Module]).

(@Maintenance Manual "Unit Module Setup" P4-17)

Label Setup



- Description of Each Label: ART (Arterial Pressure) PAP (Pulmonary Artery Pressure) CVP (Central Venous Pressure) ICP (Intra-cranial Pressure) IAP (Intra-aortic Balloon Pumping Pressure) LVP (Left Ventricular Pressure) US1 to US5: User labels (3 characters) which can be set on the "Initial Settings". (Additional Settings (1) (Composition of the Settings)
 - US3 to US5 cannot be selected for the equipment connected to DS-LANII/III.

When the BP Label is ART

By selecting [ON] for "ART Catheter Check Message" [Menu>Parameter>BP1 (ART)>Detail Setup], an alarm will be generated when the catheter is disconnected.

(@"Set the "ART Catheter Check Message"." P7-32)

- The default setting of "ART Catheter Check Message" is [OFF].
- When "ART Catheter Check Message" is set to [ON], alarm will generate when the transducers are opened to air.

When the BP Label is IAP

PDP (Peak Diastolic Pressure) of IABP can be displayed in addition to systolic, diastolic, and mean pressure. Note that Systolic Pressure (SYS) = Peak Systolic Pressure (PSP).

 Note that Systolic Pressure (SYS)=Peak Systolic Pressure (PSP) when reviewing graphic trend, data base, or when setting the alarm. · When ECG is not measured, PDP cannot be calculated.

When the BP Label is CVP

The measurement unit can be selected from "mmHg", "kPa" or "cmH₂O".

The measurement unit can be selected on the "Initial Settings" menu. The selected unit will be displayed on the BP numeric data box.

(@Maintenance Manual "Measurement Unit" P5-10)



When the BP Label is ICP

CPP (Cerebral Perfusion Pressure) can be measured.

CPP = Mean Arterial Pressure - Mean Intracranial Pressure

If the CPP value is negative, the data will not be displayed. Also, alarm cannot be set for CPP.



PCWP Measurement

When PAP is set as BP label, the mean value can be displayed as PCWP (Pulmonary Capillary Wedge Pressure).



1 Press the key for "PCWP".

 PCWP measurement screen will be displayed.On the PCWP screen, the current BP waveform and RESP waveform will be displayed.

Select the waveform scale from [20]/[50] as necessary.

3 Press the [Freeze] key.

The displayed waveform will freeze and cursor will be displayed. The cursor point indicates the current mean pressure.



4 Use the 1/4 keys to set the PCWP value.

5 Press the [Input] key after setting the PCWP value.

• The PCWP value will be displayed inside the PAP (BP label) numeric data box with the measurement time.

It will be also displayed on the trend data.

Scale Setup

· When wireless network is used, BP waveform with a scale above the set scale will not be properly transmitted. The displayed BP scale should be within the set scale.

NOTE

- · Select the full scale for displaying and printing.
- ٠ The scale selection will differ depending on the label as shown below.
- Change the scale before the freeze operation, as the waveform will be deleted if the scale is changed after the freeze operation.

								Sc	ale						
	5	10	15	20	30	40	50	75	100	150	200	250	300	mmH	lg
BP Label	1	2	3	4	5	6	8	12	16	20	24	32	40	kPa	
							20	40 cmH ₂ O							
BP1 to BP8 User Label				0			0	0	0	0	0	0	0		
ART, IAP, LVP							0	0	0	0	0	0	0		
PAP				0		0	0	0	0	0	0	0	0		
CVP		0		0	0	0	0	0	0	0	0	0	0	0	0
ICP	0	0	0	0			0	0	0	0	0	0	0		

1 Press the key for "Scale Selection", and display the scale selection window.

 $\mathbf{2}$ Select the scale from the displayed selection.

Alarm Settings

Set the BP alarm.

(@ "Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Set the upper limit in the range of 2 mmHg to 300 mmHg / 0.2 kPa to 40.0 kPa. If a value ٠ above 300 mmHg / 40.0 kPa is set, the upper alarm will turn OFF.
- ٠ Set the lower limit in the range of 0 mmHg to 295 mmHg / 0 kPa to 39.5 kPa. If a value below 0 mmHg / 0 kPa is set, the lower alarm will turn OFF.
- · Alarm will not generate until 30 seconds has passed after the zero balance or after the transducer has been opened to air.

REFERENCE

- Select ON/OFF of BP alarm and set the upper and lower alarm limit for systolic (S), diastolic (D), and mean (M) BP.
- The alarm limit should be set for each unit (mmHg/kPa).
- The adjustable increment will be according to the "BP Alarm Increment" setting. (Normal/ Small).

(@Maintenance Manual "Display/Print Setup" P5-13)

- The adjustable increment for upper and lower limit changes from 50 mmHg / 7 kPa.
- When [Auto] is set for the BP label of BP1/ART, the upper and lower limit will be automatically set to +40 mmHg / +5 kPa and -20 mmHg / -3 kPa respectively to the current value.
- When [Auto] is set for the BP label other than BP1/ART, the upper and lower limit will be automatically set to +20%, -20% respectively to the current value.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (PAINT PS-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

	"BP Alarm Increment" Setup				
	If [Normal] is selected; If [Small] is selected;				
0 mmHg to 50 mmHg	2 mmHg increment	1 mmHg increment			
50 mmHg to 300 mmHg	5 mmHg increment				
0 kPa to 7 kPa	0.2 kPa increment				
7 kPa to 40.0 kPa	0.5 kPa increment				

Detail Setup (BP Parameter)

Press the [Menu], [BP] keys to display the "BP" setup screen.

The "BP" setup screen can be also displayed by pressing the detail key (\Box) on the BP floating window.



200 ABT

Alarn NIBP ART Catheter Check Wessam

٩٢



1 Set the "Synchronized Mark/Tone". (BP1/ART)

OFF

The parameter to display the HR synchronized mark can be selected from ECG, SpO₂, and BP (BP1 or ART). If BP1 and ART are measured simultaneously, ART will be prioritized.

- ▶ [Auto]: The synchronized mark will be displayed in the priority of "ECG > SpO₂-1 > SpO₂-2 > BP".
- [ECG]: HR synchronized mark will be displayed.
- ▶ [SpO₂-1]/[SpO₂-2]: SpO₂ synchronized mark will be displayed.

× X

80

- [BP]: BP synchronized mark will be displayed.
- ▶ [OFF]: Synchronized mark will not be displayed.

\square	NOTE) ———	
•	If the correspo	nding BP (BP1/ART) is not measured, PR (BP) will be displayed as "-	".

2 Set the "Display Type".

Â CAUTION

The undisplayed BP data will not generate a BP alarm or be displayed in the tabular ٠ trend. Select the appropriate display type according to the monitoring purpose.

NOTE

- The display type of numeric data can be selected from [S/M/D]/[S/D]/[M]. The BP alarm will not be generated unless the data is displayed.
- If the BP label is CVP, IAP, PAP, ICP, the display type is fixed.

- ▶ [S/D/M]: The systolic/diastolic/mean BP value will be displayed.
- ▶ [S/D]: The systolic/diastolic BP value will be displayed.
- [M]: The mean BP value will be displayed.

 ${f 3}$ Select the "Wave Filter" from [6Hz]/[8Hz]/[12Hz]/[40Hz].

NOTE

• Select the appropriate filter from 6 Hz, 8 Hz, 12 Hz, 40 Hz. An artifact may interfere on the BP waveform depending on the combination of BP measurement circuit.

4	Set	the	"Mean	Wave".
_				

[ON]: The mean BP waveform will be displayed and <MEAN_WAVE> will be displayed inside the numeric data box.



5 Set the "Respiration Filter".

The BP waveform baseline drift caused by the respiration influence can be prevented by setting ON the respiration filter.

▶ [ON]: Respiration Filter will turn ON.

▶ [OFF]: Respiration Filter will turn OFF.

6 Select the output signal for "IBP Analog Output".

Set the "Alarm during NIBP".

- [ON]: BP alarm will generate even during NIBP measurement.
- [OFF]: BP alarm will not generate during NIBP measurement and for 30 seconds after the measurement.

8 Set the "ART Catheter Check Message".

- [ON]: When the BP label is "ART" and the catheter is disconnected, check message will be displayed.
- [OFF]: ART catheter check message will not be displayed.

- The setting is common for all BP channels. When setting is changed for BP1, the same setting will be applied for BP2 to 8.
- The default setting of "ART Catheter Check Message" is [OFF].
- When "ART Catheter Check Message" is set to [ON], alarm will generate when the transducers are opened to air.

Select ON/OFF for parameter display.

(BCG Parameter Setup" P7-6)

- When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend input will also cease.
 - If the display of waveform/numeric data labeled as BP1/ART is set to OFF, the BP pulse

BP1	(mmHg) 116 /	77 (92)
BP1	(mmHg) 116 /	77
BP1	(mmHg)	92

BΡ

rate will not be displayed.

BP Source Selection for PR_IBP

Select the BP source for the pulse rate measurement.

The PR_IBP source can be set by displaying the PR_IBP floating window, and pressing the key for "PR_IBP Source".



Selecting [Auto] will measure the pulse rate from ART or BP1.

Pulse Oximetry

This section explains the procedure and settings of SpO_2 monitoring when the SpO_2 Unit (HS-8312N/HS-8312M/DS-8007N/DS-8007M/HG-810/HG-820) manufactured by Nellcor or Masimo is used.

When using the HG-810/HG-820, it is necessary to set the SpO_2 channel manually.

(@Maintenance Manual "Unit Module Setup" P4-17)

$SpO_2\ Monitoring$

WARNING When measuring the SpO₂ of patient with high fever or peripheral circulatory insufficiency, check the sensor attachment periodically and change the attachment site. The temperature

- For the following case, accurate measurement may not be possible.
 - Patient with excessive abnormal hemoglobin (COHb, MetHb)
 - · Patient with the pigment injected to the blood
 - Patient receiving CPR treatment
 - · When a sensor is applied to a limb with NIBP cuff, arterial catheter, or intracatheter

of the attachment site will rise due to the sensor heat which may result in burn injury.

- · When measuring at site with venous pulse
- · Patient with body motion
- · Patient with small pulse
- When a patient is receiving a photodynamic therapy, measuring SpO₂ on a same site for a long duration may cause blisters from the irradiation light of the SpO₂ sensor. Make sure to periodically change the sensor attachment site.
- Do not connect unspecified sensor or cable to any I/O connector. If done so by mistake, not
 only that the equipment cannot deliver its maximum performance, the equipment may be

damaged and safety cannot be ensured.

- If any measurement seems questionable, first check the patient's vital signs by alternate means and then check the pulse co-oximeter for proper functioning.
- For the following case, accurate measurement of SpO₂ may not be possible.
 - Improper sensor application
 - Elevated levels of COHb or MetHb: High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (CO-Oximetry) of a blood sample should be performed.
 - · Elevated levels of bilirubin
 - Elevated levels of dyshemoglobin
 - · Vasospastic disease, such as Raynaud's, and peripheral vascular disease
 - Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc.
 - Hypocapnic or hypercapnic conditions
 - Severe anemia
 - Very low arterial perfusion
 - Extreme motion artifact
 - Abnormal venous pulsation or venous constriction
 - · Severe vasoconstriction or hypothermia
 - Arterial catheters and intra-aortic balloon
 - Intravascular dyes, such as indocyanine green or methylene blue
 - · Externally applied coloring and texture, such as nail polish, acrylic nails, glitter, etc.
 - Birthmark(s), tattoos, skin discolorations, moisture on skin, deformed or abnormal fingers. etc.
- Dyes or any substance containing dyes that change usual blood pigmentation may cause erroneous readings.
- The SpO₂ data should not be used as the sole basis for diagnosis or therapy decisions. It
 must be used in conjunction with clinical signs and symptoms.
- Do not use the SpO₂ data to monitor apnea condition.
- This equipment may be used during defibrillation, but this may affect the accuracy or availability of the SpO₂ parameters and measurements.
- This equipment may be used during electrocautery, but this may affect the accuracy or availability of the SpO₂ parameters and measurements.
- The SpO₂ data cannot be used for arrhythmia analysis.
- SpO₂, SpCO, SpMet, and SpHb are empirically calibrated in healthy adult volunteers with normal levels of carboxyhemoglobin (COHb) and methemoglobin (MetHb).

- If irritation such as skin reddening appears with the sensor use, change the attachment site
 or stop using the sensor.
- When attaching the sensor with tape, do not wrap the tape too tight. At the same time, check the blood flow constantly so that congestion is not generated at the peripheral site.
- Even attachment for a short duration may inhibit the blood flow and generate compression necrosis or burn injury. Also, blood flow inhibition may prevent correct measurements.
- If SpO₂ values indicate hypoxemia, a laboratory blood sample should be taken to confirm

the patient's condition.

- If the <SpO₂ Low Perfusion> message is frequently displayed, find a better perfused monitoring site. In the interim, assess the patient and, if indicated, verify oxygenation status through other means.
- Change the application site or replace the sensor and/or patient cable when a <Replace Sensor>, <Replace Cable>, <Low Signal IQ> is displayed on the monitor. These messages may indicate that patient monitoring time is exhausted on the patient cable or sensor.
- Replace the cable or sensor when a <Replace Sensor> or <Low Signal IQ> message is consistently displayed while monitoring consecutive patients after completing troubleshooting steps listed in this manual.
- If using pulse oximetry during full body irradiation, keep the sensor out of the radiation field. If the sensor is exposed to the radiation, the reading might be inaccurate or the device might read zero for the duration of the active irradiation period.
- The device must be configured to match your local power line frequency to allow for the cancelation of noise introduced by fluorescent lights and other sources.
- Variation in hemoglobin measurements may be profound and may be affected by sampling technique as well as the patient's physiological conditions. Any results exhibiting inconsistency with the patient's clinical status should be repeated and/or supplemented with additional test data. Blood samples should be analyzed by laboratory devices prior to clinical decision making to completely understand the patient's condition.
- Check the sensor attachment site constantly in every 4 hours when probes or reusable sensor are used, and at least every 8 hours when single patient use sensors are used.Be especially careful of a patient with bad perfusion. If the sensor attachment position is not changed constantly, skin irritation or skin necrosis due to compression may be developed. For the patient with bad perfusion, check the sensor attachment position at least every 2 hours.
- As skin for neonate, premature infant is immature, change the sensor attachment site more frequently depending on the condition.
- Direct sunlight to the sensor area can cause a measurement error. Place a black or dark cloth over the sensor if using in direct sunlight.
- When not measuring, unplug the relay cable and sensor from the SpO₂ connector.Otherwise, the outside light may affect to falsely display measurements.
- If "- -" is displayed for the numeric data, make sure that the sensor is properly attached.
- · Before bathing the patient, make sure to remove the sensor and equipment from the patient.

Precautions when using the Masimo Rainbow SET Sensor

WARNING

- As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.
- Do not place this equipment or accessories in any position that might cause it to fall on the patient.
- · Do not start or operate this equipment unless the setup was verified to be correct.
- Do not use this equipment during magnetic resonance imaging (MRI) or in an MRI environment.
- Do not use this equipment if it appears or is suspected to be damaged.
- Explosion hazard: Do not use this equipment in the presence of flammable anesthetics or other flammable substance in combination with air, oxygen-enriched environments, or nitrous oxide.
- · To ensure safety, avoid stacking multiple devices or placing anything on the instrument

during operation.

- To protect against injury, follow the directions below:
 - Avoid placing the equipment on surfaces with visible liquid spills.
 - Do not soak or immerse the equipment in liquids.
 - · Do not sterilize the equipment.
 - Use cleaning solutions only as instructed in this operation manual.
 - · Do not attempt to clean the equipment while monitoring patient.
- To protect from electric shock, always remove the sensor and completely disconnect this
 equipment before bathing the patient.
- If any measurement seems questionable, first check the patient's vital signs by alternate means and then check this equipment for proper functioning.
- · Inaccurate SpCO and SpMet readings can be caused by the following.
 - Improper sensor application
 - · Intravascular dyes, such as indocyanine green or methylene blue
 - Abnormal hemoglobin levels
 - Low arterial perfusion
 - · Low arterial oxygen saturation levels including altitude induced hypoxemia
 - · Elevated total bilirubin levels
 - Motion artifact
- Inaccurate SpHb and SpOC readings can be caused by the following.
 - Improper sensor application
 - · Intravascular dyes, such as indocyanine green or methylene blue
 - Externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc.
 - Elevated PaO₂ levels
 - · Elevated levels of bilirubin
 - Low arterial perfusion
 - Motion artifact
 - Low arterial oxygen saturation levels
 - Elevated carboxyhemoglobin levels
 - Elevated methemoglobin levels
 - Hemoglobinopathies and synthesis disorders such as thalassemias, Hb s, Hb c, sickle cell, etc.
 - Vasospastic disease such as Raynaud's
 - Elevated altitude
 - Peripheral vascular disease
 - Liver disease
 - EMI radiation interference
- Inaccurate SpO₂ readings can be caused by the following.
 - Improper sensor application
 - Elevated levels of COHb or MetHb: High levels of COHb or MetHb may occur with a seemingly normal SpO₂. When elevated levels of COHb or MetHb are suspected, laboratory analysis (co-oximetry) of a blood sample should be performed.

- · Intravascular dyes, such as indocyanine green or methylene blue
- · Externally applied coloring and texture such as nail polish, acrylic nails, glitter, etc.
- · Elevated levels of bilirubin
- · Elevated levels of dyshemoglobin
- · Low arterial perfusion
- Motion artifact
- Interfering Substances: Dyes or any substance containing dyes that change usual blood pigmentation may cause erroneous readings.
- This equipment is intended only as an adjunct device in patient assessment. It should not be used as the sole basis for diagnosis or therapy decisions. It must be used in conjunction with clinical signs and symptoms.
- This equipment may be used during defibrillation, but this may affect the accuracy or availability of the parameters and measurements.
- This equipment may be used during electrocautery, but this may affect the accuracy or availability of the parameters and measurements.
- · This equipment should not be used for arrhythmia analysis.
- SpCO readings may not be provided if there are low arterial saturation levels or elevated methemoglobin levels.
- SpO₂, SpCO, SpMet, and SpHb are empirically calibrated in healthy adult volunteers with normal levels of carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- Do not adjust, repair, open, disassemble, or modify this equipment or accessories. Injury to
 personnel or equipment damage could occur. Return this equipment for servicing if
 necessary..

- Do not place this equipment where the controls can be changed by the patient.
- Electrical shock and flammability hazard: Before cleaning, always turn off the equipment and disconnect from any power source.
- When patients are undergoing photodynamic therapy, they may be sensitive to light sources. Pulse oximetry may be used only under careful clinical supervision for short time periods to minimize interference with photodynamic therapy.
- Do not place this equipment on electrical equipment that may affect the operation, preventing it from working properly.
- If SpO₂ values indicate hypoxemia, a laboratory blood sample should be taken to confirm the patient's condition.
- If the <SpO₂ Low Perfusion> message is frequently displayed, find a better perfused monitoring site. In the interim, assess the patient and, if indicated, verify oxygenation status through other means.
- If using this equipment during full body irradiation, keep the sensor out of the radiation field.
 If the sensor is exposed to the radiation, the reading might be inaccurate or the equipment might read zero for the duration of the active irradiation period.
- The equipment must be configured to match your local power line frequency to allow for the cancelation of noise introduced by fluorescent lights and other sources.
- To ensure that alarm limits are appropriate for the patient being monitored, check the limits each time this equipment is used.
- Variation in hemoglobin measurements may be profound and may be affected by sampling technique as well as the patient's physiological conditions. Any results exhibiting inconsistency with the patient's clinical status should be repeated and/or supplemented with

additional test data. Blood samples should be analyzed by laboratory instruments prior to clinical decision making to completely understand the patient's condition.

- Do not submerge this equipment in any cleaning solution or attempt to sterilize by autoclave, irradiation, steam, gas, ethylene oxide or any other method. This will seriously damage this equipment.
- Electrical Shock Hazard: Carry out periodic tests to verify that leakage currents of patientapplied circuits and the system are within acceptable limits as specified by the applicable safety standards. The summation of leakage currents must be checked and in compliance with IEC 60601-1 and UL 60601-1. The system leakage current must be checked when connecting external equipment to the system. When an event such as a component drop or a spillage of blood or other liquids occurs, retest before further use. Injury to personnel could occur.
- Disposal of product Comply with local laws in the disposal of the equipment and/or its accessories.
- To minimize radio interference, other electrical equipment that emits radio frequency transmissions should not be in close proximity to this equipment.

NOTE

- A functional tester cannot be used to assess the SpO₂ accuracy.
- High-intensity extreme lights (such as pulsating strobe lights) directed on the sensor, may not allow this equipment to obtain SpO₂ readings.
- When using the maximum sensitivity setting, performance of the "Sensor Off" detection may be compromised. If the equipment is in this setting and the sensor becomes dislodged from the patient, the potential for false readings may occur due to environmental "noise" such as light, vibration, and excessive air movement.
- · Changes or modifications shall void the warranty for this equipment.

NOTE

- SpCO, SpMet, SpHb, PI, and PVI are parameters which can be measured by the Masimo unit.
- SpOC can be measured only on the DS-8007M.
- High-intensity extreme lights (such as pulsating strobe lights) directed on the sensor, may not allow this equipment to obtain SpO₂ readings.
- Do not loop the patient cabling into a tight coil or wrap around the equipment, as this can damage the patient cabling.
- Additional information specific to the Masimo sensors compatible with this equipment, including information about parameter/measurement performance during motion and low perfusion, may be found in the sensor's directions for use

Prepare an appropriate probe or sensor for the patient.

(Pulse Oximetry Measurement (Manufactured by Covidien)" P13-3)

(@"Pulse Oximetry Measurement (Manufactured by Masimo)" P13-4)

Connect the sensor to HS-8000 or module.

In Case of NellcorTM Unit:

1 Connect the DOC-10 SpO₂ Relay Cable to the SpO₂ connector on the HS-8312N or DS-8007N, HG-820. The illustration is example of connection with HS-8312.



2 Insert the sensor into the SpO₂ relay cable connector, and lock it with the transparent cover.



In Case of Masimo Unit:

- 1 Connect the SpO₂ patient cable (RD[®], LNCS[®], M-LNCSTM, Rainbow[®]) to the SpO₂ connector on the HS-8312M or DS-8007M, HG-810.
- 2 Connect the patient cable and the sensor.

Face the metallic side of the sensor upward and align the logo with that of the patient cable. Then, insert the sensor connector to the patient cable until a click sound is heard.

The SpO₂ patient cables (RD[®], LNCS[®], M-LNCSTM, Rainbow[®]) are for Masimo SET sensor, Rainbow SET sensor only. Connect them only to the HS-8312M or DS-8007M, HG-810. Otherwise, the equipment will not properly function.

NOTE

- Pull the connector slowly to ensure it is securely connected.
- · If necessary, fixate the cable to the patient.

3 Attach the sensor to the patient.

l 🗘 CAUTION

• If the nail is rough, dirty, or manicured, accurate measurement will not be possible. Change the finger or clean the nail before attaching the probe and sensor.

Probe Type Sensor

1 As shown below, the probe cable should be on the nail side.

- 1 Light Emitting Part
- 2 Light Receiving part



 $2\,$ Adjust the sensor so that the light-

emitting part (on cable side) is over the nail, or as instructed per the related sensor instruction manual.



3 Press the probe lightly so that the finger and the rubber cover are appressed. This is to stabilize the probe, and to avoid ambient light.



Single-use Type

- 1 Clean the attachment site with alcohol, etc., and allow to dry.
- 2 Align the light emitting element and light receiving element of the sensor with the measuring site in between when attaching the sensor to patient.
 - 1 Light Emitting Element
 - 2 Light Receiving Element



3 Secure the cable with surgical tape so that the sensor does not come off when the cable is pulled.





Attachment to the toe

Attachment to the finger

 $\mathbf{4}$ Verify that the SpO₂ measurement and SpO₂ waveform are displayed on the home display.



SpCO, SpMet, SpHb, SpOC Measurement (Masimo)

This section explains the SpCO, SpMet, SpHb, SpOC measurement procedure when using the HS-8312M or DS-8007M, HG-810.



REFERENCE

• SpCO, SpMet, SpHb, SpOC measurements are optional function.

Select the Rainbow sensor for the patient.

(@"Pulse Oximetry Measurement (Manufactured by Masimo)" P13-4)

The measurement procedure is the same with that of the SpO₂. Verify that the SpCO, SpMet, SpHb, SpOC value is displayed on the monitor. (PSpO2 Monitoring" P7-33)

Precautions about the Masimo Sensors and Cables

A technology called X-Cal for patient safety and reinforcement of efficiency in a clinical site is implemented for Masimo sensors and cables.

X-Cal is designed to address the following three common factors that can impact measurement accuracy and patient safety due to reliability risks.

- 1 Imitation Masimo sensors and cables
- 2 Cables and sensors used far beyond their expected life
- 3 Third-party reprocessed pulse oximetry sensors

If a sensor or cable that does not support X-Cal is used with an X-Cal enabled device, SpO₂ measurement will not be available.

Even if Masimo sensors or specified sensors and cables are used, SpO₂ measurement may not be available if the sensors and cables are used beyond their expected life.

About the Expected Life of Sensors and Cables

- The X-Cal function automatically monitors the Active Monitoring Time (actual time of SpO₂ monitoring) for each sensor and cable. If the sensors and cables are used beyond the expected life, the message, <Replace Cable> or <Replace Sensor> will be displayed.
- The measurement will not cease until it is completed even if the cable or sensor has reached end of life during the measurement.
- When a measurement with cable or sensor that has reached end of life is suspended for certain amount of time, and resumed with the same cable or sensor, a message to replace the sensor or cable will be displayed.
- The sensor or cable that has reached end of life needs to be replaced before resuming monitoring.
- The following table shows the expected life of cable and sensor. The indication of usage hours per day (24 hours/12 hours/8 hours) are also shown.

Sensors or Cables	Expected Life	When monitoring 24 hours/day	When monitoring 12 hours/day	When monitoring 8 hours/day
Single Patient Use SpO ₂ "L" Sensor with replaceable tape	336 hours	14 days	28 days	42 days
Single Patient Use SpO ₂ Sensor	168 hours	7 days	14 days	21 days
Reusable SpO ₂ Sensor (DCI, DCIP, YI, TF-I, DBI)	8,760 hours	12 months	2 years	3 years
Patient Cable	17,280 hours	24 months	4 years	6 years

Active Monitoring Time (actual time of monitoring)

SpO₂ Parameter Setup (Nellcor)

This section explains the measurement procedure when using the HS-8312N or DS-8007N, HG-820. Press the [Menu], [SpO₂] keys to display the "SpO₂" setup screen.



When Using the HS-8312N

Select the waveform size from [x1/4]/ [x1/2]/ [x1]/ [x2]/ [x4]. (shown on right)

m Z Select the label from [None]/ [Auto]/ [RH]/ [LH]/ [RF]/ [LF]/ [OT].

 When [Auto] is selected, the label will be automatically assigned depending on the SpO₂ unit type and channel number.
 Nellcor 1ch: N1, 2ch: N2
 Masimo 1ch: M1, 2ch: M2



3 Set the SpO₂ alarm.

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Whether to use the second alarm function and its threshold selection should be based on the patient's clinical indication portent and medical evaluation.
- Set the upper limit in the range of 51%SpO₂ to 100%SpO₂. If a value above 100%SpO₂ is set, the upper alarm will turn OFF.
- Set the lower limit in the range of 50%SpO₂ to 99%SpO₂. If a value below 50%SpO₂ is set, the lower alarm will turn OFF.

REFERENCE

- Also, when the SpO₂ value is unstable around the lower alarm limit, the frequently generated alarm can be corrected by setting the second alarm function.
 (P6-2)
- When [Auto] is set, the upper limit will be turned OFF and the lower limit will be set to 90%SpO₂.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Baintenance Manual "Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.
- The upper/lower limit can be set in 1%SpO2 increment.
- Indicates the current measurement value.
- The following delay occurs for the SpO2 alarm depending on the patient classification and

second alarm setting. (For Nellcor)

	Second Alarm	Patient Cla	assification
	Setup	Adult/Child	Neonate
SpO ₂ Alarm Condition Delay	For all settings	About 7 sec. to 9 sec.	About 7 sec. to 9 sec.
SpO ₂ Alarm Signal Generation Delay	OFF	About 5 sec.	0 sec.
	10	About 5 sec. to 7 sec.	About 5 sec. to 7 sec.
	25	About 11 sec. to 13 sec.	About 11 sec. to 13 sec.
	50	About 19 sec. to 22 sec.	About 19 sec. to 22 sec.
	100	About 36 sec. to 38 sec.	About 36 sec. to 38 sec.

4 Set the Ext SpO_2 alarm.

(@"Alarm Limit Setup for Each Parameter" P6-10)

NOTE

 Set the lower limit in the range of 50%SpO₂ to 98%SpO₂. If a value below 50%SpO₂ is set, the lower alarm will turn OFF.

REFERENCE

- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Plantenance Manual "Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.
- When [Auto] is set, the lower limit will be set to "SpO₂ lower limit 10%SpO₂".
- The lower limit can be set in 1%SpO₂ increment.
- **b** indicates the current measurement value.
- The following delay occurs for the Ext SpO₂ alarm depending on the patient classification and second alarm setting. (For Nellcor)

	Patient Classification		
	Adult/Child	Neonate	
SpO ₂ Alarm Condition Delay	About 7 sec. to 9 sec.	About 7 sec. to 9 sec.	
SpO ₂ Alarm Signal Generation Delay	About 5 sec.	0 sec.	

5 Set the PR alarm.

(@ "Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Set the upper limit in the range of 22 bpm to 300 bpm. The upper limit alarm will become OFF if the value exceeds 300 bpm.
- Set the lower limit in the range of 20 bpm to 295 bpm. If a value below 20 bpm is set, the lower alarm will turn OFF.

REFERENCE

 To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Plantenance Manual "Alarm Related Setup" P5-4)

- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.
- When [Auto] is set, the upper and lower limit will be automatically set to +40 bpm and -40 bpm to the current value respectively.
- The upper and lower limit can be set in 5 bpm increments. It can be set in 1 bpm increment for 25 bpm and below.
- The following delay occurs for the PR alarm depending on the patient classification. (For Nellcor)
 - PR Alarm Condition Delay: <Adult/Child/Neonate> About 5 sec. to 6 sec.
 - PR Alarm Signal Generation Delay: <Adult/Child> About 5 sec., <Neonate> 0 sec.

6 Set the "Alarm during NIBP".

NOTE

- During the NIBP measurement, the cuff inflation restricts the blood flow which disables the correct detection of the SpO₂ and PR, and may generate an improper alarm.
- Selecting [OFF] for "Alarm during NIBP" will not generate the SpO₂, Ext SpO₂, PR, SpCO (Masimo only), SpMet (Masimo only), SpHb (Masimo only) alarm until the NIBP measurement is complete.

REFERENCE

- This setup can be used when the SpO₂ sensor and the NIBP cuff is placed on the same limb for measurement.
- [ON]: Alarm will be generated even during NIBP measurement.
- ▶ [OFF]: SpO₂/PR alarm will not be generated during NIBP measurement.

8 Select ON/OFF for parameter display.

(@"ECG Parameter Setup" P7-6)

- When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend input will also cease.
- When the waveform and numeric data display is set to OFF, the pulse rate measured by SpO₂ will not be displayed either.

REFERENCE

When SpO₂ sensor is attached to the patient with the SpO₂ display set to OFF, and SpO₂ is measured for 10 seconds, the pulse wave and numeric data will be automatically displayed.

RR_SpO2 Parameter Setup (Nellcor)

This section explains the RR_SpO₂ measurement procedure when using the DS-8007N.

- The RR_SpO₂ can be measured only when using the Nellcor Respiratory Sensor.
- · For details, contact your nearest service representative.

1 Prepare the sensor.

(@"Pulse Oximetry Measurement (Manufactured by Covidien)" P13-3)

The measurement procedure is the same with that of the SpO₂. Verify that the RR_SpO₂ value is displayed on the monitor. (@ "SpO2 Monitoring" P7-33)

SpO₂ Parameter Setup (Masimo)

This section explains the measurement procedure when using the HS-8312M or DS-8007M, HG-810. Press the [Menu], $[SpO_2]$ keys to display the "SpO₂" setup screen.

- REFERENCE
- This setting is available when using the HS-8312M or DS-8007M, HG-810. PVI, SpCO, SpMet, SpHb, SpOC measurements are optional function.



SpO₂ averaging duration setting. (For Masimo)

	SpO ₂ Averaging	Patient Classification	
		Adult/Child	Neonate
SpO ₂ Alarm Condition Delay	For all settings	About 7 sec. to 9 sec.	About 7 sec. to 9 sec.
SpO ₂ Alarm Signal Generation Delay	For all settings	About 5 sec.	0 sec.

4 Set the Ext SpO_2 alarm.

(P7-43) (P7-43)

5 Set the PR alarm.

(@"SpO2 Parameter Setup (Nellcor)" P7-43)

REFERENCE

- The following delay occurs for the PR alarm depending on the patient classification. (For Masimo)
 - PR Alarm Condition Delay: <Adult/Child> About 8 sec. to 10 sec. <Neonate> About 7 sec. to 9 sec.
 - PR Alarm Signal Generation Delay: <Adult/Child> About 5 sec., <Neonate> 0 sec.

6 Set the "Alarm during NIBP".

Set the "Synchronized Mark/Tone".

8 Select the averaging duration from [2-4 sec.]/[4-6 sec.]/[8 sec.]/[10 sec.]/[12 sec.]/[14 sec.]/[16 sec.].

WARNING

 Be careful when setting the "SpO₂ Averaging" duration as the SpO₂ alarm is based on the displayed SpO₂ value which is averaged from the duration set in "SpO₂ Averaging". The alarm occurrence time will be affected or may not occur for the transient value of SpO₂ depending on the set duration.

 ${f 9}$ Select the pulse detection sensitivity from [High]/[Normal]/[APOD].

• If [High] is selected for pulse sensitivity, probe-off detection will become somewhat inaccurate.

NOTE

- To improve the low perfusion condition, or to perform fast tracking when the SpO₂ value changes abruptly, select [High].
- For standard use, select [Normal].
- If there is a high possibility of sensor getting disconnected, select [APOD].
- · When the pulse detection sensitivity is set to [High], performance of the "Sensor Off"

detection may be compromised. If the equipment is in this setting and the sensor becomes dislodged from the patient, the potential for false readings may occur due to environmental "noise" such as light, vibration, and excessive air movement.

10 Set the "FAST SAT".

NOTE

- To pick up the abrupt change of the value sooner, and to take advantage of the qualities of FAST SAT mode, SpO₂ averaging time will be fixed as [2-4 sec.] when FAST SAT is set ON.
- ▶ [ON]: Abrupt change of the SpO₂ value can be monitored.
- ▶ [OFF]: FAST SAT mode will turn OFF.

11 Set the "PI (Perfusion Index) Display".

NOTE

- The perfusion index is calculated by pulsatile signal divided by apulsatile signal times 100, and indicates patient's circulation condition at the monitoring site.
- The perfusion index assists clinicians in determining optimal placement of the SpO₂ sensor. This parameter is also useful as a troubleshooting tool by helping a clinician rule out whether a questionable value may be due to low perfusion and/or a low signal to noise condition. An added benefit is that changes in perfusion can be an indicator to the clinician of important changes in the patient's physiological status.
- ▶ [ON]: PI will be displayed.



▶ [OFF]: PI will not be displayed.



REFERENCE

- Perfusion Index (PI) is a relative assessment of the pulse strength at the monitoring site. It is a ratio of the pulsatile and the non-pulsatile blood flow at the monitoring site. It can be used to find the most appropriate sensor application site by finding the site with the highest PI. Perfusion Index (PI) is displayed in the range from 0.02% to 20%.
- Pleth Variability Index (PVI) is an index of the change in PI that occurs during the respiratory cycle. It is calculated by measuring the changes in PI over a time interval where one or more complete respiratory cycles have occurred. Pleth Variability Index (PVI) is displayed in the range from 0% to 100%.

12 Select [ON]/[OFF] for "Signal IQ Wave".



• The signal IQ wave cannot be printed.

REFERENCE

· The signal IQ wave indicates the signal confidence and pulse beat. The vertical length

indicates the signal confidence. A low vertical line indicates a lower signal confidence.



13 Select ON/OFF for parameter display.

(Pr-43) (

14 Set the SpCO alarm.

Press the [], [Sp*], [SpCO] keys to display the SpCO alarm setup screen.



- Set the upper limit in the range of 1%SpCO to 40%SpCO. If a value above 40%SpCO is set, the upper alarm will turn OFF.
- The lower limit cannot be set.
- The automatic alarm cannot be set.

15 Set the SpMet alarm.

Press the [SpMet] key to display the SpMet alarm setup screen. Set the alarm in the same procedure as SpCO.



- Set the upper limit in the range of 1%SpMet to 15%SpMet. If a value above 15%SpMet is set, the upper alarm will turn OFF.
- The lower limit cannot be set.
- The automatic alarm cannot be set.

16 Set the SpHb measurement condition. Press the [SpHb] key to display the SpHb setup screen.



- 1 Select the SpHb averaging duration from [Short] / [Medium] / [Long].
- 2 Set the SpHb alarm.

- Set the upper limit in between 2.0 g/dL to 24.5 g/dL. The upper limit alarm will turn OFF if the value above 24.5 g/dL is set.
- Set the lower limit in between 1.0 g/dL to 24.0 g/dL. The lower limit alarm will turn OFF if the value below 1.0 g/dL is set.
- The automatic alarm cannot be set.

Non-Invasive Blood Pressure

The procedure of NIBP measurement and measurement condition setup are explained.

- · For the following situation, measurements will be terminated.
 - When the measurement time has exceeded 160 seconds for adult and child, 80 seconds for neonate.
 - When the inflation value has exceeded 300 mmHg for adult, 210 mmHg for child, and 150 mmHg for neonate.
- If used with the incorrect patient classification, it will not only cause erroneous measurement, but the inflating level for the adult may be applied to child or neonate causing dangerous situation to the patient.
- · The following factors may affect the NIBP value.
 - · Body motion, arrhythmia, convulsion, low pulse pressure, slow pulse
 - · Continuous noise such as cardiac massage
 - · Noise from the electrosurgical instrument

Lineup of Cuffs

REFERENCE

- According to the AHA (American Heart Association) guideline, the appropriate cuff width is 40% of the arm circumference.
 - Select the appropriate cuff from the following selections.
 - For other usable cuffs, refer to the section on "Optional Accessories".
 - (P"Non-Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)" P13-2)

NIBP Monitoring

WARNING

 Before the NIBP measurement, make sure the patient classification ([Adult]/[Child]/ [Neonate]) is properly selected on the "Admit/Discharge" menu. Otherwise, correct measurement cannot be performed, and congestion or other injury may result.

- Correct NIBP measurement cannot be performed if oxygenator is used or if the pulse is difficult to detect.
- Pay attention when measuring the NIBP of patient with bleeding disorders or hypercoagulation. The cuff inflation constricting the arm may cause petechia or circulatory failure with blood clot.
- Do not apply the cuff to the arm or thigh where vein is secured. The blood may backflow causing the chemical injection to cease.
- Properly arrange the cuff and air hose.
- Check the condition of cuff-applied part on the patient during measurement so that the blood circulation will not be blocked over long period of time by the squashed or bent cuff hose.
- Check the patient's condition constantly while measuring over a long period of time with interval of 2.5 minutes or less. Also, periodically check the blood circulation while performing periodic measurement over a long period of time. Congestion or rash may occur at the measuring site.
- Make sure to check the patient's condition constantly when repeatedly using continuous measurement as it may cause dysfunction of patient's circulation.
- When the cuff is not applied to the patient, pay attention not to leave the cuff unattended. If
 periodic or continuous measurement is set, the cuff will automatically inflate and may cause
 the rubber bag inside the cuff to burst. When not performing the NIBP measurement, set the
 NIBP measurement interval OFF and disconnect the air hose from the NIBP connector.
- · The following factors may affect the NIBP value.
 - · Body motion, arrhythmia, convulsion
 - · Continuous noise such as cardiac massage
 - Periodic electromagnetic noise
- If the cuff inflation may adversely affect the patient's blood flow or wound, attach the cuff to an appropriate position under physician's instruction.
- Do not apply the NIBP cuff to the arm of the mastectomized side. It may cause swelling or

other circulatory failure.

- It is not intended for measuring the NIBP of pregnant patient, including pre-eclamptic. It may cause incorrect NIBP measurement.
- Pay attention when measuring the NIBP of pregnant (including pre-eclamptic) patient. It may affect the NIBP value.

NOTE

• When the [NIBP Start/Stop] key is pressed or when the NIBP measurement interval is changed, the standby mode will be canceled and the NIBP periodic measurement will start.

1 Select the appropriate cuff type for the patient.

(@"Non-Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)" P13-2)

▲ CAUTION

- Select the appropriate cuff size which best fits the arm circumference. If the cuff size is inappropriate, it may cause measurement error.
- Do not use a cuff which is worn out. The cuff may burst during inflation.

 $\mathbf{2}$ Connect the cuff to the air hose.



 ${f 3}$ Connect the air hose to the NIBP connector on the HS-8000 or module.



• Make sure that the cuff hose connection is secure. If there is any air leakage, correct NIBP measurement cannot be performed. (NOTE)

• The neonate cuff should be connected to air hose for neonate. Other cuffs should be connected to air hose for general use.

The HS-8000 automatically determines the patient classification (neonate or adult/child) according to the connected air hose. If the air hose is not connected to the cuff connector, the measurement will not start.

4 Apply cuff to the patient.



NOTE

- Position the ARTERY T mark over the artery on the patient's arm and wrap the cuff around.
- One or two fingers should just fit in between the cuff and arm.

REFERENCE

 Align the cuff height and heart position to eliminate an error caused by the blood weight. It is most appropriate to measure with the patient lying down and arms naturally extended.



5 Press the [NIBP Start/Stop] key (user key or fixed key).



- Cuff inflation and measurement will start.
- Upon completion, the measured value will be displayed inside the NIBP numeric data box. The measurement can be also started by pressing the [NIBP Start/Stop] key on the HS-8000. The blue LED will light during the measurement.

After the measurement, the LED will turn OFF, a beep tone will generate for 1 second and the measurement result will be displayed on the monitor.

REFERENCE

- About the Oscillometric Method
- The oscillometric method measures the blood pressure by detecting the pulse oscillation change by the cuff pressure. The cuff connects to the NIBP connector via the air hose. The air pressure inside the cuff is converted to voltage by the pressure sensor, converted to digital signal (A/D conversion), and transmitted to the CPU. The measurement is performed with the following process.
 - The cuff inflates to the set value and inhibits the arterial blood flow at the measured site.
 - The cuff gradually deflates.
 - The arterial blood flow of the patient will return when the cuff pressure is decreased sufficiently.
 - The oscillation (pulse signal) caused by the restricted blood circulation is transmitted to the pressure sensor via the air hose, and converted to an electric signal.
 - From the pulse signal and cuff pressure detected at the pressure measurement circuit, the systolic, diastolic, average blood pressure and pulse rate will be measured at the CPU.
- The systolic, diastolic, mean blood pressure will be displayed on the monitor. The measurement will start with the following factor.
 - · When the [NIBP Start/Stop] key (fixed key or user key) is pressed.
 - At the selected measurement interval.
 - · For fixed amount of time after the NIBP Cont. key (user key) is pressed. (Max. 15 min.)
 - If "NIBP Measurement at Alarm Occurrence" is set ON, and the set parameter generates an alarm.
 - When the change in patient's circulation condition is detected from the time difference of ECG and SpO₂ waveform.

Inflation Mode Setup

The maximum inflation value and measurement duration needs to be changed according to the patient classification. The inflation mode will automatically change according to the patient classification setting. Set the appropriate patient classification on "Admit/Discharge" menu or "Detail Setup" menu under NIBP parameter setup.

The NIBP measurement on this equipment is provided with forced exhaust system for safety purpose. When the maximum inflation value is reached or when the fixed measurement duration is exceeded, the system will automatically start to exhaust. The maximum inflation value, maximum measurement duration, initial inflation value, measurement range, and alarm limit range for this exhaust system is set according to the patient classification setting.

Patient Classification	Target Inflation Value	Maximum Inflation Value	Maximum Measurement Duration
Adult	100 mmHg to 290 mmHg (Default: 180 mmHg)	300 mmHg	160 sec.
Child	100 mmHg to 200 mmHg (Default: 140 mmHg)	210 mmHg	160 sec.
Neonate	100 mmHg to 140 mmHg (Default: 110 mmHg)	150 mmHg	80 sec.
NIBP Auto Mode Setup

Non-invasive blood pressure can be measured automatically at selected time intervals.

If continuous measurement is started during the NIBP auto mode, the auto mode will automatically resume when the continuous measurement completes.



Press the [NIBP Auto Mode] key on the home display.

> The "NIBP Auto Mode" window will be displayed.



 ${f 2}$ Select the measurement interval from the displayed selection.

- When [1min] is selected, the 1-minute interval measurement will start from the time the selection is made.
- The 1-minute interval measurement will automatically stop after 12 minutes (maximum of 15 minutes when re-measured), and 2.5-minutes interval measurement will start.
- The continuous mode will continuously measure for 12 minutes (maximum of 15 minutes when re-measured). When the measurement completes, 2.5 minute interval measurement will start.
- When "Auto Mode with Start/Stop Key" is set to [ON], the auto mode measurement needs to be started manually.
- When using the continuous mode or Lumbar mode for measurement, make sure that the setting is according to the intended purpose.
 (P7-56)
- The Lumbar mode is recommended for use during spinal anesthesia. It should be used with sufficient safety measures.

NOTE

- 1-minute interval measurement cannot be stopped by pressing the [NIBP Start/Stop] key (fixed key or user key). To stop the 1-minute interval measurement, select [OFF] or other interval on "NIBP Auto Mode" window.
- When the NIBP auto mode interval is [Cont.]/[1min]/[2min]/[2.5min]/[5min]/[Lumbar Mode], NIBP measurement cannot be started from the central monitor.

- When the [NIBP Start/Stop] key is pressed or when the NIBP measurement interval is changed, the standby mode will be canceled and the NIBP periodic measurement will start.
- > The measurement will automatically start at selected interval.
- The selected interval will be displayed inside the numeric data box.



REFERENCE

- · Select [OFF] if not performing the auto mode measurement.

(57	INIDP I	Parameter	Setup	P7-59)	

• When [60min]/[120min] is selected for the measurement interval, the measurement will start 5 minutes before the set time. If outputting the data to PC or other external device using the PC communication function of this system, an error may be generated to the NIBP measurement time depending on the input interval of the external device. This system outputs the data at completion of NIBP measurement, and if the external device inputs the data at 60 minutes interval, 60 minutes time lag will occur. By starting the measurement 5 minutes early, this time lag between the external device can be minimized.

[Menu > Initial Settings > Meas. > Other]



• On the "Initial Settings", whether or not to backup the NIBP measurement interval at discharge/power ON can be selected. (OFF/Backup/OFF→2.5min./OFF→5min.)

About the Lumbar Mode

The Lumbar mode is intended for use during spinal anesthesia. The Lumbar mode performs the measurement as follows.



If [Lumbar] is selected when the measurement is not performed, the first measurement will start.

If [Lumbar] is selected during the measurement, the current measurement will be counted as the first measurement. The second measurement will start after 1 minute, and after 7 times of 2-minute interval measurement, the Lumbar mode will end. The Lumbar mode can be manually stopped by selecting other interval or selecting [Lumbar] again. When the Lumbar mode ends, 5-minute interval measurement will automatically start.

- Pressing the [NIBP Start/Stop] key during measurement will only stop the measurement and not the Lumbar mode. To stop the Lumbar mode, select other interval or select [Lumbar] again.
- The manual measurement can be performed in between the Lumbar mode measurement. The Lumbar mode measurement will not start if the manual measurement is still in progress when the next Lumbar mode measurement time arrives.

Oscillation Graph Display



The description of the oscillation graph is as follows.

The horizontal axis shows the cuff pressure, and vertical axis shows the pulse amplitude with reference to maximum pulse amplitude.

The bar graph shown at left indicates the size of maximum pulse amplitude compared with the reference value. For example, if the maximum pulse amplitude is 1/2 of the reference value, the bar graph will be half filled in.

- 1 Bar Graph
- 2 DIA Value
- 3 SYS Value



Dyna Alert Function Status

The Dyna Alert function is a technology to prevent accidents which may occur by sudden BP change during the nonmeasured duration by estimating the variation of circulatory dynamics.

This function is available for the HS-8312N, or DS-8007N with the Nellcor SpO₂ module.

When [ON] is selected for "Dyna Alert", NIBP measurement will automatically start when the Dyna Alert estimated value exceeds the alarm limit. The function will activate with the following condition.

(P"Dyna Alert" P7-60)

- Patient Classification: Adult (20 kg or above)
- Cuff Applied Site: Upper Arm
- SpO₂ Sensor Attachment Site: Fingertip
- NIBP Measurement Interval: 5 minutes to 60 minutes

- When the SpO₂ sensor is applied to the toe or forehead, the circulatory dynamics variation monitoring by the Dyna Alert may not properly function.
- The circulatory dynamics variation monitoring by the Dyna Alert is effective only on the HS-8312N, or DS-8007N with the Nellcor SpO₂ module.

In the NIBP numeric data box, the following mark and message indicating the status of the Dyna Alert function will be displayed.



D.Alert Color of Mark	Message	Status	Dyna Alert Function Status ^{*1}
Gray	DA Setup: OFF	Dyna Alert (DA) is set to OFF.	Disable
	Patient: Child	NIBP measurement is performed on child.	Disable
	Patient: Neonate	NIBP measurement is performed on neonate.	Disable
	Pacemaker: ON	Pacemaker setting is set to ON.	Disable
	Interv.: <5min.	NIBP interval is set to Cont., 1min, 2min, or 2.5min.	Suspended
	Interv.: >60min.	NIBP interval is set to 120min.	Suspended
	Interv.: OFF	NIBP interval is set to OFF.	Suspended
	Measuring BP ^{*2}	Invasive blood pressure is measured.	Suspended
Yellow	Measure NIBP	Initialization of Dyna Alert is complete, and the NIBP measurement has not been performed since the power is turned ON.	Suspended
	Poor ECG Signal	ECG signal failure due to lead-off, noise, etc.	Disable
	Poor PTG Signal	PTG (Photoplethysmograph) signal failure due to sensor off, noise, severe low perfusion, etc.	Disable
	DA-NIBP Suspended	Within 2.5 minutes from previous Dyna Alert NIBP measurement.	Suspended
	Measuring NIBP	NIBP measurement other than Dyna Alert is in progress.	Disable
	Initializing	Waiting for stable signal after starting Dyna Alert.	Disable
Green	PTG Low Perfusion	PTG amplitude is 200 unit or above, and below 800 unit.	Enable
	Mon. BP Var.	Dyna Alert is properly monitoring circulatory dynamics variation.	Enable
Pink	Measuring DA-NIBP	Dyna Alert NIBP measurement is in progress.	Disable

*1:

*2:

Disable: Circulatory dynamics variation is not monitored.

Suspended: Circulatory dynamics variation is monitored. But the display suspends the measurement when NIBP measurement is requested. When the suspending factor is resolved, the measurement will resume as quickly as possible.

Enable: Circulatory dynamics variation is monitored. The display control software responds to NIBP measurement request as quickly as possible.

"Measuring BP" indicates the status when IBP (BP1 or ART) measurement is possible and can be displayed on the monitor.

- When using the Dyna Alert function, be aware of these risks and do not increase the NIBP interval time by relying only on the Dyna Alert function.
- After the Dyna Alert NIBP measurement, the next Dyna Alert NIBP measurement cannot be performed for 2.5 minutes.
- The Dyna Alert will not properly function for the following cases.
 - If peripheral circulatory insufficiency or very low BP is developed.
 - If highly-frequent arrhythmia is generated.
 - If an oxygenator is used.
 - · If a large noise from body movement or electric surgery equipment is interfering.
 - If autonomic nerve or circulatory dynamics is largely affected by medication.

NIBP Parameter Setup

Press the [Menu], [NIBP] keys to display the "NIBP" setup screen.



5 ECG RESP NIBP BP SpO2 TEMP 1 Ø NIBP Auto Hode NIBP S 0 asure at Alam Quick Neasur 滋 × 180 Sight Inflat: Cancel Error WAF ON 80 ON End Tone ٩Þ

1 NIBP Alarm

(@ "Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Set the upper limit in the range of 15 mmHg to 300 mmHg / 2.0 kPa to 40.0 kPa. If a value above 300 mmHg / 40.0 kPa is set, the upper alarm will turn OFF.
- Set the lower limit in the range of 10 mmHg to 295 mmHg / 1.5 kPa to 39.5 kPa. If a value below 10 mmHg / 1.5 kPa is set, the lower alarm will turn OFF.

REFERENCE

- Set ON/OFF of NIBP alarm, upper and lower alarm limits of systolic (S), diastolic (D), mean (M) NIBP.
- When [Auto] is set, the upper and lower limit will be automatically set to +40 mmHg / +5 kPa and -20 mmHg / -3 kPa respectively to the current value.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.
- The alarm limit should be set for each unit (mmHg/kPa).
- The upper/lower limit can be set in 5 mmHg / 0.5 kPa increment.

2_{NIBP} Auto Mode</sub>

NIBP measurement will be performed automatically at selected time intervals.

(PT-55) (

Patient Classification

The patient classification setting is linked with that on the "Admit/Discharge" screen. The inflation value and measurement duration will differ according to the patient classification setting.

(@"Inflation Mode Setup" P7-54)

WARNING

- The patient classification selection influences the precision of the QRS detection and NIBP measurement. Make sure the correct selection is made.
- The NIBP air hose corresponded to the set patient classification must be used to perform NIBP measurement. However, if the patient classification is child, NIBP air hose for adult can be used.

4 Dyna Alert

[ON]: Dyna Alert function will turn ON when HS-8312N or DS-8007N is used.





- When the PTG (SpO₂) sensor is applied to the toe or forehead, the circulatory dynamics variation monitoring by the Dyna Alert may not properly function.
- The circulatory dynamics variation monitoring by the Dyna Alert is effective only on the

HS-8312N or DS-8007N with the Nellcor SpO₂ module.

REFERENCE

- About the Dyna Alert:
 - Using a cuff allows to measure the blood pressure noninvasively, but on the other hand, there is a demerit of not being able to perform the measurement continuously.Therefore, there is always a risk of sudden blood pressure change in between the periodic measurements.

5 Oscillograph

[ON]: Oscillation graph will be displayed inside the numeric data box.

[Oscill. Print] key will be also displayed.

[Oscill. Print]: Oscillation graph will be output on the HR-800 Recorder Unit.

[OFF]: Oscillation graph will not be displayed.

[Real Time]: Oscillation graph will be updated during the measurement.

NOTE

 The oscillation graph can be displayed when the NIBP numeric data box size is 2-box size or larger, and "Oscillograph" is set to [ON] on the "NIBP" setup screen.

O PR Display

[ON]: PR will be displayed.

NOTE

 PR will be displayed only. It will not generate alarm, or be displayed for the tabular trend.

NIBP Erase Time

NIBP data will be erased after the set duration (60min/120min).

OMeasure at Alarm

NIBP measurement will start at alarm generation.

Select [ON] for "NIBP Measurement at Alarm Occurrence", and select the alarm factor to start the NIBP measurement. Multiple parameters can be selected.

CAUTION

• If the NIBP measurement has not been performed since the power was turned ON, NIBP measurement at alarm occurrence will not be performed.

Quick Measurement

[ON]: NIBP measurement will be performed in duration of about 20 seconds to 25 seconds in case of adult patient.

NOTE

• The quick measurement can be performed only if the patient classification is adult or child. For neonate, normal measurement will be performed regardless of this setting.

10 Sight Inflation

[ON]: Sight inflation function will turn ON.

The inflation target level will be automatically estimated during the inflation, and starts to deflate after the target level is reached.

If [ON] is selected for "Sight Inflation", the target inflation value will be increased in case such as sudden increase of blood pressure to prevent the re-inflation.

[OFF]: Sight inflation function will turn OFF.

It will inflate to the target level set according to the previous measurement result.

NOTE

- The sight inflation function can be used only during the NIBP auto mode measurement.
- The sight inflation function cannot be used when the patient classification is "Neonate".
- The sight inflation function cannot be used when performing the 1-minute interval measurement or continuous measurement.
- When performing manual measurement/measurement at alarm occurrence, it will inflate to the set target inflation value regardless of the sight inflation setting.

11 Mean BP (MAP) Display

[ON]: Mean BP (MAP) value will be displayed.



 If the mean BP (MAP) value is not displayed, the mean BP (MAP) alarm will not be generated.

12 End Tone

[ON]: A buzzer tone will be generated when the NIBP measurement completes.

13 User Interval

The interval is fixed as "Lumbar Mode".

(Bout the Lumbar Mode" P7-56)

14 Auto Mode with Start/Stop Key

NIBP measurement will be performed automatically at selected time intervals.

- [OFF]: When the power is turned ON, NIBP auto mode will resume even after the patient is discharged regardless of whether the next patient is admitted or not.
- ▶ [ON]: When the power is turned ON, NIBP auto mode will resume by starting a manual measurement for the newly admitted patient. Until the NIBP auto mode is resumed or the interval is changed, "Standby" will be displayed inside the NIBP numeric data box.

NOTE

- If the power OFF duration was within 30 seconds, the NIBP auto mode will resume at power ON even when the above setting is [ON].
- When the [NIBP Start/Stop] key is pressed or when the NIBP measurement interval is changed, the standby mode will be canceled and the NIBP periodic measurement will start.

15 Time Display

The time for the NIBP measurement will be displayed.

- [Elapsed]: The elapsed time from the previous NIBP measurement will be displayed.
- [Meas.]: The NIBP measured time will be displayed.

16 Periodic Measurement Starting Time

The starting time of periodic measurement can be set.

- [Time]: The periodic measurement will start from the integral multiple of the selected interval starting from 0min.
- [Meas.]: The periodic measurement will start from the actual starting time.

	Measurement time when [Time] is selected:	Measurement time when [Meas.] is selected:
When the interval is	15:11:15	15:11:15
[15min.] and the	15:15:00	15:26:15
measurement is started on	15:30:00	15:41:15
15:11:15	15:45:00	15:56:15
When the interval is	16:00:00	16:26:15
changed to [30min.] on	16:30:00	16:56:15
15:58	17:00:00	17:26:15

17 Target Inflation Value

The window to set the target inflation value will be displayed.

Set the target inflation value using the up/down keys. The indication of target inflation value is SYS + 40 mmHg for adult/child, SYS + 30 mmHg for neonate.



Patient Classification	Target Inflation Value	Maximum Inflation Value	Maximum Measurement Duration
Adult	100 mmHg to 290 mmHg (Default: 180 mmHg)	300 mmHg	160 sec.
Child	100 mmHg to 200 mmHg (Default: 140 mmHg)	210 mmHg	160 sec.
Neonate	100 mmHg to 140 mmHg (Default: 110 mmHg)	150 mmHg	80 sec.

For the following case, the target inflation value will be automatically set to the default value of each patient classification.

- When the patient is discharged
- When the patient classification is changed

- When the "Sight Inflation" is [OFF], the target inflation value after the first periodic measurement will be automatically set based on the previous measurement value.
- When the "Sight Inflation" is [ON], the target inflation value after the first periodic

18 Cancel Error

By pressing [Cancel Error], the measurement error can be canceled.

NOTE

 Make sure that the NIBP measurement can be properly performed after solving the cause of the NIBP system error message. If the message still remains, equipment failure can be considered.

(@"Non-Invasive Blood Pressure" P11-37)

Temperature

This section explains the measurement procedure and measurement condition setup of temperature (T1 to T8).

TEMP Monitoring

1 Select the appropriate probe for the patient.

- Before the measurement, make sure that the specified probe/relay cable is used. If unspecified probe/relay cable is used, measurement error may occur.
- Stop using the probe if it is damaged..

NOTE

• 700 series temperature probe cannot be used.

 $\mathbf{2}$ Connect the probe to HS-8000 or module.

REFERENCE

- The HS-8000 or module utilizes multiparameter amplifier input method which allows monitoring of 2 channels of temperature through the 2ch temperature relay cable (CJO-P01T-DA**) connected to the HS-8000 or module connector. The measurement is also possible using the HM-800/HM-801 Multi Module inserted to the input box.
- 1 Connect the 2ch temperature relay cable (CJO-P01T-DA**) to the multiparameter connector on the HS-8000.

2 Connect the temperature probe to the 2ch temperature relay cable.



3 Attach the probe to the patient.

In Case of Body Surface Probe 409B:

1 Attach the probe to the body surface, and secure with surgical tape.



· The probe location shown above is an example. Adjust the probe location according to the patient's condition.

In Case of Rectal Temperature Probe 401, 402:

- 1 Clean/Disinfect/Sterilize the probe according to the guidelines provided with the probe product.
- 2 Insert the probe into the rectum about 3 cm to 7 cm deep.
- 3 Secure the probe to inner thigh with surgical tape.





4 Check that the temperature is displayed.

- 1 Press the [Home] key on the user key.
- 2 Verify that the measured data is displayed on the home display. If the measured data is not displayed during the 1 channel temperature measurement, the temperature probe may be connected to incorrect channel. Connect the probe to the correct channel and verify that the

measured data is displayed.



TEMP Parameter Setup

Channel Settings of BP, TEMP

The default settings of BP and TEMP channels are as follows.



When the DS-8007 is used, the temperature measured on temperature jack is fixed as T1, T2.

NOTE

The channel settings can be changed on the "Multi-amplifier" screen ([Initial Settings > System > Unit Module]).

(@Maintenance Manual "Unit Module Setup" P4-17)

Parameter Setup

Press the [Menu], [TEMP] keys to display the "TEMP" setup screen.





Select the label from [Tx] to [US7].



• When the parameter display is set to OFF, the alarm generation and tabular/graphic trend input will also cease.



Cardiac Output and Blood Temperature

When thermodilution catheter is used to measure the cardiac output, the blood temperature (Tb) can be monitored. The CO measurement can be performed using the multiparameter connector on the HS-8000 or module. The measurement is also possible using the HM-800/HM-801 Multi Module inserted to the input box. (Cardiac Output (CO)" P8-45)

Connecting the Super Unit

1 Select the catheter relay cable.

NOTE

· The usable catheter relay cable depends on the injectate temperature measurement method.Select the appropriate cable according to the used measurement method.

Injectate Temperature Measurement Method	Catheter Relay Cable
0°C/24°C Temperature	CJO-P01C-C2.4
Flow-through Sensor	CJO-P01C-F2.4
In-line Sensor	CJO-P01C-L2.4
Injectate Temperature Probe	CJO-P01C-T2.4

2 Connect the catheter relay cable to the multiparameter connector on the HS-8000/DS-8007 or HM-800/HM-801 Multi Module, and connect the catheter to the catheter relay cable.

Example of In-line Sensor





Cardiac Output Measurement Algorithm

Cardiac output is measured using the thermodilution method.

Thermodilution Method

The thermodilution catheter is inserted from the vein through the right atrium, right ventricle, and pulmonary artery. From the side hole near the catheter tip, injectate is injected quickly to the right atrium. At this time, the heart contraction and heat diffusion mixes the injectate with blood, and causes blood temperature fall. Variable initiated by these effects are measured as time function at the pulmonary artery, and the following thermodilution curve can be drawn.

Cardiac output is calculated by applying this to the Stewart-Hamilton formula shown below.



 $CO = 60 \cdot Vi \cdot \frac{Si \cdot Ci}{Sb \cdot Cb} \cdot \frac{Ct(Tb-Ti)}{S} = CC \cdot \frac{Tb-Ti}{S}$

S

- CO : Cardiac Output [L/min]
- Vi : Injectate Volume [L]
- Tb : Blood Temperature [°C]
- Ti : Injectate Temperature [°C]
- Ct : Correction coefficient for injectate temperature rise inside catheter
- 60 : seconds
 - : Area of thermodilution curve $\int_{0}^{\infty} \Delta Tb(t) dt[^{\circ}C sec]$

 $\Delta Tb(t)$: Temperature change of Tb after "t" seconds. [°C]

- CC : Catheter Constant (Computation Constant: CC value)
- Si : Specific Gravity of Injectate [g/cm³]
- Sb : Specific Gravity of Blood [g/cm³]
- Ci : Specific Heat of Injectate [cal/(g/°C)]
- Cb : Specific Heat of Blood [cal/(g/°C)]

As shown above, cardiac output is directly proportional to the Injectate Volume (Vi) and the difference between Blood Temperature and Injectate Temperature (Tb - Ti), and is inversely proportional to the area of the thermodilution curve (S).

Hematocrit Value

Hematocrit value of 45%, (Si*Ci)/(Sb*Cb) = 1.08 is programmed for this equipment.

(NOTE

• If the hematocrit value is different, an error may be caused in cardiac output measurement.

Blood Temperature Alarm Setup

Press the [TEMP], [Tb] keys.
(Press To Display the Parameter Setup Screen P7-1)

> The alarm setup menu will be displayed.

NOTE



2 Select ON/OFF of blood temperature alarm and set the upper and lower alarm limits. (@ "Alarm Limit Setup for Each Parameter" P6-10)

• Set the upper limit in the range of 31.0°C to 45.0°C. If a value above 45.0°C is set, the

upper alarm will turn OFF.

 Set the lower limit in the range of 30.0°C to 44.0°C. If a value below 30.0°C is set, the lower alarm will turn OFF.

REFERENCE

- The upper and lower limit can be set in 0.5°C increments.
- When [Auto] is set, the upper and lower limit will be automatically set to +2.0 °C and -2 °C to the current value respectively.
- [Auto] key will be displayed only when [Enable] is set for "Auto Alarm Setup" under "Initial Settings".

CO₂ Concentration (Mainstream Method)

This section explains about the CO_2 concentration measurement procedure and measurement condition setup when using the Philips Capnostat 5 (Mainstream Method, Gas Unit I/F HPD-800/HPD-810/HPD-820).

 When the Multigas Unit (MGU-800/MGU-810 series) and HPD-800/HPD-810/HPD-820 are simultaneously used, the CO₂ concentration measurement will be performed by the equipment selected for the "CO₂ Source Priority" under ([Menu] > "Parameter" [CO₂]).

Patient Application and Display

By using the HPD-800/HPD-810/HPD-820 CO₂ Gas Unit I/F, CO₂ measurement by the Philips Capnostat 5 (Mainstream Method) can be performed.

Connect the HPD-800/HPD-810/HPD-820 Gas Unit I/F to the connector.

The connection procedure differs depending on the gas unit type.

In Case of HPD-800/HPD-810

Connect to the AUX connector of the HS-8000, DS-8007 or HM-801.

In Case of HPD-820

Connect to the CO₂ I/F connector on the DS-8007.

 \mathbf{Z} Connect the CO₂ sensor (Capnostat5) to the CO₂ connector on the HPD-800/HPD-810/HPD-820.



- The CO₂ sensor will automatically begin warming up. The CO₂ sensor requires a warming up process to achieve stable operating temperature. Warm up process will require minimum of 2 minutes.
- ▶ During the warm up period, <CO₂ Warm Up> message will be displayed on the monitor.
- > When the warm up completes, the message will disappear.

NOTE

•When using more than one HM-801, only one AUX connector on the HM-801 can be used for measurement of CO₂, BIS (with HBX-800).

•When the CO₂ is measured on both the HS-8000/DS-8007 and HM-801, the measurement of the HS-8000/DS-8007 will be prioritized.

3 Prepare an airway adapter suitable for the patient.

•The disposable airway adapter should be opened just before use.

•Do not reuse the disposable airway adapter.Do not disassemble, clean, disinfect, or sterilize it.

NOTE

•There are 4 types of airway adapters.Select the appropriate adapter according to the used endo-tracheal tube size and operating environment.

<u>Airway Adapter (Adult) 7007</u> For patients using an endo-tracheal tube more than, or equal to 4.0 mm in diameter. Reusable Type
<u>Airway Adapter (Neonate) 7053</u> For patients using an endo-tracheal tube less than, or equal to 4.0 mm in diameter. Reusable Type
<u>Airway Adapter (Disposable, Adult) 6063</u> For patients using an endo-tracheal tube more than, or equal to 4.0 mm in diameter. Single-Use Type
<u>Airway Adapter (Disposable, Neonate) 6312</u> For patients using an endo-tracheal tube less than, or equal to 4.0 mm in diameter. Single-Use Type

4 Verify that the warm up is complete, and attach the CO₂ sensor to the airway adapter until a click sound is heard.

- 1 Capnostat 5 CO₂ Sensor
- 2 Window
- 3 Airway Adapter
- A: Thick Side
- B: Thin Side



\land CAUTION

•The airway adapter should be attached with the thicker side facing to the patient. If attached oppositely, it may damage the CO₂ sensor or airway adapter.

5 Perform the setting for the O₂ compensation, N₂O compensation, anesthetic gas compensation, atmospheric pressure Set these items each time the condition changes. (GP "CO2 Parameter Setup" P7-74)

O Press the [Menu], [CO₂] ("Parameter"), [Calibrate Airway Adapter] keys to calibrate the airway adapter.

- Calibration will start.
- > During calibration, <Zeroing> will be displayed.
- ▶ Upon completion of calibration, a tone will be generated and <Cal. complete> will be displayed.
- ▶ If the calibration fails, an error tone will be generated and <Cal. error> will be displayed.

NOTE

• The airway adapter calibration must be performed before connecting to the respiration circuit.

The airway adapter calibration should be also performed for the following case.

- When the airway adapter is replaced.
- When <Zero the CO₂ Adapter> or <Check airway adapter.> is displayed.
- A clean airway adapter must be used. If reusing an airway adapter, clean and air-dry it. Then, wipe the window with a swab, and sterilize (EOG, etc.) before use.
- During the calibration, the measurement data will be displayed as "---". The measurement data during calibration may be included in the trend data causing discontinuity.
- Calibration cannot be performed if respiration is detected within 20 seconds before calibration. In such case, wait for 20 seconds and perform the calibration again.
- · When <Cal. error> is displayed, perform the airway adapter calibration again.
- The EtCO₂ accuracy is checked according to the test method of ISO 80601-2-55: 2011 (Medical electrical equipment-Part 2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors).
- The waveform sampling rate is 100 Hz.
- Quantitative effects of humidity and condensation: Full accuracy specifications will be maintained for all non-condensing humidity levels.
- The CO₂ measurement accuracy is tested at 35°C.
- The respiration rate accuracy was verified by using a solenoid test setup to deliver a square wave of known CO₂ concentration to the device. 5% and 10% CO₂ concentrations were used and respiration rate was varied over the range of the device. Pass/Fail criteria was comparison of the respiratory rate output from the sensor to the frequency of the square wave. EtCO₂ measurements at those rates were compared to the CO₂ readings under static flow conditions.

Verify that the airway adapter calibration is properly completed, disconnect the CO₂ sensor from the airway adapter temporarily, and attach the airway adapter to the patient's respiration circuit.

 $\boldsymbol{\delta}$ Connect the CO₂ sensor to the airway adapter.



 $\mathbf{9}$ Verify that the CO₂ waveform, EtCO₂ value, InspCO₂ value are displayed.



CO₂ Parameter Setup

Press the [Menu], [CO₂] keys to display the "CO₂" setup screen.



Calibrate Airway Adapter

The airway adapter will be calibrated. (@ "Patient Application and Display" P7-71)

2_{Scale}

Select from [0-50]/[0-100] if the measurement unit is mmHg, and from [0-4]/[0-8]/[0-10] if the unit is kPa or %.

$\mathbf{3}_{\text{EtCO}_2}$ (End-tidal CO₂)

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The EtCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit in the range of 3 mmHg to 100 mmHg/0.3 kPa to 13.3 kPa/0.3% to 13.3%.

Setting a value above 100 mmHg/3.3 kPa/13.3% will turn OFF the alarm.

• Set the lower limit in the range of 1 mmHg to 98 mmHg/0.1 kPa to 13.1 kPa/0.1% to 13.1%.

Setting a value below 1 mmHg/0.1 kPa/0.1% will turn OFF the alarm.

• When Capnostat 5 is used, EtCO₂ alarm will not generate unless 2 or more respiration is detected after completion of the airway adapter calibration.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper/lower limit can be set in 1 mmHg/0.1 kPa/0.1% increment.
- When [Auto] is set, the upper and lower limit will be automatically set to +10 mmHg / +1.3 kPa / +1.3%, and -10 mmHg / -1.3 kPa / -1.3% respectively to the current value.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

4 InspCO₂ (Inspired CO₂)

(@"Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The InspCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit in the range of 1 mmHg to 4 mmHg/0.1 kPa to 0.4 kPa/0.1% to 0.4%.

Setting a value equal to or above 4 mmHg/0.4 kPa/0.4% will turn OFF the alarm.

• When Capnostat 5 is used, InspCO₂ alarm will not generate unless 2 or more respiration is detected after completion of the airway adapter calibration.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper limit can be set in 1 mmHg/0.1 kPa/0.1% increments. There is no lower limit.
- When [Auto] is set, the upper limit will be set to 3 mmHg / 0.4 kPa / 0.4%.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Alarm Related Setup" P5-4)

• When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

$\mathbf{5}_{EtCO_2}$ Peak Duration

[10sec]/[20sec]: Maximum EtCO₂ value for the selected duration will be displayed.

[OFF]: EtCO₂ value for each respiration will be displayed.

NOTE

- As the EtCO₂ value display is updated each second, EtCO₂ value for each respiration cannot be displayed if respiration rate is 60 Bpm and above.
- For the InspCO₂ value, minimum value of 20 seconds will be displayed regardless of the setting.

$\mathbf{6}_{O_2}$ Compensation

By entering the used O_2 concentration value, compensation can be made to display more accurate value. Enter the O_2 compensation value on the " O_2 " screen, and press the [Set] key.

O ₂ Comp.			(\mathbf{X})
×02 G	DTP. + ANI	esthet ic	Comp. ≤ 100%
7	8	9	(0 - 100%)
4	5	6]
1	2	3	Input
0		С	Cancel

NOTE

• The value cannot be changed if the total value of O₂ compensation and anesthetic agent compensation exceeds 100%. In such case, change the O₂ compensation value after changing the anesthetic agent compensation value.

 7_{N_2O} Compensation

NOTE

• If N₂O is present in the respiration circuit, the CO₂ value tends to be displayed higher than the actual value. By setting the N₂O compensation to [ON], this can be adjusted.

8 Anesthetic Agent Compensation

By entering the used anesthetic agent concentration value, compensation can be made to display more accurate value.

Enter the anesthetic compensation value on the "Agent" screen, and press the [Set] key.

Anesthetic Comp. 🗙				
		Π	Π	\prod
*	0 2 Co	np. + An	esthet ic	Conp. \leq 100%
	7	8	9	(0.0 - 20.03)
4	4	5	6	
[1	2	3	Input
	0		С	Cancel

NOTE

• The value cannot be changed if the total value of O₂ compensation and anesthetic

agent compensation exceeds 100%. In such case, change the anesthetic agent compensation value after changing the O_2 compensation value.

9 Atmospheric Pressure

By entering the atmospheric pressure, the pressure difference will be compensated and allows more accurate measurement.

Enter the atmospheric pressure value on the "Atmos. Pressure" screen, and press the [Set] key.



10_{CO2} Source Priority

When MGU-800/MGU-810 and HS-8000/DS-8007/HM-801 are simultaneously used, the CO₂ source to prioritize the measurement can be set.

- ▶ [MGU-800]: CO₂ value measured by the MGU-800/810 Multigas Unit will be prioritized.
- ▶ [HS-8000]: CO₂ value measured by the HS-8000/DS-8007/HM-801 will be prioritized.

NOTE

 When the HS-8000/DS-8007 and HM-801 are simultaneously used, the CO₂ measurement of the HS-8000/DS-8007 will be prioritized.

11 Display ON/OFF

(@"ECG Parameter Setup" P7-6)

- When the waveform and numeric data display is set to OFF, the alarm generation and tabular trend input will also cease.
- When the waveform and numeric data display is set to OFF, the respiration rate measured by CO₂ will not be displayed either.

REFERENCE

• During "Display OFF" condition, the display will automatically return if the filter line is applied to the patient and more than 2 respirations are detected in 30 seconds.

CO2 Concentration (Sidestream Method)

The HCP-800/HCP-810/HCP-820 is a CO_2 Gas Unit which measures CO_2 concentration. The HCP-800/HCP-810/ HCP-820 CO_2 Gas Unit incorporates Microstream technology of Covidien for EtCO₂ (End-tidal CO₂ concentration) and InspCO₂ (Inspiratory CO₂ concentration) measurement. This section explains about the procedure and setup of the CO₂ concentration measurement of the HCP-800/HCP-810/HCP-820.

WARNING

- When using a sampling line for intubated patients with a closed suction system, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.
- Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.
- Do not cut or remove any part of the sampling line. It could lead to erroneous readings.
- If too much moisture enters the sampling line (i.e., from ambient humidity or breathing of unusually humid air), <Check Sample Line> will appear in the message area. Replace the sampling line when this message appears.
- Carefully route the filter line to reduce the possibility of patient entanglement or strangulation.
- Do not lift the HCP-800/HCP-810/HCP-820 by the filter line, as the filter line could disconnect from the equipment, causing the equipment to fall on the patient.
- CO₂ readings and respiration rate can be affected by sensor application, ambient environment, and patient conditions.

- When the multigas unit (MGU-800/MGU-810 series) and HCP-800/HCP-810/HCP-820 are simultaneously used, the CO₂ concentration measurement will be performed by the equipment selected for the "CO₂ Source Priority" under ([Menu] > "Parameter" [CO₂]).
- The Microstream EtCO₂ sampling lines are designed for single patient use, and are not to be reused. Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line as this can cause damage to the monitor.
- Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.
- Before use, carefully read the Directions for Use for the Microstream EtCO₂ sampling line.
- Use only the Microstream EtCO₂ sampling line to ensure proper function of the monitor.
- During nebulization or suction for intubated patient, remove the sampling line from the HCP-800/HCP-810/HCP-820 to avoid moisture buildup and sampling line occlusion.
- Replace the sampling line according to hospital protocol or when a blockage is indicated on the equipment. Excessive patient secretions or a buildup of liquids in the airway tube may occlude the sampling line, requiring more frequent replacement.

NOTE

- When connecting a sampling line to the HCP-800/HCP-810/HCP-820, screw the sampling line clockwise into the connector firmly to avoid inaccurate measurement which may be caused by gas leak from the connection point.
- When <Check Sample Line> appears on the screen indicating that the filter line connected to the HCP-800/HCP-810/HCP-820 is blocked, the CO₂ pump will stop pumping the patient's breath to the monitor. In such case, follow the instructions in the "Troubleshooting"

section of this manual. First, disconnect and reconnect the filter line. If the message still appears, disconnect and replace the filter line. Once a working filter line is attached, the pump will automatically resume operation.

- After connecting the CO₂ sampling line to the HCP-800/HCP-810/HCP-820 and patient, check that CO₂ values appear on the monitor display.
- The EtCO₂ accuracy is checked according to the test method of ISO 80601-2-55: 2011 (Medical electrical equipment - Part 2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors).
- The waveform sampling rate is 20 samples per second.
- The minimum value and maximum value of the CO₂ waveform are used for the InspCO₂ value and EtCO₂ value respectively.
- When using the HCP module with a ventilator, under high over pressures close to 10 kPa (100 cmH₂O), the module may enter into a blockage mode in order to protect the module from damage.
- The respiration rate test simulates breaths for use in respiration rate measurement with a system which uses a tank of N2 (representing no CO₂ for inhalation) and a tank of CO₂ (of the %CO₂ required for the particular test). A control board, which is triggered by a computer, uses solenoids to switch the module input between the 2 tanks of gas, creating a gas CO₂ square wave. This system can create simulated breaths over the full required range of specified respiration rates.

Patient Application and Display

The CO₂ concentration can be measured by using the HCP-800/HCP-810/HCP-820 CO₂ Gas Unit.

- NOTE
- Accurate CO₂ concentration measurement can be acquired after 40 seconds from turning the power ON.

Connect the HCP-800/HCP-810/HCP-820 CO₂ Gas Unit to the connector.

The connection procedure differs depending on the gas unit type.

In Case of HCP-800/HCP-810

Connect to the AUX connector of the HS-8000, DS-8007, or HM-801.

In Case of HCP-820

Connect to the CO₂ I/F connector on the DS-8007.

NOTE

- When using more than one HM-801, only one AUX connector on the HM-801 can be used for measurement of CO₂, BIS (with HBX-800).
- When the CO2 is measured on both the HS-8000/DS-8007 and HM-801, the measurement of the HS-8000/DS-8007 will be prioritized.

f 2 Attach the airway adapter, oral/nasal sampling line or nasal sampling line to the patient.

For intubated patient



- 1 Attach the airway adapter to respiration circuit.
- 2 Connect one end of the sampling line to the connector on the HCP-800/HCP-810/HCP-820. Verify that all the tubes are properly connected.

For patient using the nasal prong



- 1 Attach the nasal or oral/nasal patient interface of the sampling line to the patient as described in the sampling line directions for use.
- 2 Connect the sampling line to the connector on the HCP-800/HCP-810/HCP-820. Verify that all the tubes are properly connected.

3 Start the CO₂ concentration measurement.



▶ Verify that the CO₂ waveform, EtCO₂ value, InspCO₂ value are displayed.

- If the power supply is interrupted due to power failure, etc., HCP-800/HCP-810/HCP-820 will be initialized even if the power interruption was within 30 seconds.
- NOTE
- Connecting a sampling line or nasal prong to the HCP-800/HCP-810/HCP-820 will automatically start the sampling pump. To prevent the pump from deteriorating, disconnect the sampling line and nasal prong from the HCP-800/HCP-810/HCP-820

when not measuring the CO_2 concentration.

- · Set the scale, measurement unit, alarm, etc. as necessary.
- When ambient temperature or atmospheric pressure changes significantly, auto zeroing will function. During auto zeroing, "---" will be displayed inside the CO₂ numeric data box and CO₂ measurement cannot be performed.

CO₂ Parameter Setup

Press the [Menu], [CO₂] keys to display the "CO₂" setup screen.



1 Scale

Select from [0-50]/[0-100] if the measurement unit is mmHg, and from [0-4]/[0-8]/[0-10] if the unit is kPa or %.

 $\mathbf{2}_{\text{EtCO}_2}$ (End-tidal Carbon Dioxide)

(P"Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The EtCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit in the range of 3 mmHg to 100 mmHg/0.3 kPa to 13.3 kPa/0.3% to 13.3%.

Setting a value above 100 mmHg/3.3 kPa/13.3% will turn OFF the alarm.

• Set the lower limit in the range of 1 mmHg to 98 mmHg/0.1 kPa to 13.1 kPa/0.1% to 13.1%.

Setting a value below 1 mmHg/0.1 kPa/0.1% will turn OFF the alarm.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper/lower limit can be set in 1 mmHg/0.1 kPa/0.1% increment.
- When [Auto] is set, the upper and lower limit will be automatically set to +10 mmHg / +1.3 kPa / +1.3%, and -10 mmHg / -1.3 kPa / -1.3% respectively to the current value.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

3 InspCO₂ (Inspired Carbon Dioxide)

(@"Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- The InspCO₂ alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, CO₂ unit is connected, or a patient is discharged.
- Set the upper limit in the range of 1 mmHg to 4 mmHg/0.1 kPa to 0.4 kPa/0.1% to 0.4%.

Setting a value equal to or above 4 mmHg/0.4 kPa/0.4% will turn OFF the alarm.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper limit can be set in 1 mmHg/0.1 kPa/0.1% increments. There is no lower limit.
- When [Auto] is set, the upper limit will be set to 3 mmHg / 0.4 kPa / 0.4%.
- To use the auto setting, [Enable] should be selected for "Alarm Auto Setup" under "Initial Settings". (Alarm Related Setup" P5-4)
- When [Enable] is set for "Alarm Threshold Limit" under "Initial Settings", and auto setting is used, the upper and lower alarm limit will be clipped to the set limit.

4 EtCO₂ Peak Duration

[10sec]/[20sec]: Maximum EtCO₂ value, minimum InspCO₂ value for the selected duration will be displayed. [OFF]: EtCO₂ value, InspCO₂ value for each respiration will be displayed.

NOTE –

 As the EtCO₂ value display is updated each second, EtCO₂ value for each respiration cannot be displayed if respiration rate is 60 Bpm and above.

$\mathbf{5}_{\mathrm{CO}_2}$ Source Priority

When MGU-800/MGU-810 and HS-8000/DS-8007/HM-801 are simultaneously used, the CO_2 source to prioritize the measurement can be set.

- ▶ [MGU-800]: CO₂ value measured by the MGU-800/810 Multigas Unit will be prioritized.
- ▶ [HS-8000]: CO₂ value measured by the HS-8000/DS-8007/HM-801 will be prioritized.

NOTE

 When the HS-8000/DS-8007 and HM-801 are simultaneously used, the CO₂ measurement of the HS-8000/DS-8007 will be prioritized.

$\mathbf{6}_{CO_2}$ Calibration

CO₂ calibration can be performed.

(Maintenance Manual "CO2 Calibration (HCP-800/HCP-810/HCP-820)" P9-11)

Suspend CO₂

[Suspend]: The pump operation will stop, CO_2 waveform and numeric data display will disappear, and "Suspended" will be displayed inside the CO_2 numeric data box.

[Resume]: Resumes CO₂ monitoring. This key will be displayed when the measurement is suspended.

• When the measurement is suspended, the alarm generation and trend input will be also suspended.

B Display ON/OFF

(@"ECG Parameter Setup" P7-6)

- When the waveform and numeric data display is set to OFF, the alarm generation and tabular trend input will also cease.
- When the waveform and numeric data display is set to OFF, the respiration rate measured by CO₂ will not be displayed either.

REFERENCE

• During "Display OFF" condition, the display will automatically return if the filter line is applied to the patient and more than 2 respirations are detected in 30 seconds.

Multigas Unit/SPIRO

The MGU-800/810 series Multigas Unit can be connected to the DS-8400 system via U-LINK connector. (@Maintenance Manual "Connection of Multigas Unit" P1-11)

When the multigas unit is connected, monitoring conditions for CO_2 concentration, anesthetic gas concentration, O_2 concentration, and N_2O concentration, respiration (SPIRO) can be set.

The MGU-800/810 series have an internal barometer and thermistor that allow compensation for changes over a range of temperature and atmospheric pressures.

WARNING

- Make sure to use only the specified Mindray Medical Sweden AB product.
 (P"Anesthetic Gas Concentration Measurement (Manufactured by Mindray Medical Sweden AB)" P13-7)
- Be careful not to damage the water trap during operation as bacteria and/or mucus may contaminate the MGU-800/810 series.
- The airway adapter, sampling line, flow sensor are disposable products that are intended for single patient use only. Do not reuse them on other patients as it may cause cross-infection.
- Do not use the MGU-800/810 series with the flammable anesthetic agents.
- To protect the hospital staffs from unnecessary anesthetic agent, it is strongly
 recommended to connect the exhaust hole to the gas exhaust system in the hospital.
- · The sampling line may get clogged by internal condensation.

- When the multigas unit (MGU-800/MGU-810 series) and HPD-800/HPD-810/HPD-820, HCP-800/HCP-810/HCP-820 are simultaneously used, the CO₂ concentration measurement will be performed by the equipment selected for the "CO₂ Source Priority" under ([Menu] > "Parameter" [CO₂]).
- The MGU-800/MGU-810 series require warm up of about 10 minutes to correctly measure the data.

- If the power supply is interrupted due to power failure, etc., MGU-800/810 series multigas unit will initialize and enter into warm-up mode even if the power interruption is within 30 seconds.
- · Zero Calibration:

The zero calibration will automatically start when the MGU-800/810 series multigas unit is connected.

After the warm-up completes, zero calibration will be performed every 4 hours during stable operation.

During warm-up, zero calibration interval will become shorter than during normal operation. During zero calibration, measurement data will not be updated. Calibration gas is not required during zero calibration.

- Make sure the sampling line and flow sensor is securely connected to prevent any leakage.
- An environment with alcoholic vapor may adversely affect the measurement readings.
- CO₂, N₂O or anesthetic agent in the atmosphere around the MGU-800/810 series may adversely affect the measurement readings.
- SPIRO and ventilator cannot be used simultaneously.
- NOTE
- The MGU-800/810 series uses a fixed correction of 11hPa (22°C@40% RH) to compensate for the influence of water vapor in the gas sample, when converting the gas readings to ATPD. An increase in the ambient H₂0 partial pressure to 30 hPa (28°C@80% RH or 33°C@60% RH) will cause a general error for all gases of only -2% REL.
- The EtCO₂ accuracy is checked according to the test method of ISO 80601-2-55: 2011 (Medical electrical equipment-Part 2-55: Particular requirements for the basic safety and essential performance of respiratory gas monitors) under the condition of 60 Bpm and below with I:E ratio of 1:1.
- When the RR exceeds 60 Bpm, the EtCO₂ accuracy cannot be specified. (Depends on the I:E ratio.)
- The data sampling rate is 25 Hz.
- The minimum value and maximum value of the CO₂ waveform are used for the InspCO₂ value and EtCO₂ value respectively.
- For the gas measurement data, "0" will be displayed if the value becomes below the following threshold for 3 seconds or more. (Full Accuracy/during warm-up) CO₂: 0.1/0.3[vol%]
 N₂O: 3/3[vol%]
 O₂: 0/0[vol%]
 Volatile Anesthetic: 0.15/0.3[vol%]
- The ability to properly resolve end-tidal values can be measured by using the set-up described in ISO 80601-2-55: 2011, figure 201. 101.
 In short, the method consists of sampling gas from two different sources connected to an electrically controlled pneumatic valve to permit rapid switching between the two sources. During the test, the valve is set to switch gas source at a number of frequencies (simulating the range of specified breath rates) and for each frequency, the end-tidal value presented by the gas analyzer is noted. From a diagram of end-tidal value over frequency, the frequency at which the gas analyzer is no longer able to resolve end-tidal values according to specification is identified.

Connecting to the Respiration Circuit

□ Multigas Concentration Measurement (MGU-800 Series)

WARNING

- Do not use adult/pediatric type water traps and/or sampling lines with neonates to avoid high sampling flow.
- Connect only DRYLINE gas sampling lines to the water trap. Note that there may be other compatible tubing present, e.g. IV-lines.
- Do not use DRYLINE neonatal sampling lines (blue luer lock nuts) with DRYLINE adult water traps as this could result in incorrect measurement data.
- Do not use DRYLINE adult sampling lines (colorless luer lock nuts) with DRYLINE neonatal water traps as this could result in incorrect measurement data.

NOTE

If [Adult] or [Child] is selected as patient classification on the "Admit/Discharge" screen, install the DRYLINE Adult/Child Water Trap (60-13100-0).
 If [Neonate] is selected as patient classification on the "Admit/Discharge" screen, install the DRYLINE Neonatal Water Trap (60-13200-0).
 If the used water trap and the set patient classification does not match, <GAS Check Water Trap Class> will be displayed.

Install the DRYLINE Water Trap (Adult/Child: 60-13100-00, Neonate: 60-13200-00) aligning the lugs with the corresponding holes in the receptacle and pushing gently into place. (See below.)

Make sure that both barbs on the lugs are fully engaged by pulling the water trap, which should be firmly seated.



Connect the DRYLINE Airway Adapter (Straight: 60-14100-00, or Elbow: 60-14200-00) to the patient breathing system.

3 Remove the protective cap from the airway adapter and connect it to the sampling line (for adult/child: 60-15200-00, for neonate: 60-15300-00).

4 Connect the other side of the sampling line to the inhale port of the water trap. When the water trap is half full, empty the water trap's reservoir.

(@Maintenance Manual "Water Trap (Multigas Unit)" P8-5)

WARNING

• The contents of the water trap should be handled as a potential infection hazard.

□ Multigas Concentration/Spirometry Measurement (MGU-810 Series)

WARNING

• Only combine the SPIRIT Flow Sensors and DRYLINE Water Traps as described in the table below. Other combinations might lead to incorrect measurements.

Patient Category	Patient Classification Selection on "Admit/ Discharge" Screen	SPIRIT Flow Sensor	DRYLINE Water Tap
Adult	Adult	Adult (60-16100-00)	Adult/Child (60-13100-00)
Child	Child	Child (60-16200-00)	Neonate (60-13200-00)
Neonate	Neonate	Child (60-16200-00)	Neonate (60-13200-00)

- · Use the adult flow sensor for a patient whose tidal volume is above 150 mL.
- · Use the pediatric flow sensor for a patient whose tidal volume is below 300 mL.
- Make sure to use the correct flow sensor depending on the patient conditions, adult or pediatric and the tidal volume.
- · Do not confuse the gas sampling line with other compatible tubing, e.g. IV-lines.

1 Install the DRYLINE Water Trap.

- Connect the end (for adult: 22/15 mm, pediatric: 15 mm) of the flow sensor, marked **+ n** to the patient tracheal tube or similar.
- **3** Connect the end of the flow sensor to the patient breathing system. For best results, a heat and moisture exchanger (HME) or similar should be put between the flow sensor and the breathing system.

4 Connect the pressure line of the flow sensor to the flow sensor connector on the MGU-810.

5 Connect the gas sampling line of the flow sensor (for adult: colorless, for pediatric: blue) to the gas inlet of the water trap. When the water trap is half full, empty the water trap's reservoir.
(Adultion of the flow sensor (Multigas Unit) P8-5)

WARNING

• The contents of the water trap should be handled as a potential infection hazard.

To prevent accumulation of condensed fluid, the flow sensor shall be always be positioned a few degrees off the horizontal level towards the ventilator side. For the same reason, the pressure tubes shall exit the flow sensor upwards.

The pressure tubes should be routed in such a way that a water lock is formed by a section of tubing being positioned lower than the flow sensor connector on the MGU-810.

8 A patient breathing system leakage test shall be performed according to the recommendations of the ventilator manufacturer.

/1\ CAUTION

The adult flow sensor dead space is 6.9 mL and the flow resistance is 1.8 cmH₂O at 60 L/ min.

The pediatric flow sensor dead space is 0.75 mL and the flow resistance is 0.9 cmH₂O at 10 L/min.

Adjust ventilation accordingly.

- To prevent condensation, the patient breathing circuit, flow sensor and pressure tubing should not be directly exposed to cooling equipment such as fans or cooling blankets.
- ٠ Leakage of gas from the patient breathing system may occur if the pressure or gas sampling lines are not connected to the MGU-810.
- The pressure tube and gas sampling lines of the flow sensor should always be routed from the patient circuit to the MGU-810 such a way as to avoid kinking.
- Flow sensors that have suffered damage to sensor head, tubing or tubing connector must not be used.
- If liquid has entered the pressure tubes, it can be removed by gently tapping or shaking the flow sensor.

CO₂ Measurement Unit Setup



- alarm limit, measurement unit and scale are common for all the units and modules.
- When a measurement unit is changed, make sure to set the alarm condition for that unit. Set the alarm for each measurement unit.

Press the [Menu], [Initial Settings], [Meas.], [Unit] keys.

> The "Unit" menu will be displayed.

Menu 🗲 Initial Set	tings > Meas.
User Label	Unit Other
Explanation	Area (F)
C02	nalig
BP	nniig
CVP	mHg/kPa
TEMP	ť
ST	an
Height/Weight	Cn/kg

Press the [mmHg]/[kPa]/[%] key.

> The data of currently set measurement unit will be displayed on the graphic/tabular trend.



GAS Display during Undetected Breath

The gas data display when a respiration is not detected can be selected from [None] (bar display) or [Insp. Only] (displays only the inspiratory data).

1 Press the [Menu], [Initial Settings], [Meas.], [Other] keys.

> The "Other" menu will be displayed.

Menu > Initial Settings > Meas.	0
User Label Unit Other	\sum
Explanation Area	2
MIBP Start OFF HR/PR Source ECD/Sp02/BP Smin. early Gas Display Gas Display	
Under extend None	
WAP Calc. (ART. NIBP) Have Haudfacturer BIOSENS	
Arrhythmia Analysis Filter Disp Haveform ARGON	
EDWARDS	
Synchronized BarX/Tome ECG ECG	

2 Press the [None]/[Insp. Only] key.

[None]: When a respiration is not detected, inspiratory and expiratory data will become invalid and bar marks will be displayed instead.

[Insp. Only] : When a respiration is not detected, only the inspiratory data will become valid and bar marks will be displayed for expiratory data.

NOTE

 When [Insp. Only] is selected for "GAS Display during Undetected Breath" and if only inspiratory data is displayed, inspiratory and expiratory data display on the central monitor will become invalid.

• When [Insp. Only] is selected for "GAS Display during Undetected Breath" and if only inspiratory data is displayed, the GAS alarm will not be generated.

Multigas Unit Data Setup (Multigas Concentration/Spirometry)

Press the [Menu], [GAS] "Parameter" keys.

▶ The Multigas setup screen will be displayed.







GAS_MAC Screen (MGU-800 series)



GAS_SPIRO Screen (MGU-810 series)



 $\mathbf{2}$ Set the CO₂ waveform scale.

Select from [0-50]/[0-100] if the measurement unit is mmHg, and from [0-4]/[0-8]/[0-10] if the unit is kPa or %.



GAS_AGT Screen (MGU-800 series)



GAS_RESP Screen (MGU-800 series)



 $\mathbf{3}$ Set the O₂ waveform scale.

Select from [18-30]/[18-60]/[18-100]/[0-30]/[0-60]/[0-100].

4 Set the alarm.

(P"Alarm Limit Setup for Each Parameter" P6-10)

NOTE

 The following alarm will not generate unless 2 or more respiration is detected within 30 seconds after the power is turned ON, monitoring is resumed, multigas unit is connected, or a patient is discharged.

EtCO₂ Alarm



- Set the upper limit in the range of 3 mmHg to 100 mmHg, 0.3 kPa to 13.3 kPa, 0.3% to 13.3%. Setting a value above 100 mmHg, 3.3 kPa, 13.3% will turn OFF the alarm.
- Set the lower limit in the range of 1 mmHg to 98 mmHg, 0.1 kPa to 13.1 kPa, 0.1% to 13.1%. Setting a value below 1 mmHg, 0.1 kPa, 0.1% will turn OFF the alarm.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper/lower limit can be set in 1 mmHg, 0.1 kPa, 0.1% increment.

InspCO₂ Alarm

NOTE

• Set the upper limit in the range of 1 mmHg to 4 mmHg, 0.1 kPa to 0.4 kPa, 0.1% to 0.4%. Setting a value above 4 mmHg / 0.4 kPa / 0.4% will turn OFF the alarm.

REFERENCE

- The alarm limit should be set for each unit (mmHg/kPa/%).
- The upper/lower limit can be set in 1 mmHg, 0.1 kPa, 0.1% increment.

ExpO₂ Alarm

NOTE

- Set the upper limit in the range of 18% to 100%. The alarm will turn OFF if a value above 100% is set.
- Set the lower limit in the range of 18% to 100%. The alarm will turn OFF if a value below 18% is set.

REFERENCE

• The upper/lower limit can be set in 2% increment.

$InspO_2$ Alarm

NOTE

• Set the upper limit in the range of 18% to 100%.
The alarm will turn OFF if a value above 100% is set.

• Set the lower limit in the range of 18% to 100%. The alarm will turn OFF if a value below 18% is set.

REFERENCE

• The upper/lower limit can be set in 2% increment.

ExpN₂O/ InspN₂O Alarm

 NOTE
 Set the upper and lower limit in between 0 to 100%. The upper limit and lower limit will turn OFF if a value above 100% and below 0% is set respectively.

REFERENCE

• The upper/lower limit can be set in 2% increment.

AGT-E/AGT-I Alarm (MGU-810 series)

\frown		NOTE
		NOTE
	٠	The adjustable range of the upper limit differs depending on the anesthetic gas label. ISO, HAL, ENF: 0.5% to 6.0%
		SEV: 0.5 to 8.0%
		DES: 0.5% to 18.0%
		The alarm will turn OFF if a value above the range is set.
	•	The adjustable range of the lower limit differs depending on the anesthetic gas label.
		ISO, HAL, ENF: 0.5% to 6.0%
		SEV: 0.5% to 8.0%
		DES: 0.5% to 18.0%
		The alarm will turn OFF if a value below the range is set.

REFERENCE

• The upper/lower limit can be set in 0.5% increment.

MAC Alarm

NOTE
Set the upper limit in the range of 0.1 to 9.9. The upper limit alarm will turn OFF if a value below 9.9 is set.

REFERENCE

• The upper limit can be set in 0.1 increments.

RR/Apnea Alarm

NOTE

• Set the upper limit of RR alarm in the range of 10 Bpm to 60 Bpm. If a value above 150 Bpm is set, the upper alarm will turn OFF.

Set the upper limit of apnea alarm in the range of 10 sec. to 60 sec. If a value above 60 sec. is set, the upper alarm will turn OFF.

 Set the lower limit of RR alarm in the range of 5 Bpm to 145 Bpm. If a value below 5 Bpm is set, the lower alarm will turn OFF.

REFERENCE

 The adjustable increment for RR alarm depends on the patient classification and "RR Alarm Increment" setting. (Initial Settings>User I/F).

	Alarm Increment (Initial Settings > User I/F)		
	Normal	Small	
Adult	5 Bpm increment	1 Bpm increment	
Child/Neonate	2 Bpm increment	1 Bpm increment	

· The apnea alarm can be set in 1 second increment.

ExpMV/PEAK/PEEP Alarm

NOTE

- Set the upper/lower limit of ExpMV alarm in the range of 2.0 L/minute to 20 L/minute for Adult, 0.5 L/minute to 5.0 L/minute for Child/Neonate.
- Set the upper/lower limit of PEAK alarm in the range of 8 cmH₂O to 100 cmH₂O.
- Set the upper/lower limit of PEEP alarm in the range of 2 cmH₂O to 50 cmH₂O.

REFERENCE

· The upper/lower limit can be set as followings. ExpMV alarm can be set in 0.5 L/minute increment. PEAK/PEEP alarm can be set in 1 cmH₂O increment.

5 Perform the zero calibration.Press the [Zero Cal.] key to start the zero calibration.

NOTE

· While performing the zero calibration, the baseline waveform is displayed.

REFERENCE

· On the patient monitor, a zeroing (zero calibration) of the multigas unit is periodically performed, but it can also be performed manually when necessary.

6 Set the "Flow Rate" (sampling flow rate for the multigas unit).

The selectable "Flow Rate" value differs depending on the type of used water trap (adult/child or neonate) and sampling line.

- ▶ When using a water trap for adult/child, select from [120]/[150]/[200].
- ▶ When using a water trap for neonate, select from [70]/[100]/[120].

NOTE

- · If the used water trap and the set patient classification does not match, <GAS Check Water Trap Class> will be displayed.
- If <GAS Pump Regulating> is displayed, the gas sampling flow rate may be insufficient. Check the sample line for any blockage or bent. If the message is still displayed, adjust the flow rate.
- ٠ Select the appropriate water trap, sampling line, or flow sensor from 2 types according to

the patient classification.

- User water trap and sampling line for MGU-800, water trap and flow sensor for MGU-810.
- Refer to "Chapter 13 Accessories" for the usable water trap, sampling line, or flow sensor. (@"Anesthetic Gas Concentration Measurement (Manufactured by Mindray Medical Sweden AB)" P13-7)

Set the "Wave Clip".

If the gas waveform amplitude exceeds the waveform display area, whether or not to clip the exceeded part can be selected.

- [ON]: The exceeded part of the waveform will be displayed in straight line at the upper or lower scale limit.
- [OFF]: The whole part of the waveform will be displayed even if it exceeds the scale. However, the exceeded part may not be displayed depending on the sweep speed of the waveform displayed above or below the gas waveform.



Set the "CO₂ Source Priority".

When MGU-800/MGU-810 and HS-8000/DS-8007/HM-801 are simultaneously used, the CO₂ source to prioritize the measurement can be set.

- ▶ [MGU-800]: CO₂ value measured by the MGU-800/810 Multigas Unit will be prioritized.
- ▶ [HS-8000]: CO₂ value measured by the HS-8000/DS-8007/HM-801 will be prioritized.

 When the HS-8000/DS-8007 and HM-801 are simultaneously used, the CO2 measurement of the HS-8000/DS-8007 will be prioritized.

Select the agent gas label from [Auto]/[ISO]/[SEV]/[HAL]/[ENF]/[DES].

• [Auto]: The label will be automatically set according to the detected anesthetic gas.

10 When the MGU-810 series is used, set the respiratory waveform scale.

Select ON/OFF for parameter display.

(@"ECG Parameter Setup" P7-6)

CAUTION

- When the waveform and numeric data display is set to OFF, the alarm generation and tabular/graphic trend input will also cease.
- When the waveform and numeric data display is set to OFF, the respiration rate measured by the multigas unit will not be displayed either.

MAC Display

The MAC value can be displayed in the numeric data display area.

- NOTE
- The MAC value will be displayed only if [ON] is set for "MAC Value". Perform the setting if necessary.

1 Press the [Menu], [GAS] "Parameter", [MAC] keys.

▶ The MAC value setup screen will be displayed.

2 Select ON/OFF for "MAC Value".

- [ON]: The MAC value will be displayed in the numeric data display area.
- ▶ [OFF]: The MAC value will not be displayed in the numeric data display area.
- ➤ To change the displayed default value, enter the value using the numeric keys, and press the key for the corresponding constant.

The MAC value is calculated from the following formula.

 $MAC = \frac{ExN_2O}{x(N_2O)} + \frac{ExPAGT}{x(PAGT)} + \frac{ExSAGT}{x(SAGT)}$ •Ex N₂O: Expired N₂O (%) •Ex PAGT: Expired Primary Agent (%) •Ex SAGT: Expired Secondary Agent (%) •X (N₂O): N₂O Constant •X (PAGT): Primary Agent Constant

•X (SAGT): Secondary Agent Constant

BIS Data (HBX-800 with BISx)

This section explains about the BIS measurement and setup procedure when using the BISx with the BIS I/F Unit, HBX-800.

WARNING

- Clinical judgment should always be used when interpreting BIS in conjunction with other available clinical signs. Reliance on BIS alone for intra-operative anesthetic management is not recommended.
- BIS values should also be interpreted cautiously with certain anesthetic combinations, such as those relying primarily on either ketamine or nitrous oxide/narcotics to produce unconsciousness.
- As with any monitored parameter, artifacts and poor signal quality may lead to inappropriate BIS values. Potential artifacts may be caused by poor skin contact (high impedance), muscle activity or rigidity, head and body motion, sustained eye movements, improper sensor placement and unusual or excessive electrical interference.

- Generally, the BIS value decreases with the decrease of brain activity. When a patient is in hypothermia state during cardiac bypass surgery, the suppression of brain wave will cause the BIS value to decrease.
- Pay attention when artifact interferes or signal quality decreases, as it may cause incorrect BIS measurement.
- · Pay attention when AC disturbing signal interferes during Filter OFF condition, as it may



cause incorrect BIS measurement.

- Pay attention when a pacemaker pulse is displayed in the brain wave, as it may cause incorrect BIS measurement.
- The BIS value tends to increase with the EMG interference. The patient's shivering during recovery from anesthesia increases the EMG and may case the BIS value to increase.
- When attaching the BIS sensor, lightly apply pressure to the electrode part for about 5 seconds to decrease the electrode impedance.

Preparation for Monitoring

By connecting the BISx module using the HBX-800 BIS I/F Unit, BIS data can be monitored.



- When using more than one HM-801, only one AUX connector on the HM-801 can be used for measurement of CO₂, BIS (with HBX-800).
- When the BIS data is measured on both the HS-8000/DS-8007 and HM-801, the measurement of the HS-8000/DS-8007 will be prioritized.

1 Select the appropriate sensor for the patient.

Connect the HBX-800 to the AUX connector on one of the following equipments, and connect the BISx to the serial communication connector on the HBX-800.

*HS-8000 Super Unit *DS-8007 *HM-801



3 Attach the BIS sensor to the patient.

When the system detects the sensor, "Sensor Check" window will be displayed, and impedance for all the electrodes will be automatically measured.

REFERENCE

• Pressing the [Sensor Check] key will also start the sensor check process.

> The measured results will be displayed on the "Sensor Check" window.



- ▶ In this display, the impedance value for each electrode, in kilo ohms, appears on the screen along with its status.
- <PASS>: An electrode passes if the impedance for that electrode is less than 7.5 kilo ohms, and the ground electrode (electrode #2) is less than 30 kilo ohms.
- <HIGH>: The impedance value is above 7.5 kilo ohms.

As long as the combined impedance of electrodes #1 and #3 and the combined impedance of electrodes #1 and #4 are less than 15 kilo ohms, and the ground electrode is less than 30 kilo ohms, the sensor check will be considered successful.

- <LEAD OFF>: The electrode is detached from the patient.
- <NOISE>: The signal from the electrode is outside the measurable range.

NOTE

During the sensor check process, EEG waveform will become unstable.

5 If the impedance for all the electrodes are within variable range, <Sensor Check Passed> will be displayed on the "Sensor Check" window.

 $\mathbf{6}$ Press the \mathbf{x} key on the "Sensor Check" window to end the sensor check process.

▶ BIS measurement will automatically start when the "Sensor Check" window is closed.

NOTE

If the "Sensor Check" window is closed before <Sensor Check Passed> is displayed,
 <BIS Perform "Sensor Check"> will be displayed. Press the [Sensor Check] key and start the sensor check again.

BIS Setup

Press the [Menu], [BIS] keys to display the "BIS" setup screen.



7 Scale

Select the EEG waveform scale from [±25]/[±50]/[±100]/[±250].

2 Alarm

Select ON/OFF of BIS alarm and set the alarm limits.

3 Short Trend 2nd Parameter

- Select the second parameter for short trend from [SR]/[EMG]/[SQI].
- Selecting [OFF] will not display the second parameter for short trend.

4 Continuous Impedance Check

Select whether or not to perform continuous impedance check.
 If [ON] is selected, the check process will continue until it passes.
 Select [OFF] if it affects other measurements.

- The conductive parts of sensors and connectors should not contact other conductive parts, including earth.
- To reduce the hazard of burns in the high-frequency surgical neutral electrode connection, the BIS sensor should not be located between the surgical site and the electrosurgical unit return electrode.
- The BIS sensor must not be located between defibrillator pads when a defibrillator is used on a patient connected to the patient monitor.
- Continuous impedance checking may need to be disabled if the 1 nA, 128 Hz impedance check signal interferes with other equipment (e.g., evoked potential monitors).
- Considerations when using Electro-Convulsive Therapy (ECT) equipment during BIS monitoring: Place ECT electrodes as far as possible from the BIS sensor to minimize the effect of interference. Certain ECT equipment may interfere with the proper function of the BIS monitoring system.

NOTE

- During the continuous impedance check, the following impedance will be measured.
- A) Combined Impedance of Signal Electrode and Reference Electrode This check process will not affect the EEG waveform. If the impedance value is within the allowable range, the check result will not be notified.
- B) Impedance of Ground Electrode This check process will be performed every 10 minutes. During this process, <Ground Check in Progress> will be displayed, as artifact interferes to the EEG waveform.

5 Smoothing Rate

Select from [10 sec.] / [15 sec.]/ [30 sec.].

EEG Filter

Select from [ON]/[OFF].

BIS Data (A-2000/A-3000)

This section explains about the BIS setup procedure when using the A-2000 BIS Monitor or A-3000 BIS Vista (Covidien).

On the BIS setup screen, the second parameter to be displayed on the short trend can be selected. The first parameter is fixed to BIS value.



Press the [Menu], [BIS] ("Parameter") keys to display the BIS setup screen.



Short Trend 2nd Parameter

- Select the second parameter for short trend from [SR]/[EMG]/[SQI].
- Selecting [OFF] will not display the second parameter.

2_{Trend E}

▶ Trend E screen will be displayed.

Ventilator

By connecting a ventilator, numeric data and waveform measured by the ventilator can be displayed on the DS-8400 System.

(Maintenance Manual "Ventilator Connection" P4-3)

By assigning [P-V/F-V] to numeric data box, P-V (pressure-volume) loop/F-V (flow-volume) loop can be also displayed.



This section explains about the AWP/AWF/AWV scale setup procedure and P-V/F-V screen operation.

AWP/AWF/AWV Scale Setup

Press the [Menu], [Ext. Device], ("Parameter"), [VENT] key to display the "VENT" screen. The ventilator measurement will be displayed, and AWF / AWP / AWV / P-V, F-V scale can be set.



REFERENCE

• The scale setup window can be also displayed by pressing the scale on the waveform display area or [Scale] on the user key.

Press the key for "AWF Scale", and set the AWF scale.



 $\mathbf 2$ Press the key for "AWP Scale", and set the AWP scale.



Press the key for "AWV Scale", and set the AWV (Airway Volume) waveform scale.



Press the key for "P-V, F-V Scale", and set the P-V (Pressure-Volume) scale.

- ▶ Pressure: Select from [10]/[20]/[30]/[50]/[120] (cmH₂O).
- ▶ Volume: Select from [250]/[500]/[750]/[1000] (mL).

Press the key for "P-V, F-V Scale", and set the F-V (Flow-Volume) scale.

- ▶ Flow: Select from [±20]/ [±50]/ [±180] (L/min) .
- Volume: Select from [250]/[500]/[750]/[1000] (mL).

Scale	(X)
Pressure(cnH20)	
10 20 30 50	120
Volume(mL)	
250 500 750 1000	
Flov(L/nin)	
(± 20) (± 50) (± 180)	

P-V/F-V Loop Display

The ventilator data can be displayed in P-V/F-V loop for review.

• For PURITAN-BENNETT ventilator, P-V loop and F-V loop cannot be displayed or printed.

Press the P-V/F-V numeric data box.

▶ The P-V/F-V review screen will be displayed.



- ▶ P-V (pressure-volume) loop/F-V (flow-volume) loop is sampled each 60ms and displayed for each respiration. The beginning of the loop is displayed in cyan, and the rest of the loop is displayed in white.
- ▶ For the P-V loop, the horizontal axis shows AWP (unit: cmH₂O), and vertical axis shows volume (unit: mL).
- ▶ For the F-V loop, the horizontal axis shows volume (unit: mL), and vertical axis shows AWF (unit: L/min).

2 Press the [Scale] key to set the P-V/F-V scale.

Select the scale from the displayed scale selection window.

 ${f 3}$ To stop the loop drawing, press the [Freeze] key.

- The loop drawing will stop.
- ▶ Press the [Freeze] key again to resume the waveform trace.

4 A control loop can be registered to see the change in P-V/F-V loop.

- Press the [Regist] key to store the displayed P-V/F-V loop as a control loop.
- Press the [Review] key to display the registered control loop. The control loop 1 will be displayed in yellow, and control loop 2 will be displayed in green.

FLOW-i Data

The FLOW-i can be connected to the serial port, status port of the DS-8400 system or to the HP-800. (@Maintenance Manual "Connection with the FLOW-i" P4-13)

When the FLOW-i is connected, monitoring conditions for CO_2 concentration, anesthetic gas concentration, O_2 concentration, N_2O concentration, and respiration can be set.



When the numeric data acquired from FLOW-i is displayed, the following alarms cannot be set. Also, these alarms will not generate.
 InspCO₂/EtCO₂, InspO₂/ExpO₂, InspN₂O/ExpN₂O, InspAgent/ExpAgent, MAC, ExpMV, PEAK, PEEP

- The FLOW-i and MGU-800/810 cannot be used simultaneously.
- The FLOW-i and ventilator cannot be used simultaneously.

CO2 Measurement Unit Setup

- (NOTE)
- The CO2 measurement unit is not linked between the FLOW-i and this equipment.
- When the FLOW-i is connected, CO₂ alarm cannot be set. Also, the alarm will not generate.

Press the [Menu], [Initial Settings], [Meas.], [Unit] keys.

> The "Unit" menu will be displayed.



Press the [mmHg]/[kPa]/[%] key.

> The data of currently set measurement unit will be displayed on the graphic/tabular trend.



FLOW-i Setup

1 Press the [Menu], [Anes.] "Parameter" keys.

> The anesthesia setup menu will be displayed.





FLOW-i_O₂ Setup



FLOW-i_N₂O Setup



FLOW-i_MAC Setup







FLOW-i_AGT Setup



FLOW-i_RESP Setup

Select from [0-50]/[0-100] if the measurement unit is mmHg, and from [0-4]/[0-8]/[0-10] if the unit is kPa or %.

 $\mathbf{3}$ Set the O₂ waveform scale.

▶ Select from [18-30]/[18-60]/[18-100]/[0-30]/[0-60]/[0-100].

4 Set the scale for anesthetic gas concentration.

▶ Select from [0-4]/[0-8]/[0-16].

5 Set the "Wave Clip".

If the gas waveform amplitude exceeds the waveform display area, whether or not to clip the exceeded part can be selected.

- [ON]: The exceeded part of the waveform will be displayed in straight line at the upper or lower scale limit.
- [OFF]: The whole part of the waveform will be displayed even if it exceeds the scale. However, the exceeded part may not be displayed depending on the sweep speed of the waveform displayed above or below the gas waveform.

6 Set the "CO₂ Source Priority".

When the FLOW-i and HS-8000/DS-8007/HM-801 are simultaneously used, the CO₂ source to prioritize the measurement can be set.

- ▶ [Anesthesia]: CO₂ value measured by the FLOW-i will be prioritized.
- ▶ [HS-8000]: CO₂ value measured by the HS-8000/DS-8007/HM-801 will be prioritized.

NOTE

• When the HS-8000/DS-8007 and HM-801 are simultaneously used, the CO₂ measurement of the HS-8000/DS-8007 will be prioritized.

Set the RR/APNEA alarm.

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

• Only the RR/APNEA alarm can be set. The following alarms cannot be set. Also, these alarms will not generate.

InspCO₂/EtCO₂, InspO₂/ExpO₂, InspN₂O/ExpN₂O, InspAgent/ExpAgent, MAC, ExpMV, PEAK, PEEP

NOTE

- Set the upper limit of RR alarm in the range of 10 Bpm to 60 Bpm. If a value above 150 Bpm is set, the upper alarm will turn OFF.
 Set the upper limit of apnea alarm in the range of 10 sec. to 60 sec. If a value above 60
- sec. is set, the upper alarm will turn OFF.Set the lower limit of RR alarm in the range of 5 Bpm to 145 Bpm. If a value below 5 Bpm
- Set the lower limit of RR alarm in the range of 5 Bpm to 145 Bpm. If a value below is set, the lower alarm will turn OFF.

REFERENCE

• The adjustable increment for RR alarm depends on the patient classification and "RR Alarm Increment" setting. (Initial Settings>User I/F).

	Alarm Increment (Initial Settings > User		
	Normal	Small	
Adult	5 Bpm increment	1 Bpm increment	
Child/Neonate	2 Bpm increment	1 Bpm increment	

• The apnea alarm can be set in 1 second increment.

 $\mathbf{8}$ Set the respiration waveform scale.

MAC Display

The MAC value can be displayed in the numeric data display area.

NOTE

The MAC value will be displayed only if [ON] is set for "MAC Value". Perform the setting if necessary.

1 Press the [Menu], [GAS] "Parameter", [MAC] keys.

▶ The MAC value setup screen will be displayed.





2 Select ON/OFF for "MAC Value".

- [ON]: The MAC value will be displayed in the numeric data display area.
- ▶ [OFF]: The MAC value will not be displayed in the numeric data display area.

96

38

21/ 16

30/ 0

Å [Put_] [122

60

7.5 20 PEAK 18 18

129/* 82

Uventilator Data Display and Setup

By connecting the FLOW-i, the numeric data and waveform measured by the ventilator can be displayed. By assigning [P-V/F-V] to numeric data box, P-V (pressure-volume) loop/F-V (flow-volume) loop can be also displayed.

This section explains about the AWP/AWF/AWV scale setup procedure and P-V/F-V screen operation.

Press the [Menu], [Anes.], ("Parameter"), [VENT] key to display the ventilator screen.

The ventilator measurement will be displayed, and AWF / AWP / AWV / P-V, F-V scale can be set.



REFERENCE

• The scale setup window can be also displayed by pressing the scale on the waveform display area or [Scale] on the user key.

1 Set the AWF scale.

1 Press the key for "AWF Scale".

- The scale selection for AWF (airway flow) waveform will be displayed.
- 2 Select from [±5]/ [±10]/ [±20]/ [±50]/ [±180] (L/min).

2_{AWP Scale}

- 1 Press the key for "AWP Scale".
 - The scale selection for AWP (airway pressure) waveform will be displayed.
- 2 Select from [10]/[20]/[30]/[50]/[120] (cmH₂O).

3AWV Scale

- 1 Press the key for "AWV Scale".
 - The scale selection for AWV (airway volume) waveform will be displayed.
- 2 Select from [50]/[250]/[500]/[1000]/[3000] (mL).







4 Set the P-V Scale.

- 1 Press the key for "P-V, F-V Scale".
 - The scale selection for P-V (pressure-volume) loop will be displayed.
- 2 Pressure: Select from [10]/[20]/[30]/[50]/[120] (cmH₂O).
- **3** Volume: Select from [250]/[500]/[750]/[1000] (mL).

5 Set the F-V Scale.

- 1 Press the key for "P-V, F-V Scale".
 - ▶ The scale selection for F-V (flow-volume) loop will be displayed.
- 2 Flow: Select from [±20]/ [±50]/ [±180] (L/min).
- 3 Volume: Select from [250]/[500]/[750]/[1000] (mL).

P-V/F-V Loop Review

The ventilator data of the FLOW-i can be reviewed in P-V/F-V loop display.

Press the P-V/F-V numeric data box on the home display.

▶ The P-V/F-V review screen will be displayed.



- ▶ P-V (pressure-volume) loop/F-V (flow-volume) loop is sampled each 60ms and displayed for each respiration. The beginning of the loop is displayed in cyan, and the rest of the loop is displayed in white.
- ▶ For the P-V loop, the horizontal axis shows AWP (unit: cmH₂O), and vertical axis shows volume (unit: mL).
- ▶ For the F-V loop, the horizontal axis shows volume (unit: mL), and vertical axis shows AWF (unit: L/min).

 $\mathbf{2}$ Press the [Scale] key to set the P-V/F-V scale.

Select the scale from the displayed scale selection window.

 ${f 3}$ To stop the loop drawing, press the [Freeze] key.

- ▶ The loop drawing will stop.
- Press the [Freeze] key again to resume the waveform trace.

4 A control loop can be registered to see the change in P-V/F-V loop.

- > Press the [Regist] key to store the displayed P-V/F-V loop as a control loop.
- Press the [Review] key to display the registered control loop. The control loop 1 will be displayed in yellow, and control loop 2 will be displayed in green.

Scale
Pressure(onH20)
Volume(mL)
250 500 750 1000
Flow(L/min)
$[\pm 20] [\pm 50] [\pm 180]$

SvO₂/CCO Data

The DS-8400 System can display the monitoring data of oximeter/CCO measurement device, Vigilance, Vigilance CEDV, Vigilance II, Vigileo, EV-100 (Edwards Lifescience) or the hemodynamic monitoring device, PiCCO2, PulsioFlex (PULSION Medical Systems).

(@Maintenance Manual "SvO2/CCO Monitor Connection" P4-6)

On the SvO₂/CCO data screen, the displayed numeric data can be switched.



Display Example for ICO Mode

STAT Mode: When the Vigilance is in CCO mode, STAT mode display can be set ON or OFF. Index Display: When the Vigilance is in CCO mode, Index Display can be set ON or OFF.

When the Vigilance is in ICO mode, the 6 latest data of ICO (Intermittent Cardiac Output) and ICI (Intermittent Cardiac Index) will be displayed.

STAT Mode / Index Display

Press the [Menu], [Ext. Device] ("Parameter") keys.

> The Vigilance screen will be displayed.



- Index Display: When Vigilance is in CCO mode, Index display can be set ON or OFF.
- When the Vigilance is in ICO mode, the 6 latest data of ICO (Intermittent Cardiac Output) and ICI (Intermittent Cardiac Index) will be displayed.

2 Select [ON]/[OFF] for "STAT Mode" and "Index Display".

▶ STAT Mode [OFF], Index Display [OFF]: SvO₂ (or ScvO₂), CCO, EDV, BT will be displayed inside the SvO₂+CO numeric data box.



STAT Mode [OFF], Index Display [ON]: CCI and EDVI will be displayed instead of CCO and EDV.



▶ STAT Mode [ON], Index Display [OFF]: CCO_STAT and EDV_STAT will be displayed instead of CCO and EDV.



► STAT Mode [ON], Index Display [ON]: CCI_STAT and EDVI_STAT will be displayed instead of CCO and EDV.



NOTE

ON/OFF of STAT mode can be changed only when Vigilance is connected.

INVOS Data

By connecting the INVOS 5100C Cerebral Oximeter (Covidien), regional cerebral oxygen saturation (rSO₂) can be monitored non-invasively on the DS-8400 System. (@Maintenance Manual "Connecting to the INVOS" P4-12)

On the INVOS screen (shown on right), the channel can be changed for each INVOS data.

Lt-rSO₂/Rt-rSO₂ data of the selected channel will be displayed inside the INVOS numeric data box.

Menu > Parameter > Ext. Device	ڪ[
SVO2 /CCO VENT INVOS	t
Explanation Area	
Lt-rS02 ch1	
Rt-rS02 ch2	
S1-rS0z ah3	
S2-rS0z ali4	

Channel Number Setup for INVOS Data

In the INVOS numeric data box, measurement data of $Lt-rSO_2/Rt-rSO_2$ will be displayed. On the INVOS screen, the channel for $Lt-rSO_2/Rt-rSO_2$ data can be selected.

1 Press the [Menu], [Ext. Device] ("Parameter"), [INVOS] keys.

• The INVOS screen will be displayed.

2 Press the [ch*] key for the INVOS label ("Lt-rSO₂" / "Rt-rSO₂" / "S1-rSO₂" / "S2-rSO₂") to set the channel.

• The dropdown list will be displayed.

3 Select the channel from [ch1]/[ch2]/[ch3]/[ch4].

Stopwatch

The stopwatch function can be used by setting the [Stopwatch] key on the numeric data box or on the user key.

1 Press the [Stopwatch] key on the numeric data box or on the user key.

• The "Stopwatch" window will be displayed.

STOPWATCH		X
TIMER1	00:00:00	START STOP RESET
TIMER2	00:00:00	START STOP RESET
	Label	

Label Setup

1 Press the [Label] key on the "Stopwatch" window.

• The stopwatch label setup window will be displayed.



2 Enter 8 characters using alphanumeric keypad.



1 Press the [Start]/[Stop]/[Reset] key on the "Stopwatch" window.

- [Start]: The stopwatch will start.
- [Stop]: The stopwatch will suspend/resume.
- ▶ [Reset]: The stopwatch will reset to "00:00:00". If pressed during stopwatch operation, counting will resume from "00:00:00".

NOTE

- If the discharge procedure is performed during stopwatch operation, the counting will stop and will be reset to "00:00:00".
- The stopwatch counting will continue even when the monitoring is suspended.

Multiparameter Connector Setup for BP, TEMP, CO Measurement

On the Super Unit and Multi Module, multiparameter connectors are provided. The quantity of multiparameter connectors are as follows.

Multiparameter Connectors	Super Unit	Multiparameter Connectors	Super Unit
<u>3 ports</u>		2 ports	
TEMPx6 (maximum) BPx6 (maximum) COx1 (maximum)	HS-8312N, HS-8312M	TEMPx6 (maximum) [*] BPx4 (maximum) COx1 (maximum)	DS-8007N, DS-8007M

Multiparameter Connectors	Multi Module	Multiparameter Connectors	Multi Module
2 ports		<u>1 ports</u>	
TEMPx4 (maximum) BPx4 (maximum) COx1 (maximum)	um)	TEMPx2 (maximum) BPx2 (maximum) COx1 (maximum)	HM-801

*: TEMPx2 are fixed jacks.

By using the multiparameter connector, any combination of BP, TEMP and CO measurement can be performed according to the monitoring purpose.

By using the 2ch TEMP relay cable, 2ch BP relay cable, or 2ch BP conversion cable, 2 channels of temperature and BP can be monitored through one multiparameter connector.

By using the Multi Module with the Input Box, up to 8 channels of BP, 8 channels of TEMP and 1 channel of CO can be measured.

The multiparameter connector setup can be performed on the "Initial Settings" menu.

(Maintenance Manual "Unit Module Setup" P4-17)

Gr HS-8312N, HS-8312M

3 Ports	Blood Pressure	Temperature	СО
Blood Pressure	6ch (3ch)	-	-
Blood Pressure			
Blood Pressure			
Blood Pressure		2ch	
Blood Pressure	4ch (2ch)		-
Temperature			
Blood Pressure	2ch (1ch)	4ch	-
Temperature			
Temperature			
Temperature			
Temperature	- 6ch	-	
Temperature			
Blood Pressure			
Temperature	2ch (1ch)	2ch	1ch
CO			
Blood Pressure	4ch (2ch)		
Blood Pressure		_	1ch
CO			

Combination of BP, TEMP, CO Channels

Combination of BP, TEMP, CO Channels

3 Ports	Blood Pressure	Temperature	СО
Temperature			
Temperature	-	4ch	1ch
CO			

The numbers in parenthesis shows the channels when using the 1ch BP conversion cable.

Group DS-8007N, DS-8007M

2 Ports	Blood Pressure	Temperature	CO	
Blood Pressure	4ch	_		
Blood Pressure	(2ch)	-	-	
Blood Pressure	2ch	2ch	_	
Temperature	(1ch)	2011		
Temperature	_	4ch	_	
Temperature		1011		
Blood Pressure	2ch	_	1ch	
СО	(1ch)		1011	
Temperature	_	2ch	1ch	
CO		2011	1011	

Combination of BP, TEMP, CO Channels

The numbers in parenthesis shows the channels when using the 1ch BP conversion cable.

On the DS-8007, there are 2 temperature connectors in addition to the multiparameter connectors.

Multiparameter Connector Setup

It is necessary to manually set the measuring parameter for each multiparameter connector.

Example:

To assign BP5 to multiparameter connector 1 for the HM-800 Multi Module inserted to the built-in slot:



- 1 Press the [Menu], [Initial Settings], [System], [Unit Module], [Multiamplifier] keys.
- 2 Select the multiparameter connector location. The selected location will be displayed in blue.

3 Assign the parameter to the selected location. In this case, select [5] for "BP". The parameter will be assigned to the selected connector.

- The same parameter cannot be set to more than one connectors.
- By setting [OFF] for the unused connector, the parameter will become selectable on another connector.
- If the parameter assigned to the multiparameter connector and the connected relay cable does not match, the connector location will be displayed in red and the connected relay cable type will be displayed.

NOTE

 The temperature measured on the temperature connectors on the DS-8007 are fixed as T1, T2.

Chapter 8 Review Function

Common Operation	8-1
Arrhythmia Analysis	8-2
Arrhythmia Definition	8-2
Arrhythmia Alarm Setup	
Arrhythmia Learn	8-5
Graphic Trend	8-6
Graphic Trend Setup	
Description for Each Parameter	8-9
Short Trend	8-10
Tabular Trend	
To Display/Print the Tabular Trend	
The Description of the Display	
Parameter Setup for Tabular Trend	
Recall	
To Display the Recall Waveform	
To Display/Print the Enlarged Recall Waveform	
Recall Setup	8-19
OCRG	
Alarm History	8-21
Alarm History Setup	8-21
Description for Each Item	
Zoom Wave	8-23
ST Measurement	
To Display/Print the ST Measurement	8-24
Reference Waveform Setup	8-25
ST Alarm Setup	
12-Lead Analysis	
12-Lead ECG Display	
12-Lead Analysis Setup	
12-Lead ECG Analysis	
12-Lead Analyzed Result Output Example	
Full Disclosure Waveform (Optional Function)	
To Format the CF Card	
Waveform Setup	
Description of the Full Disclosure Waveform Display	
To Search by Time	
Hemodynamics	
Calculation Data	
To Display/Print the Hemodynamics Data	
New Input of Hemodynamics Calculation	
To Edit the Hemodynamics Input Data	
Lung Function	
Calculation Data	
To Display/Print the Lung Function Data	
New Input of Lung Function Calculation	
To Edit the Lung Function Input Data	
Cardiac Output (CO)	
To Display the CO Measurement Screen	
Cardiac Output Setup	
CO Measurement	
To Edit the CO Measurement Result	
Drug Calculation	ŏ-51

Other Bed Display	8-54
Other Bed Display/Alarm	

Chapter 8 Review Function

Common Operation



The common operations for all the review screens are explained below.

1 Time Bar

• Changing the time span, scrolling the time, displaying the latest data can be performed.



- 1 The display can be switched in 24 hours interval.
- 2 Pressing the $\boxed{\cdot }$ $\boxed{}$ / $\boxed{}$ key will move the cursor to the alarm generated time.
- 3 The time zone for the whole data is shown. \blacklozenge indicates the alarm occurrence point. The lower row shows the time zone for the displayed data. Pressing the time bar will display the data at pressed time.
- **4** Indicates the displayed time range with the bar length. Dragging the slider to the right will display newer data, and dragging it to the left will display older data.
- 5 Pressing the Latest will display the latest data.
- 6 Pressing the \mathbb{H} / \mathbb{H} will switch the display by page.
- 7 Pressing the [] / [] will switch the display by 1 data/block each.



2 Displays other review data at the same time.

▶ With the displayed date/time, the review data display can be switched. Other review data (graphic/tabular trend) can be displayed without moving the current cursor time.

Arrhythmia Analysis

This section explains about the arrhythmia analysis.

Arrhythmia Definition



The arrhythmia detection is performed by learning the normal waveform of the patient, and determines the VPC by comparing the waveform (QRS pattern) and R-R interval for each heartbeat.

The parameters such as QRS amplitude, QRS width, QRS polarity, RR interval are compared with the normal waveform to extract the abnormal QRS.

Then, the QRS with suspected VPC is pattern matched. The noise and VPC are distinguished to determine the VPC, and generates the arrhythmia alarm.

WARNING

• Objective and constant arrhythmia detection is possible through the fixed algorithm incorporated in this monitor.

However, excessive waveform morphology change, motion artifact, or the inability to determine the waveform pattern may cause an error, or fail to make adequate detection. Therefore, physicians should make final decisions by closely checking the data obtained by manual printing, alarm printing and recall waveform.

• For stable arrhythmia detection and ECG monitoring, verify proper electrode placement, lead selection, and ECG waveform size. If necessary, turn ON the AC filter. If not properly selected, it may cause erroneous detection.

QRS Classification

Each QRS will be classified to the following pattern.

N (Normal)	Normal QRS beat
V (VPC)	Ventricular extrasystole
P (Pacing Beat)	Pacing beat
F (Fusion Beat)	Fusion beat of pacing and spontaneous beat
S (SVPC)	Supraventricular extrasystole
? (Undetermined Beat)	Learning arrhythmia, or unmatched beat

Arrhythmia Type

With the QRS judgment, the following types of arrhythmia alarm will be generated.

Arrhythmia	Description	Detection Criteria
Asystole	ON 3 sec. to 10 sec., 1 sec. increments	Cardiac arrest is detected for more than preprogrammed time.
VF	ON	A random, rapid electrical activity of the heart is detected.
VT (Ventricular Tachycardia)	ON	9 or more continuous VPC beats are detected.*1
Slow VT	ON, OFF	9 or more continuous VPC beats are detected. ^{*2}
Run (Consecutive VPC)	ON, OFF 2 beats to 8 beats, 1 beat increments	Continuous VPC exceeding the preprogrammed value (2 beats to 8 beats) is detected. ^{*3}
Couplet (Couplet VPC)	ON, OFF	2 continuous VPC beats are detected.
Pause	ON, OFF 1.5 sec. to 5.0 sec., 0.5 sec. increments	Cardiac arrest exceeding the preprogrammed duration is detected.
Bigeminy (Ventricular Bigeminy)	ON, OFF	QRS pattern of V-x-V-x-V-x is detected.*4
Trigeminy	ON, OFF	QRS pattern of x-x-V-x-x-V is detected.*4
Frequent (Frequent VPC)	ON, OFF 1 bpm to 50 bpm, 1 beat increments	VPC exceeding the preprogrammed value is detected within 1 minute.
Tachy(Tachycardia)	ON, OFF	The upper HR alarm limit is exceeded.
Brady (Bradycardia)	ON, OFF	The lower HR alarm limit is exceeded.
Ext Tachy (Extreme Tachycardia)	ON, OFF 22 bpm to 300 bpm, 22 bpm to 60 bpm, 1 beat increments 60 bpm to 300 bpm, 5 beat increments	The upper alarm limit of extreme tachycardia is exceeded.
Ext Brady (Extreme Bradycardia)	ON, OFF 20 bpm to 295bpm 20 bpm to 60 bpm, 1 beat increments 60 bpm to 295 bpm, 5 beat increments	The lower alarm limit of extreme bradycardia is exceeded.
R on T (R on T VPC)	ON, OFF 200 ms to 600 ms, 8 ms increments	VPC is detected within the preprogrammed RR interval (200 ms to 600 ms).
Multiform (Multiform VPC)	ON, OFF	2 different forms of VPC beats are detected within 4 minutes.
Vent Rhythm (Ventricular Rhythm)	ON, OFF	Continuous VPC beats with HR below the set value for "HR Lower Limit for Run" (0 bpm to 100 bpm), and same or above value of the set beats for Run (2 beats to 8 beats) are detected.
SVT (Supraventricular Tachycardia)	ON, OFF 2 beats to 10 beats, 1 beat increments	Continuous SVPC exceeding the preprogrammed value (2 beats to 10 beats) is detected.
Irregular RR (Irregular RR Interval)	ON, OFF 10% to 20%, 5% increments	RR interval variability exceeding the preprogrammed value (10% to 20%) is detected.
Prolonged RR (Prolonged RR Interval)	ON, OFF	RR interval of 1.75 times longer than the normal RR interval is detected.
Pacer Not Capture (Non- Capture)	ON, OFF 80 ms to 480 ms, 8 ms increments	HR is not detected from the pacing pulse within the set duration.
Pacer Not Pacing (Oversensing)	ON, OFF 20 bpm to 200 bpm, 20 bpm to 150 bpm, 5 beat increments 150 bpm to 200 bpm, 10 beat increments	Pacing pulse and HR are not detected during the set instant HR.
Triplet (Triplet VPC)	ON, OFF	3 continuous VPC beats are detected.
S Frequent (Frequent SVPC)	ON, OFF 1 bpm to 50 bpm, 1 beat increments	SVPC exceeding the preprogrammed value is detected within 1 minute.

S Couplet (Couplet SVPC)	ON, OFF	2 continuous SVPC beats are detected.
VPC (Ventricular Extrasystole)	ON, OFF	VPC is detected.
SVPC (Supraventricular Extrasystole)	ON, OFF	SVPC is detected.

*1: HR of 140 bpm/120 bpm and above

*2: HR of 100 bpm to 140 bpm or 100 bpm to 120 bpm

*3: HR of same or above the set value of "HR Lower Limit for RUN" (0 bpm to 100 bpm)

*4: x indicates N, P, F, ?.

Arrhythmia Alarm Setup

Arrhythmia alarm setup procedure is explained below.

ON/OFF of arrhythmia alarm and arrhythmia detection level can be set.

When the measured value exceeds the set arrhythmia detection level, arrhythmia alarm will generate.

Arrhythmia Detection Level Setting

Item	Description
Asystole	3 sec. to 10 sec.
Run	2 beats to 8 beats
Pause	1.5 sec. to 5 sec.
Frequent	1 bpm to 50 bpm
Ext Tachy	22 beats to 300 beats
Ext Brady	20 beats to 295 beats

Item	Description
R on T	200 ms to 600 ms
SVT	2 beats to 10 beats
Irregular RR	10, 15, 20%
S Frequent	1 beat to 50 beats
Pacer Not Capture	80 ms to 480 ms
Pacer Not Pacing	20 bpm to 200 bpm

1 Press the [Menu], [Arrhy.] ("Alarm") key.

The arrhythmia alarm setup screen will be displayed.

$\mathbf{2}$ Set the detection level.

Set using the dropdown list, numeric keys, or displayed key selection.

3 Select ON/OFF for the alarm.

- ▶ [ON]: Alarm will generate.
- ▶ [OFF]: Alarm will not generate.

Menu > Alarm Basic Circ. Resp./Gas Explanation Area	Arrhy. ST List Detail Setup	
Asystole ON 5 sec.	Tachy ON	
VF ON	Brady ON	
VT (HR > 120bpm) ON	Run (HR > 40 bpm) ON 3 beats	
Ext Tachy OFF 150 bpm	Pause & OFF 3.0 sec.	Detail Setup
Ext Brady OFF 30 bpm	Triplet & OFF	
SLOW VT DN	Couplet AFF	•••

Arrhythmia Alarm Detail Setup

On the "Detail Setup" of arrhythmia alarm, HR Lower Limit for VT, RUN, and SVT can be set .

Press the [Menu], [Arrhy.] ("Alarm"), [Detail Setup] key.

> The "Detail Setup" window for arrhythmia alarm will be displayed.

 $\mathbf{2}$ Set the "HR Lower Limit for VT".

- ▶ Set the VT analyzing condition for the arrhythmia analysis. VT alarm will generate if the HR is same or above the set value (120 bpm/140 bpm). Slow VT alarm will generate when the HR is below the set value.
- Select from [120] or [140] (bpm).





3 Set the "HR Lower Limit for Run".

- > Set the Run analyzing condition for the arrhythmia analysis. Run alarm will generate if the HR is same or above the set value.
- ▶ Press the || Press the || | | keys for "HR Lower Limit for Run" to set the HR in the range from 0 bpm to 100 bpm.

4 Set the "HR Lower Limit for SVT".

- > Set the SVT analyzing condition for the arrhythmia analysis. SVT alarm will generate if the HR is same or above the set value.
- ▶ Press the ▲/▼ keys for "HR Lower Limit for SVT" to set the HR in the range from 100 bpm to 250 bpm.

Arrhythmia Learn

Learning the normal ECG largely affects the accuracy of arrhythmia analysis.

If any error occurs in arrhythmia detection and QRS judgment, performing arrhythmia learning will recover the original analyzing accuracy.

Arrhythmia learning will be performed for about 20 beats for the normal ECG, but it may take longer if the heartbeat is unstable.

During the arrhythmia learning procedure, arrhythmia alarm other than Asystole, VF, Pause, Tachy, Brady will not generate.

Press the [Menu], [ECG] "Parameter" keys.

Or, press the HR numeric data box, and press $(\overline{s_1})$.

The ECG setup screen will be displayed.



09:27

60

LEARN

 $\mathbf{2}$ Press the [Learn] key while displayed in gray.

- The key will change to blue.
- Arrhythmia learning will start.
- During arrhythmia learning, a message will be displayed.

NOTE

- If [Used] is selected for "Pacemaker", the [Learn] key will not change to blue and <LEARN> will not be displayed, but the learning process will be performed.
- Pressing the key while arrhythmia learning is in process will not stop the process.

Graphic Trend

This section explains the graphic trend function and printing procedure. If the numeric data is displayed on the home display, 24 hours of data will be automatically stored and displayed as trend data.

Graphic Trend Setup

1 Press the [Menu], [Trend] ("Data Review") keys.

Or, press the [Graphic Trend] key on the user key area.

The graphic trend will be displayed.



2 graphs are displayed on each page, and graphic trend of 4 parameters can be displayed simultaneously on each graph.

2 Changing the displayed time, scrolling the time, updating the data (

3 Set the parameter, display type, scale.

1 Press the scale area for each parameter, and display the scale selection window.



2 Press the key for "Parameter Selection", and select the parameter.

	e key loi 1 alameter Selection	i, and select the para	ameter.
		Parameter	X
	H Hodule Sv02/CCO Vent.	Other	OFF
	HR ST(I) ST(I)	ST(III) ST(aVR) ST(aVL) ST(aVF)	·)
	VPC ST(V) ST(V2)	\$1(¥3) \$1(¥4) \$1(¥5) \$1(¥6)	
	WIBP BP1 BP2	PR_IBP T1 T2	
	SP02 PI		
	NOTE)		
	 The selected parameter v 	vill be also registered	for the trend group.
	 The apnea duration will b 	e stored when it exce	eds the upper alarm threshold level. I
	lower than the alarm thre		
3 Select the	e scale.		
4 Press the	e key for "Display Selection", a	and select the display	, type
			.,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,
4 Move the cur	rsor		
,		_	
7 Pressing	the center part of	will display the tree	nd data at the cursor position.
2 The curso	or will move to left and right by	y dragging	
3 Press 🗲	/ 🔼 to adjust the cursor po	sition.	
	ata display at cursor position v		erased after fixed duration
	_	-	
4 Press ⊞	to display the 10-minute tre	nd data before and af	ter the cursor position.
5 Press 🔍	to return the display to the p	revious time range.	
	2		
5 Set the displa	ay range.		
	EFERENCE)		
•	The displayed data is compr	essed as follows dep	ending on the display interval.
	VPC: Maximum value within		
	APNEA: Maximum value wit		
	Other than above: Latest va		
	•	•	er with minimum resolution of 1 minute
	is displayed, one mark will b	e displayed for the 12	2-minute (720-second) data.

 If the display resolution is higher than the minimum resolution of the data, the same data is repeated to match the display resolution.
 Refer to the following table for resolution. The data resolution differs according to the parameter.

Display Resolution					
	Minimum Resolution				
Time Span	Line Display		Mark Display	(Mark: Small)	
	10 sec. Sample	30 sec. Sample	10 sec. Sample	30 sec. Sample	
10 min.	10 sec.	30 sec.	10 sec.	30 sec.	
1 hours	10 sec.	30 sec.	30 sec.	30 sec.	
2 hours	10 sec.	30 sec.	60 sec.	60 sec.	
4 hours	20 sec.	60 sec.	120 sec.	120 sec.	
8 hours	40 sec.	120 sec.	240 sec.	240 sec.	
12 hours	60 sec.	120 sec.	360 sec.	360 sec.	

Display Resolution

Display Resolution

	Minimum Resolution			
Time Span	Line Display		Mark Display	(Mark: Small)
	10 sec. Sample	30 sec. Sample	10 sec. Sample	30 sec. Sample
16 hours	80 sec.	240 sec.	480 sec.	480 sec.
24 hours	120 sec.	240 sec.	720 sec.	720 sec.

Data Resolution

Minimum Resolution	Parameter	
10 sec.	HR, ST, SpO ₂ , PR_SpO ₂ , BP1, BP2	
30 sec.	Other than above (Excluding NIBP [*])	

* Actual measured data will be displayed for NIBP.

• Press the [TREND-x] key to change the trend group. Maximum of 5 groups with 8 parameters each can be registered, and can be selected according to the monitoring purpose.

- 1 Select the group.
- 2 To change the name of trend group, press the [Change Name] key. (@"To Enter Characters" P3-25)
- Group TREN EMG Off OFF Soi BP1 NIBP EtCO2 nspag1 BIS Off OFF TREND-F OFF OFF OFF OFF HR BP1 T1 NIBP ST(II) BR_GAS Sp02 EtC02 OFF OFF OFF OFF Change Node Nam
 - Settor
 3

 1
 Time Bar
 Zeh
 Base common Black
 Wark
 Shall
 4

 2
 Image: State of the state of

- **7**Perform the setup for trend display.
 - 1 Time Bar
 - Select the time bar display interval from [4h]/[8h]/ [12h]/[16h]/[20h]/[24h]/[36h]/[48h].
 - 2 Alarm Display Selection

Select the alarm display status. If the alarm for the selected arrhythmia, parameter is generated during the displayed time range, it will be

indicated in red at the alarm status display area.

- ▶ [Trend Parameters]: The displayed trend parameters will be selected.
- [Select All]: All parameters including arrhythmia will be selected.
- ▶ [Cancel All]: All selections will be canceled.
- [Select All Arrhythmia]: All arrhythmia will be selected.
- Each parameter key: Each time the key is pressed, selected/unselected status will change.
- **3** Background Color
 - Select the background color of the graphic trend from [White]/[Black]/[Gray].
- 4 Mark
 - Select the mark size on the graphic trend from [Small]/[Big].

8 Press the [Print] key.

▶ To print the trend data, press the [Print] key, select the parameter, and press the [Enter] key.

Description for Each Parameter

Numeric Data	Description	Scale	Unit
HR	HR	100, 200, 300	bpm
VPC	VPC Counts	20, 50, 100	-
ST (I, II, III, aVR, aVL,		±0.2, ±0.5, ±1.0, ±2.0	mV
aVF, V1 to V6)	ST Level	±2, ±5, ±10, ±20	mm
SpO 2-1, SpO2-2	SpO ₂ Value	0 to 100, 50 to 100, 80 to 100	%SpO ₂
PR_SpO ₂ -1, PR_SpO ₂ -2	SpO ₂ Pulse Rate	100, 200, 300	bpm
RR_SpO ₂	SpO ₂ Respiration Rate	pO ₂ Respiration Rate 50, 100, 150	
NIRD	NIBP Value (SYS / DIA)	100, 150, 200, 300	mmHg
NIBP		16, 20, 24, 40	kPa
		20, 50, 100, 150, 200, 300	mmHg
BP1~8	Blood Pressure (Systolic / Mean / Diastolic)	4, 8, 16, 20, 24, 40	kPa
		20, 40	cmH ₂ O
222		20, 50, 100, 150, 200, 300	mmHg
PDP	Peak Diastolic Pressure of IABP	4, 8, 16, 20, 24, 40	kPa
		20, 50, 100, 150, 200, 300	mmHg
CPP	Cerebral Perfusion Pressure	4, 8, 16, 20, 24, 40	kPa
		20, 50, 100, 150, 200, 300	mmHg
PAP	Pulmonary Artery Pressure	4, 8, 16, 20, 24, 40	kPa
PR_IBP	BP Pulse Rate	100, 200, 300	bpm
T1 to 8	Temperature	20.0 to 45.0, 30.0 to 40.0	°C
Tb	Blood Temperature (Cardiac Output Measurement)	20.0 to 45.0, 30.0 to 40.0	°C
ΔTEMP-A to D	Temperature Difference	±10.0, ±25.0	°C
RR_IMP	Impedance Respiration Rate	50, 100, 150	Bpm
Apnea	Apnea Duration (Impedance, CO ₂ , Ventilator)	15, 30	s (second)
EtCO ₂ , InspCO ₂ ^{*1}	Gas Unit CO ₂ Concentration	50, 100	mmHg
$E(CO_2, IIISpCO_2)$		4, 8, 10	kPa, %
ExpO ₂ , InspO ₂ ^{*1}	Gas Unit O ₂ Concentration	50, 100	%
ExpN ₂ O, InspN ₂ O ^{*1}	Gas Unit N ₂ O Concentration	50, 100	%
RR_GAS ^{*1}	Gas Unit Respiration Rate	50, 100, 150	Bpm
ΔO ₂ ^{*1}	ΔΟ2	3, 6, 9	%
ExpAGT, InspAGT ^{*1}	Gas Unit Agent Concentration	4, 8, 10	%
MAC ^{*1}	Minimal Alveolar Concentration	5, 10	-
BIS	Bispectral Index (BIS Monitor Measurement)	25, 50, 75, 100	-
SR	Suppression Ratio (BIS Monitor Measurement)	25, 50, 75, 100	%
EMG ^{*1}	Electromyography (BIS Monitor Measurement)	30 to 80	dB
SQI	Signal Quality Index (BIS Monitor Measurement)	0 to 100	%

Numeric Data	Description	Scale	Unit
SvO2 ^{*2}	Mixed Venous Oxygen Saturation	0 to 100, 50 to 100, 80 to 100	%
ScvO2 ^{*2}	Central Venous Oxygen Saturation	0 to 100, 50 to 100, 80 to 100	%
CCO ^{*2}	Continuous Cardiac Output	6, 12, 20	L/min
CCI ^{*2}	Continuous Cardiac Index	6, 12, 20	L/min/m ²
BT ^{*2}	Blood Temperature (SvO ₂ /CCO Monitor)	20 to 45, 30 to 40	°C
RR_VENT	Ventilator Respiration Rate	50, 100, 150	Bpm
SpCO (1, 2)	Carboxyhemoglobin Concentration	20, 40, 100	%SpCO
SpMet (1, 2)	Methemoglobin Concentration	10, 15, 100	%SpMet
SpHb (1, 2)	Total Hemoglobin Concentration	10 to 20, 0 to 25	g/dL
PI(1,2)	Perfusion Index	10, 20	%
PI(1,2)	Pleth Variability Index	30, 60, 100	%
ExpMV ^{*1}	Expiratory Minute Ventilation Volume	6.0, 12.0, 20.0	L/min
PEAK ^{*1}	Peak Airway Pressure	10, 20, 50, 100	cmH ₂ O
PEEP ^{*1}	Peak End Expiratory Pressure	10, 20, 50, 100	cmH ₂ O
Lt-rSO2 ^{*2}			
Rt-rSO ₂ ^{*2}	Regional Cerebral Oxygen Saturation	20 to 100	%
S1-rSO2 ^{*2}			70
S2-rSO2 ^{*2}			

*1: When the FLOW-i Anesthesia Delivery System is used, the measurement by the FLOW-i will be used.

*2: The external device parameters to be displayed on the graphic trend/tabular trend needs to be selected in advance on the "Trend Data Setup" window ([Data Review>Graphic Trend or Tabular Trend] or [Initial Settings>External Device>Main Unit/HP-800]).

NOTE

 The apnea duration will be stored when it exceeds the upper alarm threshold level. If lower than the alarm threshold level, it will be stored as "0 (zero)".

Short Trend

The trend data can be displayed on the home display.

As the alarm occurrence point on the graph is displayed in red, the alarm data of up to 3 hours (*) can be verified on the home display.

(*: In case when the short trend data resolution is set to [30 sec.].)

Pressing the short trend of an alarm generated parameter will display the recall screen.

The short trend can be displayed for each display layout. When 12-lead layout is displayed, ST value of each lead can be displayed in short trend.



The short trend display can be turned ON or OFF using the [Short Trend ON/OFF] user key. (@"User Key Selection" P10-15)

NOTE

- When the cursor function or reference line function of the short trend is enabled, the alarm condition cannot be displayed.
- · When the short trend of multiple parameters are displayed overlapped, only the parameter

8-11

displayed on top will be displayed in red at alarm occurrence.

Selecting the Parameters to be Displayed

The parameters to be displayed can be changed on the "Display Config." menu. (@"Display Configuration" P10-1)

Also, by setting the auto display configuration, the short trend parameters can be changed automatically according to the displayed waveforms and numeric data. (> Maintenance Manual "Display/Print Setup" P5-13)

Maximum of 4 parameters can be displayed overlapped in the same short trend display area. (shown on right)

Changing the Trend Scale and Display Duration

The short trend scale will be displayed on the right or left side of the short trend.

The displayed scale will be in accordance with the scale set on the "Trend" screen.

For the following parameters, the short trend scale can be synchronized with the corresponding waveform scale by selecting [Waveform] for "Short Trend Scale" under [Menu>Display Config.>Detail Setup].

BP, PEAK, TV, CO₂, O₂, Agent

The short trend display width can be selected from 7 levels by pressing the waveform display area.

Also, by setting the "Data Resolution" (5 sec./10 sec./30 sec.) under [Display Config.] > [Detail Setup], maximum display duration (30 min./1 hr./3 hr.) can be changed. The display width can be selected from 7 levels.

Changing the Display for Each Parameter

The graph type and display order can be changed for each parameter.

By pressing the short trend scale area, "Short Trend Setup" window (shown on right) will be displayed.

"Display Selection"

Select the graph type.

- For example, there are following graph types.
 - Line 📈
 - Filled in with black color from the baseline
 - Filled in with black color between S-D (For BP)
 - Filled in with black color from the top
- [OFF]: Graph will not be displayed.

The displayable graph types will differ depending on the parameter.

♦"Display Order"

When the parameters are displayed overlapped (ex. short trend overlap, BP overlap), the display order can be selected.

- [Front]: The display will be on the front side.
- [Back]: The display will be on the back side.

Displaying the Reference Line

For the short trend of the following parameters, reference lines can be displayed.

- HR (Upper/Lower Limit)
- ST (Upper/Lower Limit) *Only for the ECG1 lead
- BP1 to 4 (Upper/Lower Limit) *S/D/M can be selected for each limit.
- NIBP (Upper/Lower Limit) *S/D/M can be selected for each limit.
- EtCO₂ (Upper/Lower Limit)

	Short Trend Setup	X			
HR PR_Sp02 PR_Sp02H2 PR_IBP VPC ST(I)					
Display Selection	Reference Line OFF Lover Limit	Reference Line Upper Linit OFF			
Display Order Front	40	Alarm Linit			
Back	▼ Current Value	▼ Current Yalue			





- SpO₂ (Lower Limit)
- BIS (Upper/Lower Limit)

The data within the reference lines (including the parameters without the reference line display) will be displayed with lower brightness.

The data outside the reference lines will be displayed with higher brightness.

The reference lines can be displayed by selecting [Enable] for "Reference Line Function". (Menu>Display Config.>Detail Setup)

However, it cannot be displayed for the overlapped short trend. And, when the reference line function is enabled, the function to display the alarm occurrence point on the graph in red cannot be used.

When [Enable] is set for "Reference Line Function", ON/ OFF and upper/lower limit of reference line display can be selected on the "Short Trend Setup" window for each parameter.

The "Short Trend Setup" window can be displayed by pressing the short trend scale area.



By displaying a cursor, the numeric data and review data at cursor position can be displayed.

The cursor can be displayed by selecting [Enable] for "Cursor Function". (Menu>Display Config.>Detail Setup)

Pressing the short trend display area will display the cursor at the last displayed position (time). If the last displayed position is

cleared by scrolling, the cursor will be displayed at the latest data position.

The cursor can be moved by dragging or pressing the short trend display area.



(However, zoom wave can be displayed only when the full disclosure waveform function is enabled.)

The cursor cannot be displayed for the overlapped short trend. And, when the cursor function is enabled, the function to highlight the alarm generated data cannot be used.

When the cursor function is enabled, the function to enlarge/reduce the short trend display area cannot be used. During the cursor display, the short trend data will not be updated. When the cursor is not used for 10 seconds or when other window is displayed, the cursor will be automatically cleared.




Tabular Trend

This section explains the tabular trend function and printing procedure.

If the numeric data is displayed on the home display, 24 hours of data will be automatically stored and displayed in 10 seconds/30 seconds interval.

To Display/Print the Tabular Trend

Press the [Menu], [Tabular Trend] ("Data Review") keys. Or, press the [Tabular Trend] key on the user key area.

• The tabular trend will be displayed.



2 Changing the displayed time, scrolling the time, updating the data ((Common Operation "P8-1)

3 The list will be scrolled up or down to display other parameters.

4 Select the display interval.

[NIBP]: The tabular trend display interval will be according to the NIBP measurement time.

NOTE
If the display resolution is higher than the minimum resolution of the data, the same data is repeated to match the display resolution.

- The data resolution differs according to the parameter.
- 24 hours of data will be stored regardless of the time bar display range.
- The minimum resolution of NIBP is 20 seconds. If multiple measurements have been performed within 20 seconds, the latest value will be displayed.

Data Resolution

Minimum Resolution	Parameter
10 sec.	HR, ST, SpO ₂ , PR_SpO ₂ , BP1, BP2
30 sec.	Other than above

5 Press the [List-x] key to change the tabular trend group. Maximum of 6 different groups of parameters can be registered according to the monitoring purpose.

- 1 Select a group from [A]/[B]/[C]/[D]/[E]/[F].
- 2 To change the name of trend group, press the [Change Name] key.
 - (@"To Enter Characters" P3-25)

6 Set the parameters for the tabular trend.

(Parameter Setup for Tabular Trend" P8-15)

Press the [Print]/[Print (All)] key.

- [Print]: The currently displayed tabular trend will be printed.
- [Print (All)]: All data for 12 parameters (which fits in 1 page) will be printed.

 $oldsymbol{\delta}$ The displayed data time can be scrolled by dragging the data display area.

The Description of the Display

For the data when the measurement was not performed (before admittance) or when the monitoring was suspended, the time will be displayed as " : ".

Also, if the measured data is not displayed on the home display, or BP zero balance is not performed, the data will not be displayed.

The alarm generated data will be displayed with red background.

The date column of alarm generated data will be also displayed with red background.



• The red background will be displayed for the alarm generated parameter.

The alarm display for the expiratory and inspiratory parameter such as $\rm EtCO_2$ and $\rm InspCO_2$ will be the same.

For example, if the alarm is generated for BP-S, the background color of BP1-S, BP1-M, BP1-D will be displayed in red.

[
Menu >	Data H	leview														JOD .
	Tre	nd T	abular Treod	Recall		CRG	Alarm History	٦.	Wave	Zoom Wave	ST			Full	Disc.	l(t)
			TIETM	L			motory	_ ·	BAIRH	Have					_	
	Expla	nation tre	a) (P)
[24h]	•									× 1			· 01	a F	•	Latest H
E H	_			06/03 23:00			06/10 5:00			06/10 11:00			06/10 1(10st	ວ –	_	
													X12	<u> </u>		
		0\$/10 17:35:00	17:35:10	17:35:20	17:35:30	17:35:40	17:35:50	17:38:00	17:36:10	17:36:20	17:36:30	17:38:40	17:36:50	9/	9] [] []
HR	[ben]	112	113	119	117	125	122	115	119	110	115	113	107	1	1	Alarm Review
Sp02	EX3	98	98	96	95	96	97	96	98	98	98	98	98	1	Π	
NIBP-S	[nmHg]	125	127	124	118	124	120	118	126	125	119	119	128	1		Tabular (Group)
NIBP-D	[nells]	83	79	79	11	84	79	79	80	78	80	80	83	1		(uroup)
RR-IMP	[Bpn]	22	21	19	19	20	22	22	19	22	19	21	20	1		LIST-4
T1	["0]	36.4	36.3	36.5	36.3	36.3	36.3	36.5	36.4	36.3	36.6	36.4	36.4	1		
^a Sv02	[%]	83	83	84	84	83	84	84	83	83				1		
a cco	[]	5.3	5.2	5.4	5.4	5.3	5.4	5.4	5.4	5.3				h		Setup
• 81S		58	58	58	58	58	58	59	58	58	58	58	58	9		0000
														Р		
														1		
	-													1		
L	I													1		
<u> </u>	<u> </u>													1		Print
<u> </u>	<u> </u>		-											1		
<u> </u>	<u> </u>	<u> </u>	-						<u> </u>					1	Ų.	
															9	Print
1	1	1	1	1	1		1	1	1	1	1		1 1	1		LIGHT

On the left side of the parameter, the color assigned for the corresponding parameter will be displayed.

		6	roup		(\mathbf{x})
HR VPC ST(I) ST(I) NBP-5 NBP-0 S002 PR_S002 BP1-5 BP1-0 BP1-9 BP2-S	L151-8 HR VPC ST(1) ST(1) ST(1) ST(1) ST(1) ST(1) ST(4VF) ST(4VF) ST(4VF) ST(4V) ST(V2) ST(V3) ST(V4)	LIST-C HR RR_IMP RR_GAS RR_VENT Sp02 P-PEAK P-PAUSE P-PAUSE P-HAM PEEP E-TV I-TV MV	L1ST-0 Sv02 CC0 EDV B-Tenp RVEF SV CC1 EDV1 ESV SVR Sa02 SVI	LIST-E BIS SQI EWG SR	LIST-F HR SOD-S NIBP-D NIBP-U BPI-U BPI-D BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U BPI-U AGS
					Change Vode Nane

Parameter Setup for Tabular Trend

Press the [Menu], [Tabular Trend] ("Data Review"), [Setup] keys.

• The tabular trend setup screen will be displayed.



 $\mathbf{2}$ Press the [Fix x param.] key to set the fixed parameters.

- Select from [0 param.] to [6 param.].
- ▶ The selected quantity of parameters will be always displayed on the tabular trend, and these data will be remained displayed even when scrolled.

3 Select the display location for the parameter. The selected location will be displayed with blue frame and will be displayed at the side.

- ▶ To change the location, directly press the desired location or drag the 🚺 key up or down.
- ▶ To change the displayed page, press the \blacktriangle / \blacktriangledown keys on the left.

4 Select the parameters.

- 1 Filter the data by sampling interval.
 - [OFF]: The line where [OFF] is selected will not be displayed.
 - [10 sec.]: The displayed data will be filtered in 10 seconds sampling interval.
 - [All]: All data will be displayed.
- 2 Select the category and displaying page.
 - [H Module]/[SvO₂/CCO]/[Vent.]/[Anes.]/[Other]: The parameters for the corresponding category will be displayed.
 - \blacktriangleright \checkmark / \blacksquare : The displaying page for the parameters can be selected.

H Module/Anes.	HR, VPC, ST, SpO ₂ -1, PR_SpO ₂ -1, SpO ₂ -2, PR_SpO ₂ -2, NIBP, BP1 to 8, PR-IBP, PDP, PCWP, CPP, T1 to 8, Tb, CO, EtCO ₂ , InspCO ₂ , RR-GAS, RR-IMP, RR-VENT, APNEA, O ₂ , N ₂ O, Agent, E-TV, I-TV, E-MV, I-MV, P-PEAK, P-PAUSE, PEEP, P-MEAN, RES, COMP, TV 1sec, I/E RATIO, PI, PVI, SpCO, SpMet, SpHb
SvO ₂ /CCO	SvO ₂ , ScvO ₂ , SaO ₂ , O ₂ EI, B-Temp, CCO, CCO-STAT, CCI, CCI-STAT, DO ₂ , RVEF, RVEF- STAT, VO ₂ , SV, SV-STAT, SVI, SVI-STAT, SVR, SVRI, SVV, EDV, EDV-STAT, EDVI, EDVI- STAT, MAP, ESV, ESVI, CFI, iCO, iCI, iSV, iSVI, iSVR, iSVRI, GEDV, GEDI, GEF, EVLW, ELWI, PVPI, ITBV, ITBI, VO ₂ e, VO ₂ I, VO ₂ Ie, iB-Temp, SQI, MAP, CVP, HR, PR, SpO ₂ , iMAP, iCVP, iAvgPR, PO ₂ I, HGB, dPmx, CO CAL

Parameters for each Category

Parameters for each Category

Ventilator	E-TV, I-TV, MV, SMV, P-PEAK, P-PAUSE, PEEP, P-MEAN, E-RES, I-RES, COMP, FiO ₂ , P-MIN, S-COMP, D-COMP, S-RR, I/E RATIO, RES, VTCO ₂ , etCO ₂ , VCO ₂ , Flowee, Ti, Ti/Ttot, PEEPtot, Elastance, Cdyn, D-Chara, Leakage, S-Mve//Mve, Tc, WOBvent, WOBpat, CPAP, P0.1, Edipeak, Edmin, SBI, VT/PBW
Anesthesia Delivery System	Flowee, Ti, Ti/Ttot, Sup.Air, SupO ₂ , SupN ₂ O
Other	BIS, SQI, EMG, SR, SEF, TOTPOW, IMP, Lt-rSO ₂ , Rt-rSO ₂ , S1-rSO ₂ , S2-rSO ₂ , tcpO ₂ , tcpCO ₂

NOTE

- The apnea duration will be stored when it exceeds the upper alarm threshold level. If lower than the alarm threshold level, it will be stored as "0 (zero)".
- The external device parameters to be displayed on the graphic trend/tabular trend needs to be selected in advance on the "Trend Data Setup" window ([Data Review>Graphic Trend or Tabular Trend] or [Initial Settings>External Device>Main Unit/HP-800]).
- The measurement unit of tcpO₂, tcpCO₂ can be set on the TCM4 or TCM5 FLEX. When the measurement unit is changed, the tabular trend data of tcpO₂ and tcpCO₂ on the bedside monitor will be deleted.

REFERENCE

- [H Module] is a generic term for HS-8000, HM-800/HM-801, HP-800, HG-810/HG-820.
- When the FLOW-i Anesthesia Delivery System is connected, the display will change to [H Module/Anes.].

5 Select the parameter to be displayed for the selected location.

▶ The blue frame will move to one row below.

O Set the time bar.

Select the time bar display interval from [4h]/[8h]/[12h]/[16h]/[20h]/[24h]/[36h]/[48h].

Recall

This section explains about the recall function and the setup procedure.

To Display the Recall Waveform

- 1 Date/Time at Alarm Occurrence
- 2 Recall Factor
- 3 Recall Waveform (Compressed: 12 sec.)



When the alarm for the specified recall factor occurs, waveforms (max. 2 waveforms/12 seconds) and numeric data for each recall factor will be stored up to 200 data. On the display selection menu, the data to be displayed can be selected from the stored recall data. 5 compressed recall waveforms will be displayed. Pressing the waveform area will display the enlarged waveform.

If the recall data exceeds 200, the data will be erased from the oldest one.

^{4 🔶} Mark

The recall waveform will be acquired from the point prior to alarm occurrence so that alarm-generated point will be displayed at 7 to 8 seconds point on the 12-seconds recall waveform. I mark indicates the alarm generated point.

7 Press the [Menu], [Recall] ("Data Review") keys. Or, press the [Recall] key on the user key area.

- ▶ Recall screen will be displayed.
- ▶ 5 compressed waveforms (12 sec. per each waveform) will be displayed.
- > The alarm occurrence time, the recall factor occurred at the same time, and the compressed waveform of recall waveform 1 will be displayed.



 $oldsymbol{Z}$ Changing the displayed time, scrolling the time, updating the data (@"Common Operation" P8-1)

3 Press the [Display Selection] key, and set the recall display.

1 Select the quantity of waveforms to be displayed.

- 2 Select the recall factor.
 - > The key will turn blue when pressed to indicate that it is selected as the recall factor.
 - ▶ [Select All]: All parameters including arrhythmia will be selected.
 - ▶ [Select All Arrhythmia]: All arrhythmia will be selected.
 - ▶ [Cancel All]: All selections will be canceled.

NOTE

The "Display Selection" setting will be also applied to the recall list display on the numeric data display area.

(@"Extended Function (Recall List)" P3-14)

4 Set the storing condition for recall data. (recall Setup" P8-19)



5 Deleting All Recall Waveform

- 1 Press the [Delete Sel.] key.
- 2 Select the parameters to delete. For the selected parameter, "x" will be displayed. To select all displayed waveforms, press the [Select All] key.

To cancel the selection, select again the parameter with "x" mark. "x" mark will be cleared indicating that



it has been removed from the deleting parameter selection.

3 Press [Delete]>[Delete OK] keys to delete the parameters with "x" mark.

To Display/Print the Enlarged Recall Waveform

On the enlarged recall waveform display, the recall waveform will be displayed in 25mm/s and by using the cursor, the data before and after the alarm occurrence can be checked.

1 Press the waveform display area on the recall screen.



> The enlarged recall waveform will be displayed.



1 Shifts the recall waveform display.

2 Recall Waveform

The waveform can be dragged to left and right.

3 Printing the Recall Waveform

The displayed enlarged waveform and numeric data will be printed. The output printer can be selected on the "Manual Printing" setup.

(@"Printing Setup" P9-1)

4 Deleting the Recall Waveform

The displayed recall waveform will be deleted.

5 The waveform can be scrolled by dragging the waveform area to left and right.

Recall Setup

The storing condition at alarm occurrence can be set for the recall function. The recall waveform and recall factor (numeric data, arrhythmia) can be selected.

1 Press the [Setup] key on the recall screen.

(@"To Display the Recall Waveform" P8-16)

> The "Setup" window will be displayed.



2 Select the time bar display interval from [4h]/[8h]/[12h]/[16h]/[20h]/[24h]/ [36h]/[48h].

3 Select the recall waveform. Up to 2 waveforms can be selected for recall waveform. (shown on right)



4 Select the recall factor.

(To Display the Recall Waveform" P8-16)

NOTE

• The recall waveform will start with the following delay time tracing back from the alarm occurrence.

	Adult	Child	Neo	nate
	Addit	Onid	Numeric Data Alarm	Arrhythmia Alarm
Delay Time	12 sec.	12 sec.	8 sec.	12 sec.

• For the parameters measured on the multigas unit, the delay time is 8 seconds.

OCRG

This section explains about the OCRG display.

On the OCRG display, compressed respiration waveform, HR trend and SpO_2 trend are displayed simultaneously. The trend scale is fixed as follows.

- HR: 0 bpm to 300 bpm
- SpO₂: 70%SpO₂ to 100%SpO₂

1 Press the [Menu], [OCRG] ("Data Review") keys.

OCRG screen will be displayed.



2 Display Duration

Select from [8min]/[16min].

3 Time Bar

Select the time bar display interval from [4h]/[8h]/[12h]/[16h]/[20h]/[24h]/ [36h]/[48h].

Respiration Waveform

Select from [Impedance]/[CO2]].

5 Respiration Waveform Size



Setup

Tine Bar

 \times

24h

Select the waveform size for compressed respiration waveform.

Respiration Waveform	Size, Scale	
Impedance RESP	[x1/4]/[x1/2]/[x1]/[x2]/[x4]	
CO ₂	[50]/[100] (unit: mmHg)	
002	[4]/[8]/[10] (unit: % or kPa)	

6Printing

The currently displayed trend and compressed waveform on the OCRG screen will be printed.

Alarm History

This section explains the alarm history function and printing procedure.

The alarm generation of numeric data, arrhythmia, equipment status and change in alarm settings can be stored as alarm history. Maximum of 1599 data can be stored.

NOTE

The alarm history cannot be deleted manually. When 1600 data is exceeded, the data will be deleted from the oldest one.

Alarm History Setup

1 Press the [Menu], [Alarm History] ("Data Review") keys.

> The alarm history screen will be displayed.



 $\mathbf 2$ Changing the displayed time, scrolling the time, updating the data (@"Common Operation" P8-1)

 $\mathbf{3}$ Set the alarm history display.

- 1 Select the time bar display interval from [4h]/[8h]/ [12h]/[16h]/[20h]/[24h]/[36h]/[48h].
- 2 Select the alarm level to be displayed. The selected item will be displayed in blue.
- 3 Select the alarm type to be displayed. The selected item will be displayed in blue.



> The currently displayed alarm history will be printed.



5 The displayed data can be scrolled by dragging the display area up and down.

Description for Each Item

The descriptions of each item are as follows.

Item	Details
Time	The alarm generated time or alarm setting changed time will be displayed.
Code	The code related to alarm generation or alarm setting change will be displayed in hexadecimal.
Factor	The factor for alarm generation and alarm setting change will be displayed.
	In case of numeric data/arrhythmia alarm, the numeric data and alarm setting at alarm generation will be also displayed.
	In case of equipment status alarm, a detailed code may be also displayed.
	In case of alarm setting change, the changed value will be also displayed.
Duration (sec.)	The duration of numeric data/arrhythmia/equipment status alarm generation, alarm suspend, monitor suspend, night mode will be displayed in seconds. The maximum displayable value is 99999 sec. It will not be displayed for the alarm setting change.

Print Output Example

BED-013 2011/06/16 20	:47 FUKUDA DENSHI ID:12841	SEX:	AGE : 39	ADULT		ALARM HISTORY 1/2
TIME CODE 11/06/16 20:46:49 2091 11/06/16 20:46:05 4001 11/06/16 20:46:05 4001 11/06/16 20:46:05 3001 11/06/16 20:46:05 3001 11/06/16 20:46:05 3202 11/06/16 20:46:05 3202 11/06/16 20:46:05 3002 11/06/16 20:46:05 3002 11/06/16 20:46:05 3002 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/16 20:46:05 3007 11/06/1	1 Printer Busy 1 Printer Busy 2 Tachy Setting Changed 3 RR (GAS) Lower Limit Changed 2 RR (VENT) Lower Limit Changed 3 RR (IMP) Lower Limit Changed 3 RR (GAS) Upper Limit Changed 2 RR (VENT) Upper Limit Changed 5 RP (IMP) Upper Limit Changed 6 RP (IMP) Upper Limit Changed 7 RT (PCT) Limit Changed 7 RT (PCT) Limit Changed 7 RT (PCT) Limit Changed 7 RT (IMP) Upper Limit Changed		119 120 5 5 30 30 15 30 120		DURA. 5 5	N
BED-013 2011/06/16 20:	:47 FUKUDA DENSHI ID:12841	SEX:	AGE : 3 9	ADULT		ALARM HISTORY 2/2
TIME CODE 11/06/16 20:45:15 3A00			190		DURA.	
11/06/16 20:45:15 3001 11/06/16 20:45:12 080(11/06/16 20:45:12 080(11/06/16 20:45:12 0001 11/06/16 20:45:09 3A0(1 HR Upper Limit Changed D. TACHY 1. Upper HR		190 60 > 60 > 50	50 50	3 3	H
OP050-0	DITD LOT No. 4920 🐺 FUKUDA DENSHI CO.,					

Zoom Wave

This section explains about the "Zoom Wave" window. (When using the optional CFast card) Maximum of 6 waveforms (9.8 seconds each) can be displayed.

The "Zoom Wave" window can be also displayed by pressing the waveform area on the "Full Disc. Wave" window. If the optional CFast card is not used, the latest enlarged recall waveform will be displayed.

1 Press the [Menu], [Zoom Wave] ("Waveform Review") key.

▶ The "Zoom Wave" window will be displayed.



 $\mathbf{2}$ The waveform of previous/next alarm event will be displayed.

Changing the displayed time, scrolling the time, updating the data (PC Common Operation P8-1)

4 The numeric data of the displayed time will be displayed.

5 Switch the waveform to display.

[Limb]: Limb lead ECG waveform will be displayed. [Chest]: Chest lead ECG waveform will be displayed. [User Selection]: The waveform selected at procedure 6 will be displayed.

6 When [User Selection] is set at procedure 5, select the waveforms to be displayed.

The size/scale of the displayed waveform will change.

The currently displayed waveform will be printed.

On the HR-800, 12 seconds of waveform will be printed. The printing range starts from 1 second before the left end of the enlarged waveform.

On the laser printer, 10 seconds of waveform will be printed. The printing range starts from the left end of the enlarged waveform.

9 The waveform can be scrolled by dragging the waveform area to left and right.

ST Measurement

This section explains about the ST measurement and ST alarm function.

To Display/Print the ST Measurement

On the ST display, ECG for the selected time duration (10 sec./1 min./5 min./10 min.) will be displayed overlapped in 1 block.

If 3-lead cable is used, maximum of 8 hours of ST waveform will be displayed.

- NOTE
- If 3-lead cable is used, the measurement will be performed for only the displayed leads.
- · For the following case, ST level will not be displayed.
 - · When learning arrhythmia.
 - · When the lead is off.
 - · When the reference waveform is not set.
 - When "N" or "S" is not detected for QRS within 30 seconds.

Press the [Menu], [ST] ("Waveform Review") key. Or, press the [ST] key on the user key area.

ST screen will be displayed.



2 Changing the displayed time, scrolling the time, updating the data (

3 Changing the Displayed Waveform Size

Select from [x1/4]/[x1/2]/[x1]/[x2]/[x4]. The same waveform size will be applied to all the leads. The selected size will not be applied to the ECG waveform on the home display.

Changing the Displayed Block Duration

The "Setup" window will be displayed and "Slide Show" (1 sec./5 sec./10 sec./20 sec./30 sec.) can be selected.

REFERENCE

- When 3-lead cable is used, 36 blocks will be displayed. When 4, 5, 10-lead cable is used, 3 blocks for each lead will be displayed.
- The duration of each block can be selected from [10 sec.]/[1 min.]/[5 min.]/[10 min.]. For the selections other than [10 sec.], the overlapped waveform for the selected

duration will be displayed.

5 Printing

The currently displayed ST waveform will be printed.

Reference Waveform Setup

The ST reference waveform will be automatically set after learning the arrhythmia. The reference waveform can be updated manually.

1 Press the [Menu], [ST] ("Alarm") key.

> The ST alarm setup screen will be displayed.



 $\mathbf{2}$ Update the ST reference waveform. Press the [Update Ref. Wave] key.

- For the lead which the electrode is detached, the reference waveform cannot be set. Check if the electrode is correctly attached, and perform the setup again.
- 16 beats average of the ECG judged as normal QRS by arrhythmia analysis will be set as the reference waveform.
- > While updating the reference waveform, the [Update Ref. Wave] key will be displayed in blue.
- The updated time of the reference waveform will be displayed.

NOTE

- While learning arrhythmia, or if VPC is present, it will take more than 16 beats to set the reference waveform.
- When the electrode quantity is changed, the reference waveform will be automatically updated.
- In case such as when the patient is discharged, the reference waveform will be automatically set.

3 Set the reference point and measurement point.

- 1 Slide the reference point to right and left using the i key.
- 2 Slide the measurement point to right and left using the $\overbrace{\scriptstyle ext{inst}}$ key.

NOTE

- Set the reference point in the range of –240 ms to 0 ms in increments of 10 ms from the peak of QRS to the P wave direction.
- Set the measurement point in the range of 0 ms to 560 ms in increments of 10 ms from the peak of QRS to the T wave direction.

ST Alarm Setup

Set the ST upper limit and lower limit for the reference waveform.

1 Press the [Menu], [ST] ("Alarm") key.

> The ST alarm setup screen will be displayed.



2 Select [ON]/[OFF] for "ST All Alarm" .

▶ [OFF]: Alarms will not generate even if the alarm for each lead is set to ON.

3 Select the lead to set the alarm limit.

> The selected lead will be displayed large at the right.

4 Select [ON]/[OFF] of ST alarm for each lead.

5 Set the upper and lower alarm limit.

(Alarm Limit Setup for Each Parameter" P6-10)

NOTE

- Set the upper limit in the range from -18 mm to +20 mm/-1.8 mV to +2.0 mV. If a value above +20 mm/+2.0 mV is set, the upper alarm will turn OFF.
- Set the lower limit in the range from 20mm to +18mm/ 2.0mV to +1.8mV. If a value below -20mm/-2.0mV is set, the lower alarm will turn OFF.

REFERENCE

• The upper and lower limit can be set in 1 mm/0.1 mV increments.

12-Lead Analysis

This section explains about the 12-lead analysis function. By using the 10-electrode cable, 12-lead ECG can be displayed, analyzed, stored, and printed. Maximum of 10 analyzed results can be saved.

WARNING

- The 12-Lead ECG analysis function is designed to acquire and interpret ECG data from a resting, supine patient. If ECG signals from moving or shaking patients are acquired, erroneous 12-lead interpretation result may occur. Always ensure that the patient is kept motionless during 12-lead ECG signal acquisition and analysis.
- The 12-lead ECG analysis function is intended for use with adult and pediatric patients.
- All computerized ECG analysis results should be reviewed by a physician before making decision of the patient treatment.

- Interpretation and Minnesota codes given by this equipment do not instruct the physician as to the kind and degree of cardiac disease. Accordingly, four value judgments are given for ECG waveform and although "abnormal" indicates a large possibility of organic cardiac disease, there are cases where no cardiac disease exists despite an abnormal ECG (that is, an abnormal ECG may be caused by something other than heart). On the other hand, care should be taken in the event any preclinical coronary arteriosclerosis could be present despite a normal ECG interpretation. Therefore, for a proper diagnosis, the ECG should be integrated with other interpretations.
- ECG Recording by the Mason-Likar System
 The 12-lead ECG recorded with the torso placement of the limb leads (Mason-Likar 12-lead system) may differ from the standard 12-lead ECG. Moreover, waveforms may differ somewhat also in a supine position and a standing position (sitting position).

 We recommend to carry out the recording of the ECG by taking into consideration the waveform differences according to electrode positions or postures.
- About the ECG Analysis Program

The ECG analysis program is intended to analyze the standard 12-lead ECG waveforms. Therefore, the analyzed result for waveforms recorded with the torso placement of the limb leads (Mason-Likar 12-lead system) may differ from that of the standard 12-lead ECG waveforms.

- When a pacemaker is used, select [Used] for "Pacemaker" under "Admit/Discharge" menu.
- The threshold values for classification of 12-lead ECG interpretation and Minnesota code are set by age and sex as follows.
 - 1. Male and Female of ages 19 years old and above
 - 2. Male of age 12 through 18 years old
 - 3. Female of age 12 through 18 years old
 - 4. Male and Female of ages 3 through 11 years old
 - 5. Male and Female of ages below 2 years old
- If no patient information (i.e. Default: "Class.": [Adult], "Sex": undetermined) has been entered, the system algorithm will handle the patient as a "35 years old male".
- Before the analysis, make sure the patient classification ([Adult] / [Child]) is properly selected. If [Neonate] is selected, the 12-lead ECG analysis will not function.
- Enter the age of patient if known. If no age information (i.e. Default: [0]) has been entered, the system algorithm will handle the patient as "35 years old".
- Enter the sex of patient if known. If no sex information (i.e. Default: undetermined) has been entered, the system algorithm will handle the patient as "Male".
- If the patient classification is set as [Child] and no age (i.e. Default: [0]) has been entered,

the system algorithm will handle the patient as "less than 2 years old."

NOTE

Electrode Placement for 12-Lead ECG Analysis
 When acquiring 12-lead ECG signals, it is recommended to place the limb electrodes anywhere along the arms and legs. (P "Electrode Placement" P7-2)
 If it is difficult, use the Mason-Likar 12-lead system. To reduce the waveform differences from the standard 12-lead, it is recommended that the torso placement of the RA and LA electrodes be near as possible to each arm, in the infraclavicular fossae, within the area unaffected by myoelectricity.

12-Lead ECG Display

Press the [Menu], [12-Lead] ("Waveform Review") key.

▶ The 12-lead screen will be displayed.



2 Changing the displayed time, scrolling the time, updating the data (PCommon Operation P8-1)

 $\mathbf{3}$ The real-time waveforms are displayed.

The 12-lead analysis will be performed based on the displayed waveforms.

REFERENCE

 A pacemaker pulse will not be displayed on the 12-lead analysis screen even if [ON] is set for "Pacemaker Pulse".

4 The [Chest Lead]/[Limb Lead] keys will switch the display between chest lead and limb lead.

5_{Setup}

The 12-lead waveform size, filter, analysis method can be set.

(12-Lead Analysis Setup" P8-29)

5 Printing

- The currently displayed waveform can be printed.
- The output printer will be according to the setting made for "12-Lead Waveform" ([Bedside]/[Laser]) under [Manual Printing>Printer Sel. (Graphic Printing)].
 (@"Manual Printing (Other Setup)" P9-6)

12-Lead Analysis Setup

Press the [Menu], [12-Lead] ("Waveform Review"), [Setup] key.

• The 12-lead analysis setup screen will be displayed.

ZECG Analysis

The timing to read the waveform for ECG analysis can be set.

- [Real Time]: The waveform of 10 seconds after the [Start Analyze] key is pressed will be analyzed.
- ▶ [Review]: The waveform of 10 seconds before the [Start Analyze] key is pressed will be analyzed.

3 Waveform Size

The waveform size for the real-time waveform displayed on the 12-lead screen can be set.

- Limb Lead: The waveform size for the limb lead can be changed.
- Chest Lead: The waveform size for the chest lead can be changed.

4 Filter

The setup for the AC Filter, EMG Filter, Drift Filter can be performed.

- AC Filter: If AC noise is present, select [ON]/ [OFF] for "AC Filter". If [ON] is selected, cut-off frequency will be 75 Hz.
- ▶ EMG Filter: If EMG noise is present, select [Strong (25Hz)]/ [Weak (35Hz)]/ [OFF].
- > Drift Filter: If base line drift is present, select [Strong (0.50Hz)]/ [Weak (0.25Hz)]/ [OFF].

- A baseline or notch will be generated on the ECG waveform (display, print, recall) during the filter setting (up to about 2.4 seconds).
- This equipment complies to the distortion test of IEC 60601-2-25 when all the filters are set to OFF. The frequency characteristic is 0.05 Hz to 150 Hz when all the filters are set to OFF.]
- When a pacemaker is used, the baseline fluctuation becomes large and it may be difficult to read the electrocardiogram. Set the drift filter by checking the electrocardiogram.

Background Color

The background color for the 12-lead display can be set.

- [White]: Similar display with the electrocardiograph. Background Color: White Grid Color: Orange Waveform Color: Black (Fixed)
- [Black]: Standard color
 Background Color: Black
 Grid Color: Gray
 Waveform Color: Green (Fixed)



6 Time Bar

Select the time bar display interval from [4h]/[8h]/[12h]/[16h]/[20h]/[24h]/[36h]/[48h].

12-Lead ECG Analysis

Press the [Menu], [12-Lead] ("Waveform Review"), [Start Analyze] key.

- When the analysis completes, the analyzed result will be displayed. For the analyzed result, dominant waveform and analyzed result will be displayed.
- Abnormal region will be indicated by highlight display.



2 Analyzed Time

The analyzed time will be displayed. During the analysis, [Start Analyze] key will change to [In Progress]. The analysis can be suspended by pressing the [In Progress] key.

3 Dominant Waveform

- ► The reference waveform used for the analysis will be displayed. The dominant waveform is the waveform at the point of ♥ mark on the rhythm waveform.
- On the analyzed result, the abnormal lead with the highest grade finding will be highlighted in red.
- ▶ The dominant waveform display can be switched by pressing the [Chest Lead]/[Limb Lead] keys.

4 Filter Information

The filter used for analysis will be displayed. The filter display can be selected from frequency or type (AC, MF_ST, etc.). (PMaintenance Manual "Display/Print Setup" P5-13)



▶ For the analyzed result, overall judgment, numeric data, finding will be displayed.



- 1 Overall Judgment: The highest grade judgment will be displayed.
- 2 Numeric Data: Main numeric data used for ECG analysis will be displayed. The abnormal numeric data with the highest grade finding will be highlighted in red.
- **3** Finding: The findings by the ECG analysis will be displayed. These will be classified by colors according to the grade specified for each finding.

Grade 6: Red

Grade 4: Blue

Grade 2, 0: Black

The highest grade finding will be highlighted in color specified for each abnormality level.

6 Panorama Display

▶ By pressing the [Panorama] key, overall judgment, finding, abnormal site will be indicated by heart illustration.



- 1 Overall Judgment: The highest grade judgment will be displayed.
- 2 Finding: The ECG analysis finding of highest grade will be displayed.
- 3 Abnormal Site: The finding indicated at 2 will be displayed by a heart illustration.
- During the panorama display, [Panorama] key will change to [Numeric].
 By pressing the [Numeric] key, the analyzed result display will change to numeric data display.

Analyze Real Time Waveform

Press the [Real Time] key to return to the 12-lead analyzed result screen. Press the [Start Analyze] key on the 12-lead analyzed result screen.

B Display Analyzed Waveform

Press the [Analyzed Wave] key to display the analyzed waveform.

- ▶ [Chest Lead]: Chest lead (V1 to V6) waveform will be displayed.
- ▶ [Limb Lead]: Limb lead (I to aVF) waveform will be displayed."

9 Deleting the Analyzed Result

- > Press the [Delete] key to delete the displayed analyzed result.
- [Delete OK] will delete the displayed analyzed result data.
- Press [Cancel] to cancel the delete process.

12-Lead Analyzed Result Output Example

Press the [Print] key on the analyzed result screen or analyzed waveform screen. There are following 2 types of analyzed result printing.

Displayed key when [Print] key is pressed	Printer Selection for Graphic Printing		Key Display	Remarks
Waveform Report	12-Lead Waveform	Bedside	Yes	Standard 12-lead waveform printing
		Laser	Yes	Prints the analyzed waveform.
Panorama Report	12-Lead Analysis Result	lt Bedside No		Panorama Report
		Laser	Yes	Displayed only when [Laser] is set as the printer for graphic printing.
Analyzed Report	12-Lead Analysis Result	Bedside	Yes	Standard analyzed result printing
		Laser	Yes	Prints the waveform and analyzed result. [6Wavesx2 (1 page)]/[6Wavesx2 (2 pages)]/[3Wavesx4+Rhy. (1 page)] can be selected as print format.

NOTE

• If no patient information has been entered, "Adult", "35 years old", and "Male" will be printed.

- If the patient classification is set as "Child", and no age and sex information have been entered, "Child", "2 years old", and "Male" will be printed.
- The output printer will be according to the setting made for "12L Analysis Result" ([Bedside]/ [Laser]). [Manual Printing>Printer Sel. (Graphic Printing)]
 (P⁻ "Manual Printing (Other Setup)" P9-6)

Printed Data

The following basic data will be printed.

HR	Heart rate obtained by basic arrhythmia measurement
QRS Interval	QRS interval of basic waveform measurement. Average value of measurements of leads I to V6. The equipotential part (I wave) at the beginning of QRS and the equipotential part (K wave) at the end of QRS are not included in QRS interval.
R-R Interval	R-R interval of basic waveform measurement. Average value calculated from all the heartbeats first, and then recalculated from the R-R interval within ±25% of that value.
P-R Interval	P-R interval of basic waveform measurement. Average value of measurements of leads I to V6.
QT Interval	QT interval of basic arrhythmia measurement. Average value of measurements of leads I to V6.
QTc Interval	QTc interval of basic arrhythmia measurement. This value is calculated from the following equation: $QTc = \frac{Average waveform QT time}{\sqrt{Average R-R time of arrhythmia (sec.)}}$
QRS Axis	QRS axis of basic arrhythmia measurement. This value is calculated from the following equation, where, I, II, and III are the sums of the maximum amplitude values (signed) of Q, R, S, R', and S' waves from each lead. $Axis (°) = Tan^{-1} \left(\frac{\sqrt{3}(II + III)}{2I + II - III} \right)$
R V5/V6	Maximum amplitude of R wave or R' wave of lead V5 or lead V6. Lead V5 > Lead V6: RV5 Lead V5 = Lead V6: RV6
SV1	Maximum (absolute) value of Q, S, or S' wave of lead V1.
R+S	Sum of the amplitudes of RV5/RV6 and SV1.
ST	Amplitude from the baseline. Measurement position: End of QRS complex + (QT/10) sec.

Printing on the Bedside Monitor Printer

- ▶ When [Beside] is set for "12L Analysis Result" under [Manual Printing>Printer Sel.("Graphic Printing")], pressing the [Print] key will display the [Waveform Report]/[Analyzed Report] keys.
- > The following is the output example when [Analyzed Report] key is pressed.

BED-000 2011/03/04 20 Perro	D:10 ID: SEX:M A	GE:35 ADULT							
	l_		_l_	l_	l_	-1	l	l_	l
						In	JA		
↓ <u> </u>									
25mm/s WANU	JAL 12LEAD ANALYSIS	REC		022004000000			FILTER	0.05-150Hz	PACE OFF
				8ED-000 Demo	2011/03/04 20 D AT: 1	:10 ID: SEX:	M AGE:35 ADULT		
Lave x1		$\sim\sim\sim$	\sim			$\rightarrow \sim$		$\neg \neg \neg$	
I	/		/			/	/	·····	^ i
				. 1 .				. 1 .	. /
•	DS-8500-V03-01 (4	0220) -52	(COHERENT)	25	=/s MANU/	AL 12LEAD ANALY	SIS REC		02200400
									BED D e
									De
		~	~~ <u>}</u> ~			$-\gamma$		$-\gamma$	
1 V2 x1		$-\gamma$				γ			/r
~ V3 x1 ~			_h	_h		-1		_h	
00004009991A		FILTER: 0. 05	-150Hz	PACE OFF	ĥ	DS-8500-V03-	01 (#0220) -52	[COHERENT]	· · · · · · · · · · · · · · · · · · ·
000 2011/03/04 20:10	ID: SEX:M AGE:3	5 ADULT							
VBEAT: 1 V4 x1		-l-				l_			
			1						
		-dr-							
- V6 x1		-l-	-1	-nh-	li	-			
25mm/s MANUAL 1	2LEAD ANALYSIS REC			02200400000004009	991A		FILTER:0.05	-150Hz F	ACE OFF
860-840 2011/06/23 11:3 1D: Si	N:JN EX:N AGE:15 ADULT	101:Within Norma	[Normal range ECI 11 Limits	0	HR : P-R : QRS : QT : QT : QT : QT : QT : QT : QT : QT	605m \$1 999m \$1 176m \$1		[Wisnesota Code] 9-4-1	
Height B. Ocn						609911 51 999611 51 809611 51 809611 51 80550 51 30559 51	(4VH) O. Om (4VL) O. 2M (4VF) O. OM		
Height B.Ocm Weight D.Obg BNI B.Obg/m ² HR SObpas					AXIS: RYS: SY1: R+S:	44 deg. 51 44 deg. 51 1. 18 m¥ 51 0. 59 m¥ 51 1. 77 m¥ 51	(I) O. Om (II) O. Om (III) O. Om (IVR) O. Om (VI) O. Om	Connest:	
					5¥1 = R+S =	0. 59n¥ \$1 1. 77n¥ \$1	(¥5) O. Onn (¥6) O. Onn		
DT No. 4620 🐨 FUKUDA 🍤 👯	\$1474501 (1922) - SZ						-XUnconfirmed repor	IX Doct OP050-01TD	LOT No. 4920

> Pressing the [Waveform Report] key will print the analyzed waveform in a standard format.

Laser Printer Output

- When [Laser] is set for "12-Lead Analysis Result" under [Manual Printing>Printer Sel.("Graphic Printing")], pressing the [Print] key will display [Waveform Report]/[Analyzed Report]/[Panorama Report] keys.
- Pressing the [Analyzed Report] will print the analyzed result in a format set for "12-Lead Analysis Format" under [Menu>Setup>Manual Printing]. ("Manual Printing (12-Lead)" P9-4)



- Pressing the [Waveform Report] will print the waveform in a format set for "12-Lead Waveform Format" under [Menu>Setup>Manual Printing]. ([P"Manual Printing (12-Lead)" P9-4)
- > The following is the output example when [Panorama Report] key is pressed.



NOTE

 To print out the 12-lead analysis panorama report in color, use a laser printer with the page description language in LIPS IV. If a printer with other page description language is used, the printout will be in black and white.

Full Disclosure Waveform (Optional Function)

By using the optional CFast card, 240 hours of full disclosure waveform data can be saved.

Up to six waveforms can be displayed. The alarm event and time will be also saved which allows to search the waveform by each factor.

- Use only the specified CFast card.
- Turn OFF the power before removing the CFast card.
- Check that the CFast card indicator is not lit in orange when turning OFF the power.When using the CFast card, use the standby switch.
- The CFast card can be used only on the equipment where it was formatted.
- The CFast card formatted for the central monitor full disclosure waveform data cannot be used on the DS-8400 System.

NOTE

- When the full disclosure waveform data exceeds the capacity of the CFast card, the data will be deleted from the old one.
- To delete the full disclosure waveform data, perform the discharge procedure.
 (P5-6)

To Format the CF Card

REFERENCE

•To save the full disclosure waveform, the CFast card needs to be formatted for the full disclosure waveform. (Advantage Maintenance Manual "Formatting the Full Disclosure Waveform Card" P3-4)

Waveform Setup

The displaying/printing waveform quantity and type of storing waveform, display duration (sec.) per line for the full disclosure waveform can be preprogrammed.

- **1** Press the [Menu], [Full Disc.] ("Waveform Review"), [Setup] key.
 - > The "Setup" window for full disclosure waveform will be displayed.



2 Select the quantity of displaying/printing waveforms from [1]/[2]/[3]/[4]/[5]/[6].

NOTE
The maximum waveform quantity that can be printed differs depending on the output printer.

3 Press the key for "Waveform". (shown on right)

4 Set the display duration per line from [10 sec.]/[30 sec.]/[1 min.].

- **5** Select the slide show interval from [1 sec.]/[5 sec.]/[10 sec.]/ [30 sec.].
- **6** Select the time bar display interval from [4h]/[8h]/[12h]/[16h]/[20h]/ [24h]/[36h]/[48h].



Description of the Full Disclosure Waveform Display



▶ The "Time Search" window will be displayed.

2 Enter the searching date/time using the numeric keys and press the [Search] key.

- Searching will start.
- The searched waveform will be displayed on the full disclosure waveform display.



Hemodynamics

This section explains the procedure for hemodynamics calculation and printing.

- NOTE
- When the DS-5700 system central monitor is connected via DS-LAN II network, the system will function as follows.
 - If the equipment is connected to DS-LAN, and [ON] is selected for "Synchronize Hemodynamic Data with the Central Monitor", the latest 5 hemodynamic data will be enabled and other hemodynamic data will be deleted. For the latest 5 data, the edited hemodynamic data will synchronize between this monitor and the central monitor.
 - If the equipment is connected to DS-LAN, and [OFF] is selected for "Synchronize Hemodynamic Data with the Central Monitor", the latest 5 data will be transmitted to the central monitor, but the data will not synchronize between this monitor and the central monitor. The data edited on the central monitor will be deleted. The data edited on this monitor will be transmitted to the central monitor.

Calculation Data

Data	Item	Formula
BSA	Body Surface Area (m ²)	h ^{0.725} xw ^{0.425} x71.84x10 ⁻⁴ (Dubois Formula)
СІ	Cardiac Index (L/min/m ²)	CO BSA
SV	Stroke Volume (mL/beat)	CO x 1000 HR
SVI	Stroke Volume Index (mL/beat/m ²)	SV BSA
SVR	Systemic Vascular Resistance (dynes·sec·cm ⁻⁵)	(MAP - CVP) x 79.90 CO
SVRI	Systemic Vascular Resistance Index (dynes·sec·cm ⁻⁵ •m ²)	SVRxBSA
PVR	Pulmonary Vascular Resistance (dyn·sec·cm ⁻⁵)	(MPAP-PCWP)x79.90 CO
PVRI	Pulmonary Vascular Resistance Index (dyn·sec·cm ⁻⁵ •m ²)	PVRxBSA
LVW	Left Ventricular Work (kg·m)	COx(MAP-PCWP)x0.0136
LVWI	Left Ventricular Work Index (kg·m ²)	LVW BSA
LVSW	Left Ventricular Stroke Work (g·m)	SVx(MAP-PCWP)x0.0136
LVSWI	Left Ventricular Stroke Work Index (g·m/m ²)	LVSW BSA
RVW	Right Ventricular Work (kg·m)	COx(MPAP-CVP)x0.0136
RVWI	Right Ventricular Work Index (kg•m/m ²)	RVW BSA
RVSW	Right Ventricular Stroke Work (g⋅m)	SVx(MPAP-CVP)x0.0136
RVSWI	Right Ventricular Stroke Work Index (g·m/m ²)	RVSW BSA

NOTE

 The blood pressure unit for hemodynamics is "mmHg". If the unit is "kPa" or "cmH₂O", it will be converted to "mmHg" when calculating.

To Display/Print the Hemodynamics Data

10 hemodynamic data can be viewed in list format.

Press the [Menu], [Hemodynamics] ("Calculation") keys.

 The hemodynamics screen will be displayed.



2 [Index Disp] key

The display will alternately switch between "BSA, SV, SVR, PVR, LVW, LVSW, RVW, RVSW" and "CI, SVI, SVRI, PVRI, LVWI, LVSWI, RVWI, RVSWI".

3 [Print] key

The currently displayed hemodynamic data will be printed.

New Input of Hemodynamics Calculation

The hemodynamics calculation can be performed using the newly entered data.

The data can be entered manually using the numeric keys or automatically using the current data.

1 Press the [Menu], [Hemodynamics] ("Calculation"), [New Regist.] keys.

▶ The "Edit" window will be displayed.

			Edi	t			(\mathbf{X})	
	2	000/01/21	15:09	9:01				
Input Data		n] [ks).0 68. (-D PAP-	0 [bpn] 0 60 -S PAP-# s] [nnHs]		ART-S [nnHg] 116 CVP [nnHg] 6	ART-M [nnHs] 92 PCWP [nnHs] 8	Latest Data	-2
BSA [fi]	CI [L/min/ĥ]	SV [nL/beat]	SVR [dyn-sec-cfi]	PVR [dyn•sec•c#]	PVRI [dyn•sec•cn+n]	LVW [kg•m]	Regist	-3
1.78	2.80	83	1374	127	227	5.7		
LVWI [ks·n/n]	LVS# [s•n]	LVS₩I [s•n/m̂]	RVW [ks·n]	RV₩I [ks•n/m̂]	R¥S₩ [s•m]	RVSWI [s•n/ĥ]	Cancel	
3.2	95	53	0.68	0.38	11.3	6.3	Delete	

- The current time will be displayed at the upper area.
- Unmeasured data will be left blank.

 $\mathbf{2}$ Enter the calculation data.

- Press the [Latest Data] key to display the measured data.
- Press the key for the editing data to display the numeric keys. Edit the data using the numeric keys, and press the [Set] key.

> The edited data will be displayed in blue.

(NOTE	· · ·
(, ,

• If the height, weight, BSA is changed on the "Admit/Discharge" screen, average CI will be recalculated. However, the hemodynamic result will not be recalculated with the new average CI.

Input	Data
mpm	Duiu

Item (Unit)	Editing Range
Height (cm)	0 cm to 300 cm
Weight (kg)	0 kg to 350 kg
Body Surface Area (m ²)	0 m to 9.99 m ²
Cardiac Output (L/min)	0.00 L/min to 20.00 L/min
Heart Rate (bpm)	0 bpm to 350 bpm
Systolic Arterial Pressure(mmHg / kPa)	0 mmHg to 350 mmHg / 0 kPa to 46.6 kPa
Mean Arterial Pressure (mmHg / kPa)	0 mmHg to 350 mmHg / 0 kPa to 46.6 kPa
Diastolic Arterial Pressure(mmHg / kPa)	0 mmHg to 350 mmHg / 0 kPa to 46.6 kPa
Systolic Pulmonary Artery Pressure (mmHg / kPa)	0 mmHg to 100 mmHg / 0 kPa to 13.3 kPa
Mean Pulmonary Artery Pressure (mmHg / kPa)	0 mmHg to 100 mmHg / 0 kPa to 13.3 kPa
Diastolic Pulmonary Artery Pressure (mmHg / kPa)	0 mmHg to 100 mmHg / 0 kPa to 13.3 kPa
Central Venous Pressure (mmHg / kPa)	0 mmHg to 100 mmHg / 0 kPa to 13.3 kPa
Pulmonary Capillary Wedge Pressure (mmHg / kPa)	0 mmHg to 100 mmHg / 0 kPa to 13.3 kPa
	Height (cm) Weight (kg) Body Surface Area (m ²) Cardiac Output (L/min) Heart Rate (bpm) Systolic Arterial Pressure (mmHg / kPa) Mean Arterial Pressure (mmHg / kPa) Diastolic Arterial Pressure (mmHg / kPa) Systolic Pulmonary Artery Pressure (mmHg / kPa) Mean Pulmonary Artery Pressure (mmHg / kPa) Diastolic Pulmonary Artery Pressure (mmHg / kPa) Central Venous Pressure (mmHg / kPa) Central Venous Pressure (mmHg / kPa) Pulmonary Capillary Wedge Pressure

3 Press the [Regist.]/[Cancel] key.

- > [Regist.]: The calculation will be performed using the newly entered data, and the entered data and calculation result will be registered on the list.
- [Cancel]: The entered data will be deleted.

REFERENCE

- The calculation result will not be displayed if sufficient data is not entered.
- Maximum of 10 data can be registered. If exceeded, the oldest data will be deleted.
- The edited data will be also displayed in blue on the list.

To Edit the Hemodynamics Input Data

The entered data which has been already calculated can be edited or deleted.

1 Press the [Menu], [Hemodynamics] ("Calculation"), and then the date/time display area for the data to edit.

- The "Edit" window will be displayed.
- **2** Edit the data.

(@"New Input of Hemodynamics Calculation" P8-40)

Register the edited data. (@"New Input of Hemodynamics Calculation" P8-40)



4 Delete the data.

A confirmation message will be displayed. To delete the data, press the [OK] key.

Lung Function

This section explains the procedure for lung function calculation and printing.

Calculation Data

Data	Item	Formula
BSA	Body Surface Area (m ²)	h ^{0.725} xw ^{0.425} x71.84x10 ⁻⁴
CaO ₂	Arterial Oxygen Content (mL/dL)	CaO ₂ =1.34xHbxSaO ₂ +0.003xPaO ₂
$C\bar{v}O_2$	Mixed Venous Oxygen Content (mL/dL)	$C\bar{v}O_2$ =1.34xHbxS $\bar{v}O_2$ +0.003xP $\bar{v}O_2$
a-vDO ₂	Arteriovenous Oxygen Content Difference (vol %)	a-vDO ₂ =CaO ₂ -Cv̄O ₂
DO ₂	Oxygen Transport(mL/min)	DO ₂ =CaO ₂ xCOx10
DO ₂ I	Oxygen Transport Index(mL/min/m ²)	DO ₂ I=CaO ₂ xClx10
ΫO ₂	Oxygen Consumption(mL/min)	VO₂=a-vDO₂xCOx10
ΫO ₂ Ι	Oxygen Consumption Index(mL/min/m ²)	VO₂I=a-vDO₂xCIx10
O ₂ ER	Oxygen Extraction Rate (%)	$O_2 ER = (CaO_2 - C\bar{v}O_2)/CaO_2 x100$
		AaDO ₂ =P _A O ₂ -PaO ₂
AaDO ₂	Alveolar-Arterial Oxygen Difference (Torr)	$P_AO_2=P_1O_2-(P_ACO_2/R)x(1-F_1O_2x(1-R))$ R:Respiration Quotient (0.8 for this equipment) $P_1O_2=(P_B-47)xF_1O_2$
ḋ _s ∕ḋ _t	Shunt Rate (%)	Q _s /Q _t =(CćO ₂ -CaO ₂)/(CćO ₂ -C⊽O ₂) CćO ₂ =1.34xHb+0.003xP _A O ₂

REFERENCE

• The blood pressure unit for lung function calculation is "mmHg". If the unit is other than

"mmHg", it will be converted to "mmHg" when calculating.

To Display/Print the Lung Function Data

256 lung function data can be viewed in list format.

1 Press the [Menu], [Lung Function] ("Calculation") keys.

▶ The lung function list will be displayed.

2 [Index Disp] key

The display of BSA, CaO_2 , CvO_2 , $a-vDO_2$, DO_2 , VO_2 , O_2ER , $AaDO_2$, Qs/Qt will alternately switch with that of CI, DO_2I , VO_2I .

5 [Print] key

The currently displayed lung function data will be printed.



New Input of Lung Function Calculation

The lung function calculation can be performed using the newly entered data.

The data can be entered manually using the numeric keys or automatically using the current data.

Press the [Menu], [Lung Function] ("Calculation"), [New Regist.] keys.

▶ The "Edit" window will be displayed.



 $\mathbf{2}$ Enter the calculation data.

- ▶ Press the [Latest Data] key to display the entered data of "HEIGHT", "WEIGHT", "CO".
- Press the key for the editing data to display the numeric keys. Edit the data using the numeric keys, and press the [Set] key.
- The edited data will be displayed in blue.

NOTE

 If the height, weight, BSA is changed on the "Admit/Discharge" screen, average CI will be recalculated. However, the lung function calculation result will not be recalculated with

the new average CI.

Input Data

1	
Data	Item (Unit)
HEIGHT	Height (cm)
WEIGHT	Weight (kg)
BSA	Body Surface Area (m ²)
СО	Cardiac Output (L/min)
FiO ₂	Fraction of Inspiratory Oxygen (%)
P _B	Atmospheric Pressure (mmHg)
PaCO ₂	Partial Pressure of Arterial Carbon Dioxide (mmHg)
Hb	Hemoglobin Concentration (g/dL)
PaO ₂	Partial Pressure of Arterial Oxygen (mmHg)
SaO ₂	Arterial Oxygen Saturation (%)
$P_{\bar{V}}O_2$	Partial Pressure of Mixed Venous Oxygen (mmHg)
S _V O ₂	Mixed Venous Oxygen Saturation (%)

3 Press the [Regist.]/[Cancel] key.

- [Regist.]: The calculation will be performed using the newly entered data, and the entered data and calculation result will be registered on the list.
- [Cancel]: The entered data will be deleted.

REFERENCE

- The calculation result will not be displayed if sufficient data is not entered.
- Maximum of 256 data can be registered. If exceeded, the oldest data will be deleted.
- · The edited data will be also displayed in blue on the list.

To Edit the Lung Function Input Data

The entered data which has been already calculated can be edited or deleted.

Press the [Menu], [Lung Function] ("Calculation"), and then the date/time display area for the data to edit.

▶ The "Edit" window will be displayed.

```
2 Edit the data.
```

(☞ "New Input of Lung Function Calculation" P8-43)

Register the lung function list. (☞ "New Input of Lung Function Calculation" P8-43)



4 Delete the data.

(P8-43) (

Cardiac Output (CO)

This section explains about the cardiac output measurement using the thermodilution method, setup procedure for catheter type, etc., and procedure for editing the measurement result.

To Display the CO Measurement Screen

Press the [Menu], [CO] ("Calculation") keys.

Or, press the [CO] key on the user key area.

- > The CO measurement screen will be displayed.
- ➤ The message according to the status will be displayed, and if "READY" is displayed, the measurement can be started.

(Cardiac Output Message" P11-21)



The Description of the CO Measurement Screen

- 1 Result Status
- 2 Thermodilution Curve
- 3 Time Scale
- 4 Cardiac Index (CI)
- 5 Blood Temperature
- 6 Injectate Temperature
- 7 Status Message
- 8 Cardiac Output (CO)



Cardiac Output Setup

Before measuring the cardiac output, set the measurement condition such as ON/OFF of auto start, time scale for thermodilution curve, injection condition, etc.

1 Press the [Menu], [CO] ("Calculation"), [Setup] keys.

▶ The "Setup" window will be displayed.



2 Set ON/OFF of "Auto Start" .

- 1 Press the key for "Auto Start".
 - The dropdown list will be displayed.
- 2 Select from [ON] or [OFF].
 - [ON]: The measurement will automatically start when the injectate is injected.
 - [OFF]: The measurement will start by pressing the [Start] key.

REFERENCE

Even when [ON] is selected, the measurement can be manually started by pressing the [Start] key.

3 Set the time scale.

1 Press the key for "Time Scale".

- ▶ The dropdown list will be displayed.
- 2 Select from [30 sec.]/[60 sec.].

4 Set the computation constant.

- 1 Press the key for "CC".
 - The dropdown list will be displayed.
- 2 Select from [Auto Input]/[Manual Input].
 - [Auto Input]: The computation constant will be automatically set according to the catheter size and the injection volume.
 - [Manual Input]: The computation constant for the used catheter can be manually input with the numeric keys.

Auto Input of CC Value



1 Select the catheter manufacturer from [BIOSENS]/[ARGON]/[EDWARDS].

REFERENCE

- ARGON: Argon Medical Devices Japan, K.K. (formerly Becton, Dickinson and Company)
- The manufacturer name can be changed on "Catheter Manufacturer for CC Input" setting (Menu>Initial Settings>Meas.>Other).

2 Select the "Catheter Size (F) from [5]/ [6]/ [7]/ [7.5].

3 Select the "Injectate Volume (mL) " from [3]/ [5]/ [10].

> When the above items are selected, the computation constant will be automatically set.

When the CJ0-P01C-C2.4 Catheter Relay Cable is used:

1 Select the "Injectate Temperature" from [Ice]/ [Room].

- ▶ [Ice]: The measurement will be performed at 0°C.
- ▶ [Room]: The measurement will be performed at room temperature (24°C).

4 Press the [Set]/[Cancel] key.

• [Set]: CC value will be finalized.

NOTE

- If the CC value does not correspond to the used catheter, or to use the previous CC value, press the [Cancel] key, and enter the value manually.
- To automatically enter the computation constant, the catheter relay cable needs to be connected.

Manual Input of CC Value



1 Select the catheter manufacturer from [BIOSENS]/[ARGON]/[EDWARDS].

 $\mathbf{2}$ Up to 3 types of CC value can be programmed for each manufacturer.

- ▶ If previously entered value is present, press the key for "History".
- ▶ If the previously entered value is not present, enter the CC value using the numeric keys.
- **3** Set the "Injectate Temperature".

(alue" P8-47)

4 Press the [Set]/[Cancel] key.

• [Set]: CC value will be finalized.

CO Measurement

1 Press the [Menu], [CO] ("Calculation") keys.

• The CO measurement screen will be displayed.



▶ The displayed message will change from "WAIT" to "READY".

NOTE

• While "WAIT" is displayed, the measurement cannot be started. Wait until "READY" is
displayed.

 ${f Z}$ Verify that "READY" is displayed, and press the [Start] key.

Pressing the key will generate a sound.

3 Inject as soon as the sound generates.

> When the measurement is complete, CO and CI value will be displayed.

REFERENCE

• If "Auto Start" is ON, the measurement will automatically start at injection by detecting the blood temperature.

4 Press the [Print] key.

> The displayed thermodilution curve, CO, CI value will be printed.

NOTE

- When "WAIT" message is continuously displayed, verify that catheter relay cable is properly connected to the cardiac output module, and thermodilution catheter is securely connected.
- Before injecting, check that the Ti (injectate temperature) setting is correct.
- When repeatedly performing the measurement, inject at intervals of 30–60 seconds
- The CI value will not be displayed unless height/weight or BSA value is input on the "Admit/Discharge" screen.
 (
 "Entering the Patient Information" P5-1)
- · For the following cases, measurements may be inaccurate.
 - Shunt disease, tricuspid regurgitation or pulmonic regurgitation.
 - During exercise stress
 As body temperature varies non-continuously and unevenly by exercise, constant CO value cannot be measured.
 - Excessive Arrhythmia As blood volume varies non-continuously due to arrhythmia, accurate CO value cannot be measured.

To Edit the CO Measurement Result

The average CO and average CI can be calculated by performing the CO measurement continuously and editing the measurement result.

Press the [Menu], [CO] ("Calculation") keys.

> The CO measurement screen will be displayed.

• The average CO and average CI value obtained from the measurement result will be displayed.



2 To Change the Selected Status

The selected data for the average value will be displayed in blue. Press the graph area to change the selected status.

V Mark: VPC detected during CO measurement.

*: CO value exceeding the average CO value $\pm 10\%$.

3 [Average CO Input]

The displayed average CO value will be entered to the list.

(NOTE

• If the height, weight, BSA is changed on the "Admit/Discharge" screen, average CI will be recalculated.

As the CI will not be recalculated after the hemodynamic calculation, save the average CI by hemodynamic calculation before changing the height, weight, and BSA.

4 [Delete Sel.] ([Delete])

The [Delete Sel.] key will change to [Delete] key, and the data can be deleted.

x mark will be displayed for the data to be deleted, and pressing the [Delete OK] key will delete the data.

Drug Calculation

This section explains about the drug calculation function.

The drug calculation function is a function to calculate the flow rate of drug administration to the patient. Based on the dosing rate, flow rate and dosing duration will be calculated from the weight, drug amount, diluent amount.

It is also possible to calculate the dosing rate and dosing duration from the flow rate.

REFERENCE

 Under the "Initial Settings", the drug name, and default settings for each drug (drug amount/ unit, diluent amount, dosing rate/unit) can be set. (PMaintenance Manual "Drug Calculation" P5-27)

Press the [Menu], [Drug Calc.] ("Calculation") keys.

▶ The drug calculation menu will be displayed.



▶ If the weight is entered on the "Admit/Discharge" menu, the entered weight will be displayed.



 $\mathbf{2}$ Press the key for "Drug Name".

	Drug Selection	
AMRINONE	AHINOPHYLLINE	BRETYLIUM
DOBUTAMINE	DOPAMINE	EPINEPHRINE
HEPARIN	INSULIN	ISOPROTERENOL
LIDOCAINE	NITROGLYCERIN	NITROPRUSSIDE
NOREPINEPHRINE	PHENYLEPHRINE	PROCAINAWIDE
STREPTOKINASE	(tPA	DRUG A
DRUG B	DRUG C	DRUG D
DRUG E	DRUG F	DRUG G

- > The list of registered drugs will be displayed.Select the drug to administer to the patient.
- When a drug is selected, the drug amount, diluent amount, dosing rate/unit preset for that drug under "Initial Settings" will be automatically entered.

NOTE

• The flow rate will be automatically calculated when the value for each item is updated.

• On the initial display of the drug calculation menu, the previous calculation data will be displayed. The calculation data will be cleared when the patient is discharged.

3 Enter the value for each item.

- To change the automatically entered value, press the key for each item and manually enter the value.
- ▶ The dosing rate and flow rate can be adjusted by pressing the [+], [-] keys.

NOTE

• If the selected unit for the dosing rate requires weight, the flow rate cannot be calculated if the weight is not entered.

4 Press the [Update "End by"] key to update the estimated time of completion.

• Pressing the [Update "End by"] key will recalculate the time from the pressed time and update the estimated time of completion.

Calculation Formula for Flow Rate/Dosing Rate

According to the dosing rate unit, the calculation formula from the following 10 types will be automatically selected for calculation.

Dosing Rate Unit	Flow Rate Calculation Formula
mg/min	Flow Rate (mL/hr) = Doing Rate (mg/min) x Diluent Amount (mL) x 60 Drug Amount (mg)
mg/hr	Flow Rate (mL/hr) = Dosing Rate (mg/hr) x Diluent Amount (mL) Drug Amount (mg)
mg/kg/min	Flow Rate (mL/hr) = Drug Amount (mg) Drug Amount (mg)
mg/kg/hr	Flow Rate (mL/hr) = Drug Amount (mg) Drug Amount (mg)
µg/min	Flow Rate (mL/hr) = Drug Amount (mg) x 1000 Drug Amount (mg) x 1000
µg/hr	Flow Rate (mL/hr) = Drug Amount (mg) x 1000
µg/kg/min	Flow Rate (mL/hr) = Doing Rate (µg/kg/min) x Weight (kg) x Diluent Amount (mL) x 60 Drug Amount (mg) x 1000
µg/kg/hr	Flow Rate (mL/hr) = Drug Amount (mg) x 1000 Drug Amount (mg) x 1000
units/hr	Flow Rate (mL/hr) = Drug Amount (units) Drug Amount (units)
IU/hr	Flow Rate (mL/hr) = Dosing Rate (IU/hr) x Diluent Amount (mL) Drug Amount (IU)
Dosing Rate Unit	Dosing Rate Calculation Formula
mg/min	Dosing Rate (mg/min) = Flow Rate (mL/hr) x Drug Amount (mg) Diluent Amount (mL) x 60
mg/hr	Dosing Rate (mg/hr) = Flow Rate (mL/hr) x Drug Amount (mg) Diluent Amount (mL)
mg/kg/min	Dosing Rate (mg/kg/min) = Flow Rate (mL/hr) x Drug Amount (mg) Diluent Amount (mL) x Weight (kg) x 60
mg/kg/hr	Dosing Rate (mg/kg/hr) = Flow Rate (mL/hr) x Drug Amount (mg) Diluent Amount (mL) x Weight (kg)
µg/min	Dosing Rate (µg/min) = Flow Rate (mL/hr) x Drug Amount (mg) x 1000 Diluent Amount (mL) x 60
µg/hr	Dosing Rate (µg/hr) = Flow Rate (mL/hr) x Drug Amount (mg) x 1000 Diluent Amount (mL)

Dosing Rate Unit	Dosing Rate Calculation Formula
µg/kg/min	Dosing Rate (µg/kg/min) = Flow Rate (mL/hr) x Drug Amount (mg) x 1000 Diluent Amount (mL) x Weight (kg) x 60
µg/kg/hr	Dosing Rate (µg/kg/hr) = Flow Rate (mL/hr) x Drug Amount (mg) x 1000 Diluent Amount (mL) x Weight (kg)
units/hr	Dosing Rate (units/hr) = Flow Rate (mL/hr) x Drug Amount (units) Diluent Amount (mL)
IU/hr	Dosing Rate (IU/hr) = Flow Rate (mL/hr) x Drug Amount (IU) Diluent Amount (mL)

Unit and Setting Range (Dosing Rate, Drug Amount, Diluent Amount, Flow Rate, Weight)

		Do	sing Rate	Drug Amount		
Dr	ug	Unit (Selectable)	Setting Range	Unit	Setting Range	
AMRINONE						
AMINOPHYLLIN	E					
BRETYLIUM DOBUTAMINE DOPAMINE EPINEPHRINE ISOPROTERENOL LIDOCAINE						
		mg/min,				
		mg/hr,				
		mg/kg/min, mg/kg/hr,				
		μg/min, μg/hr,		mg		
NITROGLYCER	IN	µg/kg/min,				
NITROPRUSSIC	θE	µg/kg/hr				
NOREPINEPHRINE PHENYLEPHRINE					0.01 to 1500000.00	
PROCAINAMIDE	PROCAINAMIDE tPA		0.01 to 1500000.00			
tPA						
HEPARIN		unite /h.u.		units		
		units/hr				
STREPTOKINAS	STREPTOKINASE		IU/hr			
DRUG-A to G		mg/min, mg/hr, mg/kg/min, mg/kg/hr, μg/min, μg/hr, μg/kg/min, μg/kg/hr		mg		
		units/hr	+	units	1	
		IU/hr		IU		
Diluent Amount		Fl	ow Rate	Weight		
Unit	Setting Range	Unit	Setting Range	Unit	Setting Range	
mL	1 to 1000	mL/hr	0.1 to 1000.0	kg	0.1 to 449.9	

NOTE

• The setting is not possible if it cannot be correctly calculated by the entered value.

Other Bed Display

This section explains about the function to display the waveform and numeric data and to set alarms for other bedside monitors.

The other bed alarm function generates the alarm sound for the other bed on this monitor. To use this function, wired network (DS-LAN II or DS-LAN III) connection is required.

Δ	
∕!∖	CAUTION

•	On the DS-LANII network system, maximum of 3 monitors (including the central monitor) can display the data of this monitor using the other bed display function. However, there is no restriction of numbers for the DS-7000 series central monitors and DS- 5700. These monitors will be counted as 1 monitor regardless of the numbers. Ex. 1) In case of 1 central monitor and 5 bedside monitors (A to E): The total number of monitors that can display the data of Bedside Monitor A is 3 monitors which consist of 1 central monitor and 2 out of 4 bedside monitors (B to E). Ex. 2) In case of 3 central monitors (DS-7000 series or DS-5700) and 5 bedside monitors (A to E): The total number of monitors that can display the data of Bedside Monitor A is 5 monitors which consist of 3 central monitors and 2 out of 4 bedside monitors (B to E).
٠	If the number of bedside monitors displaying the same bed exceeds the limit, the bedside monitor with smaller ID will be prioritized.
•	If monitoring 12-lead waveform on the central monitor, the total numbers of monitors that can display the same bed will be reduced by 1.
(•	NOTE This equipment cannot connect to a wired network of AU-5500N 8ch Recorder set as the administrator. Even if connected, other bed display, printing and other function cannot be used.

Other Bed Display/Alarm

The other bed display can be accessed from the menu or from the preprogrammed user key.

Also, by setting the other bed alarm to [ON], [Other Alarm] will be displayed when other bedside monitor generates an alarm.By pressing this [Other Alarm] key, the display for the other bed can be accessed.

BED-001 20-05 Tai CH6008 FUKUDA DENSHI - Adult Pacenaker	K INITIA Uther Check Electrodes	Alarm Suspend 🖄 (118sec.)	▲C Pover → 10:08 ☆ 2017/01/19
1 Press the [Menu], [Other B	ed] keys.		
2	LEBOLING AREA		
4	tarm und OFF Atarn Display Setup Area 1 1 1 1 1 1 1 1 1 1 1 1 1		
5	BED-002 BED-003 BED-004 BED-004 BED-006 BED-007 BED-008 BED-003 BED-003		
3	BED-011 BED-012 BED-013 BED-014 BED-014 BED-016 BED-017 BED-018 BED-019 BED-019		
0	BED-021 BED-022 BED-023 BED-024 BED-024 BED-026 BED-027 BED-028 BED-023 BED-023	°	
l	BED-031 BED-032 BED-033 BED-034 BED-034 BED-038 BED-037 BED-038 BED-039 BED-039		

· On the other bed selection menu, select the bed to display from maximum of 100 beds

(in case of DS-LAN III) connected to the wired network. The Room / Bed ID for the alarm generating bed will be displayed in red. The other bed alarm generating bed will be indicated by an icon \bigwedge inside the Room/Bed ID key.

$\mathbf{2}$ Select the area.

- Select the area to be displayed.
 - [All]: The beds for all the area connected to the network will be displayed.
 - [Area 1 to 5]: The beds for each area will be displayed.



3 Press the Room/Bed ID key to display the other bed.

Waveforms and numeric data for the selected bed will be displayed. If an alarm is generated for this bed, the physiological alarm / arrhythmia alarm message will be displayed.



1 Message Area

The message for the other bed will be displayed.

2 Waveform Display Area

Maximum of 6 waveforms for the DS-LAN III network, and maximum of 2 waveforms for the DS-LAN II network can be displayed.

- **3** By pressing the [Other Bed Alarm Silence] key on the other bed display, the alarm sound for the displayed bed can be silenced.
- **4** Pressing this key will switch ON/OFF of menu title display.
- 5 Numeric Data Area

The numeric data at the bottom of the screen can be switched by using the

- 6 Press the [Waveform Selection] key to select the waveforms.
 - Waveform 1 is fixed as ECG, but other waveforms can be selected. Maximum of 6 waveforms for the DS-LAN III network, and maximum of 2 waveforms for the DS-LAN II network can be displayed.
 Select the waveform from the waveform calculation window.
 - Select the waveform from the waveform selection window.
- 7 Press the [Numeric Selection] key to display [Numeric Data Selection] window. The parameters to display

on the right side of the screen can be selected.

4 Set the other bed alarm.

Press the [Alarm Display] key to change the screen to other alarm setup mode. When the mode is changed, the [Alarm Display] key will be displayed in blue. To return to the original mode, press the [Alarm Display] key again.

Select the bed to generate the other bed alarm.

- Select the Room/Bed ID for the bed to generate the alarm. The selected bed will be indicated by blue frame. To cancel the selection, press the key for the bed again.
- [Select All], [Cancel All]: Selection/cancellation for all the beds can be performed at once.
- [Enter]: The selection will be finalized.

5 Turn ON the other bed alarm.

- [ON]: Other bed alarm will be generated.
- [OFF]: Other bed alarm will not be generated.

6 Set the area.

All the beds connected to the network can be displayed, but it is also possible to divide the beds by areas, which allows to display the beds by each area.



- 1 Press the key for "Area Setup" to change the screen to area setup mode. When the mode is changed, the key for selected area will be displayed in blue. To return to the original mode, press the key again.
- 2 Select the Room/Bed ID for the bed to assign to the area. The selected bed will be indicated by blue frame. To cancel the selection, press the key for the bed again.
 - [Select All], [Cancel All]: Selection/cancellation for all the beds can be performed at once.
 - ▶ [Enter]: The selection will be finalized.
- 3 Press the key for "Area Setup" to change the screen to area setup mode.
- 4 Press the [Area Name/Color] key.



1 Select the color to distinguish the area.

A triangle mark with the selected color will be displayed at the corner of the Room/Bed ID key.

- 2 Enter the area name using the numeric keys.
- 3 Maximum of 8 characters can be set for the area name.

Chapter 9 Printing

Printing Setup	
Manual Printing (Basic)	9-1
Manual Printing (12-Lead)	9-4
Manual Printing (Other Setup)	9-6
Automatic Printing (Alarm Printing)	9-6
Automatic Printing (Periodic Printing)	9-8
Common Setup for Printing	9-9
Freeze Printing	9-10
12-lead Waveform Printing	9-10

Chapter 9 Printing

Printing Setup

This section describes the procedure for printing and recording.

For the DS-8400 System, the following type of printing/recording can be performed.

- Manual Printing
- Automatic Printing (Periodic Printing)
- Automatic Printing (Alarm Printing)
- Freeze Printing
- Graphic Printing (Trend, Tabular Trend, Recall, etc.)

REFERENCE

- The printed HR/PR data depends on the ECG/SpO₂/BP selection for "Synchronized Mark/ Tone" under [Menu>Parameter>ECG (SpO₂, BP)]. (PT-9)
- Under the following condition, the amplitude value will be printed for the ECG calibration waveform.

*[Bar (10mm)] is set for "Waveform Size Display" under [Initial Settings>User I/F>Display/ Print].

*[ON] is set for "Print Calibration" under [Manual Printing>Common]

1 Press the [Menu], [Manual Printing] or [Auto Printing] ("Basic Setup") keys.

• The manual printing or automatic printing setup screen will be displayed.

Manual Printing (Basic)

The manual printing can be set to start from the time the key is pressed, or 8 sec./16 sec. prior to the time the key is pressed.

Also, the printing can be set to automatically stop after 24 seconds, or continue to print until the "Print Start/Stop" key is pressed again.

The printer can be selected from bedside monitor printer or central monitor printer.



1 Printer

[Bedside]: Data will be printed on the HR-800 of the bedside monitor. [Central]: Data will be printed on the central monitor printer.

2 Waveform

On the "Select Wave" window, 3 waveforms can be selected for printing. The key for the selected waveform will be displayed in blue.

3 Delay Time

[None]: Printing will start from the point the [Print Start/Stop] key is pressed. [8 sec.] / [16 sec.]: Printing will start 8 sec. or 16 sec. prior from the point the [Print Start/Stop] key is pressed.



- If [None] is selected for the manual printing delay time, QRS classification symbol will not be printed. To print the QRS symbol, set the delay time to [8 sec.] or [16 sec.].
- The HR-800 can be inserted to the internal slot, or connected using the U-LINK cable. Select which HR-800 to use from [Built-in] or [U-LINK] under [Initial Settings > System > Other].

4 Print Duration

[24sec.]: Printing will automatically stop after 24 seconds.

[Cont.]: Printing will continue until the [Print Start/Stop] key is pressed again or until paper runs out.



To Start/Stop the Printing

1 Press the user key or [Print Start/Stop] key on the HR-800.

- Pressing this key during periodic printing, alarm printing, graphic printing, or recall printing will cease the printing in process.
- ▶ Inside the [Print Start/Stop] key, the output printer status for manual printing will be displayed.



Message	Description
None	Normal Operation
PAPER OUT	There is no thermal paper.
CASSETTE	Check the cassette.
CHECK?	Other abnormality is found.

Example of Manual Printing



The 21-digit number printed at the bottom of the paper indicates the settings of the equipment. At the 14th digit from the left, filter setting (AC filter, drift filter) is printed in hexadecimal number.

0	
1	
2	AC Filter ON
3	AC Filter ON
4	
5	
6	AC Filter ON
7	AC Filter ON

8		Drift Filter ON
9		Drift Filter ON
А	AC Filter ON	Drift Filter ON
В	AC Filter ON	Drift Filter ON
С		Drift Filter ON
D		
Е	AC Filter ON	Drift Filter ON
F	AC Filter ON	Drift Filter ON

Filter setting is OFF for the numbers in blank.

Manual Printing (12-Lead)



The monitoring 12-lead waveform can be printed on the bedside monitor printer. The delay time is 6 seconds. The 12-lead waveform cannot be printed on the central monitor printer.

1 Waveform Format

[Regular]: Printing will start from the limb leads. (In the order of I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6) [Reverse]: Printing will start from the chest leads. (In the order of V1, V2, V3, V4, V5, V6, I, II, III, aVR, aVL, aVF)

2 Print Calibration

[ON]: Calibration waveform will be printed. If [Bar (10mm)] is set for "Waveform Size Display" under [Initial Settings>User I/F>Display/Print], the amplitude value corresponding to the displayed waveform size will be printed.

[OFF]: Calibration waveform will not be printed.

3 Position

[Center]: Equalizes the printing width of each lead so that the waveform baseline will be at the center. The printing scale of the waveform will be also automatically adjusted.

[Proportional]: Equalizes the blank space between each lead to avoid overlapping of the waveforms. The printing scale of the waveform will be also automatically adjusted.

[OFF]: Waveform position will not be adjusted when printing.

4 Printer Auto Scale

When position adjustment is [OFF], select whether or not to automatically adjust the scale.

NOTE

• The printer scale will be adjusted in the range of x1, x1/2, x1/4. It will not be adjusted to x2 or x4 even if the amplitude is small.

[ON]: Printing scale will be automatically adjusted.

[OFF]: Printing will be performed with the displayed scale.

5 Lead Boundary

This setting will be displayed only when [Laser] is selected as the printer for "12-Lead Waveform", "12L Analysis Result".

[ON]: Lead boundary between the leads will be printed.

[OFF]: Lead boundary will not be printed.

6 12-Lead Waveform Format

When [Bedside] is set as the printer for "12-Lead Waveform", select from [3Wavesx4]/[2Wavesx6]. When [Laser] is set as the printer for "12-Lead Waveform", select from [3Wavesx4]/[3Wavesx4+Rhy.]/ [6Wavesx2]/[12Waves].

Output Example	Waveform Layout	Length of Each Waveform
3 waves x 4	First column: I,II,III Second column: aVR, aVL, aVF Third column: V1, V2, V3 Fourth column: V4, V5, V6	6 sec.
2 waves x 6	First column: I,II Second column: III, aVR Third column: aVL, aVF Fourth column: V1, V2 Fifth column: V3, V4 Sixth column: V5, V6	6 sec.

7 12-Lead Analysis Format

When [Bedside] is set as the printer for "12L Analysis Result", the format is fixed as [3Wavesx4]. When [Laser] is set as the printer for "12L Analysis Result", select from [6Wavesx2 (2 pages)]/[6Wavesx2 (1 page)]/[3Wavesx4+Rhy.].

Manual Printing (Other Setup)

Select the printer for graphic printing and recall printing.



- 1 Press the key for [Graphic Printing] to display the "Printer Selection" window.
 - [Bedside]: Data will be printed on the HR-800 of the bedside monitor.
 - [Central]: Data will be printed on the central monitor printer.
 - [Laser]: Data will be printed on the laser printer.

		Printer Sel	ection		(X)
Trend	Bedside	Zoon Wave (Recall, Full Disc.) Bedside	Hemodynamics	Bedside
Tabular Trend	Beds i de	ST	Bedside	Lung Function	Beds i de
OCRG	Beds i de	12-Lead Haveform	Bedside	CO	Beds i de
		12L Analysis Result	Bedside		
		FD Compressed Waveform	Bedside		

- NOTE
 - The HR-800 can be inserted to the internal slot, or connected using the U-LINK cable. Select which HR-800 to use from [Built-in] or [U-LINK] under [Initial Settings > System > Other].

REFERENCE

- Graphic printing is a printing performed from the data review screen such as graphic trend and tabular trend.
- To select laser printer, it is necessary to select [ON] or [DS-LAN] for "Network Printer" under [Menu > Initial Settings > External Device > Network] in advance.
 (@Maintenance Manual "Laser Printer Setup" P4-27)

2 Recall Printing

- [Graphic Printing]: Recall data will be output on the printer selected for "Graphic Printing".
- [Manual Printing]: Recall data will be output on the printer selected for "Printer" under "Basic".

Automatic Printing (Alarm Printing)

When numeric data alarm or arrhythmia alarm occurs, printing will automatically start.

(NOTE

- The alarm detection is performed each second, and if more than one alarm occurs at the same time, one data will be stored according to the alarm priority.
- Maximum of 3 alarm data can be stored. If more than 3 alarms generate, the higher priority alarm will replace the previously stored lower priority alarm. The stored data will be deleted once it is printed.
- Priority of alarm printing factor ; ASYSTOLE > VF > VT > Ext Tachy > Ext Brady > SLOW VT > TACHY > BRADY > RUN >

 $\begin{array}{l} \mathsf{HR} \left(\mathsf{HR} / \mathsf{PR_SpO}_2 / \mathsf{PR_IBP}\right) > \mathsf{APNEA} > \mathsf{BP1} \ (\text{or } \mathsf{ART}) > \mathsf{SpO}_2 > \\ \mathsf{NIBP} > \mathsf{RR} \left(\mathsf{RR_IMP} / \mathsf{RR_CO}_2 / \mathsf{RR_GAS} / \mathsf{RR_VENT}\right) > \mathsf{EtCO}_2 > \\ \mathsf{GAS} \left(\mathsf{CO}_2\text{-}\mathsf{E} / \mathsf{CO}_2\text{-}\mathsf{I} / \mathsf{AGT}\text{-}\mathsf{E} / \mathsf{AGT}\text{-}\mathsf{I} / \mathsf{O}_2\text{-}\mathsf{E} / \mathsf{O}_2\text{-}\mathsf{I} / \mathsf{N}_2\mathsf{O}\text{-}\mathsf{I}\right) > \mathsf{MAC} > \mathsf{MV} > \mathsf{PAUSE} > \\ \mathsf{COUPLET} > \mathsf{BIGEMINY} > \mathsf{TRIGEMINY} > \mathsf{FREQUENT} > \mathsf{SVT} > \mathsf{IRREGULAR} \mathsf{RR} > \\ \mathsf{PROLONGED} \mathsf{RR} > \mathsf{SFREQUENT} > \mathsf{SCOUPLET} > \mathsf{VPC} > \mathsf{SVPC} > \mathsf{NOT} \mathsf{CAPTURE} \\ > \mathsf{NOT} \mathsf{PACING} > \mathsf{BP2} > \mathsf{BP3} > \mathsf{BP4} > \mathsf{BP5} > \mathsf{BP6} > \mathsf{BP7} > \mathsf{BP8} > \mathsf{ST} > \mathsf{TEMP} > \mathsf{Tb} > \\ \mathsf{InspCO}_2 > \mathsf{SpCO} > \mathsf{SpMet} > \mathsf{SpHb} > \mathsf{RR_SpO}_2 > \mathsf{PEAK} > \mathsf{PEEP} > \mathsf{BIS} \\ \end{array}$



1 Alarm Printing

[ON]: Printing will automatically start at alarm occurrence.

[OFF]: Printing will not start at alarm occurrence.

Alarm Factor Selection

The "Factor Selection" window will be displayed.

The selected alarm factor key will be displayed in blue. The alarm OFF mark will be displayed inside the key for the parameter in alarm OFF condition.

[Select All Arrhythmia]: All arrhythmia factors will be selected. [All ON]: All alarm factors will be selected.

[All OFF]: All selections for the alarm factor will be cancelled.

	Factor Selection						
Alarm Factor	Select All Arrhythnia	Asystole	VF	ТА	Ext Tachy	Ext Brady	
	٩̈́۴	SLOW VT	Tachy	Brady	Run	Pause	
	Weas.	HR	ST	NIBP	RR	APNEA	
Select All		BP1	BP2	BP3	BP4	BP5	
Cancel All		11	12	T3	14	15	
	•°	ъ	C02	02	N2 0	Agent	

3Printer

[Bedside]: Data will be printed on the HR-800 of the bedside monitor.

[Central]: Data will be printed on the central monitor printer.

4 Print Duration

(@"Manual Printing (Basic)" P9-1)

(NOTE

- The delay time differs depending on the print duration.
- The HR-800 can be inserted to the internal slot, or connected using the U-LINK cable. Select which HR-800 to use from [Built-in] or [U-LINK] under [Initial Settings > System > Other].

			Delay Time			
Print Duration	ation Adult Child Nec		Neona	nate		
Adult	Child	Numeric Data Alarm	Arrhythmia Alarm			

12 sec.	12 sec.	12 sec.	8 sec.	12 sec.		
12 000.	8 sec. for the multigas unit alarm					
24 sec.	16 sec.	16 sec.	16 sec.	16 sec.		

5 Waveform

(@"Manual Printing (Basic)" P9-1)

[Alarm]: Prints the waveform of the alarm factor.

Automatic Printing (Periodic Printing)

The printing will be automatically performed with the selected interval.

NOTE

- If the periodic printing is interrupted due to paper out, etc., the latest periodic printing will be performed when the printing is resumed.
- QRS classification symbol will not be printed for periodic printing.



1 Periodic Printing

[ON]: Printing will automatically start at fixed interval. [OFF]: Turns OFF the periodic printing function.

2 Printer

[Bedside]: Data will be printed on the HR-800 of the bedside monitor. [Central]: Data will be printed on the central monitor printer.

3 Timer/Interval for Periodic Printing

Periodic Interval 🗙						
Tiner	00:00	12:00				
Inter.	02:00	14:00				
	04:00	16:00				
	05:00	17:00				
	07:00	19:00 20:00				
	09:00	21:00 22:00				
	11:00	23:00				



Display Example for "Timer"

Display Example for "Interval"

[Timer]: Printing will automatically start at selected time.

[Interval]: Printing will automatically start at selected interval.

- If [5 min.] is selected for [Interval], the time will be displayed in real time such as 10:00, 10:05, ...10:25. If [60 min.] is selected, it will be displayed as 10:00, 11:00, 12:00.
- 4 Print Duration

The printing will automatically stop after the selected duration.

5 Waveform

(ranual Printing (Basic)" P9-1)

Common Setup for Printing

The printing condition common for manual printing and automatic printing can be set.



Display Example for Automatic Printing

1 QRS Classification

[ON]: QRS classification symbol will be printed with the ECG waveform.

Symbol	Description
N (Normal)	Normal QRS beat
V (VPC)	Ventricular extrasystole
S (SVPC)	Supraventricular extrasystole
P (Pacing Beat)	Pacing beat
F (Fusion Beat)	Fusion beat of pacing and spontaneous beat
? (Undetermined Beat)	Learning arrhythmia, or unmatched beat

[OFF]: QRS classification symbol will not be printed.

- NOTE
 - The QRS symbol cannot be printed for manual printing if the "Delay Time" is set to [None], and for periodic printing. To print the QRS symbol, set the "Delay Time" to [8 sec.] or [16 sec.] for manual printing.
 - The "S" (QRS symbol) will be printed as "N" on the central printer.

2 Printing Speed

[25mm/s]: The printing speed will be set to 25mm/s. [50mm/s]: The printing speed will be set to 50mm/s.

3 Print NIBP Data

[ON]: Oscillation graph and NIBP data will be printed after the waveform.

[OFF]: Oscillation graph and NIBP data will not be printed.

4 Print Calibration

[Top]: Calibration waveform will be printed at the beginning of the waveform.

[Each Page]: Calibration waveform will be printed in 18.75 cm interval.

[OFF]: Calibration waveform will not be printed.

Freeze Printing

The waveform trace can be suspended and printed from 12 seconds prior to the point the waveform trace was stopped.

The waveform selected for manual printing will be printed. The print duration is 12 seconds.

To freeze the waveform display, the [Freeze] key needs to be assigned as user key.

(To Configure the Display" P10-4)

Press the [Freeze] key on the user key.

• The waveform trace will stop.

2 Press the [Print Start/Stop] key.

- The displayed waveform will be printed.
- ► Freeze printing will be output on the bedside monitor printer. The waveforms selected for manual printing will be printed.

12-lead Waveform Printing

When the display layout is "12-Lead", pressing the [12-Lead Print] key will start 12-lead waveform printing.

Select "12-Lead" for the display layout.



2 Press the [12-Lead Print] key.

Printing will start.

• The printing duration of the waveforms for each format are as follows.

	Printing Format	Printing Duration	Delay Time
When printed on the bedside monitor	3 waves x 4	6 sec.	6 sec.
printer:	2 waves x 6	0 500.	0 500.
	3Wavesx4 ^{*1}	2.5 sec.	
When printed on the laser printer:	6Wavesx2*1	5 sec.	10 sec.
when printed on the laser printer.	3 wavesx4+Rhythm*1	12.5 sec.	10 Sec.
	12 Waves*2	10 sec.	

*1 [CONTINUOUS]: The waveform output will be in the time sequence of waveform block order.

*2 [COHERENT]: The waveform output will be in the same time phase for all waveforms.

Chapter 10 System Configuration

Display Configuration	
Numeric Data Selection	
To Configure the Display	10-4
Waveform Selection	10-14
User Key Selection	10-15
Changing the Display Layout from the Home Display	10-19
Tone/Volume	10-20
Color	10-21
Brightness	10-23
Night Mode	10-24
Night Mode	

Chapter 10 System Configuration

Display Configuration

This section describes about the display configuration type and the procedure to configure the display. The monitoring display can be configured according to the monitoring purpose. There are following types of basic display layout.

Standard

- +12-Lead
- Numeric Data/Bottom

When ECG cascade or block cascade is selected, a full disclosure waveform can be displayed. The user keys can be also assigned to the numeric data area.

If extended board (optional) is equipped, up to 2 extended displays can be used. (extended display function)

Display Example



Standard (Box Layout: Right)



12-Lead (Box Layout: Right)



On this system, 9 main modes and 6 sub modes can be preprogrammed according to the monitoring purpose. By preprogramming the configuration to each mode, the display configuration setups at admittance of patient can be simplified by just selecting one of the modes.

("To Select the User Mode" P5-8)

It is recommended to program the mode in rough classification such as patient's condition, monitoring purpose (ICU or surgery), and if necessary, perform unique setup for each patient.

Numeric Data Selection

The numeric data to be displayed can be selected on the "Numeric Data Selection" window.

The parameters of the "Numeric Data Selection" window can be assigned to the numeric data box on the home display.

("Numeric Data Box Display (for each parameter)" P3-8)

Num	eric Data Selecti	on X
OFF	ST-C	Cont inuous Set up
HR/PR	BP1	
HR	BP2	
PR_Sp02	BP3	Set
PR_Sp02N2	BP4	Reselect Area
PR_IBP	BP5	•000
VPC, PACE	BP6	• ►
ST, VPC	BP7	
ST-A	BP8	
ST-B	NIBP	

Example: Page 1

The Numeric Data Box Size for Each Parameter

		Size							
Numeric Data	Width ^{*1}	W1/2		W1			W2		
	Height ^{*2}	H1	H1	H2	H3	H1	H2	H3	
HR/PR		x	0	0	0	0	0	0	
HR		х	0	0	0	0	0	0	
PR_SpO ₂		х	0	0	0	0	0	0	
PR_IBP		0	0	0	0	0	0	0	
VPC, PACE		х	0	0	0	0	0	0	
ST, VPC		х	0	0	0	0	0	0	
ST-A, ST-B, ST-C		х	х	0	0	х	0	0	
BP1 to BP8		х	0	0	0	0	0	0	
NIBP		х	0	0	0	0	0	0	
NIBP List		х	0	0	0	0	0	0	
SpO ₂		х	0	0	0	0	0	0	
SpO ₂ , PR		x	0	0	0	0	0	0	
SpCO		х	0	0	0	0	0	0	
SpMet		x	0	0	0	0	0	0	
SpHb		х	0	0	0	0	0	0	
Sp*		х	0	0	0	0	0	0	
RR_IMP, RR_CO ₂ , RR_VENT, RR_SpO ₂		0	0	0	0	0	0	0	
T1 to T8, Tb		0	0	0	0	0	0	0	
T1/T2, T3/T4, T5/T6, T7/T8		х	0	0	0	0	0	0	
ΔΤΕΜΡ-Α, ΔΤΕΜΡ-Β, ΔΤΕΜΡ-C, ΔΤΕΜΡ-D		0	0	0	0	0	0	0	
VENT		x	х	0	0	х	0	0	

	Size							
Numeric Data	Width ^{*1}	W1/2		W1			W2	
	Height ^{*2}	H1	H1	H2	H3	H1	H2	H3
P-V, F-V		х	х	0	0	х	0	0
SvO ₂ , CO		х	х	0	0	х	0	0
SvO ₂ , CO, CI		x	х	0	0	х	0	0
CO, SV, SVV		х	х	0	0	х	0	0
BIS		x	0	0	0	0	0	0
INVOS		х	0	0	0	0	0	0
CO ₂		х	0	0	0	0	0	0
0 ₂		0	0	0	0	0	0	0
N ₂ O		0	0	0	0	0	0	0
Agent		х	0	0	0	0	0	0
RR, CO ₂ , Agent, O ₂ , N ₂ O		х	х	0	0	х	0	0
CO ₂ , Agent, O ₂ , N ₂ O		x	х	0	0	х	0	0
RR, Agent, O ₂ , N ₂ O		х	х	0	0	х	0	0
Agent, O ₂ , N ₂ O		х	х	0	0	х	0	0
Agent, N ₂ O		x	0	0	0	0	0	0
GAS, SPIRO		x	х	0	0	х	0	0
SPIRO		x	x	0	0	x	0	0
НЕМО		x	x	0	0	x	0	0
HEMO-I		x	x	0	0	x	0	0
STOPWATCH		х	0	0	0	0	0	0
VENT-A		х	0	0	0	0	0	0
VENT-B		х	0	0	0	0	0	0
Hemo/etc-A		х	0	0	0	0	0	0
Hemo/etc-B		х	0	0	0	0	0	0
Extended Function-A		х	x	0	°*3	х	0	°*3

The Numeric Data Box Size for Each Parameter

*1: For LC-8016TC, W1/2 is about 30 mm, W1 is about 60 mm, W2 is about 120 mm For LC-8018TC, W1/2 is about 34mm, W1 is about 69mm, W2 is about 138mm

*2: For LC-8016TC: H1 is about 16 mm, H2 is about 32 mm, H3 is about 48 mm (H1 is the same length as waveform areax2) For LC-8018TC: H1 is about 17mm, H2 is about 36mm, H3 is about 55mm (H1 is the same length as waveform areax2)

*3: For "Extended Function-A", H6 is the maximum height.

To Configure the Display

1 Press the [Menu], [Display Config.] ("Basic Setup") keys.

- The display configuration menu will be displayed.
- 1 Layout (Changing the Layout" P10-4)
- 2 Numeric Data (@ "Changing the Displayed Numeric Data" P10-5)
- 3 Waveform ("Changing the Displayed Waveform" P10-6)
- 4 Sweep Speed (@"Sweep Speed" P10-9)
- 5 Short Trend (Short Trend Display" P10-7)
- 6 Zoom Wave (@"Enlarged Waveform Setup" P10-10)
- 7 User Key (@"User Key Setup" P10-10)
- 8 Detail Setup (@"Detail Setup" P10-11)

Changing the Layout

The layout can be changed with the following procedure.

Press [Change] for "Layout".

▶ The "Layout" window will be displayed.

 $\mathbf{2}$ Select the layout to be displayed.

When "Bottom" is selected, select the number of rows.

 $\mathbf{3}$ Select the user key location from [Right] or [Left], and select the number of columns.

4 Select the numeric data box location from [Right] or [Left].

5 The displayed parameters will be automatically located with the selected layout. Check the home display.

▶ If there are parameters which cannot be displayed due to display area, "Delete Confirmation" window will be displayed. (shown on right)

Pressing the [Set] key will set the layout with some parameters not displayed.

Pressing the [Cancel] key will return to the "Layout" window.





Delete Confirmation						
If changed to the selected layout, the following item(s) cannot be neasured.						
BP3						
BP4						
BP5			Cancel			
			Set			
L	1					

6 If not changing the layout, press the [Cancel] key.

Adjusting the Layout Automatically

The display layout can be automatically adjusted. The automatic mode can be selected from the following two types.

Type-1 (All Auto Mode)

The measured parameters will be automatically located according to the priority. The display layout remains the same. (The layout will change if there is not enough space to display all parameters.)

The display priority can be set on the "Auto Display Configuration" under "Initial Settings". (A Maintenance Manual "Display/Print Setup" P5-13)

•Type-2 (Auto Mode depending on Parameter Quantity)

The parameters will be automatically located according to the parameter quantity using the current display configuration. The display layout, numeric data location and user keys on the numeric data area remain the same.

1 Select [Type-1] or [Type-2].



2 Select [Auto] for "Layout".

NOTE

- For both [Type-1] and [Type-2], the waveform layout is equivalent to that when the [Same with Numeric] key is pressed.
- When [Auto] is selected for the display layout, the following changes are not possible. *Changing the Displayed Waveform
 - *Changing the Displayed Numeric Data
 - *Changing the short trend parameters

Changing the Displayed Numeric Data

The displayed numeric data can be changed with the following procedure.

CAUTION

 When performing the telemetry or wired network transmission, configure the display so that the numeric data corresponding to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted.

NOTE

For HR/PR data, an alarm will be generated only for the current parameter displayed in the HR/PR numeric data box.

The parameter for the HR/PR numeric data box can be selected by pressing the key for "HR/PR" on the ECG, BP, SpO₂ parameter setup window/floating window or by pressing the [HR/PR] user key.

Press the [Change] key for "Numeric Data".

• The display will change to numeric data selection mode.

> The "Numeric Data Selection" window will be displayed.



 $\mathbf{2}$ Press the numeric data display area to change the parameter.

- By pressing the selected area again, the selection will be canceled.
- > To start again from the beginning, press the [Reselect Area] key.
- Adjust the size of the selected area which is indicated by blue box.

3 Select the parameter on the "Numeric Data Selection" window.

Press the [] / [keys to switch the displayed parameters.

(@"Numeric Data Selection" P10-2)

Press [Continuous Setup] to switch to continuous setup mode.

> On the continuous setup mode, the numeric data box area can be sequentially selected.

5 Press the [Set] key.

▶ The setup will be finalized.

NOTE

 The selected parameter may not be displayed depending on the combination of the parameters and size.

Changing the Displayed Waveform

The displayed waveform can be changed with the following procedure.

• When performing the telemetry or wired network transmission, configure the display so that the numeric data corresponding to the waveform is displayed. If not, the displayed waveform or numeric data may not be transmitted.

1

Press [Change] for "Waveform".

• The display will change to waveform selection mode.

> The "Waveform Selection" window will be displayed.



 $\mathbf{2}$ Press the waveform display area to change the parameter.

- By pressing the selected area again, the selection will be canceled.
- ▶ To start again from the beginning, press the [Reselect Area] key.
- Adjust the size of the selected area which is indicated by blue box.
- 3 Select the parameter on the "Waveform Selection" window. Press the ▲ / ▼ keys to switch the displayed parameters. (☞ "Waveform Selection" P10-14)

4 Press [Continuous Setup] to switch to continuous setup mode.

- On the continuous setup mode, the waveform display area can be sequentially selected.
- **5** Press the [Set] key.
 - The setup will be finalized.

Short Trend Display

The parameters and display duration for the short trend display can be set.



NOTE

- The short trend can be displayed when the numeric data layout is "Right"/"Right&Bottom"/ "Left"/"Left&Bottom"/"Bottom".
- For the 12-lead layout, ST value of each lead will be displayed in short trend.

- **1** Press the [Change] key to set the parameters for the short trend display.
 - 1 The parameters for the current waveform display area will be displayed.
 - 2 The selected short trend parameters will be displayed.
 - 3 Select the short trend area, and assign the parameter for that area.



- 4 [Same as Numeric]: The same parameters for the currently displayed numeric data will be set as the short trend parameters.
- 5 [Same as Waveform]: The same parameters for the currently displayed waveform will be set as the short trend parameters.
- 6 Press [Continuous Setup] to switch to continuous setup mode. On the continuous setup mode, the short trend display area can be sequentially selected.

NOTE

- The [Change] key will be displayed when [User Setup] is selected for "Short Trend" (Display Config.>Detail Setup).
- [Same as Numeric], [Same as Waveform] will be applied for the displayed parameters at the point when the key is pressed. The short trend parameters will not automatically change when the displayed parameters are changed.

 $\mathbf{2}$ Select ON/OFF of short trend display.

- [ON]: Short trend will be displayed on the home display.
- ▶ [OFF]: Short trend will not be displayed on the home display.
- [Overlap]: Short trend will be displayed overlapped with the waveform.

3 When [ON] or [Overlap] is selected, select the display duration. The selectable duration differs depending on the short trend data resolution and display width (7 levels).

		Display Width (7 levels)						
		0	1	2	3	4	5	6
Data Resolution	5 sec.	Display OFF	5 min.	10 min.	15 min.	20 min.	25 min.	30 min.
	10 sec.	Display OFF	10 min.	20 min.	30 min.	40 min.	50 min.	60 min.
	30 sec.	Display OFF	30 min.	60 min.	90 min.	120 min.	150 min.	180 min.

4 Select the display duration for the short trend.

1 Press the waveform display area on the home display.



2 The trend display time will change to the time of the pressed position.



- When an alarm is generated for the recall alarm factor, recall screen will be displayed.
- When the cursor function is enabled, a cursor will be displayed. The display duration can be changed under "Short Trend".(Display Config.>Detail Setup)

Sweep Speed

1

The sweep speed can be set with the following procedure. The sweep speed can be set differently for the circulatory system waveforms (ECG, BP) and respiratory system waveforms.

Select the circulatory sweep speed from [6.25]/[12.5]/[25]/[50] (mm/s).

 ${f 2}$ Select the respiratory sweep speed from [6.25]/[12.5]/[25] (mm/s).

Enlarged Waveform Setup

By selecting [ON] for "Zoom", the displayed waveform size and sweep speed will be doubled.

Menu > Basic Setup > Display Config.	IND-MIT FUKUDA DENSHI - MMIT R [MITTAL]	45 Power → 15:32 ■ 2017/01/19
Layout Change Numeric Data Change Size Change		NO Hone
Auto Type-1 Vaveforn Change Same as Moneric	96 116 77	Aların Si Lence
Background DFF Speed Circ. Vent. mm/s] 25 6.25		Auto Node Alarm Auto Node Alarm V ett./mini
Palette Vivid Snort Change OFF 15 min.		6.20 MEAN 21 NIBP List Recall
Detail Setup		3.5 1.2 21 16 30 0 ► Fried Start St
	[™] ₩ № a	Key Lock Night Wode

- NOTE
- When the sweep speed is set to [50 mm/s], "Zoom" cannot be set to [ON].
- Scale will not be enlarged.

User Key Setup

The user key can be set with the following procedure.

1 Press the [Change] key for "User Key".

- The display will change to user key selection mode.
- > The "User Key Selection" window will be displayed.



 $\mathbf{2}$ Select the area to change the user key.

- ▶ By pressing the selected area again, the selection will be canceled.
- > To start again from the beginning, press the [Reselect Area] key.
- Adjust the size of the selected area which is indicated by blue box.

3 Select the function to assign to the user key on the "User Key Selection" window.

- NOTE
 - The displayed user key can be switched between 2 displays using the [User Key Up] and [User Key Down] keys.
 - Press the ▲/ ▼ keys to switch the user key selection. (P10-15)
4 Press [Continuous Setup] to switch to continuous setup mode.

On the continuous setup mode, the user key display area can be sequentially selected.

5 Press the [Set] key.

• The setup will be finalized.

Detail Setup

- **1** Press the key for "Detail Setup".
 - The "Detail Setup" window will be displayed.



1 Alarm Limit Display

The alarm limit can be displayed inside the numeric data box.

[Graph]: Alarm limit will be displayed in bar graph. [Numeric]: Alarm limit will be displayed in numeric format. [OFF]: Alarm limit will not be displayed.

2 At Alarm Occurrence

The numeric data display format at alarm occurrence can be selected.

[Reversed]: The numeric data will be displayed in reversed color at alarm occurrence. [3D]: The numeric data will be displayed in 3D at alarm occurrence.

3 Grid

The ECG waveform can be displayed on the grid.

[ON]: Grid will be displayed.

[Bold]: Grid will be displayed in bold format.

[OFF]: Grid will not be displayed.

REFERENCE

Short trend and grid cannot be displayed overlapped.

4 Scale

The scale can be selected from [ON]/[Bold1]/[Bold2].

5 Thickness

The thickness of the displayed waveforms can be selected from [Thin] / [Regular] / [Thick].

6 Clip

Whether or not to clip the overlapped waveforms of the neighboring display area can be selected.

7 Fill CO₂ Waveform

Whether or not to fill in the CO_2 waveform from the baseline can be selected.

8 Fill O₂ Waveform

Whether or not to fill in the O₂ waveform from the baseline can be selected.

9 Fill Agent Waveform

Whether or not to fill in the Agent waveform from the baseline can be selected.

10 BP Overlap

The overlapping BP waveforms can be set for each overlap group 1 to 3.

11 RR Overlap

The overlapping RR waveforms can be set.

12 12-Lead ST Wave

The ST waveform to be displayed for the 12-Lead layout can be set. [Ref.]: The ST reference waveform will be displayed. [Average]: The average waveform will be displayed.

13 12-Lead ST Short Trend

The display format for the ST short trend can be selected from [Plot]/[Fill]/[OFF].

14 ST/VPC/Arrhy. Alarm Display

Whether or not to display the ST value, VPC (integrated value of 1 minute), arrhythmia alarm message inside the HR numeric data box can be selected.

15 Block Cascade

The waveform combination for block cascade display can be set.

16 Graphic/Tabular Trend

Graphic trend and tabular trend will be displayed in the waveform display area.

NOTE

• When [ON] is selected for "Graphic/Tabular Trend", the waveform set to the same display area with the graphic/tabular trend will not be displayed.

17 Graphic/Tabular Trend Size

Select the display area size for graphic/tabular trend from [Big]/[Medium]/[Small].

18 Short Trend

The short trend parameters can be linked to the displayed numeric data or waveform.

[Link with Numeric]: The short trend layout will be linked to the displayed numeric data on the home display. [Link with Waveform]: The short trend layout will be linked to the displayed waveform on the home display. [User Setup]: User settings will be applied for the short trend layout.

19 Short Trend Scale

The short trend scale for the following parameters can be synchronized with the scale of trend or waveform. BP / PEAK / TV / CO_2 / O_2 / Agent

20 Display Parameter

Whether or not to display the parameter name of the displayed short trend can be set. [ON]: Displays the parameter name with the corresponding color of the parameter. [Gray]: Displays the parameter name in gray. [OFF]: Parameter name will not be displayed.

21 Reference Line Function

Whether or not to display the reference lines can be set for the following parameters.

HR, ST, BP1 to 4, NIBP, EtCO₂, SpO₂, BIS

[Enable]: The reference line function will be enabled. On the "Short Trend Setup" window (displayed when short trend scale area is pressed), ON/OFF of reference line display and reference line position can be set for each parameter.

[Disable]: The reference line function will be disabled.

The reference line function cannot be used for the overlapped short trend display.

• When [Enable] is selected, the function to highlight the alarm generated data cannot be used.

22 Cursor Function

Whether or not to display a cursor can be selected. By displaying a cursor, the measured data and review data (tabular trend/graphic trend/zoom wave) at the time of cursor position can be displayed. [Enable]: The cursor function will be enabled. However, the function to enlarge/reduce the display duration by pressing the short trend area will be disabled.

[Disable]: The cursor function will be disabled.

NOTE

NOTE

- The cursor function cannot be used for the overlapped short trend display.
- When [Enable] is selected, the function to highlight the alarm generated data cannot be used.
- The cursor will be displayed when the short trend area is pressed, and will be automatically cleared after a short while.

23 Cursor Linkage

When [Enable] is selected for "Cursor Function", the review data to be displayed can be selected from [Tabular Trend] / [Graphic Trend] / [Zoom Wave].

The zoom wave can be displayed only when the full disclosure waveform function is enabled.

24 Short Trend Overlap

Maximum of 4 parameters can be displayed overlapped in the same short trend area.

However 2 blocks of waveform area are required for each parameter. For example, to display 3 parameters in the same short trend area, 6 blocks of waveform area are required.
 ShortTrend OverLap
 Status
 X

 ShortTrend OverLap1
 ØFF
 ØFF
 ØFF

 ShortTrend OverLap2
 ØFF
 ØFF
 ØFF

 ShortTrend OverLap3
 ØFF
 ØFF
 ØFF

25 Data Resolution, Display Duration

Select the data resolution from [5 sec.]/[10 sec.]/[30 sec.]. The display duration will differ depending on the "Data Resolution" setting.

For [5sec.], maximum display duration is 30 minutes.

- For [10sec.], maximum display duration is 1 hour.
- For [30sec.], maximum display duration is 3 hours.

 ${f Z}$ Press the [Home] key to check the configured display.

NOTE

• If the numeric data box is configured at the bottom of display, user keys cannot be assigned to the numeric data box area.

- After configuring the display, make sure to verify the configured display by pressing the [Home] key.
- To maintain the configured display even after the power is turned OFF or after the discharge procedure, store the configuration to one of the user modes, or select [Backup] for "Display Configuration" under Initial Settings>User I/F>At Power ON/At Discharge.
 (P"To Select the User Mode" P5-8)

Waveform Selection

The waveform to be displayed can be selected on the "Waveform Selection" window. This section explains the details of the displayed waveforms.



1 ECG1 to ECG12

The ECG waveform of the specified channel will be displayed. Minimum of 2 blocks are required to display the ECG waveform.

2 ECG1 to ECG12 Cascade

The ECG waveforms of the specified channel will be displayed in cascade. Minimum of 2 blocks are required to display in cascade.

3 BP Overlap 1 to 3

The BP waveform (BP1 to BP8) set on "BP Overlap Setup" will be displayed. If the waveform display area is too small to display the assigned BP waveforms, it will be displayed in the priority from smaller channel numbers.

4 RR Overlap 1 to 3

The RR waveform (CO₂, O₂, Agent) set on "RR Overlap Setup" will be displayed. If the waveform display area is too small to display the assigned waveforms, it will be displayed in the priority of CO₂>O₂>Agent.

5 Block Cascade

The waveforms (2 to 6) set on the "Block Cascade Setup" will be displayed in one block.

Other than the waveforms explained above, the selected waveform on the "Waveform Selection Window" will be displayed.

User Key Selection

The user keys can be set on the "User Key Selection" window. This section explains the function for each user key.

U	ser Key Selection	
OFF	Print Start/Stop	Cont inuous Setup
Home	Monitor Suspend	
Wenu	Night Wode	
Minimize	Freeze	Set
Restore #indow	Key Lock	Reselect Area
User Key 🌲	Mode Select	•0000
Alarm Silence	Admit/Disch.	• ►
Alarm Suspend	Rapid Discharge	
NIBP Start/Stop	HR/PR	User Key 🕇
NIBP Cont.	HR/PR Source	User Key ↓

Example: Page 1

Page 1	
OFF	Blank key will be displayed.
Home	The display will return to the home display. The [Home] key is also available as fixed key.
Menu	The menu screen will be displayed. The [Menu] key is also available as fixed key.
Minimize Window	Pressing this key will minimize the currently displayed window and will be stored to the user key.
Restore Window	Pressing this key will restore the minimized window.
User Key	The first and second page of the user key area will switch. This key will be located at the same position for both first and second page.
Alarm Silence	Alarm sound will be suspended for fixed amount of time.The [Alarm Silence] key is also available as fixed key. By pressing the key for more than 3 seconds while the alarm is not generated, it will bring the system to "Alarm Sound Suspend" condition.
Alarm Suspend	Alarm (sound and display) will be suspended for fixed amount of time.
NIBP Start/Stop	NIBP measurement will start/stop.
NIBP Cont.	NIBP continuous measurement will start/stop.
Print Start/Stop	Manual printing will start/stop.
Monitor Suspend	Confirmation window to suspend monitoring will be displayed.
Night Mode	Night mode will turn ON/OFF.
Freeze	Waveform trace will freeze for fixed amount of time. Pressing the [Print Start/Stop] key while in freeze condition will print the frozen waveform. Pressing the key again will start the waveform trace again.
Key Lock	Touch key operation will turn ON/OFF. It can be used when cleaning the display panel.
Mode Selection	User mode selection screen will be displayed.
Oxygenator Mode	The home display will switch to oxygenator mode.
Admit/Discharge	Admit/Discharge screen will be displayed.
Rapid Discharge	Confirmation window to erase the data will be displayed.
HR/PR	The HR/PR numeric data box will be switched between HR and PR.

Page 2	
HR/PR Source	The parameter for HR/PR Source will be automatically selected.
BP Zero	Zero balance of BP1 to BP8 will be performed.
Leads	List of lead groups will be displayed, and selecting a lead group will display the lead selection window.2 blocks of user key area are required to assign this key. It cannot be assigned to the numeric data area.
ECG Size (All Leads)	The waveform size for all ECG leads can be changed.
Scale	The home display will change to scale selection mode.
Scale (Extended Display)	The extended display waveform size/scale setup menu will be displayed. This setting can be performed only when [OFF] is set for "Sync wave size/scale of extended display with main unit".
SpO ₂ -1 Display ON/OFF	SpO ₂ -1 display will turn ON/OFF.
SpO ₂ -2 Display ON/OFF	SpO ₂ -2 display will turn ON/OFF.
CO ₂ Display ON/OFF	CO ₂ display will turn ON/OFF.
GAS Display ON/OFF	Multigas unit data display will turn ON/OFF.
Auto Display Config.	The display will be automatically configured with the currently measured parameters.
Enlarged Display	For the standard display layout, the numeric data box width will enlarge.
Short Trend ON/OFF	Short Trend display will turn ON/OFF.
Graphic/Tabular Trend ON/OFF	Graphic/Tabular Trend display will turn ON/OFF.
Transparent Window ON/OFF	Transparent window will turn ON/OFF.
Change Palette	Palette selection window will be displayed.
Graphic Trend	The graphic trend will be displayed.
Trend (Group)	List of trend groups will be displayed, and selecting a trend group will display the graphic trend.
Tabular Trend	The tabular trend will be displayed.
Tabular Trend (Group)	List of tabular trend groups will be displayed, and selecting a trend group will display the tabular trend.

Page 3	
NIBP List	NIBP list will be displayed.
Recall	Recall screen will be displayed.
Alarm History	Alarm history will be displayed.
OCRG	OCRG screen will be displayed.
ST	ST screen will be displayed.
Cardiac Output	CO measurement screen will be displayed.
Drug Calculation	Drug Calculation screen will be displayed.
PCWP	PCWP measurement screen will be displayed. If BP labeled as PAP is not measured, this screen will not be displayed.
Hemodynamics	Hemodynamics screen will be displayed.
Lung Function	Lung Function screen will be displayed.
Full Disclosure Waveform	Full disclosure waveform will be displayed.
12-Lead Analysis	12-lead analysis screen will be displayed.
12-Lead Print	
Tone/Volume	The "Tone/Volume" menu will be displayed.
NIBP Auto Mode	NIBP Auto Mode window will be displayed.
Alarm Setup (All)	Alarm settings for all parameters will be displayed.
Alarm Setup (Basic)	Alarm settings for basic parameters will be displayed.
Manual Printing	Manual printing setup screen will be displayed.
Display Configuration	The display configuration window will be displayed.
Time / Date	The Time/Date setup window will be displayed.

Page 4	
Other Bed	Other bed screen will be displayed.
STOPWATCH	Stopwatch screen will be displayed.
Group 1 to 2	Selection list of key group 1 to 2 will be displayed.
Group 3 to 5	Selection list of key group 3 to 5 will be displayed.
Event	Event selection list will be displayed. The selected event will be stored as recall waveform.
Print (LBP) Cancel	Printing on the laser printer will be canceled.
Oxygenator Mode	The "Oxygenator Mode" menu will be displayed.
Main Mode 1(Initial)	Main mode 1 (Initial) will be set as the monitoring mode.
Main Mode 2(Hemo.)	Main mode 2 (Hemo) will be set as the monitoring mode.
Main Mode 3 (Cardiac)	Main mode 9 (Cardiac) will be set as the monitoring mode.
Main Mode 4(Local)	Main mode 4 (Local) will be set as the monitoring mode.
Main Mode 5(Full)	Main mode 5 (FullI) will be set as the monitoring mode.
Main Mode 6(Heart)	Main mode 6 (Heart) will be set as the monitoring mode.
Main Mode 7(Neo.)	Main mode 7 (Neo.) will be set as the monitoring mode.
Main Mode 8(Recovery)	Main mode 8 (Recovery) will be set as the monitoring mode.
Main Mode 9 (Cardiac)	Main mode 9 (Cardiac) will be set as the monitoring mode.
Sub Mode 1 (Induct.)	Sub Mode 1 (Induct.) will be set as the monitoring mode.
Page 5	
Sub Mode 2 (Surgery)	Sub Mode 2 (Surgery) will be set as the monitoring mode.
Sub Mode 3 (Waking)	Sub Mode 3 (Waking) will be set as the monitoring mode.
Sub Mode 4 (12-Lead)	Sub Mode 4 (12-Lead) will be set as the monitoring mode.
Sub Mode 5	Sub Mode 5 will be set as the monitoring mode.
Sub Mode 6	Sub Mode 6 will be set as the monitoring mode.

* Default user mode names are displayed inside the brackets. The mode names can be changed.

(☞ Maintenance Manual "To Program the User Mode" P5-31)

WARNING

• After changing the mode, make sure that the monitoring setting is appropriate. When the mode is changed, patient classification, alarm settings, etc. will be changed.

NOTE

 If the [Minimize] key is not set to the user key area, pressing the (Minimize) key on the window will function the same as the (X) (Close) key.

Changing the Display Layout from the Home Display

When using the LC-8016TC/LC-8018TC Display Unit, holding down the home display area will change the mode to display layout change mode.



On the display layout change mode, the following layout change can be performed.

- · Changing the displayed position of the waveform/numeric data
- Changing the size of the waveform/numeric data
- Adding the waveform/numeric data
- Deleting the waveform/numeric data

Changing the Displayed Position of the Waveform/Numeric Data

Drag the waveform/numeric data to a desired position.

The color of the dragged position will change to yellow.

By releasing the finger where the color has changed to yellow, the data will be located to a new position.

NOTE

By pressing the [UNDO] key which is displayed after changing the position, the previous operation can be canceled.

Changing the Size of the Waveform/Numeric Data

The size of the waveform/numeric data can be changed by pinch in/out operation.



- To enlarge the size, there should be enough space on the display.
- · The size can be enlarged downwards or rightwards.

Adding the Waveform/Numeric Data

The waveform/numeric data can be added by touching the free area, or by holding down the interspace of currently displayed waveform/numeric data.

Deleting the Waveform/Numeric Data

The waveform/numeric data can be deleted by flicking the corresponding display area.

Ending the Display Layout Change Mode

The display layout change mode can be ended by pressing the [Setup completed] key on the upper right of the screen.

Tone/Volume

This section explains the tone/volume setup procedure for alarm sound, HR synchronized tone, key sound, and boot/ shutdown sound. The tone/volume setup screen also allows to turn OFF the ventilator alarm sound. The volume of the sound which notifies the completion of BP zero balance and NIBP measurement can be adjusted on "Other" setting.

- NOTE
- The tone setup for the synchronized tone is effective only for HR and BP synchronized tone. The tone for SpO₂ synchronized sound will change according to the SpO₂ value. The tone will increase as the SpO₂ value increases, and vice versa.
- The start-up sound will always generate for the DS-8400 even when the "Boot/Shutdown" sound is set to minimum level.
- When the DS-8400 is activated with the DS-8007 connected and soon enters into standby mode, a shutdown sound may generate on the DS-8007.

1 Press the [Menu], [Sound] ("Basic Setup") keys.

> The "Tone/Volume" menu will be displayed.



2 Set the volume.

WARNING

 Changing the setting for "Alarm System" (Initial Settings > Alarm) will also change the alarm volume and tone setting. Make sure to check the volume and tone when the setting is changed.

- Pay attention not to set the alarm volume too low to avoid missing any important alarms. The alarm sound for ECG, SpO₂, CO₂ will be different from the test sound. The set volume will be applied but the set tone will not be applied to these parameters.
- When [Standard Tone] is set for the "Alarm System", the alarm volume and tone for the ventilator alarm and equipment status alarm will be the same with that of the vital alarm.

REFERENCE

- 1 Slide the / up or down.

- When the slider is released, $\blacktriangle / \checkmark$ will be displayed.
- 2 Press the A/V keys.
 - ▶ The volume will be adjusted.

REFERENCE

 The order of alarm priority is Urgent (H) > Caution (M) > Status (L). The volume is also set according to the alarm priority. The volume for high priority alarm cannot be set lower than the lower priority alarm, and vice versa.

3Select the tone level.

NOTE

- The tone selection is different for the synchronized tone, alarm sound, and key sound.
- For the "Sync. Tone", [Selected Tone] will generate the HR synchronized tone with the selected tone. [Sync. with SpO₂ Value] will generate the HR synchronized tone with the same tone with the SpO₂ synchronized tone. If the SpO₂ value is invalid, [Tone 2] will be applied.

4 Press the [Test] key to check the set volume/tone.

5 Set [ON]/[OFF] for ventilator alarm sound.

Color

In this section, setup procedure for the color of background, numeric data, waveform is explained.

The colors of the background, numeric data, waveform, user key can be customized.

The colors can be customized according to the various monitoring scene such as recognizable colors from a far distance or colors which will not strain your eyes by the long time monitoring.

Press the [Menu], [Color] ("Basic Setup") keys.

> The "Color" selection window will be displayed.



 $\mathbf{2}$ Press the key for "Background", and set the background color.

- The background color for the numeric data area and waveform area can be selected from three colors (black, gray, light gray).
- The background color can be also set by pressing the [Menu], [Display Config.] ("Basic Setup"), "Background" keys.
- The selected background color will be immediately reflected.

 ${f 3}$ Set the color of the numeric data and waveforms

The color can be set for each parameter. 12 colors (+white) from each palette are selectable.

- 1 Pressing the [Palette] key will display the "Palette" selection window. (shown on right)
- 2 Select the palette from [Light] / [Clear] / [Deep] / [Vivid], and press the [Set] key.
 - The color of the numeric data and waveform will change to the selected palette color.
- 3 Press ▲ ▼ to switch the page.
- 4 Press the key for the parameter to change the color.
 - ▶ The "Color" selection window will be displayed.

	Color					\times							
	1	2	3	4	5	6	7	8	9	10	11	12	White
Light													
Clear													
Deep													
Vivid													

5 Select a color.

The selected color for the parameter will be applied to the waveform, numeric data, graphic trend, and tabular trend.

4 Set the color of the user key.

- 1 Press the key for "User Key".
 - The "User Key Color" selection window will be displayed.
- 2 Press \blacksquare to switch the page.
- 3 Select the user key to change the color.
 - Pressing the key again will cancel the selection.
- **4** Select the color displayed on the right.
 - ▶ The color of the user key will change.

	User key Color	(\times)
Home	Wonitor Suspend	
Henu	Night Wode	
Winimize	Freeze	°⊾
Restore #indo#	Key Lock	
User Key 🜲	Wode Select	
Alarm Silence	Admit/Disch.	
Alarm Suspend	Rapid Discharge	
NIBP Start/Stop	HR/PR	
NIBP Cont.	HR/PR Source	
Print Start/Stop	BP Zero	



	Palette	(\mathbf{X})
Palette	Light	
	Vivid	Set

Brightness

In this section, brightness adjustment of the monitor display is explained.

• This equipment utilizes LED for the backlight. Since this LED deteriorates by the life cycle, the display may become dark, scintillate, or may not light by the long term use. In such case, contact your nearest service representative.

1 Press the [Menu], [Brightness] ("Basic Setup") keys.

• The "Brightness" menu will be displayed.



When the slider is released, ▲/▼ will be displayed.

 $\mathbf{3}$ Press the $\mathbf{A}/\mathbf{\nabla}$ keys.

• The brightness will be adjusted.

Menu > Basic Set		<u>رم</u>
Brightness	Night Mode	$\underline{(\uparrow)}$
Explanation	Area	
Brightness Fright Fight Bare		

Night Mode

This section explains about the night mode setup procedure.

The night mode is a function to preset the screen brightness and alarm volume when turning OFF the light of the ward or when the patient is asleep, etc.

The night mode can be manually set to ON, or automatically set to ON by preprogramming the time to turn ON/OFF the night mode.

Operation flow when the night mode is set to "Timer"



Operation flow when the night mode is set to [Darker] or [Dark]

1 To manually set the night mode, select [ON] for "Night Mode" or press the [Night Mode] key on the user key area.

BED-001 FUKUDA DENSHI - MANT KITAL		± Power ⊕ 13: 301//	56
	une Lata	HR 00ml Av. st I 0.2 Annu Av. (mn) 8VR	
		V ^{IIIC} 30 UU , Menu	
1991 2000			nce
Menu > Basic Setup	<u>ک</u>	Disch. 13	Zero 1975
Social Sector		COP	
Night Mode Wanual	Night Mode		arm stup stup
	OFF	NIBP(weeka) 129/°82 (* 98)	arn tory
		B.Atert NIBP List Re-	call
Detail Setup Yotune 1 Sternal Wonitor Night Wede 08 08 08 08 08 08 08 0			end
Display Darker			int L/Stop Linu
Atarn Indicator OFF		RR(i-II) 30	t Hode

• During the night mode, "Night Mode Active" message will be displayed.

NOTE

• When the timer is set, the night mode will automatically start at the set "Start Time".

$\mathbf{2}$ Cancel the night mode.

(@Maintenance Manual "Display/Print Setup" P5-13)

NOTE

- The night mode can be manually turned ON from the menu, user key, or remote control even when the night mode is set to automatically turn ON. The night mode will automatically turn OFF at the set "End Time".
- The night mode cannot be set when the ventilator alarm is generated.
- The night mode cannot be set during the battery operation.

Night Mode

The time to start and end the night mode, and the night mode display can be set.

Press the [Menu], [Night Mode] ("Basic Setup") keys.

> The "Night Mode" menu will be displayed.





 ${f Z}$ Set the "Start Time" and "End Time" for the night mode.

- [Manual]: The night mode can be turned ON or OFF manually using the user key.
- [Timer]: The night mode will automatically turned ON or OFF at the preprogrammed time. The night mode can be manually turned ON from the user key or from the remote control unit even when the [Timer] is set.

When [Timer] is selected:

1 Press the key for "Start Time".

- ▶ The "Start Time" window will be displayed. (shown on right)
- 2 Use the numeric keys to enter the time.
- 3 Press the [Set] key.

3 Set the volume.

4 Set the "End Time" with the same procedure from Step 3 to 5.



WARNING

- When selecting [Silence], pay attention not to miss any important alarm by simultaneously monitoring the bed on other monitors such as central monitor.
- [No Change]: Standard volume will be set.
- [3]: Third level from the minimum volume will be set.
- ▶ [1]: Minimum volume will be set.
- [Silence]: Sound will be silenced.

4 Set the brightness.

WARNING

- When selecting [Time Only], pay attention not to miss any important alarm by simultaneously monitoring the patient on other monitors such as central monitor.
- ▶ [No Change]: Brightness will not change
- [Dark]: 80% of the maximum brightness will be set.
- [Darker]: 50% of the maximum brightness will be set.
- [Time Only]: Only the time will be displayed. The message will disappear after 1 minute from starting the night mode.

5 Set the alarm indicator operation.

- [ON]: The alarm indicator will light even during the night mode.
- [OFF]: The alarm indicator will not light during the night mode.

6 Set the external monitor operation.

- [ON]: Displays the home display on the external monitor.
- [OFF]: Turns OFF the external monitor display.
- ▶ [OFF (Time Only)]:

If [Time Only] is selected for "Display": Displays the [Time Only] screen on the external monitor as well as the main unit.

If [No Change], [Dark] or [Darker] is selected for "Display": Turns OFF the external monitor display.

Chapter 11 Troubleshooting

Message List	11-1
Vital Alarm Message	11-1
Vital Alarm Message (DS-LAN Standard Setup)	11-4
Equipment Status Alarm Message	11-8
Numeric Data Box Message	
Ventilator Alarm Message	
Ventilator Alarm Factor	
Cardiac Output Message	
Troubleshooting	11-22
ECG	
Respiration	11-27
Invasive Blood Pressure	
SpO2 Measurement (HS-8312N, DS-8007N, HG-820)	
SpO2 Measurement (HS-8312M, DS-8007M, HG-810)	11-33
Non-Invasive Blood Pressure	
Temperature	11-41
Cardiac Output (CO)	
CO2 Measurement (HCP-800/HCP-810/HCP-820)	
CO2 Measurement (HPD-800/HPD-810/HPD-820)	11-46
Recorder Unit (HR-800)	11-48
Network Printer	
Wired Network (DS-LANII/ DS-LANII)	11-49
Telemeter (HLX-801)	11-51
Remote Control	11-52
General	
Super Unit	11-57
Data Transfer Function	
IB-8004 Input Box	11-59
Expansion Module	11-61
Extended Display Unit	11-63
Ventilator	11-63
Multigas Unit	11-64
SvO2/CCO Monitor	11-66
BIS Monitor (A-2000/A-3000)	
BIS (When HBX-800 is used)	11-68
INVOS	
FLOW-i	
PC Communication	
TCM4/TCM5 FLEX	
Magnetic Card Reader/Barcode Reader	
External Media	11-75

Chapter 11 Troubleshooting

Message List

This section lists the alarm messages for each parameter.

For the vital alarm message, there are numeric data alarm and arrhythmia alarm, and the delay time are as follows.

- Numeric Data Alarm: Adult/Child: 5 sec., Neonate: none However, for HR alarm, there is no delay time for adult/child if "HR Delay" is set to ON.
- Arrhythmia Alarm: Adult/Child/Neonate: none

Vital Alarm Message

- The arrhythmia alarm message other than Tachy, Brady, Ext Tachy, Ext Brady will continue to be displayed for 30 seconds after the alarm is resolved.
- The alarm level shown below is the standard level set by Fukuda Denshi.
- The alarm level can be changed on the "Initial Settings".

Top Priority Alarm (Alarm Level S)

This level can be selected for some parameters only when [Fukuda Tone] is selected for the "Alarm System" ("Initial Settings"). It cannot be selected when the "Alarm System" is [Melodic Tone] or [Standard Tone].

Life Threatening Alarm (Alarm Level H)

Parameter	Message
Respiration (Impedance, CO ₂ , Ventilator)	<apnea></apnea>
SpO ₂ *	<lower ext="" spo<sub="">2# Alarm></lower>
Arrhythmia	<asystole></asystole>
	<vf></vf>
	<vt></vt>
	<slow vt=""></slow>
	<tachy></tachy>
	<brady></brady>
	<ext tachy=""></ext>
	<ext brady=""></ext>

Cautionary Alarm (Alarm Level M)

Parameter	Message
HR	<lower alarm="" hr=""></lower>
	<upper alarm="" hr=""></upper>
Blood Pressure	<lower alarm="" bp#=""> or <lower (label)="" alarm="">^{*1}</lower></lower>
	<upper alarm="" bp#=""> or <upper (label)="" alarm="">^{*1}</upper></upper>

Parameter	Message
Pulse Rate	<lower alarm="" pr=""></lower>
(BP)	<upper alarm="" pr=""></upper>
SpO ₂ *	<lower spo<sub="">2 # Alarm>^{*1}</lower>
	<upper spo<sub="">2 # Alarm>^{*1}</upper>
	<lower ext="" spo<sub="">2# Alarm></lower>
	<upper ext="" spo<sub="">2# Alarm></upper>
Pulse Rate	<lower alarm="" pr=""></lower>
(SpO ₂)	<upper alarm="" pr=""></upper>
Non-Invasive Blood Pressure	<lower alarm="" nibp=""></lower>
	<upper alarm="" nibp=""></upper>
Respiration	<lower alarm="" rr=""></lower>
(Impedance, CO_2 , Gas, Ventilator, SpO_2)	<upper alarm="" rr=""></upper>
Gas ^{*2}	<lower co<sub="">2-E Alarm></lower>
	<upper co<sub="">2-E Alarm></upper>
	<upper co<sub="">2-I Alarm></upper>
	<lower o<sub="">2-E Alarm></lower>
	<upper o<sub="">2-E Alarm></upper>
	<lower o<sub="">2-I Alarm></lower>
	<upper o<sub="">2-I Alarm></upper>
	<lower n<sub="">2O-E Alarm></lower>
	<upper n<sub="">2O-E Alarm></upper>
	<lower n<sub="">2O-I Alarm></lower>
	<upper n<sub="">2O-I Alarm></upper>
	<lower (agt="" alarm="" label)-e=""></lower>
	<upper (agt="" alarm="" label)-e=""></upper>
	<lower (agt="" alarm="" label)-i=""></lower>
	<upper (agt="" alarm="" label)-i=""></upper>
SPIRO ^{*2}	<lower alarm="" mv=""></lower>
	<upper alarm="" mv=""></upper>
BIS (When HBX-800 is used)	<lower alarm="" bis=""></lower>
	<upper alarm="" bis=""></upper>
Arrhythmia	<run></run>
	<pause></pause>

*1: # indicates the label of BP, TEMP, SpO₂. For SpO₂, N1/N2/M1/M2/HR/HL/FR/FL/OT will be displayed for #.

*2: When the numeric data acquired from FLOW-i is displayed, alarm will not generate. The alarm will not generate on the central monitor either.

Treatment Needed Alarm (Alarm Level L)

Parameter	Message	
ST1 to 12	<lower alarm="" st(lead="" type)=""></lower>	
	<upper alarm="" st(lead="" type)=""></upper>	
SpCO#	<upper alarm="" spco#="">^{*1}</upper>	
SpMet#	<upper alarm="" spmet#="">^{*1}</upper>	
SpHb	<lower alarm="" sphb#=""> ^{*1}</lower>	
	<upper alarm="" sphb#="">^{*1}</upper>	
Temperature	<lower alarm="" temp#=""> or <lower (label)="" alarm="">^{*1}</lower></lower>	
(TEMP1 to 8)	<upper alarm="" temp#=""> or <upper (label)="" alarm="">^{*1}</upper></upper>	
Blood Temperature	<upper alarm="" tb=""></upper>	
	<lower alarm="" tb=""></lower>	
Arrhythmia	<couplet></couplet>	
	<bigeminy></bigeminy>	
	<trigeminy></trigeminy>	
	<frequent></frequent>	
	<triplet></triplet>	
	<r on="" t=""></r>	
	<multiform></multiform>	
	<vent. rhythm=""></vent.>	
	<svt></svt>	
	<irregular rr=""></irregular>	
	<prolonged rr=""></prolonged>	
	<s frequent=""></s>	
	<s couplet=""></s>	
	<vpc></vpc>	
	<svpc></svpc>	
	<pacer capture="" not=""></pacer>	
	<pacer not="" pacing=""></pacer>	
SPIRO*2	<upper alarm="" peak=""></upper>	
	<lower alarm="" peak=""></lower>	
	<upper alarm="" peep=""></upper>	
	<lower alarm="" peep=""></lower>	

*1: # indicates the channel number of BP, TEMP, SpCO, SpMet, SpHb.

For SpCO, SpMet, SpHb, N1/N2/M1/M2/HR/HL/FR/FL/OT will be displayed for #.

*2: When the numeric data acquired from FLOW-i is displayed, alarm will not generate. The alarm will not generate on the central monitor either.

□Notification Alarm

Parameter	Message
All Alarm	<alarm (xxx="" sec.)="" suspend=""></alarm>
Alarm Sound Suspend	<alarm (xxx="" min.)="" silence=""></alarm>
Arrhythmia	<learn></learn>
	<arrhy. off=""></arrhy.>
Oxygenator Mode	<all alarm="" off=""></all>

NOTE

- (xxx sec) of the <Alarm Suspend (xxx sec)> message indicates the remaining time of alarm suspended duration.
- (xxx min.) of the <Alarm Silence (xxx min.)> message indicates the remaining time of alarm sound suspended duration.
- The <ARRHY OFF> message will be displayed when the Slow_VT, Tachy, Brady, Ext Tachy, Ext Brady and HR alarm is OFF.

Vital Alarm Message (DS-LAN Standard Setup)

WARNING

• The SpO₂ respiration measurement function is not intended for use as an APNEA monitor.

- The arrhythmia alarm message other than Tachy, Brady, Ext Tachy, Ext Brady will continue to be displayed for 30 seconds after the alarm is resolved.
- The alarm level shown below is the standard level set by Fukuda Denshi.

Top Priority Alarm (Alarm Level S)

This level can be selected for some parameters only when [Fukuda Tone] is selected for the "Alarm System" ("Initial Settings"). It cannot be selected when the "Alarm System" is [Melodic Tone] or [Standard Tone].

Life Threatening Alarm (Alarm Level H)

Parameter	Message
HR	<lower alarm="" hr=""></lower>
	<upper alarm="" hr=""></upper>
Pulse Rate	<lower alarm="" pr=""></lower>
(SpO ₂)	<upper alarm="" pr=""></upper>
Pulse Rate	<lower alarm="" pr=""></lower>
(BP)	<upper alarm="" pr=""></upper>
SpO ₂ *	<lower spo<sub="">2 # Alarm>^{*1}</lower>
	<upper spo<sub="">2 # Alarm>^{*1}</upper>
	<upper ext="" spo<sub="">2 Alarm></upper>
	<lower ext="" spo<sub="">2 Alarm></lower>
Blood Pressure	<lower alarm="" bp1=""></lower>
	<upper alarm="" bp1=""></upper>
	<lower alarm="" art=""></lower>
	<upper alarm="" art=""></upper>
Non-Invasive Blood Pressure	<lower alarm="" nibp=""></lower>
	<upper alarm="" nibp=""></upper>
Respiration	<lower alarm="" rr=""></lower>
(Impedance, CO ₂ , GAS, Ventilator, SpO ₂)	<upper alarm="" rr=""></upper>
	<apnea></apnea>
Gas ^{*2}	<lower co<sub="">2-E Alarm></lower>
	<upper co<sub="">2-E Alarm></upper>
	<upper co<sub="">2-I Alarm></upper>
	<lower o<sub="">2-E Alarm></lower>
	<upper o<sub="">2-E Alarm></upper>
	<lower o<sub="">2-I Alarm></lower>
	<upper o<sub="">2-I Alarm></upper>
	<lower n<sub="">2O-E Alarm></lower>
	<upper n<sub="">2O-E Alarm></upper>
	<lower n<sub="">2O-I Alarm></lower>
	<upper n<sub="">2O-I Alarm></upper>
	<lower (agt="" alarm="" label)-e=""></lower>
	<upper (agt="" alarm="" label)-e=""></upper>
	<lower (agt="" alarm="" label)-i=""></lower>
	 <upper (agt="" alarm="" label)-i=""></upper>
Arrhythmia	<asystole></asystole>
	VF>
	<vt></vt>
	<slow vt=""></slow>
	<tachy></tachy>
	<brady></brady>
	<run></run>
	<ext tachy=""></ext>
	=x radify

Parameter	Message
	<ext brady=""></ext>

*1: For SpO₂, N1/N2/M1/M2/HR/HL/FR/FL/OT will be displayed for #.

*2: When the numeric data acquired from FLOW-i is displayed, alarm will not generate. The alarm will not generate on the central monitor either.

Cautionary Alarm (Alarm Level M)

Parameter	Message
Blood Pressure	<lower 8="" alarm="" bp2="" to=""> or <lower (label="" alarm="" art)="" other="" than="">^{*1}</lower></lower>
	<upper 8="" alarm="" bp2="" to=""> or <upper (label="" alarm="" art)="" other="" than="">^{*1}</upper></upper>
ST1 to 12	<lower alarm="" st(lead="" type)=""></lower>
	<upper alarm="" st(lead="" type)=""></upper>
SpCO#	<upper alarm="" spco#="">^{*1}</upper>
SpMet#	<upper alarm="" spmet#="">^{*1}</upper>
SpHb	<lower alarm="" sphb#=""> *1</lower>
	<upper alarm="" sphb#="">^{*1}</upper>
TEMP (TEMP1 to 8)	<upper alarm="" temp#=""> or <upper (label)="" alarm="">^{*1}</upper></upper>
	<lower alarm="" temp#=""> or <lower (label)="" alarm="">^{*1}</lower></lower>
Blood Temperature	<upper alarm="" tb=""></upper>
	<lower alarm="" tb=""></lower>
BIS	<upper alarm="" bis=""></upper>
	<lower alarm="" bis=""></lower>
MV ^{*2}	<upper alarm="" mv=""></upper>
	<lower alarm="" mv=""></lower>
PEAK ^{*2}	<upper alarm="" peak=""></upper>
	<lower alarm="" peak=""></lower>
PEEP ^{*2}	<upper alarm="" peep=""></upper>
	<lower alarm="" peep=""></lower>
Arrhythmia	<pause></pause>
	<couplet></couplet>
	<bigeminy></bigeminy>
	<trigeminy></trigeminy>
	<frequent></frequent>
	<triplet></triplet>
	<r on="" t=""></r>
	<multiform></multiform>
	<vent. rhythm=""></vent.>
	<svt></svt>
	Irregular RR>
	<prolonged rr=""></prolonged>
	<s frequent=""></s>
	<s couplet=""></s>
	<vpc></vpc>
	<svpc></svpc>

Parameter	Message	
	<pacer capture="" not=""></pacer>	
	<pacer not="" pacing=""></pacer>	

*1: # indicates the channel number of BP, TEMP, SpCO, SpMet, SpHb. For SpCO, SpMet, SpHb, N1/N2/M1/M2/HR/HL/FR/FL/OT will be displayed for #.

*2: When the numeric data acquired from FLOW-i is displayed, alarm will not generate. The alarm will not generate on the central monitor either.

Notification Alarm

Parameter	Message
All Alarm	<alarm (xxx="" sec.)="" suspend=""></alarm>
Alarm Sound Suspend	<alarm (xxx="" min.)="" silence=""></alarm>
Arrhythmia	<learn></learn>
	<arrhy. off=""></arrhy.>
Oxygenator Mode	<all alarm="" off=""></all>

NOTE

- (xxx sec) of the <Alarm Suspend (xxx sec)> message indicates the remaining time of alarm suspended duration.
- (xxx min.) of the <Alarm Silence (xxx min.)> message indicates the remaining time of alarm sound suspended duration.
- The <ARRHY OFF> message will be displayed when the Slow_VT, Tachy, Brady, Ext Tachy, Ext Brady and HR alarm is OFF.

Equipment Status Alarm Message

Top Priority Alarm (Alarm Level S)

Item	Message	Delay Time (sec.)
Ventilator	<vent. alarm=""></vent.>	1
	<vent comm=""></vent>	1

Life Threatening Alarm (Alarm Level H)

Item	Message	Delay Time (sec.)
Main Unit	<dsc-8400 failure=""></dsc-8400>	10
	<dsc-8400 failure="" speaker=""></dsc-8400>	10
	<charge battery.="" the=""></charge>	10
Super Unit	<super failure="" unit=""></super>	3
	<ecg error="" unit=""></ecg>	5
	<super failure="" multiamp.="" unit=""></super>	3
	<super failure="" transducer="" unit="" voltage=""></super>	3
	<nibp (xxx-xxx)="" error="" meas.="">*1</nibp>	10 or 3
	<gas f="" failure="" i="" unit=""></gas>	3
	<super spo<sub="" unit="">2 Failure></super>	5 or 1
	<charge battery.="" the=""></charge>	10
Blood Pressure	<transducer failure="" voltage=""></transducer>	3
	<check art="" catheter.="" the=""></check>	1
GAS (MGU-800/MGU-810)	<gas failure="" unit=""></gas>	1
SPIRO	<spiro error="" unit=""></spiro>	1
BIS (When HBX-800 is used)	<bisx failure=""></bisx>	3
	<bisx incompatible=""></bisx>	3

*1: # indicates an error code.

Cautionary Alarm (Alarm Level M)

Item	Message	Delay Time (sec.)
NIBP	<nibp (###-##)="" failed.="" meas.="">^{*1}</nibp>	1
CO ₂ (HCP-800/HCP-810/	<co<sub>2 Check Sample Line></co<sub>	1
HCP-820)	<co<sub>2 Check Exhaust Port></co<sub>	1
	<co<sub>2 Unit Failure></co<sub>	1
CO ₂ (HCP-800/HCP-810/ HCP-820)	<co<sub>2 Cal. Required></co<sub>	1
Capnostat 5 CO ₂ (Gas Unit I/ F and Mainstream Module)	<co<sub>2 Sensor Failure></co<sub>	1
GAS (MGU-800)	<gas check="" class="" trap="" water=""></gas>	1
	<gas off="" pump=""></gas>	1
	<gas check="" line="" sample=""></gas>	1
	<gas failed="" zeroing=""></gas>	1
	<gas replace="" trap="" water=""></gas>	1

Item	Message	Delay Time (sec.)
	<gas check="" conn.="" trap="" water=""></gas>	1
	<gas check="" conn.=""></gas>	1
SPIRO (MGU-810)	<spiro check="" class="" flowsensor=""></spiro>	1
BIS (When HBX-800 is used)	<check bis="" check="" perform="" sensor="" sensor,=""></check>	3
Main Unit	<dsc-8400 battery="" check="" short-term=""></dsc-8400>	10
	<dsc-8400 battery="" check="" long-term=""></dsc-8400>	10
Super Unit	<super check="" conn.="" unit=""></super>	3
	<super of="" operating="" out="" range="" temp.="" unit=""></super>	3
	<super analog="" unadjusted="" unit=""></super>	3
Input Box	<ib-8000-# check="" conn.="">*2</ib-8000-#>	3
	<ib-8000-# failure="">^{*2}</ib-8000-#>	3
Display Unit	<display failure="" unit=""></display>	3
Module	<ib# failure="" module="" slot#="">^{*3}</ib#>	3
	<ib# analog="" slot#="" unadjusted="">^{*3}</ib#>	3
Monitor Suspend	<monitor suspend="" time-out=""></monitor>	1
Full Disclosure Waveform	<failed card.="" cf="" disclosure="" full="" the="" to="" write=""></failed>	1

*1: On "Initial Settings" menu, the alarm level can be selected from Level M, L, N (Notification). (Default: Level M) If [Alarm Silence] key is pressed during Level M, L alarm generation, the alarm level will change to Level N (notification).

indicates an error code.

*2: # indicates the input box number.

*3: # indicates the input box number, and the slot number of input box.

Treatment Needed Alarm (Alarm Level L)

Item	Message	Delay Time (sec.)
ECG	<check #)="" #,="" (#,="" electrodes="">^{*1}</check>	3
	<ecg attachment.="" check="" electrodes=""></ecg>	3
	<cannot analyze=""></cannot>	1
	<ecg detection="" error="" pacing=""></ecg>	1
	<ecg artifact=""></ecg>	3
	<ecg 5="" are="" electrodes="" only="" used.=""></ecg>	1
Impedance	<rr exceeded.="" is="" meas.="" range=""></rr>	3
	<cva detected=""></cva>	Adult, Child: 20, Neonate: 10
SpO ₂ (Masimo Unit)	<spo<sub>2- # Check Sensor Attach.>*2</spo<sub>	3
	<spo<sub>2- # Replace Sensor>^{*2}</spo<sub>	1
	<spo<sub>2- # Low Perfusion>^{*2, *3}</spo<sub>	1
	<spo<sub>2- # Pulse Search>*2</spo<sub>	1
	<spo<sub>2- # Noise Interference>^{*2}</spo<sub>	1
	<spo<sub>2- # Check Sensor>^{*2}</spo<sub>	1
	<spo<sub>2- # Check Sensor Conn.>*2</spo<sub>	1
	<spo<sub>2- # Replace Cable>^{*2}</spo<sub>	3
	<spo<sub>2- # Check Cable>^{*2}</spo<sub>	3
	<spo<sub>2- # Disconnected>*2</spo<sub>	3

Item	Message	Delay Time (sec.)
	<spo<sub>2- # only mode>^{*2}</spo<sub>	1
	<spo<sub>2- # Check Cable, Sensor>^{*2}</spo<sub>	1
SpO ₂ (Nellcor Unit)	<spo<sub>2- # Check Sensor Attach.>*2</spo<sub>	3
	<spo<sub>2- # Replace Sensor>^{*2}</spo<sub>	1
	<spo<sub>2- # No Pulse Detected>*2</spo<sub>	1
Blood Pressure	<bp #="" off="" transducer="">^{*4*9}</bp>	5
Temperature	<t ##="" sensor="" unknown="">^{*5}</t>	3
Non-Invasive Blood Pressure	<check cuff,="" hose="" nibp="">^{*6}</check>	3
	<nibp air="" check="" hose="" patient="" type,=""></nibp>	3
Capnostat 5 CO ₂ (Gas Unit I/ F and Mainstream Module)	<check co<sub="">2 Airway Adapter></check>	1
SPIRO (MGU-810)	<spiro check="" flow="" sensor=""></spiro>	1
BIS (When HBX-800 is used)	<replace bis="" sensor=""></replace>	3
	<bis sensor="" usage=""> 24hrs.>^{*2}</bis>	3
	<bis disconnected="" sensor=""></bis>	1
	<bis check="" high="" impedance,="" sensor=""></bis>	3
	<bis check="" lead="" off,="" sensor=""></bis>	3
	<bis 15%="" <="" sqi="">^{*2}</bis>	3
	<bisx disconnected=""></bisx>	3
Connector Off	<ecg disconnected=""></ecg>	3
	<bp #="" disconnected="">*4</bp>	3
	<spo<sub>2- # Disconnected>^{*2}</spo<sub>	3
	<t ##="" disconnected="">^{*5}</t>	3
	<co disconnected=""></co>	3
	<co<sub>2 Disconnected></co<sub>	3
Main Unit	<dsc-8400 check="" unit=""></dsc-8400>	10
	<dsc-8400 of="" operating="" out="" range="" temp.=""></dsc-8400>	10
	<charge battery.="" the=""></charge>	10
	<reinstall battery.="" the=""></reinstall>	5
	<fan failure=""></fan>	3
Super Unit	<super card="" check="" sd="" unit=""></super>	3
	<super check="" dip-sw="" unit=""></super>	3
	<super failure="" temp="" unit=""></super>	3
	<super data="" failed.="" transfer="" unit=""></super>	3
Input Box	<ib-8000-# failure="">^{*7}</ib-8000-#>	3
	<ib-8000-# of="" operating="" out="" range="" temp.="">^{*7}</ib-8000-#>	3
Display Unit	<check display="" unit=""></check>	3
	<display of="" operating="" out="" range="" temp.="" unit=""></display>	3
Module	<ib# check="" module="" slot#="">^{*8}</ib#>	3
	<ib# of="" operating="" out="" range="" slot#="" temp.="">^{*8}</ib#>	3
	<ib# failure="" module="" slot#="">^{*8}</ib#>	3
	<ib# disconnected="" module="" slot#="">^{*8}</ib#>	3
	<ib# failure="" slot#="" temp="" unit="">^{*8}</ib#>	3

Item	Message	Delay Time (sec.)
Built-in Slot	<built-in check="" module="" slot=""></built-in>	3
	<built-in of="" operating="" out="" range="" slot="" temp.=""></built-in>	3
	<built-in failure="" module="" slot=""></built-in>	3
	<built-in disconnected="" module="" slot=""></built-in>	3
	<built-in failure="" slot="" temp="" unit=""></built-in>	3
Check Connection, Check Reception, Interference	<check svo<sub="">2/CCO Monitor Conn.></check>	1
	<check bis="" conn.=""></check>	1
	<check conn.="" invos=""></check>	1
	<check conn.="" flow-i=""></check>	1
	<check conn.="" printer=""></check>	3
	<chk comm="" ds-lan=""></chk>	3
	<check conn.="" hlx=""></check>	3
	<check comm="" printer=""></check>	1
	<check conn.="" tcm=""></check>	1
Full Disclosure Waveform	<wrong card="" cf="" disclosure.="" for="" full=""></wrong>	1
	<failed card.="" cf="" disclosure="" from="" full="" read="" the="" to=""></failed>	1
	<check card="" cf="" disclosure.="" for="" full=""></check>	1

*1: # indicates an electrode type.

*2: # indicates the label of SpO₂.

*3: On "Initial Settings" menu, the alarm level can be selected from Level L/N. (Default: Level L)

*4: # indicates the label of BP.

*5: # indicates the label of TEMP.

*6: On "Initial Settings" menu, the alarm level can be selected from Level M/L/N. (Default: Level L)

If [Alarm Silence] key is pressed during Level M/L alarm generation, the alarm level will change to Level N (notification).

*7: # indicates the input box number.

*8: # indicates the input box number, and the slot number of input box.

*9: On "Initial Settings" menu, the alarm level can be selected from Level M/L. (Default: Level L)

|--|

 <NIBP meas. failed>, <Check NIBP cuff, hose>, <Connector Off>, <ECG Only 5 electrodes are used.>, <Check xx Conn.>, <Check xx Comm.>, <SPIRO Check Flow Sensor> alarms will be canceled when [Alarm Silence] key is pressed. Pay attention not to cancel the important alarm.

□Notification Alarm

Item	Message	Delay Time (sec.)
Operation	<waveform (xxsec.)="" frozen="">^{*1}</waveform>	1
	<key (xxsec.)="" locked="">^{*1}</key>	1
	<night active="" mode=""></night>	1
	<oxygenator mode=""></oxygenator>	1
ECG	<ecg amplitude="" low=""></ecg>	3
	<ecg artifact=""></ecg>	3
	<ecg emg="" interference=""></ecg>	3
	<check electrodes="">*7</check>	3
Blood Pressure	<bp #="" required="" zeroing="">^{*2}</bp>	1
Temperature	<t #="" sensor="" unknown="">^{*3}</t>	1
SpO ₂ (Masimo Unit)	<spo<sub>2- # Demo Mode>^{*4}</spo<sub>	1
	<spo<sub>2- # Zeroing>^{*4}</spo<sub>	1
	<spo<sub>2 Check Sensor Attach.>*7</spo<sub>	3
	<spo<sub>2 Cable Near Expiration></spo<sub>	3
	<spo<sub>2 Sensor Near Expiration></spo<sub>	3
SpO ₂ (Nellcor Unit)	<spo<sub>2- # Motion Artifact>^{*4}</spo<sub>	1
	<spo<sub>2 Check Sensor Attach.>^{*7}</spo<sub>	3
Capnostat 5 CO ₂ (Gas Unit I/	<co<sub>2 Warming Up></co<sub>	1
F and Mainstream Module)	<zero co<sub="" the="">2 Adapter></zero>	1
	<unknown co<sub="">2 Sensor></unknown>	1
CO ₂ (HCP-800/HCP-810/ HCP-820)	<co<sub>2 Suspended></co<sub>	1
	<co<sub>2 Zeroing></co<sub>	1
GAS (MGU-800/MGU-810)	<gas up="" warm=""></gas>	1
	<gas zeroing=""></gas>	1
GAS (MGU-800/MGU-810)	<gas pump="" regulating=""></gas>	1
	<gas agents="" mixed="">^{*5}</gas>	1
	<gas cal.="" required.="" zero=""></gas>	1
	<gas cal.="" required.=""></gas>	1
SPIRO (MGU-810)	<spiro up="" warm=""></spiro>	1
, , , , , , , , , , , , , , , , , , ,	<spiro active="" calibration=""></spiro>	1
	<spiro zeroing=""></spiro>	1
BIS (When HBX-800 is used)	<bis expired="" sensor="">^{*5}</bis>	3
	Sensor Check in Progress>	3
	Second Check in Progress	3
	<bis noise=""></bis>	3
	<pre><bis "sensor="" check"="" perform=""></bis></pre>	3
	SQI < 50%>	3
	<pre><bis demo="" sensor=""></bis></pre>	3
Non-Invasive Blood Pressure	<initializing nibp=""></initializing>	3
Recorder Unit	<check printer="">^{*6}</check>	3
	<pre></pre> <pre><</pre>	3

ltem	Message	Delay Time (sec.)
	<printer busy="">^{*6}</printer>	1
	<check cassette="">*6</check>	3
Central Printer	<check (central)="" paper="">^{*6}</check>	3
	<check cassette="">^{*6}</check>	3
	<printer (central)="" busy="">^{*6}</printer>	1
	<check central="" printer="">*6</check>	3
Central Printer	<central check="" connection="" printer=""></central>	1
(Laser Printer)	<central check="" printer="" setting=""></central>	1
	<check central="" id=""></check>	1
	<chk comm="" ds-lan=""></chk>	1
Main Unit	<pre><dsc-8400 check="" rotary="" sw=""></dsc-8400></pre>	1
	<pre><dsc-8400 check="" dipsw=""></dsc-8400></pre>	1
System Configuration	<check config.="" equip.=""></check>	1
	<some are="" display="" displayed="" due="" layout="" not="" parameters="" setting.="" the="" to=""></some>	3
Check Connection, Check Reception, Interference	<check conn.="" system=""></check>	3
Data Transfer	<uploading></uploading>	1
	<upload standby=""></upload>	1

*1: ## indicates the remaining time.

*2: # indicates the channel number of BP.

*3: # indicates the channel number of TEMP.

*4: # indicates the label of SpO_2 .

*5: On "Initial Settings" menu, the alarm level can be selected from Level M, L, N (Notification). (Default : Notification)

*6: The alarm generation can be inhibited depending on the setting.

*7: Displayed when lead-off or sensor-off condition remains after the power is turned ON, monitoring is resumed, or a patient is discharged.

Numeric Data Box Message

HR

Message
<unit failure=""></unit>
<upper alarm="" hr=""></upper>
<lower alarm="" hr=""></lower>
<lower alarm="" st=""></lower>
<upper alarm="" st=""></upper>
<cannot analyze=""></cannot>
<check electrodes=""></check>
<check attachment.="" electrodes=""></check>
<pacing detection="" error=""></pacing>
<only 5="" are="" electrodes="" used.=""></only>
<out of="" range=""></out>
<low amplitude=""></low>
<noise interference=""></noise>
<artifact></artifact>

∎st

Message	
<lower alarm="" st=""></lower>	
<upper alarm="" st=""></upper>	

BP1 to 8

Level H for BP1 and ART, Level M for other label

Message
<lower alarm="" bp=""></lower>
<upper alarm="" bp=""></upper>
<zero required=""></zero>
<check catheter.="" the=""></check>
<out of="" range=""></out>

Pulse Rate (BP Source)

Message
<upper alarm="" pr=""> (BP)</upper>
<lower alarm="" pr=""> (BP)</lower>
<check catheter.="" the=""></check>
<out of="" range=""></out>

If <NIBP Meas. Error> is displayed, the message can be cancelled by pressing [Cancel Error] on the NIBP setup screen, or [NIBP Start/Stop] key (user key).

If the same message is repeatedly displayed, a failure of the equipment can be considered. Cease the measurement and contact your nearest service representative.

(@"<NIBP Unit Error (E**-**)> is displayed on the main unit." P11-40)

Message
<nibp error="" meas.=""></nibp>
<upper alarm="" nibp=""></upper>
<lower alarm="" nibp=""></lower>
<measurement failed.=""></measurement>
<check cuff,="" hose="" nibp=""></check>
<check air="" hose="" patient="" type,=""></check>
<initializing></initializing>
<out of="" range=""></out>

□SpO₂ (Nellcor Model)

Message
<unit failure=""></unit>
<ext spo<sub="">2 Alarm></ext>
<lower spo<sub="">2 Alarm></lower>
<upper spo<sub="">2 Alarm></upper>
<replace sensor=""></replace>
<check attach.="" sensor=""></check>
<no detected="" pulse=""></no>
<motion artifact=""></motion>
<pulse search=""></pulse>

□SpO₂/SpCO/SpMet/SpHb (Masimo Model)

Message
<ext spo<sub="">2 Alarm></ext>
<lower spo<sub="">2 Alarm></lower>
<upper spo<sub="">2 Alarm></upper>
<upper alarm="" spco=""></upper>
<upper alarm="" spmet=""></upper>
<lower alarm="" sphb=""></lower>
<upper alarm="" sphb=""></upper>
<replace sensor=""></replace>
<check attach.="" sensor=""></check>
<low confidence=""></low>
<pulse search=""></pulse>
<noise interference=""></noise>
<check sensor=""></check>
<replace cable=""></replace>

Message
<check cable=""></check>
<check conn.="" sensor=""></check>
<zeroing sensor=""></zeroing>
<spo<sub>2 only mode></spo<sub>
<low iq="" signal=""></low>
<low confidence=""></low>

RR (SpO₂: Nellcor Model)

Message
<unit failure=""></unit>
<rr interference=""></rr>
<unable calculate="" to=""></unable>
<calculating></calculating>
<outside range=""></outside>
<out of="" range=""></out>
<upper rr="" spo<sub="">2 Alarm></upper>
<lower rr="" spo<sub="">2 Alarm></lower>

PR-SpO₂

Message
<upper alarm="" pr=""> (SpO₂)</upper>
<lower alarm="" pr=""> (SpO₂)</lower>
<out of="" range=""></out>

TEMP1 to 8

Message
<upper alarm="" temp=""></upper>
<lower alarm="" temp=""></lower>
<temp failure="" unit=""></temp>
<unknown sensor=""></unknown>
<out of="" range=""></out>

Пть

Message
<lower alarm="" tb=""></lower>
<upper alarm="" tb=""></upper>
<out of="" range=""></out>

RR (Impedance)

Message
<apnea alarm=""></apnea>
<upper alarm="" rr=""></upper>
<lower alarm="" rr=""></lower>

Message
<cva detected=""></cva>
<rr exceeded.="" is="" meas.="" range=""></rr>
<out of="" range=""></out>
<suspended></suspended>

RR (Ventilator)

Message
<apnea alarm=""></apnea>
<upper alarm="" rr=""></upper>
<lower alarm="" rr=""></lower>

RR (Gas)

Message
<apnea alarm=""></apnea>
<upper alarm="" rr=""></upper>
<lower alarm="" rr=""></lower>
<out of="" range=""></out>

 \Box CO₂ (When Gas Unit I/F HPD-800/HPD-810/HPD-820 and Capnostat5 are used)

Message
<upper co<sub="">2-E Alarm></upper>
<lower co<sub="">2-E Alarm></lower>
<upper co<sub="">2-I Alarm></upper>
<check adapter.="" airway=""></check>
<zeroing></zeroing>
<gas up="" warm=""></gas>
<zero co<sub="">2 Adapter></zero>
<unknown sensor=""></unknown>
<out of="" range=""></out>

CO₂ (HCP-800/HCP-810/HCP-820)

Message
<initializing></initializing>
<check line="" sample=""></check>
<zeroing></zeroing>
<check exhaust="" port="" the=""></check>
<perform calibration.=""></perform>
<gas f="" failure="" i="" unit=""></gas>
<out of="" range=""></out>
<upper co<sub="">2-E></upper>
<lower co<sub="">2-E></lower>
<upper co<sub="">2-I></upper>

Gas (When MGU-800/810 is used)

Message
<upper co<sub="">2-E Alarm></upper>
<lower co<sub="">2-E Alarm></lower>
<upper co<sub="">2-I Alarm></upper>
<upper o<sub="">2-E Alarm></upper>
<lower o<sub="">2-E Alarm></lower>
<upper o<sub="">2-I Alarm></upper>
<lower o<sub="">2-I Alarm></lower>
<upper n<sub="">2O-E Alarm></upper>
<lower n<sub="">2O-E Alarm></lower>
<upper n<sub="">2O-I Alarm></upper>
<lower n<sub="">2O-I Alarm></lower>
<upper agt-e="" alarm="">[*]</upper>
<lower agt-e="" alarm="">*</lower>
<upper agt-i="" alarm="">[*]</upper>
<lower agt-i="" alarm="">*</lower>
<upper alarm="" mac=""></upper>
<upper alarm="" rr=""></upper>
<lower alarm="" rr=""></lower>
<apnea alarm=""></apnea>
<gas check="" trap="" water=""></gas>
<gas check="" conn.="" trap="" water=""></gas>
<gas off="" pump=""></gas>
<gas pump="" regulating=""></gas>
<gas check="" line="" sample=""></gas>
<gas failed="" zeroing=""></gas>
<gas failure="" unit=""></gas>
<gas up="" warm=""></gas>
<gas zeroing=""></gas>
<gas agents="" mixed=""></gas>
<gas cal.="" required.="" zero=""></gas>
<gas cal.="" required.=""></gas>
<gas replace="" trap="" water=""></gas>
<out (co<sub="" of="" range="">2)></out>
<out (rr_co<sub="" of="" range="">2)></out>
<out (n<sub="" of="" range"="">2O)></out>
<out (o<sub="" of="" range="">2)></out>
<out (agent)="" of="" range="">*</out>

*: The selected or detected label will be displayed for the agent label.
SPIRO (MGU-810)

Message
<spiro up="" warm=""></spiro>
<spiro check="" class="" flowsensor=""></spiro>
<spiro check="" flow="" sensor=""></spiro>
<spiro active="" calibration=""></spiro>
<spiro zeroing=""></spiro>
<spiro error="" unit=""></spiro>
<out (tv)="" of="" range=""></out>
<out (mv)="" of="" range=""></out>
<out (press)="" of="" range=""></out>
<upper rr=""></upper>
<lower rr=""></lower>
<apnea></apnea>
<upper mv=""></upper>
<lower mv=""></lower>
<upper peak=""></upper>
<lower peak=""></lower>
<upper peep=""></upper>
<lower peep=""></lower>

BIS (When HBX-800 is used)

Message
<upper alarm="" bis=""></upper>
<lower alarm="" bis=""></lower>
<check sensor=""></check>
<expired sensor=""></expired>
<invalid sensor=""></invalid>
<sensor many="" too="" uses=""></sensor>
<sensor usage=""> 24hrs.></sensor>
<check conn.="" sensor=""></check>
<sensor check="" in="" progress=""></sensor>
<ground check="" in="" progress=""></ground>
<high impedance=""></high>
<artifact></artifact>
<lead off=""></lead>
<bis "sensor="" check"="" perform=""></bis>
<sqi 15%="" <=""></sqi>
<sqi 50%="" <=""></sqi>
<artifacts></artifacts>
<bisx failure=""></bisx>
<bisx incompatible=""></bisx>

Ventilator Alarm Message

Top Priority Alarm (Alarm Level S)

	Item	Message
Ì	Ventilator	<vent. alarm=""></vent.>
	Ventilator	<vent comm=""></vent>

WARNING

- When the VELIA, ASTRAL, VS ULTRA ventilator is connected, and the ventilator power is turned OFF, alarm will not generate on the DS-8400. If the connection cable is disconnected from the ventilator, <Vent. Alarm> will generate, but <VENT_COMM> alarm will not generate on the DS-8400.
- The ventilator alarm sound is set to OFF (factory default).
- The alarm sound can be turned ON on the "Tone/Volume" menu. (Tone/Volume P10-20)

Ventilator Alarm Factor

- For the ventilators other than Servo ventilators, ventilator alarm factor will not be notified to the central monitor.
- Depending on the central monitor type and software version, ventilator alarm factor may not be displayed. For details of the central monitor type and software version, refer to your nearest service representative.
- The ventilator alarm factors are displayed only on the central monitor. These will not be displayed on the bedside monitor.

Displayed Alarm Message	Remarks
VENT AWP	Airway Pressure Alarm
VENT MV	Minute Ventilation Alarm
VENT APNEA	Apnea Alarm
VENT CONT. HP	Continuous High Pressure Alarm
Upper VENT_FiO2	FiO ₂ Upper Limit Alarm
Lower VENT_FiO2	FiO ₂ Lower Limit Alarm
Upper VENT_CO ₂	EtCO ₂ Upper Limit Alarm
Lower VENT_CO ₂	EtCO ₂ Lower Limit Alarm
Upper VENT_RR	RR Upper Limit Alarm
Lower VENT_RR	RR Lower Limit Alarm
VENT_PEEP	PEEP Low Alarm
VENT_COMM	Power OFF, cable disconnected, standby condition, etc.
VENT_URGENT	Other high level alarm
Ventilator	Other ventilator alarm

Cardiac Output Message

Status Message

Message	Details
WAIT	Preparing for measurement. It will be also displayed when catheter relay cable is not connected to the CO module, or when thermodilution catheter is not connected.
READY	Ready to start the measurement.
BUSY	In process of measurement.
END	Measurement is completed.

Result Status

The result status will be displayed for 30 seconds after completion of measurement.

Message	Details
со_ок	CO is correctly measured.
UPPER_FAULT	Measurement error
	After the injection, the blood temperature is out of the measurement range.
	The thermistor connector and relay cable are not securely connected.
	The sensor or relay cable is defective.
PEAK_FAULT	Measurement error
	The peak of the thermodilution curve can not be detected.
	The thermistor connector and relay cable are not securely connected.
	The sensor or relay cable is defective.
LOWER_FAULT	Measurement error
	The blood temperature has not returned to stable condition after the measurement.
	The thermistor connector and relay cable are not securely connected.
	The sensor or relay cable is defective.
SENSOR_ERROR	Measurement error
	The thermistor connector and relay cable are not securely connected.
	The sensor or relay cable is defective.
OVER RANGE	Measurement error
	The CO value is out of the calculation range.

Troubleshooting

This section explains the troubleshooting for each case.

ECG

Check Electrodes> or <LEAD OFF> is displayed.

Cause 1

The electrode is detached, or is not making good electrical contact with the skin.

Solution

Check if the electrodes are properly attached.

Replace the electrodes.

Make sure that the lead cable or relay cable is not defective (wire break, etc.).

(Before Attaching the Electrodes" P7-1)

(Bectrode Placement" P7-2)

Cause 2

Even though the "3-lead Override" (ECG Parameter Setup) is set to [ON], electrodes other than LA, RA, LL are connected.

Solution Set the "3-lead Override" to [OFF]. Or, detach the electrodes other than LA, RA, LL.

ECG Low Amplitude> is displayed.

Cause 1

The ECG amplitude is 0.25 mV or below for the waveform size of x1, x1/2, x1/4, and 0.15 mV or below for the waveform size of x2, x4.

Solution

Change the electrode site, or select a lead with higher QRS amplitude.

NOTE

Using 4-electrode or 5-electrode instead of 3-electrode allows more accurate QRS detection.

Cause 2

The electrode contact is poor.

Electrical blanket or other noise source is near the patient.

Solution

Attach the electrodes firmly. Or, replace the electrodes.

- If the lead cable or relay cable is defective (wire break, etc.), replace it.
- If any noise source is near the patient, move it away from the patient as far as possible.

Cause 3

Even though the "3-lead Override" (ECG Parameter Setup) is set to [ON], electrodes other than LA, RA, LL are

connected.

Solution

Set the "3-lead Override" to [OFF].Or, detach the electrodes other than LA, RA, LL.

ECG Artifact> is displayed.

<u>Cause 1</u>

The electrode contact is poor. Electrical blanket or other noise source is near the patient.

Solution

Attach the electrodes firmly.

- If the lead cable or relay cable is defective (wire break, etc.), replace it.
- If any noise source is near the patient, move it away from the patient as far as possible.

Cause 2

EMG is interfering.

Solution

- Change the electrode site to a location where the myoelectricity will be less likely to interfere.
- Select ESIS for the filter mode.

\land CAUTION

 Selecting ESIS for the filter mode will decrease the QRS amplitude and may result in not counting the heart rate.

Cause 3

Even though the "3-lead Override" (ECG Parameter Setup) is set to [ON], electrodes other than LA, RA, LL are connected.

Solution

Set the "3-lead Override" to [OFF].Or, detach the electrodes other than LA, RA, LL.

The ECG waveform is in the baseline position.

The lead-off condition may have occurred by the following causes.

Cause 1

Electrode is detached.

Solution

Place the electrodes again. If the electrode contact is poor, replace the electrode. (@"Before Attaching the Electrodes" P7-1) (@"Electrode Placement" P7-2)

Cause 2

The lead cable is disconnected from the electrode terminal.

Solution

Securely connect the lead cable.

REFERENCE

• If the error persists, wire break of the lead cable or relay cable can be considered. Contact your nearest service representative.

Check Electrodes Attachment> is displayed.

Cause 1

The electrode contact with the skin is poor. There is substantial contact resistance between the electrodes.

Solution

Replace all the electrodes.Make sure to use the same type of electrodes . (@"Before Attaching the Electrodes" P7-1) (@"Electrode Placement" P7-2)

Cause 2

Even though the "3-lead Override" (ECG Parameter Setup) is set to [ON], electrodes other than LA, RA, LL are connected.

Solution

Set the "3-lead Override" to [OFF].Or, detach the electrodes other than LA, RA, LL.

ECG Unit Error> is displayed.

<u>Cause</u>

A communication error has occurred between the ECG measuring unit.

Solution

A failure of the ECG unit can be considered. Contact your nearest service representative.

The measurement data is displayed as "xxx".

<u>Cause</u>

The heart rate is outside the measurement range. Solution

- Check if the electrodes are properly attached.
 (P"Before Attaching the Electrodes" P7-1)
 (P7-2)
- Replace the electrode, or check the lead cable and relay cable.

Heart rate is not counted. Heart rate is low.

<u>Cause</u>

The ECG waveform amplitude is below the QRS detection level (0.3 mV).

Solution 1

Change the electrode site, or select a lead with higher QRS amplitude.

- Using 4-electrode or 5-electrode instead of 3-electrode allows more accurate QRS detection.
- Also, if large amount of noise is interfering, the noise may be erroneously detected as QRS.Change the electrode site and increase the ECG amplitude.

Solution 2

Increase the displayed waveform size. By increasing the waveform size, small QRS wave will become detectable. However, noise may be also detected.

Heart rate is not counted, and <LEAD OFF> is displayed.

<u>Cause 1</u>

The electrode of the displayed lead type is detached, or is not making good electrical contact with the skin. Solution

- Check if the electrodes are properly attached.
 (P"Before Attaching the Electrodes" P7-1)
 (P7-2)
- Replace the electrode, or check the lead cable and relay cable.

Cause 2

Even though the "3-lead Override" (ECG Parameter Setup) is set to [ON], electrodes other than LA, RA, LL are connected.

Solution

Set the "3-lead Override" to [OFF].Or, detach the electrodes other than LA, RA, LL.

Artificial pacemaker pulse is not displayed.

Cause 1

[Not Used] is selected for "Pacemaker" on the "Admit/Discharge" menu.

Solution

Select [Used] for "Pacemaker".

Cause 2

"Pacemaker Pulse" is set to [OFF] (ECG Parameter Setup).

Solution

Select [ON] for "Pacemaker Pulse" .

Cause 3

The electrode attachment site is not appropriate.

Solution

Check the electrode attachment site. (@"Before Attaching the Electrodes" P7-1) (@"Electrode Placement" P7-2)

ECG Pacing detection error> is displayed.

<u>Cause</u>

The pacemaker pulse is detected 16 pulses or more per second.

Solution 1

- Check if the electrodes are properly attached.
 ("Before Attaching the Electrodes" P7-1)
 ("Electrode Placement" P7-2)
- Replace the electrode, or check the lead cable and relay cable.
- If any noise source is near the patient, move it away from the patient as far as possible.

Solution 2

If the patient is not using a pacemaker, select [Not Used] for "Pacemaker"("Admit/Discharge").

ECG Disconnected> is displayed.

<u>Cause</u>

While monitoring the ECG, the relay cable was unplugged.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the ECG relay cable. The message will disappear, and the alarm will be silenced.

Cannot analyze> is displayed.

<u>Cause</u>

"Suspend Arrhy, Analysis during Noise Interference" ("Initial Settings") is set to ON, and arrhythmia analysis is suspended for more than 30 seconds due to continuous noise or EMG interference.

Solution

Check the electrode attachment, and remove the noise source.

- Check the electrode attachment, lead cable and relay cable.
- If the electrode, lead cable, or relay cable is defective, replace them.
- If any noise source is near the patient, move it away from the patient as far as possible. If EMG is interfering, change the electrode site to a location where EMG will less likely to interfere.

Arrhythmia cannot be detected, or is judged as "?".

Cause 1

The amplitude of ECG1 or ECG2 is below the QRS detection level (250 µV and below).

Solution

Change the electrode site, or select a lead with higher QRS amplitude for both ECG1 and ECG2. When the electrode site is changed, perform the arrhythmia learn process.

Cause 2

The shapes of normal heartbeat and arrhythmia are similar.

Solution

Change the electrode site or select a lead which shows a clear difference between a normal heartbeat and arrhythmia. When the electrode site is changed, perform the arrhythmia learn process.

Cause 3

Noise is interfering with the ECG.

Solution

Check the electrode attachment, and remove the noise source.

- Check the electrode attachment, lead cable and relay cable.
- If the electrode, lead cable, or relay cable is defective, replace them.
- If any noise source is near the patient, move it away from the patient as far as possible. If EMG is interfering, change the electrode site to a location where EMG will less likely to interfere.

Respiration

CVA detected> message is displayed.

<u>Cause</u>

Heartbeat is interfering and superimposed on the respiration waveform.

Solution

Place the electrode as shown below where the heartbeat will be less likely to interfere.

Or, select a lead where the heartbeat will be less likely to interfere.



Q<RR meas. range is exceeded.> message is displayed.

Cause 1

Electrode is detached.

Solution

Reattach the electrode. If the electrode contact is poor, replace the electrode. (@"Before Attaching the Electrodes" P7-1) (@"Electrode Placement" P7-2)

Cause 2

The electrode contact impedance is high.

Solution 1

Reattach the electrode. If the electrode contact is poor, replace the electrode. (
"Before Attaching the Electrodes" P7-1)
(
"Electrode Placement" P7-2)

Solution 2

Change the lead for respiration measurement.

"0" is displayed for respiration rate, or apnea alarm is generated.

Cause

The amplitude of the respiration waveform is too low.

Solution 1

Change the electrode site, or select a lead with higher QRS amplitude.

Solution 2

Increase the displayed waveform size.

The respiration waveform and respiration rate is not displayed.

Cause 1

The electrosurgery-proof type ECG relay cable is used.

Solution

The impedance respiration can not be measured if the electrosurgery-proof type ECG relay cable is used. If not

using during electrosurgery, use the standard ECG relay cable.

Cause 2

The impedance respiration measurement is ceased.

Solution

Select [ON] for "Impedance Measurement" on "Admit/Discharge" or "RESP" setup screen.

 If the pacemaker with the minute ventilation measuring function is used, turn OFF the impedance respiration measurement. Otherwise, both the pacemaker and this monitor will not be able to perform accurate measurement.

The measurement data is displayed as "xxx".

<u>Cause</u>

The respiration rate is outside the measurement range. Solution

- Check if the electrodes are properly attached.
 ("Before Attaching the Electrodes" P7-1)
 ("Electrode Placement" P7-2)
- Replace the electrode, or check the lead cable.
- Change the lead for respiration measurement.

The lead for respiration measurement cannot be changed.

<u>Cause</u>

HLX is used.

Solution

- If HLX is set, the lead will be fixed to [II].
- If the respiration amplitude for lead II is small, check the electrode attachment. (@"Before Attaching the Electrodes" P7-1)
 - (P"Electrode Placement" P7-2)

Invasive Blood Pressure

The PDP value is displayed as "---".

<u>Cause</u>

The BP measured by the HM-800/HM-801Multi Module is labeled as [IAP].

Solution

PDP will not be calculated if the BP measured by the HM-800/HM-801 is labeled as [IAP]. When using the HM-800/HM-801, do not set the BP label to IAP. When monitoring PDP, set the BP label to [IAP].

□<BP* Transducer OFF> is displayed.

<u>Cause</u>

The BP (1 to 8) transducer is not connected.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced.

Solution 2

Connect the transducer.

Solution 3

The BP relay cable or transducer may be defective. Replace the BP relay cable or transducer.

□<BP* Zero Required> is displayed.

<u>Cause</u>

The BP zero balance has not been performed since the power is turned ON.

Solution

Open the three-way valve of the transducer to air and perform zero balance.

The measurement data is displayed as "---".

<u>Cause</u>

The BP zero balance has not been performed since the power is turned ON.

Solution

Open the three-way valve of the transducer to air and perform zero balance.

BP value and waveform are not displayed properly.

<u>Cause</u>

The BP zero-balance is unstable.

Solution 1

Open the three-way valve of the transducer to air and perform zero balance.

Solution 2

Disconnect the BP transducer from the BP relay cable, and check if there is any abnormality on the connector terminal. Make sure that there is no distortion nor substance, such as blood, medicament, attached which may cause contact failure.

If any abnormality is found, replace the BP transducer or BP relay cable.

The measurement data is displayed as "xxx".

<u>Cause</u>

The BP value is outside the measurement range.

Solution

Perform BP zero balance again.

Check if the measurement data is within the measurement range. Check the BP relay cable and BP transducer.

□<BP# Disconnected> is displayed.

<u>Cause</u>

While monitoring the blood pressure, BP relay cable was disconnected from the 2ch BP conversion cable.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the BP interface cable or 2ch BP conversion cable. The message will disappear, and the alarm will be silenced.

The zero balance process fails.

<u>Cause</u>

The three-way valve may not be opened to air, or artifact is present due to movements, etc.

Solution

Check if the three-way valve is opened to air. Verify that <Zero ready> is displayed on the parameter setup screen, or <READY> is displayed on the user key before starting the zero balance.

Contract Contract

Cause 1

The BP relay cable or transducer is defective.

Solution

Replace the BP relay cable or transducer.

Cause 2

The hardware failure has occurred.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

Check the ART catheter.> is displayed.

Cause 1

During the measurement, ART catheter was disconnected.

Solution

Connect the ART catheter securely. Make sure that the ART catheter is not loose.

Cause 2

The BP relay cable or transducer is defective.

Solution

Replace the BP relay cable or transducer.

SpO₂ Measurement (HS-8312N, DS-8007N, HG-820)

\Box < SpO₂ Check Sensor Attach.> is displayed.

Cause

The sensor is detached from the patient. Solution 1 Check if the sensor is properly attached to the patient. Solution 2 Check that the light emitting and receiving parts of the sensor LED are aligned.

\Box < SpO₂ Pulse Search > is displayed.

Cause 1

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly.

Solution

Check that the light emitting and receiving parts of the sensor LED are aligned.

Cause 2

The sensor has not been attached long enough to obtain stable measurement.

Solution

After the sensor attachment, wait and see for about one minute until the waveform stabilizes.

\Box < SpO₂ No Pulse Detected> is displayed.

<u>Cause</u>

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly. Solution

Check that the light emitting and receiving parts of the sensor LED are aligned. Avoid the sensor from exposure to ambient light.

\Box < SpO₂ Motion Artifact> is displayed.

<u>Cause</u>

There is excessive body motion from the patient.

Solution

Relocate the sensor to which body motion will have less influence.

The pulse waveform is not displayed, or interrupted.

Situation: <SpO2 Check Sensor Attach.> is displayed.

Cause 1

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly.

Solution

Check that the light emitting and receiving parts of the sensor LED are aligned.

The sensor is defective. Solution Replace the sensor.

Cause 3

 $\ensuremath{\mathsf{SpO}_2}$ sensor is not firmly connected to the connector.

Solution

Make sure the SpO₂ sensor is firmly connected.

Cause 4

Sensor is exposed to light.

Solution

Place a black or dark cloth over the sensor to avoid direct sunlight. When not using the sensor for measurement, avoid placing the sensor in light or unplug the sensor from the connector.

\Box SpO₂ value is unstable.

<u>Cause 1</u>

There is excessive body motion from the patient which disables correct measurement.

Solution 1

Have the patient lie still.

Solution 2

Relocate the sensor, or change the sensor to which the body motion will have less influence.

Cause 2

The probe size is not appropriate.

Solution

Select a probe size which is appropriate for the patient.

Cause 3

Sensor is exposed to light.

Solution

Place a black or dark cloth over the sensor to avoid direct sunlight.

\Box < Super Unit SpO₂ Failure > is displayed.

<u>Cause 1</u>

The sensor is defective. Solution Replace the sensor.

Cause 2

Communication error has occurred with the SpO2 unit.

Solution

A defective cable or SpO_2 unit failure can be considered. Contact your nearest service representative.

The system was started with the sensor and cable connected.

Solution

Disconnect the SpO_2 cable and sensor from this equipment, and press the standby switch to enter into standby mode. Then, press the standby switch again to cancel the standby mode, and when the monitoring screen is displayed, connect the cable and sensor.

 \Box < SpO₂ Replace Sensor> is displayed.

<u>Cause 1</u>

The sensor is not connected securely.

Solution

Connect the sensor securely.

<u>Cause 2</u> The sensor is defective. Solution

Replace the sensor.

Cause 3

A wrong sensor is used. Solution Replace the sensor. For details of the usable sensors, refer to your nearest service representative.

 \Box < SpO₂ Disconnected> is displayed.

<u>Cause</u>

The SpO₂ relay cable is disconnected during SpO₂ monitoring.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the SpO₂ relay cable. The message will disappear, and the alarm will be silenced.

SpO2 Measurement (HS-8312M, DS-8007M, HG-810)

\Box <SpO₂ Replace Sensor> is displayed.

<u>Cause 1</u> The sensor is not connected securely. Solution Connect the sensor securely. <u>Cause 2</u>

The sensor is defective. Solution Replace the sensor.

A wrong sensor is used.

Solution

Replace the sensor.

(@"Pulse Oximetry Measurement (Manufactured by Masimo)" P13-4)

<u>Cause 4</u>

The sensor is used beyond its expected life.

Solution

Replace the sensor.

(NOTE

- The X-Cal function automatically monitors the Active Monitoring Time (actual time of SpO₂ monitoring) for each sensor and cable.
- Even if the sensor is used beyond its expected life, the measurement will not cease unless the power is turned OFF, sensor is disconnected from the cable, cable is disconnected from the monitor, or the sensor is reattached.
- When a measurement with a sensor that has reached its end of life is suspended for certain amount of time, and resumed with the same sensor, a message to replace the sensor will be displayed.
- Depending on the equipment, some sensors may not be recognized.

\Box < SpO₂ Check Sensor Attach.> is displayed.

Cause 1

The sensor is detached from the patient.

Solution 1

Check if the sensor is properly attached to the patient.

Solution 2

Check that the light emitting and receiving parts of the sensor LED are aligned.

Cause 2

The sensor is exposed to too much ambient light. The detecting part of the sensor is not covered appropriately.

Solution 1

Turn down or turn off the light.

Solution 2

Avoid the sensor from exposure to ambient light.

Solution 3

Relocate the sensor position.

\Box < SpO₂ Low Perfusion> is displayed.

<u>Cause</u>

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly.

Solution

Check that the light emitting and receiving parts of the sensor LED are aligned.

Confidence> is displayed.

<u>Cause</u>

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly.

Solution

Check that the light emitting and receiving parts of the sensor LED are aligned.

 \Box < SpO₂ Pulse Search > is displayed.

Cause 1

The amplitude of the pulse waveform is too low, or the sensor is not positioned correctly.

Solution

Check that the light emitting and receiving parts of the sensor LED are aligned.

Cause 2

The sensor has not been attached long enough to obtain stable measurement.

Solution

After the sensor attachment, wait and see for about one minute until the waveform stabilizes.

\Box < SpO₂ Noise Interference> is displayed.

<u>Cause</u>

External signal or energy is interfering with the measurement.

Solution

Remove the external interference or apply ambient shielding.

 \Box <SpO₂ Check Sensor>, <SpO₂ Replace Cable>, or <SpO₂ Check Cable> is displayed.

Cause 1

Unrecognizable sensor is connected.

A wrong patient cable is used.

When attached to the patient, the sensor was exposed to high-intensity light which lead to false recognition.

Solution

Reattach the SpO₂ sensor and patient cable.

Replace with a Fukuda Denshi specified patient cable and sensor.

(@"Pulse Oximetry Measurement (Manufactured by Masimo)" P13-4)

Cause 2

The cable is used beyond its expected life.

Solution

Replace the patient cable.

NOTE

- The X-Cal function automatically monitors the Active Monitoring Time (actual time of SpO₂ monitoring) for each sensor and cable.
- Even if the cable is used beyond its expected life, the measurement will not cease unless the power is turned OFF or the cable is reconnected.
- When a measurement with a cable that has reached its end of life is suspended for certain amount of time, and resumed with the same cable, a message to replace the cable will be displayed.

· Depending on the equipment, some cable may not be recognized.

 \Box < Super Unit SpO₂ Failure > is displayed.

<u>Cause</u>

Communication error has occurred with the SpO2 unit.

Solution

A defective cable or SpO2 unit failure can be considered. Contact your nearest service representative.

\Box < SpO₂ Disconnected > is displayed.

<u>Cause</u>

The SpO₂ relay cable is disconnected during SpO₂ monitoring.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the SpO₂ relay cable. The message will disappear, and the alarm will be silenced.

\Box < SpO₂ only mode> is displayed.

<u>Cause</u>

When the Rainbow sensor is used, SpCO, SpMet or SpHb parameter cannot be measured.

Solution 1

Remove the sensor from the patient's finger, and then reattach it.

Solution 2

Remove the sensor or patient cable from the HS-8000 or module, and then reconnect it to the SpO₂ connector.

Comparison Comparis

<u>Cause</u>

There is excessive body motion, or sensor attached position is not appropriate.

Solution 1

Check that the light emitting and receiving parts of the sensor LED are aligned.

Solution 2

Relocate the sensor to which body motion will have less influence.

PVI, SpCO, SpMet, SpHb, SpOC cannot be measured.

Cause 1

PVI, SpCO, SpMet, SpHb, SpOC measurements are optional functions.

Solution

It is necessary to add these as the measuring parameters. For details, contact your nearest service representative.

Cause 2

The used sensor cannot measure the PVI, SpCO, SpMet, SpHb, SpOC.

Solution

Use the sensor which can measure the PVI, SpCO, SpMet, SpHb, SpOC. For details, contact your nearest service representative.

Non-Invasive Blood Pressure

The cuff is not inflated although the pump is operating.

<u>Cause 1</u>

The air hose is not firmly connected, and the air is leaking.

Solution

Check if the air hose is properly connected.

Cause 2

The cuff size does not match the selected patient type.

Solution

Use the cuff with correct size for the selected patient type.

The pump is not operating.

<u>Cause</u>

The air hose is disconnected from the NIBP connector. Solution

Check if the air hose is properly connected.

The measurement data is displayed as "---".

Cause 1

The measurement accuracy is not reliable due to body motion artifact.

Solution

During the measurement, have the patient stay still.

Cause 2

The pulse is too small to acquire reliable measurement accuracy.

Solution

Check if the cuff application is proper, and if the cuff size corresponds with the selected patient type.

Cause 3

The air hose is disconnected.

Solution

Check if the air hose is tightly connected, and then measure again. If the same message is displayed again, air leakage inside the HS-8000/DS-8007 can be considered. Contact your nearest service representative.

Check NIBP cuff, hose> is displayed.

Cause 1

The connection between the cuff and air hose or the air hose and NIBP connector is loose or disconnected. Solution

If the connection is loose or disconnected, securely connect it and perform the measurement again.

If the same message is displayed again, internal air leakage can be considered. Cease the measurement, and contact your nearest service representative.

The cuff is compressed.

Solution

Make sure that the cuff is not subjected to compression, and measure with the patient arms extended as much as possible.

If the same message is repeatedly displayed, air system may be clogged. Cease the measurement, and contact your nearest service representative.

Cause 3

The cuff size is not suitable for the patient.

Solution

Check that the cuff size is appropriate for the patient, and that the cuff is properly attached, and measure again.

Cause 4

The cuff size and the patient classification setting do not match.

Solution

Make sure that the appropriate cuff size is used according to the patient classification setting.

 \Box <NIBP measurement failed (Cxx-xx)> is displayed.

Error code condition (phenomenon, or situation) and its cause are indicated below.

C02-00 When "Quick Measurement" is [OFF], the data could not be measured.

Cause 1

The blood pressure may not be correctly measured due to the patient's condition.

Solution

Check the patient's condition, and measure again.

Cause 2

The cuff application has become loose.

Solution

Check that the cuff size is appropriate for the patient, and then measure again after attaching the cuff properly.

C02-01 When "Quick Measurement" is [ON], the data could not be measured.

Cause 1

The blood pressure may not be correctly measured due to the patient's condition.

Solution

Check the patient's condition, set "Quick Measurement" to OFF, and measure again.

Cause 2

The cuff application has become loose.

Solution

Check that the cuff size is appropriate for the patient, and then measure again after attaching the cuff properly.

C02-02 The air hose was disconnected from the NIBP connector during the measurement.

<u>Cause</u>

The air hose was disconnected from the NIBP connector during the measurement.

Solution

Connect the air hose to the NIBP connector, and then measure again.

C03-xx The exhaust ventilation has ceased, or the target deflation speed was not achieved.

Cause 1

During measurement, an artifact such as body motion may have interfered.

Solution

Keep the patient still as much as possible, and measure while the patient is not moving. When performing the measurement during surgery, avoid artifact caused by the surgery.

Cause 2

During the measurement, air hose was bent or occluded by the compression.

Solution

Make sure that the air hose is not bent or compressed before the measurement.

If the error persists and C03-xx error is frequently displayed, contact your nearest service representative and notify the error code.

C04-xx The cuff inflation was insufficient for the patient's blood pressure.

<u>Cause</u>

The blood pressure has significantly increased from the previous measurement.

Solution

Check the cuff application and size and perform the manual measurement.

C06-xx The pulse signal detected during the measurement was unstable.

During the measurement, the patient has trembled or moved.

Solution

Keep the patient still as much as possible, and measure while the patient is not trembling or moving.

Cause 2

Arrhythmia has frequently occurred during the measurement.

Solution

If arrhythmia occurs many times, correct measurement cannot be performed. Measure when arrhythmia is not frequently occurring.

C07-00 The measurement time has exceeded the allowable time.

<u>Cause</u>

Measurement is automatically repeated due to body motion or insufficient inflation.

Solution

Check the cuff application and size, and measure while keeping the patient still as much as possible.

Cause 1

C08-00 The detected PR value was abnormal.

<u>Cause</u>

The patient has trembled or moved.

Solution

Keep the patient still as much as possible, and measure while the patient is not moving.

C09-00 The inflation value has exceeded the allowable maximum value.

<u>Cause</u>

The cuff was subjected to compression.

Solution

Make sure that the cuff is not subjected to compression, and measure with the patient arms extended as much as possible.

C10-xx The detected pulse amplitude was abnormal.

<u>Cause</u>

The cuff size is not suitable for the patient.

Solution

Check that the cuff size is appropriate for the patient, and that the cuff is properly attached, and measure again.

The time of measurement disappears and the numeric data is displayed as " - - - ".

<u>Cause</u>

The preprogrammed time to clear the NIBP data has elapsed.

Solution

The "NIBP Erase Time" can be selected from [60 min.], [120 min.], and after the set duration, the NIBP data will be displayed as "---".

Select the appropriate time which best fits the monitoring purpose.

The NIBP periodic measurement is ceased.

Cause

<NIBP Meas. Error (Exx-xx)> is displayed during the measurement.

Solution

When <NIBP Meas. Error (Exx-xx)> is displayed, the NIBP periodic measurement will be canceled. To resume the measurement, press the [NIBP Start/Stop] key and check that the measurement is properly performed.

\square <NIBP Unit Error (E^{**}-^{**})> is displayed on the main unit.

<u>Cause</u>

An error has occurred on the NIBP unit.

E08-01: Communication Error (Sub CPU)

E08-02: WatchDog Timeout

E08-03: Pressure Offset Error

E08-04: Pressure Comparison Error

E08-05: Sub CPU Power Supply Failure

E08-06: Pressure Sensor 2 Power Supply Failure

E08-07: Pressure Sensor 1 A/D Reference Power Voltage Failure

E08-08: Rapid Exhaust Error

E08-09: Air Hose Identification Error E09-A: Exceeded Maximum Cuff Pressure E09-B: Inflation Timeout E09-C: Quick Mode Timeout E09-D: Measurement started during the long pause E09-E: Measurement Timeout E09-F: Main CPU Pressure Data Transmission Timeout E09-G: Pressure Sensor 1 +5V Power Supply Failure E09-H: Zero Calibration Timeout E09-I: ROM Test Error E09-J: RAM Test Error E09-L: Clock Transmission Ceased E09-M: Communication Failure at Power ON E09-N: Pressure Comparison Error E09-O: Maximum Inflation Timeout E09-Q: Measurement was started before zero calibration E09-R: Zeroing Error E09-S: WatchDog Timeout E09-T: +5V Digital Power Supply Failure E09-U: Main CPU Power Supply Failure E09-V: Pump Control Signal Failure E09-W: Quick Exhaust Valve Control Signal Failure E09-X: Sub CPU Constant Exhaust Valve Control Signal Failure E09-Y: Main CPU Constant Exhaust Valve Control Signal Failure

Solution 1

These errors can be cleared by pressing the [Cancel Error] on the NIBP setup menu or [NIBP Start/Stop] key (fixed key or user key). If the same message is repeatedly displayed, a failure of the equipment can be considered. Cease the measurement, and contact your nearest service representative.

Solution 2

When <NIBP Unit Error (Exx-xx)> is displayed, make sure that the congestion is not generated, and remove the cuff if necessary.

Temperature

T* Unknown Sensor> is displayed.

Cause 1

700 series temperature probe is used.

Solution

Use the 400 series temperature probe for measurement.

Cause 2

There is a contact failure of the temperature probe.

Solution

Check if the temperature probe is properly inserted.

The measurement data is displayed as "xxx".

<u>Cause</u>

The temperature measurement is outside the measurement range.

Solution

Check if the temperature probe is properly inserted. Replace the temperature probe, or check the temperature probe.

\Box <T* Disconnected> is displayed.

<u>Cause</u>

While monitoring the temperature, the temperature probe was unplugged.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the temperature probe. The message will disappear, and the alarm will be silenced.

□<Super Unit TEMP Unit Failure> is displayed.

<u>Cause</u>

An error was detected on the temperature unit.

Solution

A unit failure can be considered. Cease the measurement, and contact your nearest service representative.

Cardiac Output (CO)

When measured consecutively, the measurement value varies. (±10% or more)

Cause 1

The injection method is not appropriate. Solution Inject within 1 to 3 seconds.

Cause 2

Injection temperature is not appropriate.

Solution

If iced injectate is used, pay attention not to warm the injector with hands.

Cause 3

The thermistor location is not appropriate. Solution Reposition the thermistor.

<u>Cause 4</u> Arrhythmia event has occurred during the measurement. Solution Wait until the patient has stable heart rhythm.

There was patient's body movement during the measurement.

Solution

Have the patient stay still during the measurement.

Cause 6

The patient's hemodynamics changed during the measurement. Solution Wait until the patient has stable hemodynamics.

Abnormal measurement value is displayed.

<u>Cause</u> The catheter size, injectate volume, catheter constant (CC) is not correct. Solution Set the proper condition, CC value for the used catheter.

The blood temperature (Tb), injectate temperature (Ti) is not displayed.

Cause
The catheter is not properly connected.
Solution
Securely connect the catheter.

The thermodilution curve is deformed.

<u>Cause</u> The injection is not smooth, steady motion. Solution Inject promptly within 1 to 3 seconds.

The baseline of the thermodilution curve is displaced to the minus side. <LOWER FAULT> is displayed.

<u>Cause</u>

The blood temperature has not returned to stable condition after the measurement.



The thermodilution curve did not return to the cut off point soon enough. The temperature must return to a point that is 30% of the peak value within 30 seconds (or 60 seconds depending on the setup).

Solution

If performing continuous measurement, wait for 30 to 60 seconds and check that "Ready" is displayed before

performing the next measurement.

The thermodilution curve is low. <PEAK FAULT> is displayed.

<u>Cause</u>

The peak of the thermodilution curve can not be detected.



After the measurement is started, the peak of the thermodilution curve was not determined within 22 seconds (when the time scale is "30 sec") or 45 seconds (when the time scale is "60 sec").

Solution

The thermistor may be contacting the pulmonary artery wall. Reposition the thermistor and measure again.

□<UPPER FAULT> message is displayed.

<u>Cause</u>

After the injection, the blood temperature is out of the measurement range.



After the measurement is started, the change in blood temperature is less than $0.1^{\circ}C / 1.8^{\circ}F$ for more than 15 seconds (when the time scale is "30 sec") or 30 seconds (when the time scale is "60 sec").

Solution

Use the iced injectate, and measure again.

OVER RANGE> is displayed.

<u>Cause</u>

The CO value is out of the calculation range.

Solution

The area of the thermodilution curve is too large to calculate. Start the measurement again.

The measurement is interrupted, and the error message, <UPPER_FAULT>, <PEAK_FAULT>, <LOWER FAULT>, <SENSOR ERROR> is displayed.

<u>Cause 1</u>

The thermistor connector and relay cable is not securely connected.

Solution

Correct measurement cannot be performed unless the thermistor connector and relay cable is securely connected. Check the connection and perform the measurement again.

Cause 2

The sensor or relay cable is defective.

Solution

If the sensor or cable is defective, measurement can not be performed. Replace the sensor or cable and perform the measurement again.

CO Disconnected> message is displayed.

<u>Cause</u>

The catheter relay cable was disconnected while monitoring the cardiac output.

Solution 1

To cease monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced. Solution 2

To continue monitoring, plug in the catheter relay cable. This will clear the message and silence the alarm.

CO₂ Measurement (HCP-800/HCP-810/HCP-820)

\Box < CO₂ Check Sample Line> is displayed.

<u>Cause 1</u> The sampling tube is clogged. Solution Replace the sampling tube. <u>Cause 2</u> The sampling line is bent or pinched.

Solution

Make sure that the sampling line is properly allocated.

Initializing> displayed inside the numeric data box does not disappear.

<u>Cause</u>

An error has occurred during the initialization at power ON.

Solution

Reconnect the cable of HCP-800/HCP-810/HCP-820 and reboot.

If the message is still displayed, CO₂ unit failure can be considered. Contact your nearest service representative.

\Box < CO₂ Unit Error> is displayed.

<u>Cause</u>

Communication error has occurred with the CO_2 unit.

Solution

A cable disconnection or CO₂ unit failure can be considered. Contact your nearest service representative.

There is substantial measurement error.

Cause 1

20 minutes have not yet elapsed since the power is turned ON.

Solution

For 20 minutes from turning ON the power, there will be a substantial measurement error.

Cause 2

The CO_2 calibration value is not appropriate. Solution Perform the CO_2 calibration again.

\Box <CO₂ Disconnected> is displayed.

<u>Cause</u>

When the filter line is disconnected during CO2 monitoring, this message will be displayed.

Solution 1

To cease monitoring, press the [Alarm Silence] key to clear the message and silence the alarm.

Solution 2

To continue monitoring, plug in the filter line. This will clear the message and silence the alarm.

Waveforms and measurement data are not displayed.

<u>Cause</u>

The gas module is used at the same time.

Solution

The HCP-800/HCP-810/HCP-820 and gas module cannot be used at the same time. If used, CO_2 measurement by the gas module will be prioritized.

CO₂ Measurement (HPD-800/HPD-810/HPD-820)

\Box < CO₂ Sensor Failure > is displayed.

Cause 1

The CO₂ sensor temperature has increased above 40°C/104°F.

Solution

Remove any heat generating source around the sensor.

The CO₂ sensor is malfunctioning.

Solution 1

Replace the CO₂ sensor.

Solution 2

If the error persists, the failure of HPD-800/HPD-810/HPD-820 can be considered. Stop using the unit and contact our service representative.

 \Box <Zero the CO₂ Adapter> is displayed.

<u>Cause</u>

The CO_2 sensor is not zero balanced. Solution Perform the zero calibration of the sensor. ($rac{reg}$ "CO2 Concentration (Mainstream Method)" P7-71)

 \Box < Check CO₂ Airway Adapter> is displayed.

Cause 1

The airway adapter is unclean.

Solution

A clean airway adapter must be used. If reusing an airway adapter, clean and air-dry it. Then, wipe the window with swab, and sterilize (EOG, etc.) before use.

Cause 2

The airway adapter is disconnected from the sensor.

Solution 1

Securely connect the airway adapter to the sensor.

Solution 2

If error persists, perform the airway adapter calibration again.

 \Box <Unknown CO₂ Sensor> is displayed.

\Box < CO₂ Disconnected> is displayed.

<u>Cause</u>

When the cable is disconnected during CO_2 monitoring, this message will be displayed.

Solution 1

To cease monitoring, press the [Alarm Silence] key to clear the message and silence the alarm.

Solution 2

To continue monitoring, plug in the cable. This will clear the message and silence the alarm.

Recorder Unit (HR-800)

Check Paper> is displayed, and printing cannot be performed. The power supply LED on the HR-800 is lit in orange.

The power supply LED on the HR-600 is in in orange.

<PAPER OUT> is displayed inside the [Print Start/Stop] user key.

<u>Cause</u>

There is no paper in the printer.

Solution

Set the paper in the paper holder.

Check Cassette> is displayed, and printing cannot be performed. The power supply LED on the HR-800 is lit in orange.

<CASSETTE> is displayed inside the [Print Start/Stop] user key.

<u>Cause</u> The paper holder is open. Solution Firmly close the paper holder.

Although the paper is fed, printing is not performed.

<u>Cause</u>

The paper is not correctly installed. The front and backside of the paper is set oppositely.

Solution

Set the paper in the paper holder so that the logo, FUKUDA DENSHI CO., LTD appears on the upper surface.

The second and third waveforms are not printed for manual printing or alarm printing.

<u>Cause</u>

The second and third waveforms are not set on the printing setup screen.

Solution

Set the second and third waveform on the corresponding printing setup screen.

The power supply LED on the HR-800 is lit in orange, and [Print Start/Stop] key does not function.

Cause 1 The U-LINK setting is incorrect. Solution Press the [Initial Settings]> [External Device]> [U-LINK] keys. If HR-800 is connected via MGU-800, select [MGU-800]. If HR-800 is not connected via MGU-800, select [OFF]. Cause 2 The HR-800 setting is incorrect. Solution Press the [Menu]>[Initial Settings]>[System]>[Other]>[HR-800] keys. If the HR-800 is installed to the internal slot, select [Built-in].

If connected to U-LINK, select [U-LINK].

Check Printer> is displayed and printing cannot be performed. The power supply LED on the HR-800 is lit in orange. <CHECK?> is displayed inside the [Print Start/Stop] user key.

<u>Cause 1</u> The paper is jammed. Solution Open the paper holder and properly set the paper.

Cause 2

The thermal head temperature has increased or other failure exists.

Solution

Damage to the thermal head or other failure can be considered. Contact your nearest service representative.

Network Printer

Central Printer Check Connection> is displayed and printing cannot be performed.

Cause

The central monitor selected as the output destination is not connected to the printer.

Solution

Check the printer setting on the central monitor, and make sure the communication with the printer is established.

Central Printer Check Setting> is displayed and printing cannot be performed.

<u>Cause</u>

The central monitor selected as the output destination does not support the network printing function. Or, the printer setting is set to [OFF] on the central monitor selected as the output destination.

Solution

Use the DS-7700/DS-7700W system with the software version from V06, and set the printer setting to [ON].

□<Check Central ID> is displayed and printing cannot be performed.

<u>Cause</u>

The central monitor selected as the output destination does not support the network printing function.

Solution

Select the central monitor which supports the network printing function.

Wired Network (DS-LANII/ DS-LANIII)

The data is not displayed on the central monitor.

Cause 1

The DS-LAN setup is not correct.

Solution

Make sure that the DS-LAN Setup (DS-LANII/DS-LANIII) for all bedside monitors and central monitors in the same network are the same. If the DS-LAN setting is changed, make sure to restart the system.

A central monitor which is not compatible is used.

Solution

The following central monitors can not be used on the DS-LANIII network.

- DS-5700
- DS-5800N/NX/NX^{MB}
- DS-7600/7600W with software version V05 and prior

When using these central monitors, all monitors in the same network should be set to DS-LANII.

Cause 3

Inappropriate HUB is used.

Solution

For the DS-LANII network, use the repeater HUB. For the DS-LANIII network, use the switching HUB.

Cause 4

The bed ID is duplicated in the same network.

Solution

If bedside monitors with the same bed ID exist in the same network, communication is not possible. Make sure to set a unique bed ID for each bedside monitor.

Cause 5

An equipment not specified by Fukuda Denshi is connected to the network.

Solution

Do not connect PC, printer, or other unspecified equipment to the DS-LAN network.

Cause 6

The DS-LAN cable is not properly connected.

Solution

The DS-LAN connection will be performed by our service representative. Contact our service representative.

The CO₂ waveform is not displayed on the central monitor although the CO₂ numeric data is displayed.

Cause 1

[Impedance] is selected for "RR/APNEA Alarm Source" on the RESP setup screen.

Cause 2

[Ventilator] is selected for "RR/APNEA Alarm Source" on the RESP setup screen.

Solution

Select [CO₂] for "RR/APNEA Alarm Source" on the RESP setup screen. In this case, RR and apnea alarm will be generated based on CO₂ measurement.

The impedance respiration waveform is not displayed on the central monitor although the RR numeric data is displayed.

Cause 1

[CO2] is selected for "RR/APNEA Alarm Source" on the RESP setup screen.

 $\label{eq:constraint} \ensuremath{\left[\mathsf{Ventilator} \right]}\xspace is selected for "\mathsf{RR}/\mathsf{APNEA}\xspace Alarm Source" on the \mathsf{RESP}\xspace setup screen.$

Solution

Select [Impedance] for "RR/APNEA Alarm Source" on the RESP setup screen.

\square	NOTE
•	The impedance respiration waveform will not be displayed if [CO ₂] is set for "RR/APNEA
	Alarm Source". AWF, AWP waveform will be displayed.

- The CO₂ waveform will not be displayed if [Impedance] is set for "RR/APNEA Alarm Source". AWF, AWP waveform will be displayed.
- The CO₂ waveform and impedance waveform will not be displayed if [Vent.] is set for "RR/APNEA Alarm Source" .

Check DS-LAN Comm> is displayed.

Cause 1

The LAN cable is loose, or contact failure has occurred. The power of the central monitor has been turned OFF. Solution

Check the LAN connection on both the main unit and wall side. Disconnect and connect it again to make sure that it is firmly connected

Check the LAN connection on the central monitor. Disconnect and connect it again to make sure that it is firmly connected.

Turn ON the power of the central monitor.

Telemeter (HLX-801)

The data cannot be received at the telemetry center.

Cause 1

The channel ID or group ID is not corresponded with the telemetry receiver.

Solution

Set the correct channel ID and group ID.

Cause 2

The HLX connection setting is incorrect.

Solution

Press the [Menu]>[Initial Settings]>[System]>[Other]>[HLX Connection] keys.

If connected to internal port, select [Internal Port].

If connected to COM port, select [COM Port].

After the setting, make sure that the set channel is displayed as "CHxxxx" at the upper left of the screen.

The impedance respiration waveform cannot be received at the telemetry center.

Cause 1

[CO₂] is selected for "RR/APNEA Alarm Source" on the RESP setup menu.

[Ventilator] is selected for "RR/APNEA Alarm Source" on the RESP setup screen.

Solution

Select [Impedance] for "RR/APNEA Alarm Source" on the RESP setup menu.

The BP waveform of 100 mmHg and above cannot be properly received.

<u>Cause</u>

The BP waveform and scale are not the same.

Solution

When the BP waveform is above 100 mmHg, set the BP scale above 100 mmHg.

Check HLX Conn.> is displayed.

<u>Cause</u>

The connection with the HLX is interrupted.

Solution

Check the connection between the HLX and DSC-8410. Check if [HLX] is set for the corresponding port under [Initial Settings] > [External Device] > [Main Unit HP-800].

□<HLX Ver.> is displayed.

<u>Cause</u>

Installation has failed.

Solution

Check the software version of the HLX.

If "HLX-801 V99-99" is displayed, perform the installation again.

If the software version of "HLX-501 V01-09" or older is displayed, contact your nearest service representative.

Remote Control

The remote control does not function.

<u>Cause 1</u> The remote control bed ID is not correct.

Solution

Set the correct remote control ID.

Cause 2

The section number is not correct.

Solution

Set the correct section number.

The remote control does not properly function.

<u>Cause</u>

The remote control setting on the monitor does not correspond to the function key on the remote control unit.

Solution

Make sure the remote control setting on the monitor and the function key on the remote control unit is corresponded.

General

Even though the numeric data displayed on the extended display unit or central monitor is exceeding the alarm limit, alarm does not generate.

<u>Cause</u>

The parameters not displayed on the display unit (LC-8016TC/8018TC) are displayed on the central monitor/ extended display unit as [All Data] is selected for "Numeric Data External Output" under [Initial Settings] > [System] > [Other].

Solution 1

For the parameters which requires alarm monitoring on the extended display unit/central monitor, make sure to display those on the display unit (LC-8016TC/8018TC).

Solution 2

For the extended display unit/central monitor, if monitoring is necessary for only the parameters displayed on the display unit (LX-8016TC/8018TC), select [Displayed Data] for "Numeric Data External Output" under [Initial Settings] > [Cystem] > [Other].

Nothing is displayed on the screen, and the power supply LED is not lit.

Cause 1

The display unit is not properly attached to the main unit.

Solution

Connect the apparatus correctly. (@Maintenance Manual "System Construction" P1-2)

Cause 2

The main unit or display unit is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

The power supply LED and charge LED lights alternately.

Cause 1

The display unit is not properly attached to the main unit.

Solution

Connect them correctly. (@Maintenance Manual "System Construction" P1-2)

Cause 2

The main unit or display unit is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

The data is initialized each time the power is turned ON.

Cause 1

The internal switch setting is incorrect.

Solution

The internal switch setting needs to be changed. Contact your nearest service representative.

Cause 2

The battery for the backup memory is depleted.

Solution

The battery needs to be replaced. Contact your nearest service representative.

The display is dark, or cannot be seen clearly.

<u>Cause 1</u> The night mode is set. Solution Cancel the night mode.

Cause 2

The service life of the LCD backlight has expired.

Solution

The LCD unit needs to be replaced. Contact your nearest service representative.

 The display unit utilizes LED for the backlight. Since this LED deteriorates by the life cycle, the display may become dark, scintillate, or may not light by the long term use. In such case, contact your nearest service representative.

The system does not start although the power supply cable is connected.

Cause 1

Incorrect CF card is inserted.

Solution

Remove the CF card, turn OFF the power, and turn ON the power again.

Cause 2

The main unit or display unit is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

The clock is often delayed.

<u>Cause</u>

The battery for the backup memory is depleted.

Solution

Check if the time is delayed when the power is turned OFF.

The battery needs to be replaced. Contact your nearest service representative.
The touch panel key does not function properly.

Cause 1

Due to change in installation environment, the detecting location is misaligned.

Solution

Calibration is required. Contact your nearest service representative.

• The calibration will be performed by our service representative.Users should not perform this procedure as incorrect calibration may cause malfunction of the equipment.

Cause 2

The LCD unit is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□<DSC-8410 Failure> or <DSC-8410 Check Unit> is displayed.

<u>Cause</u>

The hardware failure has occurred.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

Content of Content

<u>Cause</u>

The main unit or display unit is used outside the specified environment condition.

Solution

Use the equipment in the specified environment condition (10°C to 40°C).

□<Display Unit Failure>, <Check Display Unit> is displayed.

<u>Cause</u>

The display unit failure has occurred.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□<DSC-8410 Check Rotary SW> is displayed.

<u>Cause</u>

The rotary switch setting is incorrect.

Solution

If the rotary switch is not set to "0", the equipment will not function properly. Immediately turn OFF the power and cease the operation. Contact your nearest service representative. Control Con

<u>Cause</u>

Rebooting of the system is required.

Solution

Power cycle the system. If the same message is repeatedly displayed, turn OFF the power and contact your nearest service representative.

□<DSC-8410 Check Short-Term Battery>, <DSC-8410 Check Long-Term Battery> is displayed.

<u>Cause</u>

The battery is depleted or malfunctioning.

Solution

The battery needs to be replaced. Contact your nearest service representative.

Some parameters are not displayed due to the display layout setting.> is displayed.

Cause 1

The measured parameter is not set to be displayed.

Solution

On the "Display Config." setting, select the measured parameter to be displayed.

Cause 2

During auto display configuration, the quantity of measured parameters exceeded the displayable parameters. Solution

If there are parameters which measurements are not actually performed, please disconnect their probes/cables.

Check Equip. Config.> is displayed.

Cause 1

The "Multiamplifier" setting does not correspond to the connected cable.

Solution

Check the "Multiamplifier" setting (Initial Settings>System>Unit Module>Multiamplifier), and make sure that the setting corresponds to the connected cable.

Cause 2

On the "External Device" setting, the set external device is duplicated.

Solution

Check the "External Device" setting, and make sure that the selected external device is not duplicated. The external devices other than Vigilance, INVOS, BIS cannot be duplicated. The combinations of FLOW-i and MGU-800/MGU-810, FLOW-i and ventilator are not possible.

Check Module-LAN Comm> is displayed.

<u>Cause</u>

The connection of the module-LAN connector on the main unit or external device is not secured.

Solution

Securely connect the cable to the module-LAN connector.

Securely connect the cable to the external device.

If the error persists, contact your nearest service representative.

Super Unit

The system does not start although the power is turned ON.

The power supply LED on the Super Unit does not light in green. <Super Unit Check Conn.> is displayed.

Cause 1

The power cable of the main unit is not connected.

Solution

Turn OFF the power and connect the power cable.

Cause 2

The module connection cable is not properly connected.

Solution

Turn OFF the power, and make sure that the module connection cable is securely connected. If the knob is loose, tighten it securely.

Cause 3

The fuse inside the Super Unit has blown out.

Solution

Immediately turn off the power and stop using the device. Contact your nearest service representative.

Cause 4

The Super Unit is not properly connected to the adapter (HSA-80/HSA-81 or DSA-82).

Solution

Insert the Super Unit into the HSA-80/HSA-81 or DSA-82 until a click sound is heard.

Super Unit Out of Operating Temp. Range> is displayed.

<u>Cause</u>

The temperature inside the Super Unit has exceeded the operating temperature range.

Solution

The operation cannot be guaranteed.Immediately turn off the power and stop using the device. Contact your nearest service representative.

Super Unit Analog Unadjusted> is displayed.

Cause

One of ECG, respiration, or BP is not adjusted.

Solution

Parameter cannot be measured properly in this situation. Contact your nearest service representative.

□<Super Unit Check DIP-SW> is displayed.

<u>Cause</u>

The DIP switch setting has been changed.

Solution

Contact your nearest service representative.

□<Super Unit Check SD Card> is displayed.

<u>Cause</u>

The SD Card is defective or the Super Unit is malfunctioning. Solution

Contact your nearest service representative.

Data Transfer Function

The patient name is flashing.

<u>Cause</u>

This is a normal operation which indicates the data updating process.

An error occurs during the data update process.

<u>Cause</u>

The HS-8000 is disconnected during the data update process.

Solution 1

Do not disconnect the HS-8000 during the data update process. If the same error persists, refer to your nearest service representative.

Solution 2

If the error occurs during the write process on the DS-8400 System, start again from the read process on the original patient monitor. If the same error persists, refer to your nearest service representative.

Cause 2

The module connection cable is not properly connected.

Solution

Turn OFF the power, and make sure that the module connection cable is securely connected. Reconnect the cable if necessary. If the knob is loose, tighten it securely.

When the HS-8000/DS-8007 is connected, the alarm sound is suspended.

<u>Cause</u>

This is a normal operation. To not suspend the alarm sound, set the alarm sound suspend function OFF.

The recall data cannot be transferred.

Cause 1

The SD card is not inserted to the Super Unit. Solution

Insert the SD card to the Super Unit, and format it.

<u>Cause 2</u> The SD card is not formatted. Solution Format the SD card.

The upload process does not start.

<u>Cause 1</u>

[Transport] is not selected for "Data Transfer" under [Initial Settings > System > Other].

Solution

Check if [Transport] is selected for "Data Transfer" under [Initial Settings > System > Other]. If not, select [Transport].

Cause 2

On the central monitor, "Data Transfer" function is set to [OFF].

Solution

Check the setting on the central monitor.

NOTE
For the software version and model type of the central monitor compatible to data transfer function, refer to your nearest service representative.

The data cannot be transferred from the DS-8007 to DS-8400.

<u>Cause</u>

CFast Card is not inserted to the DS-8400.

Solution

Use the CFast Card to save the full disclosure waveform data.

Alarm settings, parameter settings are not transferred from the transport monitor. The confirmation window to apply the alarm settings, parameter settings of the transport monitor is not displayed.

<u>Cause</u>

The "Data for Transfer" setting is set to [OFF].

Solution

Check if the "Data for Transfer" setting is set to [ON]. If not, set it to [ON].

IB-8004 Input Box

The system does not start although the power of the DS-8400 is turned ON. <IB-8000-* Check Conn.> is displayed.

Cause 1

The power cable of the main unit is not connected.

Solution

Turn OFF the power and connect the power cable.

Cause 2

The module connection cable is not properly connected.

Solution

Turn OFF the power, and make sure that the module connection cable is securely connected.

□<IB*Slot* Module Disconnected> is displayed.

Cause 1

Infrared communication port is unclean.

Solution

Remove the expansion module, clean the infrared communication port, and insert the expansion module again.

Cause 2

The slot used in the IB-8004 is malfunctioning.

Solution

Insert the expansion module into another slot. If it operates properly, the IB-8004 needs to be repaired. Cease the measurement, and contact your nearest service representative.

Cause 3

IB-8004 or expansion module is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

The LAN-ID indicator keeps flashing.

Cause 1

The LAN ID setting is not correct.

Solution

Check the LAN-ID setting ID and make sure to set the correct LAN-ID.

Cause 2

The IB-8004 is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

Power is not supplied to the expansion module.

Cause 1

The expansion module is not properly connected.

Solution

Remove and reinsert the expansion module.

Cause 2

The slot used in the IB-8004 is malfunctioning.

Solution

Insert the expansion module into another slot. If it operates properly, the IB-8004 needs to be repaired. Cease the measurement, and contact your nearest service representative.

Cause 3

IB-8004 or expansion module is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□<IB-8000-* Check Unit>, <IB-8000-* Out of Operating Temp. Range> is displayed.

<u>Cause</u>

The IB-8004 is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□<IB-8000-* Failure> is displayed.

<u>Cause</u>

This message may be occasionally displayed when the expansion module is removed/inserted.

Solution

If the message automatically disappears, there is no problem on the equipment. If the message is repeatedly displayed, immediately turn OFF the power and cease the operation as the failure of the IB-8004 can be considered.

Expansion Module

The system does not start although the power of the DS-8400 is turned ON.

Cause 1

The power cable of the main unit is not connected.

Solution

Turn OFF the power, and make sure that the module connection cable is securely connected.

Cause 2

The module connection cable is not properly connected.

Solution

Turn OFF the power, and make sure that the module connection cable is securely connected.

Cause 3

The standby switch of the display unit is set to OFF.

Solution

Turn ON the standby switch on the display unit.

Check Conn.> is displayed.

Cause 1

Infrared communication port is unclean.

Solution

Remove the expansion module, clean the infrared communication port, and insert the expansion module again.

Cause 2

The slot used in the IB-8004 is malfunctioning.

Solution

Insert the expansion module into another slot. If it operates properly, the IB-8004 needs to be repaired. Cease the measurement, and contact your nearest service representative.

IB-8004 or expansion module is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

Power is not supplied to the expansion module.

Cause 1

The expansion module is not properly connected.

Solution

Disconnect and connect the expansion module.

Cause 2

The slot used in the IB-8004 is malfunctioning.

Solution

Insert the expansion module into another slot. If it operates properly, the IB-8004 needs to be repaired. Cease the measurement, and contact your nearest service representative.

Cause 3

IB-8004 or expansion module is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□ The following error messages related to the expansion module are displayed on the monitor. <IB* Slot* Module Failure>, <IB* Slot* Analog Unadjusted>,

<IB* Slot* Check Module>, <IB* Slot* Out of Operating Temp. Range>

<u>Cause</u>

The module connected to the IB-8004 slot is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

□<IB* Slot* Module Disconnected> is displayed.

Cause 1

The expansion module is not properly connected.

Solution

Disconnect and connect the expansion module.

Cause 2

The slot used in the IB-8004 is malfunctioning.

Solution

Insert the expansion module into another slot. If it operates properly, the IB-8004 needs to be repaired. Cease the measurement, and contact your nearest service representative.

Cause 3

IB-8004 or expansion module is malfunctioning.

Solution

Immediately turn OFF the power and cease the operation. Contact your nearest service representative.

Extended Display Unit

Nothing is displayed on the extended display unit. The same display with the main unit is displayed.

Cause

The video cable of the extended display unit is connected to the external monitor connector of the main unit.

Solution

Connect the cable to the extended display unit connector on the main unit.

The touch panel does not function on the extended display unit.

<u>Cause</u>

The serial communication cable is not connected.

Solution

Connect the serial communication cable of the extended display unit to the extended serial connector (COM A, COM B) of the main unit.

Ventilator

 \Box <Vent. Alarm> is displayed.

<u>Cause</u>

The following alarm has generated on the ventilator.

- Parameter alarm such as AWP, MV, FiO₂
- Technical alarm such as battery replacement of the ventilator

Solution

Check the alarm cause of the ventilator, and take appropriate action.

Vent. Offline> is displayed.

<VENT COMM> is displayed on the monitor and the ventilator.

Cause 1

The cable between the DS-8400 System and the ventilator is disconnected or not securely connected. Solution

Make sure the cable is properly connected.

Cause 2

The power of the ventilator is turned OFF.

Solution

Turn ON the power of the ventilator.

Cause 3

The ventilator is in standby mode. Solution Start the ventilation on the ventilator.

The network setting of the monitor does not match with the ventilator.

Solution

Make sure that the network setting of the connecting equipments are as follows.

SV-300/SERVO-i/SERVO-s

• No network setting.

SERVO-U/n/air

• No network setting.

VELIA, ASTRAL, VS ULTRA

• No network setting.

PB-740/760/840

- Baud Rate: 9600 bps
- Parity Bit: None
- Stop Bit: 1
- Data Bit: 8

Evita4/2dura/XL

- Communication Protocol: Medibus
- Baud Rate: 19200 bps
- Parity Bit: Even
- Stop Bit: 1

Multigas Unit

GAS Unit Failure> is displayed.

<u>Cause</u>

A hardware failure was detected on the gas unit.

Solution

Contact your nearest service representative.

GAS Check Sample Line> is displayed.

<u>Cause</u>

The sampling line or water trap is completely occluded.

The moisture inside the sampling line is drawn towards the water trap to be removed.

Solution 1

Check if the sampling line is occluded. Remove the occlusion if found.

Solution 2

Replace the sampling line, water trap.

GAS Check Water Trap> is displayed.

<u>Cause 1</u>

The water trap of the gas unit is not inserted, or not properly attached.

Solution

Insert the water trap. Make sure the water trap is properly connected.

Cause 2

Water trap is partly clogged or damaged.

Solution

Replace the water trap.

GAS Check Water Trap Class> is displayed.

<u>Cause</u>

The patient classification is not corresponded to the used water trap and the sampling tube.

Solution

Make sure the patient classification is corresponded to the used water trap and the sampling tube. When the patient classification is "Adult" or "Child", make sure to use the water trap and sampling line intended for adult/pediatric.

When the patient classification is neonate, make sure to use the water trap and sampling line intended for neonate.

GAS Mixed Agents Detection> is displayed.

Cause

More than one halogenated anesthetic gas exists.

Solution 1

Make sure that multiple anesthetic gases are not used. Make sure that the anesthetic gas carburetor setting is correct.

Solution 2

If the problem persists, contact our service representative.

GAS Zeroing Failed> is displayed.

<u>Cause</u>

The zero calibration process has not been properly completed.

Solution

Perform the manual zero calibration again.

□<SPIRO Unit Failure> is displayed.

<u>Cause</u>

The hardware failure of the SPIRO unit was detected.

Solution

Contact your nearest service representative.

SPIRO Check FlowSensor Class> is displayed.

Cause 1

The flow sensor is disconnected or not securely connected.

Solution

Make sure that the flow sensor is securely connected.

Cause 2

The flow sensor is damaged.

Solution

Replace the flow sensor.

Cause 3

The used flow sensor does not correspond to the patient classification setting on the monitor.

Solution

Make sure that the used flow sensor corresponds to the patient classification setting. When the patient classification is adult, use the flow sensor intended for adult. When the patient classification is "Child" or "Neonate", use the flow sensor intended for pediatric.

□<SPIRO Check Flow Sensor> is displayed.

<u>Cause</u>

If the flow sensor is disconnected during multigas monitoring, the message will be displayed.

Solution 1

To cease multigas monitoring, press the [Alarm Silence] key. The message will disappear, and the alarm will be silenced.

Solution 2

To continue monitoring, plug in the flow sensor. The message will disappear, and the alarm will be silenced.

SvO₂/CCO Monitor

The numeric data is not displayed.

Cause 1

The cable is not properly connected.

Solution

Connect the following cable securely.

SvO ₂ /CCO Monitor	Connection Cable		
	For Status II Connector	For Serial Connector	
Vigilance	CJ-406RI-70Vigi (x1)	CJO-04RS4	
Vigilance CEDV	CJ-406RI-70Vigi (x1)	CJO-04RS4	
Vigilance II	CJ-402RI-70SVi (x1)	CJ-502	
Vigileo	CJ-402RI-70SVi (x1)	CJ-502	
EV1000	CJ-406RI-70Vigi (x1)	CJO-04RS4	
PiCCO2	CJO-19RS5 (x1)	CJO-18RS5	
PulsioFlex	-	CJ-725 ^{*1}	

*1: To connect the PulsioFlex PC4000, USB to RS-232C connector cable (UC-232A, ATEN) specified by PULSION Medical Systems is required.

Cause 2

The "External Device" setting is not correct.

Solution

Select [Vigilance/Vigileo], [PiCCO] or [PulsioFlex] for the port function on the "External Device" setup screen.

Cause 3

The measurement data is not displayed on the corresponding external device.

Solution

The measurement data of SvO_2 , CO, etc. will not be displayed on the monitor unless the data is displayed on the used external device. Check if the data is displayed on the used external device.

Cause 4

The CCO is not measured.

Solution

The CCO/CCI data will be displayed on this equipment only when CCO is measured on each external device.

Cause 5

The network setting of the monitor does not match with each external device.

Solution

The network setting of the monitor is fixed to the default setting of each external device and cannot be changed. Make sure that the network setting of the connecting equipment is in default setting.

In Case of Vigilance/Vigileo:

Make sure that the network is set as follows.

For procedure to check the Vigilance/Vigileo network setting, refer to the operation manual for the Vigilance/Vigileo.

- Device: IFM Out
- Baud Rate: 19200 bps
- Parity Bit: None
- Stop Bit: 1
- Data Bit: 8
- Flow Control: 2 sec.

In Case of PiCCO

Make sure that the network is set as follows.

For procedure to check the PiCCO network setting, refer to the operation manual for the PiCCO2.

• RS-232C Protocol: PiCCO2 V3.0

In Case of PulsioFlex:

Check if the network is set as follows.

For procedure to check the PulsioFlex network setting, refer to the operation manual for the PulsioFlex.

• RS232C protocol: PulsioFlex V1.0

Cause 6

The software version of Vigilance does not correspond.

Solution

If the Vigilance without the STAT function is connected, the STAT data will not be displayed. Check the software

version of the Vigilance.

Cause 7

The software version of PiCCO does not correspond.

Solution

The compatible version of PiCCO2 is from V3.0. Check the version of the PiCCO2.

Cause 8

The software version of PulsioFlex does not correspond.

Solution

Check the software version of PulsioFlex. The compatible version is PulsioFlex V1.0.

BIS Monitor (A-2000/A-3000)

The numeric data is not displayed.

Cause 1

If the SQI value is lower than 15, BIS data and SR data will not be displayed.

Solution

Refer to the BIS monitor operation manual and set the SQI value above 15.

Cause 2

The communication setting of the BIS monitor is incorrect.

Solution

ASCII should be set to communicate with this system.

Make sure that ASCII is set on the BIS monitor communication setting.

Refer to the BIS monitor operation manual for procedures.

Check BIS Conn.> is displayed.

<u>Cause</u>

The cable is not properly connected. The connection cable is disconnected.

Solution

Securely connect the connection cable to the serial or status connector of the main unit or the Status II connector of the HP-800.

BIS (When HBX-800 is used)

□<BISx Disconnected> is displayed.

<u>Cause 1</u> The BISx is disconnected. Solution Verify all cable connections and connect the BISx correctly.

The BISx cable is defective.

Solution

Check the cable including the connector part, and replace the cable if necessary.

Cause 3

The BISx is defective. Solution Replace the BISx.

□<BIS High Impedance, Check Sensor> is displayed.

Cause 1

The sensor is not fully in contact with patient's skin.

Solution

Attach the electrode firmly to patient's skin.

Cause 2

The sensor application is incorrect.

Solution

Read instructions on sensor package to ensure correct sensor placement.

Cause 3

The patient interface cable (PIC cable) is defective.

Solution

Replace the patient interface cable (PIC cable).

<u>Cause 4</u> The BISx is defective. Solution Replace the BISx.

□<BIS Sensor Disconnected> is displayed.

Cause 1

The sensor is disconnected. Solution Connect the sensor.

Cause 2

Poor or contaminated connection between the sensor and patient interface cable (PIC cable).

Solution

Clean the connection part, and connect them properly.

Cause 3

The patient interface cable (PIC cable) is disconnected.

Solution

Connect the patient interface cable (PIC cable) correctly.

The patient interface cable (PIC cable) is defective.

Solution

Replace the patient interface cable (PIC cable).

<u>Cause 5</u> The BISx is defective. Solution Replace the BISx.

□<BIS Perform "Sensor Check"> is displayed.

<u>Cause 1</u>

At least one element of sensor has too high impedance, and "Sensor Check" window is closed before sensor check completes.

Solution

Press the "Sensor Check" key to start the sensor check process and ensure that <PASS> is displayed.

Cause 2

The sensor application is incorrect.

Solution

Read instructions on sensor package to ensure correct sensor placement.

Cause 3

The sensor is not properly connected.

Solution

Verify that the sensor is properly connected.

Cause 4

The patient interface cable (PIC cable) is defective.

Solution

Replace the patient interface cable (PIC cable).

Cause 5

The BISx is defective. Solution Replace the BISx.

Artifacts> is displayed.

Situation: The signal quality is less than half of the level desirable for optimal monitoring conditions.

(NOTE

 This message may occur as the results of artifact (non-EEG signal) such as those generated from motion (patient movement or eye blinks) or the presence of electrocautery, warming blankets, or other devices.

Cause 1

Artifact, such as those generated by motion or eyeblinks, is causing loss of EEG recognition.

Solution

If <Artifacts> appears on the display, attempt to identify and eliminate artifact source.

EMG bar indicates electrical activity that may be interfering with EEG recognition.

Solution

If EMG bar is illuminated, attempt to determine and eliminate cause.

Cause 3

The patient interface cable (PIC cable) is defective.

Solution

Make sure that sensor check passes. If not, replace the patient interface cable (PIC cable).

Cause 4

The BISx is defective. Solution Replace the BISx.

\Box < BIS SQI < 15% > is displayed.

Situation: The signal quality is too low to accurately calculate a BIS value.

The BIS value and other trend variables that are adversely affected by artifact are not displayed.

NOTE

 This message may occur as the results of artifact (non-EEG signal) such as those generated from motion (patient movement or eye blinks) or the presence of electrocautery, warming blankets, or other devices.

Cause 1

Artifact, such as those generated by motion or eyeblinks, is causing loss of EEG recognition. Solution

If <BIS SQI < 15%> appears on the display, attempt to identify and eliminate artifact source.

Cause 2

EMG bar indicates electrical activity that may be interfering with EEG recognition.

Solution

If EMG bar is illuminated, attempt to determine and eliminate cause.

Cause 3

The patient interface cable (PIC cable) is defective.

Solution

Make sure that sensor check passes. If not, replace the patient interface cable (PIC cable).

Cause 4

The BISx is defective. Solution Replace the BISx.

□<BISx Incompatible> is displayed.

<u>Cause</u>

The sensor is not compatible with the monitor configuration.

Solution

Replace the BISx.

Check BIS Sensor, Perform Sensor Check> is displayed.

<u>Cause</u>

Problem is detected relating to sensor ground element, or sensor is using too much current. Solution 1 Disconnect and examine sensor connection, clean any contamination, then perform "Sensor Check". Solution 2 Replace the sensor if necessary, then perform "Sensor Check". Solution 3 Replace the patient interface cable (PIC cable), then perform "Sensor Check". Solution 4 Replace the BISx, then perform "Sensor Check".

□<Replace BIS Sensor, Too Many Uses>, <Replace BIS Sensor, Invalid Sensor> is displayed.

Cause 1

Sensor has been connected and disconnected too many times.

Solution

Replace the sensor.

<u>Cause 2</u> The sensor is invalid. Solution Replace the sensor.

□<Sensor Usage > 24hrs.> is displayed.

<u>Cause</u>

The sensor was attached to the system for more than 24 hours. Solution Replace the sensor.

□<BISx Failure> is displayed.

<u>Cause</u> The BISx is defective. Solution Replace the BISx, then perform "Sensor Check".

The power indicator on the HBX-800 is lit in red.

<u>Cause</u>

The HBX-800 is defective.

Solution

Cease using the equipment and contact your nearest service representative to repair the equipment.

INVOS

The numeric data is not displayed.<Check INVOS Conn.> is displayed.

<u>Cause</u>

The cable is not properly connected. The connection cable is disconnected.

Solution

Securely connect the connection cable to the serial or status connector of the main unit or the Status II connector of the HP-800.

FLOW-i

The numeric data is not displayed.<Check FLOW-i Conn.> is displayed.

Cause 1

The cable is not properly connected. The connection cable is disconnected.

Solution

Securely connect the connection cable to the serial or status connector of the main unit or the Status II connector of the HP-800.

Cause 2

The FLOW-i is in standby mode.

Solution

The numeric data will be displayed when the measurement is started on the FLOW-i.

Cause 2

The software version of the FLOW-i is not compatible with the DS-8400.

Solution

The compatible software version of DS-8400 is from 01-01. The compatible software version of FLOW-i is system software version 02 and 03 (FCI Protocol version 0004 and 0005 respectively).

PC Communication

Check System Conn.> is displayed.

Cause 1

The cable is disconnected or not properly connected. The power is not supplied to the communication port.

Solution

Connect the cable securely. Check if the power is supplied to the communication port by checking the communication indicator.

Cause 2

Communication with the PC is not performed. The communication is ceased.

Solution

Resume the communication with the PC. The communication time out period is about 1 minute.

TCM4/TCM5 FLEX

The numeric data is not displayed.

<u>Cause 1</u>

The cable is not properly connected.

Solution

Connect the following cable securely.

Transcutaneous Blood Gas	Connection Cable		
Monitor	For Status II Connector For Serial Connecto		
TCM4	-	CJ-726 (straight) ^{*1}	
TCM5 FLEX	-	CJ-725 (cross) ^{*2}	

*1: To connect the TCM4, the cable specified by Radiometer Medical ApS is required. The communication will be enabled by connecting the CJ-726 and the cable specified by Radiometer Medical ApS.

*2: To connect the TCM5 FLEX and CJ-725, D-sub 9-pin male to male gender changer (inch screw) is required.

<u>Cause 2</u>

The "External Device" setting is not correct.

Solution

Select [TCM4/TCM5] for the port function on the "External Device" setup screen.

Cause 3

The measurement data is not displayed on the corresponding external device.

Solution

The measurement data of $tcpO_2$, $tcpCO_2$ will not be displayed on the monitor unless the data is displayed on the used external device. Check if the data is displayed on the used external device.

Cause 4

The network setting of the monitor does not match with each external device.

Solution

The network setting of the monitor cannot be changed. Make sure that the network setting of the connecting device is as follows. For details of the network setting on the TCM4 and TCM5 FLEX, refer to Radiometer Medical ApS.

In case of TCM4:

RS-232C Protocol: Monlink

In case of TCM5 FLEX:

• RS-232C Protocol: Monlink2.0

<u>Cause 5</u>

The software version of TCM4 or TCM5 FLEX does not correspond.

Solution

Check the software version of TCM4 or TCM5 FLEX. For details of the network setting on the TCM4 and TCM5 FLEX, refer to Radiometer Medical ApS.

TCM4: Version 3.04 TCM5 FLEX: Version 1.18

Check TCM Conn.> is displayed.

<u>Cause 1</u>

The cable is disconnected, or not securely connected.

Solution

Connect the cable correctly.

Cause 2

The power of the external device has been turned OFF.

Solution

Turn ON the power of the external device.

Cause 3

The TCM series device other than TCM4, TCM5 FLEX is connected.

Solution

Only TCM4, TCM5 FLEX can be connected. Check the model type of the external device.

Magnetic Card Reader/Barcode Reader

The magnetic card reader or barcode reader does not function.

<u>Cause</u>

The conversion cable (CJ-756) is not connected.

Solution

If the magnetic card reader or barcode reader is connected directly to the serial port on this equipment without the conversion cable, it will not function. Make sure to use the conversion cable.

External Media

CF/CFast Card Slot: There is no card in the slot.> is displayed.

<u>Cause</u>

CF/CFast card is not inserted or not correctly set in the CF/CFast card slot.

Solution

Set the CF/CFast card into the CF/CFast card slot.

Control Con

Cause 1

There is no data on the CF card.

Solution

Check if the CF card is readable. Or, check if the data is present on the CF card. Pressing "Yes" will not start reading the compatible data.<Card access error.> will be displayed.

Error is detected during the read process.

Solution

The data may not be correctly written on the CF card. Format the card again on the used equipment and try the write/read process again. Pressing "Yes" will not start reading the compatible data.

□<Card access error.> is displayed.

Cause 1

There is not enough capacity on the CF card to write the data.

Solution

Check the remaining card capacity.

Format the card again on the used equipment and try the write/read process again.

Cause 2

Error is detected during the write process.

Solution

Make sure that the CF card is properly inserted and try the write process again.

Format the card again on the used equipment and try the write/read process again.

Cause 3

Unspecified CF card is used. Solution Use the specified CF card.

□No data on the CF card.

<u>Cause</u>

There is no data on the CF card.

Solution

Check if the CF card is readable. Or, check if the data is present on the CF card.

Wrong CF card for full disclosure.>, <Failed to read full disclosure from the CF card.>

<u>Cause</u>

Specified memory card is not used. The card is not formatted. The data stored in the card is damaged. The card has been already used on another equipment.

Solution 1

Use the recommended memory card.

Disconnect and connect the full disclosure waveform card again to make sure that it is properly inserted. Format the card on the used equipment. (All previous data will be deleted.)

Solution 2

If the error persists, contact your nearest service representative.

12-1
12-2
12-6
12-10
12-15

Contents

Chapter 12 Setup Item/Default Value

This section lists selection, default setting, and backup status for each setup item. The following indicates the selection, default setting and backup status for each setup item.

Patient Admit / Discharge

Item	Description	Default	At Power ON	At Discharge
Mode Selection	Main Mode 1 to 9, Sub Mode 1 to 6, Extended Display1 1 to 3, Extended Display2 1 to 3	Main Mode 1	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/ Discharge].	
ID	Numeric, Alphabet, Symbol (20 characters)	Blank	Backup	Initialize
Patient Name	Numeric, Alphabet, Symbol (16 characters)	Blank	Backup	Initialize
Patient Classification	Adult, Child, Neonate	Adult	Depends on the "Patient Classification" setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Sex	Male, Female	No selection	Backup	Initialize
Team	Red, Orange, Yellow, Yellow-green, Green, Light Blue, Blue, Purple	Red	Backup	Initialize
Birth Date	Birth Date	Blank	Backup	Initialize
Age	0 year to 150 years or 0 day to 999 days	0 year	Backup	Initialize
Height	0.0 cm to 300.0 cm	0.0 cm	Backup	Initialize
Weight	0.0 kg to 350.0 kg	0.0 kg	Backup	Initialize
BSA	0.00 m ² to 9.99 m ²	0.00 m ²	Backup	Initialize
Blood Type	A, B, O, AB Rh +/-	Blank	Backup	Initialize
Pacemaker	Used, Not used	Not Used	Depends on the setting under [Initial Settings>User I/F>Power ON/Discharge].	
Impedance Measurement	ON, OFF	ON	Depends on the setting under [Initial Settings>User I/F>Power ON/Discharge].	
Admit Date	Year, Month, Day	Blank	Backup	Initialize

Item	Description	Default	At Power ON	At Discharge
System Alarm	Suspend, ON	Suspend	-	-
HR*2	ON, OFF 20 bpm to 300 bpm 5 bpm increments	ON 40 bpm to 120 bpm	Depends on the "Main Mod setting under [Initial Settings>U I/F>Power ON/Discharge]. If "Main Mode" setting is [Back Depends on the "Alarm" setti under [Power ON/Discharge	
PR_IBP ^{*2} PR_SpO ₂ ^{*2}	ON, OFF 20 bpm to 300 bpm 5 bpm increments	OFF OFF to OFF		
Asystole	ON 3 sec. to 10 sec. 1 sec. increments	ON 5 sec.		
VF	ON	ON		
VT	ON	ON		
Slow_VT	ON, OFF	ON	1	
Run	ON, OFF 2 beats to 8 beats 1 beat increments	ON 3 beats		
Couplet	ON, OFF	OFF		
Pause	ON, OFF 1.5 sec. to 5 sec. 0.5 sec. increments	OFF 3.0 sec.		
BIGEMINY	ON, OFF	OFF	1	
TRIGEMINY	ON, OFF	OFF	1	
FREQUENT	ON, OFF 1 bpm to 50 bpm 1 bpm increments	OFF, 10 bpm	_	
Tachy	ON, OFF	ON	1	
Brady	ON, OFF	ON		
Ext Tachy ^{*2}	ON, OFF 22 bpm to 300 bpm 5 bpm increments	OFF 150 bpm		
Ext Brady ^{*2}	ON, OFF 20 bpm to 295 bpm 5 bpm increments	OFF 30 bpm		
Triplet	ON, OFF	OFF	1	
R on T	ON, OFF 200 ms to 600 ms 8 ms increments	OFF 320ms		
Multiform	ON, OFF	OFF	1	
Vent Rhythm	ON, OFF	OFF	1	
SVT	ON, OFF 2 beats to 10 beats 1 beat increments	OFF 6 beats		
Irregular RR	ON, OFF 10% to 20% 5% increments	OFF 10%		
Prolonged RR	ON, OFF	OFF	1	
S FREQUENT	ON, OFF 1 bpm to 50 bpm 1 bpm increments	OFF 10 bpm		
S Couplet	ON, OFF	OFF	1	
VPC	ON, OFF	OFF	1	
SVPC	ON, OFF	OFF	1	
Pacer not Capture	ON, OFF 80 ms to 480 ms 8 ms increments	OFF 320 ms	1	

Item	Description	Default	At Power ON	At Discharge
Pacer not Pacing	ON, OFF 20 bpm to 200 bpm 5 bpm increments	OFF 50 bpm		
HR Lower Limit for VT	120 bpm, 140 bpm	120	Depends on th setting under [Init	e "Main Mode"
HR Lower Limit for Run ^{*3}	0 bpm to 100 bpm 10 bpm increments	40 bpm	I/F>Power O If "Main Mode" se	N/Discharge]. etting is [Backup]; e "Alarm" setting
HR Lower Limit for SVT	100 bpm to 250 bpm 10 bpm increments	150 bpm	under [Power (ON/Discharge].
ST1 to ST12(mm) ^{*1}	ST All Alarm ON, OFF Individual Alarm ON, OFF ±20 mm 1 mm increments	ST All Alarm OFF Individual Alarm OFF OFF to OFF		
ST1 to ST12(mV) ^{*1}	ST All Alarm ON, OFF Individual Alarm ON, OFF ±2.00mV0.1mV increments	ST All Alarm OFF Individual Alarm OFF OFF to OFF		
BP1 (mmHg) ^{*4}	ON, OFF 0 mmHg to 300 mmHg 5 mmHg increments	ON SYS: 80 to 180 DIA: OFF to OFF MEAN: OFF to OFF		
BP1 (kPa) ^{*4}	ON, OFF 0 kPa to 40.0 kPa 0.5 kPa increments	ON SYS: 10.0 to 24.0 DIA: OFF to OFF MEAN: OFF to OFF		
BP2 to BP8 (mmHg) ^{*4}	ON, OFF 0 mmHg to 300 mmHg 5 mmHg increments	OFF SYS: OFF to OFF DIA: OFF to OFF MEAN: OFF to OFF		
BP2 to BP8 (kPa) ^{*4}	ON, OFF 0 kPa to 40.0 kPa 0.5 kPa increments	OFF SYS: OFF to OFF DIA: OFF to OFF MEAN: OFF to OFF		

^{*1}: The same setting applies for "mm" and "mV".

^{*2}: For HR, Ext Tachy, Ext Brady, 60 bpm or lower can be set in 1 bpm increments. For PR_SpO₂, 25 bpm or lower can be set in 1 bpm increments.

 $^{\star3}\!\!:$ "HR Lower Limit for Run" can be set in 5 bpm increments for 50 bpm and above.

^{*4}: For BP, 50 mmHg/7.0 kPa or lower can be set in 2 mmHg/0.2 kPa increments.

Item	Description	Default	At Power ON At Discharge
CVP (mmHg) (kPa) ^{*2}	ON, OFF 0 mmHg to 300 mmHg 5 mmHg increments 0 kPa to 40 kPa 0.5 kPa increments	OFF SYS: OFF to OFF DIA: OFF to OFF MEAN: OFF to OFF	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/Discharge]. If "Main Mode" setting is [Backup]; Depends on the "Alarm" setting under [Power ON/Discharge].
CVP (cmH ₂ O)	ON, OFF 0 cmH ₂ O to 40 cmH ₂ O 1 cmH ₂ O increments	OFF SYS: OFF to OFF DIA: OFF to OFF MEAN: OFF to OFF	
RR_IMP ^{*3} RR_VENT ^{*3} RR_CO ₂ ^{*3}	ON, OFF 5 Bpm to 150 Bpm 5 Bpm increments 5 Bpm to 150 Bpm (Neonate) 2 Bpm increments	ON 5 Bpm to 30 Bpm	
Apnea	ON, OFF 10 sec. to 60 sec. 1 sec. increments	ON 15 sec.	
SpO ₂	ON, OFF 50% SpO ₂ to 100%SpO ₂ 1%SpO ₂ increments	ON 90%SpO ₂ to OFF	
EXT SpO ₂	ON, OFF 50-90%SpO ₂ 1%SpO ₂ increments	ON 80%SpO ₂	
SpCO	ON, OFF 1%SpCO to 40%SpCO 1%SpCO increments	OFF	
SpMet	ON, OFF 1%SpMet to 15%SpMet 1%SpMet increments	OFF	
SpHb	ON, OFF 1.0 g/dL to 24.5 g/dL 0.1 g/dL increments	OFF	
NIBP (mmHg)	ON, OFF 10 mmHg to 300 mmHg 5 mmHg increments	ON SYS: 80 to 180 DIA: OFF to OFF MEAN: OFF to OFF	
NIBP (kPa)	ON, OFF 1.5 kPa to 40.0 kPa 0.5 kPa increments	ON SYS: 10.0 to 24.0 DIA: OFF to OFF MEAN: OFF to OFF	
TEMP1 to TEMP8 (°C)	ON, OFF 30-45°C 0.5°C increments	OFF, OFF to OFF	
Tb (°C)	ON, OFF 30-45°C 0.5°C increments	OFF, OFF to OFF	
CO ₂ -E (mmHg) ^{*1}	ON, OFF 1 mmHg to 100 mmHg 1 mmHg increments	OFF	Depends on the "Main Mode" setting under [Setup>Initial Settings>User I/F>Power ON/
CO ₂ -E (kPa) ^{*1}	ON, OFF 0.1 kPa to 13.3 kPa 0.1 kPa increments	OFF	Discharge].
CO ₂ -E (%) ^{*1}	ON, OFF 0.1% to 13.3% 0.1% increments	OFF	
CO ₂ -I (mmHg) ^{*1}	ON, OFF 1 mmHg to 4 mmHg 1 mmHg increments	OFF	
CO ₂ -I (kPa) ^{*1}	ON, OFF 0.1 kPa to 0.4 kPa 0.1 kPa increments	OFF	
CO ₂ -I (%) ^{*1}	ON, OFF 0.1% to 0.4% 0.1% increments	OFF	
O ₂ -E (%) ^{*1}	ON, OFF 18-100%	OFF	1
O ₂ -I (%) ^{*1}	ON, OFF 18-100%	OFF	
N ₂ O-E (%) ^{*1}	ON, OFF 0-100%	OFF	
N ₂ O-I (%) ^{*1}	ON, OFF 0-100%	OFF	1

ľ	tem	Description	Default	At Power ON	At Discharge
ISO-E (%), HAL-E (%), ENF-E (%) ^{*1}		ON, OFF 0.5-6.0%	OFF	Depends on th setting under Initi	
ISO-I (%) , HAL-I (%) , ENF-I (%) ^{*1}		ON, OFF 0.5-6.0%	OFF	I/F>Power OI If "Main Mode" se	N/Discharge].
SEV-E (%) *1		ON, OFF 0.5-8.0%	OFF	Depends on the under [Power ("Alarm" setting
SEV-I (%) *1		ON, OFF 0.5-8.0%	OFF		JN/Dischargej.
DES-E (%) *1		ON, OFF 0.5-18.0%	OFF		
DES-I (%) *1		ON, OFF 0.5-18.0%	OFF		
MAC ^{*1}		ON, OFF 0.1 to 9.9	OFF		
PEAK ^{*1}		ON, OFF 8-100cmH ₂ O	OFF		
PEEP ^{*1}		ON, OFF 2 cmH ₂ O to 50 cmH ₂ O	OFF		
MV-E ^{*1}		Adult: ON, OFF 0.5 L/min to 20 L/min Child, Neonate: ON, OFF 0.5 L/min to 5 L/min	OFF		
BIS (When HBX-	800 is used)	ON, OFF 1 to 99 increments of 1	ON 40 to OFF		
Alarm Settings (Setup)	Alarm Suspend Time	1 min., 2 min.	2 min.	Depends on th setting under [Init I/F>Power OI	ial Settings>User
	Alarm Silence Time	1 min., 2 min.	2 min.		
	Alarm Sound Suspend	ON, OFF	ON		
	Alarm Sound Suspend Time	[1min.]/[2min.]/[5min.]/[10min.]/ [30min.]/[60min.]/[90min.]/ [20min.]/[240min.]/[360min.]	60 min.		
	Status Alarm Control Status Alarm Control	Link to alarm silence time, Link to each new occurrence	Link to each new occurrence		
	Alarm Limit Display	Graph, Numeric, OFF	Graph		

*1: When the numeric data acquired from FLOW-i is displayed, alarms cannot be set. Also, these alarms will not generate.

*2: For CVP, 50 mmHg/7.0 kPa or lower can be set in 2 mmHg/0.2 kPa increments.

*3: For RR, 1 Bpm increments may be applied depending on the "RR Alarm Increment" settings. (@Maintenance Manual "User I/F" P5-13)

NOTE

 By selecting [Backup] for "Power ON" and "Discharge" under [Initial Settings>User I/F >Power ON/Discharge], the settings will be retained at "Power ON" and "Discharge" respectively. If [Initialize] is selected, the settings will be initialized at "Power ON" and "Discharge".

Parameter

ECG

Item	Description	Default	At Power ON	At Discharge
Leads	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6	ECG1: II ECG2: aVR ECG3: I ECG4: III ECG5: aVL ECG6: avF ECG7: V1 ECG8: V2 ECG9: V3 ECG10: V4 ECG11: V5 ECG12: V6		e setting under ngs>User I/ I/Discharge].
Size	Auto, x1/4, x1/2, x1, x2, x4	ECG1 to ECG12 x1	Backup	Initialize
Filter Mode	Monitor, Diagnosis, ESIS	Monitor	Backup	Backup
Synchronized Mark/Tone	ECG, SpO ₂ -1, SpO ₂ -2, BP, Auto, OFF	Auto	Backup	Backup
Pacemaker	*Same with "Patient Admit/Discharge" section.			
Pacemaker Pulse	ON, OFF	OFF	Backup	Backup
Pace Pulse Mask Time	Auto, 10ms, 20ms, 40ms, OFF	Auto	Backup	Initialize
HR Average	Average, Instant	Average	Backup	Backup
HR Delay	ON, OFF	OFF	Backup	Backup
Drift Filter	ON, OFF	OFF	Backup	Backup
AC Filter	ON, OFF	ON	Backup	Backup
Auto Lead	ON, OFF	OFF	Backup	Backup
3-lead Override	ON, OFF	OFF	Backup	Backup
ST/VPC/Arrhy. Alarm Display	ON, OFF	ON	Backup	Backup
ECG Analog Output	Disp. Lead, Selected Lead	Disp. Lead	Backup	Backup
ECG Waveform Display during Lead-OFF	ON, OFF	OFF	Backup	Backup
Noise Detection	ON, OFF	OFF	Backup	Backup
Chest Lead-OFF	Enable, Disable	Enable	Backup	Backup

RESP

Item	Description	Default	At Power ON	At Discharge
Size	x1/4, x1/2, x1, x2, x4	x1	Backup	Initialize
RR Synchronized Mark	ON, OFF	ON	Backup	Backup
RR/APNEA Alarm Source	Auto, Impedance, Ventilator, CO ₂ /GAS	Auto	Backup	Backup
CVA Detect	ON, OFF	OFF	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Impedance Measurement	*Same with "Patient Admit/Discharge" section.			
Impedance Detection Lead	1, 11	11	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	

Item	Description	Default	At Power ON	At Discharge
Size	x1/4, x1/2, x1, x2, x4	x1	Backup	Initialize
Synchronized Mark/Tone	*Same with ECG setting.			
Alarm during NIBP	ON, OFF	ON	Backup	Backup
Label	None/Auto/RH/LH/RF/LF/OT	None	Backup	Backup

 SpO_2 (NellcorTM)

Item	Description	Default	At Power ON	At Discharge
Second Alarm	OFF, 10, 25, 50, 100	OFF	Backup	Backup

SpO₂ (Masimo Unit)

Item	Description	Default	At Power ON	At Discharge
SpO ₂ Averaging	2-4 sec, 4-6 sec, 8 sec, 10 sec, 12 sec, 14 sec, 16 sec	8 sec.	Backup	Backup
Pulse Sensitivity	High, Normal, APOD	Normal	Backup	Backup
FAST SAT	ON, OFF	OFF	Backup	Backup
Perfusion Index	ON, OFF	ON	Backup	Backup
Signal IQ Wave	ON, OFF	OFF	Backup	Backup
SpHb Averaging	Short, Medium, Long	Medium	Backup	Backup

NIBP

Item	Description	Default	At Power ON	At Discharge
Patient Classification	*Same with "Patient Admit/Discharge" section.	·		
Quick Measurement	ON, OFF	ON	Depends on the [Initial Settin F>Power ON	ngs>User I/
NIBP Auto Mode	Cont., 1min, 2min, 2.5min, 5min, 10min, 15 min, 20min, 30min, 60min, 120min, Lumbar Mode, OFF	OFF	Depends on the [Initial Settin F>Power ON	ngs>User I/
Dyna Alert	ON, OFF	ON	Backup	Backup
Sight Inflation	ON, OFF	OFF	Backup	Backup
Oscillograph	ON, OFF	OFF	Backup	Backup
MAP	ON, OFF	ON	Backup	Backup
PR Display	ON, OFF	OFF	Backup	Backup
End Tone	ON, OFF	ON	Backup	Backup
NIBP Erase Time	60 min., 120 min.	120 min.	Backup	Backup
User Interval	Lumbar Mode	Lumbar Mode	Backup	Backup
Measure at Alarm	ON, OFF	OFF	Backup	Backup
	Asystole, VF, VT, Slow VT, Run, Bigeminy, Trigeminy, Pause, Couplet, Tachy, Brady, Frequent	No Selection	Backup	Backup
	HR, ST, RR, APNEA, SpO ₂ , BP1, BP2, BP3, BP4, BP5, BP6, BP7, BP8, T1, T2, T3, T4, T5, T6, T7, T8, Tb, CO ₂ , O ₂ , N ₂ O, AGENT, SpCO, SpMet, SpHb, MV, PEEP, PEAK	No Selection	Backup	Backup
Auto Mode with Start/ Stop key	ON, OFF	ON	Backup	Backup
Time Display	Elapsed, Meas.	Elapsed	Backup	Backup

Item	Description	Default	At Power ON	At Discharge
Periodic Measurement Starting Time	Time, Meas.	Time	Backup	Backup
Target Inflation Value	Adult: 100 mmHg to 290 mmHg Child: 100 mmHg to 200 mmHg Neonate: 100 mmHg to 140 mmHg	Adult: 180 mmHg Child: 140 mmHg Neonate: 110 mmHg	Backup	Initialize

BP1 to 8

Item	Description	Default	At Power ON	At Discharge	
Scale [*]	20, 50, 75, 100, 150, 200, 250, 300 mmHg	200 mmHg 50 mmHg (BP2)		ne setting under ings>User I/ N/Discharge].	
	4, 8, 12, 16, 20, 24, 32, 40 kPa	24 kPa 8 kPa (BP2)	-		
Label	BP*, ART, PAP, CVP, ICP, IAP, LVP, US1 to US5	BP* (BP1 to BP8)	Backup	Backup	
Synchronized Mark/Tone	*Same with ECG setting.				
Display Type	S/M/D, S/D, M	S/M/D	Backup	Backup	
Wave Filter	6, 8, 12, 40 Hz	12Hz	Backup	Backup	
Mean Wave	ON, OFF	OFF	Backup	Backup	
Respiration Filter	ON, OFF	OFF	Backup	Backup	
Alarm during NIBP	ON, OFF	ON	Backup	Backup	

*: The scale selection will differ depending on the label.

TEMP1 to TEMP8

Item	Description	Default	At Power ON	At Discharge
Label	T#, Tsk, Tre, Tes, Tco, US1 to US7	T* (T1 to T8)	Backup	Backup

ΔTEMP-A to TEMP-D

Item	Description	Default	At Power ON	At Discharge
ΔTemp-A	(T1-T8) to (T1-T8)	T1 to T2	Backup	Backup
ΔTemp-B	(T1-T8) to (T1-T8)	T3 to T4	Backup	Backup
ΔTemp-C	(T1-T8) to (T1-T8)	T5 to T6	Backup	Backup
ΔTemp-D	(T1-T8) to (T1-T8)	Т7-Т8	Backup	Backup

CO₂ (Capnostat 5/HPD-800/HPD-810/HPD-820)

Item	Description	Default	At Power ON	At Discharge
Scale	0-50, 0-100 mmHg	0-50	Depends on the	e setting under
	0-4, 0-8, 0-10 kPa	0-4	[Initial Settin F>Power ON	ngs>User I/
	0-4, 0-8, 0-10%	0-4		/Discharge].
EtCO ₂ Peak Duration	10 sec, 20 sec, OFF	10 sec.	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
CO ₂ Source Priority	MGU-800, HS-8000	HS-8000	Backup	Backup
O ₂ Compensation	0-100%	21%	Backup	Backup
N ₂ O Compensation	ON, OFF	OFF	Backup	Backup
Anesthetic Compensation	0.0-20.0%	0.0%	Backup	Backup
Atmospheric Pressure	400-850 mmHg	760 mmHg	Backup	Backup

CO2 (COVIDIEN/HCP-800/HCP-810/HCP-820)

Item	Description	Default	At Power ON	At Discharge
Scale	0-50, 0-100 mmHg	0-50	Depends on the	e setting under
	0-4, 0-8, 0-10 kPa	0-4	[Initial Setti	ngs>User I/
	0-4, 0-8, 0-10%	0-4	F>Power ON	I/Discharge].
EtCO ₂ Peak Duration	10 sec, 20 sec, OFF	10 sec.	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
CO ₂ Source Priority	MGU-800, HS-8000	HS-8000	Backup	Backup

SPIRO, Ventilator, FLOW-i

Item	Description	Default	At Power ON At Discharge	
AWP Scale	10, 20, 30, 50, 120 cmH ₂ O	50 cmH ₂ O	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
AWF Scale	5, 10, 20, 50, 180 L/min	50 L/min	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
AWV Scale	50, 250, 500, 1000, 3000 mL	500 mL	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	

Cardiac Output (CO)

Item	Description	Default	At Power ON	At Discharge
Auto Start	ON, OFF	ON	Backup	Backup
Time Scale	30 sec., 60 sec.	30 sec.	Backup	Backup

Multigas Concentration, FLOW-i

Item	Description	Default	At Power ON	At Discharge
GAS_CO ₂ Scale [*]	50, 100 mmHg	0-50 mmHg		
	4, 8, 10kPa	4 kPa	Backup	Backup
	4, 8, 10%	4%		
GAS_O ₂ Scale [*]	18-30, 18-60, 18-100, 0-30, 0-60, 0-100%	18-30%	Backup	Backup
Agent Selection	ISO, HAL, ENF, SEV, DES, Auto	Auto	Backup	Backup
Agent Scale [*]	0-4, 0-8, 0-16%	4%	Backup	Backup
Flow Rate (When adult/child water trap is used.)	120, 150, 200 ml/min	200 ml/min	Backup	Backup
Flow Rate (When neonate water trap is used.)	70, 100, 120 ml/min	120 ml/min	Backup	Backup
Wave Clip [*]	ON, OFF	ON	Backup	Backup
CO ₂ Source Priority	MGU-800, HS-8000	MGU-800	Backup	Backup
	Anesthesia, HS-8000 [*]	Anesthesia	Backup	Backup

*: This setting is enabled when FLOW-i is connected.

BIS (A-2000/A-3000)

Item	Description	Default	At Power ON	At Discharge
Short Trend 2nd Parameter	SR, EMG, SQI, OFF	SR	Backup	Backup

BIS (When HBX-800 is used)

Item	Description	Default	At Power ON	At Discharge
Scale	EEG1, EEG2	±50µV	Backup	Backup
Short Trend 2nd Parameter	SR, EMG, SQI, OFF	SR	Backup	Backup
Continuous Impedance Check	ON, OFF	OFF	Backup	Backup
Smoothing Rate	10, 15, 30 sec.	15 sec.	Backup	Backup
EEG Filter	ON, OFF	ON	Backup	Backup

Stopwatch

Item	Description	Default	At Power ON	At Discharge
Label 1	9 alphanumeric characters	TIMER1	Backup	Backup
Label 2		TIMER2	Backup	Backup

Data Review

Graphic Trend

Item		Description	Default	At Power ON	At Discharge
Trend A	HR, ST (I to V6), SpO ₂ , PR_SpO ₂ , VPC, NIBP, BP1 to 8, PR_IBP, PDP, CPP, TEMP1 to 8, Tb, ΔTEMP-A to D, RR_IMP, APNEA, EtCO ₂ , InspCO ₂ , RR_GAS, ExpN ₂ O, InspN ₂ O, O ₂ , ExpAGT, InspAGT, MAC, BIS, SR, EMG, SQI, SvO ₂ , ScvO ₂ , CCO, CCI, BT, RR_VENT, RR_SpO ₂ , ExpO ₂ , Insp O ₂ , PI, PVI, SpCO, SpMet, SpHb, SpOC, PEAK, PEEP, ExpMV		Upper Row: HR, NIBP Lower Row: SpO ₂ , TEMP1, RR_IMP	Backup	Backup
Trend B			Upper Row: HR, BP1, TEMP1, NIBP Lower Row: SpO ₂ , EtCO ₂ , ST (II) , RR_CO ₂	Backup	Backup
Trend C			Upper Row: HR, TEMP1, BP1, NIBP Lower Row: SpO ₂ , InspO ₂ , EtCO ₂ , InspAGT	Backup	Backup
Trend D	-		N/A	Backup	Backup
Trend E			Upper Row: EMG, SQI Lower Row: BIS, SR	Backup	Backup
Time Range	10min, 1h, 2h,	4h, 8h, 12h, 16h, 24h	4 hours	Backup	Backup
Display Selection	$ \blacksquare, \bigtriangledown, \blacktriangle, \blacksquare, \blacksquare, +, \varkappa, $ $ \blacksquare, \blacksquare, \frown, \varkappa, $ $ \blacksquare, \blacksquare, \blacksquare, \blacksquare, \blacksquare, \blacksquare$				
Scale, Display Selection	HR, PR_SpO ₂ , PR_IBP	100, 200, 300bpm	300bpm	Backup	Backup
	ST (V to V6)	± 0.2, ± 0.5,± 1.0, ± 2.0mV ±2.0, ±5.0, ± 10.0, ±20.0mm	± 0.5mV± 5.0mm	Backup	Backup
	VPC	20, 50, 100 beats	20 beats	Backup	Backup
	BP1 to BP8	20, 50, 100, 150, 200, 300mmHg 4, 8, 16, 20, 24, 40kPa	200mmHg 24kPa	Backup	Backup

Graphic Trend

Item		Description	Default	At Power ON	At Discharge
Scale, Display Selection	PDP, CPP	20, 50, 100, 150, 200, 300mmHg 4, 8, 16, 20, 24, 40kPa	200mmHg 24kPa	Backup	Backup
	NIBP	100, 150, 200, 300mmHg 16, 20, 24, 40kPa	200mmHg 24kPa	Backup	Backup
	TEMP1 to TEMP8,	20.0-45.0, 30.0-40.0°C	30.0-40.0°C ■	Backup	Backup
	Tb	20.0-45.0, 30.0-40.0°C	20.0-45.0°C	Backup	Backup
	SpO ₂	0-100, 50-100, 80-100%	80-100%	Backup	Backup
	SpCO	0-20, 0-40, 0-100%	0-20%	Backup	Backup
	SpMet	0-10, 0-15, 0-100%	0-10%	Backup	Backup
	SpHb	10 to 20, 0 to 25 g/dL	10 to 20 g/dL	Backup	Backup
	SpOC	10 to 26, 0 to 36 mL/dL	10 to 26 mL/dL	Backup	Backup
	RR_IMP, RR_VENT, RR_GAS, RR_SpO ₂	50, 100, 150Bpm	50Bpm	Backup	Backup
	Apnea	15 sec., 30 sec.	15 sec.	Backup	Backup
	CO ₂	50, 100 mmHg 4.0, 8.0, 10.0kPa 4.0, 8.0, 10.0%	50mmHg 4.0 kPa 4.0%▲	Backup	Backup
	0 ₂	50, 100%	100%	Backup	Backup
	ΔΟ ₂	3.0, 6.0, 9.0%	3%	Backup	Backup
	N ₂ O	50, 100%	100%	Backup	Backup
	Agent	4.0, 8.0, 10.0%	8%	Backup	Backup
	PI	0-10, 0-20%	0-10%	Backup	Backup
	PVI	0-30, 0-60, 0-100%	0-30%	Backup	Backup
	PEAK	0-10, 0-20, 0-50, 0-100 cmH ₂ O	0-20 cmH ₂ O	Backup	Backup
	PEEP	0-10, 0-20, 0-50, 0-100 cmH ₂ O	0-20 cmH ₂ O	Backup	Backup
	MV	0.0-6.0, 0.0-12.0, 0.0 to 20.0.0L/min	0.0 to 12.0.0L/min	Backup	Backup
	SvO ₂ , ScvO ₂	0-100, 50-100, 80-100%	0-100%	Backup	Backup
	CCO	6, 12, 20.0L/min	6.0L/min	Backup	Backup
	CCI	6.0, 12.0, 20.0L/min/m ²	6.0L/min/m ²	Backup	Backup
	BT	20.0-45.0, 30.0-40.0°C	20.0-45.0°C	Backup	Backup

Graphic Trend

Item	Description		Default	At Power ON	At Discharge
	BIS	25, 50, 75, 100	100	Backup	Backup
	SR	25, 50, 75, 100%	100%	Backup	Backup
	SQI	0-100%	100%	Backup	Backup
	EMG	30-80dB	30-80dB	Backup	Backup
	Lt-rSO ₂	20-100	20-100	Backup	Backup
	Rt-rSO ₂	20-100	20-100	Backup	Backup
	S1-rSO ₂	20-100	20-100	Backup	Backup
	S2-rSO ₂	20-100	20-100	Backup	Backup

Tabular Trend

Item	Description	Default	At Power ON	At Discharge
Interval	10sec., 30sec., 1min., 2min., 2.5min., 5min., 10min., 15min., 30min., 60min., NIBP	5 min.	Backup	Backup
Group	A to F	A	Backup	Backup
Fixed Parameters	0 to 6 param.	0 param.	Backup	Backup
List Setup	[H Module] OFF, HR, VPC, ST (I to V6), SpO ₂ , PR_SpO ₂ , NIBP-S/D/M, BP1 to 8- S/D/M, PR_IBP, PDP, PCWP, CPP, TEMP1 to 8, Tb, CO, EtCO ₂ , InspCO ₂ , RR_GAS, RR_IMP, RR_VENT, APNEA, O ₂ -E, O ₂ -I, N ₂ O-E, N ₂ O-I, AGT-E, AGT-I, AGT2-E, AGT2-I, PI, PVI, SpCO, SpMet, SpHb, SpOC, E-TV, I-TV, E-MV, I-MV, P-PEAK, P-PAUSE, PEEO, P-MEAN, RES, COMP, TV1sec, I/E RATIO			
	[SvO ₂ /CCO] SvO ₂ , ScvO ₂ , SaO ₂ , O ₂ EI, B-Temp, CCO, CCO-STAT, CCI, CCI-STAT, DO ₂ , RVEF, RVEF-STAT, VO ₂ , SV, SV-STAT, SVI, SVI-STAT, SVR, SVRI, SVV, EDV, EDV-STAT, EDVI, EDVI-STAT, MAP, ESV, ESVI, dPmx, CO CAL, OFF			
	[Ventilator] E-TV, I-TV, MV, SMV, P-PEAK, P-PAUSE, PEEP, P-MEAN, E-RES, I-RES, COMP, FiO ₂ , P-MIN, S-COMP, D-COMP, S-RR, I/E RATIO, RES			
	[Other] BIS, SQI, EMG, SR, SEF, TOTPOW, IMP, Lt-rSO ₂ , Rt-rSO ₂ , S1-rSO ₂ , S2-rSO ₂ , tcpO ₂ , tcpCO ₂			
	Group A	HR, VPC, ST (I), ST (II), NIBP-S, NIBP-D, SpO ₂ , PR_SpO ₂ , BP1-S, BP1-D, BP1-M, BP2-S, BP2- D, BP2-M, EtCO ₂ , RR_CO ₂ , RR_IMP, APNEA, TEMP1, TEMP2	Backup	Backup
	Group B	HR, VPC, ST(I) to ST(V6)	Backup	Backup
Tabular Trend

Item	Description	Default	At Power ON	At Discharge
List Setup	Group C	HR, RR_IMP, RR_GAS, RR_VENT, SpO ₂ , P-PEAK, P-PAUSE, P-MEAN, PEEP, E-TV, I-TV, MV, E-RES, I-RES, COMP, O ₂ -I, EtCO ₂ , APNEA	Backup	Backup
	Group D	SvO ₂ , CCO, EDV, B-Temp, RVEF, SV, CCI, EDVI, ESV, SVR, SaO ₂ , SVI, ESVI, SVRI, CCO_STAT, EDV_STAT	Backup	Backup
	Group E	BIS, SQI, EMG, SR	Backup	Backup
	Group F	HR, SpO ₂ , NIBP-S, NIBP-D, NIBP-M, BP1-S, BP1-D, BP1- M, RR_GAS, EtCO ₂ , O ₂ -I, AGT-I	Backup	Backup
Filtering (Sampling Interval)	10sec., All	All	Initialize	Initialize

OCRG

Item	Description	Default	At Power ON	At Discharge
Display Duration	8, 16 min	8 min.	Backup	Backup
Waveform	Impedance, CO ₂	Impedance	Backup	Backup
Respiration Waveform Size (Impedance)	x 1/4, x1/2, x1, x2, x4	x1	Backup	Backup
Respiration Waveform Size (CO ₂)	50, 100 mmHg	50 mmHg	Backup	Backup

Recall

Item	Description	Default	At Power ON	At Discharge
Waveform	ECG1, ECG2, BP1 to 8, SpO ₂ , RESP, CO ₂ , GAS_CO ₂ , EEG1 to 8	ECG1, ECG2	Backup	Backup
Recall Factor	Asystole, VF, VT, Slow VT, Run, Bigeminy, Trigeminy, Pause, Couplet, Tachy, Brady, Frequent, ExtTachy, ExtBrady, RR IREG, Prolong RR, R ON T, TRIPLET, MLTIFORM, VENT RHYTHM, NOT CAPT, NOT PACING, S COUPLET, VPC, SVT, SVPC, S FREQUENT, HR, ST, NIBP, RR, APNEA, SpO ₂ , PR, BP1 to 8, TEMP1 to 8, Tb, CO ₂ , O ₂ , N ₂ O, AGENT, SpCO, SpMet, SpHb, PEAK, PEEP, MV	All ON	Backup	Backup
List	14 waves	14 waves	Backup	Backup

Recall

Item	Description	Default	At Power ON	At Discharge
Recall Display Selection	Asystole, VF, VT, Slow VT, Run, Bigeminy, Trigeminy, Pause, Couplet, Tachy, Brady, Frequent, ExtTachy, ExtBrady, RR IREG, Prolong RR, R ON T, TRIPLET, MLTIFORM, VENT RHYTHM, NOT CAPT, NOT PACING, S COUPLET, VPC, SVT, SVPC, S FREQUENT, HR, ST, NIBP, RR, APNEA, SpO ₂ , PR, BP1 to 8, TEMP1 to 8, Tb, CO ₂ , O ₂ , N ₂ O, AGENT, Event 1 to 8, SpCO, SpMet, SpHb, PEAK, PEEP, MV	All ON	Backup	Backup

ST Measurement

Item	Description	Default	At Power ON	At Discharge
Measurement Point	0 ms to 560 ms	120 ms	Depends on the "Main Mode" setting	Initialize
Reference Point	0 ms to -240 ms	-80 ms	under [Initial Settings>User I/F>Power ON/ Discharge].	Initialize
ST Waveform Size	x1/4, x1/2, x1, x2, x4	x1	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/	
Slide Show Interval	1, 5, 10, 20, 30 sec.	5 sec.	Disch	arge].
ST Waveform Interval	10 sec., 1 min., 5 min., 10 min.	10 sec.		

NOTE

 The graphic trend, tabular trend, alarm history will be saved even after the power is turned OFF.

• The recall, ST, OCRG data will be saved even after about one hour of standby mode (AC power or battery operation).

Item		Description	Default	At Power ON	At Discharge
ECG Analysis		Real Time, Review	Real Time	Depends on the "Main	Initialize
Limb Lead Siz	e	x1/4, x1/2, x1, x2, x4	x1	Mode" setting	
Chest Lead Si	ze	x1/4, x1/2, x1, x2, x4	x1	under [Initial Settings>User //F>Power ON/	
Filter	AC Filter	ON, OFF	OFF	Discharge].	
	EMG Filter	OFF, Strong (25Hz), Weak (35Hz)	OFF		
	Drift Filter OFF, Strong (0.50Hz), Weak (0.25Hz)		Strong (0.50 Hz)		
Background Color		White, Black	Black	setting un Settings>User	e "Main Mode" der [Initial I/F>Power ON/ arge].

12-lead Display

Basic Setup

Item		Description	Default	At Power ON At Discharge
Vital	Urgent	Volume: 11 levels	4	Depends on the "Main Mode" setting under [Setup>Initial
Alarm Sound		Tone: 5 types [*]	1	Settings>User I/F>Power ON/ Discharge].
	Caution	Volume: 11 levels	4	Dischargej.
		Tone: 5 types [*]	1	
	Status	Volume: 11 levels	4	
		Tone: 4 types [*]	1	
Ventilator	ON/OFF		OFF	
Alarm Sound	Volume: 11 l	evels	4	
	Tone: 1 type		1	
Status Alarm	Urgent	Volume: 11 levels	4	
Control Alarm Sound		Tone: 1 type [*]	1	
	Caution	Volume: 11 levels	4	
		Tone: 1 type [*]	1	
	Status	Volume: 11 levels	4	
		Tone: 1 type [*]	1	
Sync. Tone	Volume: 11 levels		2	
	Tone: 5 types		1	
	Sync. Tone: Selected Tone, Sync. with SpO ₂ Value		Selected Tone	
Key Sound	Volume: 11 I	evels	4	
	Tone: 3 type	S	1	
Other Bed Alarm	Volume: 11 levels		4	
	Tone: 1 type		1	
Boot/Shutdown	Volume: 11 l	evels	2	
Sound	Tone: 3 types		1	
Other	Volume: 11 I	evels	4	
	Tone: 1 type		1	

* When [Fukuda Tone] is selected for "Alarm System", the tone can be selected from 8 levels.

Display Configuration

Item	Description	Default	At Power ON	At Discharge
Layout	Right (1 column), Right (1 column) + Bottom, Right (2 columns), Right (2 columns) + Bottom, Left (1 column), Left (1 column) + Bottom, Left (2 columns), Left (2 columns) + Bottom, Bottom (2 rows to 6 rows)	Bottom (5 rows)	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Auto Display Config.	Туре-1, Туре-2	Type-1	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Background Color	Refer to the Color Setup.			
Palette	Refer to the Color Setup.			

Display Configuration

Item	Description	Default	At Power ON At Discharge	
Numeric Data	OFF, HR/PR, HR, PR_IBP, VPC/PACE, ST/VPC, ST-A to C, BP1 to 8, NIBP, NIBP LIST, SpO ₂ -1, SpO ₂ -1/PR_SpO ₂ -1, PR_SpO ₂ -1, RR_IMP, RR_CO ₂ , RR_VENT, RR_SpO ₂ , TEMP1 to 8, TEMP1/2, TEMP3/4, TEMP5/6, TEMP7/8, SpO ₂ -2, SpO ₂ -2/PR_SpO ₂ -2, PR_SpO ₂ -2 Δ TEMP-A to D, VENT, P-V F-V, SvO ₂ CO, BIS, CO ₂ , O ₂ , N ₂ O, Agent, RR/CO ₂ /Agent/O ₂ /N ₂ O, CO ₂ /Agent/O ₂ /N ₂ O, RR/Agent/O ₂ /N ₂ O, Agent/O ₂ /N ₂ O, Agent/N ₂ O, HEMO, HEMO-I, STOPWATCH, SpCO, SpMet, SpHb, GAS/SPIRO, SPIRO, VENT-A, VENT-B, Hemo/etc-A, Hemo/etc-B, Extended Function-A	HR, SpO ₂ -1, NIBP, BP1, RR_IMP, CO ₂	Depends on the setting unde [Initial Settings>User I/ F>Power ON/Discharge].	
Waveform	OFF, ECG1 to ECG12, ECG1 Cascade to ECG12 Cascade, BP1 to BP8, BP Overlap 1 to BP Overlap 3, SpO ₂ -1, SpO ₂ -2, RESP, AWF, AWP, AWV, CO ₂ , O ₂ , Agent, Block Cascade, RR Overlap 1 to 3, EEG1 to 8	ECG1, SpO ₂ -1, BP1, RESP, CO ₂	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Enlarged Waveform	ON, OFF	OFF		
Sweep Speed	Circ.: 6.25, 12.5, 25, 50 Vent.: 6.25, 12.5, 25	Circ.: 25 Vent.: 6.25	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
Short Graphic Trend	Short Trend Selection ON, OFF, Overlap Display Length: 0, 5, 10, 15, 20, 25, 30 min.	OFF 15 min.	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	
User Key	OFF, Home, Menu, Minimize, User Key, Main 1 to 9, Extended Display 1 Mode 1 to 3, Extended Display 2 Mode 1 to 3, Sub Mode 1 to 6, BP1 to 6 Scale, Initialize Scale, Alarm Silence, Alarm Suspend, NIBP Start/Stop, NIBP Cont., Print Start/Stop, Monitor Suspend, Night Mode, Freeze, Key Lock, Mode Select., Oxygenator Mode, Admit/Discharge, Rapid Discharge, NIBP Start/Stop, HR/PR, HR/PR Source, NIBP Cont., BP Zero, Lead, ECG Size (All Leads), Monitor Suspend, Scale, Scale (Extended Display), SpO ₂ Display ON/OFF, CO ₂ Display ON/OFF, GAS Display ON/OFF, Suspend CO ₂ , Auto Display Config., Enlarged Display, Short Trend ON/OFF, Transparent Window ON/OFF, Change Palette, Graphic Trend, Trend (Group), Tabular Trend, Tabular Trend (Group), NIBP List, Recall, OCRG, ST, Cardiac Output, Hemodynamics, Lung Function, Full Disc. Wave, 12-Lead Analysis, Tone/Volume, NIBP Auto Mode, Alarm Setup (Basic, All), Manual Printing, Display Config., Time/Date, Stopwatch, Group 1, Group 2, Group 3, Group 4, Group 5, Event, Print (LBP) Cancel, Oxygenator Mode	Home, Menu, Alarm Silence, Admit/Disch., BP Zero, NIBP Start/Stop, NIBP Auto Mode, Alarm Setup (All) NIBP Cont., Alarm History, NIBP List, Recall, Graphic/Tabular Trend, Print Start/Stop, Key Lock, Night Mode	Depends on the setting under [Initial Settings>User I/ F>Power ON/Discharge].	

Item	Desc	cription	Default	At Power ON	At Discharge	
Detail Setup (Numeric Data)	Alarm Limit Display	Graph, Numeric, OFF	Graph		e setting under ngs>User I/	
	At Alarm Occurrence	Reversed, 3D	Reversed	F>Power ON	F>Power ON/Discharge].	
Detail Setup	Grid	ON, OFF, Bold	ON			
(Wave)	Scale	ON, Bold1, Bold2	ON			
	Thickness	Thin, Regular, Thick	Regular			
	Clip	ON, OFF	ON			
	CO ₂ Wave Fill	ON, OFF	ON			
	O ₂ Wave Fill	ON, OFF	OFF	Depends on th	e setting under	
	Agent Wave Fill	ON, OFF	OFF	[Initial Setti	ngs>User I/	
	BP Overlap 1	BP1 to 8	BP1 to 4	F>Power Or	I/Discharge].	
	BP Overlap 2, 3	-	N/A			
	RR Overlap 1	CO ₂ , O ₂ , Agent	CO ₂ , O ₂ , Agent			
	BP Overlap 2, 3	-	N/A			
	12-Lead ST Short Trend	OFF, Fill, Plot	Fill			
	ST/VPC/Arrhy. Alarm Display	ON, OFF	ON	Backup	Backup	
	Block Cascade	Waveform Quantity: 2 to 6 Displayed Waveform: OFF, ECG1 to 12, BP1 to 8, SpO ₂ , RESP, AWF, AWP, CO ₂ , O ₂ , Agent	Waveform Quantity: 2 Displayed Waveforms: ECG1, ECG2	Depends on th Initial Setti	e setting under ngs>User I/	
	Waveforms	Standard, Extended	Standard		I/Discharge].	
	Graphic/Tabular Trend	ON, OFF	OFF			
	Graphic/Tabular Trend Size	Big, Medium, Small	Small			
Detail Setup (Short Trend)	Short Trend	Link with Numeric, Link with Waveform, User Setup	Link with Numeric			
	Short Trend Scale	Trend, Waveform	Trend			
	Display Parameter	ON, Gray, OFF	OFF			
	Reference Line Function	Enable, Disable	Disable	Depends on th	•	
	Cursor Function	Enable, Disable	Disable		ngs>User I/ I/Discharge].	
	Cursor Linkage	Tabular Trend, Graphic Trend, Zoom Wave	Tabular Trend			
	Short Trend Overlap 1, 2, 3		OFF, OFF, OFF, OFF			
			_			

Display Configuration

NOTE

Data Resolution

 By selecting [Backup] for "Power ON" and "Discharge" under [Initial Settings>User I/ F>Power ON/Discharge], the display configuration settings will be retained at "Power ON"

5 sec., 10 sec., 30 sec. 5 sec.

and "Discharge" respectively. If [Initialize] is selected, the settings will be initialized at "Power ON" and "Discharge".

Manual Printing

	Item	Description	Default	At Power ON At Discharge
Basic	Printer	Bedside, Central	Bedside	Depends on the "Main Mode" setting under [Initial
	Waveform	ECG1, ECG2, ECG3, BP1 to 8, SpO ₂ , RESP, CO ₂ , O ₂ , Agent, AWF, AWP	ECG1	Settings>User I/F>Power ON/ Discharge].
	Print Duration	24 sec., Cont.	24 sec.	
	Delay Time	None, 8sec., 16 sec.	8 sec.	
12-lead	12-Lead Waveform Format (Bedside)	3 wavesx4, 2 wavesx6	3 waves x 4	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/ Discharge].
	12-Lead Waveform Format (Laser)	3 wavesx4, 3 wavesx4+Rhy., 6 wavesx2, 12 waves	3 waves x 4	
	12-Lead Analysis Format (Bedside)	3 waves x 4	3 waves x 4 (fixed)	
	12-Lead Analysis Format (Laser)	6 wavesx2 (2 pages), 6 wavesx2 (1 page), 3 wavesx4+Rhythm	6 wavesx2 (2 pages)	
	Position	Center, Proportional, OFF	Proportional	
	Wave Format	Regular, Reverse	Regular	
	Printer Auto Scale	ON, OFF	ON	
	Print Calibration	ON, OFF	ON	
	Lead Boundary	ON, OFF	ON	
Other	Graphic Trend	Bedside, Central, Laser	Bedside	Depends on the "Main Mode"
Setup: Graphic	Tabular Trend	Bedside, Central, Laser	Bedside	setting under [Initial Settings>User I/F>Power ON/
Printing	OCRG	Bedside, Laser	Bedside	. Discharge].
	Zoom Wave (Recall, Full Disc.)	Bedside, Central, Laser	Bedside	
	ST	Bedside, Central, Laser	Bedside	
	12-Lead Waveform	Bedside, Laser	Bedside	
	12-Lead Analysis Result	Bedside, Laser	Bedside	
	Full Disc. Compressed Wave	Bedside, Laser	Bedside	
	Hemodynamics	Bedside, Central, Laser	Bedside	
	Lung Function	Bedside, Central, Laser	Bedside	
	СО	Bedside, Central, Laser	Bedside	
Other Setu	p: Recall Printing	Graphic Printing, Manual Printing	Graphic Printing	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/ Discharge].

Auto Printing

Item		Description	Default	At Power ON	At Discharge
Alarm	Printing	ON, OFF	OFF	Depends on the "Main Mode" setting under [Initial	
Printing	Factor	Alarm for each arrhythmia, parameter	All	Settings>User	I/F>Power ON/
	Printer	Bedside, Central	Bedside	Discharge].	
Waveform ECG1, ECG2, ECG3, BP1 to 8, SpO ₂ , RESP, CO ₂ , O ₂ , Agent, AWF, AWP, Alarm		ECG1, Alarm Factor			
	Print Duration	12 sec., 24 sec.	12 sec.		
Periodic	Periodic Printing	ON, OFF	OFF	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/ Discharge].	
Printing	Printer	Bedside, Central	Bedside		
	Waveform	ECG1, ECG2, ECG3, BP1 to 8, SpO ₂ , RESP, CO ₂ , O ₂ , Agent, AWF, AWP	ECG1		
	Periodic Interval	Interval, Timer	Timer		
	Interval	1, 2, 3, 5, 10, 15, 20, 30, 60, 120 min.	120 min.		
	Timer	0:00 to 23:00 (1:00 interval)	None	-	
	Print Duration	6, 12, 24 sec.	12 sec.	1	

Common Setup for Printing

Item	Description	Default	At Power ON At Discha	
QRS Classification	ON, OFF	ON	Depends on the setting une	
Speed	50 mm/S, 25 mm/S	25 mm/S	Settings>User I/F>Power ON/ Discharge].	
Print Calibration	Top, Each Page, OFF	OFF	Disch	aigoj.
Print NIBP Data	ON, OFF	OFF		

Other Setup

	Item	Description	Default	At Power ON	At Discharge
Night	Mode	Manual, Timer	Manual		ne "Main Mode" Inder [Initial
Mode	Start Time	00:00 to 23:59	Start Time: 21:00	Settings>User	I/F>Power ON/ arge].
	End Time	00:00 to 23:59	End Time: 07:00	Disci	laigej.
	Volume	No Change, 3, 1, 0	1		
	Display	No Change, Dark, Darker, Time Only	Darker		
	Alarm Indicator	ON, OFF	OFF		
	External Monitor Display during Night Mode	ON, OFF, OFF (Time Only)	ON		
Color	Background Color (Meas.) Background Color (Wave)	Black, Gray, Light Gray	Numeric Data: Black Waveform: Black	Depends on the "Main Mode" setting under [Initial Settings>User I/F>Power ON/ Discharge].	
	Palette	Light, Clear, Deep, Vivid	Vivid		
	HR	12 colors + White	6	-	
	ST		6		
	VPC		White		
	PACE	1	White		
	NIBP		8	1	

	Item	Description	Default	At Power ON	At Discharge
	SpO ₂ (Ch1, Ch2)		4		
	SpCO (Ch1, Ch2)		4		
	SpMet (Ch1, Ch2)		4		
	SpHb (Ch1, Ch2)		4		
	CO ₂		8		
	RESP		White		
	BP1, ART		1		
	PAP		4		
	CVP		8		
	ICP		7		
	IAP		12		
	LVP		2		
	US1 to US5 (BP)		White		
	BP2		8		
	BP3		4		
	BP4		6		
	BP5		2		
	BP6		12		
	BP7		9		
	BP8		7		
	TEMP1 to 8, Tb		2		
	Tsk, Tre, Tes, Tco, US1 to US7		2		
	AWF		6		
	AWP		4		
	AWV		8		
	BIS		2		
	INVOS		White		
	SvO ₂ +CO		White		
	Stopwatch		White		
Brightness	Brightness	7 levels	Тор	setting un Settings>User	e "Main Mode" der [Initial I/F>Power ON/ arge].
Stopwatch	1	9 alphanumeric characters	TIMER1	Backup	Backup
Label	2	1	TIMER2	Backup	Backup

Chapter 13 Accessories

Accessories	13-1
Optional Accessories	13-1
ECG, Impedance Respiration Measurement (Manufactured by	
Fukuda Denshi)	.13-1
Invasive Blood Pressure Measurement (Manufactured by Fukud	la
Denshi)	.13-2
Non-Invasive Blood Pressure Measurement (Manufactured by	
Fukuda Denshi)	.13-2
Temperature Measurement (Manufactured by Fukuda Denshi)	.13-3
Pulse Oximetry Measurement (Manufactured by Covidien)	.13-3
Pulse Oximetry Measurement (Manufactured by Masimo)	.13-4
CO Measurement (Manufactured by Fukuda Denshi)	.13-5
CO2 Concentration Measurement (Manufactured by Philips)	.13-5
CO2 Concentration Measurement (Manufactured by Covidien).	.13-6
Anesthetic Gas Concentration Measurement (Manufactured by	
Mindray Medical Sweden AB)	.13-7
BIS Measurement (Manufactured by Covidien)	.13-7
Others (Manufactured by Fukuda Denshi)	.13-8

Chapter 13 Accessories

Accessories

This section lists the accessories for the main unit.

- Use only the spare parts specified for this equipment. Otherwise, proper function cannot be executed.
- · Specifications are subject to change without prior notice to improve the quality of products.
- Power Supply Cable: CS-34 (3 m)
- DS-8400 System Operation Manual (This Manual)
- DS-8400 System Maintenance Manual
- Double Washer Sems Screw M4x30: Q'ty 4 (For fall-prevention bracket)
- Fall-Prevention Bracket
- Cable Lock

Optional Accessories

The following products are available as optional accessories for the DS-8400 System. Purchase them as required.

- Use only the accessories specified for this system. Otherwise, proper function cannot be executed.
- For quality improvement, specifications are subject to change without prior notice.

ECG, Impedance Respiration Measurement (Manufactured by Fukuda Denshi)

Item	Model Type	Note
ECG Lead Cable	CMO-07FT-3NAB	3-electrode AAMI, clip type
ECG Lead Cable	CMO-07FT-4NAB	4-electrode AAMI, clip type
ECG Lead Cable	CMO-07FT-5NAB	5-electrode AAMI, clip type
ECG Relay Cable	CIO-07CTP-3NA	3-electrode AAMI, standard type
ECG Relay Cable	CIO-07CTP-4NA	4-electrode AAMI, standard type
ECG Relay Cable	CIO-07CTP-5NA	5-electrode AAMI, standard type
ECG Lead Patient Cable	CMO-07FTP-10NAB	10-electrode AAMI, clip, type, standard type

Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)

Item	Model Type	Note
BP Relay Cable	CJO-P01B-SA3.6	1 channel, 3.6m For use with Argon Medical Devices CDX III/Press Disposable Pressure Transducers
BP Relay Cable	CJO-P01B-SB3.6	1 channel, 3.6m For use with Argon Medical Devices DTX Plus Disposable Pressure Transducers
2ch BP Relay Cable	CJO-P01B-DA0.8	2 channels, 0.8m For use with Argon Medical Devices CDX III/Press Disposable Pressure Transducers
2ch BP Relay Cable	CJO-P01B-DA4.3	2 channels, 4.3m For use with Argon Medical Devices CDX III/Press Disposable Pressure Transducers
2ch BP Relay Cable	CJO-P01B-DB0.8	2 channels, 0.8m For use with Argon Medical Devices DTX Plus Disposable Pressure Transducers
2ch BP Relay Cable	CJO-P01B-DB4.3	2 channels, 4.3m For use with Argon Medical Devices DTX Plus Disposable Pressure Transducers
2ch BP Conversion Cable	CJO-P01B-DJ0.5	2 channel-1 channel Conversion Relay Cable

REFERENCE

Argon Medical Devices: Former Becton Dickinson

)•

Non-Invasive Blood Pressure Measurement (Manufactured by Fukuda Denshi)

Item	Model Type	Note
Infant Cuff	CUF-8501	Width 8 cm to 13 cm
Pediatric Cuff	CUF-8502	Width 12 cm to 19 cm
Adult Cuff (Small)	CUF-8503	Width 17 cm to 25 cm
Adult Cuff (Medium)	CUF-8504	Width 23 cm to 33 cm
Adult Cuff (Large)	CUF-8505	Width 31 cm to 40 cm
Adult Cuff (Thigh)	CUF-8506	Width 38 cm to 50 cm
Air Hose (1.5m) General	OA-80APR1.5	For Rectus Connector Type
Air Hose (3.5m) General	OA-80APR3.5	For Rectus Connector Type
Air Hose (1.5m) Neonate	OA-80NE1.5	For SunTech Medical Neonatal Soft Disposable BP Cuff
Air Hose (3.5m) Neonate	OA-80NE3.5	For SunTech Medical Neonatal Soft Disposable BP Cuff

Temperature Measurement (Manufactured by Fukuda Denshi)

Item	Model Type	Q'ty	Note	
2ch Temperature Relay Cable	CJO-P01T-DA0.5	1	0.5m Use with YSI400 compatible probe	
2ch Temperature Relay Cable	CJO-P01T-DA4.0	1	4m Use with YSI400 compatible probe	
NOTE 700 series temperature probe cannot be used.				

Pulse Oximetry Measurement (Manufactured by Covidien)

ltem	Model Type	Note
DURASENSOR	DS-100A	Reusable For adult finger (weight of 40kg and over)
OxiMax	MAX-N	Single-Patient-Use For neonate foot/adult finger (Neonate: weight of less than 3kg, Adult: weight of 40kg and over)
OxiMax	MAX-I	Single-Patient-Use For infant toe (weight of 3 to 20kg)
OxiMax	MAX-P	Single-Patient-Use For pediatric finger (weight of 10 to 50kg)
OxiMax	MAX-A	Single-Patient-Use For adult finger (weight of 30kg and over)
OxiMax	MAX-R	Single-Patient-Use For adult nose (weight of 50kg and over)
OxiMax	MAX-FAST	Single-Patient-Use For adult/pediatric forehead (weight of 10kg and over)
SpO ₂ Relay Cable	DOC-10	3m

NOTE

• There are various types of sensors available. For details, refer to your nearest service representative.

RR_SpO₂ Measurement

Item	Model Type	Remarks
Nellcor Respiratory Sensor	10068119	For adult weighing 30 kg and above, usable only with DS-8007N

Pulse Oximetry Measurement (Manufactured by Masimo)

□SpO ₂ , PR,	PI, PVI	Measurement
-------------------------	---------	-------------

Item	Model Type	Note
Masimo SET Sensor	LNCS DCI	Reusable Sensor for Adult
Masimo SET Sensor	LNCS Adtx	Adhesive Sensor for Adult
Masimo SET Sensor	LNCS Pdtx	Adhesive Sensor for Pediatric
Masimo SET Sensor	LNCS Neo-L	Adhesive Sensor (L-Shape) for Neonate
Masimo SET Sensor	LNCS Inf-L	Adhesive Sensor (L-Shape) for Infant
Masimo SET Sensor	LNCS NeoPt-L	Adhesive Sensor (L-Shape) for Premature Neonate
Masimo RD SET Sensor	RD SET DCI	Reusable Sensor for Adult
Masimo RD SET Sensor	RD SET Adt	Adhesive Sensor for Adult
Masimo RD SET Sensor	RD SET Pdt	Adhesive Sensor for Pediatric
Masimo RD SET Sensor	RD SET Inf	Adhesive Sensor for Infant
Masimo RD SET Sensor	RD SET Neo	Adhesive Sensor for Neonate
Masimo RD SET Sensor	RD SET NeoPt	Adhesive Sensor for Premature Neonate
LNCS Patient Cable	Red LNC-04	For LNCS sensor, 1.2m
LNCS Patient Cable	Red LNC-10	For LNCS sensor, 3.0m
LNCS Patient Cable	Red LNC-14	For LNCS sensor, 4.2m
RD Patient Cable	RD SET MD20-1.5	For RD SET sensor, 0.5m
RD Patient Cable	RD SET MD20-05	For RD SET sensor, 1.5m
RD Patient Cable	RD SET MD20-12	For RD SET sensor, 3.7m

□SpO₂, PR, PI, PVI, SpMet, SpCO Measurement

Item	Model Type	Note
Masimo Rainbow Sensor	Rainbow DCI-dc3	Reusable Direct Connect Sensor for Adult (0.9m)
Masimo Rainbow Sensor	Rainbow DCI-dc8	Reusable Direct Connect Sensor for Adult (2.4m)
Masimo Rainbow Sensor	Rainbow DCI-dc12	Reusable Direct Connect Sensor for Adult (3.6m)
Masimo Rainbow Sensor	Rainbow R25	Adhesive Sensor for Adult
Masimo Rainbow Sensor	Rainbow R25-L	Adhesive Sensor (L-Shape) for Adult/Neonate
Masimo Rainbow Sensor	Rainbow R20	Adhesive Sensor for Pediatric
Masimo Rainbow Sensor	Rainbow R20-L	Adhesive Sensor (L-Shape) for Pediatric/Infant
Rainbow RC Patient Cable	Rainbow RC-1	For Rainbow Sensor, 0.3m
Rainbow RC Patient Cable	Rainbow RC-4	For Rainbow Sensor, 1.2m
Rainbow RC Patient Cable	Rainbow RC-12	For Rainbow Sensor, 3.6 m
RD Rainbow Patient Cable	RD Rainbow SET MD20-1.5	For RD SET sensor, 0.5m
RD Rainbow Patient Cable	RD Rainbow SET MD20-05	For RD SET sensor, 1.5m
RD Rainbow Patient Cable	RD Rainbow SET MD20-12	For RD SET sensor, 3.7m

□SpO₂, PR, PI, PVI, SpMet, SpHb, SpOC Measurement

Item	Model Type	Note
Masimo Rainbow ReSposable Sensor System (For Adult)	Rainbow ReSposable R2-25	ReSposable Sensor Cable (For Adult) x1 ReSposable Sensor (Adhesive Tape for Adult) x10

Item	Model Type	Note
Masimo Rainbow ReSposable Sensor System (For Child)	Rainbow ReSposable R2-20	ReSposable Sensor Cable (For Child) x1 ReSposable Sensor (Adhesive Tape for Child) x10
Masimo Rainbow ReSposable Sensor Tape (For Adult)	Rainbow ReSposable R2-25a	To be used with ReSposable sensor (adhesive tape for adult), ReSposable sensor cable (for adult), 25 per box
Masimo Rainbow ReSposable Sensor Tape (For Child)	Rainbow ReSposable R2-20a	To be used with ReSposable sensor (adhesive tape for child), ReSposable sensor cable (for child), 25 per box
Masimo Rainbow ReSposable Sensor Cable (For Adult)	Rainbow ReSposable R2-25r	To be used with ReSposable sensor tape (for adult), 5 per box
Masimo Rainbow ReSposable Sensor Cable (For Child)	Rainbow ReSposable R2-20r	To be used with ReSposable sensor tape (for child), 5 per box
Rainbow RC Patient Cable	Rainbow RC-1	For Rainbow Sensor, 0.3m
Rainbow RC Patient Cable	Rainbow RC-4	For Rainbow Sensor, 1.2m
Rainbow RC Patient Cable	Rainbow RC-12	For Rainbow Sensor, 3.6 m
RD Rainbow Patient Cable	RD Rainbow SET MD20-1.5	For RD SET sensor, 0.5m
RD Rainbow Patient Cable	RD Rainbow SET MD20-05	For RD SET sensor, 1.5m
RD Rainbow Patient Cable	RD Rainbow SET MD20-12	For RD SET sensor, 3.7m

NOTE

• SpCO and SpHb cannot be measured at the same time for all the sensors.

NOTE

• There are various types of sensors available.For details, contact your nearest service representative.

CO Measurement (Manufactured by Fukuda Denshi)

Item	Model Type	Note
Catheter Relay Cable	CJO-P01C-C2.4	
Flow-through Sensor Relay Cable	CJO-P01C-F2.4	
In-line Sensor Relay Cable	CJO-P01C-L2.4	
Injectate Probe Relay Cable	CJO-P01C-T2.4	

CO₂ Concentration Measurement (Manufactured by Philips)

\Box For HPD-800/HPD-810/HPD-820 Gas Unit I/F with Capnostat 5 CO₂ Sensor

Item	Model Type	Note
Capnostat 5 CO ₂ Sensor	1015928	
Single-Patient Use Adult Airway Adapter	6063-00	Single patient use, for ET tube sizes > 4.0 mm (10 per box)
Single-Patient Use Neonatal Airway Adapter	6312-00	Single patient use, for ET tube sizes = < 4.0 mm (10 per box)

Reusable Adult Airway Adapter 7007-00 Reusable, for ET tube sizes 7007-01 7007-00: 10 per box, 7007 Reusable Neonatal Airway Adapter 7053-00	
Pousshie Neonatel Aimuny Adoptor 7052.00 Pousshie for ET tube size	
Reusable Neonatal Airway Adapter7053-00Reusable, for ET tube sizes7053-01(7053-00: 10 per box, 7053)	

 There are various types of sampling device available. For details, refer to our service representative.

CO2 Concentration Measurement (Manufactured by Covidien)

□For HCP-800/HCP-810/HCP-820 CO₂ Gas Unit

Sampling Devices

Item	Model Type	Note	
Intubated EtCO ₂			
Filter Line H Set (Adult/Pediatric)	XS04624	For long term use	
Filter Line H Set (Infant/Neonate)	006324	For long term use	
Vital Line H Set (Adult/Pediatric)	010787	For long term use	
Vital Line H Set (Infant/Neonate)	010807	For long term use	
Non-Intubated EtCO ₂			
Smart CapnoLine Plus (Adult/Intermediate)	009818	For oral nasal, short term use	
Smart CapnoLine Plus O ₂ (Adult/Intermediate)	009822	For oral nasal, short term use	
Smart CapnoLine (Pediatric)	007266	For oral nasal, short term use	
Smart CapnoLine H Plus O ₂ (Adult/Intermediate)	010433	For oral nasal, long term use	
Smart CapnoLine H (Pediatric)	010581	For oral nasal, long term use	
Smart CapnoLine H/O ₂ (Pediatric)	010582	For oral nasal, long term use	
CapnoLine H (Adult)	008177	For nasal, long term use	
CapnoLine H (Pediatric)	008178	For nasal, long term use	
CapnoLine H (Infant/Neonate)	008179	For nasal, long term use	
Smart CapnoLine H/O ₂ (Adult)	008180	For nasal, long term use	
CapnoLine H/O ₂ (Pediatric)	008181	For nasal, long term use	

*Packaged in 25 units unless otherwise specified.

NOTE

• There are various types of sampling device available. For details, refer to our service representative.

Anesthetic Gas Concentration Measurement (Manufactured by Mindray Medical Sweden AB)

For MGU-800/810 Series, Artema Model

Sampling Devices

Item	Model Type	Note
DRYLINE Water Trap, Adult	60-13100-00	Non-sterile
DRYLINE Water Trap, Neonate	60-13200-00	Non-sterile
DRYLINE Airway Adapter, Straight	60-14100-00	Non-sterile, disposable
DRYLINE Airway Adapter, Elbow	60-14200-00	Non-sterile, disposable
DRYLINE Sampling Line, Adult	60-15200-00	Non-sterile, 2.5m, disposable
DRYLINE Sampling Line, Neonate	60-15300-00	Non-sterile, 2.5m, disposable
SPIRIT Flow sensor, Adult	60-16100-00	For MGU-810 series, single-use only
SPIRIT Flow sensor, Pediatric	60-16200-00	For MGU-810 series, single-use only

BIS Measurement (Manufactured by Covidien)

Item	Model Type	Remarks
BISx	186-0195-SF	SW 1.13
Patient Interface Cable	186-0107	
BIS Extended Use Sensor	186-0160	
BIS Pediatric Sensor	186-0200	
BIS Quatro Sensor	186-0106	

- Avoid liquid ingress to the patient interface cable (PIC). Contact of fluids with the PIC sensor connector can interfere with PIC performance.
- To minimize the risk of patient strangulation, the patient interface cable (PIC) must be carefully placed and secured.
- When installing the BISx, it should not be closely attached to the patient. Secure it on the bedside rail or pole using a clip.
- BIS sensor is disposable. Do not reuse it.
- Do not reuse the BIS sensor to other patients. It may cause cross-infection.
- The duration for one usage should be within 24 hours.

Others (Manufactured by Fukuda Denshi)

Item	Model Type	Remarks
AC Power Cable	CS-24	
Ground Cable	CE-12	
Ground Cable	CE-01A	
Display Unit Connection Cable	CJ-731B	2.5m
Display Unit Connection Cable	CJ-731C	6m
Display Unit Connection Cable	CJ-731D	10 m
Remote Control Unit	CF-820	
Printing paper	OP050-01TDR	10 per box
Ethernet Branch Cable	CJ-522A	Length 1m (For DS-LAN)
Ethernet Branch Cable	CJ-522B	Length 2m (For DS-LAN)
Ethernet Branch Cable	CJ-522C	Length 4m (For DS-LAN)
Ethernet Branch Cable	CJ-522D	Length 10m (For DS-LAN)
Ethernet Branch Cable	CJ-522E	Length 20 m (For DS-LAN)
RS-232C Cable	CJ-725	Cross Cable with Core
CF Card	FCF-128	128 MB, For data transfer (for DS-8400)
CF Card	FCF-1000	1 GB, For data transfer (for DS-8400)
SD Card	SD-1G	1 GB, For data transfer (for HS-8000)
SD Card	SD-8G	8 GB, For data transfer (for HS-8000)
SD Card	SD-16G	16 GB, For full disclosure waveform
SD Card	FSD-8GA	8 GB, For full disclosure waveform
CFast Card	FCS-64G	64 GB, For full disclosure waveform
Telemetry Transmitter Module	HLX-801 (FA) / HLX-801 (G)	
Expansion Board	CC-84	
HS Adapter for DS-8400	HSA-81	
DS-8007 Adapter	DSA-82	
Module Connection Cable	CJO-08SS0.3	module-LAN Cable 0.3m
Module Connection Cable	CJO-08SS1.5	module-LAN Cable 1.5m
Module Connection Cable	CJO-08SS3.5	module-LAN Cable 3.5m
Module Connection Cable	CJO-08SS5	module-LAN Cable 5m
Module Connection Cable	CJO-08SS10	module-LAN Cable 10m
Unit Connection Cable	CJO-09SS0.3	U-Link Cable 0.3m
Unit Connection Cable	CJO-09SS1.5	U-Link Cable 1.5m
Unit Connection Cable	CJO-09SS5	U-Link Cable 5 m
AUX Connection Cable (0.65m)	CJO-15RR0.65	Relay cable for HCP-810/HPD-810/HBX-800 Connects to the HS-8000
AUX Connection Cable (1.5m)	CJO-15RR1.5]
AUX Connection Cable (3 m)	CJO-15RR3	
AUX Connection Cable LEM (0.36)	CJO-25TR0.36	Relay cable for HCP-810/HPD-810/HBX-800 Connects to the DS-8007/HM-801
AUX Connection Cable LEM (0.65)	CJO-25TR0.65	7
AUX Connection Cable LEM (1.5)	CJO-25TR1.5	7

Item	Model Type	Remarks
AUX Connection Cable LEM (2.7)	CJO-25TR2.7	
IB Clamp Base	OAO-51A	For attaching the IB-8004 to the pole
HS Fixing Base	OAO-52A	For fixing the HSA-80, HR-800, IB-8004
HS Attachment Spacer	OAO-46A	For attaching the HSA-80 to IB-8004
HS Rail Clamp	OAO-48A	For attaching the HSA-80 to a medical rail
HS Suspension Base	OAO-49A	For suspending the HSA-80, HR-800, IB-8004
Cover Panel	OAO-45A	For extending the OAO-49A
HS Pole Clamp	OAO-50A	For attaching the HSA-80 to a pole
GCX Attachment for Monitor	OAO-70A	For attaching to the GCX arm
Gas Unit/External Output Box Mounting Bracket	OAO-100A	For attaching the gas unit/external output box to the HSA-81
Lower Trolley Unit for Monitor	OTO-13	
Upper Trolley Unit for Monitor	OAO-8400	
Storage Box for Trolley	OAO-91A	
Mouse	-	Use the product recommended by Fukuda Denshi.
Relay Cable (Straight)	CJ-726	Connects the extended display unit and the serial connector (COM A/B) on the main unit (CC-84).
External Output Box	CJO-C01Q-SJ0.3	For HS-8000 series
MGU Calibration Gas	60-12001-00	Calibration gas for MGU-800/810 series
Gas Regulator	60-12000-00	For MGU-800/810 series
Exhaust Tube	60-12120-00	For MGU-800/810 series
Lithium-ion Battery	BTO-005	

External Equipment Connection Cable

Equipment	Model Type	Remarks
SV-300	CJ-401RI-70SV3	For Status II Connector
SERVO-i / SERVO-s/ SERVO-U/ SERVO-n/ SERVO-air	CJ-402RI-70SVi	For Status II Connector
SERVO-U/n/air (Alarm Detection Only)	CJO-27DJ2	For Status II Connector
PB 740/760/840	CJ-403RI-70PB	For Status II Connector
Evita (XL, 4, dura)	CJ-402RI-70SVi	For Status II Connector
VELIA, ASTRAL	CJO-23DR2	For Status II Connector
VS ULTRA	CJO-24DR2	For Status II Connector
Vigilance, Vigilance CEDV, EV1000	CJ-406RI-70Vigi	For Status II Connector
	CJO-04RS4	For Serial Connector
Vigilance II, Vigileo	CJ-402RI-70SVi	For Status II Connector
	CJ-502	For Serial Connector
BIS	CJ-407-RI-70BIS	For Status II Connector
	CJO-03RS4	For Serial Connector
INVOS 5000C	CJ-406RI-70Vigi	For Status II Connector
	CJO-04RS4	For Serial Connector
PiCCO2	CJO-18RS5	For Serial Connector
	CJO-19RS5	For Status II Connector
FLOW-i	CJ-502	For Serial Connector
	CJ-402RI-70SVi	For Status II Connector
Magnetic Card Reader/Barcode Reader	CJ-756	For Serial Connector
PulsioFlex PC4000	CJ-725	For Serial Connector ^{*1}
TCM4	CJ-726	For Serial Connector
TCM5 FLEX	CJ-725	For Serial Connector

*1: To connect the PulsioFlex PC4000, USB to RS-232C connector cable (UC-232A, ATEN) specified by PULSION Medical Systems is required.

Chapter 14 Specification

Specification	14-1
Main Unit: DSC-8400 Series	14-1
Transport Monitor: DS-8007 and DSA-82	14-2
Display Unit: LC-8018TC/LC-8016TC	14-3
Super Unit: HS-8312N/8312M and HSA-80/HSA-81	
Expansion Unit: MGU-800/810 Series and HR-800	14-5
Expansion Module: HM-800/801, HP-800, HG-810/820	14-6
Gas Unit I/F: HPD-800/HPD-810 and CO2 Gas Unit: HCF	P-800/HCP-
810	14-7
Gas Unit I/F: HPD-820 and CO2 Gas Unit: HCP-820	14-8
Input Box: IB-8004	14-9
BISx I/F Unit: HBX-800	14-10
Performance	14-11
Measurement Unit for Each Parameter	14-21
Alarm Limit Range for Each Parameter	14-26
About the SpO2 Clinical Test	
•	-

Chapter 14 Specification

Specification

This section states the specification of this equipment.

NOTE

- The DS-8400 system comforms to waterproof/dustproof rating of IPX1 with the following system configuration.
 - Main Unit: DSC-8410
 - Display Unit: LC-8018TC/LC-8016TC
 - Super Unit: HS-8312N/HS-8312M
 - + HS Adapter: HSA-81
 - Multi Module: HM-801
 - Recorder Unit: HR-800
 - Input Box: IB-8004
 - Multi Module: HM-800
 - Multiport Module: HP-800
 - SpO₂ Module: HG-810
 - SpO₂ Module: HG-820
 - Gas Unit I/F: HPD-810 and CO₂ Gas Unit: HCP-810
 - Multigas Unit: MGU-811P

Main Unit: DSC-8400 Series

Size

360(W) x 310(H) x 255(D) mm (not including the protrusion)

Weight

7.5 kg (not including the accessory)

Environmental Conditions

Operating Temperature	10°C to 40°C
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C
Transport/Storage Humidity	10% to 95% (40°C, non-condensing) However, for the CF-820 IR Remote Control Unit, the following condition applies. 10% to 90% (38°C, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	100-240 V AC
Frequency	50 Hz / 60 Hz
Power Consumption	120 VA
Battery for Operating the Equipment	
Rated Voltage	14.4V
Rated Capacity	4100 mAh
Operation Time	1.5 hours or more (at 25°C, DSC-8410, LC-8018TC, power saving mode, DS-8007 standard measurement, continuous NIBP measurement of 15 min.)
Charging Time	2.5 hours (rapid charge-standby), 5 hours (normal charge-operation)
Usable Life	
6 years	According to self-certification (@Maintenance Manual "Periodic Replacement" P7-1)

Transport Monitor: DS-8007 and DSA-82

Size	
DS-8007	200 (W) mm x 108 (D) mm x 185 (H) mm (not including the protrusion)
DSA-82	224(W) x 76(D) x 185(H) mm (not including the protrusion)
Weight	
DS-8007	2.4 kg (not including the accessory)
DSA-82	0.7 kg (not including the accessory)
Environmental Conditions	
Operating Temperature	10°C to 40°C/50°F to 104°F
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C/14°F to 140°F
Transport/Storage Humidity	10% to 95% (40°C/104°F, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Degree of protection against electric shock	ECG/RESP, SpO ₂ ,SpCO*,SpMet*,SpHb*,TEMP,BP,CO: Type CF Applied Part NIBP, CO ₂ Concentration: Type BF Applied Part *DS-8007M only
Operation Mode	Continuous Operating Equipment
Waterproof/Dustproof	DS-8007 Main Unit: IP32 Only when temperature connector cover, USB memory slot cover, CO ₂ I/F connector cover, button cover, battery cover are attached.
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	AC 100-240 V (When DSA-81 is used) DSA-82: DC18V (Supplied from DSC-8410 main unit)
Usable Life	
6 years	According to self-certification. (@Maintenance Manual "Periodic Replacement" P7-1)

Display Unit: LC-8018TC/LC-8016TC

Size	
LC-8018TC	475 (W) mm x 307 (H) mm x 62.5 (D) mm (not including the protrusion)
LC-8016TC	410 (W) mm x 265 (H) mm x 62.5 (D) mm (not including the protrusion)
Weight	
LC-8018TC	4.5 kg (not including the accessory)
LC-8016TC	3.5 kg (not including the accessory)
Environmental Conditions	
Operating Temperature	10°C to 40°C
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C
Transport/Storage Humidity	10% to 95% (40°C, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Degree of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	DC 18 V (Supplied from DSC-8410 main unit)
Usable Life	
6 years	According to self-certification

Super Unit: HS-8312N/8312M and HSA-80/HSA-81

Size		
HS-8312N/8312M	85 (W) mm x 100 (H) mm x 200 (D) mm (not including the protrusion)	
HSA-80	85(W) mm x 68(H) mm x 188(D) mm (not including the protrusion)	
HSA-81	108(W) mm x 85(H) mm x 216(D) mm (not including the protrusion)	
Weight		
HS-8312N/8312M	1.2 kg (not including the accessory)	
HSA-80	0.2 kg (not including the accessory)	
HSA-81	0.25 kg (not including the accessory)	
Environmental Conditions		
Environmental Conditions Operating Temperature	10°C to 40°C	
	10°C to 40°C 30% to 85% (non-condensing)	
Operating Temperature		
Operating Temperature Operating Humidity Operating Atmospheric	30% to 85% (non-condensing)	
Operating Temperature Operating Humidity Operating Atmospheric Pressure Transport/Storage	30% to 85% (non-condensing) 70 kPa to 106 kPa	

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Degree of protection against electric shock	ECG /RESP (Impedance), SpO ₂ , SpCO [*] , SpMet [*] , SpHb [*] , TEMP, BP, CO: Type CF Applied Part NIBP: Type BF Applied Part *: For HS-8312M only
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	HSA-80: DC 18 V (Supplied from DSC-8410 Main Unit) HSA-81: DC 18 V (Supplied from DSC-8410 Main Unit) HS-8000 Series: DC 12 V (Supplied from DSC-8410 Main Unit via HSA-80 or HSA-81)
Usable Life	
6 years	According to self-certification (@Maintenance Manual "Periodic Replacement" P7-1)

Expansion Unit: MGU-800/810 Series and HR-800

Size			
MGU-801P	125 (W) mm x 110 (H) mm x 200 (D) mm (not including the protrusion)		
MGU-811P	125 (W) mm x 108.5 (H) mm x 200 (D) mm (not including the protrusion)		
HR-800	87 (W) mm x 108.5 (H) mm x 100 (D) mm (not including the protrusion)		
Weight			
AGO ₂ Gas Unit	MGU-801P	1.8 kg (not including the accessory)	
	MGU-811P	1.8 kg (not including the accessory)	
HR-800	0.44 kg (not including the accessory)		
Environmental Conditions			
Operating Temperature	MGU Series 10°C to 35°C		
	HR Series 10°C to 40°C		
Operating Humidity	30% to 85% (non-condensing)		
Operating Atmospheric Pressure	70 kPa to 106 kPa		
Transport/Storage Temperature	-10 °C to 60°C		
Transport/Storage Humidity	10% to 95% (40°C, non-condensing)		
Storage Atmospheric Pressure	70 kPa to 106 kF	^p a	

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Degree of protection against electric shock	Respiration Gas (MGU-800/810): Type BF Applied Part
Protection against Ignition of Flammable Gas	Not provided
Voltage	MGU-800/810 Series: DC 18 V (Supplied from DSC-8410 Main Unit)
	HR-800: DC 18 V (Supplied via DSC-8410 Main Unit or MGU-800/810 series)
Usable Life	
6 years	According to self-certification

Expansion Module: HM-800/801, HP-800, HG-810/820

Size

40 (W) mm x 100 (H) mm x 130 (D) mm (not including the protrusion)

Weight

•	
HM-800/HM-801	0.5 kg (not including the accessory)
HP-800	0.5 kg (not including the accessory)
HG-810/HG-820	0.5 kg (not including the accessory)
Environmental Conditions	
Operating Temperature	10°C to 40°C
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa

Transport/Storage Temperature	-10°C to 60°C
Transport/Storage Humidity	10% to 95% (40°C, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Degree of protection against electric shock	TEMP, BP, CO (HM-800/HM-801): Type CF Applied Part BIS: Type BF Applied Part (When connected to BISx)
	SpO ₂ , SpCO, SpMet, SpHb (HG-810): Type CF Applied Part
	SpO ₂ (HG-820) : Type CF Applied Part
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	DC 12 V (Supplied from the DSC-8410 Main Unit via IB-8004 Input Box)
Usable Life	
6 years	According to self-certification

Gas Unit I/F: HPD-800/HPD-810 and CO_2 Gas Unit: HCP-800/HCP-810

Size

HPD-800/HPD-810/HCP-800/	36(W) mm x 91(H) mm x 87(D) mm (not including the protrusion)
HCP-810	

Weight

HPD-800	0.4 kg (not including the accessory)
HPD-810	0.18 kg (not including the accessory)
HCP-800	0.28 kg (not including the accessory)
HCP-810	0.22 kg (not including the accessory)

Environmental Conditions

Operating Temperature	10°C to 40°C
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C
Transport/Storage Humidity	10% to 95% (40°C, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Degree of protection against electric shock	CO ₂ : Type BF Applied Part
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	HCP-800/HCP-810: DC12V HPD-800/HPD-810: DC5V/12V (Supplied from the DSC-8410 Main Unit via AUX connector on the HS-8000 series or HM-801)
Usable Life	
6 years	According to self-certification (@Maintenance Manual "Periodic Replacement" P7-1)

Gas Unit I/F: HPD-820 and CO_2 Gas Unit: HCP-820

Size	
------	--

0.20	
HCP-820/HPD-820	120 (W) mm x 53 (D) mm x 80 (H) mm (not including the protrusion)
Weight	
HCP-820	0.3 kg (not including the accessory)
HPD-820	0.2 kg (not including the accessory)
Environmental Conditions	
Operating Temperature	10°C to 40°C/50°F to 104°F
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C/14°F to 140°F
Transport/Storage Humidity	10% to 95% (40°C/104°F, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety- Electromagnetic Compatibility -Requirements and Tests)
Type of protection against electric shock	Class I Equipment (DS-8007 System)/Internally Powered Equipment (DS-8007 System)
Degree of protection against electric shock	CO ₂ : Type BF Applied Part
Protection against Ignition of Flammable Gas	Not provided
Waterproof/Dustproof	IPX0
Power Supply	
Voltage	HCP-820: DC 12 V HPD-820: DC 12 V / 5 V (Via DS-8007 Main Unit)
Usable Life	
6 years	According to self-certification. (@Maintenance Manual "Periodic Replacement" P7-1)

Input Box: IB-8004

Size	
IB-8004	180 (W) mm x 137.5 (H) mm x 160 (D) mm (not including the protrusion)
Weight	
IB-8004	1.3 kg (not including the accessory)
Environmental Conditions	
Operating Temperature	10°C to 40°C
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C
Transport/Storage Humidity	10% to 95% (40°C, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa

Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical Electrical Equipment- Part 1-2: General Requirements for Safety – Electromagnetic Compatibility – Requirements and Tests)
Type of protection against electric shock	Class I Equipment (During AC power operation) Internally Powered Equipment (During battery operation)
Protection against Ignition of Flammable Gas	Not provided
Power Supply	
Voltage	DC 18 V (Supplied from DSC-8410 Main Unit)
Usable Life	
6 years	According to self-certification

BISx I/F Unit: HBX-800

Size

36 (W) mm x 87 (D) mm x 91 (H) mm (not including the protrusion)

Weight

0.2 kg (not including the accessory)

Environmental Conditions	
Operating Temperature	10°C to 40°C/50°F to 104°F
Operating Humidity	30% to 85% (non-condensing)
Operating Atmospheric Pressure	70 kPa to 106 kPa
Transport/Storage Temperature	-10°C to 60°C/14°F to 140°F
Transport/Storage Humidity	10% to 95% (40°C/104°F, non-condensing)
Storage Atmospheric Pressure	70 kPa to 106 kPa
Safety	
General Standard	IEC 60601-1:2005+A1:2012 (Medical electrical equipment - Part 1: General requirements for basic safety and essential performance)
EMC Standard	IEC 60601-1-2: 2007 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests)
Type of protection against electric shock	Class I equipment (DS-8400 System)/Internally Powered Equipment (DS-8400 System)
Degree of protection against electric shock	BIS: Type BF Applied Part (When connected to BISx)
Protection against Ignition of Flammable Gas	Not provided
Waterproof/Dustproof	IPX0
Voltage	DC 12 V / 5 V

Usable Life

6 years

According to self-certification

Performance

This section states the performance of the DS-8400 system. The EMC essential performance is indicated with X.

Display Panel	
Display Device	18.5 inch TFT Color LCD (LC-8018TC)
	15.6 inch TFT Color LCD (LC-8016TC)
Resolution	18.5 inch: 1366 pixel × 768 pixel, refresh frequency 60 Hz
	15.6 inch: 1366 pixel × 768 pixel, refresh frequency 60 Hz
Function Control	Touch Screen Method
Waveform Trace	Stationary Trace
Sweep Speed	ECG/SpO ₂ /BP/EEG (6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s) RESP/ CO ₂ /O ₂ /AG (6.25 mm/s, 12.5 mm/s, 25 mm/s)
Operation	
Touch Panel	Capacitive Touch Panel
Fixed Keys	3 keys (NIBP Start/Stop, NIBP Auto Mode, Alarm Silence)
Sound Pressure	
Alarm Sound (Standard Tone)	Maximum: 75.0 dB, Minimum: 45.0 dB
HR Synchronized Tone	Maximum: 86.0 dB, Minimum: 36.0 dB
SpO ₂ Synchronized Tone	Maximum: 83.0 dB, Minimum: 39.0 dB
Clock Accuracy	
	±2 min. per year (25°C)
Alarm	
Alarm Function	For each alarm level, the respective alarm sound generates, and the alarm indicator flashes.
Alarm Indicator	Visual check is possible from 4m distance.
Alarm Display	Visual check is possible from 1 m distance.

ECG

Х

Lead Type	Wired 3, 4, 5, 10-electrode
Frequency Characteristic (HS-8000)	150Hz/40Hz/15Hz (3, 4, 5, 10-electrode)
Input impedance	2.5 M Ω or above
Maximum Input Voltage	10 mVp-p
Polarization Voltage	±825 mV or above
Common Mode Rejection Ratio	90 dB or above
HR Measurement Range	Adult: 0, 12 bpm to 300 bpm Neonate: 0, 30 bpm to 300 bpm
HR Measurement Accuracy	±3 bpm
HR Display Response Time	Adult/Child: 6 sec., Neonate: 3 sec.
Instant HR	Calculated each second based on the latest RR interval.
Waveform Size Selection	1/4, 1/2, 1, 2, 4
Accuracy of Input Signal Reproduction	Overall system error and frequency response is set using method A, B, C, and D.
Defibrillation Proof	Provided:
Lead-off Detection Current	100nA and below
Heart rate meter accuracy and response to irregular rhythm	80 bpm Ventricular Bigeminy : 80 bpm
	60 bpm Ventricular Bigeminy : 60 bpm

120 bpm Ventricular Bigeminy : 120 bpm

ላ

, L

90 bpm Bidirectional Systoles : 90 bpm

Response time of heart rate meter to change in heart rate HR change from 80 bpm to 120 bpm: Range 4.7 sec. to 5.1 sec., Average 4.8 sec.

HR change from 80 bpm to 40 bpm: Range 5.0 sec. to 5.5 sec., Average 5.3 sec. Time to ALARM for tachycardia

Ventricular Tachycardia 1 mVpp, 206 bpm: Range 7.3 sec. to 8.1 sec., Average 7.6 sec.

Ventricular Tachycardia 2 mVpp, 206 bpm: Range 7.4 sec. to 8.2 sec., Average 7.7 sec.

Ventricular Tachycardia 0.5 mVpp, 206 bpm: Range 8.5 sec. to 9.4 sec., Average 8.9 sec.

Ventricular Tachycardia 2 mVpp, 195 bpm: Range 5.0 sec. to 5.4 sec., Average 5.2 sec.

MMMMMM

Ventricular Tachycardia 4 mVpp, 195 bpm: Range 4.1 sec. to 5.8 sec., Average 5.0 sec.

Ventricular Tachycardia 1 mVpp, 195 bpm: Range 6.3 sec. to 8.0 sec., Average 7.0 sec.

Active Noise Suppression RL Drive Maximum 10.8 mV

5 μV/LSB and below 100 μs and below

1.2 mV T-wave can be removed when tested according to IEC 60601-2-27.

3.2 sec, 0.3 sec, 0.1 sec (time constant can be changed)
a) Pacemaker Pulse without Over/Undershoot Capable to reject pulses of pulse width 0.1 ms to 2 ms, amplitude ±2 mV to ±700 mV

 b) Pacemaker Pulse with Over/Undershoot Rejection is not possible.

c) Pacer Pulse Detector Rejection of Fast ECG Signals Slew Rate 3.2V/S

Analog Front End: 8000 samples/s/channel Digital Signal Processing: 500 samples/s/channel and above (without skew)

Resolution Skew

Sampling Rate

Tall T-wave Rejection

Transient Characteristic

Rejection of Pacemaker

Capability

Pulse

Respiration (HS-8000)

Method	Impedance Method	
Frequency Characteristic	1.5 Hz (adult, child) / 2.5 Hz (neonate)	
Current	100 μA and below (at 66.65 kHz±5%)	
Measurement Range	0, 4 Bpm to 150 Bpm	
Measurement Accuracy	±3 Bpm	
Respiration (DS-8007)		
Method	Impedance Method	
Frequency Characteristic	1.5 Hz (adult, child) / 2.5 Hz (neonate)	
Current	100 µA and below (at 33.3 kHz±5%)	
Measurement Range	0, 4 Bpm to 150 Bpm	
Measurement Accuracy	±3 Bpm	

	Temperature	
	Measurement Method	Thermistor Method
	Probe	400 only
	Measurement Range	0°C to 45°C
	Measurement Accuracy	±0.2°C at 25°C to 45°C Outside above range ±0.4°C
	No. of Channels	Maximum 8 channels
	Temperature Delay Time (From temperature probe to monitor display)	10 sec. or less (Not including the time constant of temperature probe.)
	Operating Mode	Direct Mode
	SpO ₂ (Arterial Oxygen Satura	ation)
	Measurement Value Update Rate	1 sec.
	Nellcor Unit	
	Measurement Method	2 Wavelength Pulse Wave Method Wavelength: Approx. 660 nm (red light) Approx. 900 nm (infrared light) Output: 15 mW and below
	Measurement Range	1%SpO ₂ to 100%SpO ₂
	Resolution	1%SpO ₂
*	Measurement Accuracy	Adult: $\pm 3\%$ SpO ₂ when 70%SpO ₂ to 100%SpO ₂ (When DS-100A is used) Neonate: $\pm 2\%$ SpO ₂ when 70%SpO ₂ to 100%SpO ₂
	PR Measurement Range	20 bpm to 250 bpm
	PR Accuracy	± 3 bpm when 20 bpm to 250 bpm
	Measurement Response Time	6 sec. to 7 sec. (averaging duration)
	Respiration Rate (Pulse Wave Analysis)	
	Display Range	4 Bpm to 40 Bpm
×	RR Measurement Accuracy	Mean Error: Within ±1 Bpm
		Mean Square Deviation: Below 3 Bpm
	(<u>NOTE</u>	

• RR_SpO₂ measurement (pulse wave analysis) is optional function of the DS-8007.

Masimo Unit

Measurement Method	2 Wavelength Pulse Wave Method Masimo LNOP/LNCS Sensor Wavelength: Approx. 660 nm (red light) Approx. 905 nm (infrared light)	
	Output: 15 mW and below Masimo Rainbow Sensor Wavelength: 12 different wavelengths are used within the range of 620 nm to 1270 nm Output: 25 mW and below	
	SpO ₂	
---	-------------------------------	--
	Measurement Range	1%SpO ₂ to 100%SpO ₂
	Resolution	1%SpO ₂
Ж	Measurement Accuracy	Adult: ±2%SpO ₂ when 70%SpO ₂ to 100%SpO ₂ Neonate: ±3%SpO ₂ when 70%SpO ₂ to 100%SpO ₂
	SpCO	
	Measurement Range	0%SpCO to 99%SpCO
	Resolution	1%SpCO
	Measurement Accuracy	±3%SpCO (SpCO: 1%SpCO to 40%SpCO)
	SpMet	
	Measurement Range	0%SpMet to 99.9%SpMet
	Resolution	0.1%SpMet
	Measurement Accuracy	±1%SpMet (SpMet: 1%SpMet to 15%SpMet)
	SpHb	
	Measurement Range	0 g/dL to 25.0 g/dL
	Resolution	0.1 g/dL
	Measurement Accuracy	±1 g/dL (SpHb: 8 g/dL to 17 g/dL)
	PI (Perfusion Index)	
	Measurement Range	0.02% to 20% (disposable sensor), $0.15%$ to 20% (reusable sensor)
	Minimum Display Unit	0.01%
	PVI (Pleth Variability Index)	
	Measurement Range	0 to 100%
	Calculation Time	15 sec.
	SpOC	
	Measurement Range	0 ml/dL to 35 ml/dL
	Minimum Display Unit	0.1 ml/dL
	Pulse Rate	
	Measurement Range	26 bpm to 239 bpm
Х	Measurement Accuracy	± 3 bpm when 26 bpm to 239 bpm (without body motion)
	Measurement Response Time	7 levels 2 to 4 sec., 4 to 6 sec., 8 sec., 10 sec., 12 sec., 14 sec., 16 sec. (averaging duration)
	(NOTE	

The SpO₂ measurement accuracy is determined based on the values of the root-mean-square (rms) difference between SpO₂ readings of the pulse oximeter equipment and values of SaO₂ determined with a CO-oximeter, by healthy adult volunteers. The pulse oximeter equipment measurements are statistically distributed; ±2% measurement accuracy means that only about two-thirds of pulse oximeter equipment measurements can be expected to fall within ±2% of the value measured by a CO-oximeter.

• PVI, SpCO, SpMet, SpHb, SpOC measurements are optional functions.

BP

Transducer Sensitivity	5 µV / V / mmHg
Measurement Range	-50 mmHg to 300 mmHg
Frequency Characteristic	DC 6 Hz / 8Hz / 12Hz / 40Hz
Measurement Accuracy	Within $\pm 2\%$ or ± 1 mmHg of full scale, whichever is greater
Zero Balance Range	Within ±150 mmHg
PR Measurement Range	Adult: 12 bpm to 300 bpm Neonate: 30 bpm to 300 bpm
PR Accuracy	Within ± 3% or 1bpm, whichever is greater
No. of Channels	Maximum 8 channels

NIBP (Non-Invasive Blood Pressure) (AAMI SP10: 2002+A1:2003+A2:2006+(R)2008 Manual, electronic or automated sphygmomanometers) (ISO81060-2:2013 Non-invasive sphygmomanometers -- Part 2: Clinical investigation of automated measurement type)

	Measurement Method	Oscillometric Method
	Pressure Measurement Range	0 mmHg to 300 mmHg
	Resolution	1 mmHg
	Pressure accuracy	±3 mmHg
	BP Measurement Error according to th	e Clinical Performance Test
	Mean Error	Within ±5 mmHg
	Standard Deviation of Error	8 mmHg and below
	Error of Cuff Pressure Display	Within ±3 mmHg
Ж	Measurement Error (including simulator)	±10 mmHg
	PR Measurement Range	40 bpm to 240 bpm
	PR Accuracy	±5%
	Deflation Speed	5±1 mmHg/sec. (Quick Measurement OFF) 10±2 mmHg/sec. (Quick Measurement ON)
	Safety Mechanism	Adult: 300 mmHg and below Child: 210 mmHg and below Neonate: 150 mmHg and below

Chapter 14 Specification

CO₂ (Carbon Dioxide Concentration)

Philips Capnostat 5 (Gas Unit I/F and Mainstream Module) Measurement Method Infra-Red Solid-State Method, Mainstream Method 0 mmHg to 150 mmHg Measurement Range Measurement Accuracy 0 mmHg to 40 mmHg: ±2 mmHg 41 mmHg to 70 mmHg: ±5% 71 mmHg to 100 mmHg: ±8% 101 mmHg to 150 mmHg: ±10% CO₂ value error compensation when interference gas is present 0 mmHg to 40 mmHg: Additional error of ±1mmHg 41 mmHg to 70 mmHg: Additional error of ±2.5% 71 mmHg to 100 mmHg: Additional error of ±4% 101 mmHg to 150 mmHg: Additional error of ±5% These are maximum error only if compensation of atmospheric pressure, O2, N2O, anesthetic agent are properly performed. **RR** Measurement Range 0 Bpm to 150 Bpm **RR** Measurement Accuracy ±1 Bpm **Rise Time** 60 ms and below Covidien Unit Measurement Method Infra-Red Solid-State Method, Microstream Method Measurement Range 0 mmHg to 99 mmHg 0 mmHg to 38 mmHg: ±2 mmHg X Measurement Accuracy 39 mmHg to 99 mmHg: ± { 0.05 x displayed value +0.08x (displayed value - 39 mmHg) } (RR: 80 Bpm and below) The larger of ± 4 mmHg or ±12% : (RR: over 80 Bpm) Variation of Measurement ±2 mmHg (Within 6 hours after power ON) Accuracy CO₂ measurement accuracy when interference gas is present 0 mmHg to 38 mmHg: ± (2 mmHg + 0.04 x displayed value) 39 mmHg to 99 mmHg: ± { 0.09 x displayed value + 0.08 x (displayed value - 39 mmHg) } **RR** Measurement Range 0 Bpm to 150 Bpm **RR** Measurement Accuracy 0 Bpm to 70 Bpm: ±1 Bpm 71 Bpm to 120 Bpm: ±2 Bpm 121 Bpm to 150 Bpm: ±3 Bpm Flow Rate 50 mL/min +15, -7.5 mL/min. System Response Time 4.2 sec. Delay Time 4.0 sec. **Rise Time** 0.2 sec. CO Measurement Method Thermodilution Method Measurement Range 0.1 L/min to 20.0 L/min Measurement Range and Accuracy **Blood Temperature** ±0.3°C at 17°C to 45°C/±0.5°F at 63°F to 113°F Injectate Temperature ±0.5°C at -1°C to 35°C/±0.9°F at 30°F to 86°F BIS BISx (Covidien) Bispectral Index (BIS)

Measurement Range	0% to 100%
Resolution	1%
Signal Quality Index (SQI)	
Measurement Range	0% to 100%
Resolution	0.1%
EMG	
Measurement Range	25 dB to 100 dB
Bar Graph Display Range	30 dB to 55 dB
Graphic Trend Display Range	30 dB to 80 dB
Resolution	1 dB
Suppression Ratio (SR)	
Measurement Range	0% to 100%
Spectral Edge Frequency	
Measurement Range	0.5 Hz to 30 Hz
Total Power (TOTPOW)	
Measurement Range	40 dB to 100 dB
Waveform Display Scale	±25µV, ±50µV, ±100µV, ±250µV ±10%
Frequency Characteristic	
Filter ON	2.0 Hz to 70 Hz, AC Filter ON (50 Hz or 60 Hz)
Filter OFF	0.25 Hz to 100 Hz, AC Filter OFF
Anesthetic Agent Concentration	MGU-800/MGU-810
Sidestream Method	
Measurement Method	CO_2 , N ₂ O, Volatile Anesthetic: Infra-Red Solid-State Method O_2 : Paramagnetic Method
Warm-Up Time	Multigas Module ISO Accuracy: 45 sec. Full Accuracy: 10 sec.
Auto Zeroing	Multigas Module ISO Accuracy: 30 sec. Full Accuracy: 4 hours
Measurement Range	
CO ₂ :	0 to 10.0% (0 to 76mmHg, 0 to 10kPa)
N ₂ O:	0 to 100%
O ₂ :	0 to 100%
AG Halothane:	0 to 5%
AG Enflurane:	0 to 5%
AG Isoflurane:	0 to 5%
AG Sevoflurane:	0 to 8%
AG Desflurane:	0 to 18%
RR:	0, 2 to 100Bpm
Measurement Accuracy	
CO ₂ :	0 to 1[vol%]: ±0.1[vol%] 1 to 5[vol%]: ±0.2[vol%] 5 to 7[vol%]: ±0.3[vol%] 7 to 10[vol%]: ±0.5[vol%]

N ₂ O:	0 to 20[vol%]: ±2[vol%] 20 to 100[vol%]: ±3[vol%]
0 ₂ :	MGU-801P/MGU-811P 0 to 25[vol%]: ±1[vol%] 25 to 80[vol%]: ±2[vol%] 80 to 100[vol%]: ±3[vol%]
Volatile Anesthetic	Halothane, enflurane, and isoflurane
	0 to 1[vol%]: ±0.15[vol%] 1 to 5[vol%]: ±0.2[vol%]
	Sevoflurane 0 to 1[vol%]: ±0.15[vol%] 1 to 5[vol%]: ±0.2[vol%] 5 to 8[vol%]: ±0.4[vol%]
	Desflurane 0 to 1[vol%]: ±0.15[vol%] 1 to 5[vol%]: ±0.2[vol%] 5 to 10[vol%]: ±0.4[vol%] 10 to 15[vol%]: ±0.6[vol%] 15 to 18[vol%]: ±1.0[vol%]
Respiration Rate	±1 bpm when below 60 bpm
Respiration Detection	Changes with CO ₂ level in 1[vol%].

Interference from other gases

		Interference to Measurement Data [vol%]				
Interference Gas or Vapor	CO ₂	N ₂ O	0 ₂	Volatile Anesthetic		
CO2*1*2	-	0.1	0.2	0		
N ₂ O *1*2	0.1	-	0.2	0.1		
02 ^{*1*2}	0.1	0.1	-	0.1		
Volatile Anesthetic ^{*1*2}	0.1	0.1	1.0	Secondary 0.1 (Average)		
<100% Xenon	0.1	0	0.5	0		
<50% Helium	0.1	0	0.5	0		
Metered dose inhaler propellants	Not specified.	Not specified.	0.5	Not specified.		
<0.1% Ethanol	0	0	0.5	0		
Saturated Isopropanol Vapor	0.1	0	0.5	0		
<1% Acetone	0.1	0.1	0.5	0		
<1% Methane	0.1	0.1	0.5	0		

*1: This is the maximum influence within the gas level of specified measurement accuracy. The total influence will not exceed 5% of the gas level.

*2: For CO_2 , N_2O , O_2 , the influence from mixed agent is the same as that from single agent.

Threshold

 Volatile
 Primary 0.15[vol%] (Full Accuracy)

 Anesthetic
 0.4[vol%] (during warm-up)

 For halothane, add 0.1[vol%] to above value.

 Secondary 0.3[vol%] (Full Accuracy)

 0.5[vol%] (during warm-up)

 If primary agent is larger than 10[vol%], 5% of primary gas level.

 (10% for isoflurane)

 For halothane, add 0.1[vol%] to above value.

 70 mL/min to 200 mL/min

 ±10 mL/min or ±10%, whichever is greater

Flow Rate

14-19

Delay Time	4	s (When genuine accessory is used)
Rise Time (When	genuine access	sory is used)
CO ₂	2	50 ms (Fall Time 200 ms)
N ₂ O	2	50ms
O ₂	A	t flow rate of 200 mL/min: 500 ms (15% to 21%) , 700 ms (21% to 60%)
	A	t flow rate of 120 mL/min: 600 ms (15% to 21%) , 800 ms (21% to 60%)
Halotl	nane, Isoflurane	, Sevoflurane, Desflurane, Enflurane 300 ms
Enflur	ane 3	50 ms
DRYLINE Water	. А	imptying interval (half full, worst case) dult/Child: 17 hours @ 200 mL/min, 37°C, 100% RH leonate: 20 hours @ 120 mL/min, 37°C, 100% RH
Spirometry (MG	U-810 series)	
AWP [cmH ₂ O]		
Meas	urement Range:	-20 cmH ₂ O to 100 cmH ₂ O (Adult, Pediatric*)
		±1 cmH ₂ O (Adult, Pediatric*)
AWF (both direct.)[L/min]	
Meas	urement Range:	1.5 L/min to 100 L/min (Adult), 0.25 L/min to 25 L/min (Pediatric*)
Tidal Volume (ins	p. and exp.) [mL]
Meas	urement Range:	: 150 mL to 2000 mL (Adult), 15 mL to 300 mL (Pediatric*)
	Accuracy:	±6% or 30 mL, whichever is greater (Adult), ±6% or 4 mL, whichever is greater (Pediatric*)
Minute Ventilatior	Volume (insp. a	and exp.)[L/min]
Meas	urement Range:	2 L/min to 20 L/min (Adult), 0.5 L/min to 5 L/min (Pediatric*)
Compliance [mL/o		
Meas	urement Range:	4 mL/cmH ₂ O to 100 mL/cmH ₂ O (Adult), 1 mL/cmH ₂ O to 100 mL/cmH ₂ O (Pediatric*)
Airway Resistanc	e [cmH ₂ O/L/s]	
Meas	urement Range:	0 cmH ₂ O/L/s to 40 cmH ₂ O/L/s (Adult, Pediatric*)
Peak, Plateau, PE	EP, and Mean	Pressure [cmH ₂ O]
Meas	urement Range:	-20 cmH ₂ O to 100 cmH ₂ O (Adult, Pediatric*)
I:E Ratio		
Meas	urement Range:	: 1:4.5 to 2:1 (Adult, Pediatric*)
Conditions of Use	for Stated Accu	Jracy
Respir	ation Rate (RR):	4 Bpm to 35 Bpm (Adult), 4 Bpm to 50 Bpm (Pediatric*)
	I:E Ratio:	1:4.5 to 2:1 (Adult, Pediatric*)
	ntubation Tube:	5.5 mm to 10 mm (Adult), 3 mm to 6 mm (Pediatric*)
*Including neonat	э.	
Printing (Record	ler Unit)	
Printing Speed		50 mm/s, 25 mm/s (Error: within ±5%)
Resolution		Head Direction: 8 dots/mm Feed Direction: 40 lines/mm (at printing speed of 25mm/s)
Printing Wavefor	ms 3	3 waveforms
Printing Type	١	Naveform, List, Graphic
Detection	F	Paper out, printhead temperature

Input Box (IB-8004)	
Connectable Units	Maximum 1 unit
Number of Slots	Maximum 4 slots
Analog Waveform Output	
Output Voltage	ECG Output 1 V/mV (fixed), BP Output 1 V/100 mmHg (fixed)
Output Voltage Accuracy	within ±10% (Both ECG and BP output)
Analog Output Frequency Range	ECG Output: 0.5 Hz to 40 Hz
	BP Output: DC to 40 Hz
Delay Time	35 ms and below (ECG waveform) 35 ms and below (BP waveform: when 40 Hz is set for waveform filter)
Output Impedance	100Ω±10%
Load Impedance	1kΩ to ∞
Pacemaker Pulse	None
QRS Synchronization Outpu	ıt
Output Waveform	Square Wave (Positive/negative logic can be selected.)
Output Voltage	+4.3 V to +5.0 V (High Level) +0.3 V and below (Low Level)
Synchronized Signal Width	100 ms / 60 ms / 20 ms (Selectable)
Delay Time	35 ms and below (when the "Filter" setting is [Monitor] or [Diag.])
Output Impedance	Open Collector Output (with +5 V 500 Ω pull-up resistor)
	me of analog waveform output and QRS synchronization output depends on the

- The delay time of analog waveform output and QRS synchronization output depends on tr filter setting and the input waveform type. For details, refer to your nearest service representative.
- The QRS synchronized signal is not intended to be used as synchronized signal for defibrillator. When using the QRS synchronized signal, refer to your nearest service representative.

Measurement Unit for Each Parameter

The measurement units for this equipment are as follows.

Description	Parameter	Display	Unit	Default
Heart Rate / Pulse Rate	ECG	HR	bpm (beats per minute)	
	Blood Pressure	PR_IBP	bpm	
	SpO ₂	PR_SpO ₂	bpm	
ST Level	ECG	ST	mm, mv	mm
VPC	ECG	VPC	beat/minute	
		PACE	beat/minute	
Respiration Rate	Impedance	RR_IMP	Bpm (breaths per minute)	
	Ventilator	RR_VENT	Bpm	
	CO ₂	RR_GAS	Bpm	
	SpO ₂	RR_SpO ₂	Bpm	

Description	Parameter	Display	Unit	Default
Apnea	Impedance	APNEA	s (second)	
	CO ₂	APNEA	s (second)	
	Ventilator	APNEA	s (second)	
Blood Pressure	Blood Pressure	BP	mmHg, kPa cmH ₂ O (CVP only)	mmHg
Non-Invasive Blood Pressure	Non-Invasive Blood Pressure	NIBP	mmHg, kPa	mmHg
Arterial Oxygen Saturation	SpO ₂	SpO ₂	%	
Perfusion Index	Perfusion Index	PI	%	
	Pleth Variability Index	PVI	%	
Carboxyhemoglobin Concentration	SpCO	SpCO	%	
Methemoglobin Concentration	SpMet	SpMet	%	
Total Hemoglobin	SpHb	SpHb	g/dL	
Arterial Oxygen Saturation	SpOC	SpOC	mL/dL	
Temperature	Temperature	TEMP	°C	
End Tidal CO ₂ Concentration	CO ₂	EtCO ₂	mmHg, kPa, %	mmHg
Inspiratory CO ₂ Concentration	CO ₂	InspCO ₂	mmHg, kPa, %	mmHg
со	СО	СО	L/minute	
Blood Temperature	Blood Temperature	Tb	°C	
Injectate Temperature	Injectate Temperature	Ti	°C	
Airway Flow	Airway Flow	AWF	L/minute	
Airway Pressure	Airway Pressure	AWP	cmH ₂ O	
Ventilatory Volume	Ventilatory Volume	AWV	mL	
Tidal Volume	Expiratory Tidal Volume	E-TV	mL	
	Inspiratory Tidal Volume	I-TV	mL	
	Ventilatory Volume per second	TV/1Sec	%	
Minute Ventilation	Minute Ventilation Volume	MV	L/minute	
Volume	Spontaneous Minute Volume	SMV	L/minute	
Compliance	Compliance	COMP	mL/cmH ₂ O	
Airway Resistance	Expiratory Resistance	E-RES	cmH ₂ O/L/sec	
	Inspiratory Resistance	I-RES	cmH ₂ O/L/sec	
Airway Pressure	Mean Airway Pressure	MEAN	cmH ₂ O	
	Peak Airway Pressure	PEAK	cmH ₂ O	
	Pause Airway Pressure	Pause	cmH ₂ O	
	Plateau Pressure	PLATEAU	cmH ₂ O	
Peak End Expiratory Pressure	Peak End Expiratory Pressure	PEEP	cmH ₂ O	
Fraction of Inspiratory Oxygen	Fraction of Inspiratory Oxygen	FIO ₂	%	

Description	Parameter	Display	Unit	Default
	Mixed Venous Oxygen Saturation	SvO ₂	%	
	Central Venous Oxygen Saturation	ScvO ₂	%	
	Arterial Oxygen Saturation	SaO ₂	%	
	Oxygen Uptake Index	O ₂ EI	%	
	Oxygen Transport	DO ₂	mL/minute	
	Oxygen Consumption	VO ₂	mL/minute	
	Stroke Volume	SV	mL/beat	
	Stroke Volume (STAT Mode)	SV_STAT	mL	
	Stroke Volume Index	SVI	mL/m ²	
	Stroke Volume Index (STAT Mode)	SVI_STAT	mL/m ²	
	HR	HR	bpm (beats per minute)	
	Mean Arterial Pressure	MAP	mmHg	
	Central Venous Pressure	CVP	mmHg	
/igilance Data Vigilance	Continuous Cardiac Output	ссо	L/minute	
Vigilance CEDV	Continuous Cardiac Output (STAT Mode)	CCO_STAT	L/minute	
Vigilance II Vigileo	Continuous Cardiac Index	CCI	L/minute/m ²	
0	Continuous Cardiac Index (STAT Mode)	CCI_STAT	L/minute/m ²	
	Systemic Vascular Resistance	SVR	dyn-sec-cm ⁻⁵	
	Systemic Vascular Resistance Index	SVRI	(dyn-sec-cm ⁻⁵ -m ²)	
	Blood Temperature	ВТ	°C, °F	°F
	Ejection Fraction	RVEF	%	
	Ejection Fraction (STAT Mode)	RVEF_STAT	%	
	End-Diastolic Volume	EDV	mL	
	End-Diastolic Volume (STAT Mode)	EDV_STAT	mL	
	End-Diastolic Volume Index	EDVI	mL/m ²	
	End-Diastolic Volume Index (STAT Mode)	EDVI_STAT	mL/m ²	
	End-Systolic Volume	ESV	mL	
	End-Systolic Volume Index	ESVI	mL/m ²	
	Stroke Volume Variance	SVV	%	
				1

Description	Parameter	Display	Unit	Default
Multigas Unit	End-tidal Carbon Dioxide	CO ₂ -E	mmHg, kPa, %	mmHg
	Inspired Carbon Dioxide	CO ₂ -I	mmHg, kPa, %	mmHg
	End Tidal Oxygen	0 ₂ -E	%	
	Fraction of Inspiratory Oxygen	O ₂ -I	%	
	Expired Nitrous Oxide	N ₂ O-E	%	
	Inspired Nitrous Oxide	N ₂ O-I	%	
	End Tidal Anesthetic Gas	AGT-E	%	
	Inspired Anesthetic Gas	AGT-I	%	

Description	Parameter	Display	Unit	Default
BIS Data	Bispectral Index	BIS	(no unit)	
	Signal Quality Index	SQI	%	
	Electromyograph	EMG	dB	
	Suppression Ratio	SR	%	
	Spectral Edge Frequency	SEF	Hz	
	Total Power	TOTPOW	dB	
	Impedance	IMP	Kohms	
INVOS 5100C Monitor Data	Regional Cerebral Oxygen Saturation (Left)	Lt-rSO ₂	%	
	Regional Cerebral Oxygen Saturation (Right)	Rt-rSO ₂	%	

Description	Parameter	Display	Unit	Default
PiCCO Data	Pulse Contour Cardiac Output	CCO	L/min	
	Pulse Contour Cardiac Output Index	CCI	L/min/m ²	
	Stroke Volume	SV	mL	
	Stroke Volume Index	SVI	mL/m ²	
	Stroke Volume Index	SVV	%	
	Systemic Vascular Resistance	SVR	dyn x s x cm ⁻⁵	
	Systemic Vascular Resistance Index	SVRI	dyn x s x cm ⁻⁵ x m ²	
	Central Venous Oxygen Saturation	ScvO ₂	%	
	Oxygen Delivery	DO ₂	ml/min	
	Oxygen Consumption	VO ₂	ml/min	

Description	Parameter	Display	Unit	Default
	Pulse Contour Cardiac Output	ссо	L/minute	
	Pulse Contour Cardiac Output Index	ССІ	L/minute/m ²	
	Stroke Volume	SV	mL/beat	
	Stroke Volume Index	SVI	mL/m ²	
	Stroke Volume Variation	SVV	%	
	Systemic Vascular Resistance	SVR	dyn x sec x cm ⁻⁵	
	Systemic Vascular Resistance Index	SVRI	dyn x sec x cm ⁻⁵ x m ²	
	Central Venous Oxygen Saturation	ScvO ₂	%	
PulsioFlex Data	Oxygen Delivery	DO ₂	mL/minute	
	Oxygen Consumption	VO ₂	mL/minute	
	index of Left Ventricular Contractility	dPmx	mmHg/sec	
	Calibrated Cardiac Output	CO CAL	L/min	
	Heart Rate	HR	bpm (beats per minute)	
	Mean Arterial Pressure	MAP	mmHg	
	Global End-Diastolic Volume	GEDV	mL	
	Global End-Diastolic Volume Index	GEDI	mL/m ²	
	Extravascular Lung Water	EVLW	mL	
	Extravascular Lung Water Index	ELWI	mL/kg	
	Pulmonary Vascular Permeability Index	PVPI		
	Global Ejection Fraction	GEF	%	
	Cardiac Function Index	CFI	1/min	
	Blood Temperature	BT	°C, °F	°F
	Oxygen Delivery Index	DO ₂ I	mL O ₂ /min/m2	
	Oxygen Consumption Index	VO ₂ I	mL O ₂ /min/m2	

Description	Parameter	Display	Unit	Default
TCM4, TCM5 FLEX Data	Transcutaneous Oxygen Partial Pressure	tcpO ₂	mmHg, kPa	*
	Transcutaneous Carbon Dioxide Partial Pressure	tcpCO ₂	mmHg, kPa	*

*: The measurement unit of $tcpO_2$, $tcpCO_2$ can be set on the TCM4 or TCM5 FLEX. When the measurement unit is changed, the tabular trend data of $tcpO_2$ and $tcpCO_2$ on the bedside monitor will be deleted.

Alarm Limit Range for Each Parameter

	Adjustable Range	
Item	Lower Limit Upper Lin	nit [Auto] Setting *
	Adjustable Increments	
HR	20 bpm to 295 bpm 22 bpm to 300 b	pm Upper: current value +40 bpm
	25 bpm and below: 1 bpm increments 25 bpm and above: 5 bpm increments	Lower: current value –40 bpm
ST	-2.0 mV to +1.8 mV -1.8 mV to +2.0	mV
12-Lead ST	0.1 mV increments	Upper: current value +0.2 mV (+2 mm)
	-20 mm to +18 mm -18 mm to +20 m	nm Lower: current value -0.2 mV (-2 mm)
	1 mm increments	
Ext Tachy	- 22 bpm to 300 b	pm
	50 bpm and below: 1 bpm increments 50 bpm and above: 5 bpm increments	HR Lower Limit +10 bpm
Ext Brady	20 bpm to 295 bpm -	
	50 bpm and below: 1 bpm increments 50 bpm and above: 5 bpm increments	HR Lower Limit - 10 bpm
RR (Adult)	5 Bpm to 145 Bpm 10 Bpm to 150 B	3pm
	5 Bpm increments	Upper: current value +20 Bpm
RR (Child/Neonate)	0 Bpm to 148 Bpm 4 Bpm to 150 Bp	
	2 Bpm increments	
RR_SpO ₂ (Adult)	5 Bpm to 30 Bpm 10 Bpm to 35 Bp	om
	5 Bpm increments	
RR_SpO ₂ (Child)	6 Bpm to 32 Bpm 8 Bpm to 34 Bpr	N/A
	2 Bpm increments	
Apnea	- 10 sec. to 60 sec	
	1 sec. increments	15 sec.
BP1 to 4	0 mmHg to 295 mmHg 2 mmHg to 300 i	 mmHg
	0 mmHg to 50 mmHg: 2 mmHg increments 50 mmHg and above: 5 mmHg increments	When BP label is BP1/ART: Upper: current value +40 mmHg (+5.0
	0.0 kPa to 39.5 kPa 0.2 kPa to 40.0 k	kPa kPa)
	0 kPa to 7.0 kPa: 0.2 kPa increments 7.0 kPa and above: 0.5 kPa increments	Lower: current value -20 mmHg (-3.0 kPa) When BP label is other than BP1/ART: Upper: current value +20%
CVP	0.0 cmH ₂ O to 38 cmH ₂ O 2 cmH ₂ O to 40 c	cmH ₂ O Lower: current value -20%
	1 cmH ₂ O increments	
NIBP	10 mmHg to 295 mmHg 15 mmHg to 300) mmHg
	5 mmHg increments	Upper: current value +40 mmHg (+5.0
	1.5 kPa to 39.5 kPa 2.0 kPa to 40.0 k	kPa kPa) kPa Lower: current value -20 mmHg (-3.0 kPa)
	0.5 kPa increments	
SpO ₂	50%SpO ₂ to 99%SpO ₂ 51%SpO ₂ to 100	0%SpO ₂ Upper: OFF
	1%SpO ₂ increments	Lower: 90%SpO ₂
Ext SpO ₂	50%SpO ₂ to 98%SpO ₂ -	Upper: OFF
	1%SpO ₂ increments	Lower: 90%SpO ₂

The alarm limit can be set in the following range.

	Adjustable Range			
Item	Lower Limit	Upper Limit	[Auto] Setting *	
	Adjustab	le Increments		
EtCO ₂	1 mmHg to 98 mmHg	3 mmHg to 100 mmHg		
	1 mmHg increments		1	
	0.1 kPa to 13.1 kPa	0.3 kPa to 15.0 kPa	 Upper: current value +10 mmHg (+1.3 kPa / +1.3%) 	
	0.1 kPa increments		Lower: current value -10 mmHg (-1.3 kPa / -1.3%)	
	0.1% to 13.1%	0.3% to 15.0%	(-1.5 KF a / -1.5 /0)	
	0.1% increments			
InspCO ₂	-	1 mmHg to 4 mmHg		
	1 mmHg increments			
	-	0.1 kPa to 3.0 kPa	2 mm l m (0.2 kDz (0.20())	
	0.1 kPa increments		– 3 mmHg (0.3 kPa / 0.3%)	
	-	0.1% to 3.0%		
	0.1% increments			
TEMP	30.0°C to 44.0°C	31.0°C to 45.0°C		
	0.5°C increments		Upper: current value +2.0°C/+4.0°F	
	86.0°F to 111.0°F	88.0°F to 113.0°F	Lower: current value -2.0°C/-4.0°F	
	1.0°F increments			
SpCO	-	1%SpCO to 40%SpCO	_ N/A	
	1%SpCO increments			
SpMet	-	1%SpMet to 15%SpMet	– N/A	
	1%SpMet increments			
SpHb	1.0 g/dL to 24.0 g/dL	2.0 g/dL to 24.5 g/dL	– N/A	
	0.1 g/dL increments			
BIS	1 to 98	2 to 99	– N/A	
	increments of 1			

*: If the value exceeds the adjustable range, the limit within the range will be set.

The automatic setup will not be performed for the turned OFF limit.

About the SpO₂ Clinical Test

Covidien Unit

The SpO_2 and pulse rate measurement accuracy have been validated for each range by testing on healthy adult male and female volunteers against a laboratory CO-Oximeter.

The SpO₂ accuracy has been validated for the range from 70% to 100% by testing on healthy adult male and female volunteers (age: 19 to 48 years old) with light to dark skin pigmentation. The standard deviation is $\pm 2\%$ which encompasses 68% of the population.

The pulse rate accuracy has been validated for the range from 70% to 100% by testing on healthy adult male and female volunteers (age: 19 to 48 years old) with light to dark skin pigmentation The standard deviation is ± 3 bpm which encompasses 68% of the population.

These clinical test data are disclosed based on the data provided from Covidien.

Masimo Unit

The SpO₂, SpCO, SpMet, and pulse rate measurement accuracy have been validated for each range by testing on healthy adult male and female volunteers against a laboratory CO-Oximeter.

 SpO_2 and SpMet accuracy have been validated by testing on 16 neonatal NICU patients ranging in age from 7 days to 135 days old and weighing between 0.5 kg to 4.25 kg. Seventy-nine (79) data samples were collected over a range of 70% to 100% SpO_2 and 0.5% to 2.5% HbMet with a resultant accuracy of 2.9% SpO_2 and 0.9% SpMet.

The SpO₂ accuracy has been validated for the range from 70% to 100% by testing on healthy adult male and female volunteers (age: 21 to 36 years old) with light to dark skin pigmentation. Without body motion, the standard deviation is $\pm 2\%$ which encompasses 68% of the population. With body motion, the standard deviation is $\pm 3\%$ which encompasses 68% of the population, frictional or contact motion of 1 cm to 2 cm, and random vibration of 1 Hz to 5 Hz were tested.

The pulse rate accuracy has been validated for the range from 70% to 100% by testing on healthy adult male and female volunteers (age: 24 to 37 years old) with light to dark skin pigmentation The standard deviation is ± 3 bpm which encompasses 68% of the population.

The SpCO accuracy has been validated for the range from 0% to 40% by testing on healthy adult male and female volunteers with light to dark skin pigmentation. The standard deviation is $\pm 3\%$ which encompasses 68% of the population.

The SpMet accuracy has been validated for the range from 0% to 15% by testing on healthy adult male and female volunteers with light to dark skin pigmentation. The standard deviation is $\pm 1\%$ which encompasses 68% of the population.

The SpHb accuracy has been validated for the range from 8 g/dL to 17 g/dL by testing on healthy adult male and female volunteers with light to dark skin pigmentation. The standard deviation is ± 1 g/dL which encompasses 68% of the population.

These clinical test data are disclosed based on the data provided from Masimo.

Index

Numerics

12 Lead	
Printing 9-4	, 9-10
12-Lead Analysis	8-28
12-Lead Display	10-1
12-Lead Print Key	9-10
1ch BP Relay Cable	7-20
2ch BP Conversion Cable	7-20

А

A	
Airway Adapter	7-72, 7-74
Airway Adapter Calibration	11-47
Alarm during NIBP	7-45
Alarm Factor	
Alarm History	8-21
Alarm Indicator	2-7
Alarm Printing	9-6
Alarm Silence	6-8, 6-9
Alarm Silence Key	2-7, 3-1
Alarm Sound Suspend	6-8
Alarm Threshold Limit	6-11
Analog Output Connector	
ARRHY OFF	6-1
Arrhythmia	
Alarm	3-17
Analysis	5-3
Detection	3-17
Learn	3-9
Arrhythmia Alarm	
Asystole	8-3
Bigeminy	8-3
Brady	
Couplet	
Frequent	
Pause	8-3
Run	8-3
Slow VT	8-3
Tachy	8-3
Trigeminy	8-3
VF	8-3
VT	8-3
Arrhythmia Alarm Type	8-3
Arrhythmia Analysis	
Arrhythmia Detection Level	
AUX Connection Cable	
AUX Connector	2-8, 2-10

В

Battery Charging LED	2-3, 2-7, 4-1, 4-2
Battery Cover	2-7
Blood Temperature	8-45
BP Zero Balance Indicator	2-9, 2-10
BP Zero Balance Key	2-9, 2-10

С

8	
Calibration	
Airway Adapter	
Print Calibration	9-4, 9-10, 12-19
Capnostat 5 CO2 Sensor	
Card Format	8-36
Cardiac Index (CI)	8-45
Cardiac Output (CO)	8-45
CVA Detect	

D

Daily Check	4-1
Data Transfer Function	
Drug Calculation	8-51
Drug Calculation Formula	8-52
DS I/F Connector	2-6, 2-8
Dyna Alert	7-57
Dyna Alert Function	7-57
Dyna Alert Message	

Е

2-8
3-7
8-2
4-7
1-1
1-1

F

12-13
7-12, 12-6
, 7-11, 12-6
7-32, 12-8
7-32, 12-8
7-77, 7-83
2-7
9-10

G

Graphic Trend 8-6

Н

HR Delay	7-11
HR Lower Limit for Run	
HR Lower Limit for SVT	8-5
HR Lower Limit for VT	8-5
HR Lower Limit for VT/RUN	6-2

I

7-107
2-9, 2-10, 2-11
8-45

	L	
Lead-Off Message		. 3-17
List of Setup Items		
Basic Setup		12-15

Μ

Magnetic Card5-4Menu Configuration1-5
Message
NIBP Meas. Error 11-15
Message List 11-1
module-LAN Connector 2-5
Monitor New Patient 4-6
Monitor Suspend 5-10
Resume 5-11
Monitor Suspend Label 5-9
Monitor Suspend Timer 5-9
Monitoring Suspend 5-9
Multiparameter 1-2
Multiparameter Connector 2-8, 2-9, 2-10, 7-20

Ν

NIBP Auto Mode Key	3-1
NIBP Connector	2-8
NIBP Inflation Mode	
NIBP Start/Stop Key	3-1
Numeric Data	
Alarm Limit Range	14-26
Numeric Data Box Area	3-2
Numeric Data Box Message	11-14

0

Operation Keys 3-	·1
-------------------	----

Ρ

2-12
5-4
2-9, 2-10, 2-11
4-1
2-3, 2-6, 2-9, 2-10, 2-11, 4-1
12-7

R

Recall Factor	8-16
Recall List	
Recall Waveform	8-16
Record Key	2-12
Recording Paper	4-12
Release Lever	2-5, 2-6, 2-9, 2-10, 2-11
Release Lock Button	
Remote Control Sensor	2-3
Result Status	
RR Synchronized Mark	

RR_SpO2 Setup	7-46
RR/APNEA Alarm Source	7-17

S

Scale
CO2 Scale 7-81
Printer Auto Scale
SD Card Slot 2-8
Sensitivity
BP Transducer 14-16
size/scale
SpCO,SpMet,SpHb,SpOC Measurement 7-41
Specification 14-1
SpO2 Connector
SpO2 Second Alarm 6-2
Standby Switch 2-3, 2-6, 3-1
STAT Mode
Status Message 8-45
Super Unit Connector 2-5, 2-6
Synchronize Hemodynamic Data with the Central Monitor 8-39

т

Tabular Trend	8-13
Temperature Connector	2-8
Thermodilution Curve	8-45
Time Bar	8-1
Time Scale	8-45
Trend Parameters	8-8
Trend Setup	8-8

U

U-LINK Connector	2-12
USB Memory Slot	2-7
User Key	3-2
User Selectable Numeric Data	

V

Ventilator Alarm Message 11-20

W

Waveform Size
Auto
ECG 7-6, 11-22, 11-24
ECG Waveform Size Selection 7-6, 12-6, 14-12
Pulse Wave7-43
RESP 7-15, 11-27, 12-13
Respiration Waveform Size Selection 12-6
SpO2 Waveform Size Selection 12-7
ST 12-14
Waveform Size/Scale 3-20

Х

X-Cal	′- 42
-------	--------------

	Υ	
Y-Piece		
	Z	
Zero Calibration		
Zeroing		11-8, 11-12

SD Logo is a trademark of SD-3C, LLC.

Microstream, FilterLine and CapnoLine are trademarks or registered trademarks of Oridion Medical 1987 Ltd. X-Cal is a trademark of Masimo Corporation.

Other company and product names used in this manual are trademarks or registered trademarks of respective companies.

FUKUDA DENSHI CO., LTD.

3-39-4 Hongo, Bunkyo-ku, Tokyo, Japan Tel: +81-3-5684-1455 Fax: +81-3-3814-1222 http://www.fukuda.com

Printed in Japan 4L011439B 201808