





# DS-8100

## Mounting Solution Guide

# DS-8100 with No Additional Units Configuration and Mounting Parts



**DS-8100N/M**  
Main Unit



**FD-0026-16**  
5" Plate Adapter  
(attach to bottom of DS-8100)

## Roll Stand Option

## Wall Mount Option



**FD-0013-61E**  
Roll Stand



**WMM-0006-04**  
Dual Utility Hook  
*\*do not need one for the roll stand*



**FD-0009-21D**  
M Series Fixed Height Arms

**OR**

**OR**

**OR**



**WS-0003-26A**  
Quad Utility Hook




**FD-0021-13**  
Variable Height




**FD-0033-01**  
Wall Mount  
Variable Height

# DS-8100 with EtCO2 Module

## Configuration and Mounting Parts




**DS-8100N/M**  
Main Unit




**CJO-15RR0.17G**  
AUX Connection Cable (0.17M)  
\*Use with OAO-86A HCP Fixing Bracket for DS-8100



**HCP-810**  
Sidestream Unit



**FD-0026-16**  
5" Plate Adapter  
(attach to bottom of DS-8100)



**OAO-86A**  
Fixing Bracket for HCP/HPD



**HPD-810**  
Mainstream Unit

### Roll Stand Option

### Wall Mount Option



**FD-0013-61E**  
Roll Stand



**WMM-0006-04**  
Dual Utility Hook  
\*do not need one for the roll stand



**FD-0009-21D**  
M Series Fixed Height Arms

**OR**

**OR**


**OR**



**WS-0003-26A**  
Quad Utility Hook




**FD-0021-13**  
Variable Height



**FD-0033-01**  
Wall Mount  
Variable Height

# DS-8100 with Recoder

## Configuration and Mounting Parts



**DS-8100N/M**  
Main Unit



**OAO-66A**  
DS-8100 Adapter



**FD-0020-16**  
Base Plate Adapter  
for OAO-66A



**HR-811**  
Printer / Expansion Port Unit  
Slave Monitor Output /  
Analog Output

### Roll Stand Option

### Wall Mount Option



**FD-0013-61E**  
Roll Stand



**WMM-0006-04**  
Dual Utility Hook  
\*do not need one for the roll stand



**FD-0009-21D**  
M Series Fixed Height Arms

**OR**

**OR**

**OR**



**WS-0003-26A**  
Quad Utility Hook




**FD-0021-13**  
Variable Height



**FD-0033-01**  
Wall Mount  
Variable Height



## Bed Mount for DS-8100



**DS-8100N/M**  
Main Unit



**OAO-85A**  
Bed Mount Attachment



\*compatible with roll stand and wall mount  
\*can not be used with the HR-811 printer attached



FUKUDA DENSHI USA, INC.  
17725-C NE 65th Street Redmond, WA 98052  
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (425) 869-2018  
www.fukudaamerica.com

GCX® is a registered trademark of GCX Corporaon.  
The information contained in this document is subject to change  
without notice. © 2022 FUKUDA DENSHI USA, INC. All rights reserved.