

DS-8100N/M General Specifications

Display		
Size	10.2 inch wide color LCD	
Resolution	1024 x 600, WSVGA	
Waveforms	Max. 14	
Sweep Speed	Circulatory	6.25, 12.5, 25, 50 mm/s
	Respiratory	6.25, 12.5, 25 mm/s
Waveform Display	Stationary Trace Mode	

Operation		
Touch Screen Method		
Jog Dial with Push Key		
5 Fixed Keys	Home, Menu, Previous Display, NIBP (Start/Stop) and Alarm Silence	

Parameters		
Waveform	ECG, IBP (max. 2ch), SpO ₂ , RESP, CO ₂	
Measurements*	HR, ST, VPC, IBP (max. 2ch), SpO ₂ , PI (Masimo®), SpCO/SpMet/SpHb/PVI (Masimo®, optional), RR, PR, APNEA, NIBP, TEMP (max. 4ch), CO ₂ , CO	
Arrhythmia Analysis (27)	ASYSTOLE, VF, VT, SVT, RUN, Tachy, Brady, Bigeminy, Frequent, Couplet, Trigeminy, PAUSE, ST, Prolonged RR, Ext. Tachy, Ext. Brady, S Frequent, S Couplet, SVPC, Irregular RR, Triplet, R on T, Multiform, Vent Rhythm, VPC, Pacer Not Capture, Pacer Not Pacing	


*Masimo® measurement is available only for Masimo® model devices.

ECG		
Measuring Range	Adult/Child	0, 12 to 300 bpm
	Neonate	0, 30 to 300 bpm
Measuring Accuracy	±3 bpm	
Size	1/4, 1/2, 1, 2 and 4	
HR Display Response Time	Adult/Child	6 sec
	Neonate	3 sec

Defibrillation Proof Provided

SpO ₂		
Measuring Method	2 Wavelength Pulse Wave	
Module	Nellcor™/Masimo® Technology	
Measuring Range	1 ~ 100%	
Measuring Accuracy	±3% (Nellcor™)/ ±2% (Masimo®)	
PR Measurement Range	20 ~ 250 bpm (Nellcor™)	
	26 ~ 239 bpm (Masimo®)	
PR Accuracy	±3 bpm	

Respiration		
Measuring Method	Impedance	
Measuring Range	0, 4 to 150 bpm	
Measuring Accuracy	±3 bpm	



Temperature		
Measuring Method	Thermistor	
Measuring Range	32°F to 113°F / 0 to 45°C	
Measuring Accuracy	±0.4°F (between 77°F and 113°F) / ±0.2°C (between 25°C and 45°C)	
Number of Channels	Max. 4	

EtCO ₂ (optional)		
Sidestream (Microstream™)		
Measuring Range	0 to 99 mmHg	
Measuring Accuracy	0 to 38 mmHg : ±2 mmHg	
	39 to 99 mmHg : ±0.05 reading + 0.08x (reading -39 mmHg)	

Mainstream (RESPIRONICS®)		
Measuring Range	0 to 150 mmHg	
Measuring Accuracy	0 to 40 mmHg : ±2 mmHg	
	41 to 70 mmHg : ±5%	
	71 to 100 mmHg : ±8%	
	101 to 150 mmHg : ±10%	

Invasive Blood Pressure		
Measuring Range	-50 to 300 mmHg	
Measuring Accuracy	±1 mmHg or within ±2%	
PR Measuring Range	Adult	12 to 300 bpm
	Neonate	30 to 300 bpm
PR Accuracy	±1 bpm or ±3%	
Number of Channels	Max. 2	

Non-Invasive Blood Pressure		
Measuring Method	Oscillometric	
Measuring Range	Adult	10 to 280 bpm
	Child	10 to 180 bpm
	Neonate	10 to 130 bpm
Static Pressure Accuracy	±3 mmHg	
PR Measuring Range	40 to 240 bpm	
PR Accuracy	±2% or ±2 bpm	

Printer (HR-811, optional)	
Number of Waveform	Max. 3
Printing Type	Thermal
Printing Speed	50, 25 mm/sec
Waveform Printed	ECG, RESP, SpO2, IBP and CO2

Full Disclosure

When optional CF card (16 GB) is used, 48 hours with a maximum of 6 waveforms

Graphic Trend

24 hours / 4 selectable groups / 46 items

Tabular Trend

24 hours / 6 selectable groups / 91 items

Dimensions (W x H x D) (not including the protrusion)	
Main Unit	11.8 x 10.4 x 3.0 inch 300 x 265 x 75 mm

Weight	
Main Unit	7.7 lbs 3.5 kg

Power	
Power Requirements	AC 100 ~ 240V, 50/60 Hz
Power Consumption	60 VA max.
Battery Operation	3 hours with optional battery BTO-008

Environmental Condition

Operating Environment	
Ambient Temperature	50 to 104°F / 10 to 40°C
Relative Humidity	30 to 85%

Transport / Storage	
Ambient Temperature	14 to 140°F / -10 to 60°C
Relative Humidity	10 to 95%

Safety

General Standard	IEC 60601-1 : 1988 + A1 : 1991 + A2 : 1995 IEC 60601-1-1 : 2000
EMC Standard	IEC 60601-1-2 : 2007
Electrical Shock Protection	Class I (AC Power Operation)
	Internal Powered Equipment (Battery Operation)
Conformity	CE Marking per 93/42/EEC Directive

RoHS Compliant

External Connection

Serial Connector (COM x 2 ports)
CF / SD Card Slots
DS-LAN Connector
AUX Connector
Status I/O Connector (Status II x 2 ports)
Telemetry Connector
Analog Output Connector (via HR-811 Printer)
VGA Output Connector (via HR-811 Printer)



FUKUDA DENSHI USA, INC.

17725-C NE 65th Street Redmond, WA 98052

Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (888) 224-7090

www.fukudaamerica.com

Nellcor™ and Microstream™ are registered trademarks of Covidien. Masimo® is a registered trademark of Masimo Corporation.

RESPIRONICS® is a registered trademark of Koninklijke Philips.

The information contained in this document is subject to change without notice.

© 2022 FUKUDA DENSHI USA, INC. All rights reserved.

Document No. USRB-0121-00