




# DS-8200


## Mounting Solution Guide

# DS-8200 with No Additional Units

## Configuration and Mounting Parts



**DS-8200**  
System



**OAO-70A**  
GCX® Attachment  
for Monitor (BS-8210)

### Roll Stand Option

### Wall Mount Option



**FD-0013-61E**  
Roll Stand



**WMM-0006-04**  
Dual Utility Hook  
*\*do not need one for the roll stand*

**OR**

**OR**

**OR**



**FD-0009-21D**  
M Series Fixed Height Arms



**WS-0003-26A**  
Quad Utility Hook



**FD-0021-13**  
Variable Height




**FD-0033-01**  
Wall Mount  
Variable Height

# DS-8200 with EtCO2 Module

## Configuration and Mounting Parts



**DS-8200**  
System



**CJO-15RR0.65G**  
AUX Connection Cable (0.65M)




**HCP-810**  
Sidestream Unit


**OR**



**HPD-810**  
Mainstream Unit



**OAO-70A**  
GCX® Attachment  
for Monitor (BS-8210)



**OAO-72A**  
Gas Unit / External Output Mounting Bracket

### Roll Stand Option



**FD-0013-61E**  
Roll Stand

### Wall Mount Option



**WMM-0006-04**  
Dual Utility Hook  
*\*do not need one for the roll stand*

**OR**



**FD-0009-21D**  
M Series Fixed Height Arms

**OR**



**WS-0003-26A**  
Quad Utility Hook

**OR**



**FD-0021-13**  
Variable Height



**FD-0033-01**  
Wall Mount  
Variable Height

# Fukuda Mounting Kit for DS-8200 with Recorder


## Configuration and Mounting Parts



**DS-8200**  
System




**OAO-70A**  
GCX® Attachment  
for Monitor (BS-8210)



**CJO-09SS1.5**  
Unit Connection Cable  
(U-LINK) (1.5M)



**HR-800**  
3CH Thermal Printer



**OAO-68A**  
Countertop for DS-8200



**OAO-74A**  
Recorder Mounting Bracket  
for OAO-68A

### Roll Stand Option



**FD-0013-61E**  
Roll Stand

### Wall Mount Option



**WMM-0006-04**  
Dual Utility Hook  
\*do not need one for the roll stand



**FD-0009-21D**  
M Series Fixed Height Arms

**OR**

**OR**

**OR**



**WS-0003-26A**  
Quad Utility Hook



**FD-0021-13**  
Variable Height




**FD-0033-01**  
Wall Mount  
Variable Height

# GCX® Mounting Kit for DS-8200 with Recorder

## Configuration and Mounting Parts




**DS-8200**  
System




**CJO-09SS1.5**  
Unit Connection Cable  
(U-LINK) (1.5M)



**HR-800**  
3CH Thermal Printer



**OAO-70A**  
GCX® Attachment  
for Monitor (BS-8210)



**WMM-0006-01C**  
Down Post (6") for  
M/VHM Mount

\*Single Stage Fixing Post for HR-800  
Recorder Only



**FD-0028-81**  
Mounting Tray  
for HR-800



**FD-0027-83**  
Single-Sided  
Downpost-tray  
Adapter

## Wall Mount Option



**FD-0009-21D**  
M Series Fixed  
Height Arms

**OR**



**FD-0021-13**  
Variable Height

**OR**



**FD-0033-01**  
Wall Mount  
Variable Height




## DS-8200 with Telemetry Transmitter



**DS-8200**  
System



**HLX-801(600MHz)**  
**HLX-801(G) (1.4GHz)**  
Telemetry Transmitter Module



**OAT-02A**  
Installation Cover

### Bed Mount Kit for DS-8200



**DS-8200**  
System



**OAO-69A**  
Bed Mount Attachment



\*compatible with roll stand and wall mount

### VESA Attachment for LC-8210



**LC-8210**  
Display Unit



**CJO-16SS1.5G (1.5M)**  
**CJO-16SS3G (3.0M)**  
**CJO-16SS5G (5.0M)**  
Display Unit Extension Cable



**OAO-71A**  
VESA Attachment  
for LC-8210



FUKUDA DENSHI USA, INC.  
17725-C NE 65th Street Redmond, WA 98052  
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (425) 869-2018  
www.fukudaamerica.com

GCX® is a registered trademark of GCX Corporation.  
The information contained in this document is subject to change without notice. © 2022 FUKUDA DENSHI USA, INC. All rights reserved.