

DS-8000 Series Patient Monitor Accessory List





How to order

If you would like to order accessories email/call Customer Support at **(800) 365-6668** or **CustomerService@fukudaamerica.com** and list the following:



Customer number



Ship to address




P.O. number

x2

Quantity

#

Item number



DS-8000 Series Accessory List

TABLE OF CONTENTS

NEW

DS-8007 Transport Monitor	2	DS-8200 Transport Monitor	23
Configuration	3	Configuration	24
ECG	4	ECG	25
NIBP	4	NIBP	25
NIBP	5	NIBP	26
IBP	5	IBP	26
SpO2	6	SpO2	27
Temperature	6	Temperature	27
CO	7	CO	28
EtCO2	7	EtCO2	28
BIS™	8	Print	29
Print	8	Other	29
Other	8		
NEW		DS-8500 Bedside Monitor	30
DS-8400 Bedside Monitor	9	Configuration	31
Configuration	10	ECG	32
ECG	11	NIBP	33
NIBP	12	IBP	34
IBP	13	SpO2	34
SpO2	13	Temperature	35
Temperature	14	CO	35
CO	14	EtCO2	36
EtCO2	15	BIS™	37
BIS™	16	Print	38
Print	17	Agent	38
Agent	17	Other	38
Other	17		
DS-8100 Portable Bedside Monitor	18	DS-8900 Central Station	39
ECG	19	Print	40
NIBP	19	Other	40
NIBP	20		
IBP	20		
SpO2	21		
Temperature	21		
CO	21		
EtCO2	22		
Print	22		
Other	22		

DS-8007 TRANSPORT MONITOR





DS-8007N - Nellcor™ Model
DS-8007M - Masimo® Model

ECG
12-Lead

IBP
max. 4ch

NIBP

EtCO2

BIS™

SpO2

TEMP
max. 6ch

CO

Print

Optional:

Configuration

— Standalone —



DS-8007N - Nellcor™ Model
DS-8007M - Masimo® Model

*Required to have AC Power



DSA-81 AC Unit

— To Connect to DS-8400 —

Dock on the Back



DS-8007N - Nellcor™ Model
DS-8007M - Masimo® Model


OR



DS-8400
Bedside Monitor

Module Connection Cable

*Connect to Main Unit DSC-8410



○ **CJO-08SS1.5** (1.5M)

Special Order Item

- **CJO-08SS0.3** (0.3M)
- **CJO-08SS3.5** (3.5M)
- **CJO-08SS5** (5M)



DSA-82
DS-8007 Adapter


— To Connect to DS-8500 —



DS-8500
Bedside Monitor

Module Connection Cable

*Connect to Main Unit DSC-8530



○ **CJO-08SS1.5** (1.5M)

Special Order Item

- **CJO-08SS0.3** (0.3M)
- **CJO-08SS3.5** (3.5M)
- **CJO-08SS5** (5M)



DSA-82
DS-8007 Adapter


ECG



DS-8007
Main Unit




○ **CIO-07CTP-3NA** - 3 Lead ECG Trunk Cable



○ **CMO-07FT-3NAB** - 3 Lead ECG Patient Cable




○ **CIO-07CTP-4NA** - 4 Lead ECG Trunk Cable




○ **CMO-07FT-4NAB** - 4 Lead ECG Patient Cable



○ **CIO-07CTP-5NA** - 5 Lead ECG Trunk Cable



○ **CMO-07FT-5NAB** - 5 Lead ECG Patient Cable




○ **CMO-07FT-10NAB** - 12 Lead ECG Trunk & Patient Cable for 10-electrode

NIBP



DS-8007
Main Unit


Device with
Serial #
Starting with
5040xxxx




○ **OA-80APR3.5** - 3.5M Cuff Hose

Special Order Item

○ **OA-80APR1.5** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

NIBP



DS-8007
Main Unit


Device with
Serial #
Starting with
5041xxxx




- **OA-80APS3.5-S** - 3.5M Cuff Hose

Special Order Item

- **OA-80APS1.5-S** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5-S** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5-S** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

IBP



DS-8007
Main Unit

***For transducer cables other than the listed below, please contact our Customer Service**




- **CJO-P01B-DJ0.5** - 2ch BP Conversion Cable

**need if 2 IBP cables will be used*



- **CJO-P01B-SB3.6** - 3.6M DTX Plus


Special Order Item



- **CJO-P01B-DB0.8** - 0.8M 2ch DTX Plus
- **CJO-P01B-DB4.3** - 4.3M 2ch DTX Plus


SpO₂

DS-8007
Main Unit


DS-8007N - Nellcor™ Model

- **1M1120** - Nellcor™ SpO₂ Sensor Extender

***Can only order with Equipment Orders.
Call Nellcor/Medtronic for replacements.**



- **OA-MAX-AL** - Nellcor™ Disposable SpO₂



- **OA-DS-100A** - Nellcor™ Durasensor SpO₂

*Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors

DS-8007M - Masimo® Model


Please contact Masimo® for Accessories

TEMP

DS-8007
Main Unit

**Required for 3, 4, 5 or 6 x TEMP*

**NOT Required for 1 or 2 x TEMP*



- **CJO-P01T-DA0.5** - 0.5M 2ch Temperature Relay Cable

- **CJO-P01T-DA4.0** - 4.0M 2ch Temperature Relay Cable



- **CS409** - Temperature Probe Surface

CO

DS-8007
Main Unit○ **CJO-P01C-C2.4** - Catheter Relay Cable○ **CJO-P01C-T2.4** - Injectate Probe Relay Cable○ **CJO-P01C-L2.4** - In-line Sensor Relay Cable○ **CJO-P01C-F2.4** - Flow-through Sensor Relay CableEtCO₂○ **HCP-820** Sidestream Unit


*Contact your local Medtronic Representative for sampling tube

○ Covidien Microstream™

○ **HPD-820** Mainstream Unit


*Contact your local Philips Representative

○ PHILIPS RESPIRONICS®


DS-8007
Main Unit

*cannot be used with BIS™ module HBX-800 and OAO-103A

AUX LEM Connection Cable

○ **CJO-25TR0.65** (0.65M)

*Use with OAO-103A Mounting Bracket for DS-8007

○ **CJO-25TR1.5** (1.5M)○ **CJO-25TR2.7** (2.7M)○ **HCP-810** Sidestream Unit

*Contact your local Medtronic Representative for sampling tube

○ Covidien Microstream™


○ **HPD-810** Mainstream Unit

*Contact your local Philips Representative

○ PHILIPS RESPIRONICS®


BIS™

*cannot be used with EtCO2 module HCP-810/HPD-810 and OAO-103A

AUX LEM Connection Cable

○ **CJO-25TR0.65** (0.65M)

*Use with OAO-103A Mounting Bracket for DS-8007



- **CJO-25TR1.5** (1.5M)
- **CJO-25TR2.7** (2.7M)



- **HBX-800 BISx™** Module Interface



*shown connected to HBX-800 BISx™ Module Interface

- **186-0195-SF BISx™** System


DS-8007
Main Unit

Print

DS-8007
Main Unit



- **DSA-81** AC Unit

Unit Connection Cable (U-LINK)

- **CJO-09SS1.5** (1.5M)

Special Order Item

- **CJO-09SS0.3** (0.3M)
- **CJO-09SS5** (5.0M)



- **HR-800** 3CH Thermal Printer
use with paper OP050-01TDR

*Countertop Recorder Mount OAO-107A is a Special Order Item

Accessories

HLX-801 (600MHz)
HLX-801(G) (1.4GHz)
Telemetry Transmitter
Module

*use with DSA-81 AC Unit
and OAO-111A HLX Holder Mount



CF-820
Remote Controller



BTO-008
Rechargeable
Battery



SD-16G
16GB SD
Full Disclosure Card

*for 240 Hours of Full Disclosure




FSD-8GA
8GB FSD
Full Disclosure Card

*for 120 Hours of Full Disclosure

DS-8400 BEDSIDE MONITOR





DS-8400

ECG
12-Lead

IBP
max. 8ch

NIBP

EtCO2

BIS™

SpO2
max. 2ch

TEMP
max. 6ch

CO

Print

Agent

Optional:

Configuration

DS-8007 - Portable
Bedside Monitor

DSC-8410 -Main Unit




LC-8016TC - 15.6" Display Unit
LC-8018TC - 18.5" Display Unit

Dock on the Back



DS-8007N - Nellcor™ Model
DS-8007M - Masimo® Model

OR



HSA-81 - HS Adapter




HS-8312N - Nellcor™ HS Module
HS-8312M - Masimo® HS Module

Separate Connection

Module Connection Cable

*Connect to Main Unit DSC-8410



○ **CJO-08SS1.5** (1.5M)



DSA-82
DS-8007 Adapter



DS-8007N
Nellcor™ Model
DS-8007M
Masimo® Model



HSA-80
HS- Adapter



HS-8312N
Nellcor™ HS Module
HS-8312M
Masimo® HS Module

***Can fit one module in the DSC-8410 unit - requires input box for more than one**

Special Order Item

- **CJO-08SS0.3** (0.3M)
- **CJO-08SS3.5** (3.5M)
- **CJO-08SS5** (5M)



IB-8004
Input Box



HG-820
Nellcor™
SpO2 Module



HG-810
Masimo®
SpO2 Module



HM-800
Multi-Module



HP-800
Multi Port
Module



HM-801
Multi-Module

ECG




DS-8007
Main Unit



HS-8312
HS Module




○ **CIO-07CTP-3NA** - 3 Lead ECG Trunk Cable



○ **CMO-07FT-3NAB** - 3 Lead ECG Patient Cable




○ **CIO-07CTP-4NA** - 4 Lead ECG Trunk Cable




○ **CMO-07FT-4NAB** - 4 Lead ECG Patient Cable



○ **CIO-07CTP-5NA** - 5 Lead ECG Trunk Cable



○ **CMO-07FT-5NAB** - 5 Lead ECG Patient Cable



○ **CMO-07FT-10NAB** - 12 Lead ECG Trunk & Patient Cable for 10-electrode

NIBP




DS-8007
Main Unit



HS-8312
HS Module


Device with
Serial #
Starting with
5040xxxx




- **OA-80APR3.5** - 3.5M Cuff Hose

Special Order Item

- **OA-80APR1.5** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5** - 3.5M Neonate Cuff Hose

Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

NIBP




DS-8007
Main Unit



HS-8312
HS Module


Device with
Serial #
Starting with
5041xxxx



- **OA-80APS3.5-S** - 3.5M Cuff Hose


Special Order Item

- **OA-80APS1.5-S** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)

Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs



- **OA-80NE1.5-S** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5-S** - 3.5M Neonate Cuff Hose

IBP



DS-8007
Main Unit



HS-8312
HS Module




HM-800
Multi-Module



HM-801
Multi-Module

***For transducer cables other than the listed below, please contact our Customer Service**




- **CJO-P01B-DJ0.5** - 2ch BP Conversion Cable

**need if 2 IBP cables will be used*



- **CJO-P01B-SB3.6** - 3.6M DTX Plus

Special Order Item



- **CJO-P01B-DB0.8** - 0.8M 2ch DTX Plus
- **CJO-P01B-DB4.3** - 4.3M 2ch DTX Plus

SpO2



DS-8007
Main Unit




HG-820
Nellcor™ SpO2 Module




HG-810
Masimo® SpO2 Module

DS-8007N or HG-820 - Nellcor™ Model



- **1M1120** - Nellcor™ SpO2 Sensor Extender

***Can only order with Equipment Orders.
Call Nellcor/Medtronic for replacements.**



- **OA-MAX-AL** - Nellcor™ Disposable SpO2



- **OA-DS-100A** - Nellcor™ Durasensor SpO2

**Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors*

DS-8007M or HG-810 - Masimo® Model

Please contact Masimo® for Accessories

TEMP



DS-8007
Main Unit



HS-8312
HS Module



HM-800
Multi-Module




HM-801
Multi-Module

**Required to measure TEMP with HS-8312, HM-800 and HM-801*

**NOT Required for 1 or 2 x TEMP with DS-8007*

**Required for 3, 4, 5 or 6 x TEMP with DS-8007*



○ **CJO-P01T-DA0.5** - 0.5M 2ch Temperature Relay Cable

○ **CJO-P01T-DA4.0** - 4.0M 2ch Temperature Relay Cable



○ **CS409** - Temperature Probe Surface

CO



DS-8007
Main Unit



HS-8312
HS Module



HM-800
Multi-Module



HM-801
Multi-Module



○ **CJO-P01C-C2.4** - Catheter Relay Cable



○ **CJO-P01C-T2.4** - Injectate Probe Relay Cable



○ **CJO-P01C-L2.4** - In-line Sensor Relay Cable



○ **CJO-P01C-F2.4** - Flow-through Sensor Relay Cable

EtCO₂

HS-8312
HS Module

AUX Connection Cable




- **CJO-15RR0.17G** (0.17M)
- **CJO-15RR0.65G** (0.65M)
- **CJO-15RR1.5G** (1.5M)

Special Order Item

- **CJO-15RR3G** (3M)



- **HCP-810** Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™



- **HPD-810** Mainstream Unit



*Contact your local Philips Representative

- PHILIPS RESPIRONICS®



- **HCP-820** Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™




- **HPD-820** Mainstream Unit



*Contact your local Philips Representative

- PHILIPS RESPIRONICS®




DS-8007 Main Unit

AUX LEM Connection Cable



- **CJO-25TR0.65** (0.65M)


*Use with OAO-103A Mounting Bracket for DS-8007



- **CJO-25TR1.5** (1.5M)
- **CJO-25TR2.7** (2.7M)



- **HCP-810** Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™



- **HPD-810** Mainstream Unit



*Contact your local Philips Representative


- PHILIPS RESPIRONICS®

*cannot be used with BIS™ module HBX-800 and OAO-103A

BIS™

HM-801
Multi-Module

AUX LEM Connection Cable



- **CJO-25TR0.65** (0.65M)
- **CJO-25TR1.5** (1.5M)
- **CJO-25TR2.7** (2.7M)



○ **HBX-800 BISx™** Module
Interface




*shown connected to HBX-800 BISx™ Module Interface


○ **186-0195-SF BISx™** System

*cannot be used with EtCO2 module HCP-810/HPD-810 and OAO-103A

AUX LEM Connection Cable



- **CJO-25TR0.65** (0.65M)
- *Use with OAO-103A Mounting Bracket for DS-8007



- **CJO-25TR1.5** (1.5M)
- **CJO-25TR2.7** (2.7M)



○ **HBX-800 BISx™** Module
Interface




*shown connected to HBX-800 BISx™ Module Interface

○ **186-0195-SF BISx™** System



DS-8007 Main Unit

AUX Connection Cable



- **CJO-15RR0.17G** (0.17M)
- **CJO-15RR0.65G** (0.65M)
- **CJO-15RR1.5G** (1.5M)

Special Order Item

- **CJO-15RR3G** (3M)



○ **HBX-800 BISx™** Module
Interface



*shown connected to HBX-800 BISx™ Module Interface

○ **186-0195-SF BISx™** System



HS-8312
HS Module


Print



DS-8400

***The HR-800 can also be built in directly into the DSC-8410**

Unit Connection Cable (U-LINK)



○ CJO-09SS1.5 (1.5M)

Special Order Item

- CJO-09SS0.3 (0.3M)
- CJO-09SS5 (5.0M)




○ **HR-800** 3CH Thermal Printer
use with paper OP050-01TDR

Agent



DS-8400

Unit Connection Cable (U-LINK)



○ CJO-09SS1.5 (1.5M)

Special Order Item

- CJO-09SS0.3 (0.3M)
- CJO-09SS5 (5.0M)



○ **MGU-801P** AGO2 Gas Unit

*Contact your local Mindray Representative for Water Trap and Sampling Tube

Accessories



**HLX-801 (600MHz)
HLX-801(G) (1.4GHz)**
Telemetry Transmitter
Module



CF-820
Remote Controller



BTO-005
Rechargeable Battery



FCS-64G
64GB CFast
Full Disclosure Card
*for 240 Hours of Full Disclosure



**FCF-128
(1Z7333B)**
128MB CF Card
Flash Memory Card
*for Data Transfer

DS-8100 PORTABLE BEDSIDE MONITOR





DS-8100N - Nellcor™ Model

DS-8100M - Masimo® Model

ECG

IBP
max. 2ch

NIBP

SpO2

TEMP
max. 4ch

CO

Optional:

EtCO2

Print


ECG



DS-8100
Main Unit




CIO-07CTP-3NA - 3 Lead ECG Trunk Cable



○ **CMO-07FT-3NAB** - 3 Lead ECG Patient Cable




CIO-07CTP-4NA - 4 Lead ECG Trunk Cable



○ **CMO-07FT-4NAB** - 4 Lead ECG Patient Cable



CIO-07CTP-5NA - 5 Lead ECG Trunk Cable




○ **CMO-07FT-5NAB** - 5 Lead ECG Patient Cable

NIBP



DS-8100
Main Unit


**Device with
Serial #
Starting with
5040xxxx**




○ **OA-80APR3.5** - 3.5M Cuff Hose

Special Order Item

○ **OA-80APR1.5** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

NIBP



DS-8100
Main Unit


Device with
Serial #
Starting with
5041xxxx




- **OA-80APS3.5-S** - 3.5M Cuff Hose

Special Order Item

- **OA-80APS1.5-S** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5-S** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5-S** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

IBP



DS-8100
Main Unit

***For transducer cables other than the listed below, please contact our Customer Service**




- **CJO-P01B-DJ0.5** - 2ch BP Conversion Cable

**need if 2 IBP cables will be used*




- **CJO-P01B-SB3.6** - 3.6M DTX Plus

Special Order Item




- **CJO-P01B-DB0.8** - 0.8M 2ch DTX Plus
- **CJO-P01B-DB4.3** - 4.3M 2ch DTX Plus

SpO₂DS-8100
Main Unit**DS-8100N - Nellcor™ Model**

- **1M1120** - Nellcor™ SpO₂ Sensor Extender

***Can only order with Equipment Orders.
Call Nellcor/Medtronic for replacements.**



- **OA-MAX-AL** - Nellcor™ Disposable SpO₂



- **OA-DS-100A** - Nellcor™ Durasensor SpO₂

*Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors


DS-8100M - Masimo® Model

Please contact Masimo® for Accessories

TEMPDS-8100
Main Unit

**Required for 3 or 4 x TEMP*

**NOT Required for 1 or 2 x TEMP*



- **CJO-P01T-DA0.5** - 0.5M 2ch Temperature Relay Cable
- **CJO-P01T-DA4.0** - 4.0M 2ch Temperature Relay Cable



- **CS409** - Temperature Probe Surface

CODS-8100
Main Unit

- **CJO-P01C-C2.4** - Catheter Relay Cable



- **CJO-P01C-T2.4** - Injectate Probe Relay Cable



- **CJO-P01C-L2.4** - In-line Sensor Relay Cable




- **CJO-P01C-F2.4** - Flow-through Sensor Relay Cable

EtCO₂

DS-8100 Main Unit


AUX Connection Cable



○ **CJO-15RR0.17G** (0.17M)

*Use with OAO-86A HCP Fixing Bracket for DS-8100

OAO-86A



○ **CJO-15RR0.65G** (0.65M)


○ **CJO-15RR1.5G** (1.5M)

Special Order Item

○ **CJO-15RR3G** (3M)



○ **HCP-810** Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

○ Covidien Microstream™



○ **HPD-810** Mainstream Unit



*Contact your local Philips Representative

○ PHILIPS RESPIRONICS®

Print



DS-8100 Main Unit



○ **HR-811** Printer/ Expansion Port Unit
Slave Monitor Output/Analog Output
use with paper OP050-01TDR

Accessories



HLX-801 (600MHz)
HLX-801(G) (1.4GHz)
Telemetry Transmitter
Module



CF-820
Remote Controller



BTO-008
Rechargeable Battery



FCF-16GA
(9B0102420)
16GB Full Disclosure Card

*for 48 Hours of Full Disclosure



SD-8G
8GB SD Card
*to store patient data



SD-1G
1GB SD Card
*to store patient data



FCF-128
(1Z7333B)
128MB CF Card
Flash Memory Card
*for transferring monitor settings

DS-8200 TRANSPORT MONITOR



**DS-8200****ECG**

12-Lead

IBP

max. 6ch

NIBP**SpO2****TEMP**

max. 6ch

CO

Optional:

EtCO2**Print**

Configuration

**LC-8210**
Display Unit**BS-8210**
Base Unit**HS-8312N** - Nellcor™ HS Module
HS-8312M - Masimo® HS Module**HSB-80**
HS Adapter


ECG



DS-8200
(HS-8312 HS Module)




○ **CIO-07CTP-3NA** - 3 Lead ECG Trunk Cable



○ **CMO-07FT-3NAB** - 3 Lead ECG Patient Cable




○ **CIO-07CTP-4NA** - 4 Lead ECG Trunk Cable




○ **CMO-07FT-4NAB** - 4 Lead ECG Patient Cable



○ **CIO-07CTP-5NA** - 5 Lead ECG Trunk Cable



○ **CMO-07FT-5NAB** - 5 Lead ECG Patient Cable




○ **CMO-07FT-10NAB** - 12 Lead ECG Trunk & Patient Cable for 10-electrode

NIBP



DS-8200
(HS-8312 HS Module)


Device with
Serial #
Starting with
5040xxxx




○ **OA-80APR3.5** - 3.5M Cuff Hose

Special Order Item

○ **OA-80APR1.5** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

NIBP



DS-8200
(HS-8312 HS Module)


Device with
Serial #
Starting with
5041xxxx




- **OA-80APS3.5-S** - 3.5M Cuff Hose

Special Order Item

- **OA-80APS1.5-S** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5-S** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5-S** - 3.5M Neonate Cuff Hose


Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

IBP



DS-8200
(HS-8312 HS Module)

***For transducer cables other than the listed below, please contact our Customer Service**




- **CJO-P01B-DJ0.5** - 2ch BP Conversion Cable

**need if 2 IBP cables will be used*



- **CJO-P01B-SB3.6** - 3.6M DTX Plus


Special Order Item



- **CJO-P01B-DB0.8** - 0.8M 2ch DTX Plus
- **CJO-P01B-DB4.3** - 4.3M 2ch DTX Plus


SpO₂

DS-8200
(HS-8312 HS Module)


HS-8312N - Nellcor™ Model

- **1M1120** - Nellcor™ SpO₂ Sensor Extender

***Can only order with Equipment Orders.
Call Nellcor/Medtronic for replacements.**



- **OA-MAX-AL** - Nellcor™ Disposable SpO₂



- **OA-DS-100A** - Nellcor™ Durasensor SpO₂


*Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors

HS-8312M - Masimo® Model

Please contact Masimo® for Accessories

TEMP

DS-8200
(HS-8312 HS Module)



- **CJO-P01T-DA0.5** - 0.5M 2ch Temperature Relay Cable
- **CJO-P01T-DA4.0** - 4.0M 2ch Temperature Relay Cable



- **CS409** - Temperature Probe Surface

CO



DS-8200
(HS-8312 HS Module)



○ **CJO-P01C-C2.4** - Catheter Relay Cable



○ **CJO-P01C-T2.4** - Injectate Probe Relay Cable



○ **CJO-P01C-L2.4** - In-line Sensor Relay Cable




○ **CJO-P01C-F2.4** - Flow-through Sensor Relay Cable

EtCO₂

DS-8200
(HS-8312 HS Module)

AUX Connection Cable



○ **CJO-15RR0.65G** (0.65M)


○ **CJO-15RR1.5G** (1.5M)

Special Order Item

○ **CJO-15RR3G** (3M)



○ **HCP-810** Sidestream Unit



*Contact your local Medtronic Representative
for sampling tube

○ **Covidien Microstream™**




○ **HPD-810** Mainstream Unit



*Contact your local Philips Representative

○ **PHILIPS RESPIRONICS®**



DAO-72A - Gas Unit/External Output
Mounting Bracket


*2 brackets per package

Print



DS-8200

Unit Connection Cable (U-LINK)



○ **CJO-09SS1.5** (1.5M)

Special Order Item

- **CJO-09SS0.3** (0.3M)
- **CJO-09SS5** (5.0M)




○ **HR-800** 3CH Thermal Printer
use with paper OP050-01TDR

Accessories



**HLX-801 (600MHz)
HLX-801(G) (1.4GHz)**
Telemetry Transmitter
Module



OAT-02A
Installation Cover



CF-820
Remote Controller



BTO-008
Rechargeable Battery

*MAX 2 x BTO-008 supported for 5 hour battery power operation
- 1 x BTO-008 for HSB-80 and 1 x BTO-008 for BS-8210



**FCF-16GA
(9B0102420)**
16GB Full Disclosure Card
*for 48 Hours of Full Disclosure



**FCF-128
(1Z7333B)**
128MB CF Card
Flash Memory Card
*for transferring monitor settings



SD-8G
8GB SD Card
*to store patient data



SD-1G
1GB SD Card
*to store patient data

DS-8500 BEDSIDE MONITOR





DS-8500

Configuration

ECG
12-Lead

IBP
max. 8ch

NIBP

SpO2
max. 2ch

TEMP
max. 8ch

CO

Optional:

EtCO2

BIS™


Print

Agent



LC-8015TC - 15" Display Unit
LC-8019TC - 19" Display Unit

DSC-8530 - Main Unit




OAO-44A - Stand
(optional)



Module Connection Cable

*Connect to Main Unit DSC-8530



○ **CJO-08SS1.5** (1.5M)



DSA-82
DS-8007 Adapter



DS-8007N
Nellcor™ Model

DS-8007M
Masimo® Model



HSA-80
HS- Adapter



HS-8312N
Nellcor™ HS Module

HS-8312M
Masimo® HS Module

Special Order Item

- **CJO-08SS0.3** (0.3M)
- **CJO-08SS3.5** (3.5M)
- **CJO-08SS5** (5M)



IB-8004
Input Box



HG-820
Nellcor™
SpO2 Module



HG-810
Masimo®
SpO2 Module



HM-800
Multi-Module



HP-800
Multi Port
Module



HM-801
Multi-Module

ECG




DS-8007
Main Unit



HS-8312
HS Module




○ **CIO-07CTP-3NA** - 3 Lead ECG Trunk Cable



○ **CMO-07FT-3NAB** - 3 Lead ECG Patient Cable




○ **CIO-07CTP-4NA** - 4 Lead ECG Trunk Cable




○ **CMO-07FT-4NAB** - 4 Lead ECG Patient Cable



○ **CIO-07CTP-5NA** - 5 Lead ECG Trunk Cable



○ **CMO-07FT-5NAB** - 5 Lead ECG Patient Cable



○ **CMO-07FT-10NAB** - 12 Lead ECG Trunk & Patient Cable for 10-electrode


NIBP

DS-8007
Main Unit



HS-8312
HS Module


Device with
Serial #
Starting with
5040xxxx




- **OA-80APR3.5** - 3.5M Cuff Hose

Special Order Item

- **OA-80APR1.5** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)



- **OA-80NE1.5** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5** - 3.5M Neonate Cuff Hose

Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs


NIBP

DS-8007
Main Unit



HS-8312
HS Module


Device with
Serial #
Starting with
5041xxxx




- **OA-80APS3.5-S** - 3.5M Cuff Hose

Special Order Item

- **OA-80APS1.5-S** - 1.5M Cuff Hose



- **98-0600-40** - Infant Cuff (8-13cm)
- **98-0600-41** - Pediatric Cuff (12-19cm)
- **CUF-8503** - Adult Cuff (Small) (17-25cm)
- **CUF-8504** - Adult Cuff (Medium) (23-33cm)
- **CUF-8505** - Adult Cuff (Large) (31-40cm)
- **CUF-8506** - Adult Cuff (Thigh) (38-50cm)
- **CUF-KIT-ADULT** - Adult Cuff Kit (Small, Medium, Large)
- **CUF-KIT-FULL** - Full Cuff Kit (Infant, Pediatric, Small, Medium, Large, Thigh)




- **OA-80NE1.5-S** - 1.5M Neonate Cuff Hose
- **OA-80NE3.5-S** - 3.5M Neonate Cuff Hose

Use with Trimline™ Neonatal Tempa-Kuff® BP Cuffs

IBPDS-8007
Main UnitHS-8312
HS ModuleHM-800
Multi-ModuleHM-801
Multi-Module

***For transducer cables other than the listed below, please contact our Customer Service**




○ **CJO-P01B-DJ0.5** - 2ch BP Conversion Cable

**need if 2 IBP cables will be used*



○ **CJO-P01B-SB3.6** - 3.6M DTX Plus

Special Order Item




○ **CJO-P01B-DB0.8** - 0.8M 2ch DTX Plus

○ **CJO-P01B-DB4.3** - 4.3M 2ch DTX Plus


SpO2DS-8007
Main UnitHS-8312
HS ModuleHG-820
Nellcor™ SpO2 ModuleHG-810
Masimo® SpO2 Module

DS-8007N, HS-8312N or HG-820 - Nellcor™ Model




○ **1M1120** - Nellcor™ SpO2 Sensor Extender

***Can only order with Equipment Orders.
Call Nellcor/Medtronic for replacements.**



○ **OA-MAX-AL** - Nellcor™ Disposable SpO2



○ **OA-DS-100A** - Nellcor™ Durasensor SpO2

**Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors*

DS-8007M, HS-8312M or HG-810 - Masimo® Model

Please contact Masimo® for Accessories

TEMP



DS-8007
Main Unit



HS-8312
HS Module



HM-800
Multi-Module




HM-801
Multi-Module

**Required for 3, 4, 5 or 6 x TEMP*

**Required to measure TEMP with HS-8312 or HM-800*

**NOT Required to measure 1 or 2 X TEMP with DS-8007*



○ **CJO-P01T-DA0.5** - 0.5M 2ch Temperature Relay Cable

○ **CJO-P01T-DA4.0** - 4.0M 2ch Temperature Relay Cable



○ **CS409** - Temperature Probe Surface

CO



DS-8007
Main Unit



HS-8312
HS Module



HM-800
Multi-Module



HM-801
Multi-Module




○ **CJO-P01C-C2.4** - Catheter Relay Cable



○ **CJO-P01C-T2.4** - Injectate Probe Relay Cable



○ **CJO-P01C-L2.4** - In-line Sensor Relay Cable




○ **CJO-P01C-F2.4** - Flow-through Sensor Relay Cable

EtCO₂

HS-8312
HS Module

AUX Connection Cable



- CJO-15RR0.65G (0.65M)
- CJO-15RR1.5G (1.5M)

Special Order Item

- CJO-15RR3G (3M)



- HCP-810 Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™



- HPD-810 Mainstream Unit



*Contact your local Philips Representative

- PHILIPS RESPIRONICS®



- HCP-820 Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™




- HPD-820 Mainstream Unit



*Contact your local Philips Representative

- PHILIPS RESPIRONICS®



DS-8007 Main Unit


*cannot be used with BIST™ module HBX-800 and OAO-103A

AUX LEM Connection Cable



- CJO-25TR0.65 (0.65M)

*Use with OAO-103A Mounting Bracket for DS-8007



- CJO-25TR1.5 (1.5M)
- CJO-25TR2.7 (2.7M)



- HCP-810 Sidestream Unit



*Contact your local Medtronic Representative for sampling tube

- Covidien Microstream™







- HPD-810 Mainstream Unit



*Contact your local Philips Representative

- PHILIPS RESPIRONICS®


BIS™  HM-801 Multi-Module	AUX LEM Connection Cable  ○ CJO-25TR0.65 (0.65M) ○ CJO-25TR1.5 (1.5M) ○ CJO-25TR2.7 (2.7M)	 ○ HBX-800 BISx™ Module Interface	 *shown connected to HBX-800 BISx™ Module Interface ○ 186-0195-SF BISx™ System
 DS-8007 Main Unit	AUX LEM Connection Cable  ○ CJO-25TR0.65 (0.65M) *Use with OAO-103A Mounting Bracket for DS-8007  OAO-103A	 ○ HBX-800 BISx™ Module Interface	 *shown connected to HBX-800 BISx™ Module Interface ○ 186-0195-SF BISx™ System
 HS-8312 HS Module	AUX Connection Cable  ○ CJO-15RR0.17G (0.17M) ○ CJO-15RR0.65G (0.65M) ○ CJO-15RR1.5G (1.5M) Special Order Item ○ CJO-15RR3G (3M)	 ○ HBX-800 BISx™ Module Interface	 *shown connected to HBX-800 BISx™ Module Interface ○ 186-0195-SF BISx™ System

Print



DS-8500

Unit Connection Cable (U-LINK)



- CJO-09SS1.5 (1.5M)


Special Order Item

- CJO-09SS0.3 (0.3M)
- CJO-09SS5 (5.0M)



- **HR-800** 3CH Thermal Printer
use with paper OP050-01TDR

Generic Ethernet Cable




- **HP** Network Printer

Agent



DS-8500

Unit Connection Cable (U-LINK)



- CJO-09SS1.5 (1.5M)

Special Order Item

- CJO-09SS0.3 (0.3M)
- CJO-09SS5 (5.0M)



- **MGU-801P** AGO2 Gas Unit


*Contact your local Mindray Representative for Water Trap and Sampling Tube

Accessories



HLX-801 (600MHz)
Telemetry Transmitter
Module

Special Order Item



OAT-8185A
Mounting Bracket



CF-820
Remote Controller



FCF-16GA
(9B0102420)
16GB Full Disclosure Card
*for 48 Hours of Full Disclosure



SD-8G
8GB SD Card
*to store patient data



SD-1G
1GB SD Card
*to store patient data



FCF-128
(1Z7333B)
128MB CF Card
Flash Memory Card
*for transferring monitor settings

DS-8900

CENTRAL STATION






DS-8900 - Central Station Main Unit



LC-8026T - External Display

Print

Unit Connection Cable (U-LINK)



○ **CJO-09SS1.5** (1.5M)

Special Order Item


○ **CJO-09SS0.3** (0.3M)

○ **CJO-09SS5** (5.0M)



○ **HR-800** 3CH Thermal Printer
use with paper OP050-01TDR

Ethernet Cable




○ **CJO-17SS2.5** (2.5M)




○ **HP** Network Printer

Accessories



OP050-01TDR
Thermal Roll Paper



CJO-17SS2.5 (2.5M)
Ethernet Cable



FUKUDA DENSHI USA, INC.
17725-C NE 65th Street Redmond, WA 98052
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (425) 869-2018
www.fukudaamerica.com

BIS™, BISx™, Nellcor™ and Microstream™ are registered trademarks of Medtronic. Masimo® is a registered trademark of Masimo.
PHILIPS RESPIRONICS® is a registered trademark of Koninklijke Philips. Trimline™ and Tempa-Kuff® are registered trademarks of Welch Allyn.

The information contained in this document is subject to change without notice.

© 2023 FUKUDA DENSHI USA, INC. All rights reserved.