

DS-1200 Bedside Monitor Accessory List



FUKUDA

Accessible Healthcare for Everyone

DS-1200 Bedside Monitor Main Configuration

- ECG
- NIBP
- SpO2
- TEMP
- IBP
- CO
- EtCO2
- BIS
- GAS AGENT



Model	SpO2 Unit	Recorder	Telemetry
DS-1200MR	Masimo®	Yes	
DS-1200MRT	Masimo®	Yes	Yes
DS-1200NR	Nellcor™	Yes	
DS-1200NRT	Nellcor™	Yes	Yes



DS-1200MR
Masimo® SpO2

DS-1200MRT
Masimo® SpO2 /
Telemetry



DS-1200NR
Nellcor™ SpO2

DS-1200NRT
Nellcor™ SpO2 /
Telemetry

Additional Modules



HC-110
EtCO2 Oridion
Sidestream
Module
*1 Slot Module



HC-120
EtCO2 Respiration
Mainstream
Module
*1 Slot Module



MG-110
Multigas Unit
CO2/N2O/AGENT/O2
*2 Slot Module



HP-800
Serial Port
Module
*1 Slot Module



HM-800
Multi-Parameter
Module
*1 Slot Module



MG-120
Multigas Unit
CO2/N2O/AGENT/
O2/SPIROMETRY
*3 Slot Module



HG-820
Nellcor™
SpO2 Module
*1 Slot Module




HG-810
Masimo®
SpO2 Module
*1 Slot Module

ECG




CIO-07CTP-3NA
3 Lead ECG Trunk Cable



CMO-07FT-3NAB
3 Lead ECG Patient Cable




CIO-07CTP-4NA
4 Lead ECG Trunk Cable




CMO-07FT-4NAB
4 Lead ECG Patient Cable



CIO-07CTP-5NA
5 Lead ECG Trunk Cable




CMO-07FT-5NAB
5 Lead ECG Patient Cable



CMO-07FT-10NAB
12 Lead ECG Trunk & Patient Cable for 10-electrode

NIBP




OA-80APS3.5-S
3.5M Cuff Hose

Special Order Item

OA-80APS1.5-S - 1.5M Cuff Hose



- CUF-1000-XS** - Adult Extra Small Cuff (4.7" - 6.7")
- CUF-1000-S** - Adult Small Cuff (6.7" - 8.7")
- CUF-1000-M** - Adult Medium Cuff (8.7" - 12.6")
- CUF-1000-L** - Adult Large Cuff (12.6" - 16.5")
- CUF-1000-XL** - Adult Extra Large Cuff (16.5" - 19.7")



OA-80NE1.5-S
1.5M Neonate Cuff Hose


OA-80NE3.5-S
3.5M Neonate Cuff Hose

Use with SunTech® Neonatal Soft Disposable BP Cuff


SpO2



1M1120
Nellcor™ SpO2 Sensor Extender



OA-MAX-AL
Nellcor™ Disposable SpO2



OA-DS-100A
Nellcor™ Durasensor SpO2

*Only sold with equipment orders. Not available for individual sale.


*Contact your local Nellcor/Medtronic Representative for forehead, ankle, child and other disposable sensors

DS-1200M Series

- OA-R8-MASIMO** - Masimo Starter Kit - Reusable
- OA-D8-MASIMO** - Masimo Starter Kit - Disposable

*Only sold with equipment orders. Please contact Masimo for accessories.

TEMP



CJO-P01T-DA0.5
0.5M 2ch Temperature Relay Cable

CJO-P01T-DA4.0
4.0M 2ch Temperature Relay Cable



CS409
Temperature Probe Surface

IBP



CJO-P01B-DJ0.5
2ch BP Conversion Cable



CJO-P01B-SB3.6
3.6M DTX Plus

FDBETWPX - Baxter Edward Interface Cable

**For use with Baxter Edwards TruWave PX-Series Transducer*

FDICUTRP4 - ICU Medical Transpac® IV Transducer Cable

**For use with ICU Medical Transpac® 4 Transducer*

CO



CJO-P01C-C2.4
Catheter Relay Cable



CJO-P01C-T2.4
Injectate Probe Relay Cable



CJO-P01C-L2.4
In-line Sensor Relay Cable




CJO-P01C-F2.4
Flow-through Sensor Relay Cable

EtCO2



HC-110
EtCO2 Oridion
Sidestream Module

*1 Slot Module



Covidien
Microstream™

**Contact your local
Medtronic Representative
for sampling tube*



HC-120
EtCO2 Respironics
Mainstream Module

*1 Slot Module



Philips
Respironics®

**Contact your local
Philips Representative*

BIS



CJO-29TR0.4
AUX Lem Communication Cable



186-0195-SF
BISx System



GAS AGENT



MG-110
Multigas Unit
CO₂/N₂O/AGENT/O₂

*2 Slot Module



MG-120
Multigas Unit
CO₂/N₂O/AGENT/O₂/SPIROMETRY

*3 Slot Module

Water Traps

115-058733-00 - Adult DRYLINE™ II Water Trap (Box of 10)

115-043024-00 - Adult/Pediatric DRYLINE™ II Water Trap (Box of 10)

115-043025-00 - Neonate DRYLINE™ II Water Trap (Box of 10)

Sampling Lines

115-043017-00 - Adult/Pediatric DRYLINE™ Sampling Line (2.5M) (Box of 25)

115-043018-00 - Neonate DRYLINE™ Sampling Line (2.5M) (Box of 25)

Airway Adapters

115-043021-00 - Adult/Pediatric DRYLINE™ Airway Adapter (Elbow) (Box of 25)

115-043020-00 - Adult/Pediatric DRYLINE™ Airway Adapter (Straight) (Box of 25)

Other



BTO-005
Rechargeable Battery



OP050-02TDR-BX **OP050-02TDR-CS**
Thermal Paper (Box) Thermal Paper (Case)

**This paper is different from the 8000 series paper*



FSD-32G
Full Disclosure SD Card
**10 Days (240 hours) of Full Disclosure*

Mounting Solution



OAO-70A
GCX Attachment for Monitor



FD-0009-21D
Wall Mount Fixed Height



FD-0026-01
Wall Mount Variable Height



WS-0003-25A
Quad Utility Hook



WMM-0006-04
Dual Utility Hook



RS-0025-02
Roll Stand



OTO-16S
Fukuda Roll Stand
**DS-1200 not pictured*



OAO-1011A
Balancer
**Need when using OTO-16S*



OTO-16L
Fukuda Roll Stand
**DS-1200 not pictured*



Accessible Healthcare for Everyone

<http://www.fukudaamerica.com>



FUKUDA DENSHI USA, INC.
17725-C NE 65th Street Redmond, WA 98052
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (888) 224-7090

BIS™, BISx™, Medtronic®, Nellcor™ are registered trademarks of Medtronic. Masimo® is a registered trademark of Masimo Corporation. PHILIPS RESPIRONICS® is a registered trademark of Koninklijke Philips. SunTech® is a registered trademark of SunTech Medical, Inc. DRYLINE™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Transpac® is a registered trademark of ICU Medical, Inc.

The information contained in this document is subject to change without notice.
© 2025 FUKUDA DENSHI USA, INC. All rights reserved. Document No. USRB-0273-01