

DS-1800 Central Station Accessory List



FUKUDA

Accessible Healthcare for Everyone

DS-1800 Central Station Main Configuration



Model	Number of Telemetry Receiving Beds	Printer	Extended Display
DS-1800L	0		
DS-1800LR	0	Yes	
DS-1800LRE	0	Yes	Yes
DS-1812	12		
DS-1812R	12	Yes	
DS-1812RE	12	Yes	Yes

Accessories



BTO-005
Rechargeable Battery



OP050-02TDR-BX
Thermal Paper (Box)

OP050-02TDR-CS
Thermal Paper (Case)

**This paper is different from the 8000 series paper
The image is of the DS-1200




FSD-64G
Full Disclosure SD Card
**14 Days (336 hours) of Full Disclosure*

Recorder Unit

**The DS-1800 has a built-in recorder*



**To externally connect a recorder unit*




CJO-09SS1.5
Unit Connection Cable
(U-LINK) (1.5M)



HR-800
Recorder Unit

Laser Printer

Direct Connection Example




Z09124
CAT5E - Patch Cable



HP® Network Printer


Shared Connection Example



Z09124
CAT5E - Patch Cable



Z89116
8 Port Ethernet Hub



Z09124
CAT5E - Patch Cable



HP® Network Printer

Touchscreen Monitor



ELO-KIT-YCABLE
Interface Cable



CJ-726
Interface Cable



**Image only. Not the actual product.*

DISPLAY-27ELO
Elo® Display
Medical Grade - LCD Monitor
27-inch



**Image only. Not the actual product.*

DISPLAY-22ELO
Elo® Display
Medical Grade - LCD Monitor
22-inch

Mouse and Keyboard

Mouse

**Image only. Not the actual product.*



USB Wired Mouse

**Please contact your Sales Representative for more information.*

With USB 2.0 Interface (Type A)

Compliant with IEC 62368-1, CISPR 32, CISPR 35 or equivalent standard

When one mouse is connected to the USB connector, both display units can be controlled by the same mouse.

When two mice are connected to the USB connector, each display unit can be operated separately with each mouse.

Keyboard

**Image only. Not the actual product.*



USB Wired Keyboard

**Please contact your Sales Representative for more information.*

With USB 2.0 Interface (Type A)

Compliant with IEC 62368-1, CISPR 32, CISPR 35 or equivalent standard

The keyboard can be used when entering patient information.

NOTE

The keyboard can be used only when the touch panel keyboard is displayed (patient name, monitor suspend setup, etc.)

The keyboard cannot be used for the display other than above (password input, etc.)



Accessible Healthcare for Everyone

<http://www.fukudaamerica.com>



FUKUDA DENSHI USA, INC.
17725-C NE 65th Street Redmond, WA 98052
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (888) 224-7090

HP® is a registered trademark of HP Development Company, LP. Elo® is a registered trademark of Elo Touch Solutions, Inc.

The information contained in this document is subject to change without notice.
© 2023 FUKUDA DENSHI USA, INC. All rights reserved. Document No. USRB-0271-00