

Specification

[illegible]

FUKUDA DENSHI USA, INC.
17725-C NE 65th Street Redmond, WA 98052
Toll Free: (800) 365-6668 / Local: (425) 881-7737 / Fax: (888) 224-7090

Medtronic® and Nellcor™ are registered trademarks of Medtronic. Masimo® is a registered trademark of Masimo Corporation.
Exergen® is a registered trademark of Exergen Corporation.
The information contained in this document is subject to change without notice.
© 2026 FUKUDA DENSHI USA, INC. All rights reserved. Document No. USRB-0269-01

<http://www.fukudaamerica.com>



DYNASCOPE

Bedside Monitor

BDS-1001 System



Medical Innovation for a Brighter Tomorrow

Portrait Style Monitor

Vertically oriented monitor with a tilted display which facilitates easy viewing for the healthcare professionals.


Space Saving Compact Size

Approx. 10-inch display that can display both waveform and numerics




User Friendly Screen Configuration

The screen is designed to display only the necessary information and allow for a more flexible screen layout. Can be configured to display only measurement values or waveforms.




Patient Scoring


SI (Shock Index) - Index that serves as a guide to bleeding and/or hemorrhagic shock



EWS (Early Warning Score) - Indicators that can assist in the early detection of physiological changes in patients



RPP (Rate Pressure Product) - Displays an index for evaluating the state of cardiac function



NEW

Seamless Temperature Data Integration

Body temperature data can be automatically captured and sent directly to the EMR.*

- Streamlined workflow from measurement to documentation
- Reduced manual entry and transcription errors
- Centralized management of vital signs, including body temperature

*Optional accessories, software updates, and EMR integration setup may be required.



Long Duration Waveform

Can be used as a patient data management tool in an environment without a central station or departmental system, and also with optional SD card (SD-16G) can back up post-discharge data for up to 14 days.

CO2 Gas Unit

Optional CO2 gas unit, HCP-110, available for CO2 measurement.




Optional HCP-110 Attached

Discharge		2023/ 8/ 1 12:57
Admit Date/Time		2023/ 8/ 1 12:56
ID	NAME	Admit Date/Time
ID: FUKUDA20232	Name: FUKUDA2	Admit Date/Time: 2023/ 8/ 1 12:57
ID: FUKUDA20233	Name: FUKUDA3	Admit Date/Time: 2023/ 8/ 1 12:57
ID: FUKUDA20234	Name: FUKUDA4	Admit Date/Time: 2023/ 8/ 1 12:56
ID: FUKUDA20235	Name: FUKUDA5	Admit Date/Time: 2023/ 8/ 1 12:56

Inflation NIBP Measurement Method

Can measure blood pressure gently and painlessly even for hemodynamically unstable patients with slow pulses, small pulse wave amplitudes, and low blood pressure.



Fukuda Original Cuffs

